

NICODOM IR Drugs, 1463 spectra

The infrared spectral library „NICODOM IR Drugs“ is an advanced collection of 1463 FTIR (ATR) digital spectra of drugs and related compounds.

Typical information about the sample contains:

Code
Name
Synonyms (up to 10 different synonyms)
CAS Number
Molecular weight
Chemical Formula
Manufacturer

Example of included information:

Code: DG0001

Name: Acebutolol Hydrochloride

CAS: 34381-68-5

Formula: C₁₈H₂₉ClN₂O₄

Molecular weight: 372.891,

Synonyms: N-[3-acetyl-4-(2-hydroxy-3-[isopropylamino]propoxy)phenyl]]-butanamide

Manufacturer: Sigma

The library includes spectra of analytical drug standards, street drugs and other helping substances.

This FTIR spectral library can be purchased separately or as a part of cost effective library package “NICODOM IR Toxic/Hazardous Package, 7313 spectra”.

The Nicodom FTIR spectra of toxic and hazardous compounds were collected using the FTIR Spectrometer Nexus 670™ (Thermo) and the Miracle™ single bounce ATR accessory (Pike Technologies) equipped by Diamond, ZnSe, or Si crystals.

The Nicodom spectral library of polymers features resolution of 4cm⁻¹ in the spectral range 600-4000 cm⁻¹, measurement time 1 minute, the instrument was purged by dried air.

The complete information is available in the full library version (not available in the demo list of compounds)

The last library update has been done 2011, the Nicodom IR Libraries include many new substances which came to market during the last decade.

This FTIR spectra library is available as digital searchable library compatible with your spectroscopic software. Also it can be delivered with SEARCH software or in *.pdf format.

To download the list of spectra, to download a free demo library compatible with your software, to check the format compatibility as well as for pricelist, ordering info and information about all NICODOM IR/NIR Libraries and other NICODOM products please visit our webpage.

List of Spectra - Chemical Names

Acebutolol Hydrochloride
Acyclovir
Ajmaline
Amlodipin Besylate
Amlodipine Maleate
Amoxicillin
Ampicillin
Ampicillin sodium
Antazoline sulfate
Atorvastatin Calcium
Atropine
Aviptadil
Bicalutamide
Boldenone 17-acetate
Boldenone undecylenate
Butobarbital
Captopril
Carbenicillin
Carprofen
Catechin
Cefaclor
Cialis
Cortisone Acetate
Diclofenac potassium
(+)-alpha-Dihydratetabenazine
(+/-)-N-Ethyl-3,4-methylenedioxyamphetamine.HCl
(R)-(+)-Desmethyisibutramine HCl
1-(1-Phenylcyclohexyl)morpholine .HCl
1-(1-Phenylcyclohexyl)pyrrolidine .HCl
1-[1-(2-Thienyl)cyclohexyl]morpholine (TCM)
1-[1-(2-Thienyl)cyclohexyl]piperidine (TCP)
17a-Epitestosterone
17b-Hydroxy-5b-androstan-3-one
17-Methyltestosterone
19-Nortestosterone decanoate
1-Phenyl-1-Butanol
2-(2,4,5-Trichlorophenoxy)Propionic Acid
2,3-Dimethylbenzoic Acid
2,4,5-Trimethoxyamphetamine
2,4,6-Trimethoxyamphetamine
2,5-Dimethoxy-4-ethylamphetamine
2,5-Dimethoxy-4-methylamphetamine.HCl; STP.HCl
2,5-Dimethoxy-4-methylamphetamine; STP
2,5-Dimethoxyamphetamine .HCl
2C-I
2-Chloroprocaine .HCl
2-Methoxy-4,5-methylenedioxyamphetamine.HCl
3,4,5-Trimethoxyamphetamine
3,4-Dimethylbenzoic Acid
3,4-Methylenedioxyamphetamine; MDA
3alpha-hydroxytibolone
3-Chloropropylamine .HCl
4-Aminobenzoic acid; PABA
4-androsten-17beta-ol-3-one enanthate

NICODOM IR Drugs, 1463 spectra

Copyright © NICODOM 2006-2012

NICODOM Ltd., Hlavni 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>
ftir@ftir.cz, ftirspectra@gmail.com

List of Spectra Page 2 from 53

4-Hydroxybutyric acid, sodium salt
4-Chlorotestosterone Acetate
4-Methoxyamphetamine.HCl
5-Pregnen-3b-Ol-20-On
6-Monoacetylmorphine
7-(Beta-Hydroxyethyl)Theophylline
7-(Beta-Hydroxypropyl)Theophylline
8-Chlortheophylline
2-Aminotridecane Undecylenate
2-Chloro-4-nitrophenol
4-Methylumbelliferone
5-Aminosalicylic acid
5-Sulfosalicylic acid
Acetaminophen
Acetazolamide
Acetophenetidin
Acetylsalicylic acid
Acriflavine
Acrinol
Adiphenine HCl
Allantoin
Allobarbitol
Allylcyclopentenyl Barbituric Acid
Alphaprodine HCl
Alphenal
Alprazolam
Ambroxol .HCl
Aminoglutethimide
Aminophylline
Aminopyrine
Amitriptyline
Amobarbital Sodium
Amodiaquine Hydrochloride
Amphetamine .HCl
Androsterone
Anileridine HCl
Anisotropine methylbromide
Antazoline mesylate
Apomorphine HCl
Aprobarbital
Ascorbic Acid:Dextrose 1:1
Ascorbic Acid:Sucrose 1:1
Atabrine hydrochloride
Atenolol
a-Tocopheryl acetate hydrate
Atropine sulfate
Azidamfenicol
Bacitracin Zinc
Barbital
Barbital Sodium
Barbituric Acid
Benactyzine .HCl
Benzocaine
Benzoylecgonine
Benzquinamide
Benztropine Mesylate

NICODOM IR Drugs, 1463 spectra

Copyright © NICODOM 2006-2012

NICODOM Ltd., Hlavni 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

ftir@ftir.cz, ftirspectra@gmail.com

List of Spectra Page 3 from 53

Betamethasone
Bethanechol
Bisabolol
Bismuth Tribromophenate
Bromazepam
Bromisovale
Bromodiphenhydramine HCl
Bufotenine
Bufotenine monooxalate
Buprenorphine HCl
Butabarbital
Butacaine sulfate
Butalbital
Butamirate citrate
Butethamate citrate
Caffeine
Caffeine Citrate
Caffeine Sodium Benzoate
Caffeine:Aspirin 1:1
Caffeine:Ephedrine 1:1
Caffeine:Paracetamol 1:1
Caffeine:Sucrose 1:1
Calcium DL-pantothenate
Calcium Dobesilate
Calcium D-pantothenate, Vitamin B5
Calcium Lactate
Cannabicyclol
Cannabinol
Carbamazepine
Carbinoxamine maleate
Carbromal
Cefadroxil
Cefaloridine
Cefixime
Cefotaxime sodium
Cefpodoxim Proxetil
Cefpodoxime
Cefuroxime Sodium
Cinchocaine.HCl
Cinnarizine
Ciprofloxacin
Citalopram Hydrobromide
Clavulanate Potassium
Clobazam
Clomipramine
Clonazepam
Clopamide
Clorazepate dipotassium
Clortermine HCl
Cloxacillin
Clozapine
Cocaine
Cocaine.HCl
Cocaine:Caffeine 1:1
Cocaine:Fructose 1:1
Cocaine:Lactose 1:1

Cocaine:Procaine 1:1
Cocaine:Sucrose 1:1
Codein Phosphate
Codeine
Codeine .HCl
Codeine Phosphate
Colecalciferol (Vitamin D3)
Collodium
Condylina, Podophyllotoxin
Cyanocobalamine (Vitamin B12)
Cyclazocine
Cyclomethycaine sulfate
Cyclopentamine .HCl
Cycrimine .HCl
Cyproheptadine .HCl
D-(-)-Fructose
D-(+)-Lactose Powder
d-2-Bromolysergic acid diethylamide
D-Alanine
D-Arginine
D-Arginine hydrochloride
D-Ascorbic acid
D-Aspartic acid
D-a-Tocopherol type V
d-Biotin
Debrisoquine sulfate
Dehydroisoandrosterone
Delorazepam
Demoxepam
Deserpidine
Desipramine .HCl
Desmethyldiazepam
Dexamethasone Acetate
Dexbrompheniramine maleate
Dexpanthenol
Dextromethorphan
Dextrorphan-D-Tartrate
D-Glucose
Diaveridine
Diazepam
Dibucaine
Diclofenac Sodium
Diethylpropion .HCl
Diethylstilbestrol
Diethyltryptamine .HCl
Digitoxin
Digoxin
Dihydrocodeine Bitartrate
Dihydromorphine
Dihydroxypropyladenine (DHPA)
Diiodohydroxyquin
Diltiazem Hydrochloride
Dimethindene Maleate
Diphenhydramine HCl
Diphenidol HCl
Diphenoxylate HCl

Diphenylpyraline
Dipropyltryptamine HCl
Dipyron
D-Isoascorbic Acid
Disulfiram
Dithiaden
Dithranole
DL-6,8-Thioctic Acid
DL-Alanine
dl-Amphetamine sulfate
DL-Asparagine monohydrate
DL-Aspartic acid
DL-a-Tocopherol acetate, Vitamin E acetate
DL-Cystathionine
DL-Dithiothreitol
DL-Ephedrine Hydrochloride
dl-Epinephrine
dl-Epinephrine .HCl
D-Leucine
DL-Methionine
DL-Phenylalanine
DL-Tryptophan
Donepezil HCl
Dopamine
Dosulepine
Doxepin HCl
Doxycycline
Doxylamine Succinate
D-Phenylalanine
Droperidol
D-Tryptophan
D-Tyrosine
D-Valine
Dyclonine HCl
Dyphylline
Ecgonine
Ecgonine methyl ester
Ephedrine Sulfate
Ephedrine:Dextrose 1:1
Ephedrine:Fructose 1:1
Ephedrine:Lactose 1:1
Epiandrosterone
Epinephrine Bitartrate
Ergocalciferol Vitamin D2
Ergonovine Maleate
Ergosterol, Provitamin D
Ergotamine Tartrate
Esculin
Ethacridine Lactate Monohydrate
Ethinamate
Ethionamide
Ethosuximide
Ethyl Biscoumacetate
Ethylmorphine .HCl
Ethylnorepinephrine hydrochloride
Etiocholanolone

NICODOM IR Drugs, 1463 spectra

Copyright © NICODOM 2006-2012

NICODOM Ltd., Hlavni 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

ftir@ftir.cz, ftirspectra@gmail.com

List of Spectra Page 6 from 53

Etorphine
Famotidine
Fenbufen
Fenfluramine .HCl
Fenofibrate
Fenofibric Acid
Fenoprofen
Fenoterol Hydrobromide
Fenpiverinium bromide
Fentanyl Citrate
Fluconazole
Flunitrazepam
Fluphenazine HCl
Flurazepam.2HCl
Flurazepam.HCl
Flurbiprofen
Folic Acid (Vitamin B11)
Furazolidone
Furosemide
Glipizide
Glutethimide
Glybenclamide
Guaiacol Glyceryl Ether
Halazepam
Haloperidol
Harmaline
Heparine 500000 Units
Heroin
Heroin HCl
Heroin:Caffeine:Paracetamol 1:1:1
Heroin:Dextrose 1:1
Heroin:Fructose 1:1
Heroin:Inositol 1:1
Heroin:Lactose 1:1
Heroin:Procaine 1:1
Heroin:Sucrose 1:1
Hesperetin
Hexobarbital
Hexocyclium Methylsulfate
Histamine
Homatropine hydrobromide
Hydralazine HCl
Hydrocodone Bitartrate
Hydrocortisone Acetate
Hydrochlorothiazide
Hydromorphone
Hydromorphone HCl
Hydroxyamphetamine hydrobromide
Hydroxyzine Dihydrochloride
Chenodeoxycholic Acid
Chloral betaine
Chloral hydrate
Chloramphenicol
Chlordiazepoxide
Chlordiazepoxide.HCl
Chlorhexidine Digluconate, 20% solution

Chlorothiazide
Chlorpheniramine maleate
Chlorphentermine .HCl
Chlorpromazine Chloride
Chlorpropamide
Chlorprothixene
Chlorthalidone
Chlorzoxazone
Cholesterol
Ibogaine .HCl
Ibuprofen
Ibuprofen Sodium
Imipramine .HCl
Indapamide
Indole
Indomethacin
Ioxitalamic Acid
Irbesartan
Isobarbituric Acid
Isocarboxazid
Isoniazid
Isopropamide Iodide
Isosorbide 5-Mononitrate
Isoxicam
Itraconazole
Kebuzone
Ketamine HCl
Ketoconazole
Ketoprofen
L-(+)-Lysine monohydrochloride
l-a-Acetylmethadol .HCl
l-a-Acetyl-N-normethadol .HCl
L-Alanine
Lamotrigine
Lansoprazole
L-Arginine
L-Arginine HCl
L-Asparagine
L-Asparagine *H2O
L-Aspartic Acid
L-Cysteine
L-Cysteine HCl
L-Cystine
L-Ephedrine.HCl
Levetiracetam
Levodopa; 3-Hydroxy-L-tyrosine
Levofloxacin
Levopropoxyphene napsylate
Levorphanol Tartrate
L-Histidine
L-Histidine Hydrochloride Monohydrate
L-Hydroxyproline
Lidocaine
Lidocaine .HCl
Lincomycin Hcl
L-Isoleucine

L-Leucine
L-Lysine acetate
L-Lysine Monohydrochloride
L-Methionine
L-Methyldopa
Loperamide HCl
Lorazepam
Lormetazepam
L-Ornithine
L-Ornithine hydrochloride
Losartan Potassium
L-Phenylalanine
L-Proline
LSD
L-Serine
L-Threonine
L-Tryptophan
L-Tyrosine
L-Valine
Lysergic acid
Maprotiline
Mazindol
MDMA
Mebutamate
Meclizine HCl
Meclofenamic Acid
Mecloqualone
Medazepam
Mefenamic acid
Meloxicam
Menadione, Vitamin K3
Mepenzolate Br
Meperidine .HCl
Mephexalone
Mephentermine
Mephentoin
Mephobarbital
Meprobamate
Mescaline .HCl
Mescaline hemisulfate
Methacycline.HCl
Methadone .HCl
Methamphetamine.HCl
Methapyrilene .HCl
Methaqualone
Metharbital
Methenamine Mandelate
Methimazole
Methixene.HCl
Methocarbamol
Methohexital
Methscopolamine Bromide
Methylclothiazide
Methylergonovine Maleate
Methylphenidate HCl
Methyprylon

Methysergide
Metoclopramide Hydrochloride
Metolazone
Metoprolol Tartrate
Metronidazole
Montelukast Na
Morphine
Morphine.HCl
Moxifloxacin HCl
Myo-Inositol
N,N-Diethyl-1-phenylcyclohexylamine HCl
N,N-Dimethyl-1-phenylcyclohexylamine HCl
N-Acetyl-L-cysteine
N-Acetyl-L-tyrosine
Naftidrofuryl oxalate
N-Allylmorphine
Naloxone
Nandrolone propionate
Naphazoline HCl
Naphazoline nitrate
Naproxen
N-Ethyl-1-phenylcyclohexylamine .HCl
N-ethyl-3-piperidyl Benzilate
Nicardipine .HCl
Nicotinamide
Nicotine Hydrogen Tartrate
Nicotinic acid
Nifedipine
Nitrazepam
Nitrofurantoin
N-Methyl-3,4-methylenedioxyamphetamine
Nonanoic acid vanillylamide
Noroxymorphine .HCl
Nortriptyline HCl
N-Propyl-3,4-methylenedioxyamphetamine
N-Propylamphetamine .HCl
Ofloxacin
Olanzapine
Opium powder
Ornidazole
Oxacillin
Oxazolam
Oxcarbazepine
Oxiconazole Nitrate
Oxybuprocaine Hydrochloride
Oxycodone
Oxycodone HCl
Oxymetazoline HCl
Oxymorphone
Oxyphenbutazone
Oxyphencyclimine HCl
Oxyphenonium Bromide
P-Aminophenol
Pantoprazole Sodium Salt
Papaverine hydrochloride
Paraquat

NICODOM IR Drugs, 1463 spectra

Copyright © NICODOM 2006-2012

NICODOM Ltd., Hlavni 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

ftir@ftir.cz, ftirspectra@gmail.com

List of Spectra Page 10 from 53

Paroxetine.HCl
Pemoline
Pentazocine
Pentazocine HCl
Pentobarbital Sodium
Pentoxifylline
Phenacetamide
Phenazocine HBr
Phenazone
Phencyclidine .HCl
Phendimetrazine Bitartrate
Phendimetrazine Hydrogen Tartrate
Pheniramine maleate
Phenmetrazine HCl
Phenobarbital
Phenobarbital Sodium
Phenoxybenzamine HCl
Phentermine
Phentermine HCl
Phenyl Salicylate
Phenylbutazone
Phenylephrine hydrochloride
Phenylpropanolamine HCl
Phenyltoloxamine Citrate
Phenytoin (5,5-Diphenylhydantoin)
Pholcodine
Phthalylsulfathiazol
Phytonadione, Vitamin K1
Pilocarpine
Pilocarpine Hydrochloride
Piperidolate
Piracetam
Piroxicam
Pitofenone
Pitofenone Hydrochloride
Podophyllotoxin
Potassium guaiacol sulfonate
Prazepam
Prazosin Hydrochloride
Pridinol Methanesulfonate Anhydrous
Prilocaine .HCl
Primidone
Probenecid
Procainamide
Procaine
Prochloroperazine Dimaleate
Promazine .HCl
Promethazine .HCl
Propafenone
Propantheline Bromide
Proparacaine HCl
Propoxyphene .HCl
Propranolol Hydrochloride
Propyphenazone
Protriptyline HCl
Pseudoephedrine

Pseudoephedrine.HCl
Psilocin
Psilocybin
Pyridoxal HCl
Pyridoxamine Dihydrochloride
Pyridoxine hydrochloride
Pyridoxine, Vitamine B6
Pyrilamine maleate
Quetiapine
Quinidine
Quinidine Gluconate
Quinidine Hydrochloride, Monohydrate
Quinidine Sulfate
Quinidine Sulfate Dihydrate
Quinine Hydrochloride
Quinine hydrochloride dihydrate
Quinine Sulfate
Quinine Sulfate Dihydrate
Ramipril
Ranitidine Hydrochloride
Reserpine
Retinol palmitate
Retinyl propionate
Riboflavine; Vitamin B2
Risedronate Sodium
Risperidone
Rosiglitazone Maleate
Roxithromycine
Rutoside
Salicylamide
Salicylic acid
Secobarbital
Sibutramine Hydrochloride
Silver Albumin Acetyl Tannate
Silver Albumin Acetyl Tannate with Borax
Sitosterol
Sodium Dichlorisocyanurate
Sodium Salicylate
Stanolone
Stigmasterol
Sucrose
Sulfadiazine
Sulfadimethoxine
Sulfamethazine
Sulfamethizole
Sulfamethoxazole
Sulfanilic acid
Sulfapyridine
Sulfathiazole
Sulpiride
Sumatriptan Succinate
Talbutal
Taurine
Temazepam
Tenoxicam
Terbinafine HCl

Testosterone
Testosterone Benzoate
Tetrabenazine
Tetracaine Hydrochloride
Tetracycline Hydrochloride
Tetraphenylethylene
Thebacon HCl
Thebaine
Theobromine
Theophylline
Thiabendazole
Thiamine hydrochloride
Thiamine, Vitamin B1
Thiamylal
Thiopental
Thioridazine
Thiothixene
Thymol
Tiapride Hydrochloride
Tiaprofenic Acid
Tibolone
Ticlopidine Hydrochloride
Tolazoline Hydrochloride
Tolbutamide
Tolfenamic Acid
Tramadol
Tranexamic Acid
Tranlycypromine
Trazodone HCl
Triamcinolone Acetonide
Triazolam
Trifluoperazine .2HCl
Trihexyphenidyl Hydrochloride
Trimecain Chloride
Trimeprazine Tartrate
Trimethadione
Trimethoprim
Trimipramine Maleate
Tripelennamine HCl
Tropacocaine HCl
Troloxerutine
Tryptanol
Tybamate
Ursodeoxycholic Acid
Valine
Verapamil .HCl
Viprynum Embonate
Vitamin B11, Folic Acid
Warfarine
Xanthinol Nicotinate
Yohimbine .HCl
Zolpidem Tartrate
Zoxazolamine
(+)-Hydroxymethyl tolterodine S-(+)- mandelate
(3,4-Dimethoxyphenyl)acetic acid
(S)-(+)-Ibuprofen

[5-[[4-[2-(5-Ethyl-2-pyridinyl)ethoxy]phenyl]methyl]-2,4-] thiazolidinedione hydrochloride
1-(2-Hydroxyethyl)-2-Hydroxymethyl-5-Nitroimidazole
1-(3-Carboxypropyl)-3,7-Dimethylxanthine
17-carboxybudesonide
2-(Methylsulfonyl)ethanol
2,6-Di-tert-butyl-(dimethylaminomethyl)phenol
2-[(1-Benzyl-4-piperidyl)methyl]-5,6-dimethoxy-2,3-dihydroinden-1-one
2-Methylhydrochlorothiazid
2-Oxazolidinone-5-Spiro-4'-[N-(2-Phenylethyl)Piperidine] Hydrochloride
3,5-Di-tert-butylcatechol
3,6-Dimethyl-2-Dimethylaminomethylphenol
3-Hydroxybenzoic acid
3-tert-butyl-4-hydroxyanisole
4'-Hydroxyphenyl-Carvedilol
5-[[4-[2-(5-Acetyl-2-Pyridinyl)Ethoxy]Phenyl]Methyl]-2,4-Thiazolidinedione
5'-Hydroxyphenyl-Carvedilol
5-Fluorouracil
5-Hydroxyindole-3-acetic acid
5-Chlorouracil
Acamprosate Calcium
Alendronate Sodium
Anastrozol
Aripiprazole
Asulam
Atomoxetine hydrochloride
Bacitracin
Berberine Chloride Hydrate
Bergapten
Betamethasone sodium phosphate
Bisoprolol Fumarate
Budesonide
Capecitabine
Capsaicin
Carvedilol
Cetirizine Hydrochloride
Clenbuterol hydrochloride
Clomiphene Citrate
Daidzein
Dehydro Aripiprazole
Desethoxytamsulosin
Dexamethasone
Dihydrocholesterol
Dihydroxydecarbamacapitabine
Diosgenin
Diosmetin
Diovan
Donepezil Hydrochloride
Doxifluridine
Dromostanolone Propionate
D-Serine
Emodin
Enalapril Maleate
Enalaprilat
Eplerenone
Eplerenone Hydroxyacid Potassium Salt
Escitalopram oxalate

NICODOM IR Drugs, 1463 spectra

Copyright © NICODOM 2006-2012

NICODOM Ltd., Hlavni 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

ftir@ftir.cz, ftirspectra@gmail.com

List of Spectra Page 14 from 53

Estriol
Eszopiclone
Exemestane
Fingolimod hydrochloride
Fingophosphate Racemate
Fluoromethalone
Gatifloxacin Sesquihydrate
Gentamycin sulfate
Gliclazide
Gliquidone
Histamine Dihydrochloride
Hydrindantin
Hydrobromothiazide
Hydrocortisone
Hydroxy Pentoxifylline
Hydroxy Pioglitazone (M-IV)
Hyperoside
Chelidonine
Chlorhexidine
Chlorhexidine acetate
Chlormequat Chloride
Cholestane-3b,5a,6b-triol
Cholesteryl stearate
Imatinib
Isopropyl unoprostone
Isotimolol
Lacidipine
Letrozole
Levocetirizine dihydrochloride
Lovacid Ammonium Salt
Lovastatin
Memantine hydrochloride
Mesterolone
Methenolone enanthate
Methylpioglitazone
Metoprolol Succinate
Mibolerone
Moclobemide
Mycophenolate Mofetil
N-alpha-(2,4-Dinitro-5-fluorophenyl)-L-valinamide
Nandrolone Base
Naringin
N-Isobutyryl-L-Cysteine
N-Methyl-gama-(1-naphthalenyloxy)-2-thiophenepropanamine
O-Demethyl Tramadol
Olanzapine Benzoate
Ondansetron
Oxandrolone
Oxymetholone
Paliperidone
Paroxdione
Paroxetine
Paroxetine Mesylate
Phenolthiadione
Pilocarpine nitrate
Polymyxin B sulfate

NICODOM IR Drugs, 1463 spectra

Copyright © NICODOM 2006-2012

NICODOM Ltd., Hlavni 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

ftir@ftir.cz, ftirspectra@gmail.com

List of Spectra Page 15 from 53

Pramipexole
Prednisolone
Prednisolone phosphate sodium
Pridinol Mesylate
Progesterone
Protopine hydrochloride
Pyrethoxypioglitazone
Pyrhydroxyapioglitazone
Pyrketopioglitazone
Quetiapine Fumarate
Quinapril Hydrochloride
R(-)-alpha-Phellandren
R(-)-Bicalutamide
Raloxifene Hydrochloride
Raloxifene-4'-glucuronide
Raloxifene-6'-D-glucuronide
Ramiprilat
Resveratrol
Ribavirin
Rivastigmine Tartrate
Ropinirole hydrochloride
Rosuvastatin calcium
S(+)-Bicalutamide
Simvacid Ammonium
Simvastatin
Somatropin
Stanozolol
Tamoxifen
Tamsulosin hydrochloride
Telmisartan
Testosterone acetate
Testosterone cypionate
Testosterone decanoate
Testosterone isocaproate
Testosterone phenylpropionate
Testosterone propionate
Tianeptine sodium salt
Timolol Maleate
Tobramycin
Tolterodine Tartrate
Trans Paroxyl
Trazodone hydrochloride
Trenbolone Acetate
Trenbolone Enathate
Triatricol
Tylosin Tartrate
Tyloxapol
Venlafaxine Base
Venlafaxine Besylate Monohydrate
Venlafaxine Hydrochloride
Venlaphenol
Zopiclone
(-)-alpha-Santonin
(-)-Bornyl acetate
(-)-Camphene
(-)-Epicatechin

(-)-Linalool
(+)-Terpinen-4-ol
(+/-)- α -Pinene
(1,2,3,4-Tetrahydro-1-naphthyl)-formylethyldiamine hcl
(E)-Anethole
(R)-(-)-Carvone
(R)-(+)-1-Phenylethyl isocyanate
1-(4-nitrophenyl)-2-amino-1,3-propanediol
1,2,4-Butanetriol
1,4-Androstadien-17 α -methyl-17 β -ol-3-one
1,4-Butanediol
1,4-Dioxane
11b,17a-Dihydroxy-6 α -methyl-1,4-pregnadiene-3,20-dione
1-Cyano-1,2,3,4-tetrahydronaphthalene
1-Ethylpiperidine
1-Hexanol
1-Methyl-L-Histidine
1-Naphtholbenzein
1-Octadecanol
1-Octanol
1-Phenyl-ethylamine
1-Thioglycerol
2,2,6,6-Tetramethyl-4-piperidone
2,2-Diphenyl-1-picrylhydrazyl
2,3,4,5,6-Pentafluorothiophenol
2',3,4',5,7- Pentahydroxyflavone
2,3,4,5-Tetrachloronitrobenzene
2,3,4,6-Tetra-O-acetyl-beta-D-glucopyranosyl isothiocyanate
2,3,4-Trinitrotoluene
2,4,5-Trinitrotoluene
2,4,6-Triaminotoluene
2,4-Diamino-6-nitrotoluene
2,5-Dimethoxytetrahydrofuran
2,6,10,15,19,23-Hexamethyltetracosan
2,6-Dichloroindophenol sodium salt
2,6-Dimethylphenol
2,6-di-t-butyl-4-methylphenol, BHT
2,3-Dimercapto-1-propanol
2-Aminoisobutyric Acid
2-Ethylhexanol
2-Hydroxyhippuric acid
2-Chloroethylamine hydrochloride
2-Thiobarbituric acid
3-(N-Methylpyrrolidino)pyridine
3-(N-Morpholino)propanesulfonic acid; MOPS
3,3',5,5'-Tetramethylbenzidine
3,3',3'',3'''-(1,2-Ethanediyldisnitrilo)tetrakis(propionic acid sodium) salt
3,3-Dimethylnaphthidine
3,3-Methylene-bis-(4-hydroxycoumarin)
3,4-Dihydroxycinnamic acid
3,4-Dinitrophenol
3,5-Dinitrophenol
3,5-Dinitrotoluene
3,7-Dimethyl-1,6-octadien-3-yl acetate
3-Acetylthiophene
3-Amino-5-methylisoxazole

NICODOM IR Drugs, 1463 spectra

Copyright © NICODOM 2006-2012

NICODOM Ltd., Hlavni 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

ftir@ftir.cz, ftirspectra@gmail.com

List of Spectra Page 17 from 53

3-Ethyl-3-methylglutaric acid
3-Methoxy-L-tyrosine monohydrate
3-Methyl-L-histidine
4-(Ethylaminomethyl)-pyridine
4,5,7-Trihydroxy-isoflavon
4,7-Dichloroquinoline
4-Amino-3,5-dinitrotoluene
4-Amino-3-hydrazino-1,2,4-triazol-5-thiol
4-Aminohippuric acid
4-Estren-7-alpha-methyl-17-beta-ol-3-one acetate
4-Fluorobenzenesulfonate Sodium Salt
4-hydroxy-3-morpholino-1,2,5-thiadiazole
4-Hydroxy-3-morpholino-1,2,5-thiadiazole-1-oxide
4-Hydroxybenzoic acid butylester
4-Chloro-17a-methyl-17b-hydroxy-1,4-androstadiene-3-one
4-Methyl hydroxybenzoate = Methylparaben
4-Methylcyclohexanone
4-methylumbelliferyl phosphate
5,7-Dihydroxy-4-methylcoumarin
5,8,11,14,17- Eicosapentaenoic acid
5-Fluoro-2-methoxy-4(1H) pyrimidinone
5-Iodouracil
5-Keto-D-gluconic acid potassium salt
6-Aminohexanoic acid
6-Azathymine
6-Azauracil
6-Methylaminopurine
6-Piperidinopurine
7,7,8,8-Tetracyanoquinodimethane
9-Aminoacridine Hydrochloride Monohydrate
Acesulfame K
Acetic Acid
Acetone
Acetonitrile
alpha-Cyclodextrin hydrate
Adenine
Adenine Sulfate
alpha-D-Glucose 1-phosphate, dipotassium salt
Adipic acid
Adonitol
Agar
alpha-Chloralose
Aloe-Emodin
Amantadine .HCl
Amaranth
Ammonium acetate
Ammonium carbonate
Ammonium ferric citrate
Ammonium hydrogen carbonate
Ammonium Phosphate
Amyl Acetate
Amylopectine
Anethole
Ascorbyl Palmitate
Aspartame
Azodicarbonamide

NICODOM IR Drugs, 1463 spectra

Copyright © NICODOM 2006-2012

NICODOM Ltd., Hlavni 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

ftir@ftir.cz, ftirspectra@gmail.com

List of Spectra Page 18 from 53

Benflumetol
Benzalacetone
Benzoic acid
Benzyl acetate
Benzyl alcohol
Benzyl benzoate
Beta Alanine
Beta-Alanine methyl ester hydrochloride
Beta-Arthemeter
beta-Carotene
Beta-cyclodextrin, sulfated sodium salt
Beta-cyclodextrine
Bicarbonate of soda
Bicuculline
Bilirubin
b-Ionone
Bismuth Gallate Basic
Boldenone
Boldenone Propionate
Borneol
Buflomedil hydrochloride
Butylated hydroxyanisole, BHA
Cacotheline
Calcium Acetate
Calcium Disodium Edetate
Calcium Hydrogenphosphate
Calcium hydroxide
Calcium Lactobionate
Calcium oxide
Calcium Phosphate Monobasic
Calcium sulfate
Camphoric acid
Caramel
Carboxy methyl cellulose
Carisoprodol
Carnauba wax
Carrageenan
Carvacrol
Casein sodium salt
Castor Oil
Ceftriaxone impurity A
Ceftriaxone sodium
Cefuroxime
Cetrimonium Tosylate
Cineole
Cinnamyl alcohol
cis-4,7,10,13,16,19-Docosahexaenoic acid
Citral
Citric Acid
Citronellal
Coccarboxylase tetrahydrate
Cocoamidopropyl Hydroxysultaine
Cochineal
Colistin Sulfate
Cyanoguanidine
Cyanuric acid

Cyclohexanol
Cyclohexylamine
Cyprofloxacin hydrochloride
Cystamine dihydrochloride
Cytidine
Cytosine
D-(-)-Arabinose
D-(-)-Ribose
D-(+)-Fucose
D-(+)-Glucosamine .HCl
D-(+)-Malic acid
D(+)-Raffinose pentahydrate
D-(+)-Xylose
Dabsyl Chloride
Dansyl Fluoride
Dansyl Hydrazine
D-Arabonic gamma-lactone
D-Cystine
Decafluorobenzophenone
Dehydrocholic acid
Deoxycholic acid; sodium salt monohydrate
Desoxyribose
Dextran 60000-80000
D-Galactose
D-Glucoheptono-1,4-lactone
D-Glucose 6-phosphate disodium
Dibutyl phthalate
Diethyl ether
Diethyl phthalate
Diethylene Glycol Dibutyl Ether
Diethylstilbestrol
Dihydrostreptomycin sesquisulfate
Dimethicone
Dimeticone
Dioctyl phthalate
Diphenyl
Dipotassium Phosphate
Dipotassium tartrate
Disodium 5-Inosinate
Disodium citrate
Disodium edetate dihydrate
Disodium EDTA
Disodium phosphate
Disodium tartrate, dihydrate
D-Isoleucine
DL-alpha-Methoxyphenylacetic acid
DL-Citronello
DL-DOPA = 3-(3,4-dihydroxyphenyl)- DL-alanine
DL-Histidine Hydrochloride
DL-Homocystine
DL-Isocitric acid trisodium
DL-Lactic acid, Sodium salt
DL-Malic acid
DL-Menthol
DL-Monosodium glutamate
DL-Proline

DL-Serine
DL-Tartaric Acid
DL-Tyrosine
DL-Valine
D-Mannitol
D-Mannoheptulose
D-Mannose
D-Rhamnose
D-Sorbitol
D-Tagatose
D-Talose
D-Tartaric Acid
D-Threonine
Dulcitol
E104(2)-Quinoline Yellow Lake
E104-Quinoline Yellow
E107-Yellow 2G
E110-FD&C Yellow 6, Sunset Yellow FCF, Orange Yellow S
E111-Orange GGN, Acid Orange 20, alpha-Naphthol Orange
E1200-Polydextrose
E122-Carmoisine, Azorubine
E126-Ponceau 6R
E127-FD&C Red 3, Erythrosine
E131-Patent Blue V
E142-Greens S
E143-Fast green FCF
E1450-Starch sodium octenyl succinate
E1518-Glyceryl triacetate, triacetin
E151-Black PNBrilliant Black BN
E155-Brown HT Chocolate brown HT
Eicosapentaenoic acid
Epichlorhydrin
Erythritol
Escin
Ethanol
Ethyl 4-Fluorobenzenesulfinate
Ethyl cellulose
Ethyl DL-mandelate
Ethyl hydroxyethyl cellulose
Ethyl maltol
Ethyl Methanesulfonate
Ethyl N-benzoyl-L-argininate hydrochloride
Ethyl palmitate
Ethylene glycol dimethyl ether
Ethylenediamine
Ethylparaben
Etofibrate
Eugenol
Farnesol
FD&C Red #40, Allura red AC
FD&C Yellow 5, Tartrazine
Fenoprofen
Ferric Hydroxide
Fluorescein
Fluoxymesterone
Formaldehyde

Formic acid
Fumaric acid
Furil
Galactaric acid
Gallic acid
g-Butyrolactone
Gelatin
Gemfibrozil
Geranyl Acetate
Glucono-delta-lactone
Glucose
Glycerin
Glycerol Stearate
Glyceryl Stearate
Glycine
Glycocholic acid
Glycyrrhizic acid ammonium salt
Glyoxylic acid monohydrate
Guaiacol
Guanidine carbonate
Guanidine Hydrochloride
Guanine
Gum arabic
Gum Tragacanth
Heparin sodium salt
Heptafluorobutyric Acid
Hexachloroplatinic (IV) acid hexahydrate
Hexanol
Hexylresorcin
Hordenine Sulfate
Hyaluronic Acid
Hydroxy propyl distarch phosphate
Hydroxyethylcellulose
Hydroxymethylcellulose
Hydroxypropyl Cellulose
Hypoxanthine
Chlorogenic acid
Chlorophyll
Chloroquine Diphosphate
Cholesteryl Acetate
Cholesteryl Chloride
Cholic acid
Choline Bitartrate
Choline chloride
Chondroitin 4-sulfate Sodium
Chondroitin 6-sulfate Sodium from Shark Cartilage
Chondroitin Sulfate
Imidazole acrylic acid
Imidazolidinyl urea
Indoprofen
Iron (III) Hydroxide
Iron (III) Oxide Hydroxide
Iron (III) oxide monohydrate, yellow
Iron (III) Oxide, Hematite
Iron (III) Oxide, Red
Iron Citrate

Iron Oxide
Iron Oxide Black
Isatin
Isoamyl alcohol
Isoamyl ether
Isomenthol
Isopropyl 4-Fluorobenzenesulfinate
Isopropylamine
Isoquinoline
Kaolin, Aluminium silicate
L(-)-Carnitine
L(-)-Sorbitol
L(+)-Glutamic acid .HCl
L(+)-Tartaric acid
Lanolin Wax
L-Arabitol
L-Arginine L-Aspartate
Lavandulol
L-Carnosine
L-Cysteine ethyl ester hydrochloride
L-Djenkolic Acid
Lecithin
Lepidine
L-Galactose
L-Glucose
L-Glutamic Acid
L-Glutamine
Limonene
Linoleic Acid
Liothyronine sodium
L-Lactic acid
L-Lactic acid, Sodium salt
L-Malic acid
L-Menthol
L-Norleucine
L-Norvaline
L-Rhamnose
L-Ribose
L-Sorbitol
Magnesium carbonate
Magnesium Citrate
Magnesium hydroxide carbonate
Magnesium chloride
Magnesium lactate
Magnesium oxide
Magnesium Stearate
Maleimide
Maltol
Mannitol
Mannose
Marfey's reagent
Melatonin
Meloxicam Sodim Salt Hydrate
Menthol
meso-2,3-Dimercaptosuccinic acid
Methocel E3

Methyl alcohol
Methyl Anthranilate
Methyl Cellulose
Methyl docosanoate
Methyl elaidate
Methyl Furoate
Methyl Linoleate
Methyl Oleate
Methyl ricinoleate
Methyl stearate
Methylcaffeine
Methyltetracosanoate
Methyltrienolone
Michler's Ketone
Miranol MTH
Mizoribine
Mofetil Mycophenolate
Mono and Diglycerides of fatty acids
Monopotassium Phosphate
Monopotassium tartrate
Monosodium Citrate Anhydrous
Monosodium phosphate
Monosodium tartrate
Morpholine
Myrcene
Myristic Acid Methylester
N-(1,2,3,4-tetrahydro-1-naphthyl)-ethylendiamine
N-(2-Aminoethyl)-1,2,3,4-tetrahydro-1-naphthamid
N-(4-Aminobenzoyl)-L-glutamic acid
N,N,N',N'-Tetramethyl-p-phenylenediamine dihydrochloride
N,N-Diethylacetamide
N,N'-Dimethylpiperazine
N-Acetyl-1-aspartyl glutamic acid
N-Acetyl-DL-Tryptophan
N-Acetylglycine
N-Acetyl-L-valine
Nandrolone phenylpropionate
n-Benzyloxycarbonyl-L-glutaminyglycine
N-Cyclohexyl-2-benzothiazolesulfenamide
Neohesperidine dihydrochalcone
Neostigmine methyl sulfate
Neryl Acetate
N-Formylmorpholine
Nitromethane
N-Methyl-D-glucamine
N-Methyl-N-(trimethylsilyl)trifluoroacetamide
N-Methylpropylamine
Noscapine hydrochloride
N-t-BOC-L-norvaline, DCHA salt
N-t-BOC-L-proline
N-t-BOC-L-valine
Oil of Peppermint, Menthol
Oil, Almond oil sweet
Oil, Caraway oil
Oil, Cinnamon Oil
Oil, Clove oil

Oil, Lavender Oil
Oil, Olive Oil
Oil, Soybean Oil
Oil, Sunflower oil
Oil, Wheat Germ Oil
Oleic acid
Orcinol
Orlistat
Orthoboric Acid
Orthophenyl phenol
Oxalate Iron dihydrate
Oxalic Acid
Oxfloxacin
o-Xylene
Palmitic acid
Palmitoleic acid
Pamam dendrimer, generation 1.5
Pamam dendrimer, generation 2
Pamam dendrimer, generation 5
p-Anisidine
Papain von Carica papaya
Papaverine
Paraffin Oil
p-Benzoquinone
Pectine
PEG-200 Glyceryl Tallowate
Pentaerythritol
Pentasodium triphosphate
Perfluorotributylamine
Phenylcyclohexane
Phytol
Picein
Picolinic acid
Pinacol
Pinene
Piperazine dihydrochloride
Piperidine
Poly(propylene glycol) 2000
Polyethylene Oxide 5000000
Polysorbate 20
Polyvidone 25, Polyvinylpolypyrrolidone
Ponceau 4R
Ponceau SX
Potassium Acetate
Potassium carbonate
Potassium carbonate, sesquihydrate
Potassium Ferrocyanide
Potassium hydrogen carbonate
Potassium lactate
Potassium nitrite
Potassium oxalate monohydrate
Potassium Sorbate
Powdered cellulose
Propyl Gallate
Propylene Glycol
Propylparaben

NICODOM IR Drugs, 1463 spectra

Copyright © NICODOM 2006-2012

NICODOM Ltd., Hlavni 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

ftir@ftir.cz, ftirspectra@gmail.com

List of Spectra Page 25 from 53

Proteinase K
Pulegon
Pyruvic acid
Quillaia extract
Quillaja Saponin
Quinhydrone
Quinine
Quinine Hydrochloride Dihydrate
Quinoline
Quinoline Yellow Lake
Raffinose
Red 10B, C.I. Food Red 12
Red 2G, Acid Red 1
Resazurine
Resorcinol
Rhodizonic acid
Riboflavin-5'-Phosphate
Riboflavine-5'-monophosphate, Na salt
Ricinoleic Acid
Ricinoleic Acid Methyleneester
Rosmarinate
Rosmarinic acid
Salbutamol
Sanguinarine Nitrate
Scopolamine hydrobromide
Seleno-DL-methionine
Silicagel
Silver Sulphadiazine
Sodium Acetate
Sodium benzoate
Sodium butyrate
Sodium carbonate
Sodium Cetylsulfate
Sodium cyclamate
Sodium formate
Sodium Gluconate
Sodium glutarate
Sodium glycocholate
Sodium Hyaluronate
Sodium Laurylsulfate
Sodium L-Malate
Sodium nitrate
Sodium nitrite
Sodium Octanoate
Sodium oleate
Sodium potassium tartrate
Sodium Tetraborate
Sodium-L-Glutamate
Sorbic Acid
Sorbitan monolaurate
Sorbitan monooleate
Sorbitan tristearate
Sotalol hydrochloride
Sparteine sulfate
Stachyose Tetrahydrate
Starch

Stearamide
Stearylstearate
Succinamide
Succinic acid
Succinimide
Sulfacetamide Sodium Salt Hydrate
Sulfaguanidine monohydrate
Sulfate Ferrous
Sunset yellow FCF
Suxamethoniumchloride
Talc
Tetrahydrofuran
Tetrasodium Edetate
Tetryzoline Hydrochloride
Thaumatococin
Thimerosal
Thiophene
Thujone
Thymidine
Tiapride hydrochloride
Titanium Dioxide
Tocopherol
Toluene
trans,trans-Farnesol
Trans-4-Hydroxy-D-proline
trans-Cinnamaldehyde
trans-Nerolidol
Trehalose dihydrate
Triammonium citrate
Tricalcium Citrate Tetrahydrate
Tricalcium phosphate
Tricaprin
Tricine
Triethyl citrate
Trifluorothymidin
Tripotassium Citrate
Tripotassium phosphate
Trisodium Citrate Dihydrate
Trisodium Citrate Dihydrate
Trisodium phosphate
Trizma Acetate
Tropic acid
Tween 40
Tween 85
Tylosin, [R-(R*,R*)]-2,3-dihydroxybutanedioate (salt)
Uracil
Urea
Uric Acid
Urotropine
Valencene
Vanillin
Vinylpyrrolidone-Vinylacetate Copolymer
Wax, Yellow Bees Wax
Xanthan gum
Xanthine
Xylazine

Xylitol
Xylometazoline Impurity A
Yeast, Furcelleran
Zein
Zinc lactate
Zinc salicylate
Zinc Stearate
Dihydroisoandrosterone
DL-Amphetamine sulfate

List of Spectra - CAS Numbers

10016-20-3
100208-62-6
100239-44-9
100-51-6
100-92-5
10094-58-3
100-97-0
100986-85-4
10098-89-2
10124-31-9
10161-33-8
10161-34-9
1020-31-1
10236-47-2
10238-21-8
10262-69-8
102-76-1
102767-28-2
10323-20-3
103577-45-3
103628-48-4
103890-78-4
10389-72-7
103-90-2
10417-94-4
10418-03-8
104344-23-2
104-46-1
104-54-1
1045-69-8
1047-25-2
104-94-9
104-98-3
10504-35-5
10540-29-1
105-87-3
106-22-9
106-23-0
106266-06-2
106-28-5
106-34-3
106-51-4
106-58-1

1066-33-7
106-89-8
107-15-3
1072-67-9
107-35-7
107724-20-9
107868-30-4
107-95-9
1082-88-8
1083-09-6
108-46-3
108-80-5
1088-11-5
108-88-3
108-91-8
108-93-0
1095-90-5
109-99-9
110-02-1
110119-84-1, 120011-70-3
110-14-5
110-15-6
110-17-8
110-44-1
110-63-4
110-71-4
110-89-4
110-91-8
11099-07-3
111-01-3
11113-50-1
111-27-3
11138-66-2
111470-99-6
111-87-5
1119-34-2
111974-69-7
111974-72-2
1120-01-0
112529-15-4
112-61-8
112-62-9
112-63-0
112-73-2
112-80-1
112809-51-5
112-92-5
1132-61-2
113299-38-0
113299-40-4
113-52-0
113-53-1
113-59-7
113-92-8
1143-38-0
114-49-8

115007-34-6
1152-76-7
115-37-7
115-38-8
1154-12-7
115-43-5
115436-72-1
115-44-6
115-63-9
115-77-5
115-95-7
1172-18-5
1176-08-5
1177-87-3
117-81-7
1183849-42-4
1185-57-5
118-55-8
118-71-8
118-76-3
119-65-3
1197-18-8
1197-21-3
120011-70-3
120014-06-4
120373-24-2
120-47-8
120511-73-1
120-51-4
120-72-9
121-00-6
121-33-5
121-57-3
121-79-9
122-09-8
122-11-2
1225-55-4
122-57-6
1229-29-4
123-30-8
123-35-3
123-51-3
123-56-8
123-77-3
123-91-1
123-94-4
124-04-9
124-10-7
124-26-5
124529-02-8
124750-99-8
1248-42-6
124-90-3
125-02-0
125-33-7
1253-84-5

125-40-6
125-52-0
1255-49-8
125-60-0
125-64-4
125-71-3
125-72-4
125-84-8
1260-17-9
12629-01-5
1264-72-8
126-52-3
12678-07-8
126-91-0
127-08-2
127-09-3
127-17-3
127-35-5
127-48-0
128-13-2
128-37-0
128794-94-5
129101-54-8
129-17-9, 20262-76-4
129-20-4
129318-43-0
129-51-1
129722-12-9
129722-25-4
130018-87-0
130-40-5
1305-62-0
1305-78-8
130-89-2
1309-33-7
1309-37-1
1309-48-4
130-95-0
131-01-1
13103-34-9
13138-48-2
1317-60-8
1317-61-9
1321-14-8
132-17-2
1321-73-9
132-20-7
132539-06-1
132679-61-9
1330-43-4
1332-37-2
1332-58-7
133-37-9
1338-39-2
1338-43-8
134-20-3

1343-98-2
134523-03-8
13463-67-7
134-71-4
134-72-5
134-80-5
13492-01-8
135-07-9
13521-83-0
135-23-9
136434-34-9
136-47-0
136-77-6
137-08-6
13710-19-5
13717-00-5
13734-41-3
137-58-6
137-66-6
137862-53-4
138-09-0
138-15-8
138-32-9
13838-08-9
138402-11-6
138729-47-2
138786-67-1
138-86-3
139-05-9
139-33-3
13943-58-3
14007-64-8
140-11-4
1405-20-5
1405-41-0
1405-54-5
1405-87-4
1405-89-6
1406-18-4
1406-65-1
141-12-8
141-22-0
141-24-2
141433-60-5
14152-97-7
141-53-7
1422-07-7
142227-49-4
142227-51-8
1424-00-6
142-47-2
142-64-3
143-19-1
14371-10-9
143-71-5
143-98-6

144-02-5
144-14-9
144-23-0
144-33-2
144-55-8
144598-75-4
144-62-7
144701-48-4
144-82-1
144-83-2
145-13-1
14521-96-1
145-50-6
146062-45-5
146-17-8
146-22-5
146-56-5
1468-83-3
147098-20-2
147-20-6
147-24-0
147-71-7
1477-39-0
147-85-3
14807-96-6
1480-96-2
148-72-1
148-79-8
14901-07-6
149-15-5
149-32-6
149-91-7
150-13-0
1502-95-0
150-30-1
151-06-4
151-21-3
151-73-5
151767-02-1
1518-16-7
1518-86-1
15205-66-0
152-11-4
152459-95-5
15262-86-9
15307-79-6
15307-81-0
153-18-4
15356-70-4
153-94-6
15402-79-6
154-23-4
154361-50-9
154-41-6
15507-76-3
155141-29-0

155-41-9
15572-79-9
15588-95-1
156-34-3
156-54-7
15676-16-1
15687-27-1
157-06-2
15761-39-4
15790-07-5
15879-93-3
16039-53-5
16051-77-7
1622-61-3
1622-62-4
162359-56-0
16283-36-6
1637-73-6
1639-60-7
16423-68-0
16698-26-3
16773-42-5
1707-14-8
171596-29-5
174264-50-7
1750-12-5
17560-51-9
17629-30-0
1783-96-6
1786-81-8
180200-66-2
1808-12-4
18109-81-4
1812-30-2
1820-81-1
182507-22-8
1837-57-6
18472-51-0
18497-13-7
18559-94-9
18624-44-7
186286-86-8
1867-66-9
18917-93-6
1898-66-4
18996-35-5
19-11-4
191217-81-9
19237-84-4
1928-81-0
1934-21-0
1937-62-8
1944-12-3
1972-08-3
19794-93-5
1984-06-1

200630-46-2
20123-80-2
202009-64-1
20236-82-2
20283-92-5
20344-49-4
20636-41-3
20702-77-6
20830-75-5
2088-71-3
2107-76-8
21187-98-4
21500-98-1
2152-34-3
215929-30-9
21602-66-4
217797-14-3
21829-25-4
21898-19-1
2192-20-3
219861-08-2
22004-32-6
2201-15-2
2201-17-4
2201-19-6
2201-39-0
2201-40-3
22071-15-4
22144-60-1
2216-51-5
2218-68-0
22189-31-7
22199-08-2
22204-53-1
22232-71-9
22316-47-8
22839-47-0
23092-17-3
23141-83-5
2315-02-8
2321-07-5
23283-97-8
2338-05-8
2353-45-9
2363-59-9
23693-18-7
23828-92-4
2391-03-9
24143-17-7
24259-59-4
24359-81-7
2437-95-8
2438-10-0
2442-49-1
2444-46-4
2446-23-3

24589-78-4
24967-93-9
25013-16-5
25086-89-9
2519-30-4
25249-54-1
25301-02-4
25322-68-3
25322-69-4
25332-39-2
25378-27-2
2578-28-1
25812-30-0
25917-35-5
25956-17-6
2595-98-4
259731-40-3
2611-82-7
2645-08-1#
26570-10-5
26658-19-5
26807-65-8
26839-76-9
26921-17-5
26937-01-9
27203-92-5
27565-41-9
2778-96-3
2783-94-0
2784-73-8
2801-68-5
28179-44-4
28721-07-5
28911-01-5
2893-78-9
2894-67-9
2898-12-6
28981-97-7
29094-61-9
29122-68-7
2955-38-6
298-14-6
298-46-4
298-57-7
298-59-9
299-39-8
300-06-4
30123-17-2
30165-97-0
302-17-0
302-40-9
302-72-7
302-84-1
302-95-4
303-42-4
303-49-1

304-20-1
304-21-2
304-55-2
304-59-6
305-84-0
3087-16-9
309-36-4
3094-09-5#
3106-85-2
31121-93-4
311-89-7
312-84-5
3130-87-8
314-19-2
315-37-7
31613-93-1
31637-97-5
317-34-0
3184-13-2
31842-01-0
31879-05-7
31884-77-2
318-98-9
3196-73-4
319-78-8
3198-07-0
321-30-2
32221-81-1
327-57-1
327-97-9
328-38-1
32986-56-4
33005-95-7
33008-06-9
331-39-5
33286-22-5
33342-05-1
33375-06-3
3337-71-1
33403-97-3
33522-95-1
3368-04-5#
338-69-2
3398-22-9
340-57-8
34381-68-5
34523-28-9
34552-83-5
34552-84-6
345-78-8
3458-28-4
3458-72-8
3466-75-9
349-46-2
3505-38-2
35323-91-2

3546-41-6
35543-24-9
35602-69-8
3567-12-2
3567-66-6
3567-69-9, 53026-69-9
357-07-3
3572-80-3
359-83-1
360-70-3
36105-20-1
361-37-5
3614-69-5
3615-37-0
3615-41-6
3615-44-9
36322-90-4
36330-85-5
3671-99-6
36791-04-5
368-16-1
3704-09-4#
37191-69-8
3734-67-6
373-49-9
37353-59-6
375-22-4
378-44-9
379-79-3
3810-80-8
3858-89-7
39236-46-9
39409-82-0
39450-01-6
3963-95-9
40077-57-4
404-82-0
404-86-4
41100-52-1
41354-29-4
41444-62-6
4180-23-8
41989-07-5
42017-89-0
4221-98-1
42417-44-7
42542-10-9
426-13-1
4268-36-4
4271-30-1#
43033-72-3
43061-13-8
43200-80-2
4330-99-8
434-07-1
434-22-0

4354-76-1
4358-88-7
4360-12-7
437-74-1
438-41-5
439-14-5
4394-85-8
440-17-5
443-48-1
443-72-1
44603-18-1
446-72-0
4548-53-2
4553-89-3
4602-84-0
461-58-5
461-89-2
4630-07-3#
465-65-6
466-99-9
4685-14-7
4691-65-0
4697-36-3
470-82-6
471-15-8
474-25-9
4752-86-7
475-31-0
476-32-4
4764-17-4
478-84-2
479-18-5
479-92-5
480-16-0
481-06-1
4812-40-2
481-29-8
481-30-1
481-37-8
481-72-1
482-36-0
484-20-8
485-49-4
486-66-8
487-54-7
487-93-4
488-81-3
491-35-0
492-62-6
492-94-4
4940-11-8
49562-28-9
496-67-3
497-19-8
498-16-8
498-59-9

499-75-2
50-00-0
50-01-1
5001-51-4
50-02-2
50-03-3
50-04-4
50-06-6
50-10-2
50-11-3
50-12-4
501-36-0
50-14-6
50-23-7
50-24-8
50-27-1
502-85-2
50-33-9
50-34-0
50-36-2
50370-12-2
50-37-3
504-15-4
504-17-6
50-41-9
50-42-0
50-48-6
50-52-2
5053-08-7#
50-54-4
50-55-5
50-58-8
50-59-9
50-63-5
506-87-6
50-69-1
50-70-4
507-70-0
50-78-2
50-89-5
50924-49-7
509-60-4
509-67-1
50978-10-4
50-98-6
50-99-7
51012-32-9
51012-33-0
5104-49-4
51-06-9
51146-56-6
512-04-9
51-21-8
51-24-1
512-69-6
51274-00-1

51333-22-3
51-35-4
51-42-3
51-43-4
5144-52-5
51-45-6
51-55-8
51-56-9
515-69-5
51-57-0
51-60-5
516-06-3
51-61-6
5175-83-7
51799-32-7
518-28-5
518-82-1
519-09-5
519-37-9
520-33-2
520-34-3
520-52-5
520-53-6
521-12-0
521-18-6
521-35-7
521-78-8
522-48-5
52-26-6
52-28-8
52315-92-1
523-44-4
52417-22-8
524-36-7
52-43-7
52-49-3
526-94-3
526-99-8
527-07-1
52740-56-4
52-86-8
52-89-1
52-90-4
529-64-6
530-14-3
530-26-7
53124-00-8
53152-21-9
531-75-9
532-03-6
53-21-4
532-32-1
533-67-5
53-39-4
53-41-8
53-42-9

53-43-0
5345-01-07
535-34-2
5355-16-8
53-60-1
536-33-4
53639-82-0, 6856-31-1
537-15-5
53716-49-7
53746-45-5
537-55-3
538-02-3
53-86-1
5392-40-5
5394-83-2
53956-04-0
53994-73-3
54063-52-4
54063-53-5
54-11-5
541-15-1
54-12-6
541-59-3
54-21-7
54-31-9
543-24-8
544-01-4
54527-84-3
54-64-8
54-71-7
548-00-5
54827-17-7
54-85-3
548-73-2
5490-27-7
549-18-8
55-06-1
550-99-2
55142-85-3
55268-75-2
55-31-2
55-48-1
555-30-6
55-56-1
55589-62-3
556-02-5
556-03-6
557-04-0
557-05-1
55904-02-4
561-20-6
561-27-3
56-17-7
561-78-4
562-10-7
56238-63-2

56-29-1
5630-53-5
56392-17-7
563-96-2
56-40-6
56-41-7
564-25-0
56-45-1
56512-49-3
56-53-1
56536-96-0
56-54-2
5655-61-8
56-75-7
56-81-5
56-84-8
56-85-9
56-86-0
56-89-3
56-92-8
56-95-1
5704-04-1#
57-10-3
57109-90-7
571-22-2
57-13-6
5721-91-5
57-27-2
57-29-4
57-30-7
57-33-0
57-41-0
57-42-1
5743-26-0
57432-61-8
57-44-3
57-48-7
57-50-1
57-53-4
57-55-6
576-26-1
57-66-9
57-68-1
577-71-9
57-83-0
57-85-2
5786-21-0
57-87-4
57-88-5
5794-04-7
5794-13-8
58-08-2
58-15-1
58-18-4
581-88-4
58-20-8

58-22-0
58-25-3
58-27-5
58-28-6
58-33-3
58-38-8
584-08-7
58-46-8
5850-44-2
58-55-9
58-56-0
586-11-8#
586-60-7
587-23-5
58-74-2
5875-06-9
58-85-5
58-86-6
58-93-5
58-94-6
58-95-7
589-92-4
590-00-1
59-02-9
590-63-6
59-23-4
59277-89-3
59-30-3
59-33-6
5934-29-2
5934-55-4
593-85-1
59-43-8
59-46-1
5950-69-6
59-51-8
59-52-9
59-63-2
5965-13-9
59-66-5
59-67-6
59729-32-7
59804-37-4
5987-82-6
59-92-7
5996-14-5
59-97-2
60046-25-5
60124-79-0
60-13-9
60-13-9
60-18-4
602-29-9
60-29-7
603-00-9
60-32-2

60-33-3
603-79-2
604-35-3
60-56-0
60-80-0
608-66-2
60-93-5
609-36-9
610-25-3
6106-24-7
611-13-2
61-12-1
61177-45-5
6119-47-7
6119-70-6
61-25-6
6132-04-3
61336-70-7
6138-23-4
6151-40-2
6164-47-2
616-68-2
61-68-7
616-91-1
61-72-3
617-45-8
617-48-1
61-76-7
61-78-9
61-80-3
61869-08-7
618-85-9
619-04-5
61-90-5
6190-55-2
619-08-9
620-45-1
6209-17-2
621-71-6
6217-54-5
622-64-0
62-33-9
62-44-2
62-50-0
62571-86-2
62-57-7
625853-76-1
627-35-0
62758-13-8
6276-54-6
627-75-8
6284-40-8
628-63-7
628-97-7
62-90-8
63-12-7

631-61-8
632-20-2
632-51-9
63409-48-3
63-42-3
634-74-2
635-65-4
6359-98-4
636-54-4
636-61-3
63-68-3
637-01-4
637-23-0
6381-79-9
6381-92-6
63-84-3
63868-62-2
63-91-2
63-92-3
6393-42-6
63-98-9
6402-23-9
64024-15-3
64-02-8
64057-70-1
640-68-6
64-17-5
64-18-6
64-19-7
64211-46-7
64-43-7
644-62-2
64485-93-4
64-55-1
6459-59-2
64-75-5
64-77-7
6487-48-5
6493-05-6
6493-06-7#
6493-07-8#
65-19-0
65-22-5
65-23-6
65277-42-1
65-31-6
65-45-2
65-46-3
657-27-2
65-85-0
6591-63-5
6600-40-4
6610-42-0
66-22-8
6629-29-4
66335-38-4

665-66-7
66-76-2
66-79-5
66-84-2
67-03-8
6706-59-8
67-20-9
673-06-3
67-45-8
67-48-1
67-52-7
67-56-1
67-64-1
67861-96-5
67874-43-5
67-97-0
68-04-2
6805-41-0
68153-76-4
68-19-9
68-35-9
68424-04-4
6850-28-8
68526-83-0
6856-31-1
685-91-6
68684-55-9
6870-94-6
68-89-3
68917-82-8
68-94-0
6898-97-1
68990-67-0
69-05-6
69-09-0
69-52-3
69-53-4
69587-11-7
696-07-1
69-65-8
696-59-3
69-72-7
69-89-6
69-93-2
70-00-8
7004-03-7
70-07-5
7021-09-2#
70-26-8
70-47-3
7054-25-3
7069-42-3
7081-40-5
70851-08-0
7085-55-4
71-00-1

71125-38-7
71125-39-8
71130-06-8
71-27-2
71-30-7
71320-77-9
71396-29-7
7143-09-1
71463-34-8
71-63-6
716-61-0
71-68-1
7179-49-9
71-81-8
71963-77-4
7207-92-3
7212-44-4
72-14-0
72-17-3
72-18-4
72-19-5
7232-21-5
723-46-6
7235-40-7
72-44-6
72-63-9
72956-09-3
73-22-3
73-24-5
73-31-4
73-32-5
73-40-5
7361-61-7
73-78-9
738-70-5
74578-69-1
74610-55-2
74-79-3
7491-74-9
75-05-8
7516-48-5
75-31-0
75330-75-5
7541-49-3
75-52-5
7558-72-7
7558-73-8
7558-79-4
7558-80-7
75708-92-8
7585-39-9
7601-54-9
76-09-5
76095-16-4
76-25-5
7631-99-4

7632-00-0
76420-72-9
76-42-6
76-43-7
7643-75-6
76-57-3
76-58-4
766-09-6
76-68-6
76-73-3
76-75-5
7675-83-4
76824-35-6
76-90-4
7695-91-2
77-02-1
771-62-0
7720-78-7
77-21-4
77-26-9
77-27-0
77-28-1
77337-73-6
77-36-1
77-39-4
7757-93-9
7758-09-0
7758-11-4
7758-23-8
7758-29-4
7758-87-4
77-65-6
77-67-8
7778-18-9
7778-53-2
7778-77-0
7786-30-3
77-92-9
77-93-0
78246-49-8
78-44-4
78628-80-5
79-33-4
79350-37-1
79-81-2
79902-63-9
8000-28-0
8000-34-8
8000-95-1
8001-21-6
8001-22-7
8001-25-0
8001-79-4
8002-43-5
8004-92-0, 100208-62-6
8006-40-4

8006-54-0
8006-90-4
8006-95-9
8007-69-0
8008-60-4
8012-95-1
8015-86-9
8015-91-6
80210-62-4
80214-83-1
80223-99-0, 106463-17-6
8028-89-5
804-63-7
8047-15-2
8048-52-0
80-50-2
80-56-8
80-97-7
81-13-0
81-23-2
81-25-4
813-94-5
814-80-2
81-81-2
82186-77-4
82248-59-7
82419-36-1
82-58-6
82586-55-8
826-36-8
82640-04-8, 84449-90-1
827-52-1
82801-81-8
82-98-4
832-92-8
83380-47-6
83-46-5
83-48-7
83-67-0
83-73-8
83881-52-1
83-88-5
84057-84-1
84485-00-7
84625-61-6
846-48-0
846-49-1
846-50-4
84-66-2
84-74-2
84-80-0
848-75-9
85-18-7
853-34-9
853-39-4
855-19-6

85721-33-1
85-73-4
85-79-0
863-57-0
86386-73-4
86-42-0
866-84-2
867-56-1
868-14-4
868-59-7
86-98-6
870-24-6
870-93-9
87239-81-4
87269-97-4
87-32-1
87333-19-5
87-67-2
87-69-4
87-78-5
87-79-6
87-81-0
87-89-8
879-39-0
87-99-0
88-02-8#
88150-47-4
88-27-7
894-71-3
89-57-6
89-65-6
89-78-1
89-82-7
89-83-8
9000-01-5
9000-07-1
9000-11-7
9000-65-1
9000-69-5
9000-70-8
9001-73-4
9002-18-0
9004-34-6
9004-54-0
9004-57-3
9004-58-4
9004-61-9
9004-64-2
9004-65-3
9004-67-5
9004-70-0
90-05-1
9005-25-8
9005-46-3
9005-49-6
9005-64-5

www.ir-spectra.com

9005-66-7
9005-70-3
9006-52-4
9006-65-9
9010-66-6
9016-00-6
90-33-5
90357-06-5
9037-22-3
9041-08-1#
90-43-7
9067-32-7
90-80-2
9082-07-9#
90-82-4
90-94-8
910-31-6
91-22-5
912-60-7
91374-20-8
91446-96-7
91-56-5
915-67-3
91-75-8
91-81-6
92-13-7
921-53-9
921-60-8
92-52-4
929-77-1
93107-08-5
93-14-1
932-53-6
93376-66-0
93390-81-9
93-40-3
93413-69-5
93-72-1
94-09-7
94-13-3
94-20-2
94-26-8
95-25-0
95-33-0
95-47-6
95713-52-3
95716-98-6
959-24-0
96-27-5
963-39-3
96-48-0
965-93-5
96-81-1
96829-58-2
97-05-2
97-53-0

97-59-6
97-67-6
977-32-2
97-77-8
98418-47-4
98-84-0
98-92-0
98-98-6
99-06-9
990-73-8
99-26-3
99294-93-6
99300-78-4
99413-49-7
99-49-0
99614-02-5
996-31-6
99-76-3
999-81-5

www.ir-spectra.com