



IR Lubricants – List of Spectra

The infrared spectral library „IR Lubricants“ is a unique collection of 2018 FTIR digital spectra of lubricants and related materials.

1612 spectra were collected by the ATR technique. The spectra were collected using the FTIR Spectrometer Nexus 670™ (Thermo) and the Miracle™ single bounce ATR accessory (Pike Technologies) equipped by ZnSe crystal.

To emphasize the absorption bands of the additives especially by lubricants with low additive contents selected 469 samples were additionally collected in transmission mode in thick KBr cell (thickness 0.1mm). Totally absorbing bands of the base oil were cut off by the transmission spectra.

The spectral library features resolution of 4cm⁻¹ in the spectral range 600-4000 cm⁻¹ (400-4000 cm⁻¹ for the transmission spectra) measurement time 1 minute.

Typical information about samples include Name, Manufacturer, Category (lubricant type) and notes.

In some cases also chemical composition and CAS number is included (when available).

Complete sample information comes with the full version of the library, following list shows only part of the information.

This library is available as digital searchable library compatible with your spectroscopic software. Also it can be delivered with SEARCH software or in *.pdf format.

To download the list of spectra from this library, to download the free Nicodom IR Demo Library (50 spectra) compatible with your software, to check the format compatibility as well as for pricelist, ordering info and information about all NICODOM IR/NIR Libraries and other NICODOM products please visit our webpage.

IR Lubricants, 2081 FTIR spectra
<http://www.ir-spectra.com/>

List of Sample Types

Here are some examples of sample types included in this spectral library. Additional to many lubricant types typical automotive fluids (Antifreeze, Brake fluids, Fuels...) are also included as well as other industrial samples (Additives, Cleaners, Emulsifiers, Plasticizers, Surfactants...). Edible and Essential Oils are also part of this library. Spectra of used industrial oils show how infrared spectroscopy can be used for used oil analysis.

Additive
Anticorrosion Oil
Antifreeze
Aviation Oil
Axle Oil
Base Oil
Bearing Oil
Brake fluids
Circulation Oil
Cleaner
Compressor Oil
Edible Oils
Emulsifier
Essential Oils
Forming Oil
Fuel
Fuel Additive
Gear Oil
Glazing Agent
Grease
Grinding Oil
Hardening Oil
Heat Transfer Medium
Historical Lubricant
Hydraulic Oil
Industrial Oil
Machine Oil
Medicinal Oil
Metalworking Fluid
Motor Oil
Natural
Oil Additive
Oil for sliding guides
Plasticizer
Process Oil
Quenching Oil
Refrigeration Oil
Rolling Oil
SBP spirit
Servo Oil
Slideway Oil
Steam Cylinder Oil
Surfactant

IR Lubricants, 2081 FTIR spectra
<http://www.ir-spectra.com/>

Tractor Lubricant used
Transformer Oil
Turbine Oil
Used Oils

List of ATR Spectra

BP ENERGOL CM, Forming Oil, BP,
BP ENERSYN HT 460, Gear Oil, BP,
BP BARTRAN HV 46, Hydraulic Oil, BP,
BP BIOHYD 46, Hydraulic Oil, BP,
BP ENERGOL SHF-LT 15, Hydraulic Oil, BP,
BP MICROTREND 217 M, Metalworking Fluid, BP,
BP OLEX XK 6014, Metalworking Fluid, BP,
BP SEVORA 532 L, Metalworking Fluid, BP,
BP SEVORA 721M, Metalworking Fluid, BP,
BP ENERFLEX 74, Special Oil, BP,
BP ENERGOL WM 2, Special Oil, BP,
BP ENERPAR M002, Process Oil, BP,
BP DIELECTRIC 200T, Transformer Oil, BP,
BP ENERGOL THB 32, Turbine Oil, BP,
BP ENERGOL THB 46, Turbine Oil, BP,
BP TURBINOL 32, Turbine Oil, BP,
BP ENERGOL LPT 68, Compressor Oil, BP ,
BP ENERTHENE 2368, Special Oil, BP ,
CASTROL ICEMATIC SW 32, Compressor Oil, Castrol,
CASTROL EPX 90, Gear Oil, Castrol,
CASTROL TAF-X 75W90, Gear Oil, Castrol,
CASTROL OPTITEMP MT, Grease, Castrol,
CASTROL VISCOGEN 0, Grease, Castrol,
CASTROL LHM PLUS, Hydraulic Oil, Castrol,
CASTROL DC282, Metalworking Fluid, Castrol,
CASTROL B353, Motor Oil, Castrol,
CASTROL SYNTEC ACTIVE SAE 10W-30, Motor Oil, Castrol,
CASTROL SYNTEC BLEND SAE 20W-50, Motor Oil, Castrol,
CASTROL OPTIFLEX A460, Special Oil, Castrol,
CASTROL VISCOGEN KL 15, Special Oil, Castrol,
TUNNELWASH, Cleaner, Castrol,
ESSO RUST PROTECTION MZ 45, Anticorrosion Oil, Esso,
ESSO RUST-BAN 335, Anticorrosion Oil, Esso,
ESSO RUST-PEL 91 NAL/1996, Anticorrosion Oil, Esso,
ESSO RUST-PEL 91 NAL/1999, Anticorrosion Oil, Esso,
ESSO NUTO H32, Hydraulic Oil, Esso,
ESSO FLEXON 843, Machine Oil, Esso,
ESSO NUTO 460, Machine Oil, Esso,
ESSO FANOX AN 32, Metalworking Fluid, Esso,
ESSO KUTWELL 40, Metalworking Fluid, Esso,
ESSO LP 32, Metalworking Fluid, Esso,
ESSO EXTRA MOTOR OIL SAE 15W-40, Motor Oil, Esso,
ESSO ULTRA SAE 10W 40, Motor Oil, Esso,
ESSOLUBE XT 301 SAE 15W-40, Motor Oil, Esso,

IR Lubricants, 2081 FTIR spectra
<http://www.ir-spectra.com/>

IR Lubricants

ESSO ANTIFOAM EGL 43263, Oil Additive, Esso,
ESSO VALVES PROTECTION PROTECT PLUS, Oil Additive, Esso,
ESSO XD 80, Quenching Oil, Esso,
ESSO WALZOEL BM 200, Rolling Oil, Esso,
ESSO WALZOEL BM 71, Rolling Oil, Esso,
ESSO FLEXON 856, Special Oil, Esso,
ESSO PRIMOL 352, Special Oil, Esso,
ESSO UNIVOLT 524, Transformer Oil, Esso,
AEROSHELL FLUID 18, Aviation Oil, Shell,
AEROSHELL FLUID 31, Aviation Oil, Shell,
DONAX YB/1998, Brake fluid, Shell,
DONAX YB/1999, Brake fluid, Shell,
AVGAS 100LL, Fuel, Shell,
SHELL DONAX TA CHECK, Gear Oil, Shell,
SHELL SPIRAX MA SAE 80W, Gear Oil, Shell,
SHELL HELIX PLUS SAE 10W-40, Motor Oil, Shell,
SHELL HELIX SAE 15W-40, Motor Oil, Shell,
SHELL HELIX SUPER SAE 15W-40, Motor Oil, Shell,
SHELL SUPER 2T, Motor Oil, Shell,
SHELL DIESEL DEPRESSER SUPER CONCENTRATE, Oil Additive, Shell,
SHELL XHVI 5.2, Base Oil, Shell,
SHELL CASSIDA CHAIN OIL 150, Special Oil, Shell,
SHELL ONDINA P 68, Special Oil, Shell,
SHELL SINOL, Special Oil, Shell,
BRANOROL 32/10, Anticorrosion Oil, BRANOPAC,
CASTROL RUSTILO DWX 32, Anticorrosion Oil, Castrol,
DINITROL 3654, Anticorrosion Oil, DINOL,
DINITROL 4942, Anticorrosion Oil, DINOL,
ANTICORIT, Anticorrosion Oil, Fuchs,
ANTICORIT 6120/42E, Anticorrosion Oil, Fuchs,
ANTICORIT OHK 360, Anticorrosion Oil, Fuchs,
HOUGHTON GRIND 590, Anticorrosion Oil, Houghton Lubricor,
OLVIKOR 700, Anticorrosion Oil, Hungary,
CONTRAKOR A 40F, Anticorrosion Oil, Kluber,
CONTRAKOR FLUID H1, Anticorrosion Oil, Kluber,
KLUBERSYNTH, Anticorrosion Oil, Kluber,
LUSIN PROTECT G11, Anticorrosion Oil, Kluber,
LUSIN PROTECT G12, Anticorrosion Oil, Kluber,
LUSIN PROTECT G31, Anticorrosion Oil, Kluber,
LUSIN PROTECT G32, Anticorrosion Oil, Kluber,
MOBIL OIL TECREX, Anticorrosion Oil, Mobil,
ANTICORROSION FOR AVIATION/1998, Anticorrosion Oil,
ANTICORROSION FOR AVIATION/1999, Anticorrosion Oil,
ANTICORROSION FOR AVIATION/2002, Anticorrosion Oil,
ANTICORROSION HIGH PERFORMANCE, Anticorrosion Oil,
ANTICORROSION OIL /2005, Anticorrosion Oil,
ANTICORROSION SPECIAL, Anticorrosion Oil,
OK-40, Anticorrosion Oil,
WD-40, Anticorrosion Oil,
KONKOR 101, Anticorrosion Oil, Paramo,
KONKOR 103, Anticorrosion Oil, Paramo,
KONKOR 210, Anticorrosion Oil, Paramo,
KONKOR 268, Anticorrosion Oil, Paramo,

IR Lubricants, 2081 FTIR spectra

<http://www.ir-spectra.com/>

IR Lubricants

SHELL ENSIS OIL 30, Anticorrosion Oil, Shell,
PEKOL 80, Anticorrosion Oil, Slovnaft,
PEKOL VK, Anticorrosion Oil, Slovnaft,
ANTICORROSION OIL /1998, Anticorrosion Oil, SPANJAARD,
ANTICORROSION OIL /2002, Anticorrosion Oil, SPANJAARD,
DINITROL AV 8, Anticorrosion Oil, Sweden,
INKOR O-190, Anticorrosion Oil, Triga,
INKOR O-192, Anticorrosion Oil, Triga,
TECTYL 802 A, Anticorrosion Oil, Valvoline,
TECTYL 900, Anticorrosion Oil, Valvoline,
FRIDIOL EXTRA, Antifreeze, AGRIMEX,
COOLER ANTIFREEZE, Antifreeze, AVANTI,
EUROLUB ANTIFREEZE FOR BRAKES, Antifreeze, EUROLUB,
EUROLUB ANTIFREEZE FOR WIPER, Antifreeze, EUROLUB,
FROSTON C 2230, Antifreeze, MOTORON C7,
ANTIFREEZE FLUID, Antifreeze,
CHEMIX AKVALOR, Antifreeze,
MOTORCRAFT SUPER PLUS, Antifreeze,
SILICATE FREE LIFETIME FILL , Antifreeze,
OZ-40 LENA, Antifreeze, Russia,
TOSOL A40M, Antifreeze, Russia,
PARAFLU, Antifreeze, SELENIA,
LONG-LIFE COOLANT D, Antifreeze, Texaco,
KOELVLOEISTOF, Antifreeze, Valvoline,
FRIDEX ECO, Antifreeze, Velvana,
AEROLUB, Aviation Oil, Aerolub,
AERO 585B - A, Aviation Oil, Castrol,
AERO 585B - B, Aviation Oil, Castrol,
EXXON AVIATION OIL 100, Aviation Oil, Exxon,
EXXON AVIATION OIL ELITE 20W50, Aviation Oil, Exxon,
AEROSHELL 100, Aviation Oil, Shell,
AEROSHELL 15W50 MULTIGRADE, Aviation Oil, Shell,
AEROSHELL FLUID 12, Aviation Oil, Shell,
AEROSHELL FLUID 3, Aviation Oil, Shell,
TOTAL AERO D100, Aviation Oil, Total,
TOTAL AERO D80, Aviation Oil, Total,
OL-J32 PARAMOL , Bearing Oil, Paramo,
PARAMO OL-J10, Bearing Oil, Paramo,
PARAMO OL-J100, Bearing Oil, Paramo,
PARAMO OL-J22, Bearing Oil, Paramo,
PARAMO OL-J3, Bearing Oil, Paramo,
PARAMO OL-J32, Bearing Oil, Paramo,
PARAMO OL-J46, Bearing Oil, Paramo,
MADIT OL-32, Bearing Oil, Slovnaft,
MADIT OL-B4/2000, Bearing Oil, Slovnaft,
SYNTALAN, Brake fluid, AGRIMEX,
SYNTOL HD 205/1996, Brake fluid, CHEMOPETROL,
DOT 3, Brake fluid, Coyote,
DOT 4, Brake fluid, Coyote,
DOT 5 PLUS, Brake fluid, DENICOL,
SILICONE BRAKE FLUID 2-1143, Brake fluid, DOW CORNING,
HBF4, Brake fluid, Total,
SYNTOL HD 205 Y.1999, Brake fluid, Velvana,

IR Lubricants, 2081 FTIR spectra

<http://www.ir-spectra.com/>

SYNTOL HD 205/2001, Brake fluid, Velvana,
SYNTOL HD 260, Brake fluid, Velvana,
SYNTOL HD 265, Brake fluid, Velvana,
ROTO-INJECTFLUID, Compressor Oil, ATLAS COPCO,
EMKARATE RL22H, Compressor Oil, ICI,
COMPRESSOR OIL ISO 46, Compressor Oil, Kernite,
GARGOYLE ARCTIC OIL, Compressor Oil, Mobil,
GARGOYLE ARCTIC OIL 300, Compressor Oil, Mobil,
MOBIL RARUS 427, Compressor Oil, Mobil,
PARAMO OT-K12, Compressor Oil, Paramo,
PARAMO OT-K18, Compressor Oil, Paramo,
PARAMO OT-K28S, Compressor Oil, Paramo,
PARAMO OT-K8, Compressor Oil, Paramo,
CLAVUS OIL 15, Compressor Oil, Shell,
SHELL CLAVUS G 68, Compressor Oil, Shell,
MADIT OK-VC 150, Compressor Oil, Slovnaft,
ACACIA GUM, Emulsifier, 9000-01-5
CARBOPOL 940, Emulsifier, 9003-01-4
ETHOXYLATED ALCOHOLS, Emulsifier, 66455-14-9
HYDROMYRISTENOL N, Emulsifier,
PEG 7 GLYCERYL COCOATE CHECK, Emulsifier, 68201-46-7
POLOXAMER 105, Emulsifier, 9003-11-6
POLYSORBATE 20, Emulsifier, 9005-64-5
STEARYL ALCOHOL AND ESTERS MIXTURE, Emulsifier,
TRITON X-100, Emulsifier, 9002-93-1
TWEEN 60, POE (20) SORBITAN MONOSTEARATE, Emulsifier, 9005-67-8
TWEEN 61, POE (4) SORBITAN MONOPALMITATE, Emulsifier, 9005-67-8
PARAMO PARAMUL ERO-SB, Emulsifier, Paramo,
QUAKEROL, Emulsifier, QUAKER,
Cardamom Oil, Essential Oil, A.R.A.,
CHINESE CINNAMON OIL, Essential Oil, H. REYNAUD,
LAVENDER OIL, Essential Oil, H. REYNAUD,
ORANGE BIL, Essential Oil, H. REYNAUD,
YOYوبا BIL, Essential Oil, 61789-91-1
OIL OF SWEET ORANGE, Essential Oil, WORLEE-CHEMIE,
TRENNFIT, Forming Oil, FRANK,
TRENNFIT 3102, Forming Oil, FRANK,
FORMING OIL CHECK, Forming Oil,
FORMING OIL CHECK, Forming Oil,
FORMING-15W, Forming Oil, Techservis,
FORMING-4K, Forming Oil, Techservis,
FORMING-4KA, Forming Oil, Techservis,
FORMING-4W, Forming Oil, Techservis,
FORMING-5K, Forming Oil, Techservis,
FORMING-5KA, Forming Oil, Techservis,
FORMING-5KG, Forming Oil, Techservis,
FORMING-5KM, Forming Oil, Techservis,
FORMING-5W, Forming Oil, Techservis,
FORMING-6K, Forming Oil, Techservis,
FORMING-EKO 1, Forming Oil, Techservis,
FORMING-EKO W, Forming Oil, Techservis,
FORMING-KP EKO, Forming Oil, Techservis,
FORMING-SV, Forming Oil, Techservis,

IR Lubricants, 2081 FTIR spectra

<http://www.ir-spectra.com/>

IR Lubricants

VALVOLINE SLUSHING OIL A-51, Forming Oil, Valvoline,
CONCRETE FLUID SUPER 120, Forming Oil, WOERMANN,
CONCRETE SEPARATOR 90, Forming Oil, WOERMANN,
CONCRETE SEPARATOR B, Forming Oil, WOERMANN,
CONCRETE SEPARATOR CONCENTRATE 96, Forming Oil, WOERMANN,
CONCRETE SEPARATOR LF 101, Forming Oil, WOERMANN,
CONCRETE SEPARATOR LF U5, Forming Oil, WOERMANN,
CONCRETE SEPARATOR SV304, Forming Oil, WOERMANN,
PLANTOCLEAN 364, Forming Oil, WOERMANN,
FUEL OIL LIGHT, Fuel, CHEMOPETROL,
FUEL OIL, Fuel, Germany,
KEROSENE F34 WITHOUT ADDITIVE CI/LI FS4, Fuel,
KEROSENE PL34, Fuel,
FUEL OIL EXTRA L, Fuel, Paramo,
PETROLEUM (DIESEL) NM 22A, Fuel, Paramo,
FUEL OIL L, Fuel, Petrochema,
FUEL OIL L COLORED, Fuel, Petrochema,
GASOLINE FOR PYROLYSIS/1999, Fuel, Petrochema,
GASOLINE FOR PYROLYSIS/2001, Fuel, Petrochema,
LIGHT FUEL OIL, Fuel, Petrochema,
LIGHT FUEL OIL LOW S, Fuel, Petrochema,
PETROLEUM (DIESEL) NM-35, Fuel, Petrochema,
ZIBRO PLUS, Fuel, PVG,
UNLEADED PETROL 91 OCTAN/1995, Fuel, Slovnaft,
UNLEADED PETROL 91 OCTAN/1996, Fuel, Slovnaft,
UNLEADED PETROL 91 OCTAN/1997, Fuel, Slovnaft,
UNLEADED PETROL 91 OCTAN/1999, Fuel, Slovnaft,
UNLEADED PETROL 95 OCTAN CHECK, Fuel, Slovnaft,
UNLEADED PETROL 95 OCTAN CHECK, Fuel, Slovnaft,
UNLEADED PETROL 95 OCTAN/1995, Fuel, Slovnaft,
UNLEADED PETROL 95 OCTAN/1997, Fuel, Slovnaft,
UNLEADED PETROL 95 OCTAN/1999, Fuel, Slovnaft,
UNLEADED PETROL 98 OCTAN CHECK, Fuel, Slovnaft,
UNLEADED PETROL 98 OCTAN CHECK, Fuel, Slovnaft,
ADERCO 5000, Fuel Additive, Aderco,
KEROFLUID MIL AL-31, Fuel Additive, BASF,
KEROFLUID MIL AL-41, Fuel Additive, BASF,
KEROFLUX 3182, Fuel Additive, BASF,
KEROFLUX 3222, Fuel Additive, BASF,
KEROFLUX 3283, Fuel Additive, BASF,
KEROFLUX 5486, Fuel Additive, BASF,
KEROFLUX ES 3217, Fuel Additive, BASF,
KEROFLUX ES 3267, Fuel Additive, BASF,
KEROFLUX ES 3424, Fuel Additive, BASF,
BENADIT BA CONCETRATE, Fuel Additive, CHEMOPETROL,
CORROSION INHIBITOR DCI-4A, Fuel Additive, CHEMOPETROL,
HITEC, Fuel Additive, CHEMOPETROL,
BENADIT BA, Fuel Additive, CONTICHEM,
DIESEL FUEL CONDITIONER, Fuel Additive, GOLD EAGLE,
GASOLINE FUEL SYSTEM TREATMENT, Fuel Additive, INNOTECH,
RXP MASTER BLEND, Fuel Additive, INNOTECH,
SUPER DIESEL ADITIV, Fuel Additive, LANG,
DIESEL RUSS, Fuel Additive, Liqui Moly,

IR Lubricants, 2081 FTIR spectra

<http://www.ir-spectra.com/>

IR Lubricants

NALCO 5403, Fuel Additive, NALCO VISCO CHEMICALS,
NYCOSOL 37, Fuel Additive, NYCO,
ET, Fuel Additive, OEMETA,
TEXACO TEXAFLOW, Fuel Additive, Texaco,
ADDITIVE FOR GASOLINE, Fuel Additive, WURTH,
WYNN'S ANTIFREEZING DIESEL ADDITIVE, Fuel Additive, WYNN OIL,
ADDINOL ATF DIID, Gear Oil, Addinol,
ADDINOL GEAR CLP 150, Gear Oil, Addinol,
ADDINOL GEAR CLP 160, Gear Oil, Addinol,
ADNOC GX 85W140, Gear Oil, Adnoc-Fod,
AGIP ROTRA, Gear Oil, Agip,
AGIP ROTRA GL 80W-90, Gear Oil, Agip,
ARAL ATF 22, Gear Oil, Aral,
ARAL GEAR EP 85W-90, Gear Oil, Aral,
BOSS 80W-90, Gear Oil, BOSSHARDT CHEMIE AG,
BP GEAR OIL 80W-90, Gear Oil, BP,
BP TRACTRAN 9, Gear Oil, BP,
CASTROL MTX, Gear Oil, Castrol,
CASTROL TQ DEXRON III, Gear Oil, Castrol,
CAT GL 5 80W-90, Gear Oil, Catterpillar,
CAT TDTO 10W, Gear Oil, Catterpillar,
CAT TDTO SAE30, Gear Oil, Catterpillar,
DEAGEAR 85W-90, Gear Oil, DEA,
ESSO ATF D 21 611, Gear Oil, Esso,
TITAN RENEP, Gear Oil, Fuchs,
TOP BLEND 75W90, Gear Oil, Kernite,
TOP BLEND 85W140, Gear Oil, Kernite,
KLUBER SYNTHESO D 1000 EP, Gear Oil, Kluber,
MOBIL ATF 200, Gear Oil, Mobil,
MOBIL GEAR 626/1997, Gear Oil, Mobil,
MOBIL GEAR 626/2000, Gear Oil, Mobil,
MOBIL GEAR 627, Gear Oil, Mobil,
MOBIL GEAR 629, Gear Oil, Mobil,
MOBIL GEAR 630, Gear Oil, Mobil,
MOBIL GEAR 632, Gear Oil, Mobil,
MOBIL GEAR 634, Gear Oil, Mobil,
MOBIL GEAR 636, Gear Oil, Mobil,
MOBIL LUBE SHC 75W-90 LS, Gear Oil, Mobil,
MOBIL MOBILFLUID 424, Gear Oil, Mobil,
HYKOMOL K SAE 80W90 - 1997, Gear Oil, MOL,
HYKOMOL K SAE 80W90 - 1999, Gear Oil, MOL,
KOMAGMOL 150, Gear Oil, MOL,
OMV GEAR OIL MP 80W-85, Gear Oil, OMV,
OMV GEAR OIL SY 75W90, Gear Oil, OMV,
OMV UNIGEAR, Gear Oil, OMV,
OMV UNIGEAR S SAE 75W90, Gear Oil, OMV,
CLP 220, Gear Oil,
CLP 320, Gear Oil,
CLP 5 80W90, Gear Oil,
ESD-M2C186-A, Gear Oil,
ROTAX 80W90, Gear Oil,
SQM-2C9008-A, Gear Oil,
WSD-M2C200-B, Gear Oil,

IR Lubricants, 2081 FTIR spectra

<http://www.ir-spectra.com/>

GYROL 80W90, Gear Oil, Paramo,
 GYROL ATF DIII, Gear Oil, Paramo,
 GYROL PP80W90, Gear Oil, Paramo,
 GYROL SAE 80W-90, Gear Oil, Paramo,
 PARAMO GYROL PP80, Gear Oil, Paramo,
 PARAMO GYROL PP90, Gear Oil, Paramo,
 PARAMO OT-HP3 HYDYN, Gear Oil, Paramo,
 PARAMO PREMOL PP13, Gear Oil, Paramo,
 PARAMO PREMOL PP44, Gear Oil, Paramo,
 PARAMO PREMOL PP7, Gear Oil, Paramo,
 PENNZOIL ATF+2, Gear Oil, Pennzoil,
 PETROMIN ATF II-D, Gear Oil, PETROMIN OILS,
 PETROMIN GEAR OIL, Gear Oil, PETROMIN OILS,
 AXLE OIL 80W90, Gear Oil, Q8 Oil,
 DIFFERENTIAL OIL SAE 80W90 , Gear Oil, SAFETY-KLEEN,
 SHELL VITREA 32, Gear Oil, Shell,
 SHELL VITREA 46, Gear Oil, Shell,
 MADIT AUTOMATIC ATF/1996, Gear Oil, Slovnaft,
 MADIT AUTOMATIC ATF/1999, Gear Oil, Slovnaft,
 MADIT HYPOL PP SAE 80W/90 , Gear Oil, Slovnaft,
 MADIT OT-HP 32, Gear Oil, Slovnaft,
 MADIT PP 80 SAE 80W, Gear Oil, Slovnaft,
 MADIT PP 80/1997, Gear Oil, Slovnaft,
 MADIT PP 80/2000, Gear Oil, Slovnaft,
 MADIT PP 90/1999, Gear Oil, Slovnaft,
 MADIT PP 90/2001, Gear Oil, Slovnaft,
 SWEPCO 201, Gear Oil, SWEPCO,
 TEXACO TEXTRAN TDH, Gear Oil, Texaco,
 TRANSMISSION TM 80W90, Gear Oil, Total,
 ATF TYPE D/2002, Gear Oil, Valvoline,
 HP GL-4 80W90, Gear Oil, Valvoline,
 HP-GL-5 80W90, Gear Oil, Valvoline,
 VAT ATF TYPE D/1999, Gear Oil, VAT OIL,
 CARNAUBA WAX, Glazing Agent, 8015-86-9
 MICROCRYSTALLINE WAX CHECK, Glazing Agent, 63231-60-7
 PARAFFIN WAX, Glazing Agent, 9083-41-4
 TRIMYRISTIN , Glazing Agent, 555-45-3
 TRIPELARGONIN CHECK, Glazing Agent, 126-53-4
 BP ENERSYN MTX 320, Grease, BP,
 CASTROL SYNPLEX, Grease, Castrol,
 CASTROL OPTIMOL OPTITEMP TT1, Grease, Castrol,
 GREASE MS3, Grease, Castrol,
 LUBRICAST 852, Grease, CHEMTOOL,
 DOW CORNING PLASTIC MOLYKOTE PG 30L, Grease, DOW CORNING,
 DOW CORNING PLASTIC MOLYKOTE PG 65, Grease, DOW CORNING,
 DOW CORNING PLASTIC MOLYKOTE PG 75, Grease, DOW CORNING,
 MOLYKOTE DX , Grease, DOW CORNING,
 MOLYKOTE PG602, Grease, DOW CORNING,
 WSD-M1C226-A, Grease, Ford,
 WSD-M1C234-A, Grease, Ford,
 WSD-M1C253-A, Grease, Ford,
 GLEITMO 523, Grease, GLEITMO,
 GLEITMO 582, Grease, GLEITMO,

IR Lubricants, 2081 FTIR spectra

<http://www.ir-spectra.com/>

K-LAST SUPER, Grease, Kernite,
 KLUBER NBU 8 EP, Grease, Kluber,
 KLUBER TOPAS NB52, Grease, Kluber,
 MOBIL 22, Grease, Mobil,
 MOBIL GREASE HP 222, Grease, Mobil,
 MOBIL LUX EP 2, Grease, Mobil,
 MOBIL UNIVERSAL GREASE UF 2, Grease, Mobil,
 MOBILPLEX 44, Grease, Mobil,
 MOBILTAC 81, Grease, Mobil,
 MOBILTEMP 78, Grease, Mobil,
 MOBILTEMP SHC 32, Grease, Mobil,
 NYCO GREASE GN 37, Grease, NYCO,
 NYCO GREASE GN 43, Grease, NYCO,
 OKS 245 COPPER PASTE, Grease, OKS,
 OKS 350, Grease, OKS,
 AK2, Grease,
 AKM2, Grease,
 APIEZON N, Grease, 12704-92-6
 CV JOINT , Grease,
 GRAPHITE GREASE, Grease,
 KN2, Grease,
 LA2, Grease,
 MOGUL LV2-3, Grease,
 MOLYKA-G, Grease,
 NH2, Grease,
 SAM-1C9107-A, Grease,
 T-A4, Grease,
 T-N1, Grease,
 TPD 33-078-62, Grease,
 WSD-M1C228-A, Grease,
 WSD-M1C234-A2, Grease,
 COSMETIC GREASE, Grease, Petrochema,
 CIATIM 203, Grease, Russia,
 CIATIM 205, Grease, Russia,
 GREASE CIATIM 202, Grease, Russia,
 GREASE GOI-54P, Grease, Russia,
 GREASE M9, Grease, Russia,
 GREASE MS-70, Grease, Russia,
 GREASE NK-50, Grease, Russia,
 GREASE SU, Grease, Russia,
 OKB 122-7/1999, Grease, Russia,
 SVINCOL 01, Grease, Russia,
 AERO SHELL 22, Grease, Shell,
 SHELL ALAVANIA 3, Grease, Shell,
 SHELL RETINAX GREASE A, Grease, Shell,
 SHELL RETINAX GREASE AM, Grease, Shell,
 SWEPCO MOLY GREASE 101, Grease, SWEPCO,
 SWEPCO MOLY GREASE 105, Grease, SWEPCO,
 TEXACO MULTIFAK ALL PURPOSE EP 2, Grease, Texaco,
 TEXACO MULTIFAK EP 2, Grease, Texaco,
 TOTAL MULTIS, Grease, Total,
 TOTAL SPECIS CU, Grease, Total,
 ADDINOL HLP 100, Hydraulic Oil, Addinol,

IR Lubricants, 2081 FTIR spectra
<http://www.ir-spectra.com/>

IR Lubricants

ADDINOL HLP 22, Hydraulic Oil, Addinol,
ADNOC HYDRAULIC OIL, Hydraulic Oil, Adnoc-Fod,
BP ENERGOL HLP 46, Hydraulic Oil, BP,
CAT BIO HYDO SAE 10W (HEES), Hydraulic Oil, Catterpillar,
CAT HYDRAULIC SAE 10W, Hydraulic Oil, Catterpillar,
ESSO UNIVIS J13, Hydraulic Oil, Esso,
HYDRAX HLP 46, Hydraulic Oil, EXFOX,
EXXON HYJET IV-A PLUS, Hydraulic Oil, Exxon,
PLANTOHYD 40N, Hydraulic Oil, Fuchs,
RENOLIN HM 46, Hydraulic Oil, Fuchs,
RENOLIN HV 32, Hydraulic Oil, Fuchs,
RENOLIN OTHP 32, Hydraulic Oil, Fuchs,
VG 46, Hydraulic Oil, Fuchs,
MOBIL DTE OIL 21, Hydraulic Oil, Mobil,
MOBIL DTE OIL 25, Hydraulic Oil, Mobil,
MOBIL DTE OIL 26, Hydraulic Oil, Mobil,
MOBIL HLPD 32, Hydraulic Oil, Mobil,
MOBIL SHC 524, Hydraulic Oil, Mobil,
HIDROKOMOL HVU46, Hydraulic Oil, MOL,
HYDRAUNYCOIL FH 14, Hydraulic Oil, NYCO,
HYDRAUNYCOIL FH 15, Hydraulic Oil, NYCO,
STEOL M, Hydraulic Oil, NYCO,
OEMETA HF/AS , Hydraulic Oil, OEMETA,
HLP 46, Hydraulic Oil,
HM32, Hydraulic Oil,
HM46, Hydraulic Oil,
HM68, Hydraulic Oil,
HTG40, Hydraulic Oil,
MIL-L-24459, Hydraulic Oil,
MOGUL GREEN SAVE, Hydraulic Oil,
POWER ROOF HYDRAULIC OIL, Hydraulic Oil,
HM-32, Hydraulic Oil, Paramo,
HV32, Hydraulic Oil, Paramo,
OH-HM 46, Hydraulic Oil, Paramo,
TELEX E, Hydraulic Oil, REPSOL,
AIRCRAFT SHOCKERS OIL, Hydraulic Oil, Russia,
SHELL DONAX TA CHECK, Hydraulic Oil, Shell,
SHELL HYDROL DO 32, Hydraulic Oil, Shell,
SHELL TELLUS 29, Hydraulic Oil, Shell,
MADIT OH-HM 46/1998, Hydraulic Oil, Slovnaft,
MADIT OH-HM 46/1999, Hydraulic Oil, Slovnaft,
SUNOCO SUNVIS 800 WR-HV, Hydraulic Oil, Sunoco,
EQUIVIS ZS15, Hydraulic Oil, Total,
TOTAL AEROHYDRAULIC 520, Hydraulic Oil, Total,
TOTAL AZZOLA ZS 46, Hydraulic Oil, Total,
VAT LHM, Hydraulic Oil, VAT OIL,
HLP46, Hydraulic Oil, WINTERSHALL,
BP MACCURAT 68, Machine Oil, BP,
ESSO FEBIS K 220, Machine Oil, Esso,
ESSO MACHINERY OIL 22, Machine Oil, Esso,
ESSO SPINASSO 22, Machine Oil, Esso,
FESTO MACHINE OIL, Machine Oil, FESTO,
PAPER MACHINE OIL, Machine Oil, Fina,

IR Lubricants, 2081 FTIR spectra

<http://www.ir-spectra.com/>

OIL FOR SEWING MACHINES AND BIKES, Machine Oil, GELNICA,
 ETNA SS LIGHT OIL, Machine Oil, Mobil,
 MOBIL SHC PM 220, Machine Oil, Mobil,
 MOBIL VACTRA 2, Machine Oil, Mobil,
 MOBIL VACTRA 4, Machine Oil, Mobil,
 MOBIL VACTRA OIL 2, Machine Oil, Mobil,
 MOBIL VELOCITE OIL 6/1998, Machine Oil, Mobil,
 MOBIL VELOCITE OIL 6/2000, Machine Oil, Mobil,
 MACHINE OIL LOW TEMPERATURE, Machine Oil,
 MADOL 920, Machine Oil,
 OL-N22, Machine Oil, Paramo,
 OL-P4A, Machine Oil, Paramo,
 PARAMO OD-11/1995, Machine Oil, Paramo,
 PARAMO OD-4, Machine Oil, Paramo,
 PARAMO OD-8, Machine Oil, Paramo,
 PARAMO OL-B1/1995, Machine Oil, Paramo,
 PARAMO OL-B1/1999, Machine Oil, Paramo,
 PARAMO OL-B2, Machine Oil, Paramo,
 PARAMO OL-B4/1995, Machine Oil, Paramo,
 OD-11 DARK OIL Y.1997, Machine Oil, Petrochema,
 OD-16 DARK OIL, Machine Oil, Petrochema,
 ON 15 REFRIGERATOR OIL, Machine Oil, Petrochema,
 ON 22 REFRIGERATOR OIL, Machine Oil, Petrochema,
 ON 22V REFRIGERATOR OIL, Machine Oil, Petrochema,
 ON 32 REFRIGERATOR OIL, Machine Oil, Petrochema,
 ON 46 REFRIGERATOR OIL, Machine Oil, Petrochema,
 ON 68 REFRIGERATOR OIL, Machine Oil, Petrochema,
 TOTAL AZZOLA 22, Machine Oil, Total,
 TOTAL AZZOLA ZS 100, Machine Oil, Total,
 LIMPIOL 354, Machine Oil, WOERMANN,
 CASTROL SYNTILO R, Metalworking Fluid, Castrol,
 LUBRICUT, Metalworking Fluid, CHEMTOOL,
 CIMCLEAN 50, Metalworking Fluid, Cimcool,
 CIMSTAR 450, Metalworking Fluid, Cimcool,
 CIMTECH 10, Metalworking Fluid, Cimcool,
 AQUASOL 4100, Metalworking Fluid, CMT,
 CMT SUPERCUT TOOL, Metalworking Fluid, CMT,
 SUPERCUT 4-95.1, Metalworking Fluid, CMT,
 UNIPLUS 1-77, Metalworking Fluid, CMT,
 DICROL, Metalworking Fluid, CRODA,
 DASCOOL S14, Metalworking Fluid, DA STUART,
 STUART EXCELENE 416, Metalworking Fluid, DA STUART,
 STUART EXCELENE 425, Metalworking Fluid, DA STUART,
 ECOCUT 915, Metalworking Fluid, Fuchs,
 MOBIL GENREX 22, Metalworking Fluid, Mobil,
 MOBIL GENREX 26, Metalworking Fluid, Mobil,
 MOBIL GENREX 26B, Metalworking Fluid, Mobil,
 MOBILMET 25, Metalworking Fluid, Mobil,
 METAL-CUTTING OIL, Metalworking Fluid,
 PARAMO PARACUT 15, Metalworking Fluid, Paramo,
 PARAMO PARACUT 25, Metalworking Fluid, Paramo,
 PARAMO PARACUT 32 CF, Metalworking Fluid, Paramo,
 PARAMO PARACUT 46 CF, Metalworking Fluid, Paramo,

IR Lubricants, 2081 FTIR spectra
<http://www.ir-spectra.com/>

IR Lubricants

PARAMO PARACUT P3CF, Metalworking Fluid, Paramo,
CUTTING OIL EMULZIN D, Metalworking Fluid, Petrochema,
CUTTING OIL KATOL PP, Metalworking Fluid, Petrochema,
ECOCUT, Metalworking Fluid, Petrochema,
ECOCUT W 921, Metalworking Fluid, Petrochema,
SHELL GARIA OIL, Metalworking Fluid, Shell,
ADDINOL SUPER STAR MX SAE 15W40, Motor Oil, Addinol,
ADNOC SUPER 20W50, Motor Oil, Adnoc-Fod,
AGIP SIGMA TURBO SAE 15W40, Motor Oil, Agip,
AGIP SINT 2000 10W40, Motor Oil, Agip,
AGIP SINT TURBODIESEL SAE 10W40, Motor Oil, Agip,
AGIP SUPER DIESEL SAE 15W40, Motor Oil, Agip,
AGIP UNIVERSAL SAE 15W40 Y.1999, Motor Oil, Agip,
AGIP UNIVERSAL SAE 15W40 Y.2002, Motor Oil, Agip,
ARAL EXTRA TURBORAL 10W40, Motor Oil, Aral,
2T SUPER SELFMIX, Motor Oil, AVANTI,
AVANTEL DIESEL 15W40, Motor Oil, AVANTI,
CARRERA BASIC 10W40, Motor Oil, AVANTI,
CARRERA SYNTHEC 5W30, Motor Oil, AVANTI,
BP ENERGOL 40, Motor Oil, BP,
BP TERRAC EXTRA SAE 15W40, Motor Oil, BP,
BP VANELLUS C3 EXTRA SAE 5W40, Motor Oil, BP,
BP VANELLUS C3 MULTIGRADE SAE 5W40 Y.1999, Motor Oil, BP,
BP VANELLUS C3 MULTIGRADE SAE 5W40 Y.2000, Motor Oil, BP,
BP VISCO 5000 SAE 5W40, Motor Oil, BP,
API SG/CD, Motor Oil, BRIGGS&STRATTON,
CALTEX MOTOR OIL SAE 5W40, Motor Oil, CALTEX,
CALTEX RPM DELO 400 OIL MULTIGRADE SAE 5W40, Motor Oil, CALTEX,
RED RAM BLUE VELVET 10W40, Motor Oil, CANADIAN TURBO,
TURBO MULTI-SERVICE 15W40, Motor Oil, CANADIAN TURBO,
CHEVRON CUSTOM MULTIGRADE SAE 5W40, Motor Oil, Chevron,
CHEVRON MOTOR DELO 400, Motor Oil, Chevron,
CHEVRON SUPREME MULTIGRADE SAE 5W40, Motor Oil, Chevron,
MULTIGRADE SUPREME 10W30, Motor Oil, Chevron,
DBV LEICHTLAUF MOTOR OIL SAE 10W40, Motor Oil, DBV,
SYNTEX GRADE 10W40, Motor Oil, DENICOL,
TURBO SELECT 15W40, Motor Oil, DENICOL,
ESSO SUPER OIL X SAE 5W40, Motor Oil, Esso,
ESSO ULTRA X SAE 10W40, Motor Oil, Esso,
ESSO UNIFARM SAE 10W30, Motor Oil, Esso,
ESSO UNIFARM SAE 5W30, Motor Oil, Esso,
ESSOLUBE MHC UNIVERSAL, Motor Oil, Esso,
ESSOLUBE MHX, Motor Oil, Esso,
ESSOLUBE TDS SAE 5W40, Motor Oil, Esso,
ESSOLUBE TDX SAE 10W40, Motor Oil, Esso,
MULTI USE 15W40, Motor Oil, EUROL,
FIAT MOTOR OIL 20W50, Motor Oil, FIAT,
PANOLIN UNIVERSAL SAE 5W40, Motor Oil, GARANTOL PANOLIN,
POWER STEERING FLUID, Motor Oil, GOLD EAGLE,
SUPER DUTY SAE 15W40, Motor Oil, GULF,
TURBO DIESEL SAE 15W40, Motor Oil, GULF,
MOGUL DIESEL DTT SAE 10W40, Motor Oil, KORAMO,
MOGUL DIESEL DTT SAE 5W40, Motor Oil, KORAMO,

IR Lubricants, 2081 FTIR spectra

<http://www.ir-spectra.com/>

IR Lubricants

MOGUL SUPER SAE 15W50, Motor Oil, KORAMO,
FORMULA SUPER 15W40 HD, Motor Oil, Liqui Moly,
LUKOIL 5W40, Motor Oil, LUKOIL,
MOBIL DELVAC 1220 20W20, Motor Oil, Mobil,
MOBIL DELVAC 300 SUPER SAE 15W40, Motor Oil, Mobil,
MOBIL DELVAC XHP 15W40, Motor Oil, Mobil,
MOBIL PEGASUS 489, Motor Oil, Mobil,
MOBIL RALLY FORMULA 5W40, Motor Oil, Mobil,
MAXIMOL 99 SAE 10W40, Motor Oil, MOL,
MAXIMOL SINT SAE 5W40, Motor Oil, MOL,
SUPER DIESEL SAE 15W40, Motor Oil, MOL,
MOTUL 2000 MULTIGRADE SAE 15W40, Motor Oil, MOTUL,
OMV AUSTROTRAC 15W30, Motor Oil, OMV,
OMV AUSTROTRAC SAE 10W30, Motor Oil, OMV,
OMV ECOSYN DIESEL SAE 10W-40, Motor Oil, OMV,
OMV FULL SYN, Motor Oil, OMV,
OMV GHD 40, Motor Oil, OMV,
OMV KOREX MB, Motor Oil, OMV,
OMV START EL, Motor Oil, OMV,
OMV START F SAE 20W-50, Motor Oil, OMV,
OMV SYN COM, Motor Oil, OMV,
OMV SYN STAR 4T 5W50, Motor Oil, OMV,
OMV TRUCK LD SAE 5W40, Motor Oil, OMV,
OMV TRUCK LIGHT 5W40, Motor Oil, OMV,
OMV TRUCK M PLUS SAE 5W40, Motor Oil, OMV,
COLD ROOM MOTOR OIL 5W30, Motor Oil,
COLD ROOM MOTOR OIL SAE 10W30, Motor Oil,
ENGINE SAE 15W40, Motor Oil,
MOTOR OIL 5W30, Motor Oil,
SERVICE MOTOR OIL SAE 10W30, Motor Oil,
SERVICE MOTOR OIL SAE 10W40, Motor Oil,
SERVICE MOTOR OIL SAE 15W40, Motor Oil,
SYNTHETIC 5W50, Motor Oil,
TWO-STROKE MOTOR OIL, Motor Oil,
OA-M6ADS II 20W10, Motor Oil, Paramo,
PARAMO M2T, Motor Oil, Paramo,
PARAMO TRYSK, Motor Oil, Paramo,
PARAMO TRYSK GLOBAL SAE 10W40, Motor Oil, Paramo,
PARAMO TRYSK M6A, Motor Oil, Paramo,
PARAMO TRYSK M6AD, Motor Oil, Paramo,
PARAMO TRYSK M6ADS II PLUS, Motor Oil, Paramo,
PARAMO TRYSK M9ADV, Motor Oil, Paramo,
PARAMO TRYSK RACING, Motor Oil, Paramo,
PARAMO TRYSK SUPER 20W40, Motor Oil, Paramo,
PARAMO TRYSK SUPER SAE 15W40, Motor Oil, Paramo,
PARAMO TRYSK SUPER TURBO, Motor Oil, Paramo,
PARAMO TRYSK SUPER TURBO (M7ADSI) SAE 5W40, Motor Oil, Paramo,
TRYSK GLOBAL, Motor Oil, Paramo,
PETROMIN 2-S ENGINE OIL (PETROBOARD), Motor Oil, PETROMIN OILS,
PETROMIN DIESEL HD 40, Motor Oil, PETROMIN OILS,
PETROMIN EXTRA SUPER SAE 5W40, Motor Oil, PETROMIN OILS,
PETROMIN PREMIUM 40, Motor Oil, PETROMIN OILS,
SUPER HD II SAE 5W40 ENERGY CONSERVING, Motor Oil, PHILLIPS 66,

IR Lubricants, 2081 FTIR spectra

<http://www.ir-spectra.com/>

IR Lubricants

FORMULA RALLYE 15W40, Motor Oil, Q8 Oil,
FORMULA TOP 10W40, Motor Oil, Q8 Oil,
ELITE SAE 15W40, Motor Oil, REPSOL,
PERFORMANCE PLUS TC-A, Motor Oil, SAFETY-KLEEN,
PERFORMANCE PLUS ZR SAE 10W40, Motor Oil, SAFETY-KLEEN,
SELENIA MOTOR OIL 10W60, Motor Oil, SELENIA,
VS MAX SAE 15W40, Motor Oil, SELENIA,
SHELL HARVELLA T SAE 10W30, Motor Oil, Shell,
SHELL HELIX PLUS SAE 10W40, Motor Oil, Shell,
SHELL MYRINA SAE 15W40, Motor Oil, Shell,
SHELL RIMULA TX 20W50, Motor Oil, Shell,
SHELL RIMULA TX SAE 15W40, Motor Oil, Shell,
MADIT DIESEL 2T M9ADS, Motor Oil, Slovnaft,
MADIT DIESEL M9ADS, Motor Oil, Slovnaft,
MADIT EMOL M3 ADX CHECK, Motor Oil, Slovnaft,
MADIT EMOL M3 ADX CHECK, Motor Oil, Slovnaft,
MADIT FARMAN, Motor Oil, Slovnaft,
MADIT LTA, Motor Oil, Slovnaft,
MADIT M2T, Motor Oil, Slovnaft,
MADIT M3AD 10W30, Motor Oil, Slovnaft,
MADIT M6ADS II PLUS, Motor Oil, Slovnaft,
MADIT M7ADSIII 20W40, Motor Oil, Slovnaft,
MADIT SUPER M7AD SAE SAE 0W40, Motor Oil, Slovnaft,
MADIT TIRMAN SAE 5W40/1997, Motor Oil, Slovnaft,
MADIT TIRMAN SAE 5W40/1999, Motor Oil, Slovnaft,
MADIT TURBO PLUS M7ADSIII SAE 5W40, Motor Oil, Slovnaft,
MADIT TURBO PLUS SAE 5W40, Motor Oil, Slovnaft,
MADIT UNIOL M7ADX SAE 5W40, Motor Oil, Slovnaft,
MARINE TS, Motor Oil, Sun Oil,
SUNOCO FORZA SAE 5W40, Motor Oil, Sunoco,
SUNOCO FUTURA SAE 10W40, Motor Oil, Sunoco,
SUNOCO SUPER C MULTIGRADE SAE 15W40, Motor Oil, Sunoco,
SUNOCO SUPER DYNALUBE SAE 15W40, Motor Oil, Sunoco,
SUNOCO SUPER HPD SAE 15W40, Motor Oil, Sunoco,
SUNOCO SUPER MARINE, Motor Oil, Sunoco,
TEXACO HAVOLINE SUPERIOR GRADE SAE 10W30, Motor Oil, Texaco,
TEXACO URSA SUPER LA SAE 15W40, Motor Oil, Texaco,
TEXACO URSA SUPER SAE 20W50, Motor Oil, Texaco,
TEXACO URSA SUPER TD SAE 15W40, Motor Oil, Texaco,
TOTAL MULTAGRI TM SAE 15W30, Motor Oil, Total,
TOTAL QUARTZ 3000 SAE 5W40, Motor Oil, Total,
TOTAL QUARTZ 5000 SAE 5W40, Motor Oil, Total,
TOTAL RUBIA 4400 15W40, Motor Oil, Total,
TOTAL RUBIA TIR 6400 15W40, Motor Oil, Total,
TOTAL RUBIA XT 15W40/1996, Motor Oil, Total,
TOTAL RUBIA XT A 15W40, Motor Oil, Total,
TOTAL RUBIA XT SAE 15W40/2000, Motor Oil, Total,
ALL-CLIMATE 15W40, Motor Oil, Valvoline,
ALL-FLEET PLUS 15W40, Motor Oil, Valvoline,
DURABLEND 10W40, Motor Oil, Valvoline,
SYNPOWER 5W40, Motor Oil, Valvoline,
SYNPOWER 5W50, Motor Oil, Valvoline,
VAT SYNTECH MOTOR OIL 10W40, Motor Oil, VAT OIL,

IR Lubricants, 2081 FTIR spectra

<http://www.ir-spectra.com/>

IR Lubricants

SUPER MOTOR OIL HDC-IL TURBO 15W40, Motor Oil, VECO,
VECO TURBO DIESEL 15W40, Motor Oil, VECO,
FACTORY FILL MOTOR OIL SAE 10W30 - A, Motor Oil, VW,
FACTORY FILL MOTOR OIL SAE 10W30 - B, Motor Oil, VW,
FACTORY FILL MOTOR OIL SAE 10W30 - C, Motor Oil, VW,
FACTORY FILL MOTOR OIL SAE 10W30 - E, Motor Oil, VW,
FACTORY FILL MOTOR OIL SAE 10W30 - F, Motor Oil, VW,
LUBRIGUARD SAE 10W30, Motor Oil, WESTLAND,
RUNNING-IN OIL SAE 15W40, Motor Oil, WINTERSHALL,
WINTERSHALL REKORD SAE 15W40, Motor Oil, WINTERSHALL,
WINTERSHALL TFE SAE 10W40, Motor Oil, WINTERSHALL,
CLOVE OIL, Natural, H. REYNAUD, 8000-34-8
MARINE FISH OIL, Natural, 8016-13-5
RAPESEED OIL, Natural,
ADIBIS ADX 1664 B, Oil Additive, Adibis,
ADIBIS ADX 1667, Oil Additive, Adibis,
ADIBIS ADX 3831 H, Oil Additive, Adibis,
ADIBIS ADX 3831 R, Oil Additive, Adibis,
ADIBIS ADX 3866, Oil Additive, Adibis,
ADIBIS ADX 3891, Oil Additive, Adibis,
ADIBIS GBX 3804, Oil Additive, Adibis,
EMPICRYL 6059, Oil Additive, ALBRIGHT AND WILSON ,
EMPICRYL 6113, Oil Additive, ALBRIGHT AND WILSON ,
EMPICRYL HV 210, Oil Additive, ALBRIGHT AND WILSON ,
EMPICRYL HV 220, Oil Additive, ALBRIGHT AND WILSON ,
EMPICRYL HV 226, Oil Additive, ALBRIGHT AND WILSON ,
EMPICRYL PPT 144, Oil Additive, ALBRIGHT AND WILSON ,
EMPICRYL PPT 148, Oil Additive, ALBRIGHT AND WILSON ,
EMPICRYL PT 1544, Oil Additive, ALBRIGHT AND WILSON ,
VIPRAL - MF/150, Oil Additive, AUSCHEM,
BACTRACINE, Oil Additive, BACTRIA LABORATORIES ,
TRIAZINE TRZ, Oil Additive, BACTRIA LABORATORIES ,
BASF PIBSA, Oil Additive, BASF,
BASF PIBSA COPOLYMER, Oil Additive, BASF,
KEROPON ES 3237, Oil Additive, BASF,
SEPACORR HT, Oil Additive, BASF,
ACEPOL 512 M, Oil Additive, BECKER CHEMIE,
BECROSAN 2050 T CHECK, Oil Additive, BECKER CHEMIE,
BECROSAN 2050 T CHECK, Oil Additive, BECKER CHEMIE,
BECROSAN 2079, Oil Additive, BECKER CHEMIE,
BECROSAN 21019, Oil Additive, BECKER CHEMIE,
BECROSAN 540, Oil Additive, BECKER CHEMIE,
BIOCIDE BF-1 CONCENTRATE, Oil Additive, BECKER CHEMIE,
BIOCIDE FKK - A, Oil Additive, BECKER CHEMIE,
BIOCIDE NFK, Oil Additive, BECKER CHEMIE,
BIOCIDE - BAKZID 2/1998, Oil Additive, BODE CHEMIE,
BIOCIDE - BAKZID 2/1999, Oil Additive, BODE CHEMIE,
BIOCIDE - BAKZID 80, Oil Additive, BODE CHEMIE,
BIOCIDE - BODROB, Oil Additive, BODE CHEMIE,
BIOCIDE - CINON, Oil Additive, BODE CHEMIE,
POWER BACKFOR, Oil Additive, CD-2,
FX1 MONOMOLECULAR SURFACE MODIFIER, Oil Additive, CFA,
ORONITE C 40601 - OLOA 8805, Oil Additive, Chevron,

IR Lubricants, 2081 FTIR spectra

<http://www.ir-spectra.com/>

IR Lubricants

ORONITE C 40602 - OLOA 9725-X, Oil Additive, Chevron,
ORONITE C 40603 - OLOA 278, Oil Additive, Chevron,
ORONITE C 41805 - OLOA 4976 FP, Oil Additive, Chevron,
ORONITE C 46414 - OLOA 4900C, Oil Additive, Chevron,
ORONITE C 46415 - OLOA 4917C, Oil Additive, Chevron,
DASCO ARC, Oil Additive, DA STUART,
DASKONZ, Oil Additive, DA STUART,
DEFOAMER FULLY SYNTHETIC, Oil Additive, DA STUART,
DEFOAMER HALF SYNTHETIC, Oil Additive, DA STUART,
EXKONZ 340, Oil Additive, DA STUART,
EXKONZ 3838, Oil Additive, DA STUART,
ECA 547, Oil Additive, Esso,
ECA 6482, Oil Additive, Esso,
ESSO WYROL 4, Oil Additive, Esso,
ESSO WYROL 6, Oil Additive, Esso,
HITEC 151, Oil Additive, ETHYL PETROLEUM ADDITIVES ,
HITEC 315, Oil Additive, ETHYL PETROLEUM ADDITIVES ,
HITEC 320, Oil Additive, ETHYL PETROLEUM ADDITIVES ,
HITEC 526, Oil Additive, ETHYL PETROLEUM ADDITIVES ,
DWG OIL TREATMENT KITKANPOISTAJA , Oil Additive, Finland,
SOCCO CN 10, Oil Additive, Fuchs,
SOCCO R1, Oil Additive, Fuchs,
SOCCO R2, Oil Additive, Fuchs,
TITAN UNIMOLY, Oil Additive, Fuchs,
FUNCTIONAL HF-500, Oil Additive, FUNCTIONAL PRODUCTS,
FUNCTIONAL V 508, Oil Additive, FUNCTIONAL PRODUCTS,
FUNCTIONAL WA-24, Oil Additive, FUNCTIONAL PRODUCTS,
GLEITMO 200, Oil Additive, GLEITMO,
CERECLOR 63L, Oil Additive, ICI,
CERECLOR 65L, Oil Additive, ICI,
INNOTECH INTERNATIONAL OIL TREATMENT, Oil Additive, INNOTECH,
KOMPADIT KA, Oil Additive, ISTROCHEM,
KOMPADIT UN, Oil Additive, ISTROCHEM,
MARKER FOR FUEL OIL LIGHT, Oil Additive, KAUCUK,
NA - SUL BSN-NF, Oil Additive, KING INDUSTRIES,
K-MART OIL TREATMENT REDUCES WEAR, Oil Additive, K-MART,
LEUNASOL 1000D, Oil Additive, LEUNA WERKE ,
LEUNASOL 2000D, Oil Additive, LEUNA WERKE ,
LUBRIZOL 8098, Oil Additive, LUBRIZOL,
LUBRIZOL LCA 200, Oil Additive, LUBRIZOL,
LUBRIZOL LZ 8295 B, Oil Additive, LUBRIZOL,
MILITEC-1, Oil Additive, MILITEC,
ALCONA 10 148, Oil Additive, MINERAL CHEMIE,
ALCONA A 949, Oil Additive, MINERAL CHEMIE,
ALCONA B 0031, Oil Additive, MINERAL CHEMIE,
ALCONA B-625, Oil Additive, MINERAL CHEMIE,
ALCONA D 209, Oil Additive, MINERAL CHEMIE,
BIOCIDE FKK - B, Oil Additive, MINERAL CHEMIE,
BIOCIDE NBF 21, Oil Additive, MINERAL CHEMIE,
BIOCIDE ST-1, Oil Additive, MINERAL CHEMIE,
BIOCIDE ST-9010, Oil Additive, MINERAL CHEMIE,
SYSTEMCLEANER CB - 1 CONCENTRATE, Oil Additive, MINERAL CHEMIE,
MOBILAD G-351, Oil Additive, Mobil,

IR Lubricants, 2081 FTIR spectra

<http://www.ir-spectra.com/>

GEAR UP, Oil Additive, MOTOR UP,
7-RK-228, Oil Additive, NALCO VISCO CHEMICALS,
8-RK-344, Oil Additive, NALCO VISCO CHEMICALS,
9-RK-156, Oil Additive, NALCO VISCO CHEMICALS,
NALCO 5375, Oil Additive, NALCO VISCO CHEMICALS,
DYEGUARD GREEN M, Oil Additive,
DYEGUARD RED D, Oil Additive,
DYEGUARD RED MC 25, Oil Additive,
DYEGUARD YELLOW Q 12, Oil Additive,
GEAR OIL ADDITIVE, Oil Additive,
HYDROCOR CV 50, Oil Additive,
LOW TEMPERATURE OIL ADDITIVE, Oil Additive,
METAL CONDITIONER, Oil Additive,
POWER SUPREME, Oil Additive,
PARADYN 60 (PD 60), Oil Additive, PARAMINS ESSO,
PARANOX 390, Oil Additive, PARAMINS ESSO,
PARATONE (PTN) 8000, Oil Additive, PARAMINS ESSO,
SYNACTO 245, Oil Additive, PARAMINS ESSO,
DUBACID, Oil Additive, Petrochema,
PETROSULPHONATE CALCIUM, Oil Additive, Petrochema,
PETROSULPHONATE SODIUM, Oil Additive, Petrochema,
PRO-TEC OIL BOOSTER, Oil Additive, PRO-TEC,
VANLUBE 704B, Oil Additive, R.T. VANDEBILT,
VANLUBE 819B, Oil Additive, R.T. VANDEBILT,
VANLUBE 849, Oil Additive, R.T. VANDEBILT,
VANLUBE 848, Oil Additive, R.T. VANDEBILT ,
REWOCOROS AC 28/86, Oil Additive, REWO,
REWOCOROS AC 28/88, Oil Additive, REWO,
REWOPON B 2045, Oil Additive, REWO,
REWOPON IM OA, Oil Additive, REWO,
REWOPON RA 280, Oil Additive, REWO,
ADDITIN RC 2540, Oil Additive, RHEIN CHEMIE RHEINON ,
ADDITIN RC 3741, Oil Additive, RHEIN CHEMIE RHEINON ,
ADDITIN RC 8103, Oil Additive, RHEIN CHEMIE RHEINON ,
ADDITIN RC 9303, Oil Additive, RHEIN CHEMIE RHEINON ,
ADDITIN RC 9308, Oil Additive, RHEIN CHEMIE RHEINON ,
ADDITIN RC 9802, Oil Additive, RHEIN CHEMIE RHEINON ,
ADDITIN RC FI 1814, Oil Additive, RHEIN CHEMIE RHEINON ,
PLEXOL 1019M, Oil Additive, ROHM&HAAS,
PLEXOL 104C, Oil Additive, ROHM&HAAS,
PLEXOL 105C, Oil Additive, ROHM&HAAS,
PLEXOL 1460, Oil Additive, ROHM&HAAS,
PLEXOL 156, Oil Additive, ROHM&HAAS,
PLEXOL 618E, Oil Additive, ROHM&HAAS,
PLEXOL 702, Oil Additive, ROHM&HAAS,
PLEXOL 812, Oil Additive, ROHM&HAAS,
PLEXOL 812P, Oil Additive, ROHM&HAAS,
PLEXOL 812W, Oil Additive, ROHM&HAAS,
PRIMENE 81R, Oil Additive, ROHM&HAAS,
NEOVAL D.F.R. 625, Oil Additive, SAVEX,
DIESEL FUEL - DETERGENT ADDITIVE SAP 7500, Oil Additive, Shell,
SHELL ADDITIVES SAP-005, Oil Additive, Shell,
SHELL ADDITIVES SAP-2045B, Oil Additive, Shell,

IR Lubricants, 2081 FTIR spectra

<http://www.ir-spectra.com/>

SWEPCO 502, Oil Additive, SWEPCO,
 TEXACO TAD-921, Oil Additive, Texaco,
 TEXACO TLA-6205, Oil Additive, Texaco,
 HERMOSAL, Oil Additive, Velvana,
 PROMOL 101, Oil Additive, Velvana,
 PROMOL 110, Oil Additive, Velvana,
 PROMOL 111, Oil Additive, Velvana,
 PROMOL 131, Oil Additive, Velvana,
 PROMOL 140, Oil Additive, Velvana,
 PROMOL 141, Oil Additive, Velvana,
 PROMOL 142, Oil Additive, Velvana,
 PROMOL 143, Oil Additive, Velvana,
 PROMOL 150, Oil Additive, Velvana,
 PROMOL 160, Oil Additive, Velvana,
 PROMOL 170A, Oil Additive, Velvana,
 PROMOL 211, Oil Additive, Velvana,
 PROMOL 213, Oil Additive, Velvana,
 PROMOL 224, Oil Additive, Velvana,
 PROMOL 424, Oil Additive, Velvana,
 PROMOL 431, Oil Additive, Velvana,
 PROMOL 442, Oil Additive, Velvana,
 PROMOL 481, Oil Additive, Velvana,
 PROMOL 483, Oil Additive, Velvana,
 PROMOL 485, Oil Additive, Velvana,
 PROMOL 530, Oil Additive, Velvana,
 PROMOL 81, Oil Additive, Velvana,
 PROMOL 82, Oil Additive, Velvana,
 ADDITIVE FOR MOTOR OIL, Oil Additive, WURTH,
 COOL LUBRICANT, Oil Additive, WURTH,
 WYNN'S X-TEND SUPREME PLUS, Oil Additive, WYNN OIL,
 LMC - 70RR, Cleaner, CHEMTOOL,
 SILICONE OIL B, Lubricant, Germany,
 GLEITMO 871, Anticorrosion Oil, GLEITMO,
 P3 NEUTRAPON 5088, Cleaner, HENKEL,
 COCOAMIDOPROPYL HYDROXYSULTAINE , Surfactant, 70851-08-0
 MINERAL OIL, Base Oil, 8020-83-5
 MINERAL OIL, LIGHT AND HEAVY, Base Oil, 8042-47-5
 SHELLSOL 100/140, Military Use, Shell,
 POLY(ETHYLENE GLYCOL) 8000, Plasticizer, 25322-68-3
 POLYETHYLENE OXIDE 4000000, Plasticizer, 25322-68-3
 ESSO XD100, Quenching Oil, Esso,
 ISORAPID 277, Quenching Oil, Houghton Lubricor,
 MARQUENECH 600E, Quenching Oil, PETROFER,
 MARQUENECH 72254, Quenching Oil, PETROFER,
 ESSO SO 43, Rolling Oil, Esso,
 MOBIL GLYGOYLE 30, Rolling Oil, Mobil,
 MOBIL NP PLUS/1997, Rolling Oil, Mobil,
 MOBIL NP PLUS/1999, Rolling Oil, Mobil,
 MOBIL SHC 629, Rolling Oil, Mobil,
 MOBIL W261, Rolling Oil, Mobil,
 MOBISOL PM, Rolling Oil, Mobil,
 PARAMO B28, Rolling Oil, Paramo,
 PARAMO B31, Rolling Oil, Paramo,

IR Lubricants, 2081 FTIR spectra
<http://www.ir-spectra.com/>

IR Lubricants

PARAMO P31, Rolling Oil, Paramo,
DAMINOL, Rolling Oil, Petrochema,
ROLLING OIL K, Rolling Oil, Petrochema,
EXXON ISOPAR G, SBP spirit, Exxon,
EXXSOL D40, SBP spirit, Exxon,
BENZINE, SBP spirit,
KEROSENE PS2, SBP spirit,
N-ALCANES HEAVY, SBP spirit, Slovnaft,
N-ALCANES LIGHT, SBP spirit, Slovnaft,
ARAL EURAL CHAIN 220, Special Oil, Aral,
BP DIELECTRIC 180, Special Oil, BP,
ISOSYNTH VX 110 BF, Special Oil, DA STUART,
GLEITMO 880 BIO, Special Oil, GLEITMO,
HALLIS 225, Special Oil, HALLIS,
LUBRA K AHT, Special Oil, Kernite,
KLUBER LAMORA SUPER POLADD, Special Oil, Kluber,
KLUBER QUIETSCH-EX, Special Oil, Kluber,
KLUBER STRUCTOVIS GHD, Special Oil, Kluber,
AIR FILTER SPRAY, Special Oil, MOTUL,
OIL FOR VACUUM PUMPS, Special Oil,
TRIUMF K CHECK, Special Oil, Paramo,
TRIUMF K CHECK, Special Oil, Paramo,
BIOROL, Special Oil, Petrochema,
OD-3, Special Oil, Petrochema,
OIL FOR MEDICAL USE (GREASE), Special Oil, Petrochema,
WHITE OIL 45, Process Oil, Petrochema,
ISOGRIN 130, Special Oil, PETROFER,
PETROFER ULTRA-SAFE 620, Special Oil, PETROFER,
INA POLAR 150, Special Oil, RAFINERIJA NAFTE,
ENGINE RESTORER AND LUBRICANT WITH CLS, Special Oil, RESTORE,
SHELL CRYSTAL OIL 30, Special Oil, Shell,
LUBRANOL, Special Oil, SOPURA,
TOTAL HOUSEHOLD LUBRICANT, Special Oil, Total,
BASE OIL SN 150, Special Oil, Velvana,
STEAM CYLINDER OIL 2, Steam Cylinder Oil, Schaeffer,
ALKANOL XC, Surfactant, 151-21-3
SODIUM ALKYL POLYGLYCOL ETHER SULFATE, Surfactant,
SYNTOPON A 100, Surfactant,
VANFRE IL-1, Surfactant, 151-21-3
TOMAH PA-14 AC, Surfactant, PARAMINS ESSO,
TRF-HX, Transformer Oil, Addinol,
ITE 320, Transformer Oil, Agip,
ISOLAN 4376, Transformer Oil, Aral,
ISOLAN T, Transformer Oil, Aral,
ISOLAN TT, Transformer Oil, Aral,
T2015, Transformer Oil, DUSSEK CAMPBELL,
ESSO UNIVOLT 84, Transformer Oil, Esso,
ESSO UNIVOLT N55, Transformer Oil, Esso,
RENOLIN E7, Transformer Oil, Fuchs,
RENOLIN HM68, Transformer Oil, Fuchs,
MOBIL ECT 37, Transformer Oil, Mobil,
MOBILMET 151, Transformer Oil, Mobil,
NYTRO 10X, Transformer Oil, NYNAS,

IR Lubricants, 2081 FTIR spectra

<http://www.ir-spectra.com/>

IR Lubricants

TRANSFORMER OIL 200, Transformer Oil, Petrochema,
TRANSFORMER OIL 300, Transformer Oil, Petrochema,
DIALA BX, Transformer Oil, Shell,
DIALA DX, Transformer Oil, Shell,
DIALA GX, Transformer Oil, Shell,
REGENOL, Transformer Oil, TECHNOL,
TECHNOL 2000, Transformer Oil, TECHNOL,
TECHNOL Y3000, Transformer Oil, TECHNOL,
ESSO TURBO OIL 10F, Turbine Oil, Esso,
TURBO OIL 2389, Turbine Oil, Esso,
EXXON TURBO OIL 2380, Turbine Oil, Exxon,
EXXON TURBO OIL 25, Turbine Oil, Exxon,
MOBIL DTE BB, Turbine Oil, Mobil,
MOBIL DTE EXTRA-HEAVY, Turbine Oil, Mobil,
MOBIL DTE HEAVY-MEDIUM, Turbine Oil, Mobil,
MOBIL DTE LIGHT, Turbine Oil, Mobil,
MOBIL DTE MEDIUM, Turbine Oil, Mobil,
MOBIL JET OIL 254, Turbine Oil, Mobil,
MOBIL JET OIL II, Turbine Oil, Mobil,
MOBIL MCP 319 A, Turbine Oil, Mobil,
TURBONYCOIL 210A, Turbine Oil, NYCO,
TURBONYCOIL 525-2A, Turbine Oil, NYCO,
T-68, Turbine Oil,
TB-32, Turbine Oil,
TB-46, Turbine Oil,
PARAMO LT 68, Turbine Oil, Paramo,
TURBINE OIL 46, Turbine Oil, Petrochema,
AVIATION OIL MS-8P, Turbine Oil, Russia,
AEROSHELL TURBINE OIL 500, Turbine Oil, Shell,
AEROSHELL TURBINE OIL 555, Turbine Oil, Shell,
SHELL TURBO OIL T32, Turbine Oil, Shell,
SHELL TURBO OIL T46, Turbine Oil, Shell,
MADIT TB 46, Turbine Oil, Sloznaft,
Olive Oil, Edible, 8001-25-0
Sesame Oil, Edible, 8008-74-0
Sunflower oil, Edible, 8001-21-6
Agar, Emulsifier, 9002-18-0
Brij 35, Polyoxyethylene (35) Lauryl Ether, Emulsifier, 9002-92-0
Cremophor, Emulsifier, 61791-12-6
Dimethicone, Emulsifier, 9016-00-6
Ethoxylated mono- and diglycerides, Emulsifier,
Gelatine, Emulsifier, 9000-70-8
Glycerol, Emulsifier, 56-81-5
Glycerol esters of wood rosins, Emulsifier,
Glycerol Stearate, Emulsifier, 123-94-4
Gum Tragacanth, Emulsifier, 9000-65-1
Magnesium Stearate, Emulsifier, 557-04-0
Manganese stearate, Emulsifier, 3353-05-7#
Methyl Linoleate, Emulsifier, 112-63-0
Methyl stearate, Emulsifier, 112-61-8
Mono and Diglycerides of fatty acids, Emulsifier,
Oleyl alcohol, Emulsifier, 143-28-2
Palmitic acid, Emulsifier, 57-10-3

IR Lubricants, 2081 FTIR spectra

<http://www.ir-spectra.com/>

IR Lubricants

PEG-200 Glyceryl Tallowate, Emulsifier, 68153-76-4
Polyglycerol esters of fatty acids, Emulsifier,
Polyglycerol Polyricinoleate, Emulsifier,
Polymethacrylate, Emulsifier, 25087-26-7
Polynol A, mono and diglycerides, Emulsifier,
Polyoxyethene (40) stearate, Emulsifier, 9004-99-3
Sodium Stearyl Fumarate, Emulsifier, 4070-80-8
Sorbitan monolaurate, Emulsifier, 1338-39-2
Sorbitan monopalmitate (Span 40), Emulsifier, 26266-57-9
Sorbitan monostearate, Span 60, Emulsifier, 1338-41-6
Sorbitan Oleate, Span 80, Emulsifier, 1338-43-8
Sorbitan tristearate, Emulsifier, 26658-19-5
Stearamide, Emulsifier, 124-26-5
Stearic Acid, Emulsifier, 57-11-4
Tween 20, POE (20) sorbitan monolaurate, Emulsifier, 9005-64-5
Tween 40, POE (20) sorbitan monopalmitate, Emulsifier, 9005-66-7
Tween 65, POE (20) sorbitan tristearate, Emulsifier, 9005-71-4
Tween 80, POE (20) sorbitan monooleate, Emulsifier, 9005-65-6
Tween 85, POE (20) sorbitan trioleate, Emulsifier, 9005-70-3
Zinc Stearate, Emulsifier, 557-05-1
Anise Oil, Essential Oil, H. REYNAUD,
Almond oil sweet, Essential Oil, 8007-69-0
Bergamot Oil, Essential Oil, 8007-75-8
Coconut oil, Essential Oil, 8001-31-8
Lemon Oil Type 2, Essential Oil, 8008-56-8
Myrrh Bil Type B, Essential Oil, 8016-37-3
Peanut Oil, Essential Oil, 8002-03-7#
Pine Oil, Essential Oil, 8002-09-3#
Sassafras Oil, Essential Oil, 8006-80-2
Thyme Oil, Essential Oil, 8007-46-3
Wheat Germ Oil, Essential Oil, 8006-95-9
Eucalyptus oil, Essential Oil, SHEIDON DRUG,
Kerosene, Fuel,
Beeswax white, Glazing Agent, 8012-89-3
Beeswax yellow, Glazing Agent, 8006-40-4
Candelilla wax, Glazing Agent, 8006-44-8
Ceresin, Glazing Agent, 8001-75-0
Cetaceum, Cetyl Palmitate 15, Glazing Agent, 540-10-3
Glycerol Tristearate, Glazing Agent, 555-43-1
Trilaurin, Glazing Agent, 538-24-9
Tripalmitin, Glazing Agent, 555-44-2
Lanoline, Grease, 8006-54-0
Caraway Oil, Natural, DRAGOCOL,
Apricot seed oil, Natural, 8002-78-6
Castor oil, Natural,
Castor oil hydrogenated, Natural, 8001-78-3
Castor oil, ethoxylated, Natural, 61791-12-6
Cedar Wood Oil, Natural, 8000-27-9
Cod-Liver Oil, Natural, 8001-69-2
Cottonseed Oil, Natural, 8001-29-4
Fennel oil, Natural, 8006-84-6
Linseed Oil, Natural, 8001-26-1
Maize Oil, Natural, 8001-30-7

IR Lubricants, 2081 FTIR spectra

<http://www.ir-spectra.com/>

Soybean Oil, Natural, 8001-22-7
 Ethoxylated Lanoline 4000, Emulsifier, 61790-81-6
 Ethyl palmitate, Emulsifier, 628-97-7
 Glyceryl monooleate, Emulsifier, 25496-72-4
 Paraffin Oil, Base Oil, 8012-95-1
 Silicone oil A, Lubricant, 63148-62-9
 Talc, Lubricant, 14807-96-6
 Benzyl benzoate, Plasticizer, 120-51-4
 Dibutyl phthalate, Plasticizer, 84-74-2
 Diethyl phthalate, Plasticizer, 84-66-2
 Dioctyl phthalate, Plasticizer, 117-81-7
 Poly(ethylene glycol) 1000, Plasticizer, 25322-68-3
 Poly(ethylene glycol) 2000, Plasticizer, 25322-68-3
 Poly(ethylene glycol) 300, Plasticizer, 25322-68-3
 Triethyl citrate, Plasticizer, 77-93-0
 Sodium Dodecyl Sulfate, Surfactant, 151-21-3
 AEON 3000, Compressor Oil, Gardner Denver
 AGIP Antifreeze Econom, Antifreeze
 AGIP Arnica 104/FR, Hydraulic Oil
 AGIP ATF IID, Gear Oil, Automatic Transmission Fluid
 AGIP Exidia HG 68, Slideway Oil
 AGIP Formula LL FO 5W-30, Motor Oil, Ford, Jaguar, Mazda
 AGIP OSO 46, Hydraulic Oil
 AGIP Rotra HY 80W90, Gear Oil
 AGIP Rotra HY DB, Gear Oil
 AGIP Rotra SX/S, Gear Oil
 Amber Cup Grease, Grease, Red Giant, year 1960
 Amber Pressure Gun Grease, Grease, Red Giant, year 1960
 Anderol 3486, Compressor Oil, synthetic
 Anderol 555, Compressor Oil, synthetic
 ARAL Dural SR 5, Industrial Oil, spindle
 ARAL Mega Turboral 10W-40, Motor Oil
 ARAL Multi Turboral 15W-40, Motor Oil, Mineral
 ARAL Multielastic 15W-40, Motor Oil
 ARAL Sarol CI 100, Cutting Fluid
 ARAL Sulnit K, Cutting Fluid
 ARAL Sumurool CM 10, Cutting Fluid
 ARAL Sumurool CM 22, Cutting Fluid
 ARAL Vitam DE 22, Hydraulic Oil
 ARAL Vitam GF 22, Hydraulic Oil
 ARAL Vitam GF 32, Hydraulic Oil
 ARAL Vitam GF 46 (HLP), Hydraulic Oil
 ARAL Vitam HF 46 (HVLP), Hydraulic Oil
 ATF, Gear Oil, VW, Automatic Transmission Fluid
 Axle Grease, Grease, Red Giant, year 1960
 Azalea Absolute Oil, Essential Oil
 Bay Laurel, Essential Oil
 Bay Oil, Essential Oil
 Bihol 32 NT, Hydraulic Oil, Vasoil
 Biocon PewaS Protec 40, Oil for sliding guides, Pewas
 Bitter Fennel Oil, Essential Oil
 Bitzer BSE 170, Compressor Oil, Bitzer, ester oil, refrigerators
 Black Pepper Oil, Essential Oil

IR Lubricants, 2081 FTIR spectra

<http://www.ir-spectra.com/>

IR Lubricants

Blaser Grindex 10 Conc., Cutting Fluid, Blaser Swissslube, Synthetic
BP Bartran 46, Hydraulic Oil
BP Bartran HV 46, Hydraulic Oil
BP Bartran HV 68, Hydraulic Oil, biodegradable
BP Bartran SHF-S 46, Hydraulic Oil, Mineral
BP Diesel 15W-40, Motor Oil
BP Energear EP 80W 90, Gear Oil
BP Energear SGX 75W-90, Gear Oil, Manual Transmission
BP Energol CS HB 460, Circulation Oil
BP Energol HLP-HM 32, Hydraulic Oil
BP Energol HLP-HM 46, Hydraulic Oil
BP Energol HLP-HM 68, Hydraulic Oil
BP Energol RC 150, Compressor Oil
BP Energol THB 32, Turbine Oil
BP Gear Oil XF 4040 M, Gear Oil
BP Transcal N, Heat Transfer Medium
BP Vanellus C3 extra, Motor Oil
BP Vanellus C5 global, Motor Oil, diesel engines
BP Vanellus C6 10W40 used 90000 km, Motor Oil
BP Vanellus C7 15W40 used 66000 km, Motor Oil
BP Vanellus E6 15W40, Motor Oil
BP Vanellus E6 15W40 used 86000 km, Motor Oil
BP Vanellus E7 plus, Motor Oil, diesel engines
BP Visco 2000 15W-40, Motor Oil
BP Visco 3000, Motor Oil
BP Visco 5000 5W-40, Motor Oil
Cainnamon Leaf Oil, Essential Oil
Cajeput Oil, Essential Oil
Capital Axle Grease, Grease, Red Giant, year 1960
Cardamom Oil, Essential Oil
CASE MS 1230, Tractor Oil, Synthetic
CASTROL 291/100, Cutting Fluid
CASTROL 50501, Motor Oil, TDI
CASTROL Careclean, Anticorrosion Oil
CASTROL EBD 10 (concentrate), Additive
CASTROL Edge 5W30, Motor Oil, For modern VW and Audi Engines
CASTROL Edge Sport 0W40, Motor Oil, High Performance
CASTROL Enduron plus 5W-40, Motor Oil
CASTROL GP 20W-50, Motor Oil, Motorbikes
CASTROL GPS 10W-40, Motor Oil, synthetic motorcycle oil
CASTROL GTX 3 15W-40, Motor Oil
CASTROL GTX Powertec 10W-40, Motor Oil
CASTROL GTX Protec, Motor Oil
CASTROL Honilo 919, Cutting Fluid, Honing
CASTROL Hysol T 15, Metalworking Fluid
CASTROL LX 15W-40, Motor Oil
CASTROL Magnaglide D 220, Slideway Oil
CASTROL Magnatec 10W40, Motor Oil
CASTROL Magnatec 5W40, Motor Oil
CASTROL Optigear BM 460, Gear Oil
CASTROL Optimol Optigear BM 150, Gear Oil
CASTROL Protection 15W-40, Motor Oil
CASTROL R 15W-40, Motor Oil

IR Lubricants, 2081 FTIR spectra

<http://www.ir-spectra.com/>

IR Lubricants

CASTROL R4 Superbike, Motor Oil, Fully Synthetic Motorcycle
CASTROL SLX 0W30, Motor Oil, Long Life II
CASTROL SLX Formula 0W-30, Motor Oil
CASTROL SLX Longlife II 0W-30, Motor Oil
CASTROL SLX Longlife III 5W-30, Motor Oil, VW and Audi
CASTROL SLX Longtec 0W-30, Motor Oil
CASTROL SLX Professional Power Flow 5W30, Motor Oil, Long Life III
CASTROL Syntilo RHS, Metalworking Fluid
CASTROL Tecton 15W40, Motor Oil
CASTROL Transmax Z, Gear Oil, Fully Synthetic
CASTROL Tribol 800/220, Gear Oil, Synthetic
CASTROL TTS, Motor Oil, Fully Synthetic
CASTROL TWS Motorsport 10W60, Motor Oil
CASTROL TXT Softec 5W40, Motor Oil, Advanced Engines
CASTROL TXT Softec Plus 5W-30, Motor Oil
CASTROL Viscogen G, Industrial Oil, hollow glass machines
CASTROL Viscogen KL 23, Industrial Oil, thermal stability
CAT Deo, Motor Oil, Catterpillar
CAT Hees, Hydraulic Oil, Catterpillar, biodegradable
Cedar Wood Oil, Essential Oil
Ceramic Paste Specialty, Valvoline, spray
Cimcool Milpro, Cutting Fluid, Cimcool
Cinnamon Oil, Essential Oil
Citronella Oil From Asia, Essential Oil
Citronella Oil from USA, Essential Oil
Clary Sage Oil, Essential Oil
Classic 15W40, Motor Oil, Starline, Hightech
Clove Oil, Essential Oil
Compressor Oil EP-VDL 46, Compressor Oil , Chevron, Marine
Condaform S 450, Industrial Oil, Condat, Neat Oil
Coriander Oil, Essential Oil
Cygnus 220, Gear Oil, Equilon
Cypress Oil, Essential Oil
DEA Astron HLP 46, Hydraulic Oil
DEA Astron HLP 68, Hydraulic Oil
Delphi 10W-40, Motor Oil
Delphi Prestige Plus 5W-40, Motor Oil
Dentax G 80W90, Gear Oil, Iveco
Dexron III Automatic Ultra, Gear Oil, Starline
Diamond 10W40, Motor Oil, Starline, Synthetic
Diamond 5W40, Motor Oil, Starline, Synthetic
Diesel System Purge, Cleaner, Wynns
Dill Weed Oil, Essential Oil
Ditridecyl Phthalate
Divinol Synthogear Extra LS 75W90, Gear Oil, Zeller-Gmelin
Drum Cleaner, Aresal
Durablend 10W40, Motor Oil, Valvoline
Durixol W 25, Hardening Oil, Burgdorf
Durixol W 72, Hardening Oil, Burgdorf
Ecocut 3022, Cutting Fluid, Petrochema
Ecogear 320S, Gear Oil, Addinol
Ecotec 10W-40, Motor Oil, Motor Gold, Synthetic
Elemi Oil, Essential Oil

IR Lubricants, 2081 FTIR spectra

<http://www.ir-spectra.com/>

IR Lubricants

ELF Aviation AD 100, Aviation Oil
ELF Evolution SXR 5W30, Motor Oil, Synthetic
ELF Planetelf PAG 244, Compressor Oil, air condition cars
ELF Planetelf PAG 488, Compressor Oil
ELF Solaris 5W30, Motor Oil, Synthetic
ELF Spirit MS 5000 (concentrate), Cutting Fluid
Emkarate RL 68H, Refrigeration Oil, ICI
Emulkat UNI 100, Cutting Fluid, Triga
Engine Treatment, Additive, Slick 50
ESSO 85W 140 GX, Gear Oil
ESSO ATFD Dexron II, Gear Oil
ESSO Dortan N 33, Cutting Fluid
ESSO Estor PC 40, Motor Oil, Gas Engines
ESSO Gear Oil 85W-90, Gear Oil
ESSO Gearoil GP-D 85W-90, Gear Oil
ESSO Haerteoel XD 100, Hardening Oil
ESSO HLP 46, Hydraulic Oil
ESSO HLPD 46, Hydraulic Oil
ESSO Hydraulic Oil J 32, Hydraulic Oil
ESSO LD 0W-30, Motor Oil, diesel engines
ESSO Pruefoel 4113, Special Oil, ISO 4113, diesel pumps testing
ESSO Schneidoel DN 16, Cutting Fluid
ESSO Schneidoel DN 37, Cutting Fluid
ESSO Schneidoel HP 22, Cutting Fluid
ESSO Schneidoel LP 32, Cutting Fluid
ESSO Somentor EL 45, Rolling Oil
ESSO Spartan EP 150, Gear Oil, Industrial, extreme pressure
ESSO Spartan EP 320, Gear Oil, Industrial, extreme pressure
ESSO Teresso 32, Turbine Oil
ESSO Teresso 46, Turbine Oil
ESSO Teresso EP 46, Turbine Oil
ESSO Thermaloeel T, Heat Transfer Medium
ESSO Ultra 10W-40, Motor Oil
ESSO Ultra Turbo 10W-40, Motor Oil, diesel engines
ESSO Ultron 5W-40, Motor Oil, Fully Synthetic
ESSO Unifarm 10W-30, Tractor Oil
ESSO Unirex N3, Grease, High Temperature Bearing
ESSO Univis HVI 26, Hydraulic Oil
ESSO Univis J 13, Hydraulic Oil
ESSO Univis J 26, Hydraulic Oil
ESSO Univis N 46, Hydraulic Oil
ESSO Vacmul XLM 13, Cutting Fluid, Grinding
Essolube P-40, Motor Oil
Essolube PX 40, Motor Oil, Waste gas
Essolube XT 301, Motor Oil
Essolube XT 4, Motor Oil
Essolube XT 401, Motor Oil
Exel HLPV 46, Hydraulic Oil, Sun Oil
Extensia 10W-40, Motor Oil, Renault, Renault Trucks
EXXON Exxsol D 220/230, SBP spirit
EXXON Febis K 46, Slideway Oil
EXXON Febis K 68, Slideway Oil
EXXON NL GI-2, Grease

IR Lubricants, 2081 FTIR spectra

<http://www.ir-spectra.com/>

IR Lubricants

EXXON Nuto H 46, Hydraulic Oil
EXXON Nuto H 68, Hydraulic Oil
EXXON Spartan EP 220, Gear Oil
EXXON Torgue Fluid 56, Gear Oil, for tractors
EXXONMOBIL Primol 352, Medicinal Oil
EXXONMOBIL Teresstic T 46, Turbine Oil
Fennel Sweet Oil, Essential Oil
Fina 400 5W40, Motor Oil
Fina Caprano, Motor Oil, Diesel engines of boats
Fina Excel 10W40, Motor Oil, Semisynthetic
Fina Kappa Supra, Motor Oil
Finamatic HP, Gear Oil, Automatic Transmission Fluid
Fir Needle Oil, Essential Oil
Fir Oil, Essential Oil
First 300 10W40, Motor Oil, Fina
First 300 5W40, Motor Oil, Fina
Flagroot Oil, Essential Oil
Fomblin Y LOX, Solvay
Frankincense Oil, Essential Oil
FUCHS Afrostol KA 68, For cooling machines
FUCHS Biorol, Cutting Fluid
FUCHS Ecocut, Cutting Fluid
FUCHS Plantohyd 46 S, Hydraulic Oil
FUCHS Renep 68, Gear Oil
FUCHS Renep CGLP 68, Special Oil, industrial machines
FUCHS Reniso PAG 46, Compressor Oil, air condition cars
FUCHS Renolin B15 VG 46, Hydraulic Oil
FUCHS Renolin B46 HVI, Hydraulic Oil
FUCHS Renolin CLP 100, Gear Oil
FUCHS Renolin CLP 150, Gear Oil
FUCHS Renolin Eterna 46, Turbine Oil
FUCHS Renolin MA 32, Industrial Oil, Circulating
FUCHS Renolin MR 15 VG 46, Industrial Oil, for bearings
FUCHS Renolin MR 46 MC, Hydraulic Oil
FUCHS Renolin MR5, Hydraulic Oil
FUCHS Renolin PG 46, Hydraulic Oil
FUCHS Renolin Unisyn 150, Gear Oil, PAO based, high temperature
FUCHS Renolin Unisyn 68, Gear Oil, PAO based, high temperature
FUCHS Renolin VG 32, Hydraulic Oil
FUCHS Renolin VG 46, Hydraulic Oil
FUCHS Renolit IPR 2, Grease, Li, for valves
FUCHS Specinet, Bio Fluid
FUCHS Thermisol 251, Quenching Oil
FUCHS Thermisol QH 40, Quenching Oil
FUCHS Titan GM 2040 super, Motor Oil
FUCHS Titan Hydramot 1040 MC, Motor Oil
FUCHS Titan Supergear GL-4, Gear Oil
FUCHS Titan Superior HD 15W-40, Motor Oil
FUCHS Titan Supersyn F 5W30, Motor Oil, Synthetic
FUCHS Titan Supersyn Plus 5W40, Motor Oil, Synthetic
FUCHS Titan Supersyn SL LL plus 0W-30, Motor Oil
FUCHS Titan Trans HD 1550, Motor Oil
FUCHS Titan Truck 15W-40, Motor Oil, Diesel, Trucks

IR Lubricants, 2081 FTIR spectra

<http://www.ir-spectra.com/>

IR Lubricants

FUCHS Titan Universal HD 1540, Motor Oil, diesel engines
FUCHS Titan Universal HD 30, Motor Oil, diesel engines
FUCHS Wisura AK 3009, Industrial Oil
Gear Metro, Gear Oil, Starline, Mineral Based Multigrade
Geranium Oil, Essential Oil
Ginger Oil, Essential Oil
Glue and Sealant Remover, Loctite
Green Cup Grease, Grease, Red Giant, year 1960
Green Pressure Gun Grease, Grease, Red Giant, year 1960
GX-FE 10W40, Motor Oil, Valar, Semisynthetic
Hand Separator, Red Giant, year 1960
Heavy Castor, Red Giant, year 1960
Heavy Harvester, Red Giant, year 1960
High Performace Copper Spray, Spray, Valvoline
Houghto Safe 620, Metalworking Fluid, Houghton
Houghton Lubricor HF 122, Hydraulic Oil
Hydraulic G0040000M2, Hydraulic and Steering Oil, VW
Hydrosafe 620, Hydraulic Oil, G-C Lubricants, Fire Resistant
Chamomile Oil, Essential Oil
Jasmine Absolute Oil, Essential Oil
Juniper Berry Oil, Essential Oil
KLUBER Syntheso D 220 EP, Industrial Oil, Machinery Textile
Lemon Oil, Essential Oil
Lemon Verbena Oil, Essential Oil
Lemon Wormwood Oil, Essential Oil
Lemongrass Oil, Essential Oil
Lime Oil, Essential Oil
Limetta Oil, Essential Oil
Liqui Moly 5W-40, Motor Oil
Litsea Oil, Essential Oil
Mace Oil, Essential Oil
Madit M6AD, Motor Oil, Slovnaft
Madit SW Multi 68, Oil for sliding guides, Slovnaft
Mandarin Oil, Essential Oil
Marjoram Oil, Essential Oil
Martemp 722, Quenching Oil, Houghton Lubricor
Maxlife 15W40, Motor Oil, Valvoline, premium blend
Maxlife 5W40, Motor Oil, Valvoline, Synthetic
MOBIL 10W-40, Motor Oil
MOBIL Arctic Gargoyle 300, Compressor Oil, Marine
MOBIL Delvac Hydraulic 10W VG 46, Hydraulic Oil
MOBIL Delvac XHP 10W-40, Motor Oil
MOBIL Delvac XHP extra, Motor Oil, diesel engines
MOBIL DTE 15 M, Hydraulic Oil
MOBIL DTE 26 new, Hydraulic Oil
MOBIL DTE 26 thermal degradation, Hydraulic Oil
MOBIL DTE 832, Turbine Oil
MOBIL DTE 846, Turbine Oil
MOBIL DTE Excel 46, Hydraulic Oil
MOBIL DTE Oil 26, Hydraulic Oil, ISO-VG 68
MOBIL EAL Arctic 68, Compressor Oil, Fully Synthetic, POEs based
MOBIL Formula LD 0W-30, Motor Oil, Synthetic oil for VW
MOBIL Hydraulic Premium 46, Hydraulic Oil

IR Lubricants, 2081 FTIR spectra

<http://www.ir-spectra.com/>

IR Lubricants

MOBIL Hydraulic Premium 68, Hydraulic Oil
MOBIL Hydraulic Universal 46, Hydraulic Oil
MOBIL Chassis Grease, Grease
MOBIL Mobilfluid 316 M, Hydraulic Oil
MOBIL Mobilgear 636, Gear Oil
MOBIL Mobilmet , Cutting Fluid
MOBIL Mobilmet 426, Cutting Fluid
MOBIL Mobilmet 766, Cutting Fluid
MOBIL Mobilmet S 122, Cutting Fluid
MOBIL Mobilmet XLM, Cutting Fluid
MOBIL Mobiltherm 605, Heat Transfer Medium
MOBIL Mobiltrans HD 10W, Gear Oil
MOBIL Mornop 55, Cutting Fluid
MOBIL Pegasus 15W-40 - new, Motor Oil
MOBIL Pegasus 15W-40 after 2535 hours, oxidized, Motor Oil
MOBIL Pegasus 480, Motor Oil, Gas Engines
MOBIL Pegasus 610 new, Motor Oil, Gas Engines
MOBIL Pegasus 610, after 317 hours, bio fuel, oxidized, Motor Oil
MOBIL Pegasus 610, after 867 hours, bio fuel, oxidized, Motor Oil
MOBIL Pegasus 705, Motor Oil, Gas Engines
MOBIL Pegasus 905, Motor Oil
MOBIL Pegasus 905, after 1770 hours, oxidized, Motor Oil
MOBIL Racing 2T, Motor Oil, 2 Cycle Motorcycle Oil
MOBIL Rarus 829, Compressor Oil, Synthetic
MOBIL Rarus SHC 1025, Compressor Oil
MOBIL SHC 624, Gear Oil
MOBIL SHC Formula LD 0W30 , Motor Oil, Fully Synthetic
MOBIL Super M 10W-40, Motor Oil
MOBIL Supersyn 0W40 Protection Formula, Motor Oil, Fully Synthetic
MOBIL Supersyn 15W50, Motor Oil, Motorsport Fully Synthetic
MOBIL Vacmul A6, Cutting Fluid
MOBIL Vacmul EDM 3, Cutting Fluid
MOBIL Vacmul O3C, Cutting Fluid
MOBIL Vacmul S 5, Cutting Fluid
Mogul HLPD 46, Hydraulic Oil, Paramo
Mogul HM 46, Hydraulic Oil, Paramo
Mogul TB 32 S, Turbine Oil, Paramo
Mogul TB 46 S, Turbine Oil, Paramo
Mogul Traktol UTTO, Gear Oil, Paramo, for tractors
Mogul Traktol UTTO-EKO, Gear Oil, Paramo, for tractors
Myrrh Oil from Asia, Essential Oil
Myrrh Oil from USA, Essential Oil
Myrtle Oil, Essential Oil
Neatsfood Comp., Red Giant, year 1960
Neroli Oil, Essential Oil
No.600 SR Transmission, Gear Oil, Red Giant, year 1960
Nujol, Thermo
Nutmeg Oil from Asia, Essential Oil
Nutmeg Oil from USA, Essential Oil
Off Shore, Spray, Normfest, Fully Synthetic Ceramic Spray
OH HM 46, Hydraulic Oil, Paramo
OH-HM 32, Hydraulic Oil, Paramo
OH-HM 68, Hydraulic Oil, Paramo

IR Lubricants, 2081 FTIR spectra

<http://www.ir-spectra.com/>

OH-HV 32, Hydraulic Oil, Paramo
 OMV LD SAE 15W-40, Motor Oil
 OMV Superdraw, Motor Oil
 OMV Truck FE plus 10W-40, Motor Oil
 Orange Bitter Oil, Essential Oil
 Oregano Oil, Essential Oil
 OT-HP3 ISO VG32, Hydraulic Oil, Paramo
 Palmarosa Oil, Essential Oil
 Panolin Target 15W-40, Motor Oil
 Paramo Cut 22, Cutting Fluid, Paramo
 Paramo HM 46, Hydraulic Oil, Paramo
 Paramo OT-HP 3, Gear Oil, Paramo, Hydrodynamic Transmission
 PARAMO PP 44, Gear Oil, Paramo
 PARAMO PP 90H, Gear Oil, Paramo
 Paramol TB 46, Turbine Oil, Paramo
 Patchouli Oil, Essential Oil
 Penetrating oil + MoS₂ Lubricating, Valvoline
 PENNZOIL ATF+3, Gear Oil, Automatic Transmission Fluid
 PENNZOIL AW 46, Hydraulic Oil
 PENNZOIL GL-4 80W90, Gear Oil
 Pensoline A SAE Body 50, Pensoline, year 1960
 Pensoline B SAE Body 60, Pensoline, year 1960
 Pensoline C SAE Body 40, Pensoline, year 1960
 Pensoline E SAE Body 20, Pensoline, year 1960
 Pensoline K SAE Body 30, Pensoline, year 1960
 Pensoline Zero SAE Body 10, Pensoline, year 1960
 Peppermint Oil, Essential Oil
 Petitgrain Oil, Essential Oil
 Pfeiffer Vacuum Pump Oil, Compressor Oil
 Pfinder AP 167/22, Metalworking Fluid
 Pine Oil from USA, Essential Oil
 Pine Scotch Oil, Essential Oil
 Pneumofore PL 493, Heat Transfer Medium
 Power Steering Fluid, Additive, Slick 50
 Premium High Vacuum Pump oil, Compressor Oil, Robinair
 Primol Eko HTG 40, Hydraulic Oil, Aniveg Nat
 Q8 Mahler MA, Motor Oil, Mahler
 Ravensara Oil, Essential Oil
 Red Giant A SAE Body 50, year 1960
 Red Giant B SAE Body 60, year 1960
 Red Giant C SAE Body 40, year 1960
 Red Giant E SAE Body 20, year 1960
 Red Giant K SAE Body 30, year 1960
 Red Giant Zero SAE Body 10, year 1960
 Reflo 68A, Refrigeration Oil, Petro Canada, Food Grade Ammonia
 Reifenschaum, Berner, for glossy tyres
 Rose Geranium Oil, Essential Oil
 Rosemary Oil, Essential Oil
 Rosewood Oil, Essential Oil
 Sabesto Rost off Plus, Würth
 Sandalwood Oil, Essential Oil
 Servo E-AT/PS, Servo Oil, Ford
 SHELL Albida EP2, Grease

IR Lubricants, 2081 FTIR spectra
<http://www.ir-spectra.com/>

IR Lubricants

SHELL Clavus G 68, Compressor Oil, for refrigerator compressors
SHELL Comptela S 46, Compressor Oil
SHELL Corena D46, Compressor Oil
SHELL Diala D, Electroinsulating Oil
SHELL Diala DX, Electroinsulating Oil
SHELL Fenella, Rolling Oil
SHELL Garia 405 OHA 16CNC4, Metalworking Fluid
SHELL Gear Oil EP 80W, Gear Oil
SHELL Helix Plus 10W40, Motor Oil, for Diesel, Synthetic
SHELL Helix Plus VB 10W-40, Motor Oil, Synthetic
SHELL Helix Ultra 0W30, Motor Oil, Long Life
SHELL Helix Ultra 5W-40, Motor Oil, Synthetic
SHELL Helix Ultra AC 5W-30, Motor Oil, Synthetic
SHELL Hydrol DO 46, Hydraulic Oil
SHELL Omala 100, Gear Oil
SHELL Rimula 15 W-40 used, Motor Oil
SHELL Rimula Super 15W-40, Motor Oil
SHELL Rimula Ultra 10W40 used 52000km, Motor Oil, 52000/9000km
SHELL Rimula Ultra 10W-40, Motor Oil
SHELL Spirax AX 80W90, Axle Oil, Iveco Truck
SHELL Spirax GX 80W90, Gear Oil, Iveco Truck
SHELL Tellus 46, Hydraulic Oil
SHELL Tellus S 32, Hydraulic Oil
SHELL Tellus S-46 New, Hydraulic Oil
SHELL Tellus S-46 Used in Engine, Hydraulic Oil
SHELL Tellus S-46 Used in Press, Hydraulic Oil
SHELL Tellus T 32, Hydraulic Oil, HVLP
SHELL Thermia Oil B, Heat Transfer Medium
SHELL Turbo T 32, Turbine Oil
SHELL Turbo T 46, Turbine Oil
Siempezka PL, Fuel Oil, Slovtherov
Silicone Spray Lubricating, Spray, Valvoline
Silkal 93, Coyote
Silver Synt 10W40 SL/CF, Motor Oil, Global Oil, Belgium
Sintogrind TTK, Grinding Oil, Oelheld
Spearmint Oil, Essential Oil
Special Cement BL, Rema Tip Top
Spruce Oil, Essential Oil
Star Anise Oil, Essential Oil
Steam Cylinder, Steam Cylinder Oil, Red Giant, year 1960
Stihl Two Stroke for Saw, Motor Oil, Stihl, Two Stroke for Saw
Sultex D 220, Cutting Fluid, Equilon
SUNOCO Oil Treatment, Add to Oil, Additive
SUNOCO Sungrind Range, Cutting Fluid
SUNOCO Suniso 4GS, Compressor Oil, Refrigeration
Sweet Basil Oil, Essential Oil
Sweet Orange Oil, Essential Oil
Synpower 5W40 sample A, Motor Oil, Valvoline, Fully Synthetic
Synpower 5W40 sample B, Motor Oil, Valvoline, Fully Synthetic
Syntholube G 46, Gear Oil, Molyduval, Synthetic
Synto 75W90, Gear Oil, Starline
Tagetes Oil, Essential Oil
Tangerine Oil, Essential Oil

IR Lubricants, 2081 FTIR spectra

<http://www.ir-spectra.com/>

IR Lubricants

Tarbolene, Red Giant, year 1960
Tea Tree Oil, Essential Oil
TEXACO Alcor DD 22, Hydraulic Oil
TEXACO Capella HFC, Compressor Oil, Refrigeration POE
TEXACO Cetus S 46, Compressor Oil
TEXACO Geartex EP A 85W90, Gear Oil
TEXACO Havoline 10W40, Motor Oil
TEXACO Havoline 505.01, Motor Oil
TEXACO Havoline 5W30, Motor Oil
TEXACO Havoline 5W40, Motor Oil, Synthetic
TEXACO Havoline Extra 10W40, Motor Oil
TEXACO Meropa 100, Gear Oil
TEXACO Meropa 150, Industrial Oil
TEXACO Meropa 220, Gear Oil
TEXACO Meropa 320, Gear Oil
TEXACO Meropa 320, Used 14 days, water contents, Gear Oil
TEXACO Motex 2T SX, Motor Oil, two stoke engines
TEXACO Rando HD 32, Hydraulic Oil
TEXACO Rando HD 46, Hydraulic Oil
TEXACO Rando HD 68, Industrial Oil, pumps, compressors, bearings
TEXACO Rando HD-Z 22, Hydraulic Oil
TEXACO Rando HD-Z 46, Hydraulic Oil
TEXACO Rando HD-Z 68, Hydraulic Oil
TEXACO Regal Premium 46, Turbine Oil, Marine
TEXACO Regal Premium EP 32, Turbine Oil
TEXACO Regal R&O 32, Turbine Oil, industrial, multipurpose
TEXACO Regal R&O 46, Turbine Oil, industrial, multipurpose
TEXACO Regal R&O 68, Turbine Oil, industrial, multipurpose
TEXACO Taxaflow Super Concentrate, Additive
TEXACO Texatherm HT 22, Heat Transfer Medium
TEXACO Ursa Super LA 15W40, Motor Oil
TEXACO Ursa Super TD, Motor Oil
Thuja Oil, Essential Oil
Top Blend 80W-90, Gear Oil, Kernite, Roller Tables Gearboxes
TOTAL Azolla AF 46, Hydraulic Oil
TOTAL Azolla ZS 46, Hydraulic Oil
TOTAL Azolla ZS 68, Hydraulic Oil
TOTAL Biohydran TMP 46, Hydraulic Oil
TOTAL Carter EP 220, Gear Oil
TOTAL Cortis, Industrial Oil
TOTAL Cortis MS 460, Industrial Oil, rollings
TOTAL Equivis 46 ZS, Hydraulic Oil
TOTAL Fluide ATX, Gear Oil, automated transmission
TOTAL Hydrobio 46, Hydraulic Oil, Biodegradable Synthetic
TOTAL Nateria MJ 40, Motor Oil, Gas Engines
TOTAL Preslia 32, Turbine Oil
TOTAL Preslia 46, Turbine Oil
TOTAL Quartz 10W-40, Motor Oil
TOTAL Quartz 5000 15W40, Motor Oil
TOTAL Quartz 5W40, Motor Oil, Synthetic
TOTAL Rubia S 40, Motor Oil, diesel engines
TOTAL Seriola 1510, Heat Transfer Medium
TOTAL TM 80W-90, Gear Oil

IR Lubricants, 2081 FTIR spectra

<http://www.ir-spectra.com/>

TOTAL Torilis 2500, Process Oil
TOTAL Torilis 6200, Process Oil
TOTAL Transmission BV 75W-80W Advanced, Gear Oil
TOTAL Valona MS 7023, Cutting Fluid
Trans Premium 150, Gear Oil, Statoil
Trans Premium 220, Gear Oil, Statoil
Trans Premium 220, Gear Oil, Statoil
Transmission 1#9120541, Gear Oil, GM
Traxon 80W-90, Gear Oil, Petro Canada
Tutela Truck FE AXLE 75W90, Motor Oil, Urania, Iveco Truck
Ultrasafe, PG based
Urania FE 5W30, Motor Oil, Iveco Truck
Urania LD 15W40, Motor Oil, Iveco Truck
Valar Atea, Hydraulic Oil
Valar egida XQ 15W-40, Motor Oil
Valar HLP 46, Hydraulic Oil
Valerian Root Oil, Essential Oil
Vetiver Oil, Essential Oil
Virginian Cedar Wood Oil, Essential Oil
VW Hydraulic, Hydraulic Oil, For Central Hydraulic and Servo
Way Lubricant X 220, Metalworking Fluid, Chevron
Way Lubricant X 68, Metalworking Fluid, Chevron
WD-40, Anticorrosion Oil, WD-40, spray
White Camphor Oil, Essential Oil
White Grease, Grease, Valvoline, spray
Wintergreen Oil, Essential Oil
Wormwood Oil, Essential Oil
Ylang Ylang Extra Oil, Essential Oil
Ylang-Ylang Oil, Essential Oil
York H, Compressor Oil, York
ZHS-Servo, Hydraulic and Steering Oil, Startol

List of Transmission Spectra

ADNOC GX 85W140, Gear Oil, Adnoc-Fod,
ADNOC HYDRAULIC OIL, Hydraulic Oil, Adnoc-Fod,
ADNOC SUPER 20W50, Motor Oil, Adnoc-Fod,
AEON 3000, Compressor Oil, Gardner Denver
AGIP ATF IID, Gear Oil, Automatic Transmission Fluid
AGIP Exidia HG 68, Slideway Oil
AGIP Formula LL FO 5W-30, Motor Oil, Ford, Jaguar, Mazda
AGIP OSO 46, Hydraulic Oil
AGIP Rotra HY 80W90, Gear Oil

IR Lubricants, 2081 FTIR spectra
<http://www.ir-spectra.com/>

AGIP Rotra HY DB, Gear Oil
 ALL-CLIMATE 15W40, Motor Oil, Valvoline,
 ALL-FLEET PLUS 15W40, Motor Oil, Valvoline,
 Anderol 3486, Compressor Oil, synthetic
 ANTICORROSION OIL /1998, Anticorrosion Oil, SPANJAARD,
 ANTICORROSION OIL /2002, Anticorrosion Oil, SPANJAARD,
 API SG/CD, Motor Oil, BRIGGS&STRATTON,
 Aral Kosmol TF 46; Turbine Oil; Aral
 ARAL Mega Turboral 10W-40, Motor Oil
 ARAL Multi Turboral 15W-40, Motor Oil, Mineral
 ARAL Multielastic 15W-40, Motor Oil
 ARAL Sulnit K, Cutting Fluid
 ARAL Sumurol CM 22, Cutting Fluid
 ARAL Vitam DE 22, Hydraulic Oil
 ARAL Vitam GF 22, Hydraulic Oil
 ARAL Vitam GF 32, Hydraulic Oil
 ARAL Vitam GF 46 (HLP), Hydraulic Oil
 ATF TYPE D/2002, Gear Oil, Valvoline,
 ATF, Gear Oil, VW, Automatic Transmission Fluid
 AVANTEL DIESEL 15W40, Motor Oil, AVANTI,
 AXLE OIL 80W90, Gear Oil, Q8 Oil,
 BP Bartran 46, Hydraulic Oil
 BP BARTRAN HV 46, Hydraulic Oil, BP,
 BP Bartran HV 46, Hydraulic Oil
 BP Bartran HV 68, Hydraulic Oil, biodegradable
 BP DIELECTRIC 200T, Transformer Oil, BP,
 BP Diesel 15W-40, Motor Oil
 BP Energear EP 80W 90, Gear Oil
 BP Energol CS HB 460, Circulation Oil
 BP Energol HLP-HM 32, Hydraulic Oil
 BP Energol HLP-HM 46, Hydraulic Oil
 BP Energol HLP-HM 68, Hydraulic Oil
 BP Energol LPT 68; Hydraulic Oil; BP
 BP Energol RC 150, Compressor Oil
 BP ENERGOL THB 32, Turbine Oil, BP,
 BP Energol THB 32, Turbine Oil
 BP ENERGOL THB 46, Turbine Oil, BP,
 BP ENERGOL WM 2, Special Oil, BP,
 BP ENERPAR M002, Process Oil, BP,
 BP Gear Oil XF 4040 M, Gear Oil
 BP SEVORA 532 L, Metalworking Fluid, BP,
 BP SEVORA 721M, Metalworking Fluid, BP,
 BP Transcal N, Heat Transfer Medium
 BP TURBINOL 32, Turbine Oil, BP,
 BP Vanellus C3 extra, Motor Oil
 BP Vanellus C5 global, Motor Oil, diesel engines
 BP Vanellus C6 10W40 used 90000 km, Motor Oil
 BP Vanellus C7 15W40 used 66000 km, Motor Oil
 BP Vanellus C8 Ultima; Motor Oil; BP
 BP Vanellus E6 15W40 used 86000 km, Motor Oil
 BP Vanellus E6 15W40, Motor Oil
 BP Vanellus E7 plus, Motor Oil, diesel engines
 BP Visco 2000 15W-40, Motor Oil

IR Lubricants, 2081 FTIR spectra
<http://www.ir-spectra.com/>

IR Lubricants

BP Visco 3000, Motor Oil
BP Visco 5000 5W-40, Motor Oil
CARRERA BASIC 10W40, Motor Oil, AVANTI,
Castrol 20W-30;Motor Oil used;Castrol
Castrol 20W-30;Motor Oil;Castrol
CASTROL 50501, Motor Oil, TDI
CASTROL Edge 5W30, Motor Oil, For modern VW and Audi Engines
CASTROL Edge Sport 0W40, Motor Oil, High Performance
CASTROL Enduron plus 5W-40, Motor Oil
CASTROL EPX 90, Gear Oil, Castrol,
CASTROL GP 20W-50, Motor Oil, Motorbikes
CASTROL GPS 10W-40, Motor Oil, synthetic motorcycle oil
CASTROL GTX 3 15W-40, Motor Oil
CASTROL GTX Powertec 10W-40, Motor Oil
CASTROL GTX Protec, Motor Oil
CASTROL LX 15W-40, Motor Oil
CASTROL Magnaglide D 220, Slideway Oil
CASTROL Magnatec 10W40, Motor Oil
CASTROL Magnatec 5W40, Motor Oil
CASTROL Optigear BM 460, Gear Oil
CASTROL Optimol Optigear BM 150, Gear Oil
CASTROL Protection 15W-40, Motor Oil
CASTROL R 15W-40, Motor Oil
CASTROL SLX 0W30, Motor Oil, Long Life II
CASTROL SLX Formula 0W-30, Motor Oil
CASTROL SLX Longlife II 0W-30, Motor Oil
CASTROL SLX Longlife III 5W-30, Motor Oil, VW and Audi
Castrol SLX Longtec 0W-20;Motor Oil;Castrol
CASTROL SLX Longtec 0W-30, Motor Oil
CASTROL SLX Professional Power Flow 5W30, Motor Oil, Long Life III
CASTROL SYNTEC ACTIVE SAE 10W-30, Motor Oil, Castrol,
CASTROL SYNTEC BLEND SAE 20W-50, Motor Oil, Castrol,
CASTROL Tecton 15W40, Motor Oil
CASTROL TXT Softec 5W40, Motor Oil, Advanced Engines
CASTROL TXT Softec Plus 5W-30, Motor Oil
CAT Deo, Motor Oil, Catterpillar
Classic 15W40, Motor Oil, Starline, Hightech
Compressor Oil EP-VDL 46, Compressor Oil , Chevron, Marine
COMPRESSOR OIL ISO 46, Compressor Oil, Kernite,
CONCRETE FLUID SUPER 120, Forming Oil, WOERMANN,
CONCRETE SEPARATOR 90, Forming Oil, WOERMANN,
CONCRETE SEPARATOR CONCENTRATE 96, Forming Oil, WOERMANN,
CONCRETE SEPARATOR LF 101, Forming Oil, WOERMANN,
CONTRAKOR FLUID H1, Anticorrosion Oil, Kluber,
DEA Astron HLP 46, Hydraulic Oil
DEA Astron HLP 68, Hydraulic Oil
Delphi 10W-40, Motor Oil
Delphi Prestige Plus 5W-40, Motor Oil
Dentax G 80W90, Gear Oil, Iveco
Dexron III Automatic Ultra, Gear Oil, Starline
Diamond 10W40, Motor Oil, Starline, Synthetic
Diamond 5W40, Motor Oil, Starline, Synthetic
DURABLEND 10W40, Motor Oil, Valvoline,

IR Lubricants, 2081 FTIR spectra

<http://www.ir-spectra.com/>

IR Lubricants

Durablend 10W40, Motor Oil, Valvoline
Durixol W 25, Hardening Oil, Burgdorf
Durixol W 72, Hardening Oil, Burgdorf
Ecocut 3022, Cutting Fluid, Petrochema
ECOCUT 915, Metalworking Fluid, Fuchs,
Ecogear 320S, Gear Oil, Addinol
Ecotec 10W-40, Motor Oil, Motor Gold, Synthetic
ELF Aviation AD 100, Aviation Oil
ELF Evolution SXR 5W30, Motor Oil, Synthetic
ELF Solaris 5W30, Motor Oil, Synthetic
ELITE SAE 15W40, Motor Oil, REPSOL,
Engine Treatment, Additive, Slick 50
ESSO 85W 140 GX, Gear Oil
ESSO ATF D 21 611, Gear Oil, Esso,
ESSO ATFD Dexron II, Gear Oil
ESSO Dortan N 33, Cutting Fluid
ESSO Estor PC 40, Motor Oil, Gas Engines
ESSO EXTRA MOTOR OIL SAE 15W-40, Motor Oil, Esso,
ESSO FLEXON 843, Machine Oil, Esso,
ESSO FLEXON 856, Special Oil, Esso,
ESSO Gear Oil 85W-90, Gear Oil
ESSO Gearoil GP-D 85W-90, Gear Oil
ESSO Haerteoel XD 100, Hardening Oil
ESSO HLP 46, Hydraulic Oil
ESSO HLPD 46, Hydraulic Oil
ESSO NUTO 460, Machine Oil, Esso,
ESSO NUTO H32, Hydraulic Oil, Esso,
ESSO PRIMOL 352, Special Oil, Esso,
ESSO RUST PROTECTION MZ 45, Anticorrosion Oil, Esso,
ESSO RUST-BAN 335, Anticorrosion Oil, Esso,
ESSO Schneidoel HP 22, Cutting Fluid
ESSO Schneidoel LP 32, Cutting Fluid
ESSO Somentor EL 45, Rolling Oil
ESSO Spartan EP 150, Gear Oil, Industrial, extreme pressure
ESSO Spartan EP 320, Gear Oil, Industrial, extreme pressure
ESSO Teresso 32, Turbine Oil
ESSO Teresso 46, Turbine Oil
ESSO Teresso EP 46, Turbine Oil
ESSO Thermaloel T, Heat Transfer Medium
ESSO Ultra 10W-40, Motor Oil
ESSO ULTRA SAE 10W 40, Motor Oil, Esso,
ESSO Ultra Turbo 10W-40, Motor Oil, diesel engines
ESSO Ultron 5W-40, Motor Oil, Fully Synthetic
ESSO Unifarm 10W-30, Tractor Oil
ESSO Univis N 46, Hydraulic Oil
ESSO UNIVOLT 524, Transformer Oil, Esso,
ESSO Vacmul XLM 13, Cutting Fluid, Grinding
ESSO XD 80, Quenching Oil, Esso,
Essolube P-40, Motor Oil
Essolube PX 40, Motor Oil, Waste gas
ESSOLUBE XT 301 SAE 15W-40, Motor Oil, Esso,
Essolube XT 301, Motor Oil
Essolube XT 4, Motor Oil

IR Lubricants, 2081 FTIR spectra

<http://www.ir-spectra.com/>

IR Lubricants

Essolube XT 401, Motor Oil
Exel HLPV 46, Hydraulic Oil, Sun Oil
Extensia 10W-40, Motor Oil, Renault, Renault Trucks
EXXON AVIATION OIL 100, Aviation Oil, Exxon,
EXXON AVIATION OIL ELITE 20W50, Aviation Oil, Exxon,
EXXON Febis K 46, Slideway Oil
EXXON Febis K 68, Slideway Oil
EXXON Nuto H 46, Hydraulic Oil
EXXON Nuto H 68, Hydraulic Oil
EXXON Spartan EP 220, Gear Oil
EXXON Torgue Fluid 56, Gear Oil, for tractors
EXXONMOBIL Primol 352, Medicinal Oil
EXXONMOBIL Teresstic T 46, Turbine Oil
FESTO MACHINE OIL, Machine Oil, FESTO,
Fina 400 5W40, Motor Oil
Fina Caprano, Motor Oil, Diesel engines of boats
Fina Excel 10W40, Motor Oil, Semisynthetic
Fina Kappa Supra, Motor Oil
Finamatic HP, Gear Oil, Automatic Transmission Fluid
First 300 10W40, Motor Oil, Fina
First 300 5W40, Motor Oil, Fina
FORMULA RALLYE 15W40, Motor Oil, Q8 Oil,
FORMULA SUPER 15W40 HD, Motor Oil, Liqui Moly,
FORMULA TOP 10W40, Motor Oil, Q8 Oil,
FUCHS Ecocut, Cutting Fluid
FUCHS Renep 68, Gear Oil
FUCHS Renep CGLP 68, Special Oil, industrial machines
FUCHS Renolin B15 VG 46, Hydraulic Oil
FUCHS Renolin B46 HVI, Hydraulic Oil
FUCHS Renolin CLP 100, Gear Oil
FUCHS Renolin CLP 150, Gear Oil
FUCHS Renolin Eterna 46, Turbine Oil
FUCHS Renolin MA 32, Industrial Oil, Circulating
FUCHS Renolin MR 15 VG 46, Industrial Oil, for bearings
FUCHS Renolin MR 46 MC, Hydraulic Oil
FUCHS Renolin MR5, Hydraulic Oil
FUCHS Renolin VG 32, Hydraulic Oil
FUCHS Renolin VG 46, Hydraulic Oil
Fuchs Renolin VG32;Hydraulic oil;Fuchs
FUCHS Thermisol 251, Quenching Oil
FUCHS Thermisol QH 40, Quenching Oil
FUCHS Titan GM 2040 super, Motor Oil
Fuchs Titan GTI W-30;Motor Oil used;Fuchs
Fuchs Titan GTI W-30;Motor Oil;Fuchs
FUCHS Titan Hydramot 1040 MC, Motor Oil
FUCHS Titan Supergear GL-4, Gear Oil
FUCHS Titan Superior HD 15W-40, Motor Oil
FUCHS Titan Supersyn F 5W30, Motor Oil, Synthetic
FUCHS Titan Supersyn Plus 5W40, Motor Oil, Synthetic
FUCHS Titan Trans HD 1550, Motor Oil
FUCHS Titan Truck 15W-40, Motor Oil, Diesel, Trucks
FUCHS Titan Universal HD 1540, Motor Oil, diesel engines
FUCHS Titan Universal HD 30, Motor Oil, diesel engines

IR Lubricants, 2081 FTIR spectra

<http://www.ir-spectra.com/>

FUCHS Wisura AK 3009, Industrial Oil
 Gear Metro, Gear Oil, Starline, Mineral Based Multigrade
 GX-FE 10W40, Motor Oil, Valar, Semisynthetic
 HIDROKOMOL HVU46, Hydraulic Oil, MOL,
 HP GL-4 80W90, Gear Oil, Valvoline,
 HP-GL-5 80W90, Gear Oil, Valvoline,
 HYDRAX HLP 46, Hydraulic Oil, EXFOX,
 HYKOMOL K SAE 80W90 - 1999, Gear Oil, MOL,
 CHEVRON MOTOR DELO 400, Motor Oil, Chevron,
 KOMAGMOL 150, Gear Oil, MOL,
 LIMPIOL 354, Machine Oil, WOERMANN,
 MADIT AUTOMATIC ATF/1999, Gear Oil, Slovnaft,
 MADIT DIESEL 2T M9ADS, Motor Oil, Slovnaft,
 MADIT DIESEL M9ADS, Motor Oil, Slovnaft,
 MADIT EMOL M3 ADX CHECK, Motor Oil, Slovnaft,
 MADIT EMOL M3 ADX CHECK, Motor Oil, Slovnaft,
 MADIT FARMAN, Motor Oil, Slovnaft,
 MADIT LTA, Motor Oil, Slovnaft,
 MADIT M3AD 10W30, Motor Oil, Slovnaft,
 Madit M6AD, Motor Oil, Slovnaft
 MADIT M6ADS II PLUS, Motor Oil, Slovnaft,
 MADIT OH-HM 46/1998, Hydraulic Oil, Slovnaft,
 MADIT OH-HM 46/1999, Hydraulic Oil, Slovnaft,
 MADIT OK-VC 150, Compressor Oil, Slovnaft,
 MADIT OL-B4/2000, Bearing Oil, Slovnaft,
 MADIT PP 80 SAE 80W, Gear Oil, Slovnaft,
 MADIT PP 80/2000, Gear Oil, Slovnaft,
 MADIT PP 90/2001, Gear Oil, Slovnaft,
 Madit SW Multi 68, Oil for sliding guides, Slovnaft
 MADIT TB 46, Turbine Oil, Slovnaft,
 MARQUENECH 600E, Quenching Oil, PETROFER,
 MARQUENECH 72254, Quenching Oil, PETROFER,
 Martemp 722, Quenching Oil, Houghton Lubricor
 MAXIMOL 99 SAE 10W40, Motor Oil, MOL,
 MAXIMOL SINT SAE 5W40, Motor Oil, MOL,
 Maxlife 15W40, Motor Oil, Valvoline, premium blend
 Maxlife 5W40, Motor Oil, Valvoline, Synthetic
 MD7ADS III Super Trysk; Motor Oil used; Paramo
 MD7ADS III Super Trysk; Motor Oil; Paramo
 MOBIL 10W-40, Motor Oil
 MOBIL Delvac Hydraulic 10W VG 46, Hydraulic Oil
 Mobil Delvac MX 15W-40; Motor Oil used; Mobil
 Mobil Delvac MX 15W-40; Motor Oil; Mobil
 MOBIL Delvac XHP 10W-40, Motor Oil
 MOBIL Delvac XHP extra, Motor Oil, diesel engines
 MOBIL DTE 15 M, Hydraulic Oil
 MOBIL DTE 26 new, Hydraulic Oil
 MOBIL DTE 26 thermal degradation, Hydraulic Oil
 MOBIL DTE 832, Turbine Oil
 MOBIL DTE Excel 46, Hydraulic Oil
 MOBIL DTE Oil 26, Hydraulic Oil, ISO-VG 68
 MOBIL Hydraulic Premium 46, Hydraulic Oil
 MOBIL Hydraulic Premium 68, Hydraulic Oil

IR Lubricants, 2081 FTIR spectra
<http://www.ir-spectra.com/>

IR Lubricants

MOBIL Hydraulic Universal 46, Hydraulic Oil
MOBIL Mobilmet , Cutting Fluid
MOBIL Mobilmet 426, Cutting Fluid
MOBIL Mobilmet 766, Cutting Fluid
MOBIL Mobilmet XLM, Cutting Fluid
MOBIL Mobiltherm 605, Heat Transfer Medium
MOBIL Mobiltrans HD 10W, Gear Oil
MOBIL Pegasus 15W-40 - new, Motor Oil
MOBIL Pegasus 15W-40 after 2535 hours, oxidized, Motor Oil
MOBIL Pegasus 480, Motor Oil, Gas Engines
MOBIL Pegasus 610 new, Motor Oil, Gas Engines
MOBIL Pegasus 610, after 867 hours, bio fuel, oxidized, Motor Oil
MOBIL Pegasus 705, Motor Oil, Gas Engines
MOBIL Pegasus 905, Motor Oil
Mobil Rarus 425; Compressor Oil; Mobil
MOBIL SHC 624, Gear Oil
MOBIL SHC Formula LD 0W30 , Motor Oil, Fully Synthetic
MOBIL Super M 10W-40, Motor Oil
MOBIL Supersyn 0W40 Protection Formula, Motor Oil, Fully Synthetic
MOBIL Supersyn 15W50, Motor Oil, Motorsport Fully Synthetic
MOBIL Vacmul A6, Cutting Fluid
MOBIL Vacmul EDM 3, Cutting Fluid
MOBIL Vacmul O3C, Cutting Fluid
MOBIL Vacmul S 5, Cutting Fluid
Mogul HLPD 46, Hydraulic Oil, Paramo
Mogul HM 46, Hydraulic Oil, Paramo
MOGUL SUPER SAE 15W50, Motor Oil, KORAMO,
Mogul TB 32 S, Turbine Oil, Paramo
Mogul TB 46 S, Turbine Oil, Paramo
Mogul Traktol UTTO, Gear Oil, Paramo, for tractors
MULTI USE 15W40, Motor Oil, EUROL,
OH HM 46, Hydraulic Oil, Paramo
OH-HM 68, Hydraulic Oil, Paramo
OH-HV 32, Hydraulic Oil, Paramo
OIL FOR SEWING MACHINES AND BIKES, Machine Oil, GELNICA,
OMV AUSTROTRAC 15W30, Motor Oil, OMV,
OMV ECOSYN DIESEL SAE 10W-40, Motor Oil, OMV,
OMV LD SAE 15W-40, Motor Oil
OMV START F SAE 20W-50, Motor Oil, OMV,
OMV Superdraw, Motor Oil
OMV Truck FE plus 10W-40, Motor Oil
OMV UNIGEAR, Gear Oil, OMV,
ON 22V REFRIGERATOR OIL, Machine Oil, Petrochema,
ON 32 REFRIGERATOR OIL, Machine Oil, Petrochema,
ON 46 REFRIGERATOR OIL, Machine Oil, Petrochema,
ON 68 REFRIGERATOR OIL, Machine Oil, Petrochema,
OT-HP3 ISO VG32, Hydraulic Oil, Paramo
Panolin Target 15W-40, Motor Oil
Paramo Cut 22, Cutting Fluid, Paramo
Paramo HM 46, Hydraulic Oil, Paramo
Paramo OT-HP 3, Gear Oil, Paramo, Hydrodynamic Transmission
PARAMO PP 90H, Gear Oil, Paramo
Paramol TB 46, Turbine Oil, Paramo

IR Lubricants, 2081 FTIR spectra

<http://www.ir-spectra.com/>

IR Lubricants

PENNZOIL AW 46, Hydraulic Oil
PENNZOIL GL-4 80W90, Gear Oil
Pfeffer Vacuum Pump Oil, Compressor Oil
Pneumofore PL 493, Heat Transfer Medium
POWER STEERING FLUID, Motor Oil, GOLD EAGLE,
Power Steering Fluid, Additive, Slick 50
Premium High Vacuum Pump oil, Compressor Oil, Robinair
Q8 Mahler MA, Motor Oil, Mahler
Red Giant A SAE Body 50, year 1960
Red Giant B SAE Body 60, year 1960
Red Giant C SAE Body 40, year 1960
Red Giant K SAE Body 30, year 1960
RED RAM BLUE VELVET 10W40, Motor Oil, CANADIAN TURBO,
Reflo 68A, Refrigeration Oil, Petro Canada, Food Grade Ammonia
RENOLIN HM68, Transformer Oil, Fuchs,
ROTO-INJECTFLUID, Compressor Oil, ATLAS COPCO,
Sentron GG 40; Motor Oil used; Petro Canada
Servo E-AT/PS, Servo Oil, Ford
SHELL CASSIDA CHAIN OIL 150, Special Oil, Shell,
SHELL Comptela S 46, Compressor Oil
SHELL Corena D46, Compressor Oil
SHELL Diala DX, Electroinsulating Oil
SHELL DONAX TA CHECK, Gear Oil, Shell,
SHELL Gear Oil EP 80W, Gear Oil
SHELL Helix Plus 10W40, Motor Oil, for Diesel, Synthetic
SHELL HELIX PLUS SAE 10W-40, Motor Oil, Shell,
SHELL Helix Plus VB 10W-40, Motor Oil, Synthetic
SHELL HELIX SAE 15W-40, Motor Oil, Shell,
SHELL Helix Ultra 5W-40, Motor Oil, Synthetic
SHELL Helix Ultra AC 5W-30, Motor Oil, Synthetic
SHELL Hydrol DO 46, Hydraulic Oil
SHELL Omala 100, Gear Oil
SHELL ONDINA P 68, Special Oil, Shell,
SHELL Rimula 15 W-40 used, Motor Oil
SHELL Rimula Super 15W-40, Motor Oil
SHELL Rimula Ultra 10W40 used 52000km, Motor Oil, 52000/9000km
SHELL Rimula Ultra 10W-40, Motor Oil
SHELL Spirax AX 80W90, Axle Oil, Iveco Truck
SHELL Spirax GX 80W90, Gear Oil, Iveco Truck
SHELL SPIRAX MA SAE 80W, Gear Oil, Shell,
SHELL SUPER 2T, Motor Oil, Shell,
SHELL Tellus 46, Hydraulic Oil
SHELL Tellus S 32, Hydraulic Oil
Shell Tellus S32; Hydraulic oil; Shell
SHELL Tellus S-46 New, Hydraulic Oil
SHELL Tellus S-46 Used in Engine, Hydraulic Oil
SHELL Tellus S-46 Used in Press, Hydraulic Oil
SHELL Tellus T 32, Hydraulic Oil, HVLP
SHELL Thermia Oil B, Heat Transfer Medium
SHELL Turbo T 32, Turbine Oil
SHELL Turbo T 46, Turbine Oil
SHELL XHVI 5.2, Base Oil, Shell,
Siempezka PL, Fuel Oil, Slovtherov

IR Lubricants, 2081 FTIR spectra

<http://www.ir-spectra.com/>

IR Lubricants

Silver Synt 10W40 SL/CF, Motor Oil, Global Oil, Belgium
Stihl Two Stroke for Saw, Motor Oil, Stihl, Two Stroke for Saw
SUNOCO Sungrind Range, Cutting Fluid
SUPER DIESEL SAE 15W40, Motor Oil, MOL,
SUPER MOTOR OIL HDC-IL TURBO 15W40, Motor Oil, VECO,
SWEPCO 201, Gear Oil, SWEPCO,
Synpower 5W40 sample A, Motor Oil, Valvoline, Fully Synthetic
SYNPOWER 5W40, Motor Oil, Valvoline,
SYNPOWER 5W50, Motor Oil, Valvoline,
SYNTEX GRADE 10W40, Motor Oil, DENICOL,
TB-32; Turbine Oil used; Paramo
TB-32; Turbine Oil; Paramo
TB-46; Turbine Oil used
TB-46; Turbine Oil; Paramo
TELEX E, Hydraulic Oil, REPSOL,
TEXACO Alcor DD 22, Hydraulic Oil
TEXACO Cetus S 46, Compressor Oil
TEXACO Geartex EP A 85W90, Gear Oil
TEXACO Havoline 10W40, Motor Oil
TEXACO Havoline 505.01, Motor Oil
TEXACO Havoline 5W30, Motor Oil
TEXACO Havoline 5W40, Motor Oil, Synthetic
TEXACO Havoline Extra 10W40, Motor Oil
TEXACO Meropa 100, Gear Oil
TEXACO Meropa 150, Industrial Oil
TEXACO Meropa 220, Gear Oil
TEXACO Meropa 320, Used 14 days, water contents, Gear Oil
TEXACO Meropa 320, Gear Oil
TEXACO Rando HD 32, Hydraulic Oil
TEXACO Rando HD 46, Hydraulic Oil
TEXACO Rando HD 68, Industrial Oil, pumps, compressors, bearings
TEXACO Rando HD-Z 22, Hydraulic Oil
TEXACO Rando HD-Z 46, Hydraulic Oil
TEXACO Rando HD-Z 68, Hydraulic Oil
TEXACO Regal Premium 46, Turbine Oil, Marine
TEXACO Regal Premium EP 32, Turbine Oil
TEXACO Regal R&O 32, Turbine Oil, industrial, multipurpose
TEXACO Regal R&O 46, Turbine Oil, industrial, multipurpose
TEXACO Regal R&O 68, Turbine Oil, industrial, multipurpose
TEXACO Ursa Super LA 15W40, Motor Oil
TEXACO Ursa Super TD, Motor Oil
TOTAL Azolla AF 46, Hydraulic Oil
TOTAL Azolla ZS 46, Hydraulic Oil
TOTAL Azolla ZS 68, Hydraulic Oil
TOTAL Carter EP 220, Gear Oil
TOTAL Cortis MS 460, Industrial Oil, rollings
TOTAL Cortis, Industrial Oil
TOTAL Equivis 46 ZS, Hydraulic Oil
TOTAL Fluide ATX, Gear Oil, automated transmission
TOTAL Nateria MJ 40, Motor Oil, Gas Engines
TOTAL Preslia 32, Turbine Oil
TOTAL Preslia 46, Turbine Oil
TOTAL Quartz 10W-40, Motor Oil

IR Lubricants, 2081 FTIR spectra

<http://www.ir-spectra.com/>

IR Lubricants

TOTAL Quartz 5000 15W40, Motor Oil
TOTAL Quartz 5W40, Motor Oil, Synthetic
TOTAL RUBIA 4400 15W40, Motor Oil, Total,
TOTAL Rubia S 40, Motor Oil, diesel engines
TOTAL RUBIA TIR 6400 15W40, Motor Oil, Total,
TOTAL RUBIA XT 15W40/1996, Motor Oil, Total,
TOTAL RUBIA XT A 15W40, Motor Oil, Total,
TOTAL Seriola 1510, Heat Transfer Medium
TOTAL TM 80W-90, Gear Oil
TOTAL Torilis 2500, Process Oil
TOTAL Torilis 6200, Process Oil
Trans Premium 150, Gear Oil, Statoil
Trans Premium 220, Gear Oil, Statoil
Trans Premium 220, Gear Oil, Statoil
TRANSMISSION TM 80W90, Gear Oil, Total,
Traxon 80W-90, Gear Oil, Petro Canada
TRENNFIT 3102, Forming Oil, FRANK,
TRENNFIT, Forming Oil, FRANK,
TURBINE OIL 46, Turbine Oil, Petrochema,
TURBO MULTI-SERVICE 15W40, Motor Oil, CANADIAN TURBO,
TURBO SELECT 15W40, Motor Oil, DENICOL,
Unifarm 10W-30;Tractor Lubricant used;Esso
Unifarm 10W-30;Tractor Lubricant;Esso
Urania LD 15W40, Motor Oil, Iveco Truck
Vacmul 30A40;Cutting Fluid, used;Exxon
Vacmul 30A40;Cutting Fluid;Exxon
Valar Atea, Hydraulic Oil
Valar egida XQ 15W-40, Motor Oil
Valar HLP 46, Hydraulic Oil
VAT ATF TYPE D/1999, Gear Oil, VAT OIL,
VAT SYNTECH MOTOR OIL 10W40, Motor Oil, VAT OIL,
VECO TURBO DIESEL 15W40, Motor Oil, VECO,
VS MAX SAE 15W40, Motor Oil, SELENIA,
Way Lubricant X 220, Metalworking Fluid, Chevron
Way Lubricant X 68, Metalworking Fluid, Chevron
ZHS-Servo, Hydraulic and Steering Oil, Startol

IR Lubricants, 2081 FTIR spectra
<http://www.ir-spectra.com/>