



# **IS NIR Spectra**

## **Complete collection (6049 spectra)**

This collection includes all NIR spectra collected within this project. Those are organic and inorganic chemicals, pharmaceutical excipients, drugs, polymers, additives, coatings, dyes, pigments, toxic substances, hazardous substances like explosives, agrochemicals and other types of materials. This collection should be chosen by customers with broad range of analytical applications.

The near infrared spectral library „IS NIR Spectra“ is a unique collection of 6049 NIR (near infrared) digital spectra. The spectra have been collected in range 4200-10000 cm<sup>-1</sup> (=2.4 – 1 micrometers).

„IS NIR Spectra“ is the biggest and newest NIR library available on the market. This spectral library has been designed for users who use their NIR spectrometer not only for quantitative and chemometric analysis, but also want to identify unknown samples or just want to use the library as reference spectra database.

IS means “Intergrating Sphere” and at the same time “Identify Sample”.

Technical parameters:

The „IS NIR Spectra“ were collected using the FT-NIR Spectrometer equipped with Integrating Sphere Accessory.

Resolution: standard 4cm<sup>-1</sup>

Customized resolution: upon special request the spectra can be delivered in customized resolution to reach better correspondence with spectra collected by customer NIR instrument

Spectral range in wavenumbers units: 4200-10000 cm<sup>-1</sup>

Spectral range in wavelength units: 2.4 – 1 micrometers

Only one spectrum/compound covering the mentioned spectral range

Typical text information about the sample includes (when available):

Compound name including synonyms

CAS number

Chemical formula

Molecular weight

Additional information about the sample like manufacturer and comments.

Download a free NIR library of 30 NIR spectra to check the compatibility of you instrument with the „IS NIR Spectra“.

Samples were selected and NIR spectra collected 2010-2012.

---

**IS NIR Spectra – Complete Collection**

6049 NIR spectra

Copyright © NICODOM 2012

NICODOM Ltd., Hlavni 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

[ftirspectra@gmail.com](mailto:ftirspectra@gmail.com)

This library is available as digital searchable library compatible with your spectroscopic software. Visit the „Download“ section of our webpage to select the correct library format. The spectra are available in wavenumbers (cm-1) units or wavelength (micrometer) units. Also the library can be delivered with SEARCH software or in \*.pdf format (if you do not have a spectroscopic software).

The near infrared spectral library „IS NIR Spectra“ can be purchased as a complete collection of over 6049 spectra or as following application specific sublibraries.

Ordering Code	Library Name	Number of Spectra
ISNNCM-X	IS NIR Complete Collection	6049
IS NIR Sublibraries*:		
ISNNST-X	IS NIR Starter	1043
ISNNPH-X	IS NIR Pharma, Drugs	1340
ISNNPO-X	IS NIR Polymers, Additives	1740
ISNNHA-X	IS NIR Hazardous, Toxic	1718
ISNNDY-X	IS NIR Dyes, Pigments	1005

\* all sublibraries spectra are included in the complete collection ISNNCM-X and have the same technical parameters

---

### IS NIR Spectra – Complete Collection

6049 NIR spectra

Copyright © NICODOM 2012

NICODOM Ltd., Hlavni 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971


<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

[ftirspectra@gmail.com](mailto:ftirspectra@gmail.com)

## Example Spectra

Pictures: NIR spectrum of Sucrose (powder), C<sub>12</sub>H<sub>22</sub>O<sub>11</sub>, CAS: 57-50-1, in cm<sup>-1</sup> units (range 4200-10000 cm<sup>-1</sup>), the second spectrum in micrometers (2.4 – 1 micrometers), both resolution 4cm<sup>-1</sup>, Copyright Nicodom 2012



**IS NIR Spectra – Complete Collection**

6049 NIR spectra

Copyright © NICODOM 2012

NICODOM Ltd., Hlavní 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

[ftirspectra@gmail.com](mailto:ftirspectra@gmail.com)

## List of Spectra /Compound Names

1-Chloro-4-nitrobenzene  
4-Nitroaniline  
p-Nitrophenol  
4'-Methoxyacetophenone  
4-Methoxybenzoic acid  
p-Dimethylaminobenzaldehyde  
p-Nitrobenzyl bromide  
a-Cyclodextrin hydrate  
4-Nitrophenylhydrazine  
4'-Nitroacetophenone  
Quinoline Yellow Lake  
Terephthalic acid  
p-Dinitrobenzene  
Pentadecanoic acid  
Pyrrole-2-carboxaldehyde  
2-Thiophenecarbonitrile  
N,N-Diethylethylenediamine  
4-Picoline-N-oxide  
N,N-Diethylethanolamine  
Benzyl bromide  
Sodium Picosulfate  
Ethylbenzene  
Styrene  
2,4,6-Triaminopyrimidine  
Benzyl chloride  
1,4,7,10-tetraazacyclododecane tetrahydrochloride  
Benzylamine  
Benzonitrile  
4-Cyanopyridine  
Benzyl alcohol  
Benzaldehyde  
3-Cyanopyridine  
N-Methylcyclohexylamine  
10-{N-[3,3-bis(phosphonic acid)propyl]aminocarboxymethyl}-1,4,7,10-tetraazacyclododecane-1,4,7-triacetic acid  
Phenylhydrazine  
Thioanisole  
9,10-Bis(phenylethynyl)anthracene  
2,2-Dimethyl-1,3-dioxolane-4-methanol  
Naugard 445

---

IS NIR Spectra – Complete Collection

6049 NIR spectra

Copyright © NICODOM 2012

NICODOM Ltd., Hlavni 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

[ftirspectra@gmail.com](mailto:ftirspectra@gmail.com)

3-Hydroxybenzaldehyde  
Aranox  
Hexamethylenetetramine  
L-Lysine Monohydrochloride  
Triphenyl phosphite  
1,4,7-Oxadiazonane dihydrobromide  
Triethanolamine zirconate  
N,N'-bis(ethylphosphonic Acid Diethylester)Piperazine  
Basic Brown 1 | Bismarck Brown  
Rhodium, carbonyl(2,4-pentanedionato-O,O')(tricyclohexylphosphine)-, (SP-4-2)- (9CI)  
2,6-Dichloroquinone-4-chloroimide  
3-Phenyl-1-propyne  
10-{[(2,2-Bisphosphonoethyl)hydroxyphosphoryl]methyl}-1,4,7,10-tetraazacyclododecane-  
1,4,7-triacetic acid  
4,4'-Methylenebis(N,N-dimethylaniline)  
Nonox ZA  
Diaminodiphenylmethane  
Diphenylmethane  
Phenyl Ether  
ASM 4010  
Ferrocene, 1-carboxy-1'-(diphenylphosphinothioyl)-  
Topanol M  
4-Chloro-7-nitrobenzo-2-oxa-1,3-diazole  
Acetoacetanilide  
Phenylbiguanide  
3,5-Di-tert-butylcatechol  
1,3-Diphenylurea  
1,3-Diphenyl-2-thiourea  
Hydrazinium hydroxide  
Glybenclamide  
Ferrocene  
Maprotiline  
Tripropylamine  
Triethanolamine  
Levetiracetam  
2-(Dibutylamino)ethanol  
Trihexylamine  
endo-3-Bromo-d-camphor  
3,9-Dibromocamphor  
beta-Nitrostyrene  
6-Chloro-1-hexyne  
n-(8-quinolyl)-p-toluenesulfonamide  
2-Ethylhexyl acrylate  
D-(-)-Arabinose

---

**IS NIR Spectra – Complete Collection**

6049 NIR spectra

Copyright © NICODOM 2012

NICODOM Ltd., Hlavní 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

[ftirspectra@gmail.com](mailto:ftirspectra@gmail.com)

Methyl cinnamate  
D-(-)-Lactic acid  
N-Methylisatoic anhydride  
Trans-Stilbene  
Azobenzene  
Ethyl cinnamate  
2-Phenylethyl acetate  
2-(N-Cyclohexylamino)ethanesulfonic acid; CHES  
Dibenzylamine  
2-Phenoxyethyl isobutyrate  
b-Bromostyrene  
N-Benzylmethylamine  
5,7,7,12,14,14-hexamethylcyclam  
Ethylaniline  
Phenyl Isocyanate  
Phenyl isothiocyanate  
N-(2-Hydroxyethyl)piperazine  
2-Phenylacetamide  
Phenylacetic acid  
N,N-Dimethylbenzylamine  
Acetanilide  
1-Phenyl-2-thiourea  
Potassium thioacetate  
Lacidipine  
4'-Methylacetanilide  
Acetaminophen  
4'-Aminooxanilic acid  
Diethyl Iodomethylphosphonate  
3-Acetyl umbelliferone  
Tetrabutylammonium cyanide  
Anethole  
Cinnamyl alcohol  
Oxalate Iron dihydrate  
2-Ethylhexylamine  
Tinuvin 1130  
p-Tolunitrile  
4-Chlorobenzaldehyde  
5-Ethyl-2-methylpyridine  
4-Bromoanisole  
3-Dodecylthiophene  
p-Anisidine  
4-Cyanobenzaldehyde  
Cresyl Fast Violet  
Sokrat 6492

---

**IS NIR Spectra – Complete Collection**

6049 NIR spectra

Copyright © NICODOM 2012

NICODOM Ltd., Hlavni 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

[ftirspectra@gmail.com](mailto:ftirspectra@gmail.com)

Methyl cyanoacetate  
Ethyl bromoacetate  
Ethyl chloroacetate  
Ethyl crotonate  
Diethyl-3-oxoglutarate  
Diethyl malonate  
Ethyl butyrate  
Ethyl cyanoacetate  
N-methyldiethanolamine  
2-Oxohexamethyleneimine  
N,N'-bis-(3-Aminopropyl)Ethylenediamine  
Dibutyl maleate  
O-(N-succinimidyl)-N,N,N',N'-tetramethyluronium tetrafluoroborate  
Eriochrom Blue SE  
Carbendazim  
Butyl laurate  
DL-Citronellol  
Risperidone Hydrochloride Hemipentahydrate  
Cholesteryl Caproate  
Ethyl caprylate  
Ethyl dodecanoate  
Quinhydrone  
p-Dibromobenzene  
4-Bromochlorobenzene  
4-Bromoaniline  
4-Bromophenol  
4-Chlorotoluene  
p-Cresol  
p-Dichlorobenzene  
4-Chloroaniline  
4-Chlorophenol  
p-Toluidine  
1,4-Phenylenediamine  
p-Benzoquinone  
Dimethyl succinate  
Aminomethylphosphonic acid  
Dibutyltin diacetate  
Diethyl (2-oxopropyl)phosphonate  
Diethyl dithiophosphate, ammonium salt  
4-Vinylcyclohexene dioxide  
Epichlorhydrin  
Glycidyl methacrylate  
1,2-Dibromoethane  
1-Bromopropane

---

**IS NIR Spectra – Complete Collection**

6049 NIR spectra

Copyright © NICODOM 2012

NICODOM Ltd., Hlavni 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

[ftirspectra@gmail.com](mailto:ftirspectra@gmail.com)

1,5,8,12-tetrakis[4,6-bis(n-butyl-n-1,2,2,6,6-pentamethyl-4-piperidylamino)-1,3,5-triazin-2-yl]-1,5,8,12-tetraazadodecane  
Allyl chloride  
1,2-Dichloroethane  
2-Chloroethanol  
Allylamine  
Acrylonitrile  
Chloroacetonitrile  
Ethylenediamine  
Allyl alcohol  
Propargyl alcohol  
Chloroacetaldehyde  
Ethylene glycol  
Lead Stearate  
Bismuththiol  
2-Acetylpyrrole  
Methyl formate  
Taurine = 2-Aminoethanesulfonic acid  
4-Chlorophenylhydrazine hydrochloride  
2-Methyl-2,4-pentanediol  
Betaine  
Surfynol® 61  
DL-Thioctic acid  
2-Bromopentane  
6-Methoxy-1-tetralone  
2-Pentanone  
1,3-Butanediol  
2-Cyanoacetamide  
Crotonic Acid  
Beta Alanine  
3-Mercaptopropionic acid  
Chlorodiphenylphosphine  
1-Methoxy-2-propanol  
2-Dimethylaminoethylamine  
2-Dimethylaminoethanol  
Diethyl benzylphosphonate  
Vinyl acetate  
4-Methyl-2-pentanone  
Diisopropylamine  
Biuret  
Diisopropyl ether  
Acetic anhydride  
g-Valerolactone  
Succinic anhydride

---

**IS NIR Spectra – Complete Collection**

6049 NIR spectra

Copyright © NICODOM 2012

NICODOM Ltd., Hlavni 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

[ftirspectra@gmail.com](mailto:ftirspectra@gmail.com)

Maleic anhydride  
Propylene carbonate  
4-(4-Nitrobenzyl)pyridine  
m-Cresol  
3-Chloroaniline  
3-Chlorophenol  
m-Toluidine  
1,3-Benzenediol  
2,4-Lutidine  
2,6-Lutidine  
8,11-Dibenzyl-1,4-bis(Methylenephosphonic Acid) Cyclam  
Dimethyl malonate  
Dowanol PMA  
3,5-Dimethylphenol = Xylenol  
3,5-Dimethylaniline  
Phloroglucinol, 1,3,5-Benzenetriol  
1,3,5-Trimethylhexahydro-1,3,5-triazine  
2,4,6-Collidine  
Melamine = sym-Triaminotriazine  
2,4,6-Trihydroxy-1,3,5-triazine  
Diisobutylketone  
Toluene  
4-Methylpyridine  
Cyclohexylamine  
Cyclohexanol  
Cyclohexanone  
Phenol  
3-Picoline = 3-Methylpyridine  
3-Hydroxypyridine  
4-Methylmorpholine  
2-Bromopyridine  
2-Methylpiperidine  
2-Picoline  
2-Methylpiperazine  
2-Aminopyrimidine  
Di(ethylene glycol) distearate  
Allylacetone  
1-Valeric acid  
3-Dimethylaminopropylamine  
1-Allyl-2-thiourea  
Propyl acetate  
1,3-Dibromopropane  
1-Bromobutane  
1-Chlorobutane

---

**IS NIR Spectra – Complete Collection**

6049 NIR spectra

Copyright © NICODOM 2012

NICODOM Ltd., Hlavní 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

[ftirspectra@gmail.com](mailto:ftirspectra@gmail.com)

1-Bromo-3-chloropropane  
1,3-Diaminopropane  
Malononitrile  
3-Hydroxypropionitrile  
2-Methylaminoethanol  
2-Methoxyethylamine  
Ethylene glycol monomethyl ether  
Diethylamine  
Ethyl formate  
Pyrrole  
Tetrahydrofuran  
Furan  
Thiophene  
Sodium tauroglycocholate  
Acetylacetone  
Succinamide  
Succinic acid  
Maleic acid  
Fumaric acid  
N,N,N',N'-Tetramethylethylenediamine  
N,N'-Methylenebisacrylamide  
Isopropyl myristate  
Adwawax 280  
2-Heptanone  
Sorbic Acid  
Isoamyl formate  
1,4-Dibromobutane  
Irganox 1520  
Amylamine  
2-Butene-1,4-diol (cis and trans)  
2-Butyne-1,4-Diol  
Ethylene glycol dimethyl ether  
2-(Ethylamino)ethanol  
Ethylene glycol monoethyl ether  
Cyclohexane  
Cyclohexene  
Diethylenediamine  
Pyridine  
Piperidine  
6-Methyl-5-hepten-2-one  
Glutaric acid  
Glyceryl Stearate  
2-Ethoxyethyl acetate  
Pimelic acid

---

**IS NIR Spectra – Complete Collection**

6049 NIR spectra

Copyright © NICODOM 2012

NICODOM Ltd., Hlavni 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

[ftirspectra@gmail.com](mailto:ftirspectra@gmail.com)

Sebacic Acid  
1,5-Dibromopentane  
1-Bromohexane  
Tetrabutylammonium chloride  
1,5-Pentanediol  
n-Butyl vinyl ether  
Bis(pentamethylcyclopentadienyl)titanium dichloride  
Diethylenetriamine  
N-(2-hydroxyethyl)ethylenediamin  
Diethanolamine  
2,2'-Dichloroethyl ether  
Diethylene glycol  
2,2'-Thiodiethanol  
Adipoyl chloride  
Hydroxypivalyl hydroxypivalate  
Octane  
Trioctylamine  
1-Heptanol  
Heptaldehyde  
1,3-Butanediol diacetate  
Ethylene glycol monobutyl ether  
1,5-Cyclooctadiene  
1-Bromooctane  
Octylamine  
2,2,6,6-Tetramethyl-3,5-heptadione  
1-Octanol  
Ethylene Diglycol Monoethyl Ether  
5-Chlorovaleric acid  
Diethylene glycol dimethyl ether  
Quetiapine Fumarate  
Sodium Cetylsulfate  
1-Tetradecene  
Nonanoic acid  
Alpha,alpha-Diphenyl-L-prolinol  
2-Butoxyethyl acetate  
3-Hydroxy-2-methylpyridine  
2-Pyridinecarboxaldehyde  
Benzylcyclene  
Nonylamine  
Triethylenetetramine  
4-Acetylpyridine  
4-Dimethylaminopyridine  
Nitrocyclohexane  
2-Acetylpyridine

---

**IS NIR Spectra – Complete Collection**

6049 NIR spectra

Copyright © NICODOM 2012

NICODOM Ltd., Hlavní 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

[ftirspectra@gmail.com](mailto:ftirspectra@gmail.com)

Glycol Bis(Hydroxyethyl) Ether  
4-Bromobenzaldehyde  
1-Decanol  
N-Decyl aldehyde  
2-(2-Butoxyethoxy)ethanol  
Undecylenic acid  
Methyl Palmitate  
Dodecane  
1-Dodecene  
1-Undecanol  
4-Nitropyridine N-oxide  
10-Undecen-1-ol  
1,10-Decanediol= Decamethylene glycol  
[5-[[4-[2-(5-Ethyl-2-pyridinyl)ethoxy]phenyl]methyl]-2,4-] thiazolidinedione hydrochloride  
1-Dodecanol  
Dodecyl aldehyde  
5,6,7,8-Tetrahydro-2-naphthol  
2-(2-Hexyloxyethoxy)Ethanol  
Tetraethylenglycol  
Methyl stearate  
2,6-Dihydropyridine Ditosylate  
Dodecyl acetate  
Cetyltrimethylamine  
1-Bromotetradecane  
Diethylene Glycol Dibutyl Ether  
Oleoyl chloride  
Oleic acid  
Letrozole  
1-Bromohexadecane  
1-Octadecene  
1-Bromooctadecane  
1-Octadecanol  
2,6-Diacetylpyridine  
1-Phenyl-1-hexyne  
Octadecyl isocyanate  
Guanidine  
3,4-Dimethoxyacetophenone  
Sodium pyruvate  
3-(N-Morpholino)propanesulfonic acid; MOPS  
S(+)-Bicalutamide  
R(-)-Bicalutamide  
4,4'-Dimethyl-2,2'-dipyridyl  
Dosulepine  
3-(Cyclohexylamino)propane-1-sulfonic acid; CAPS

---

**IS NIR Spectra – Complete Collection**

6049 NIR spectra

Copyright © NICODOM 2012

NICODOM Ltd., Hlavní 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

[ftirspectra@gmail.com](mailto:ftirspectra@gmail.com)

Chlorprothixene  
Caryophyllene oxide  
PAR, 4-(2-pyridylazo) resorcinol  
Dithranole  
Scopolamine .HBr  
1-Acetyl-2-phenylhydrazine  
Triethyl 1,4,7,10-tetraazacyclododecane-1,4,7-triacetate  
Glyoxalbis(2-hydroxyanil)  
Diethyl 1,4-dihydro-2,6-dimethyl-3,5-pyridinedicarboxylate  
Diethyl 2,6-dimethyl-3,5-pyridinedicarboxylate  
Benzyltrimethylammonium dichloriodate  
2-Methyl-3-Butyne-2-ol  
Chlorotriphenyllead  
Aminomethylphosphonic acid ethyl ester  
Bromophenol Blue  
Bromocresol Purple  
Dithiaden  
Pentaerythritol  
2-Butyl-2-ethyl-1,3-propanediol  
Bis(tetramethylcyclopentadienyl)titanium dichloride  
Triphenyl phosphate  
3,7-Dimethyl-1,6-octadien-3-yl acetate  
Bis(1,2,3-trimethyl-cyclopentadienyl) titanium (IV) dichloride  
Triisopropyl phosphite  
{[Hydroxy(phenyl)phosphoryl]methyl}(methyl)amino]methyl}phosphonic acid  
Methyl dichloroacetate  
Bis(2,4,6-trichlorophenyl)ethanedioate  
1-Amino-2-naphthol-4-sulfonic acid  
Dodecyl gallate  
2,5-Dibromo-3-hexylthiophene  
Tetrachlorophthalic anhydride  
Bis(trimethylsilyl-tetramethyl-cyclopentadienyl) titanium (IV) dichloride  
Tetraphenoxysilane  
1-Naphthol-5-sulfonic acid  
Dexamethasone Acetate  
1,3-Dichloro-1,4-naphthoquinone  
Diethylphthalate  
Guanosine  
Dinitrosoresorcinol  
Cinchonine  
Cholesteryl Heptanoate (Enanthate)  
Cholesteryl Pelargonate  
N-(Hydroxymethyl)phthalimide  
Cholesteryl Caprylate

---

**IS NIR Spectra – Complete Collection**

6049 NIR spectra

Copyright © NICODOM 2012

NICODOM Ltd., Hlavni 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

[ftirspectra@gmail.com](mailto:ftirspectra@gmail.com)

2-Naphthol-6,8-disulfonic acid  
Ethanox 398  
Tris(hydroxymethyl)aminomethane hydrochloride  
Ammonium ferric citrate green  
Phenyl Salicylate  
4-Acetamidophenyl salicylate  
Dibutyltin Dilaurylmercaptide, 16,3-17,5 % Sn  
Diethyl(3-bromopropyl)phosphonate  
Ethyl salicylate  
Maltol  
Tetrachlorobenzoquinone = Chloranil  
Tribromophenol  
2-Methylheptanoic acid  
Ionox 220  
o-Toluic acid  
2-Chlorobenzoic acid  
Anthranilic acid = 2-Aminobenzoic acid  
2'-Hydroxyacetophenone  
4-(2,4-Dinitroanilino)phenol  
2-Acetylfuran  
2,4-Dinitrophenylhydrazine  
1-Amino-7-sulfonaphthalene  
Methyl pyrrole-2-carboxylate  
Methylsalicylate  
(R)-(+)-N,N-dimethyl-1-(1-naphthyl)ethylamine  
Bis(tetramethylcyclopentadienyl)zirconium dichloride  
4,4'-dimethoxybenzoin  
Cyasorb UV 9  
Diethyl 9,10-dihydro-9,10-ethanoanthracene-11,11-dicarboxylate  
2,6-Pyridinedimethanol  
Benzophenone  
1,2,3,4-Tetrahydronaphthalene  
Isoquinoline  
1-Naphthylamine-6-sulfonic acid  
Uric acid sodium salt  
o-Dianisidine, tetrazotized  
2,2'-Biquinoline= 2,2'-Diquinolyl  
o-Tolidine  
Trimethyl(phenylethynyl)tin  
2-[(1-Benzyl-4-piperidyl)methyl]-5,6-dimethoxy-2,3-dihydroinden-1-one  
2,2'-(Phenylimino)diethanol  
Anthracene  
3,4-Dimethoxybenzaldehyde  
b-Naphthoxyacetic acid

---

**IS NIR Spectra – Complete Collection**

6049 NIR spectra

Copyright © NICODOM 2012

NICODOM Ltd., Hlavni 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

[ftirspectra@gmail.com](mailto:ftirspectra@gmail.com)

2-Benzyl-4-chlorphenol  
Dibenzoylmethane  
4-Ethyl hydroxybenzoate  
Benzyl Benzoate  
Diethylene glycol dibenzoate  
(4-Nitro-benzyl)-phosphonic acid  
Piperonal  
Isosafrole  
Dimethyl Terephthalate  
N-acetylorthotoluidine  
Sodium Toluenesulfonate  
Indole  
Pneumax DM  
Cyclopentadienylmanganese tricarbonyl  
Dichloro(1,5-cyclooctadiene)platinum(II)  
Catechol  
Palladium, di-m-chlorobis[(1,2,3-h)-2-methyl-2-propen-1-yl]di-  
2,4-Dichlorophenol  
Bis(tris(acetylacetonato)titanium(IV)) hexachlorotitanate(IV)  
Cyclopentanone  
Rhodium, di-m-chlorobis[(1,2,5,6-h)-1,5-cyclooctadiene]di-  
2-Imidazolidone  
Diphenylphosphinoferrocene  
3-tert-butyl-4-hydroxyanisole  
9-Anthracenecarbonitrile  
Dibenzosuberone  
Dichloro(1,5-cyclooctadiene)palladium(II)  
2,4-Dinitrotoluene  
3-Ethoxy-4-hydroxybenzaldehyde = Ethylvanillin  
Vanillic acid  
6-Benzylaminopurine  
Triethylamine  
Palladium, dibromo[(1,2,5,6-h)-1,5-cyclooctadiene]-  
Trimethyl phosphite  
3-Aminobenzenesulfonic acid  
1,1'-Bis(diphenylphosphino)ferrocene  
Benzethonium chloride  
Sulfanilic acid  
N,N-Dimethylaniline  
9,10-Dicyanoanthracene  
Propyl Gallate  
2-Chloro-4-nitroaniline  
3'-Nitroacetophenone  
3-Nitrobenzoic acid

---

**IS NIR Spectra – Complete Collection**

6049 NIR spectra

Copyright © NICODOM 2012

NICODOM Ltd., Hlavni 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

[ftirspectra@gmail.com](mailto:ftirspectra@gmail.com)

Acid Red 118 | Rybanyl Red NB  
Basic Red 14 | Basacryl Red X 4G  
Pigment Red 176  
Reactive Brown 11  
Solvent Blue 37 | Lampronol Blue BR  
Solvent Blue 66 | Brilliant Blue Base SM  
Solvent Green 20 | Waxoline Green 5G  
FD&C Yellow 6 Aluminium Lake  
m-Nitroacetanilide  
Malonaldehyde bis(diethyl acetal)  
Pigment Orange 36  
Reactive Blue 49  
Diphenylamine  
Disperse Yellow 82  
Triethyl orthoformate  
Triethyl phosphite  
1,4,7-Tris(tert-butoxycarbonylmethyl)-1,4,7,10-tetraazacyclododecane  
Amberlite IRA 400  
Benzalacetone  
Dioctyl sebacate  
Basic Blue 41  
N,N'-Dibenzylethylenediamine diacetate  
N-Acetyl-1,4-phenylenediamine  
4-Acetamidobenzaldehyde  
Pigment Yellow 62  
Food Red 16 | Sudan Red G  
1,4-Piperazinediethanol  
N-(2-hydroxyethyl)anilin  
1-Hydroxy-2-phenoxyethane  
Cetylpyridinium chloride monohydrate  
3-(Chloromethyl) heptane  
4-Hydroxybenzaldehyde  
4-Methoxybenzaldehyde  
Diethyl succinate  
P-Aminophenol  
Hydroquinone  
Maleic Hydrazide  
Propionaldehyde  
Choline  
Diacetone Alcohol  
3-Methyl-1-butanol  
Acetylacetone  
Bis[(pentamethylcyclopentadienyl)dichloro-rhodium]  
Succinimide

---

**IS NIR Spectra – Complete Collection**

6049 NIR spectra

Copyright © NICODOM 2012

NICODOM Ltd., Hlavni 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

[ftirspectra@gmail.com](mailto:ftirspectra@gmail.com)

Propionic anhydride  
Paraldehyde  
Ethyl caproate  
Butyraldehyde  
Levulinic acid  
Dioctyl adipate  
Butyl acetate  
1,4-Dioxane  
Isoamyl acetate  
Acid Blue 193  
Glycerol Stearate  
Ethidium bromide  
2-Octanol  
Sumilizer GS  
Azelaic acid  
Adipic acid  
1,6-Diaminohexane  
Octanal  
2-(2-butoxyethoxy)ethyl acetate  
Decane  
N,N',N''-trisubstituted triazacyclononane with methyl(2-carboxyethyl)phosphinic acid  
Dodecylamine  
Stearamide  
Octadecylamine  
N-Isobutyryl-L-Cysteine  
Methanesulfonyl chloride  
2-Amino-2-methyl-1-propanol  
Losartan Potassium  
Ethyl [bis(2,2,2-trifluoroethoxy)phosphinyl]acetate  
Pitofenone Hydrochloride  
Ferrocene, 1-(diphenylphosphino)-1'-formyl-  
Thymolphthalein  
Cholestane-3b,5a,6b-triol  
Fenpiverinium bromide  
2-(1h-benzotriazole-1-yl)-1,1,3,3-tetramethyluronium tetrafluoroborate  
3-[6-(4-Methoxyphenoxy)hexyl]-thiophene  
Sephadex A 25  
Tris(hydroxymethyl)nitromethane  
Sucrose octaacetate  
N-phthaloyl-aminomethyl-(ethoxycarbonyl)phosphinic acid  
Indigo blue  
Neopentyl glycol  
Tetramethylene sulfone  
Solvent Yellow 98 | Hostasol Yellow 3G

---

**IS NIR Spectra – Complete Collection**

6049 NIR spectra

Copyright © NICODOM 2012

NICODOM Ltd., Hlavni 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

[ftirspectra@gmail.com](mailto:ftirspectra@gmail.com)

Tributyl phosphate  
Phenanthroline  
5,5-Dimethyl-1,3-cyclohexanedione  
2,4,7,9-Tetramethyl-5-Decyne-4,7-Diol  
Acetone oxime  
Potassium Acetate  
Sodium Acetate  
Titanocene dichloride  
Ferrocenecarboxylic acid  
Ethynylferrocene  
Acetylferrocene  
Ferrocene, (methoxycarbonyl)- (9CI)  
N,N-Dimethylacetamide  
Ferrocenemethanol  
1,1'-Diacetylferrocene  
Chloramine B hydrate  
Diphenyl sulfone  
Chloramine T  
3-Nitrobenzenesulfonic acid sodium salt  
b-Pinene  
N-Bromosuccinimide  
N-Chlorosuccinimide  
2,6-Di-tert-butyl-p-cresol  
2,6-Di-tert-butylphenol  
Dinitrostilbenedisulfonic acid  
Saccharin Sodium  
3,4,9,10-Perylenetetracarboxylic dianhydride  
Ferrocenecarboxamide  
1,4-Diaminoanthraquinone  
1,4,5,8-Naphthalenetetracarboxylic acid  
Vat Yellow 2 | Ostanthren Yellow GC  
2-Ethynyl-6-methoxynaphthalene  
Bis(cyclopentadienyl)zirconium dichloride  
Acid Blue 1 | Triacid Blue V  
Patent Blue V  
Bis(cyclopentadienyl) titanium bromide  
1,1'-Ferrocenedicarboxylic acid  
Pigment Yellow 191  
1,5-Diaminoanthraquinone  
2,4,6-Trinitrobenzoic acid  
Aripiprazole  
Naphthionic Acid Sodium Salt  
Alizarin S, Sodium Alizarin Sulfonate  
Methyl 4-aminobutyrate hydrochloride

---

**IS NIR Spectra – Complete Collection**

6049 NIR spectra

Copyright © NICODOM 2012

NICODOM Ltd., Hlavni 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

[ftirspectra@gmail.com](mailto:ftirspectra@gmail.com)

Bromomethyl methyl ether  
Quinine Hydrochloride  
Quinine  
3-(Trimethylsilylethynyl) thiophene  
Anthraquinone-2-sulfonic acid, sodium salt, monohydrate  
Dimethyl phthalate  
a-Naphthyl red hydrochloride  
2-Naphthylamine-4,8-disulfonic acid  
(1-Aminobutyl)phosphonic acid  
Boc-L-Tryptophan  
N-(tert-butoxycarbonyl)-L-Leucine  
Uvinul 490  
2,2,4,4 -Tetrahydroxy-Benzophone 97%  
2,4-Dihydroxybenzophenone  
2-Hydroxy-4-methoxybenzophenone  
2,2',4,4',6,6'-Hexanitrodiphenylamine  
Ammonium Picrate  
Fructosazine Deoxy  
1-Nitroso-2-naphthol  
5-Hydroxy-DL-lysine hydrochloride  
Potassium guaiacol sulphonate  
5-Hydroxyindole-3-acetic acid  
Glycerol triglycidyl ether  
Diisopropylammonium Iodide  
Acid Blue 138 | Carbolan Blue B  
2-Nitroso-1-naphthol  
Olanzapine  
Direct Orange 40  
Mercury(II) trifluoroacetate  
N-alpha-(2,4-Dinitro-5-fluorophenyl)-L-valinamide  
Ammonium tris(oxalato)ferrate(III)trihydrate  
Benzylcyclam  
Aluminium Chlorohydrate  
Pigment Green 7 | Heliogen Green L 8730  
Naphthoresorcinol  
(1,2-Cyclohexylenedinitrilo)tetraacetic acid  
2-Bromo-2-butene  
Xylenes Mixed  
2-(2-Hydroxy-5-Ethyl-phenyl) benzotriazole  
Tritolyl phosphate  
4-Aminosalicylic acid sodium salt  
Diethyl ethylmalonate  
3-Indolebutyric acid  
Stanclere 55

---

**IS NIR Spectra – Complete Collection**

6049 NIR spectra

Copyright © NICODOM 2012

NICODOM Ltd., Hlavní 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

[ftirspectra@gmail.com](mailto:ftirspectra@gmail.com)

Methyl 2-furoate  
1-Phenyl-1-propanol  
Cupric acetylacetonate  
Sodium ascorbate  
Pentanol  
Sodium decane-1-sulfonate  
8-Hydroxyquinoline sulfate  
Zinc tosylate  
Aminacrine hydrochloride  
Light Blue 2R  
Pentacarbonyl iron  
Trimethyl(2,3,4,5-tetramethyl-2,4-cyclopentadien-1-yl)silane  
Benzil  
2-Methoxybenzaldehyde  
2-Naphthol  
Sodium glutarate  
1,2-Bis(chlorodimethylsilyl)ethane  
Difluorodioxouranium  
Peptazin BAFD  
Advastab A 70  
Reinecke salt  
N-Phenyl-2-naphthylamine  
Triamminetrinitrocobalt(III); hydrate  
Sodium DL-3-Hydroxybutyrate  
1,3-Diphenyltriazene  
N-Methyl-gama-(1-naphthalenyloxy)-2-thiophenepropanamine  
3-Nitro-1-propyldiphosphonic acid teraethylester  
2-Cyclopentene-1-acetic acid  
Chromium Acetylacetonate  
Hermat ZDC  
2-Benzothiazolamine  
2-Aminobenzenethiol  
Calcium D-pantothenate, Vitamin B5  
2,4-Diaminophenol .2HCl  
Tolfenamic Acid  
N-Lauroylsarcosine sodium salt  
Bis(dimethylcarbamo-dithioato-S,S') copper  
Hermat ZDM  
2,5-Dibromoterephthalic acid  
2-Methyl-1-butanol  
Morpholinium tosylate  
N-t-BOC-L-valine  
Sodium propionate  
Cyclopentyl bromide

---

**IS NIR Spectra – Complete Collection**

6049 NIR spectra

Copyright © NICODOM 2012

NICODOM Ltd., Hlavni 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

[ftirspectra@gmail.com](mailto:ftirspectra@gmail.com)

Ascorbyl Palmitate  
Sodium L-Malate  
Tetra(Benzylmethylphoshinate) Cyclene  
L-(+)-Glutamic acid .HCl  
Poly[2-methoxy-5-(2-ethylhexyloxy)-1,4-phenylenevinylene]  
Trimethylphenylammonium chloride  
Cetrimonium Tosylate  
Azidamfenicol  
Salicin  
Chelidamic acid monohydrate  
Eszopiclone  
Aminomethanesulfonic acid  
Limonene  
N,N-Dimethyl-4-nitrosoaniline  
Ethynyl p-tolyl sulfone  
Natural Red 4 | Cochenille  
Aluminium acetate basic  
Nitrilotriacetic acid  
Fmoc-L-lysine hydrochloride  
Nickel Dibutyl dithiocarbamate  
Disodium EDTA  
4-Tolylboronic acid  
4-Phenoxyaniline  
UOP 88  
Aluminium acetylacetonate  
Bis(triphenylphosphine)palladium(II) chloride  
3,4-Dihydroxybenzaldehyde = Protocatechualdehyde  
trans-Cinnamic acid  
Benzyl acetate  
Gallotannin  
1,5-Diphenylcarbazine  
Iron(III) acetylacetonate  
Palladium bis(acetylacetonate)  
Zinc(II) acetylacetonate  
Benzyl cyanide  
N-Aminoethylpiperazine  
(1-Aminopropyl)phosphonic acid  
Bacitracin  
Bacitracin Zinc  
4-(t-Octyl)phenol  
N-Trityl Aminomethylphosphonic Acid Diethylester  
1-Hexadecylpyridinium bromide  
5,10,15,20-tetraphenyl-21h,23h-porphyrin zinc  
Cyanoacetohydrazide

---

**IS NIR Spectra – Complete Collection**

6049 NIR spectra

Copyright © NICODOM 2012

NICODOM Ltd., Hlavni 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

[ftirspectra@gmail.com](mailto:ftirspectra@gmail.com)

Ethylxanthic acid, potassium salt  
Basic Red 18 | Crypur Red GTL  
Cyclopropyltriphenylphosphonium bromide  
Diethyl adipate  
Butyl acrylate  
Monoethanolamine  
2,3,4,6-Tetra-O-acetyl-beta-D-glucopyranosyl isothiocyanate  
Sodium formate  
Ethyl acetate  
Mesityl oxide  
Malonic Acid  
Rhodanine  
2,6-Diaminopyridine  
Dibenzo-18-crown-6  
2-Thiouracil  
Ethyl acetoacetate  
Aniline hydrochloride  
Methyl-3-aminocrotonate  
Tetrakis(triphenylphosphine)palladium  
4'-Hydroxyphenyl-Carvedilol  
5'-Hydroxyphenyl-Carvedilol  
Cyclopentene  
Sodium-L-Glutamate  
Hexanoic acid  
Piperazine dihydrochloride  
Nickel Diethyldithiocarbamate  
Iminodiacetic acid  
Heptane  
Pigment Red 88 | Thioindigoid Red  
Butyl ether  
Pamam dendrimer, generation 1  
Pigment Green 36 | Heliogen Green K 9360  
Dodecanoic acid  
1-Nonanol  
1-Bromododecane  
Sodium oleate  
Zinc phthalocyanine  
Triethylene Glycol Monobutylether  
Hermet ZDK  
Tetraethylene glycol dimethyl ether  
Oleyl alcohol, Cis-9-octadecen-1-ol  
2,7-Dibromo-9H-fluoren-9-one  
Tetraphenylboron Sodium  
2,6-bis(Dibromomethyl)pyridine

---

**IS NIR Spectra – Complete Collection**

6049 NIR spectra

Copyright © NICODOM 2012

NICODOM Ltd., Hlavni 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

[ftirspectra@gmail.com](mailto:ftirspectra@gmail.com)

trans-Cinnamaldehyde  
Phenol Red; Phenolsulfonephthalein  
Zoldine MS Plus  
Tris(2-cyanoethyl)phosphine oxide  
Dextrorphan-D-Tartrate  
Sodium 5,5-diethylbarbiturate  
2,2,4-Trimethyl-1,3-pentanediol  
Magnesium Citrate  
Diethanolamine hydrochloride  
D-Glucose Monohydrate  
Disodium Citrate  
2-Iodoacetamide  
1-Methyl-4-piperidone  
Telmisartan  
5-Pregnen-3 $\beta$ -Ol-20-On  
Palladate(2-), tetrakis(thiocyanato-S)-, potassium (1:2)  
Manganese pentacarbonyl bromide  
2,5-Bis(tributylstannyl)thiophene  
1-Heptadecanol  
3,5-Diamino-1,2,4-triazole  
1,3,5-Tri-t-butylbenzene  
(4-Nitrophenyl)acetaldehyde  
trans-1,2-Cyclohexanediol  
Calcein  
Chlorotributyltin  
Bromotributyltin  
Riboflavin-5'-Phosphate  
Hermet FEK  
1-amino-2-( $\alpha$ -naphthylamino)ethane dihydrochloride  
Ethyl 5-bromovalerate  
2-(p-Iodophenyl)-3-(p-nitrophenyl)-5-phenyl-2H-tetrazolium chloride  
3-Acetylthiophene  
Chlorotris(triphenylphosphine)rhodium(I)  
Rosuvastatin calcium  
Pigment Blue 15:6 | Phthalo Blue Reddish  
Arazate  
Tinuvin 1577 FF  
Agerite Resin D  
2-Iodohippuric acid  
D-Tartaric Acid  
m-Xylylenediamine  
2,4,6-Tribromoaniline  
L-Proline  
Thiosalicylic acid

---

**IS NIR Spectra – Complete Collection**

6049 NIR spectra

Copyright © NICODOM 2012

NICODOM Ltd., Hlavni 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

[ftirspectra@gmail.com](mailto:ftirspectra@gmail.com)

Diethyldithiocarbamate sodium  
8-Hydroxyquinoline  
Chromotropic Acid  
2-Ethylpiperidine  
Tetraoctylammonium bromide  
Bismuth Subsalicylate  
5-Cyano-1-pentyne  
Antipyrilazo III  
2-Benzothiazolethiol  
Trifluoromethanesulfonic acid  
Erythritol  
1-Fluoro-2-nitrobenzene  
Hexadecyltributylphosphonium bromide  
Leptacit  
Tiron=4,5-Dihydroxy-1,3-benzenedisulfonic acid, disodium  
9,10-Diphenylanthracene  
Gallic acid  
Basic Blue | Basacryl Blue X3 GL  
4-Aminobenzoic acid; PABA  
Cyclododecene  
DL-Phenylalanine  
Hydroquinone monomethylether  
3,7,11,15-Tetramethyl-2-hexadecen-1-ol  
Sodium succinate  
Sodium lauryl sulfate  
2-Hydroxy-4-butoxybenzophenone  
Ethyleneimine  
Diethyldithiocarbamic acid diethylammonium salt  
Hydroxyamphetamine hydrobromide  
2-(Methylsulfonyl)ethanol  
Verapamil .HCl  
2-acrylamido-2-methyl-1-propanesulfonic acid  
2-Naphthyl acetate  
UOP 688  
Gadolinium acetate trihydrate  
1,2,3,4-Tetrabromobutane  
Ethyltriphenylphosphonium bromide  
Diclofenac Sodium  
Diclofenac potassium  
DL-Menthol  
2-Aminofluorene  
D-Tryptophan  
Catechin  
Potassium diphenylphosphide

---

**IS NIR Spectra – Complete Collection**

6049 NIR spectra

Copyright © NICODOM 2012

NICODOM Ltd., Hlavni 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

[ftirspectra@gmail.com](mailto:ftirspectra@gmail.com)

Hermet Zn MBT  
Cholesteryl methyl carbonate  
Rosiglitazone Maleate  
1,5-cyclooctadiene  
Nickel bis(dimethyldithiocarbamate)  
Potassium Hexathiocyanatochromate(III) tetrahydrate  
Dibutyltin-bis-2-ethylhexylmaleate  
Chlorotriphenylsilane  
Di-n-Octyltin bis(2-ethylhexyl) mercaptoacetate  
Pamam dendrimer, generation 0  
Pyridinium perchlorate  
1,4-Dibromo-2,5-bis(octyloxy)benzene  
Chloropentamethyl disilane  
Pamam dendrimer, generation 5  
Galloyanine  
[1,3-Bis(diphenylphosphino)propane]nickel(II) chloride  
p-Phenetidine  
9-Bromoanthracene  
Sodium butyrate  
2-Aminoethanethiol hydrochloride  
2,6-Pyridinedicarboxylic acid diethyl ester  
cis-1,2-Dichloroethylene  
3-Amino-1-propanol  
Cyclohexylmercaptan  
1,7-Dibenzylcyclene  
Ethyl isonicotinate  
D-Arginine  
Edetate Dipotassium Magnesium  
Chlorotris(triphenylphosphine)copper  
Methyl 4-formylbenzoate  
Phenylphosphonic Acid  
1-Ethynyl-naphthalene  
N-(tert-Butoxycarbonyl)-L-Alanine  
N-t-BOC-L-proline  
p-Toluenesulfonhydrazide  
Pigment Red 48:3  
Sunset yellow FCF  
Adamantylamine methylphosphinate  
Diethylenetriaminopentamethylenephosphonic acid  
α-Chloralose  
Aminomethyl(phenyl)phosphinic acid  
Stearate Calcium  
2-[2-hydroxy-5-(1-phenylethyl)-phenyl]benzotriazol  
2-Amino-3-Picoline

---

**IS NIR Spectra – Complete Collection**

6049 NIR spectra

Copyright © NICODOM 2012

NICODOM Ltd., Hlavni 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

[ftirspectra@gmail.com](mailto:ftirspectra@gmail.com)

Ethyl benzoylformate  
Zinc lactate  
Potassium Antimony Tartrate  
Creatinine zinc chloride  
Isosorbide 5-Mononitrate  
Dicyclohexyl-18-crown-6  
Bromosulfonazo III  
1,6-Hexanediol diglycidyl ether  
Xylenol Orange  
Dimethyl Aminomalonate Hydrochloride  
(4-bromophenylethynyl)-trimethylsilane  
Piperidylmethylphosphinic Acid  
1,7-Bis(tert-butoxycarbonylmethyl)-1,4,7,10 tetraazacyclododecane  
Ferrocene, [(4S)-4,5-dihydro-4-(1-methylethyl)-2-oxazolyl]-  
Flunitrazepam  
p-Phenylenediamine Sulfate  
Zinc salicylate  
Solvent Orange 63 | Hostasol Red GG  
Chlorodiphenylsilane  
Pamam dendrimer, generation 4  
2-(4-Hydroxyphenylazo)benzoic acid  
2-Ethoxy-1-ethoxycarbonyl-1,2-dihydroquinoline  
DL-Isocitric acid trisodium  
3-Thiophenecarbonitrile  
o-Tolylboronic acid  
Tetrabutylammonium bromide  
2,7-Dibromofluorene  
1-Ethynyl-2-nitrobenzene  
1,4,7,10-tetraazacyclododecane-1,4,7,10-Tetramethylphosphinic Acid  
Dimyristyl-thio-dipropionate  
Pigment Yellow 4  
Phenylphosphonic Acid Methyl Ester Sodium Salt  
Acid Red 88  
3,5-Bis(trifluoromethyl)phenyl isocyanate  
4-(2-Pyridylazo)resorcinol sodium salt  
Tetraethyl methylenediphosphonate  
Bathophenanthroline  
1,2-Bis(diphenylphosphino)ethane  
2',4',6'-trimethylacetophenone  
Magnesium acetate tetrahydrate  
Glycinamide hydrochloride  
(2,2'-dihydroxy)diethylamine tritosylate  
Benzamidine hydrochloride  
Diglycidyl ether of bisphenol-A

---

**IS NIR Spectra – Complete Collection**

6049 NIR spectra

Copyright © NICODOM 2012

NICODOM Ltd., Hlavni 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

[ftirspectra@gmail.com](mailto:ftirspectra@gmail.com)

mono-Methyl terephthalate  
Silver tosylate  
N,N,N',N'-tetrakis(2-pyridylmethyl)ethylenediamine  
2-Amino-3-hydroxypyridine  
1,2-cyclohexanedicarboxylic acid  
4-Phenylazophenol  
Dastib 255  
Acid Violet 49 | Acid Violet 6B  
3,3'-Dicarboxy-1,1'-dihydroxy-2,2'-bipyridinium  
Ferrocene, [[bis(1-methylethyl)amino]carbonyl]-  
1,4,8-Tris-(t-butoxycarbonyl) cyclam  
(+)-Dibenzoyl-D-tartaric acid  
Diphenylphosphinic acid  
Irganox 1330  
1-Ethynyl-3,5-dimethoxybenzene  
1-Bromopyrene  
Phenethyldimethylchlorosilane  
Benzyl trans-3a,5a,8a,10a-tetraazaperhydropyrene  
Dibenzyl phosphite  
Tetrazolium violet  
L-Mandelic Acid  
Imino-bis(methylphosphonic acid)  
2-Chloro-3-fluoropyridine  
Cresol Red  
Ferrocene, 1-carboxy-1'-(diphenylphosphino)-  
3-(Trifluoromethyl) Phenylthiourea  
Sulfonazo III  
3-(Dimethylamino)propionitrile  
tert.butyl 2,2'-(1,4,7-triazacyclononane-1,4-diyl)diacetate  
Raloxifene-6'-D-glucuronide  
4-Amino-1,8-naphthalimide  
18-Crown-6  
Powderlink 1174  
Hematoporphyrin HCl  
4-Amino-3-hydrazino-1,2,4-triazol-5-thiol  
Zirconium(IV) acetylacetonate  
Polypox R14  
Metolazone  
1,2-Diaminoanthraquinone  
Sodium iminodiacetate  
Cyclopropanecarboxylic acid  
Silquest A-1120  
Salicylidene-o-Aminophenol  
5-Bromosalicylaldehyde

---

**IS NIR Spectra – Complete Collection**

6049 NIR spectra

Copyright © NICODOM 2012

NICODOM Ltd., Hlavni 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

[ftirspectra@gmail.com](mailto:ftirspectra@gmail.com)

D(+)-Raffinose pentahydrate  
Sodium 3-mercaptopropanesulphonate  
(+)-Dibenzoyl-l-tartaric anhydride  
Hematoporphyrin dihydrochloride  
Pigment Red 264  
Trimethylsulfoxonium iodide  
Dibenzyl-cross-bridged-Cyclam  
Phenylphosphinic acid  
Butyltriphenylphosphonium bromide  
D-Aspartic acid  
Eriochrome Black T  
Dimethyl-2,5-dibromoterephthalate  
Gatifloxacin Sesquihydrate  
1,4-bis(2-(5-phenyl)-oxazolyl) benzene  
Diisopropyl phosphite  
1-Amino-2-methyl-propylphosphonic acid  
Butamirate citrate  
1-(Trimethylsilyl)imidazole  
tert-Butyldimethylsilyl chloride  
1-Iodo-2-(trimethylsilyl)acetylene  
t-butyl Dimethylsilanol  
5-Chlorouracil  
Direct Yellow 9  
Dimethyl 1,3-acetonedicarboxylate  
6,9-Diamino-2-ethoxyacridine lactate .H2O  
Diphenyl Vinyl Chlorosilane  
1,1,3-tris(2-methyl-4-hydroxy-5-tert-butylphenyl)butane  
Dioctadecyl sulfide  
Chlorhexidine Digluconate, 20% solution  
Sulfofon O 680  
S-binol, (S)-(-)-1,1'-Bi-2-naphthol  
Ammonium Benzoate  
Tetrazolium blue chloride  
4-Bromophenylacetic acid  
Topanol A  
1,10-Phenanthroline chloride monohydrate  
Phenylchloroformate  
Anthranilonitrile  
Magnesium lactate  
1-(3-Dimethylaminopropyl)-3-ethylcarbodiimide  
Ferrocene, (tricyclo[3.3.1.1.3,7]dec-1-ylmethyl)- (9CI)  
Mono-L-Valinium Nitrate  
trans-4-Phenyl-3-buten-2-one  
Palladium, di-m-chlorobis[2-[(dimethylamino-kN)methyl]phenyl-kC]di-

---

**IS NIR Spectra – Complete Collection**

6049 NIR spectra

Copyright © NICODOM 2012

NICODOM Ltd., Hlavni 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

[ftirspectra@gmail.com](mailto:ftirspectra@gmail.com)

Monosodium Citrate Anhydrous  
Morpholinoethanesulfonic acid  
Thymolphthalexon  
Chlorophosponazo  
Ammonium-D-gluconate  
2,6-Dichloriodobenzene  
Tetrabutylammonium perchlorate  
6-Piperidinopurine  
FD&C Yellow 5, Tartrazine  
Direct Red 79  
Naphthol Green B  
Tetrabutylammonium Nitrate  
2-Ethynylpyridine  
2,2'-Bipyridine-4,4'-diylbis(phosphonic acid)  
2:2',6':2''-Terpyridine-4-yl-4-phenylphosphonic Acid  
tert. Butylhydroquinone  
Cholesteryl Elaidate  
2-Dodecen-1-yl succinic anhydride  
Trazodone HCl  
1,3-Diphenylacetone p-tosylhydrazone  
Sodium Octanoate  
Methylenediphosphonic acid  
Cholesteryl Myristate  
9,9-Dioctyl-2,7-dibromofluorene  
Cholesteryl Cinnamate  
2-(2-hydroxy-5-n-butyl-phenyl) benzotriazol  
Tetraphenylphosphonium chloride hydrate  
2-(2-hydroxy-5-n-octyl-phenyl) benzotriazol  
3-Methoxy-L-tyrosine monohydrate  
Calcium Dobesilate  
1-Tetradecylamine  
Phosclere T 215  
Bromoacetaldehyde diacetyl acetal  
2,2-Dimethyl-1,3-dioxane-4,6-dione  
1-Naphthol-3,6-disulfonic acid  
2-Bromobenzonitrile  
Cyclohexanecarboxaldehyde  
2,4-Dichloro-1-naphthol  
Hydropalat 875  
Sodium diethyldithiocarbamate trihydrate  
Fluoranthene  
Cobalt (II) tosylate  
Bis(diphenylphosphino)methane  
8-Hydroxyquinoline-5-Sulfonic Acid Hydrate

---

**IS NIR Spectra – Complete Collection**

6049 NIR spectra

Copyright © NICODOM 2012

NICODOM Ltd., Hlavni 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

[ftirspectra@gmail.com](mailto:ftirspectra@gmail.com)

Trimethylsilyl Benzoate  
2,6-Diisopropyl Phenol  
Thiolane 1,1-dioxide  
Digoxin  
Lead(II) cyanamide hydrate  
5,5-Diphenyl-2-thiohydantoin  
Diphenyllead(IV) dichloride  
10-Camphorsulfonyl chloride  
Dibenzyl dimethylcyclam  
4-Aminophenylboronic acid pinacol ester  
Alizarin Safirol SE  
Busan BL1166  
Diethyl propylmalonate  
Chromium(III) acetylacetonate  
(4-Nitrophenyl)phosphonic acid  
Diphenylfulvene  
Triphenylene  
Dicetyl phosphate  
3,5-dimethylpyrrole-2-carboxaldehyde  
1,4,7,10-Tetraazacyclodecane-1,4,7,-Triacetic Acid Monolactame  
Polypox R 6 (o-cresyl glycidyl ether)  
1-Phenyl-1-Butanol  
Phenylthiohydantoic acid  
L-Menthol  
L-Canavanine sulfate  
Tramadol hydrochloride  
Tetrabutylammonium Fluoride Hydrate  
Cadmium Stearate  
Nickel Stearate  
4-Hydroxy-2,2,6,6-tetramethylpiperidin-1-oxyl  
Thiobenzamide  
Propyl-cis-decahydro-2a,4a,6a,8a-tetraazacyclopent[fg]acenaphthylene  
2-(2-hydroxy-5-Methyl-phenyl)-5- chlorobenzotriazol  
Mordant Orange 1  
4-(2-Thiazolylazo)resorcinol, TAR  
Aminoacetaldehyde dimethyl acetal  
Metipranolol  
Ferrocene, (fluorocarbonyl)-  
Sodium bis(2-methoxyethoxy)aluminumhydride  
Butyltin hydroxide oxide  
Sodium 1-heptanesulfonate  
N-Acetyl-D-Tryptophan  
2,4-Thiazolidinedione  
Phenanthridine

---

**IS NIR Spectra – Complete Collection**

6049 NIR spectra

Copyright © NICODOM 2012

NICODOM Ltd., Hlavni 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

[ftirspectra@gmail.com](mailto:ftirspectra@gmail.com)

Hydrobromothiazide  
3,5-bis(trifluoromethyl)-phenyl isothiocyanate  
Fluorescein  
Tris(N-piperidinomethyl)phosphine Oxide  
Iron Citrate  
4-(1-Hydroxyethyl)pyridine  
N-(2-aminoethyl)-N-benzylethane-1,2-diamine  
1,1-bis(4-hydroxy-3-methylphenyl)cyclohexane  
2-(2-Hydroxyphenyl)-5-chlorobenzotriazol  
p-Sulfophenylazochromotropic acid  
N,N,5-trimethyl-1,3,2-benzodioxaphosphol-2-amine  
Ambroxol .HCl  
Cholesteryl Ethyl Carbonate  
Fascat 4100  
1,3-Dicyclohexylurea  
4,7-Dibromo-2,1,3-benzothiadiazole  
Tetraoctylphosphonium Bromide  
Dastib 244  
Acid Yellow 7 | Brilliant Sulpho Flavin FF  
1,6-Heptadiyne  
1,8-Nonadiyne  
2,2,6,6-Tetramethyl-4-piperidinol  
1,2,2,6,6-Pentamethyl-4-Piperidinol  
Pyridinium toluene-p-sulfonate  
2-tert-Butyl-4-methylphenol  
Cresolphthalexon  
1,4-Butanediol diglycidyl ether  
Pigment Red 3  
Mycophenolic acid  
2-Bromo-3-thiophene carboxylic acid  
L-Fucose  
N-bis(Aminomethylphosphonic acid) Glycine  
2-(2-Hydroxy-5-methyl phenyl) Benzotriazole 97%  
Di-t-butyl dicarbonate  
Nonanoic acid vanillylamide  
Dianate Na  
9-(Chloromethyl)anthracene  
Terpin hydrate  
Cholesteryl Erucate  
2,6-Diisopropylaniline  
4-Pyrrolidinopyridine  
2-Ethylhexyl glycidyl ether  
Auramine O  
1,2,3-Tris-(2-cyanoethoxy)propane

---

**IS NIR Spectra – Complete Collection**

6049 NIR spectra

Copyright © NICODOM 2012

NICODOM Ltd., Hlavni 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

[ftirspectra@gmail.com](mailto:ftirspectra@gmail.com)

Didodecyl sulfide  
5,5',7,7'-Tetrabromoindigo  
2',5'-Dichloroacetophenone  
Solvent Yellow 56 | Oilsol Yellow DEA  
1-Azidoadamantane  
4-Hydroxydeoxybenzoin  
di-n-Octadecyl disulfide  
Butylated hydroxyanisole, BHA  
Pigment Yellow 1  
Butyltin Trilaurate  
Nonylphenol  
Mark 954  
2-Azetidinecarboxylic acid  
3-Methoxy-1-butanol  
Poly(n-isopropyl acrylamide)  
Poly(tetrahydrofuran) based polymer  
Texanol®  
3-Glycidoxypropyltrimethoxysilane  
3-Methacryloxypropyltrimethoxysilane  
Acid Red 103  
2,2'-Dihydroxyazonaphthalene-4-sulfonic acid, sodium salt  
Phosclere T 310  
Rhodium (triphenylphosphine)carbonylacetylacetonate  
TPM, Tripropylene Glycol (Mono) Methyl Ether  
3-Chloropropyltrichlorosilane  
4-(1,1,3,3-tetramethylbutyl)phenyl salicylate  
Hexabromocyclododecane  
Ethylene glycol succinate  
3-Methyl-1,3-butanediol  
Poly(sodium 4-styrenesulfonate)  
L-Alanine-15N  
2-Bromoethylamine hydrobromide  
Aquazol® 200  
3-Methyl-4-nitrophenol  
Aminoguanidine Bicarbonate  
Acid Blue 127:1 | Sandolan Walk Blue N-GLN 180  
1-Bromo-2-Methylnaphtalene  
Pinacryptol Yellow  
1-Hydroxybenzotriazole  
1-(3-Dimethylaminopropyl)-3-ethylcarbodiimide hydrochloride  
D-Talose  
Tinuvin 328  
Polyethylenimine 50% in Water  
2,6-dimethyl-3,5-pyridinedicarboxylic acid

---

**IS NIR Spectra – Complete Collection**

6049 NIR spectra

Copyright © NICODOM 2012

NICODOM Ltd., Hlavni 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

[ftirspectra@gmail.com](mailto:ftirspectra@gmail.com)

Igepal CO 520  
Acridine  
Direct Blue 1 | Chicago Blue 6B  
Direct Red 81 | Durazol Red 2B  
Ponceau 4R, Cochineal Red A, Brilliant Scarlet 4R  
1-(triphenylmethyl)pyridinium tetrafluoroborate  
Tricyclohexyl phosphine  
2,7,9-Tribromo-9H-fluorene  
Mark OTM/68  
Methylcyclopentadiene dimer  
Ethyl-cis-decahydro-2a,4a,6a,8a-tetraazacyclopent[fg]acenaphthylene  
Benzylmethyl--cis-decahydro-2a,4a,6a,8a-tetraazacyclopent[fg]acenaphthylene  
Xylene Cyanol FF  
Phosphit TNPP  
Dodecenylsuccinic Anhydride  
4,4- thiobisphenol  
Sorbitan tristearate  
trans-3a,5a,8a,10a-tetraazaperhydropyrene  
Glycidyl Neodecanoate  
Methylene Green  
Indapamide  
tert.butyl 1-methyl-2-propynyl ether  
Cinnamyl chloride  
1-Ethyl-2-pyrrolidinone  
N-Octyl pyrrolidone  
1-Lauryl-2-pyrrolidone  
2,2'-[(Diocetylstannylene)bis(thio)]diacetic acid dibutyl ester  
Direct Black 4 | Diazol Black ERN Quadruple  
Bis-(sodium sulfopropyl)-disulfide  
Benzamide  
2,1,3-Benzothiadiazole  
1,10-Phenanthroline-5,6-dione  
Diethyl diphenylmethylphosphonate  
Dibenzoyl L-tartaric acid  
Tetraphenylphosphonium bromide  
Glycine tert butyl ester hydrochloride  
DL-Dithiothreitol  
(2-Bromoethyl)trimethylammonium bromide  
3-Methyl-2-cyclopenten-1-one  
Acid Red 44 | Crystal Scarlet 6R  
Silquest A-171  
D-Galactonic acid g-lactone  
1,12-Diaminododecane  
Lithium D-(-)-lactate

---

**IS NIR Spectra – Complete Collection**

6049 NIR spectra

Copyright © NICODOM 2012

NICODOM Ltd., Hlavní 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

[ftirspectra@gmail.com](mailto:ftirspectra@gmail.com)

Pigment Red 170  
5-Nitrobenzimidazole nitrate  
2,6-bis(2-benzimidazolyl)pyridine  
1,4-Diazobicyclo(2.2.2)octane  
1-Hydroxyethylidene-1,1-diphosphonic acid  
Adamantane  
Pigment Red 4  
Ioxitalamic Acid  
Xylyl bromide  
Sodium 1-hexanesulfonate  
11-Bromoundecanoic acid  
2-Phenyl Glycine  
2-Aminobenzophenone  
DL-2-Aminobutyric acid  
Poly(2-Vinylnaphthalene)  
1h-benzotriazole-1-methanol  
Cotton Blue  
D-Arabinose  
Cyclopentane  
p-Cyanobiphenyl  
Pyrazole  
Imidazole  
1,2,4-Triazole  
Sodium Dichlorisocyanurate  
Watermelon Ketone  
Silquest Y-15078  
Alprazolam  
Z-glycine-N-succinimidylester  
Pyrazine  
4-Ethynylbiphenyl  
Atenolol  
1-Hexadecanethiol  
4-Phenylbenzotrionitril  
Potassium trifluoroacetate  
Sodium 2-methyl-2-propanethiolate  
Dowanol PnB  
Methylxlenol blue  
4-chlorophenylsalicylate  
2,4-dichlorophenylsalicylate  
N-benzyl-2-bromoacetamide  
Direct Black 17 | Azogen Black D  
Cyclen  
12-Crown-4  
Basic Yellow 40 | Maxilon Flavine 10 GFF

---

**IS NIR Spectra – Complete Collection**

6049 NIR spectra

Copyright © NICODOM 2012

NICODOM Ltd., Hlavni 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

[ftirspectra@gmail.com](mailto:ftirspectra@gmail.com)

Sumilizer TP-D  
Aminomethyldiphosphonic Acid  
N-Dimethyl Aminomethyldiphosphonic Acid  
N-Acetyl-Aminomethylbisphosphonic Acid  
Clopidol, 3,5-Dichloro-2,6-dimethyl-4-pyridinol  
Glyoxylic acid  
4-Dodecyloxy-2-hydroxybenzophenone  
Cinnarizine  
Nitrotetrazolium blue chloride  
Neotetrazolium chloride  
2,3,5-Triphenyl-2H-tetrazolium bromide  
2,4-Dichloriodobenzene  
tris(Methylenephosphonic Acid) Cyclam  
Sparteine sulfate  
Naftidrofuryl oxalate  
Triphenylgermanium bromide  
Hydropalat 680  
3,5-Dimethylisoxazole  
4-Hydroxy-3-nitrobenzaldehyde  
Ammonium Citrate Dibasic  
N-Glycyl-aminomethylphosphonic acid  
Chloral hydrate  
DL-Alanine  
Diphenylguanidine  
DL-Serine  
Deoxycholic acid; sodium salt monohydrate  
Clomipramine  
meso-2,3-Dimercaptosuccinic acid  
Sodium Potassium Tartrate  
N-Benzoyl-N-phenylhydroxylamine  
Polypox R20  
Murexide  
Brilliant Yellow  
Paraform  
Basic Orange 21  
L-Carnosine  
Homovanillic acid  
Titanium(IV) 2-ethylhexyloxyde  
Bis-[2-hydroxy-5-methyl-3-(benzotriazol-2-yl)-phenyl]-methane  
Thallium(I) oxalate  
1,3,5,7,9,11-Octaisobutyltetracyclo[7.3.3.15,11]octasiloxane-endo-3,7-diol  
1,3,5,7,9,11,14-Heptaisobutyltricyclo[7.3.3.15,11]heptasiloxane-endo-3,7,14-triol  
Diethyl (hydroxymethyl)phosphonate  
Titanium Ethoxide

---

**IS NIR Spectra – Complete Collection**

6049 NIR spectra

Copyright © NICODOM 2012

NICODOM Ltd., Hlavní 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

[ftirspectra@gmail.com](mailto:ftirspectra@gmail.com)

Isoamyl alcohol  
Aldrin  
Benzylmethylcyclam  
1,7-Dimethyl-4,11-bis(Methylphosphonic Acid) Cyclam  
Irgastab 2002  
Cholesteryl Undecylate  
Ammonium hydrogentartrate  
Tetrabutylammonium hexafluorophosphate  
Ibuprofen Sodium  
Tetrabutylammonium iodide  
Tetrabutylphosphonium Bromide  
Zinc 5,10,15,20-tetra(4-pyridyl)-21 h,23 h-porphine  
Solvent Orange 7 | Oilsol Orange KB  
Alpha-naphthyl propionate  
Bis-(4-fluoro-3-nitrophenyl)sulfone  
D-Serine  
4-phthalimidobutyric acid  
DL-Asparagine monohydrate  
Octafluoronaphthalene  
3,4-Dibromothiophene  
2-(2-hydroxy-5-n-dodecylphenyl)-benzotriazol  
3-Methyl-3-oxetanemethanol  
D-Camphorsulfonic acid  
Glyceryl monostearate  
Tris(2,4-di-tert-butylphenyl) phosphite  
Ammonium Tartrate  
Aminophylline  
Cyclohexyl isocyanate  
1,18-Dibromooctadecane  
Pigment Red 208 based  
9,9-Dioctylfluorene-2,7-diboronic acid bis(1,3-propanediol) ester  
3-Dimethylamino-1-Propanol  
L-(+)-Ornithine monohydrochloride  
Propranolol Hydrochloride  
Beta-Alanine methyl ester hydrochloride  
D-Isoleucine  
Bis(dibenzylideneacetone)palladium  
5-Chlorosalicylic acid  
Adenine Sulfate  
Direct Yellow 50  
DL-Monosodium glutamate  
Fuchsin acid  
Tetrabutylammonium hydrogen sulfate  
Nickel (II) Acetylacetonate

---

**IS NIR Spectra – Complete Collection**

6049 NIR spectra

Copyright © NICODOM 2012

NICODOM Ltd., Hlavni 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

[ftirspectra@gmail.com](mailto:ftirspectra@gmail.com)

Ferrocene, 1,1'-bis(diphenylphosphinyl)-  
Irganox MD 1024  
Chloromethyl thiocyanate  
2-Thenoyl-trifluoroacetone  
Ethanolamine  
D-Norleucine  
L-Norleucine  
Phenyl phosphate, disodium salt .2H<sub>2</sub>O  
Pivaloyl chloride  
Dimethylphosphinic acid  
D-Leucine  
2-Ketoglutaric acid  
9-Ethynylphenanthrene  
Trimethylolpropane Trimethacrylate  
Tiaprofenic Acid  
Dansyl Hydrazine  
Ferrocene, 1-[(dimethylamino)methyl]-2-(diphenylphosphino)-  
3,4-Dihydroxycinnamic acid= Caffeic acid  
N-(p-methylphenyl)-1,4-phenylenediamine  
2,5-Dimethoxy-2,5-dihydrofuran (cis+ trans)  
Diltiazem Hydrochloride  
Aminomethyldiphosphonic Acid Triammonium Salt  
Decanoic acid  
Nickel Formate  
Manganese stearate  
Direct Blue 67  
(Methylthio)dimethyl sulfoxide  
Palladium diacetate  
N,N,4-trimethyl-1,3,2-benzodioxaphosphol-2-amine  
Triclosan  
D-Alanine  
5-(2-Pyridyl)tetrazole  
Boc-Glycine N-Hydroxysuccinimide  
2,2,6,6-Tetramethyl-4-piperidone hydrochloride  
Trans-4-Hydroxy-D-proline  
Aminomethylphosphonic Acid Diphenylester  
Acebutolol Hydrochloride  
2-hydroxy-4-octyloxy-5-tert. Butyl-benzophenone  
N-(2-Hydroxyethyl)-2-pyrrolidone  
Dansyl Fluoride  
Diethyl(2-oxo-2-phenylethyl)phosphonate  
Isoxicam  
D-Mannose  
Ammonium Citrate Tribasic

---

**IS NIR Spectra – Complete Collection**

6049 NIR spectra

Copyright © NICODOM 2012

NICODOM Ltd., Hlavni 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

[ftirspectra@gmail.com](mailto:ftirspectra@gmail.com)

Dowanol DPM, Di(Propylene Glycol) Methyl Ether  
Diethylenetriamine-N'-methylenephosphonic-N,N,N'',N''-tetraacetic acid  
(+)-alpha-Dihydrotetrabenazine  
Pigment Orange 5 | Hansa Red GG  
2,7-Dimethoxynaphthalene  
Ferrocene, (tricyclo[3.3.1.1<sup>3,7</sup>]dec-1-ylcarbonyl)- (9CI)  
D-Gluconic acid, hemicalcium salt  
Bromothymol Blue  
o-Fluoroiodobenzene  
2-Fluoroaniline  
1-iodoanthraquinone  
3-Acetylpyridine  
1-Ethynyl-2,4,5-trimethylbenzene  
Pigment Orange 13 | Vynamon Orange G FW  
Mono-L-Alaninium Nitrate  
Guanidineacetic acid  
Di-L-Leucine Selenate  
Arsenazo I  
Uvistab 247  
Cholesteryl stearate  
Eriochrome Cyanine R  
Acid Red 33, CI 17200  
Carmoisine, Azorubine  
Fluorescein mercuric acetate  
Tri-L-Valine Selenate Monohydrate  
Brucine  
Trifluoromethanesulfonic anhydride  
9,10-Dihydro-9-oxa-10-phosphaphenanthrene 10-oxide  
L(-)-Camphorsulfonic acid  
1-Benzyl-4-piperidone  
D-(+)-Fucose  
L-Rhamnose  
D-Mannoheptulose  
N-butyl-2,2,6,6-tetramethyl-4-piperidinamine  
Cyclododecene  
Nile blue A  
Fenbufen  
Tris[4,4'-thiobis[6-tert-butyl-m-tolyl]] phosphite  
Tosylaziridine  
3-Fluorocatechol  
4'-Methylbiphenyl-4-carbaldehyde  
Irganox 245  
Basic Red 13 | Astrazon Pink FG  
Cetyl alcohol

---

**IS NIR Spectra – Complete Collection**

6049 NIR spectra

Copyright © NICODOM 2012

NICODOM Ltd., Hlavni 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

[ftirspectra@gmail.com](mailto:ftirspectra@gmail.com)

2-Fluorophenol  
D-Glucose 6-phosphate, disodium salt  
Ethyl 2,3-dibromopropionate  
2,2,6,6-tetramethyl-4-aminopiperidin  
Isopropyl-beta-D-thiogalactopyranoside  
2,4,6-Tri(2-pyridinyl)-1,3,5-triazine  
Triethyloxonium tetrafluoroborate  
3-Amino-1-propanesulfonic acid  
Pigment Yellow 139  
di-toluene-4-sulfimide  
5,10,15,20-tetrakis-(n-methyl-4-pyridyl)-21,23h-porphyrin tetratosylate ethylenediamine  
5-Methoxy-2-mercaptobenzimidazole  
Biscyclohexanone Oxalyldihydrazone  
Beta-cyclodextrin, sulfated sodium salt  
2-(Aminomethyl)pyridine  
3-(Aminomethyl)pyridine  
4-Pyridinemethanamine  
1,8-diaminooctane  
Red 2G, Acid Red 1  
Hydroxymethylcellulose  
Heptafluorobutyric Acid  
Polycaprolactone triol  
Ethyl 2-indolecarboxylate  
Carboxyarsenazo III  
2-(tert-butylamino)ethyl methacrylate  
Ergotamine Tartrate  
3-Bromo-1-(Trimethylsilyl)-1-Propyne  
Glycinethymol Blue  
2,6-Di-tert-butyl-4-methylpyridine  
N-(N'-Phthaloylglycyl)-Aminomethylphosphonic Acid  
FD&C Blue 1, Brilliant Blue FCF  
1-Trimethylsilyl-1-hexyne  
Tinuvin 320  
Sandostab P-EPQ  
Acid Blue 92 | Coomassie Blue RL  
Tinuvin 327 UV - 2  
1,1,1-triphenyl-N-(phenylmethylidene)methylamine  
Bis(2-hydroxypropyl)-p-toluidine  
D-Glucose, 2,3,4,5,6-pentaacetate  
Aerosol 22 N  
Tinuvin 326 UV - 1  
Pigment Red 166 | CPT - Scarlet Red  
1,1-diphenyl-2-propyn-1-ol  
Imidazolidinyl urea

---

**IS NIR Spectra – Complete Collection**

6049 NIR spectra

Copyright © NICODOM 2012

NICODOM Ltd., Hlavni 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

[ftirspectra@gmail.com](mailto:ftirspectra@gmail.com)

Sodium chloroacetate  
1,9-Nonanediol  
Jeffamine T 403  
Polypropylene glycol diglycidyl ether  
N,N'-bis(salicylidene)-1,2-phenylenediamine  
2'-Deoxycytidine .HCl  
1-Hydroxy-7-azabenzotriazole  
P-Dithiane-2,5-diol  
Bromoxlenol Blue  
N-(2-iodoethyl)-phthalimide  
Isopropyl Glycidyl Ether  
2-((Trimethylsilyl)-ethynyl)thiophene  
Palladium, bis[m-(acetato-O:O')]bis[2-[(dimethylamino-N)methyl]phenyl-C]- (9CI)  
Pigment Red 177 | Cromophtal Red A2B  
4-Chloro-1,8-naphthalic anhydride  
Pigment Red 214  
Uvinul MS 40  
ASM ZKF  
Sodium Stearyl Fumarate  
(3-phenoxypropyl)phosphonic acid  
1,10-Dibromodecane  
Pigment Yellow 5  
Memantine hydrochloride  
3,5-Di-t-butyl-4-hydroxybenzyl N,N-dimethyldithiocarbamate  
6-hydroxymethyl-pyridine-2-carboxylic acid ethyl ester  
Codeine Phosphate  
Irganox 1035  
N-Benzylethylenediamine  
Tinuvin 292  
1-Methylthymine  
Isostearyl isostearate  
Crotonaldehyde  
(E)-Anethole  
1-Phenyl-2-propyn-1-ol  
Pentaerythritol Tetrabenzoate  
5-Nitro-1,10-phenanthroline  
1-Ethynyl-4-phenoxybenzene  
Dimethyl acetylmethylphosphonate  
4-tert-butylbenzotrile  
5-t-butyl 2,4-Dihydroxy benzophenone  
1,1,1-Triphenylsilylamine  
Pigment Yellow 108 | Pyramid-Yellow medium  
Pigment Red 194 | Permanent Red dark  
Tinuvin 120

---

**IS NIR Spectra – Complete Collection**

6049 NIR spectra

Copyright © NICODOM 2012

NICODOM Ltd., Hlavni 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

[ftirspectra@gmail.com](mailto:ftirspectra@gmail.com)

3-(4-methylphenyl)-1-phenyl-prop-2-en-1-one  
1-(4-methylphenyl)-3-phenylprop-2-en-1-one  
2,2',4'-Trichloroacetophenone  
4-Methoxybenzoin  
N-(4-Aminobenzoyl)-L-glutamic acid  
Palladium, [(2,3,5,6-h)-bicyclo[2.2.1]hepta-2,5-diene]dibromo-  
Pigment Orange 68  
1-(1-naphthyl)ethylamine  
1-Iodododecane  
2,3-Butanedione  
Methylthymol Blue  
N-Methyltaurine Sodium Salt  
Amitriptylinoxide  
1,2-dichlorotetramethyldisilane  
Ajmaline  
2-Bromomethyl-1,3-dioxolane  
Acid Blue 62  
Pigment Red 168 | Scarlet Red  
N-Formylmorpholine  
Glyoxal hydrate trimer  
Potassium phthalate  
Benzoin oxime = Cupron  
Chlorophenol Red  
2-(N-Morpholino)ethanesulfonic acid; MES  
Metronidazole  
6-Methylaminopurine  
Ammonium Citrate Monobasic  
Ruthenium, (acetato-O)dicarbonyl-  
Isonitrosocamphor  
1,4-Bis(trimethylsilyl)-1,3-butadiyne  
Pigment Yellow 17 based | Hostavinyl Yellow GG DL CU 005  
Food Red 1 | Ponceau SX  
a,a,a-Trifluoro-p-tolunitrile  
4-Fluorobenzoic acid  
2,2'-[2,5-bis(octyloxy)-1,4-phenylene]bis[4,4,5,5-tetramethyl-1,3,2-Dioxaborolane  
4-Fluorobenzaldehyde  
Dicyandiamide = Cyanoguanidine  
Hydantoin  
4-Trifluoromethylmercaptophenol  
6-Azauracil  
3-Aminopyridine  
N,N-Dimethylformamide dimethyl acetal  
Naloxone  
Basic Orange 22

---

**IS NIR Spectra – Complete Collection**

6049 NIR spectra

Copyright © NICODOM 2012

NICODOM Ltd., Hlavni 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

[ftirspectra@gmail.com](mailto:ftirspectra@gmail.com)

Piperidylmethylphosphonic Acid  
Acid Green 3 | Lissamine Green G  
N-Phthaloylglycine  
Cineole  
Propiolic acid  
Diphenyl Phosphite  
Oxamide  
Piperoin  
Diethyl 4,4'- biphenyldicarboxylate  
1,1'-Diethyl-4,4'-carbocyanine iodide  
N-Morpholinomethylphosphonate  
1,1-bis(4-hydroxyphenyl)butane  
Bathocuproine  
Chenodeoxycholic Acid  
Glycocholic acid  
Ethyl oxalyl chloride  
2,2-Bis(hydroxymethyl)propionic acid  
Safranin T  
1,8-diaminonaphtalene  
Tetraphenyl Cyclopentadienon  
Morin=2',3,4',5,7- Pentahydroxyflavone  
Phthiocol, 2-Hydroxy-3-methyl-1,4-naphthoquinone  
Flavianic acid hydrate  
Neocuproine  
Ninhydrin  
Bicuculline  
9-Fluorenone  
Indole-3-carboxaldehyde  
D-glycero-D-talo-heptitol= Volemitol  
Di-L-Valine Selenate Monohydrate  
Adonitol  
Methyl Phenylpropiolate  
3,5,N,N-tetramethylaniline  
Lepidine  
Furil  
1,2-cyclohexanedione dioxime  
Methyl Red-acid  
Ethyl maltol  
Acridine Orange hydrochloride hydrate  
Pigment Red 149  
Mercury(II) trifluoromethanesulfonate  
Fenofibrate  
Hippuric acid  
Indoline

---

**IS NIR Spectra – Complete Collection**

6049 NIR spectra

Copyright © NICODOM 2012

NICODOM Ltd., Hlavni 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

[ftirspectra@gmail.com](mailto:ftirspectra@gmail.com)

Bromisovale  
Arbutin  
2-Methylsuccinic acid  
Citraconic acid  
L-Djenkolic Acid  
2,6-pyridinedicarboxylic acid  
Guanidine Hydrochloride  
Calcium Lactobionate  
Hydrocortisone Acetate  
Phenobarbital  
Benzyl chloroformate  
Diphenylacetylene  
3-(4-Hydroxyphenyl)propionic acid  
Lactic Acid  
Cycloheptanol  
Cyclooctanone  
1,1-Bis(p-chlorophenyl)-2,2,2-trichloroethane = DDT  
Cefadroxil  
2-Thiohydantoin  
1,3-Cyclohexanedione  
2-Thiobarbituric acid  
4-Aminopyridine  
2-Aminopyridine  
Pigment Yellow 109 | Isoindole Yellow  
1,3-Propanediol  
2-Oxazolidinone-5-Spiro-4'-[N-(2-Phenylethyl)Piperidine] Hydrochloride  
Quinidine Sulfate  
Suberic acid  
1,4-diazacycloheptane  
12-Hydroxydodecanoic acid  
Cefaloridine  
Heptadecanoic acid  
Chloroquine Diphosphate  
Cerotic Acid  
Dimethylamine Hydrochloride  
Terfenadine  
D-(-)-Ribose  
Guanidine Nitrate  
D-Sorbitol  
Thioacetic acid  
Alloxan hydrate  
Borneol  
Acetylsalicylic acid  
L-Ascorbic acid, Vitamin C

---

**IS NIR Spectra – Complete Collection**

6049 NIR spectra

Copyright © NICODOM 2012

NICODOM Ltd., Hlavní 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

[ftirspectra@gmail.com](mailto:ftirspectra@gmail.com)

4,7-Diazafluoren-9-one  
Thymidine  
(1S,2R)-(+)-Ephedrine .HCl  
Ammonium Pyrrolidine Dithiocarbamate  
Benzimidazole  
Trimethyl phosphate  
2,4-Dinitrophenol  
4-Chloro-1,3-phenylenediamine  
1-Butoxy-2-Propanol  
L-Hydroxyproline  
Aliquat 336  
2,3-Butanediol  
Abietic acid  
Epinephrine Bitartrate  
Bismuth(III) acetate oxide  
Epinephrine  
Histamine  
(E)-2-Nitroethenylbenzene  
Homatropine hydrobromide  
Bisabolol  
Sodium Sulfanilate  
4,4'-Dinitro-2,2'-Bipyridine N,N'-Dioxide  
Neostigmine methyl sulfate  
DL-Valine  
Dopamine  
Hematoxylin  
Rubrene  
Bismuth Tribromophenate  
Ethyl Carbamate  
Uranine K  
2,2'-Dithiobis(ethylamine)  
Manganese Lactate Dihydrate  
Ethyl 4-iodobenzoate  
Triphenylmethane  
Diosmetin  
Dehydroacetic acid  
Cholesteryl Butyrate  
1,2-Naphthoquinone-4-sulfonic acid, sodium salt  
3-Aminophthalhydrazide  
5,7-dibromo-8-hydroxyquinoline  
Solvent Blue 49 | Irisol Fast Blue BLN  
L-Lysine acetate  
Pigment Yellow 75  
9,10-Dibromoanthracene

---

**IS NIR Spectra – Complete Collection**

6049 NIR spectra

Copyright © NICODOM 2012

NICODOM Ltd., Hlavni 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

[ftirspectra@gmail.com](mailto:ftirspectra@gmail.com)

Uvinul N 35  
Orange GGN, Acid Orange 20, alpha-Naphthol Orange  
Dichloro(p-cymene)ruthenium(II) dimer  
1-Nitroso-2-naphthol-3,6-disulfonic acid, disodium salt  
2-Bromo-2-nitro-1,3-propanediol  
Ferrocene, 1-(cyanomethyl)-2-(diphenylphosphino)-  
Kinetin  
Palladium, chloro(2,4-pentanedionato-O,O')(triphenylphosphine)-, (SP-4-3)- (9CI)  
Rewophat EAK 8190  
Tinopal SFR  
Monosodium tartrate  
Galactaric acid  
Sodium Gluconate  
2-(Methylthio)ethanol  
4-(2-Pyridylazo)resorcinol disodium salt, dihydrate  
3-Trimethylsilyl-2-propyn-1-ol  
Bathophenanthrolinedisulfonic acid, Na salt  
2-Furoylchloride  
Pigment Red 146 based  
Pigment Red 144  
Pigment Yellow 95  
Tinuvin 770 DF  
Cellobiose  
L-Cysteine HCl  
L-Cysteine  
Magneson II  
2,4,6-Trichlorophenylhydrazine  
o-Tolunitrile  
t-butyl Bromoacetate  
2-Iodoanisole  
Chloromethyltriphenylphosphonium chloride  
Ethyl 2-isocyanato-3-methylbutyrate  
Methylcyclam  
Esculin  
Benzidine dihydrochloride  
Benzidine sulfate  
N,N'-Diphenylbenzidine  
Dastib 242  
Naphthalene-2-sulfonic acid; sodium salt  
Sodium benzoate  
Sulfamic acid  
5-Hexynoic acid  
Sodium Hippurate Hydrate  
N-(Bromomethyl)phthalimide

---

**IS NIR Spectra – Complete Collection**

6049 NIR spectra

Copyright © NICODOM 2012

NICODOM Ltd., Hlavni 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

[ftirspectra@gmail.com](mailto:ftirspectra@gmail.com)

3-Methoxypropylamine  
(4-Aminophenyl)phosphonic acid  
3-Nitrophenylphosphonic acid  
Aleuritic acid  
2-Amino-1,3-propanediol  
Hydrazine dihydrochloride  
5-a-Androstane-3-a-ol-17-one  
Dehydroisoandrosterone  
2-Aminobenzyl alcohol  
3-Ethyl-3-methylglutaric acid  
4,6-Dinitro-o-cresol  
Glycine Hydroxamic Acid  
Benzophenone hydrazone  
DL-Cystathionine  
Diaveridine  
5-[4-(Dimethylamino)benzylidene]rhodanine  
N-Glycyl Aminoethylphosphonic Acid  
Pridinol Mesylate  
N,N-Dimethyl-p-phenylenediamine dihydrochloride  
m-Anisidine  
Carprofen  
2,6-Dibromoquinone-4-chloroimide  
N-Acetyl-L-tyrosine  
N-Benzylideneaniline  
Dibenzylideneacetone  
Indomethacin  
n-Butyl Cinnamate  
Benzyl sulfide  
1,3-Dicyclohexylcarbodiimide  
Isobutylbenzene  
4-Pentyn-1-ol  
2-Methoxy-1-naphthaldehyde  
Furylacrylic acid  
p-Xylylenediamine  
Camphoric acid  
Ethyl-3-Bromopropanoate  
Ethyl levulinate  
Cefaclor  
Cetaceum, Cetyl Palmitate 15  
Santoflex 15  
4-Methylaniline hydrochloride  
2,5-dichlorobenzenesulfonyl chloride  
4-Iodoaniline  
4-Iodophenol

---

**IS NIR Spectra – Complete Collection**

6049 NIR spectra

Copyright © NICODOM 2012

NICODOM Ltd., Hlavni 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

[ftirspectra@gmail.com](mailto:ftirspectra@gmail.com)

Bis(pentamethylcyclopentadienyl)zirconium dichloride  
2-Bromoethanol  
Basic Yellow 28 | Anilan Golden Yellow GL  
Ammonium formate  
Isooctane  
t-Butyl Acetate  
3-(N-Methylpyrrolidino)pyridine  
Di-tert-Butyl malonate  
DL-Tryptophan  
Ethyl chloroformate  
3,3-Dimethylacrylic acid  
5-Hydroxyindole-3-acetic acid  
1,3-Phenylenediamine dihydrochloride  
1-(Chloromethyl)-1h-benzotriazole  
Aniline sulfate  
Sodium Salicylate  
DL-2-Aminoadipic acid  
2-Hydroxyethyl acetate  
1-Iodobutane  
2-Acetamidoacrylic acid  
2,6-pyridinedicarboxaldehyde  
N-Acetylglycine  
Benzyl 2-bromoacetate  
Barium acetate hydrate  
Cadmium acetate dihydrate  
Isoamyl ether  
1-Chlorohexane  
Calcium formate  
Dibutyl sulfide  
Myristic acid  
n-Dotriacontane  
Nicardipine .HCl  
Butyl Cyanoacetate  
H Acid  
Thimerosal  
Pigment Red 255 | Irgazin Scarlet DPP EK  
2-(Acetylamino)-3-phenyl-2-propenoic acid  
Hydroxylamine hydrochloride  
1,3-Diphenylisobenzofuran  
Methyl Orange  
7-Iodo-8-hydroxyquinoline-5-sulfonic acid  
Ethyl Biscoumacetate  
3,3\*,5,5\*-Tetramethylbenzidine  
4-Pyridinecarboxylic acid hydrazide

---

**IS NIR Spectra – Complete Collection**

6049 NIR spectra

Copyright © NICODOM 2012

NICODOM Ltd., Hlavni 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

[ftirspectra@gmail.com](mailto:ftirspectra@gmail.com)

Basic Violet 3 | Gentian Violet  
Chromotrope 2B  
Amitriptyline hydrochloride  
Picronic acid  
a-Naphthyl salicylate  
Pigment Red 179 | Perrindo Maroon R-6424  
Trimethyl-beta-cyclodextrin  
Isonicotinic acid  
2-Nitrobenzaldehyde  
Protoporphyrin IX  
Neutral Red=3-Amino-7-dimethylamino-2-methylphenazine  
4,4'-bipyridyl  
1,4-Bis[4-benzoyl-3-hydroxy-6-(1-phenylethyl)]phenoxybutane  
Sodium methacrylate  
Dimethyloxalate  
2,4-Dichloroaniline  
Troysan Polyphase AF 3  
Triethylene Tetraamine Tetratosylate  
Atropine sulfate  
3-Nitrophenol  
Ammonium iron(III) oxalate hydrate  
Trimesic acid  
Disodium 5-Guanylate  
L-Methyldopa  
Tripalmitin  
Trimyristin  
4-Methylamino-phenol-sulfate  
Chlorhexidine  
Aluminum Ethoxide  
DL-Tyrosine  
1,4-Cyclohexanediol  
Glycylglycine  
Magnesium Stearate  
Stearate Zinc  
Zinc Undecylenate  
Zinc Caprylate  
Magnesium Formate  
Tetrabutylammonium phosphate  
Ethylamine hydrochloride  
Pigment Yellow 93  
Pigment Yellow 94  
Dihydroxypropyladenine (DHPA)  
1,4-Bis(4-benzoyl-6-tert-octyl-3-hydroxy)phenoxybutane  
2-Methoxybenzoin

---

**IS NIR Spectra – Complete Collection**

6049 NIR spectra

Copyright © NICODOM 2012

NICODOM Ltd., Hlavni 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

[ftirspectra@gmail.com](mailto:ftirspectra@gmail.com)

3-Methoxybenzoin  
2,3-Dimethoxybenzoin  
TPNB, Tripropylene Glycol n-Butyl Ether  
4-(4-Nitrophenyl)butyric acid  
Cacotheline  
4-Aminobutyric acid  
2,2-bis(4-hydroxy -3,5-di-methyl-phenyl) propane  
Cystamine dihydrochloride  
N,N',N''-tritosyldiethylenetriamine  
3,3'-Diaminodipropylamine  
Dibasic Lead Stearate  
Cefuroxime Sodium  
Dowanol EPH  
Guaiacol phosphate  
Tibolone  
Semicarbazide hydrochloride  
Bis(tributyltin) oxide  
Thallium(I) acetate  
N-Benzyl-N,N,N-triethylammonium chloride  
2-Bromopropionyl bromide  
3-Methyl-2-butanone  
Metoprolol Tartrate  
Glycine  
L-Alanine  
Doxycycline  
L-Serine  
1,4,7-triazacyclononane-N,N',N''-triacetic acid  
Dabsyl Chloride  
Pigment Yellow 138  
Chloramphenicol  
Heptadecafluorooctanesulfonic acid tetraethylammonium salt  
Glycine methyl ester hydrochloride  
DL-Glyceraldehyde  
3-(Trimethylsilyl)Propiolic Acid  
L-Aspartic acid  
L-Glutamine  
L-Glutamic acid  
L-Cystine  
Histamine Dihydrochloride  
Benzyltrimethylammonium chloride  
Mordant Violet 39  
Basic Red 9 | Fuchsin  
4,5-Dichlorophthalic acid  
Diethyl acetylmalonate

---

**IS NIR Spectra – Complete Collection**

6049 NIR spectra

Copyright © NICODOM 2012

NICODOM Ltd., Hlavní 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

[ftirspectra@gmail.com](mailto:ftirspectra@gmail.com)

2-Methyl-6-nitroaniline = 6-Nitro-o-toluidine  
Tricine  
6-Acetamidohexanoic acid  
Hexadecyl Trimethyl Ammonium Bromide  
Octadecanoic Acid, Stearic Acid  
Urea  
1,1,1-Trichloro-2-methyl-2-propanol  
1,4-Dimethylnaphthalene  
1,5-Dimethylnaphthalene  
Pamidronate  
2,6-Dinitrophenol  
Congo Red  
Calcium Acetate  
Levulinic Acid Calcium Salt  
Barbital  
Ultramarine Blue  
D-(-)-Fructose  
Sucrose (powder), Sugar, Confectioner's  
Propylene Glycol  
Semicarbazide  
N-Acetyl-Aminomethylphosphonic Acid  
Sulfaguanidine  
Sulfamethazine  
2,3-Butanedione monoxime  
2-Bromoanisole  
Ergosterol, Provitamin D  
Cholesterol  
L-Asparagine \*H<sub>2</sub>O  
Benzoin  
Pigment Yellow 128  
o-Anisic acid  
Caffeine  
2,4'-bipyridyl  
Aminopyrine  
2,7-Dihydroxynaphthalene  
1-Amino-2-phenylpropane  
Triphenylaminomethane  
Menadione, Vitamin K3  
Promethazine hydrochloride  
2-mercaptobenzimidazole  
1-Bromo-2-Iodobenzene  
2,5-Dichlorophenol  
1-Bromo-2,4-dinitrobenzene  
Tetrabenazine

---

**IS NIR Spectra – Complete Collection**

6049 NIR spectra

Copyright © NICODOM 2012

NICODOM Ltd., Hlavni 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

[ftirspectra@gmail.com](mailto:ftirspectra@gmail.com)

2-Methyl-m-phenylene diisocyanate  
Acid Green 20 | Egacide Black HG  
Theophylline = 1,3-Dimethylxanthine  
Acid Red 4  
Lactitol  
Dioctyl-2,5-dibromoterephthalate  
Adenosine  
Acid Blue 15 | Coomasie Blue FF  
Inosine  
4-Bromobenzoic acid  
1-Bromo-4-nitrobenzene  
2-Pyridilcarbinol  
3-Chlorobenzaldehyde  
DL-10-Camphorsulfonic acid  
Papaverine  
Diphenylmercury(II)  
Pyridine 4-Phosphonic Acid  
D-(+)-Xylose  
4-Bromobenzyl bromide  
4-Bromobenzyl chloride  
a-Tocopheryl acetate hydrate  
1-Bromo-4-Iodobenzene  
4-Methylcyclohexanone  
Potassium Sorbate  
DL-a-Tocopherol, Vitamin E  
Potassium formate  
2-Methyl-2-pentanol  
Betaine hydrochloride  
Methylcyclohexanol  
Isovaleraldehyde  
Barium succinate  
Guanyurea Sulfate  
1-Acetyl-2-thiourea  
3-Bromoaniline  
3-Aminophenol  
3-Methoxybenzaldehyde  
trans-1,2-Diaminocyclohexane  
Iodobenzene  
N-(2,6-Dimethylphenylcarbonylmethyl)iminodiacetic acid  
4-Pentenoic acid  
D-Galactose  
1-Heptene  
Acyclovir  
Butyl formate

---

**IS NIR Spectra – Complete Collection**

6049 NIR spectra

Copyright © NICODOM 2012

NICODOM Ltd., Hlavni 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

[ftirspectra@gmail.com](mailto:ftirspectra@gmail.com)

Allyl sulfide  
L-Histidine Hydrochloride Monohydrate  
Methylamine Hydrochloride  
Guanidine carbonate  
Thiamine chloride  
2,3-Dimethyl-2-butanol  
Procaine  
2,2,2-Trichloroacetamide  
Citric Acid Monohydrate  
DL-Methionine  
Ferrocene, 1,1',2,2',3,3',4,4'-octamethyl-  
Tetraphenyllead  
Tetraphenyltin  
o-Cresolphthalein  
Acetazolamide  
Nicotinic acid = Niacin  
Zinc acetate dihydrate  
Ammonium oxalate-1-hydrate  
Dibutylsulfon  
DL-2-Bromopropionic acid  
2-Chloropropionic acid  
Phenylhydrazin hydrochloride  
DL-Lactic acid  
a-D-Lactose Monohydrate  
5-(3-Thienyl)-2H-1,2,3,4-Tetrazole  
3,4-Dihydroxy-L-Phenylalanine  
[(Methylimino)dimethanediy]bis(phosphonic acid)  
4-Cumylphenol  
Ethylenediaminetetraacetic acid = EDTA  
Tributylin  
4-Aminoazobenzene  
Diphenylthiocarbazone  
Methyl Yellow  
DL-Amphetamine sulfate  
1,4,7-Oxadiazonane-ditosylate  
Dowex  
Acid Brown 355  
L-Tyrosine  
Tyramine hydrochloride  
1,4,7,10-Tetraazacyclododecane-1,4,7,10-tetrakisacetic acid  
2-Mercaptoethanol  
Diethyl ether  
7-(Beta-Hydroxypropyl)Theophylline  
3-Nitrophthalic acid

---

**IS NIR Spectra – Complete Collection**

6049 NIR spectra

Copyright © NICODOM 2012

NICODOM Ltd., Hlavni 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

[ftirspectra@gmail.com](mailto:ftirspectra@gmail.com)

6-Aminohexanoic acid  
Triphenylarsine  
Linoleic Acid  
Triphenylamine  
Triphenyl Phosphine  
Triphenylantimony  
Rosolic acid  
Diphenylurethane  
Acetamide  
1-Methylindole  
2,3-Dimethylbenzoic Acid  
DL-Homocysteine thiolactone hydrochloride  
Pigment Red 2  
Cholesteryl Benzoate  
Cholesteryl Acetate  
a-Naphthoflavone  
Cupric acetate .H2O  
Dansyl chloride  
N-Alanyl-aminomethylphosphonic acid  
N-Leucyl-aminomethylphosphonic acid  
N-Alanyl-Aminoethylphosphonic Acid, SR Isomer  
N-Hydroxysuccinimide  
N-Leucyl-Aminoethylphosphonic Acid, SS Isomer  
Ethylenediamine Tetramethyl Phenylphosphinate  
p-Cresylsalicylate  
Lead(II) acetate trihydrate  
Indoxyl Acetate  
Phlorizin  
1,3,5-Trimethyl-2,4-dinitrobenzene  
Dulcitol  
Dimethyl-L-tartrate  
2-Fluoropropionic acid  
L-(-)-Xylose  
Diethyl methylmalonate  
Ethyl 2-methylacetoacetate  
Ethyl 2-chloroacetoacetate  
Ruthenium, bis(2,2'-bipyridine-N1,N1')dichloro-, hydrate (1:2)  
1,3,5-Trihydroxybenzene dihydrate  
3,5-Dinitrosalicylic acid  
Tripotassium Citrate Monohydrate  
2,3,4-Trihydroxybenzoic acid  
4-Nitrophthalic acid  
Pigment Yellow 151  
Acid Blue 90 | Supranol Cyanine G

---

**IS NIR Spectra – Complete Collection**

6049 NIR spectra

Copyright © NICODOM 2012

NICODOM Ltd., Hlavní 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

[ftirspectra@gmail.com](mailto:ftirspectra@gmail.com)

Disodium succinate hexahydrate  
Disodium tartrate, dihydrate  
Methyl 2-iodobenzoate  
2-Carbethoxycyclopentanone  
Methyl Furoate  
2-Ethyltoluene  
Solvent Brown 43 | Orasol Brown 2RL  
Cinchocaine.HCl  
Sumilizer GM-S  
Quinine Hydrochloride Dihydrate  
Quinine Sulfate Dihydrate  
1-Ethyl-2-nitrobenzene  
o-Tolidine dihydrochloride  
1,2,4-Triacetoxybenzene  
2,2,7,7-Tetramethyl-3,5-octadiyne  
Sodium Acetate Trihydrate  
Sodium Citrate Dihydrate  
Amoxicillin  
2,4-Dimethoxybenzaldehyde  
Resacetophenone oxime  
2-Naphthyl salicylate  
Trehalose dihydrate  
Benzoylhydrazine  
Benzoylacetone nitrile = beta-Oxohydrocinnamone nitrile  
Ethyl nicotinate  
Cobalt(II) acetate tetrahydrate  
Quinidine Hydrochloride, Monohydrate  
Orcinol monohydrate  
Oxalic Acid Dihydrate  
2-Iodotoluene  
1,2-Diiodobenzene  
2-Iodoaniline  
2,5-Diaminotoluene sulfate  
Saccharin Sodium Dihydrate  
Chlorohydroquinone  
2-Bromohexanoic acid  
Methyl glycinate  
Tri(2-tolyl)phosphine  
2-Pyrrolidinone  
1-Methylimidazole  
Acetylcysteine  
Cloxacillin  
Direct Blue 85  
Methylene Blue

---

**IS NIR Spectra – Complete Collection**

6049 NIR spectra

Copyright © NICODOM 2012

NICODOM Ltd., Hlavni 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

[ftirspectra@gmail.com](mailto:ftirspectra@gmail.com)

Ethyl pyruvate  
DL-Malic acid  
Oil, Yoyoba Oil  
Cocoalkylamine ethoxylate  
Solvent Red 33 | Luxol Fast Red B  
Acid Yellow 194  
Pigment Red 221 | Cromophtal Red 2B  
3-Amino-1,2,4-triazole  
Isodecylalcohol ethoxylate  
Ethyl-3-nitrobenzoate  
L-Leucine  
2-Chloro-4-nitrophenol  
Acid Violet 90  
p-Toluenesulfonic acid monohydrate  
Solvent Yellow 88 | Orasol Yellow 2GLN  
Eriochrom Red G  
Pigment Red 210  
Methyl 4-bromobenzoate  
Direct Brown 218 | Saturn Brown L2R  
4-Iodobenzoic acid  
Basic Blue 95 | Remacryl Blue B  
Direct Yellow 132 | Rybanyl Yellow N-RLT  
Solvent Yellow 89 | Orasol Yellow 2RLN  
2-Ethylhexyl 2-cyano-3,3-diphenylacrylate  
Santowhite L  
5-Bromo-m-xylene  
Tribenzylamine  
N-(p-Methylphenyl)-N-(p-nitrophenyl) amine  
Cinnamic Acid  
o-Benzyl Carbamate  
N,N-Dimethyl-p-phenylenediamine sulfate  
2-Benzyloxyethanol  
Phenethyl chloride  
4-Nitrobenzoic acid  
9-Fluorenone-4-carboxylic acid  
2-Morpholinoethanol  
Hordenine Sulfate  
Acid Red 112 | Ponceau S extra  
Solvent Red 30 | Alcovar Scarlet CSD  
4-Bromophenylhydrazine hydrochloride  
Triphenylsilylacetylene  
3-Diethylamino-1-propanol  
4-Bromobenzonitrile  
4-Chlorobenzonitrile

---

**IS NIR Spectra – Complete Collection**

6049 NIR spectra

Copyright © NICODOM 2012

NICODOM Ltd., Hlavní 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

[ftirspectra@gmail.com](mailto:ftirspectra@gmail.com)

Alpha, alpha'-Dibromo-p-xylene  
a, a'-Dichloro-p-xylene  
Terephthalonitrile  
Terephthalaldehyde  
3-Hexanol  
Methyl crotonate  
Ethyl propiolate  
Ethyl glycolate  
Phenylmercuric acetate  
Dowanol PGDA  
Diethyl fumarate  
p-Phenylenediamine dihydrochloride  
1,4-Diiodobenzene  
Dimethyl Maleate  
Dimethyl fumarate  
Pentyl propionate  
Ethyl Methanesulfonate  
2-Methyl-2-hexanol  
Aniline  
Calcium Diacetate  
Ethylurea  
Thioacetamide  
Thiourea  
Captopril  
2-Aminoisobutyric Acid  
2,5-Dimethyl-1H-pyrrole  
2,5-Dimethylfuran  
m-Phthalodinitrile  
m-Phthalaldehyde  
Bromocresol Green  
Trimethylamine N-oxide dihydrate  
1,3,5-tribromobenzene  
3,5-Dichloroaniline  
3-Methylglutaric acid  
1,3-diphenylprop-2-en-1-ol  
2-Hexanol  
1,4-Pentanediol  
1-Nitrobutane  
Dibenzoyl-L-tartaric acid Monohydrate  
3-Chloro-1-propanol  
Methyl propargyl ether  
D-Arginine hydrochloride  
Pyridine hydrochloride  
N-(Hydroxymethyl)benzamide

---

**IS NIR Spectra – Complete Collection**

6049 NIR spectra

Copyright © NICODOM 2012

NICODOM Ltd., Hlavni 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

[ftirspectra@gmail.com](mailto:ftirspectra@gmail.com)

4-Chlorobutyronitrile  
Meglumine  
2-Chlorooctane  
Amyl Acetate  
Ethyl palmitate  
1,6-Dibromohexane  
1-Bromoheptane  
1,6-Hexandiol  
Ethylene glycol diformate  
1,7-Heptanediol= Heptamethylene glycol  
1-Hexadecene  
Heptadecane  
2-(4-Pentynyloxy) tetrahydro-2h-pyran  
1-Octacosane  
Hexatriacontane  
t-Butyl cyanide  
1,4-Dibromo-2,3-butanedione  
Tributylammonium chloride  
Isopentyl dodecanoate  
(4-tolylethyl)methyldichlorosilane  
Poly(methylhydrosiloxane)  
1,2-bis(tosyloxy)ethane  
Ammonium Acetate  
Basic Red 12 | Astra Phloxine G  
D-Threonine  
1,1,3,3-Tetramethylurea  
1-Amino-2-phenyl-ethylphosphonic acid  
Rose bengal  
Tetrabromophthalic anhydride  
Fuchsine Basic  
Cholesteryl Propionate  
Primid - XL - 552  
Calcium DL-pantothenate  
D-(+)-Lactose Powder  
Methyl red Hydrochloride  
2,7-Dibromo-9,9-bis[4-(octyloxy)phenyl]-9H-Fluorene  
Pyrrole-2-carboxylic acid  
Bilirubin  
Pigment Violet 23  
Pigment Yellow 74  
Acid Yellow 11 | Egacide Yellow G 140  
Acid Yellow 25 | Egacide Yellow R  
Yellow 2G  
Acid Red 37

---

**IS NIR Spectra – Complete Collection**

6049 NIR spectra

Copyright © NICODOM 2012

NICODOM Ltd., Hlavni 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

[ftirspectra@gmail.com](mailto:ftirspectra@gmail.com)

2,4,5-Trichloroaniline  
Acid Brown 365  
Clopamide  
D-(+)-Malic acid  
2-Thiophenemethanol  
N,O-Ditosylethanolamine  
L-Methionine  
Fat Red 7B  
Aluminium Stearate  
Acid Orange 3 | Anthralan Yellow RRT  
Sulfanilamide  
1-Chloro-4-iodobenzene  
1,4-Cyclohexanedione  
2,5-Dimethylthiophene  
Potassium sodium tartrate tetrahydrate  
Disodium edetate dihydrate  
(+)-S-Carboxymethyl-L-cysteine  
DL-DOPA = 3-(3,4-dihydroxyphenyl)- DL-alanine  
Tinuvin 144  
1-Iodohexane  
Metilox  
Direct Green 26 based  
1,2-Dipalmitoyl-SN-glycero-3-phosphocholine  
L-Phenylalanine  
Triphenyltin chloride  
Tetrasodium Edetate  
D-(+)-Lactose, Monohydrate  
D-Valine  
Solvent Blue 63 | Ceres Blue GN  
Acid Blue 25 | Lissamine Blue RB  
Pigment Red 8  
Pigment Red 12 | Monolite Rubine 2R  
Pigment Red 9  
Phenylurea  
Ethanol  
Formic acid  
Tetramethylammonium bromide  
Direct Yellow 34  
2,2-bis(4-hydroxy -3-methylphenyl)butane  
9-Anthraldehyde  
Direct Black 19 | Chlorazol Black BV  
Hostavin (tm) N20  
o-Phthaldialdehyde  
4-Phenyltoluene

---

**IS NIR Spectra – Complete Collection**

6049 NIR spectra

Copyright © NICODOM 2012

NICODOM Ltd., Hlavni 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

[ftirspectra@gmail.com](mailto:ftirspectra@gmail.com)

Acid Violet 1  
Acid Red 35  
Meclofenamic Acid  
Cefotaxime sodium  
cis-Stilbene  
1,3-Dioxolane  
Tetracosane  
Iodoacetic acid  
Direct Green 28  
Pigment Red 23 | Irgalit Red RBS 2344  
Morine  
Tolbutamide  
Solvent Red 18 | Ceres Red 3R  
Acid Blue 103 | Brilliant Indo Blue 5G  
Pigment Yellow 3  
Potassium oxalate monohydrate  
Hydroxy Pentoxifylline  
1-(3-Carboxypropyl)-3,7-Dimethylxanthine  
Sodium trichloroacetate  
Indole-3-acetic acid sodium  
Irganox 1425  
Yohimbine .HCl  
Pigment Yellow 183 | Paliotol Yellow K2270 VPL 2269  
1-(2,3,4,5,6-pentafluorophenyl)ethanone  
Pyridoxine, Vitamine B6  
Nicotine Hydrogen Tartrate  
Solvent Red 3 | Fat Brown B  
Pigment Red 112  
2-(diphenylphosphino)benzenamine  
Cytidine  
2,4,6-Triisopropylbenzenesulfonyl chloride  
D-Glucuronic acid  
Acridine orange hydrochloride hydrate  
L-Lysine Monohydrochloride  
2-Chloro-3-(trifluoromethyl)pyridine  
Sodium p-toluenesulfonate  
Ferrocene, hexadecyl- (9CI)  
Benzoic acid  
N-(2-chlorobenzoyloxy-carbonyloxy)succinimide  
L-Norvaline  
N-Alanyl-Aminoethylphosphonic Acid, SS Isomer  
2-Chloronicotinonitrile  
Diethylenetriamine-N'-methylene(phenyl)phosphinic-N,N,N'',N''-tetraacetic acid

---

**IS NIR Spectra – Complete Collection**

6049 NIR spectra

Copyright © NICODOM 2012

NICODOM Ltd., Hlavni 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

[ftirspectra@gmail.com](mailto:ftirspectra@gmail.com)

Diethylenetriamine-N'-methylene(dibenzylaminomethyl)-phosphinic-N,N,N'',N''-tetraacetic acid  
N- {2-[(2-aminoethyl)(benzyl)amino]ethyl}-4-methylbenzenesulfonamide  
1-Docosanol  
Uracil  
Acetylcholine bromide  
Methyl Methanesulfonate  
Alendronic acid  
4-Nitro-1,8-naphthalic anhydride  
N-Methionyl-Aminoethylphosphonic Acid  
2-methyl-1-(trimethyl-silyloxy)-1-propene  
Mercaptoacetyltriglycine  
Amantadine .HCl  
o-Phenanthroline  
1,8-Diazabicyclo[5.4.0]undec-7-ene  
Oxacillin  
D-(+)-Glucosamine .HCl  
Chromium(III) acetate dihydrate  
2-Naphthaldehyde  
Thiamine hydrochloride (Vitamin B1)  
Tetracyanoethylene  
L-Sorbitol  
2,2'-Dimethoxybenzoin  
Pentabenzyl cyclopentadiene  
1-Phenyl-1-Propyne  
1,3-bis(diphenylphosphino)propane  
(Hexamethylbenzene)ruthenium(ii) dichloride dimer  
Diethylenetriaminepentaacetic acid  
Ferrocene, [(1E)-2-carboxyethenyl]-  
Barbituric acid  
Methyl alcohol  
Solvent Yellow 14 | Sudan I  
2-Propanol  
Acetone  
Chloroform  
Silwet L-7608  
Dimethylsulfoxide  
10- {[di(hydroxy)phosphoryl]methyl}-1,4,7,10-tetraazacyclododecane-1,4,7-triacetic acid  
10- {[ (4-nitrobenzyl)(hydroxy)phosphoryl]methyl}-1,4,7,10-tetraazacyclododecane-1,4,7-triacetic acid  
Cyasorb UV - 2908  
N-(tert-butoxycarbonyl)-L-norvaline dicyclohexylammonium  
Fluorad FC 120  
Hexamethylphosphoramide

---

**IS NIR Spectra – Complete Collection**

6049 NIR spectra

Copyright © NICODOM 2012

NICODOM Ltd., Hlavní 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

[ftirspectra@gmail.com](mailto:ftirspectra@gmail.com)

Sodium Citrate  
Tallow Oil Fatty Acid  
5-t-butyl-2-hydroxy-4-butoxybenzophenone  
Thioglycolic Acid  
N,N-Dimethylformamide  
Cocamidopropyl amine oxide  
Tetramethyl orthosilicate  
Pigment Red 233 | Pinkcolor Deep  
Cyanocobalamine (Vitamin B12)  
Diethyl vinylphosphonate  
TAPSO, 2-Hydroxy-N-(tris(hydroxymethyl)methyl)-3-aminopropanesulfonic acid  
Pigment Orange 60 | Hostaperm Orange HGL  
Igepal CO-990  
Silwet L-7500  
Petroleum ether  
Pigment Yellow 150  
Pigment Yellow 155  
Acid Blue 102 | Supranol Blue GL  
Pridinol Methanesulfonate Anhydrous  
D(+)-Galacturonic acidmonohydrate  
Diethyl bromomalonate  
N,N-Diethylacetamide  
Wingstay L  
Ethyl Lactate  
AgeRite Superlite  
Tributyltin hydride  
Dipyrone  
Oil, Caraway oil  
Silwet L-7644  
Silwet L-7200  
Silwet L-7600  
Hypoxanthine  
Propylphosphonic anhydride  
Eosin Methylene Blue, Wrights stain  
Chlorpromazine Chloride  
Tripropargylamine  
2,2',2'',2'''-[(hydroxyphosphoryl)bis(methanediylnitrilo)]tetraacetic acid  
2-Bromo-3-hexylthiophene  
N-benzylcinchonidinium chloride  
Dodecanedioic acid  
Distearyl thiodipropionate  
1-Bromononane  
2-Methylimidazole  
5,10,15,20-tetrakis(4-(trimethylammonio)-phenyl)-21h,23h-porphine tetratosylate

---

**IS NIR Spectra – Complete Collection**

6049 NIR spectra

Copyright © NICODOM 2012

NICODOM Ltd., Hlavni 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

[ftirspectra@gmail.com](mailto:ftirspectra@gmail.com)

g-Caprolactone  
Ampicillin sodium  
Ampicillin  
2-(Chloromethyl)pyridine hydrochloride  
5-Iodouracil  
4-Pyridinealdoxime  
D-Mannitol  
4-Iodoanisole  
Ferrocene, (1-oxohexadecyl)-  
Salicylic acid  
Maltose  
Pigment Red 175 | Purple-Red  
5-Nitro-2-furaldehyde  
Xanthine  
Pigment Brown 25  
2,3-Pyridinedicarboxylic anhydride  
2-Adamantanone  
Mephenoxalone  
2-Bromoacetophenone  
g-L-Glutamyl-L-cysteinylglycine  
L-Ornithine  
Phenol, 2-(2H-benzotriazol-2-yl)-4,6-bis(1-methyl-1-phenylethyl)  
2,4-Dinitrofluorobenzene  
4-tert-Butyl-4'-methoxy-dibenzoylmethane  
L-Asparagine  
5-Chloro-1,3-dimethoxybenzene  
4'-Trifluoromethylphenyl acetylene  
p-Toluenesulfonamide  
4'-Hydroxypropiophenone  
Benzyloxycarbonyldiethanolamine  
Cocoamidopropyl Hydroxysultaine  
N,N-Diisopropylethylamine  
L-Histidine  
Meloxicam  
N,N'-dibenzyl-N,N'-ditosylethylenediamine  
1-Propanol  
Cytosine  
Diethyl Dinitrobenzylmalonate  
Benzene  
Cholesteryl chloroformate  
Copper Tosylate  
Natural Red 9 | Madder Lake  
1,1,1-Trichloroethane  
3-Thiophenemethanol

---

**IS NIR Spectra – Complete Collection**

6049 NIR spectra

Copyright © NICODOM 2012

NICODOM Ltd., Hlavni 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

[ftirspectra@gmail.com](mailto:ftirspectra@gmail.com)

5-Ethynyl-1-methyl-1H-imidazole  
Ampicillin Trihydrate  
Basic Orange 12  
Lincomycin Hcl  
Sodium Glucuronate  
Pigment Yellow 127 | Permanent Yellow GRL  
Beta-Arthemeter  
Sulfathiazole  
DL-Lactic acid, Sodium salt  
L-Valine  
L-Threonine  
[1,1'-Bis(diphenylphosphino)ferrocene]dichloropalladium(II)  
5-Bromo-4-chloro-3-indolyl b-D-galactopyranoside  
Alizarin = 1,2-Dihydroxyanthraquinone  
2,2'-Dipyridyl N,N'-dioxide  
4,4'-Di-O-tert-butyl-2,2'-bypyridine  
Nonox WSO  
N-(p-methoxyphenyl)-N-(p-nitrophenyl) amine  
Deoxyribonucleic acid  
2,5-Dibromo-3-carboxythiophene  
L-Tryptophan  
Adenine  
L-Isoleucine  
Guanine  
4,4'-Dinitrobibenzyl  
N-(2-Hydroxyethyl)piperazine-N'-2-ethanesulfonic acid; HEPES  
Dibenzylphosphinic acid  
4-Chloropyridinium chloride  
4-Hydroxycinnamic acid  
4-Chlorobenzoic acid  
3,3-Diaminobenzidine tetrahydrochloride  
Dibenzyl-cis-3a,5a,8a,10a-tetraazaperhydropyrene  
Plasticizer, Tegosoft EE  
trans-1,2-Dibromocyclohexane  
Ferrocene, 1-[(dimethylamino)methyl]-2-(diphenylphosphinothioyl)-  
N,N-Diphenyl-p-phenylenediamine  
Isobutyl-2-ethyl hexanoate  
4-(4-Nitrophenylazo)resorcinol  
Calix[4]arene  
3,3'',4,4'',5,5''-Hexamethoxybenzoin  
L-Arginine  
Benzyl-N-acetyl-4,6-O-benzylidene muramic acid  
Methyl iodide  
Methylamine water

---

**IS NIR Spectra – Complete Collection**

6049 NIR spectra

Copyright © NICODOM 2012

NICODOM Ltd., Hlavní 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

[ftirspectra@gmail.com](mailto:ftirspectra@gmail.com)

Dimethylaminoborane  
Dibromomethane  
Dipropyl-cis-decahydro-2a,4a,6a,8a-tetraazacyclopent[fg]acenaphthylene  
Acetonitrile  
Dichloromethane  
Diiodomethane  
Formamide  
Isocyanic acid  
2-Bromopropane  
4-Methylmorpholine N-oxide  
2-Iodopropane  
Lovastatin  
Isobornyl methacrylate  
Acetyl chloride  
Trimethylacetamide  
Iodoform  
Nitromethane  
Pigment Orange 62  
t-Butylamine  
t-Butyl alcohol  
Dimethyl methylphosphonate  
Folic Acid (Vitamin B11)  
Methanesulfonic acid  
4-Octylpyrocatechol  
Neopentyl alcohol  
2-Methyl-2-butanol  
Acetone cyanohydrin  
2,2,2-Trifluoroethanol  
1,4,8,11-Tetraazatricyclo[9.3.1.1(4,8)]hexadecane  
Diethyl t-butylmalonate  
Diethyl-Nitrobenzylmalonate  
1,8-dimethyl-1,4,8,11-tetraazatricyclo(9.3.1.14.8)hexadecanium(2+) diiodide  
Trichloroacetic acid  
Pinacol  
Enalapril Maleate  
Pigment Yellow 147  
Pigment Yellow 185  
Diethyl phosphite  
DL-Camphor  
Dimethylacetylene dicarboxylate  
Alloxantin  
Ethyl 3-ethoxypropionate  
Methoxyflurane  
Enalaprilat

---

**IS NIR Spectra – Complete Collection**

6049 NIR spectra

Copyright © NICODOM 2012

NICODOM Ltd., Hlavni 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

[ftirspectra@gmail.com](mailto:ftirspectra@gmail.com)

2-Hexyne  
L-Arabitol  
Fumarodinitrile  
Phosphonic acid, P,P'-(3-aminopropylidene)bis-  
1-Pentadecyne  
Ethylenediamine tetra(methylenephosphonic acid)  
2',7'-Dichlorofluorescein  
3-Methyl-1,2-cyclopentanedione  
Bromothymol Blue  
Bromocresol Green  
Thymol Blue  
t-Butyl Glycidyl Ether  
Trimethylsilyl cyanide  
2-Ethynylanisole  
Famotidine  
Chlorodimethylphenylsilane  
Triphenylmethyl Chloride  
Triphenylsilyl chloride  
Hydroxytriphenyltin  
1-Bromoadamantane  
1-Aminoadamantane  
1-Adamantanol  
Benzilic acid  
1,3-Dimethylbarbituric acid  
Pyrophosphoric Acid Tris(Tetrabutylammonium) Salt  
1-phenoxy-2-propanol  
Phenolphthalein  
Indole-3-carboxylic acid  
2,3,4,5,6-Pentafluoroaniline  
4-tert-Butylphenylacetylene  
Cholesteryl Valerate  
Butobarbital  
Water  
Chlorthalidone  
Chloroxine  
2,2-Bis(4-hydroxyphenyl)butane  
1,3-dibromo-5,5-dimethylhydantoin  
Dibutyltin dilaurate  
Nonox WSL  
Nonox WSP  
Bromodiphenylmethane  
5,5-Dimethylhydantoin  
Dicyclopentadiene  
Pigment Yellow 180 | Permajet Yellow HG

---

**IS NIR Spectra – Complete Collection**

6049 NIR spectra

Copyright © NICODOM 2012

NICODOM Ltd., Hlavni 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

[ftirspectra@gmail.com](mailto:ftirspectra@gmail.com)

1,1,1-Tris(hydroxymethyl)ethane  
(1S)-(-)- $\alpha$ -Pinene  
Naugawhite  
Citric Acid  
Triethyl citrate  
(-)-Quinic acid  
o-(Phenyliminomethyl)phenol  
Tetraethylammonium hydroxide (20% solution in water)  
N-Z-ethanolamine  
1,3,5-Triacetylbenzene  
1,1,1-tris(hydroxymethyl)propane  
Benzyl-cis-3a,5a,8a,10a-tetraazaperhydropyrene  
Sulfamide  
Estabex EN  
Aminomethyldiphosphonic Acid Copper Complex  
Paroxetine.HCl  
1-Ethynyl-1-cyclohexanol  
Diethyl ethylphosphonate  
Triethyl orthoacetate  
9-Acetylanthracene  
2,2'-Azobis(isobutyronitrile)  
1,2-Dibromopropane  
Isobutylbromide  
2-Methyl-1,3-butadiene  
Isobutyronitrile  
2-Methyl-1-propanol  
1,2-Diaminopropane  
Triphenylsilane  
2-Butanone  
3-Buten-2-on  
Chloroacetone  
1-Amino-2-Propanol  
Methylglyoxal  
Trichloroethylene  
Propionamide  
Acrylamide  
Chloroacetamide  
Ethanecarboxylic acid  
Triphenylphosphine oxide  
Triphenylsilanol  
Glycolic acid  
Dinitrobenzyl-cis-decahydro-2a,4a,6a,8a-tetraazacyclopent[fg]acenaphthylene  
N-Methylacetamide  
Thiosemicarbazide

---

**IS NIR Spectra – Complete Collection**

6049 NIR spectra

Copyright © NICODOM 2012

NICODOM Ltd., Hlavní 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

[ftirspectra@gmail.com](mailto:ftirspectra@gmail.com)

Nitroethane  
Isobutanoic acid  
Flexzone 7 L  
L-Lactic acid  
1,1,2,2-Tetrachloroethane  
Cefixime  
Dithiooxamid  
Methacrylic acid  
Dichloroacetic acid  
2-Nitropropane  
Santovar A  
N-Phenylalanyl-Aminoethylphosphonic Acid  
1-Ethynyl-4-pentylbenzene  
2,2-Bis(4-hydroxy-3-methylphenyl)propane  
Stearate Calcium  
Caffeine Sodium Benzoate  
Ceresin  
Solvent Yellow 33 | Waxoline Yellow T  
Mordant Brown 90 | Korostan Brown G 140  
Direct Yellow 44  
Bisphenol A  
Basic Brown 4 | Bismarck Brown R  
Wax, Yellow Bees Wax  
Pigment Yellow 53 | Sicotan Yellow K1011  
4,4'-Dichlorodiphenyl sulfone  
4,4'-Diaminodiphenylsulfone  
N-Methyl-N-nitrosotoluene-4-sulphonamide  
Acid Brown 75 | Korostan Brown ER 40 TEK.  
Cefpodoxime  
Tamsulosin hydrochloride  
Mineral Oil, light and heavy  
Tetraethyl (N,N-Dibenzyl)aminomethyl-bis(phosphonate)  
Methyl p-toluenesulfonate  
Rosin  
Veratrine  
1,8-diamino-p-menthane  
Pinene  
2-Bromobutyric acid  
Trans-2,3-Dimethylacrylic Acid  
Dihydrocholesterol  
a-Naphthyl acid phosphate, Na salt  
Dexpanthenol  
2-Amino-1-naphthalenesulfonic acid  
Dehydrocholic acid

---

**IS NIR Spectra – Complete Collection**

6049 NIR spectra

Copyright © NICODOM 2012

NICODOM Ltd., Hlavni 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

[ftirspectra@gmail.com](mailto:ftirspectra@gmail.com)

Cholic acid  
Calcium L-Lactate  
Diethyl chlorophosphate  
Calcium Lactate  
Solvent Violet 13 | Waxoline Violet A  
Strontium oxalate  
Trimethylpyruvic acid  
Purpurin  
Quinizarine  
2-(2,5-dibromophenyl)-5-phenyl-1,3,4-Oxadiazole  
Warfarine  
Diethyl glutarate  
1,8-Naphthalic anhydride  
1-Naphthylacetic acid  
Rhodamine B  
1-(Pyridine-N-Oxide-2yl-methyl)-4,7,10-tris(carboxymethyl)-Cyclene  
Phenosafranine  
Xanthene 9-carboxylic Acid  
1-Triisopropylsilyl-1-propyne  
Pigment Yellow 194  
Oxfloxacin  
Calix(8)Aren  
1,5-Dichloroanthraquinone  
Phenylphosphonic dichloride  
Sodium toluene-4-sulphinate  
Quinapril Hydrochloride  
Raloxifene Hydrochloride  
Styphnic acid  
2-Vinylnaphthalene  
1-Naphthylamine-8-sulfonic acid  
8-Anilino-1-naphthalenesulfonic acid  
4-(Trifluoromethoxy)phenol  
1-Adamantanecarboxylic acid  
4-Nitrophenyl acetate  
Cyclododecanone  
4-Aminoantipyrene  
3-Indolepropionic acid  
Acenaphthene  
1-Indanone  
Sitosterol  
1,5-Dihydroxynaphthalene  
Theobromine  
Jeffamine M-1000  
2-Hydroxy-1,4-naphthoquinone

---

**IS NIR Spectra – Complete Collection**

6049 NIR spectra

Copyright © NICODOM 2012

NICODOM Ltd., Hlavni 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

[ftirspectra@gmail.com](mailto:ftirspectra@gmail.com)

5,7-Diiodo-8-hydroxyquinoline  
Cetirizine Hydrochloride  
Riboflavin  
o-Terphenyl  
Sibutramine Hydrochloride  
Dicyclohexyl phthalate  
Pigment Red 254 based | Irgazin DPP Red BO  
Anthraquinone  
Diethyl phthalate  
Naphthol Yellow  
Dibutylphthalate  
2,13-dimethyl-6,9-dioxa-3,12,18-triazabicyclo[12.3.1]octadeca-1(18),14,16-triene, MnCl  
complex  
Phytonadione, Vitamin K1  
8-Hydroxyquinoline-5-sulfonic acid  
Phenanthrene  
Kebuzone  
Decafluorobenzophenone  
Phthalimide  
Phthalic anhydride  
8-Quinolinesulfonic acid  
Reactive Blue 182 | Cibacron Blue F-R  
Santowhite Powder  
Butyl benzyl phthalate  
Ciprofloxacin  
Trioxocyclam  
Solvent Red 24 | Oilsol Red TAX  
PAN Indicator  
Pigment Orange 69 | Pyranthrone Orange  
Bis(pentamethylcyclopentadienyl)hafnium dichloride  
FD&C Blue 2, Indigo carmine, Indigotine  
Phensuximide  
Sodium glycocholate  
1-Hydroxy-2-naphthoic acid  
1-Cyanonaphthalene  
Mono-L-Leucinium Nitrate  
Cyclene Methylphosphonate Tetrahydrobromide  
1-Nitronaphthalene  
7-Amino-1,3-naphthalenedisulfonic acid  
Dicopper(II) citrate hemipentahydrate  
Tripotassium Citrate  
Triethyl phosphonoacetate  
Fluorene  
2-Methyl-2-thiopseudourea sulfate

---

**IS NIR Spectra – Complete Collection**

6049 NIR spectra

Copyright © NICODOM 2012

NICODOM Ltd., Hlavní 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

[ftirspectra@gmail.com](mailto:ftirspectra@gmail.com)

Carbazole  
L-Lactic acid, Sodium salt  
Sodium Tartrate  
L-Cysteine ethyl ester hydrochloride  
Dimethyl phosphonate  
BOC-p-Nitro-phenylalanine  
4,7-Dichloroquinoline  
2-Chloroethylamine hydrochloride  
2-Amino-5-naphthol-7-sulfonic acid  
DL-Homocystine  
Diethyl ethoxymethylenemalonate  
1,7-Octadecane  
1-Decene  
Isoamyl salicylate  
Cefpodoxim Proxetil  
1,4,7,10-Tetraazacyclododecane-1,4,7-triacetic acid, 10-[2-[(diphosphonomethyl)amino]-2-oxoethyl]-  
N-methylpyrrolidone  
N-Acetyl-DL-Tryptophan  
2-Chlorobenzonitrile  
Ramipril  
Benzenesulfinic acid, sodium salt  
4-Bromobenzyl alcohol  
6-Aminouracil  
Phthalide  
2,4-Dichlorobenzaldehyde  
5-Bromoisatin  
4-Methoxybenzonitrile  
Indole-3-acetic acid  
2,6-Dimethylaniline  
1-(Triisopropylsilyl)pyrrole  
2,6-Dichlorophenol  
Pyrogallol = 1,2,3-Trihydroxybenzene  
Choline Bitartrate  
L(+)-Tartaric acid  
7-tert-Butoxy-2,5-norbornadiene  
Phthalic acid, K salt  
Tetrabutylammonium fluoride trihydrate  
Mannitol  
L(-)-Sorbitol  
D-Tagatose  
Hexamethylbenzene  
Pentachlorophenol  
Inositol

---

**IS NIR Spectra – Complete Collection**

6049 NIR spectra

Copyright © NICODOM 2012

NICODOM Ltd., Hlavní 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

[ftirspectra@gmail.com](mailto:ftirspectra@gmail.com)

Bis(dimethylisopropylsilyl-tetramethyl-cyclopentadienyl) titanium (IV) dichloride  
Xylitol  
2,4,6-Trichlorophenol  
N-Vinylpyrrolidone  
3-Thiophenecarboxylic acid  
2-Furoic acid  
Tetrabutylammonium sulfamate  
2-Acetylthiophene  
2-tert-Butylphenol  
2-methyl-4-(phenylmethylene)oxazol-5(4H)-one  
2,2'-Methylenbis-(4-ethyl-6-tert.butylphenol)  
Ionox 100  
2,6-Di-tert-butyl-(dimethylaminomethyl)phenol  
2-Phenylbenzothiazole  
di-1-Benzotriazolyl Carbonate  
Tetrabenzylcyclam  
2,5-Di-tert-butylhydroquinone  
1,4-Diphenyl-1,3-butadiene  
o-Iodobenzoic acid  
o-Nitrotoluene  
2-Nitroaniline  
2-Nitrophenol  
Picryl chloride  
Picric acid  
10-[[ (2-carboxyethyl)(hydroxy)phosphoryl]methyl]-1,4,7,10-tetraazacyclododecane-1,4,7-triacetic acid  
o-Phthalic acid  
Pyromellitic acid  
1-Phenyl-3-Methyl-5-Pyralozone  
5-Bromosalicylic acid  
5-Aminosalicylic acid  
4-Bromocyclohexanol  
Ethyl(triphenyl)phosphonium chloride  
D-Isoascorbic Acid  
4-Chloro-2-isopropyl-5-methylphenol  
N,N'-bis-4{[4,7,10-tris(carboxymethyl)-1,4,7,10-tetraazacyclododecane-1-yl-methyl]hydroxymethylphenyl} Thiourea  
2-Isopropylidene-5-methylcyclohexanone  
4-Nitrophenol  
Thymol = 3-p-Cymenol  
Resacetophenone  
Basic Red 46 | Anilan Red GRL  
2-Chlorobenzaldehyde  
2-Hydroxybenzyl alcohol = Saligenin

---

**IS NIR Spectra – Complete Collection**

6049 NIR spectra

Copyright © NICODOM 2012

NICODOM Ltd., Hlavni 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

[ftirspectra@gmail.com](mailto:ftirspectra@gmail.com)

Agar powder  
Salicylaldehyde  
Brij 35, Polyoxyethylene (35) Lauryl Ether  
Triton, X-100  
Poly(vinyl methyl ether)  
o-Anisidine  
Dextrin  
Methocel K15  
Nitrocellulose E-150  
PEG 400 Monostearate  
Brij 78  
Guaiacol  
Starch  
Sodium alginate  
Tween 81, POE (5) sorbitan monooleate  
Tween 60, POE (20) Sorbitan Monopalmitate  
Tween 85, POE (20) sorbitan trioleate  
Tween 65, POE (20) sorbitan tristearate  
Polyglyceryl Oleate  
Triphenyltin(IV) acetate  
Zein  
1-Bromonaphthalene  
1-Methylnaphthalene  
Agarose  
Chitosan  
1-Chloronaphthalene  
1-Naphthol  
4-Amino-5-hydroxy-2,7-naphthalenedisulfonic acid  
Pigment Yellow 126 | Permanent Yellow DGR  
2-Phenylbutyric acid  
Dimethylcyclen  
Additin PAN  
4-Methyl-7-hydroxycoumarin  
2-Phenylphenol  
Anthrone  
Bovine albumin  
Sumilizer GA-80  
Melaminium Tartrate Monohydrate  
Mandelic acid  
2,2'-Thiobis(6-tert-butyl-p-cresol)  
Diethylenetriamine-N'-methylenephosphonic-N,N,N'',N''-tetraacetic acid Hydrochloride  
Paraffin wax  
Michler's Ketone  
Aminodiphenylmethane

---

**IS NIR Spectra – Complete Collection**

6049 NIR spectra

Copyright © NICODOM 2012

NICODOM Ltd., Hlavni 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

[ftirspectra@gmail.com](mailto:ftirspectra@gmail.com)

Benzhydrol  
Cholesteryl Chloride  
1,2-Dicyanobenzene  
Veratrole  
Decahydronaphthalene  
Naphthalene  
Noscapine hydrochloride  
4-Aminodiphenylamine-2-sulfonic acid  
Ropinirole hydrochloride  
N-Phenylanthranilic acid  
Isatin  
Amaranth  
2-Methylnaphthalene  
N-Methyl Aminomethylphenylphosphinate  
2-methylquinoline  
Coumarin  
N,N-Diethylaniline  
Meso-tetraphenylporphine  
Antazoline mesylate  
Silquest A 1100  
Pilocarpine  
Dipotassium tartrate  
L-Glucose  
1-Phenyl-3-pyrazolidone  
2,3-Dihydroxynaphthalene  
Methyl Ditosyl Diethylenediamine  
Biphenyl  
4-Bromobiphenyl  
4-Phenylphenol  
3-Hydroxy-2-naphthoic acid  
2,5-Diphenyloxazol  
Phenothiazine  
3-Hexyne  
4,4'-Diamino-diphenyl  
5-Hexyn-1-ol  
4-Acetylbiphenyl  
p-Terphenyl = 1,4-Diphenylbenzene  
2,2'-(Ethylenedioxy)diethylamine  
2-Naphthylamine-6-sulfonic acid  
2,5-Dimethoxybenzaldehyde  
3,4-Dimethoxybenzyl alcohol  
2-Methoxynaphthalene  
N,N-Diethyl-1,4-phenylenediamine  
Reactive Blue 222 based

---

**IS NIR Spectra – Complete Collection**

6049 NIR spectra

Copyright © NICODOM 2012

NICODOM Ltd., Hlavni 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

[ftirspectra@gmail.com](mailto:ftirspectra@gmail.com)

3,4-Dimethoxybenzoic acid  
2'-Acetonaphthone  
2-Naphthoic acid  
Quinaldic acid  
2-Ethyl-4-methylimidazol  
Guaifenesin  
Cyclohexylacetylene  
2-Ethoxynaphthalene  
Cyclooctene  
6-Azathymine  
Umbelliferone  
Caffeine Citrate  
(3,4-Dimethoxyphenyl)acetic acid  
Benzyloxycarbonyl-tosyl-ethanolamine  
Venlafaxine Base  
2-Aminobenzimidazole  
Santowhite CI  
2-Azacyclononanone  
1-Phenyl-1,2-ethanediol  
Methyl benzoate  
N-Methylformanilide  
trans-Perhydrotetraazaphenanthroline  
2-Amino-4-chlorophenyl phenyl ether  
2,4'-Dichloroacetophenone  
Ethyl benzoate  
Benzanilide  
Hydroxypropyl-beta-cyclodextrin  
Propylparaben  
1'-Acetonaphthone  
Butyl 4-aminobenzoate  
Undecylenamidopropyltrimonium Methosulfate  
Dibenzoyl peroxide  
Chalcone  
Benzyl nicotinate  
1,7-Dicarboxy-4,11-bis(Methylenephosphonic acid) Cyclene  
Salicylaldoxime  
Diphenylsilanediol  
2,4-Dichlorophenoxyacetic acid  
Solvent Black 45 | Savinyl Black RLSN  
Copper inhibitor 50  
N,N'-bis(salicylidene)ethylenediamine  
p-Nitro-phenylalanine  
2,4-Dihydroxybenzaldehyde  
Indene

---

**IS NIR Spectra – Complete Collection**

6049 NIR spectra

Copyright © NICODOM 2012

NICODOM Ltd., Hlavní 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

[ftirspectra@gmail.com](mailto:ftirspectra@gmail.com)

Benzotriazole  
N-Cyclohexyl-2-benzothiazolesulfenamide  
Dimethylglyoxime  
1,1'-Bis(diphenylphosphino)ferrocene-palladium(II)dichloride dichloromethane complex  
o-Xylene  
o-Cresol  
1,2-Dichlorobenzene  
Solvent Blue 26 | Alcovar Blue FLH  
2-Chloroaniline  
o-Toluidine  
1,2-Phenylenediamine  
2-Aminophenol  
o-Chlorophenol  
4'-Chlorochalcone  
4-Chlorochalcone  
3,4-Dimethylaniline  
2,6-Dichlorophenol indophenol  
2,4-Dimethylaniline  
4-Chloro-2-methylaniline  
Marfey's reagent  
Eplerenone Hydroxyacid Potassium Salt  
3-Chloro-4-methylaniline  
3,4-Dichloroaniline  
p-Xylidine  
2,4-Diaminotoluene  
Deoxyadenosine  
2,5-Dichloroaniline  
2-Hydroxy-5-methylaniline  
5-Chloro-2-hydroxyaniline  
Diethyl oxalate  
4'-Methoxychalcone  
1,2,4,5-Tetramethylbenzene  
4-Methoxychalcone  
2,4,5-Trichlorophenol  
Chloromaleic anhydride  
Deoxyguanosine  
Calix[6]arene  
1,3-Dichloro-2-propanol  
3-Chloro-1,2-propanediol  
2-Butanone oxime  
Methyl acrylate  
Cyclopentanol  
2,2'-(1,4-phenylene)bis-1,3,2-dioxaborinane  
2-(3-(2h-benzotriazol-2-yl)-4-hydroxyphenyl)ethyl methacrylate

---

**IS NIR Spectra – Complete Collection**

6049 NIR spectra

Copyright © NICODOM 2012

NICODOM Ltd., Hlavní 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

[ftirspectra@gmail.com](mailto:ftirspectra@gmail.com)

g-Butyrolactone  
2-Aminothiazole  
2-Mercaptothiazoline  
3-(2H-Benzotriazol-2-yl)-4-Hydroxy Phenetyl Alcohol  
4,4'-thiobis(2-tert-butyl-5-methylphenol)  
N-Acetyl-L-valine  
1-Chloro-2,4-dinitrobenzene  
2,4-Dinitroaniline  
5-Sulfosalicylic acid  
Bis(5-chloro-2-hydroxyphenyl)methan  
Fluorone Black  
4-Allyl-2-methoxyphenol  
Nuosept 145  
5-Ureidohydantoin  
Itaconic acid  
Igranox 1222  
L-Malic acid  
Bluo-gal  
Butyl methacrylate  
Ethylene glycol dimethacrylate  
2-Chloro-1,1-dimethoxyethane  
Tetrahydrofurfuryl alcohol  
2-Furaldehyde  
Pigment Red 122  
2-Thiophenecarboxaldehyde  
Phenyltrimethylammonium iodide  
Benzenearsonic acid  
Benzenesulfonic acid  
Naugard 431  
4-t-Butylcatechol  
2-Chloroaniline-5-sulfonic acid  
Metroprolol Succinate  
Diisopropyl Bromomethylphosphonate  
3-Nitrobenzenesulfonic acid  
4-t-Butylphenol  
Toluenesulfonyl Chloride  
Solvent Blue 50 | Irisol Fast Blue GLE  
4-Nitrobenzenesulfonyl chloride  
Phenylboronic acid  
Cumene  
Acetophenone  
a,a-Dichlorotoluene  
3-Pyridinecarboxamide  
Rhodamine 6G

---

**IS NIR Spectra – Complete Collection**

6049 NIR spectra

Copyright © NICODOM 2012

NICODOM Ltd., Hlavni 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

[ftirspectra@gmail.com](mailto:ftirspectra@gmail.com)

Nitrobenzene  
Pyrazinecarboxamide  
2-Pyrazinecarboxylic acid  
Picolinic acid  
m-Toluic acid  
3-Hydroxybenzoic acid  
1-Amino-3-nitrobenzene  
Morpholine  
Irganox 565  
(2-hydroxypropyl)-alpha-cyclodextrin  
Irganox 858  
Benzopurpurin 4B  
Bismuth Gallate Basic  
3,5-Dinitrobenzoyl chloride  
3,5-Dinitrobenzoic acid  
1,3,5-Trinitrobenzene  
Ondansetron  
Aminoguanidine sulfate  
m-Dinitrobenzene  
N,N-Diethyl-1,1,1-trimethylsilylamine  
1,1,1-Trifluoro-2-phenyl-3-butyn-2-ol  
2,4-Dibromoacetophenone  
N-Acetyl-L-aspartic acid  
4-Methyl hydroxybenzoate = Methylparaben  
Ethyl-4-nitrobenzoate  
Alpha-terpinene  
4-Isopropyltoluene, p-Cymene  
4'-Bromoacetophenone  
4'-Chloroacetophenone  
Benzyltrimethylcyclam  
p-Aminoacetophenone  
4-Hydroxyacetophenon  
p-Toluic acid  
4-Hydroxybenzoic acid  
N,N-Dimethyl-p-toluidine  
p-Nitrotoluene  
1,1,1,3,3,3-Hexamethyldisilazane  
(+)-Hydroxymethyl tolterodine S-(+)- mandelate  
(Carbonato) bis(ethylenediamine)cobalt(III) Chloride  
[(Dicarboxyamino)methyl]phosphinic acid  
1,8-bis (Methylenephosphonic acid monoethyl ester) Cyclam  
1,8-bis[ethoxy(hydroxy)phosphorylmethyl]-1,4,8,11-tetraazacyclotetradecane  
17-carboxybudesonide  
1-Anthraquinone-N-aminophenylphosphonate

---

**IS NIR Spectra – Complete Collection**

6049 NIR spectra

Copyright © NICODOM 2012

NICODOM Ltd., Hlavní 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

[ftirspectra@gmail.com](mailto:ftirspectra@gmail.com)

1-Methyl-2-dimethylaminomethylpyrrole  
1-Octyl-4-(thien-3-yl)-1,2,3-triazole  
2,13-Dimethyl-6,9-dioxa-3,12,18-triazabicyclo[12.3.1]octadeca-1(18),14,16-triene  
2,2',3,3'-Tetramethoxybenzoin  
2,4,6-Triaminopyrimidinium Malonate Monohydrate  
2,4,6-Triaminopyrimidinium Succinate Dihydrate  
2-Aminopyrimidine-Succinic Acid  
2-Aminopyrimidinium (1+) Hydrogenoxalate Monohydrate  
2-Aminopyrimidinium (1+) Perchlorate 1:2  
2-Aminopyrimidinium (1+) Perchlorate 2:1  
2-ethynyl- $\alpha,\alpha,\alpha$ -trifluorotoluene  
2-Hydroxy-4-ethoxybutylformiate-benzophenon dimer  
2-Hydroxy-4-ethoxyethylformiate-benzophenon dimer  
2-Hydroxy-4-ethoxyethylformiate-benzophenone dimer  
2-Hydroxy-4-ethoxymethylformiate-benzophenone dimer  
2-Methyl-4,6-bis(hydroxy benzyl N,N-dimethyldithiocarbamate)phenol  
2-Methylhydrochlorothiazid  
2-Naphthylamine-3,6-disulfonic acid  
3,5-Dialkyl-2-hydroxybenzyl N,N-dimethyldithiocarbamate  
3-Amino-1,2,4-triazolinium (1+) hydrogen L-tartrate  
3-Propenyltriphenylphosphonium bromide  
4,5-Dimethoxy-1,3,2-benzodioxaphosphole  
4-Amino-1,2,4-triazole-1-ium hydrogen oxalate  
4-Amino-1,2,4-triazole-adipic acid (2:1)  
4-Amino-1,2,4-triazole-succinic acid (1:1)  
4-Fluorobenzenesulfonate Sodium Salt  
4-Nitrobenzylmethylphosphonate Ethyl Ester  
4-tert-Octyl phenylsalicylate  
5,7-Dimethyl 12,14-Diphenyl Cyclam  
Algol red  
Allyl-cis-decahydro-2a,4a,6a,8a-tetraazacyclopent[fg]acenaphthylene  
Aquapentaamminecobalt(III) iodide  
Benzylaminomethylphosphonic Acid Diphenyl Ester Protected  
Betacare  
Biguanidinium (1+) hydrogen succinate  
Biguanidinium (1+) Oxalate Hemihydrate  
Biguanidinium (1+) Phosphite Trihydrate  
Biguanidinium (2+) Carbonate Monohydrate  
Biguanidinium (2+) Nitrate  
Biguanidinium (2+) Phosphite Monohydrate  
Biguanidinium (2+) Sulfate Monohydrate  
Biguanidinium (2+) Tartrate  
Bis(benzyl-tetramethyl-cyclopentadienyl) titanium (IV) dichloride  
Bis(fluorophenyl-tetramethyl-cyclopentadienyl) titanium (IV) dichloride

---

**IS NIR Spectra – Complete Collection**

6049 NIR spectra

Copyright © NICODOM 2012

NICODOM Ltd., Hlavni 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

[ftirspectra@gmail.com](mailto:ftirspectra@gmail.com)

Bis(melaminium) Phosphite Tetrahydrate  
Bis(phenyl-trimethyl-cyclopentadienyl) titanium (IV) dichloride  
Boc-Tyr DCHA  
Carbaethopendecinate Bromide  
Cetyl Cholesteroyl Carbonate  
Copper(II) glycinate  
Cyclam Methylbenzylphosphinate  
Desethoxytamsulosin  
Diacetone-L-sorbose  
Dibenzylaminomethyl Piperidylmethylphosphinic Acid  
Dibenzylaminomethylphosphonic acid ethyl ester  
Dibenzyltrifluoroacetyl cyclam  
Dihydroisoandrosterone  
Diiodofluorescein sodium salt  
Ethoxyphenoldimer  
Ethyl 4-Fluorobenzenesulfinate  
Fingophosphate Racemate  
Flavine Acide Brill.  
Cholesteryl 3-Phenylpropionate  
Isopropyl 4-Fluorobenzenesulfinate  
Isopropylxanthogenate Nickel  
Linear Alpha Olefin C 14  
Linear Alpha Olefin C 18 isomerized  
Linear Alpha Olefin C 24-28  
Malachite green (Sulfonate)  
Melaminium Dihydrogenphosphate  
Melaminium Tartrate 2.5 Hydrate  
Methylpioglitazone  
Mono-beta-Alaninium Nitrate  
Morpholine Methylphenylphosphinate  
Morpholine Methyl-t-butylphosphinate  
N-(4-Hexyloxybenzylidene)-p-Toluidine  
N-(4-Octyloxybenzylidene)-p-Toluidine  
N-(chlorosulfonyl)-2-Pyrrolicarboxamide  
N,N-Dimethyl-4-propyl-1,3,2-benzodioxaphosphol-2-amine  
N-Glycyl-aminomethylphosphonic acid ethyl ester  
N-Isopropyl-N'-(p-isopropylphenyl)-1,4-phenylenediamine  
Nitrobenzylhydroxomethylphosphinic Acid  
N-Leucyl-Aminoethylphosphonic Acid, SR Isomer  
N-Phenylalanyl-aminomethylphosphonic acid  
Octacyanomolybdate(IV) potassium dihydrate  
O-Demethyl Tramadol  
Ostion (HPLC sorbent)  
Palladon

---

**IS NIR Spectra – Complete Collection**

6049 NIR spectra

Copyright © NICODOM 2012

NICODOM Ltd., Hlavni 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

[ftirspectra@gmail.com](mailto:ftirspectra@gmail.com)

p-Azobenzene sulfonic acid  
Pigment Orange 59 | Paliotol Orange  
Piperazine N,N-bis(Aminomethylphosphonic Acid) Diethyl Ester  
Piperidine Methylhydroxomethylphosphinate  
Piperidine Methylphenylphosphinate  
p-Lobelanine  
Polynol A, mono and diglycerides of fatty acid  
Pyrethoxyioglitazone  
Reactive Red 238 | Cibacron Red C-R  
Salicylaldehyde semicarbazone  
Silquest A-174 Silane  
Simvacid Ammonium  
Sodium Alkyl Aryl Sulfonate  
TAA OL-methacrylate  
t-Butyl Gamma Cyclodextrin  
Tetraformylcyclen  
Tinopal 2B Whitener  
Topanol CA  
Tridecanamine Undecylenate  
Triphenylantimony(V) hydroxide  
Triphenylantimony(V) chloride  
True Red E, Echt Rot E  
Venlafaxine Besylate Monohydrate  
Venlaphenol  
Zinc hexacyanoferrate(II) trihydrate  
Potassium Ferricyanide  
Phosphate Ammonium Dibasic  
Persulfate Sodium  
Iron (III) Oxide, Hematite  
Amitriptyline  
Sulfate Copper Pentahydrate  
Benzocaine  
Ibuprofen  
Nitrofurantoin  
Ornidazole  
Oxiconazole Nitrate  
Papaverine hydrochloride  
Phenobarbital Sodium  
Phenylephrine hydrochloride  
Phthalylsulfathiazol  
Pilocarpine Hydrochloride  
Piroxicam  
Prazosin Hydrochloride  
Propafenone

---

**IS NIR Spectra – Complete Collection**

6049 NIR spectra

Copyright © NICODOM 2012

NICODOM Ltd., Hlavni 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

[ftirspectra@gmail.com](mailto:ftirspectra@gmail.com)

Propyphenazone  
Ranitidine Hydrochloride  
Reserpine  
Roxithromycine  
Rutoside  
Silver Albumin Acetyl Tannate  
Silver Albumin Acetyl Tannate with Borax  
Sulpiride  
Tetracycline Hydrochloride  
Thioridazine  
Triamcinolone Acetonide  
Trimecain Chloride  
Trimethoprim  
Trolox  
Ursodeoxycholic Acid  
Viprynum Embonate  
Xanthinol Nicotinate  
Glucono-delta-lactone  
Polyvinylidene fluoride  
Potassium tetracyanocadmate  
Ammonium thiocyanate  
Nickel (II) cyanide  
Dibutyltin oxide  
Polyalkylbenzene  
Cobaltous acetate  
Cobalt (II) Oxalate  
Nitrate Cerium (III) Ammonium  
2-(2,5-Dibromophenyl)-5-pyrid-2-yl-1,3,4-Oxadiazole  
Magnesium Perchlorate  
Nickel Ammonium Sulfate Hexahydrate  
Sulfate Cerium (III) hydrate  
Cadmium Acetate  
Nitrate Cadmium tetrahydrate  
Lead Acetate  
Manganese Perchlorate  
Sulfate Cadmium  
Nitrate Copper  
Sodium m-arsenite  
Cesium selenate  
Lithium Selenate  
Nitrate Thallium  
Barium selenite  
Barium perchlorate  
Cadmium selenate

---

**IS NIR Spectra – Complete Collection**

6049 NIR spectra

Copyright © NICODOM 2012

NICODOM Ltd., Hlavni 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

[ftirspectra@gmail.com](mailto:ftirspectra@gmail.com)

Copper(II) Acetate  
Magnesium selenate hexahydrate  
Nitrate Magnesium  
Ammonium hydrogenselenite  
Zinc selenate  
Nitrate Mercurous Dihydrate  
Nitrate Cadmium  
Copper(ii) selenite dihydrate  
Chromium (III) hydroxide n-hydrate  
Manganese selenate  
Cobalt Selenite  
Carbonate Cobalt (II) basic  
Sulfate Nickel(II) Heptahydrate  
Nickel Selenate  
Lead(II) Hydrogenselenite  
Telluric acid  
Bismuth Phosphate  
Sodium stannate trihydrate  
Phosphorous acid  
Periodic Acid  
Cerium(III) Oxalate Nonahydrate  
Nitrate Uranyl(VI)  
Nitrate Lanthanum(III)  
Nitrate Cerium (III) hydrate  
Nitrate Samarium  
Nitrate Gadolinium hexahydrate  
Nitrate Erbium(III) pentahydrate  
Nitrate Neodymium Hexahydrate  
Nitrate Yttrium hexahydrate  
Nitrate Lutetium hydrate  
Nitrate Praseodymium(III) hexahydrate  
Nitrate Holmium(III) pentahydrate  
Nitrate Thulium (III)  
Nitrate Ytterbium pentahydrate  
Nitrate Dysprosium  
Nitrate Terbium(III) pentahydrate  
Nitrate Europium(III) hydrate  
Carbonate Cadmium  
Perchlorate Cadmium Hexahydrate  
Perchlorate Zinc Hexahydrate  
Perchlorate barium trihydrate  
Magnesium Perchlorate Hexahydrate  
Nitrate Strontium  
Lithium Perchlorate Trihydrate

---

**IS NIR Spectra – Complete Collection**

6049 NIR spectra

Copyright © NICODOM 2012

NICODOM Ltd., Hlavni 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

[ftirspectra@gmail.com](mailto:ftirspectra@gmail.com)

Sodium Perchlorate Monohydrate  
Silver acetate  
Copper Trifluoromethanesulfonate  
Copper Perchlorate Hexahydrate  
Sulfate Silver  
Sulfate Beryllium  
Sodium selenite  
Tellurium (IV) Oxide  
Ammonium chromate  
Ammonium dichromate  
Potassium dichromate  
Manganese Perchlorate Hexahydrate  
Periodic Acid Dihydrate  
Iodic acid  
Cobalt (II) Perchlorate Hexahydrate  
Sulfate Cobalt heptahydrate  
Nitrate Cobalt Hexahydrate  
Tris-ethylenediamine-cobalt(III) Bromide Dihydrate  
Cobalt(III) acetylacetonate  
Nickel(II) Hydroxide Carbonate  
Nickel Perchlorate Hexahydrate  
Nickel Acetate Tetrahydrate  
Bis(triphenylphosphine)nickel(II)chloride  
Nitrate Mercuric Monohydrate  
Sulfate Mercurous(I)  
Mercuric (II) Thiocyanate  
Sulfate Mercuric(II)  
Sulfate Thallium  
Nitrate Gallium  
Thallium(i) 2,4-pentanedionate  
Carbonate Thallium  
Carbonate Lead (II)  
Nitrate Zirconium (IV)  
Nitrate Zirconyl  
Sulfate Zirconium tetrahydrate  
Hafnium Oxide  
Sulfate Europium  
Ytterbium Oxalate  
Nitrate Cerium(III) Hexahydrate  
Lanthanum (III) Amidosulfate  
Cerium Sulfate Tetrahydrate  
Molybdenum hexacarbonyl  
Rhenium pentacarbonyl bromide  
Dirhenium decacarbonyl

---

**IS NIR Spectra – Complete Collection**

6049 NIR spectra

Copyright © NICODOM 2012

NICODOM Ltd., Hlavni 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

[ftirspectra@gmail.com](mailto:ftirspectra@gmail.com)

Tungsten hexacarbonyl  
Sulfate Cerium  
Lanthanum Oxide  
Nitrate Lanthanum hexahydrate  
Nitrate Yttrium pentahydrate  
Nitrate Uranyl(VI) Hexahydrate  
Uranyl acetate dihydrate  
Sodium cyanate  
Samarium Oxide  
Tetraphenylarsonium chloride  
(Methoxymethyl)triphenylphosphonium chloride  
Triruthenium dodecacarbonyl  
Diiron nonacarbonyl  
Mixture of alkylated phenols  
Mixture of cumylated phenols  
Nitrate Chromic  
Chromium (III) potassium trisoxalate  
Cadmium tosylate  
Chromic acetate  
Sodium Peroxide  
Sulfate Chromic  
Potassium chromate  
Chromium (III) potassium sulfate dodecahydrate  
Boron Trifluoride, etherate  
Tetrafluoroboric acid in HF  
Tetramethyl ammonium chloride  
Phosphinic acid  
Sulfuryl chloride  
Chromium hexacarbonyl  
Pentaamminecobalt(III) carbonate  
Pentaamminecobalt(III) Iodide Carbonate  
Pentaamminecobalt(III) Carbonate Nitrate  
Pentaammineiodidecobalt(III) iodide  
Pentaamminecobalt(III) carbonate iodide  
Potassium Cobalt(II) Sulfate  
Lanthanum Diphosphate  
Cyanuric chloride  
2,4-Di-T,-butylphenol based additive  
Lead Sulfate Tribasic Modified  
Nitrate Mercurous  
Beryllium acetate, basic  
Uranyl Acetate  
Sulfate Uranyl  
Zinc bis(acetato-O)dioxouranate

---

**IS NIR Spectra – Complete Collection**

6049 NIR spectra

Copyright © NICODOM 2012

NICODOM Ltd., Hlavni 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

[ftirspectra@gmail.com](mailto:ftirspectra@gmail.com)

Ammonium Diuranate  
Sodium Metaarsenate  
Sodium Arsenate Hydrate  
Sulfate Beryllium Tetrahydrate  
Sodium borohydride  
Hydrazine monohydrochloride  
Hydrazine sulfate  
Pentaamminechlorocobalt(III) chloride  
Cobalt(II) selenite dihydrate  
Magnesium Hydrogen Selenite Tetrahydrate  
Sulfate Nickel Hexahydrate  
Copper(II) Diammonium Chloride Dihydrate  
Diantimony Trioxide  
Chromic (VI) Trioxide  
Perrhenic Acid  
Sulfate Nickel  
Carbonate Nickel(II) (basic) hydrate  
Nitrate Platinum (II)  
Carbonate Nickel  
Potassium Dicyanoargentate  
Cuprous (I) Cyanide  
Sulfate Stannous  
Selenious acid  
Cadmium(II) Hydroxide with Carbonate  
Carbonate Nickel(II) Hydroxide Tetrahydrate  
Nickel(II) Hydroxide  
Molybdenum(V) Tetrachloride Oxide  
Lithium bromide  
Phosphonitrilic chloride trimer  
Bromocyan  
Lead(II) nitrate  
Nitrate Zirconyl Dihydrate  
Trinitrotoluene (TNT)  
Semtex 1A  
Octogen (HMX)  
Hexogen (RDX)  
Pentaerythritol tetranitrate (Pentrite)  
DMDNB  
Tetryl  
Dynamite  
Nitrocellulose #1 for explosives  
Pentritol Np 10T  
Bees Wax for explosives  
Paraffin liquid for explosives

---

**IS NIR Spectra – Complete Collection**

6049 NIR spectra

Copyright © NICODOM 2012

NICODOM Ltd., Hlavní 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

[ftirspectra@gmail.com](mailto:ftirspectra@gmail.com)

Permonex V19  
Permon 10  
Polonit V  
Permon G3  
Ostravit C  
Perunit HG  
A-IX-1  
Danubit 3  
Perunit 22  
Perunit 20  
Danubit 1  
Emsit 20  
Perunit 28  
Slavit V  
Permon DAP P  
DM-1  
Semtex S30  
DAP-1  
Np10  
Semtex 10SE  
Semtex 1H  
T/H/E  
DAP-2  
Harmonit AD  
A-IX-2  
Synthesit V18  
Obrysit  
Permon DAP M  
Emsit I  
PI U-EP 14  
Permon Extra 11  
Semtex 2  
Danubit Geofex 2  
Carbodanubit  
2,5-Dinitrophenol  
2-Amino-4-nitrotoluene  
4-Amino-3-nitrotoluene  
2-Amino-6-nitrotoluene  
Perchlorate Sodium  
Perchlorate Ammonium  
Nitrate Barium  
Di-isononyl Phthalate  
Acrobat MZ  
Actara 25WG

---

**IS NIR Spectra – Complete Collection**

6049 NIR spectra

Copyright © NICODOM 2012

NICODOM Ltd., Hlavni 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

[ftirspectra@gmail.com](mailto:ftirspectra@gmail.com)

Afalon 45SC  
Alette 80 WP  
Alister  
Ally 20 SX  
Applaud 25 WP  
Atlantis WG  
Atonik  
Attribut S670  
Aurora 40 WG  
Axial 100 EC  
Bancol 50 WP  
Banvel 480S  
Basagran Forte  
Basta 15  
Baycor 25 WP  
Bentazone Na 600  
Betanal Expert  
Bladex  
Bladex 50 SC  
Bofix  
Boxer 800 EC  
Breakthru  
Breakthru S240  
Butisan 400 SC  
Butisan Star  
Calaris 400 SC  
Callisto 480 SC  
Calypso 480 SC  
Casoar  
Chisel 75 WG  
Command 4EC  
Delan 700 WDG  
Dicopur  
Discus  
Dithane DG  
Dithane M 45  
Dominator  
Dual Gold 960 EC  
Dursban 10 G  
Equip  
Escort  
Esteron 60  
Euparen Multi  
Falcon 460 EC

ir-spectra.com IS NIR Spectra

---

**IS NIR Spectra – Complete Collection**

6049 NIR spectra

Copyright © NICODOM 2012

NICODOM Ltd., Hlavni 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

[ftirspectra@gmail.com](mailto:ftirspectra@gmail.com)

Fertimag Mgn 8-8  
Folpan 80WG  
Frontier 900 EC  
Fusilade Super  
Gallant Super  
Gardoprim Plus Gold 500 CS  
Genamin T200 BM  
Gesagard 80WP  
Gesaprim 500 FW  
Gesaprim 90 WDG  
Glean 75 PX  
Glean 75 WG  
Glyphogan 480 SL  
Goal 2E  
Granstar 50 SX  
Granstar 75WG  
Grid  
Grodyl Plus  
Guardian Safe Max  
Harmony Extra  
Hattrick  
Husar  
Husar Liquid  
Juwel  
Karathane LC  
Kompal 365 SE  
Lentagran WP  
Lintur 70 WG  
Lintur 75 WG  
Logran 75 WG  
Maister  
Maister Liquid  
Marshal 250 EC  
Merlin 480 SC  
Merlin 750 WG  
Merpan 80 WG  
MesuroI  
Mikado SC 300  
Milagro 040 CS  
Monitor 75 WG  
Mospilan 20 SP  
Mythos 30 SC  
Nissorun 10 WP  
Novozir Mn 80

ir-spectra.com IS NIR Spectra

---

**IS NIR Spectra – Complete Collection**

6049 NIR spectra

Copyright © NICODOM 2012

NICODOM Ltd., Hlavni 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

[ftirspectra@gmail.com](mailto:ftirspectra@gmail.com)

Nurelle D  
Omite 30 W  
Omite 570 EW  
Pantera 40 EC  
Pinoxaden 100 EC  
Pirimor 50 WG  
Pirimor DP  
Polyram WG  
Previcur 607 SL  
Puma Extra  
Racer 25 EC  
Refine  
Reldan 40 EC  
Rider 400 EC  
Ridomil Gold Mz 68 WP  
Roundup Bioactive  
Roundup Classic  
Rovral Flo  
Rubigan 12 EC  
Safari 50 DF  
Saprol  
Sekator  
Sencor 70 WG  
Solar  
Sporgon 50 WP  
Stabilan 750 SL  
Stemat 200 EC  
Stomp 400 SC  
Sumithion Super  
Syllit 65 WP  
Synbetan P  
Talstar 10 EC  
Tango Super  
Tatoo  
Teldor 500 SC  
Teridox 500 EC  
Tiara 500 SC  
Tilt 250 EC  
Titus 25 WG  
Tomigan 250 EC  
Topik 080 EC  
Touchdown  
Trophy  
U-46-D-Fluid

ir-spectra.com IS NIR Spectra

---

**IS NIR Spectra – Complete Collection**

6049 NIR spectra

Copyright © NICODOM 2012

NICODOM Ltd., Hlavni 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

[ftirspectra@gmail.com](mailto:ftirspectra@gmail.com)

U-46-M-Fluid  
Zato 50 WG  
Zolone 35 EC  
Silver Nitrate  
Ammonium Chloride  
Phosphate Calcium Monobasic  
Sulfate Sodium Monobasic monohydrate  
Dithionite Sodium  
Stannous Chloride  
Tricalcium Orthophosphate  
Food Brown 3 | Chocolate brown HT  
Potassium nitrite  
Phosphate Ammonium Monobasic  
Potassium Nitrate  
Hexadecanoic Acid, Palmitic acid  
Raloxifene-4'-glucuronide  
Doxifluridine  
Bisoprolol Fumarate  
Budesonide  
Capecitabine  
Eplerenone  
Fingolimod hydrochloride  
Hydroxy Pioglitazone (M-IV)  
Imatinib  
Levocetirizine dihydrochloride  
1-(2-Hydroxyethyl)-2-Hydroxymethyl-5-Nitroimidazole  
Alendronate Sodium  
Cialis  
Paliperidone  
Dehydro Aripiprazole  
Atomoxetine hydrochloride  
Dihydroxydecarbamapecitabine  
Donepezil Hydrochloride  
Escitalopram oxalate  
Exemestane  
Gliclazide  
5-[[4-[2-(5-Acetyl-2-Pyridinyl)Ethoxy]Phenyl]Methyl]-2,4-Thiazolidinedione  
Olanzapine Benzoate  
Phenolthiadione  
Pramipexole  
Pyrhydroxy pioglitazone  
Pyrketopioglitazone  
Rivastigmine Tartrate  
Zopiclone

---

**IS NIR Spectra – Complete Collection**

6049 NIR spectra

Copyright © NICODOM 2012

NICODOM Ltd., Hlavní 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

[ftirspectra@gmail.com](mailto:ftirspectra@gmail.com)

Acamprosate Calcium  
Amlodipin Besylate  
Amlodipine Maleate  
Carvedilol  
Gliquidone  
Lovacid Ammonium Salt  
Moclobemide  
Paroxdione  
Paroxetine  
Paroxetine Mesylate  
Ramiprilat  
Ribavirin  
Simvastatin  
Tolterodine Tartrate  
Trazodone hydrochloride  
Venlafaxine Hydrochloride  
Hexachloroplatinic (IV) acid hexahydrate  
Atropine  
Fenoterol Hydrobromide  
Ketoprofen  
Metoclopramide Hydrochloride  
Nifedipine  
Piracetam  
Tiapride Hydrochloride  
Tolazoline Hydrochloride  
Trihexyphenidyl Hydrochloride  
3,4-Dimethylbenzoic Acid  
Flurbiprofen  
Ethosuximide  
Citalopram Hydrobromide  
Fenoprofen  
Irbesartan  
Hydroxyzine Dihydrochloride  
Pentoxifylline  
Naproxen  
8-Chlorotheophylline  
Glipizide  
(R)-(+)-Desmethylsibutramine HCl  
Lansoprazole  
Donepezil HCl  
Levofloxacin  
Clavulanate Potassium  
Sulfadimethoxine  
2-(2,4,5-Trichlorophenoxy)Propionic Acid

---

**IS NIR Spectra – Complete Collection**

6049 NIR spectra

Copyright © NICODOM 2012

NICODOM Ltd., Hlavní 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

[ftirspectra@gmail.com](mailto:ftirspectra@gmail.com)

Ticlopidine Hydrochloride  
Heparine 500000 Units  
Bicalutamide  
Itraconazole  
Fluconazole  
Lamotrigine  
Hesperetin  
3alpha-hydroxytibolone  
Zolpidem Tartrate  
Risedronate Sodium  
Tenoxicam  
Sulfamethizole  
Quetiapine  
Sumatriptan Succinate  
Moxifloxacin HCl  
Ketoconazole  
Phenytoin (5,5-Diphenylhydantoin)  
Tranexamic Acid  
Ammonium Bromide  
Phenol-formaldehyde resin  
Cellulose Acetate (39.8% acetyl)  
Aluminium Hydroxide  
Polyvinylbutyral  
Mowiol® 4/98, PVA  
Pigment Green 14 | Cadmium Green Dark  
Pigment Yellow 37:1 | Cadmium Yellow no. 6 medium  
Sulfate Aluminum Ammonium Dodecahydrate  
Quillaia extract  
Thiosulfate Sodium  
Iron (III) Hydroxide  
Phosphate Potassium Dibasic  
Carnauba wax  
Sulfate Sodium Monobasic  
Sodium Silicate  
Persulfate Ammonium  
Persulfate Potassium  
Sulfate Ammonium  
Metabisulfite Potassium  
Sodium Metasilicate  
Phosphate Sodium Tribasic Dodecahydrate  
Diphosphate Tetrasodium  
Phosphate Potassium Monobasic  
Sulfate Aluminium  
Sulfate Ferrous

---

**IS NIR Spectra – Complete Collection**

6049 NIR spectra

Copyright © NICODOM 2012

NICODOM Ltd., Hlavni 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

[ftirspectra@gmail.com](mailto:ftirspectra@gmail.com)

Sulfate Magnesium Heptahydrate  
Magnesium Chloride  
Sulfate Calcium Dihydrate  
Calcium Hydroxide  
Calcium Oxide  
Aluminosilicate  
Iron (III) Oxide Hydroxide  
Sulfate Potassium  
Phosphate Potassium Tribasic n-hydrate  
Poly(propylene glycol) 200000  
Beta-cyclodextrine  
Pentasodium triphosphate  
Hydroxy propyl distarch phosphate  
Magnesium bromide diethyl etherate  
Food Red 8 | Ponceau 6R  
Bicarbonate of soda  
Ammonium carbonate  
Magnesium carbonate  
Potassium Ferrocyanide  
Sodium cyclamate  
Thaumatococcus  
Caramel  
Phenoxy resin  
Phosphate Sodium Dibasic Dihydrate  
Titanium Dioxide  
Potassium carbonate, sesquihydrate  
Clove oil  
Silica Gel, 70-230 mesh  
Magnesium Carbonate Basic Hydrate  
Iron (III) Oxide monohydrate, yellow  
Mono and Diglycerides of fatty acids  
Phosphate Calcium Dibasic, Monetite  
Quinoline Yellow  
Potassium chloride  
Neohesperidine dihydrochalcone  
Oil, Cinnamon Oil  
Paragum 501 (Polyacrylate)  
Food Black 1 | Brilliant Black BN  
Food Green 4, Acid Green 50 | Green S  
Sulfate Sodium Anhydrous  
Solvent Yellow 94 | Fluorescein  
Quillaja Saponin  
Sulfate Zinc Heptahydrate  
Acesulfame K

---

**IS NIR Spectra – Complete Collection**

6049 NIR spectra

Copyright © NICODOM 2012

NICODOM Ltd., Hlavni 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

[ftirspectra@gmail.com](mailto:ftirspectra@gmail.com)

Orthoboric Acid  
Sulfate Magnesium  
Oil of Peppermint, Menthol  
Sodium Perborate\*4H2O  
Polyethylene Oxide 4000000  
Gelatine  
Polyoxymethylene  
Acacia gum  
Xanthan gum  
Carrageenan  
Starch sodium octenyl succinate  
2,4,6-tris (dimethylaminomethyl) phenol  
Soya Bean Oil  
Castor Oil  
Amylopectine  
Cellulose Microcrystalline (Avicel)  
Talc  
Carboxymethylcellulose (Carmellose)  
Bentonite  
Phosphate Calcium Tribasic  
Gum Tragacanth  
Kaolin, Aluminum silicate dihydrate  
Lanolin Anhydrous  
Phosphate Calcium Dibasic  
Carbonate Zinc Basic  
Casein sodium salt  
Pectine  
Dioctadecyl 3,3'-thiodipropionate  
Uvinul 4050 H  
Tinuvin 622  
Chimasorb 944  
Silicate - pigment  
Hostavin N30 (HALS)  
Chimasorb 81, 2-hydroxy-4-n-octyloxybenzophenone (UV-531)  
Cyasorb UV 3346  
Glycerides of natural fatty acids  
Escorene ultra UL 00723, Copolymer EVA type  
Evatane 1005 VN 4, Copolymer EVA type  
Lotryl 35MA08  
Polypropylene+Vistalon 719, 1:1  
Polypropylene with rubber content  
Polypropylene, isotactic #1  
Poly(propylene:ethylene), Novolen 2900NC  
Polypropylene-graft-maleic anhydride

---

**IS NIR Spectra – Complete Collection**

6049 NIR spectra

Copyright © NICODOM 2012

NICODOM Ltd., Hlavni 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

[ftirspectra@gmail.com](mailto:ftirspectra@gmail.com)

Poly(ethylene:propylene), Daplen KSR 4525  
Ethylene/Vinylacetate copolymer with Aluminium Oxide pigment  
Poly(N-vinyl carbazole), Luvican  
Dutral TER 9046, Copolymer EPDM type  
Poly(ethylene:propylene), Keltan 740  
Buna CB HX 529 SIC, Polymer PB type  
Poly(butadiene:styrene), Solprene 1322  
Poly(butadiene:styrene), Cariflex 1712  
Irganox 1010 + Irganox 168  
Hostanox 03  
Bis (2-hydroxyethyl) tallow amine  
Glycerides of fatty acids C12-C22  
Millad 3988  
Leucopur EGM  
Loxamid E, Stearoamides mixture  
Polyvinylchloride  
Nylon 6 #2  
Chlorinated paraffin  
Poly(dimethylsiloxane) #2  
Methyl methacrylate  
Diurethane Dimethacrylate  
Triethylene glycol dimethacrylate  
1,4-butanediol dimethacrylate  
Cyclac G 101, Copolymer ABS type  
Poly(acrylonitrile:ethylene:propylene:styrene) #1  
Polystyrene high impact #1  
Polystyrene general purpose, Krasten 127/9001  
HIPS/PP/PA6 alloy  
HIPS/HDPE alloy #2  
Poly[(2,3-epoxypropyl methacrylate):(ethylene dimethacrylate)]  
Poly(ethylene dimethacrylate)  
Poly[(2,3-epoxypropyl methacrylate):styrene:(ethylene dimethacrylate)]  
Reilex 425, Poly[(divinyl benzene):(vinyl pyridine)]  
Poly[(divinyl benzene):styrene]  
Poly[styrene-co-acrylonitrile-co-(divinyl benzene)]  
Poly(acrylonitrile:divinylbenzene)  
Cellulose Modified  
Poly[(acrylic acid):(2-hydroxyethyl methacrylate)]  
Phenol resin #1  
Phenolic Resin  
Melamine-urea-formaldehyde resin #1  
Melamine-urea-formaldehyde resin #2  
Poly(methyl methacrylate)  
Poly(styrene:n-butyl methacrylate)

---

**IS NIR Spectra – Complete Collection**

6049 NIR spectra

Copyright © NICODOM 2012

NICODOM Ltd., Hlavni 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

[ftirspectra@gmail.com](mailto:ftirspectra@gmail.com)

Poly(methyl methacrylate):(butyl methacrylate)  
Kibisan 127, Copolymer SAN type  
Poly(acrylonitrile:styrene), PD Alfa 468  
Poly(butadiene::styrene)  
Poly(butadiene) (55% cis & trans, 45% vinyl)  
Polystyrene, dicarboxy terminated  
Poly(styrene::butadiene) graft  
Poly(styrene::butadiene) ABA block (28% styrene), Kraton D 1101  
Poly(butadiene), phenyl terminated  
Polyether polyol  
Color masterbatch Polypropylene + yellow pigment #1  
Teflon  
Poly(vinyl alcohol)  
Poly(vinyl isobutyl ether)  
Miravithen D23EA, Copolymer EVA type  
Levapren 8338, Copolymer EVA type  
Poly(acrylonitrile:butadiene)  
Poly(styrene::isoprene) ABA block (14% styrene)  
Poly(alpha-methylstyrene-block-butadiene)  
Poly(methyl methacrylate) m.w. 30000  
Poly(acrylamide) #1  
Polyamide 6 #2  
Nylon 11  
Nylon 610  
Ultramid 6A, Poly[caprolactam:(hexamethylene adipamide)]  
Poly(ethylene terephthalate) #1  
Bisphenol-A-polycarbonate, 20% glass reinforced  
Poly(ethylene oxide)  
Phenoxy resin  
Hydroxybutyl Methyl Cellulose  
Cellulose Acetate  
Lignine  
Casein  
Poly(vinylidene fluoride) #1  
Poly(ethylene:tetrafluoroethylene)  
Poly(ethylene:propylene:diene) (50% ethylene, 4% diene)  
Poly(ethylene:propylene:diene) (50% ethylene, 8% diene)  
Color masterbatch Polypropylene + brown pigment  
Color masterbatch Polypropylene + white pigment  
Poly(cyclohexylene terephthalate)  
Nylon 612  
Ultramid C-35  
Nylon 12 #1  
Nylon 66

---

**IS NIR Spectra – Complete Collection**

6049 NIR spectra

Copyright © NICODOM 2012

NICODOM Ltd., Hlavni 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

[ftirspectra@gmail.com](mailto:ftirspectra@gmail.com)

Poly(ethylene terephthalate)  
Bisphenol-A-polycarbonate  
Sodium hydrogensulfite hydrate  
Lead(II) oxide  
Ammonium hydrogenfluoride  
Beryllium oxalate  
Copper(II) oxide  
Ferric sulfate nonahydrate  
Magnesium thiosulfate hexahydrate  
Cadmium chloride hemipentahydrate  
Poly(diallylphthalate) #1  
Aluminium chloride hexahydrate  
Ammonium magnesium chloride  
Barium oxalate  
Barium chloride dihydrate  
Nickel(II) chloride hexahydrate  
Cobalt(II) chloride hexahydrate  
Potassium perchlorate  
Lead(II) phosphate  
Barium borate  
Sulfate Aluminium Potassium Dodecahydrate  
Tetrahydro-1,4-oxazine  
Sodium tetrathioantimonate  
Potassium hydroxide  
Sodium hydroxide  
Lead(II) thiocyanate  
Barium thiocyanate  
Polyethylene Oxide 900000  
Di(ethylene glycol) diethyl ether  
Potassium hydrogentartrate  
Cyclopropylamine  
Sodium nitroferricyanide(III) dihydrate  
Sodium borohydride  
Pebax  
Vespel  
Ultem 2300  
PEEK  
Torlon 4203  
Fluorosilicone Rubber  
Polysulfone Resin  
Turcite X  
Basic Yellow 62 | Crypur Yellow 5GL  
Basic Orange 50 | Crypur Orange GL  
Crypur Blue

---

**IS NIR Spectra – Complete Collection**

6049 NIR spectra

Copyright © NICODOM 2012

NICODOM Ltd., Hlavni 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

[ftirspectra@gmail.com](mailto:ftirspectra@gmail.com)

Basic Blue 102 | Crypur Blue GL  
Blerit GRL  
Cationic Dye | Anilan Yellow 2RL  
Cationic Dye | Anilan Orange GRL  
Cationic Dye | Anilan Red GTL  
Cationic Dye | Anilan Blue BRL  
Cationic Dye | Anilan Blue GL  
Cationic Dye | Anilan Brown GRL  
Cationic Dye | Anilan Black GRN  
Direct Red 76 | Saturn Scarlet LGG 200  
Direct Brown 115 | Saturn Brown L2G  
Direct Orange 26  
Direct Green 97  
Direct Black 80  
Direct Orange 37  
Direct Red 75  
Direct Brown 112  
Gold RL Aluminium  
Blue LR Aluminium  
Mediosol Navy Blue R  
Mediosol Black AB  
Mediosol Blue L  
Mediosol Black G  
Solvent Yellow 16 | Sudan Yellow 3G  
Sudan Brown BBN  
Korostan Khaki G  
Acid Brown 235 | Korostan Brown DGR  
Acid Brown 214 | Korostan Brown RT  
Velour Black B  
Velour Black BD  
Velour Black AN  
Basic Red 37  
Acid Violet 48 based | Triacid Violet 4BN  
Acid Red 73 | Acid Brilliant Croceine MOO  
Acid Green 25 | Alizarine Brilliant Green G  
Alizarine Blue  
Acid Red 1 based  
Acid Yellow 252  
Acid Red 257  
Acid Violet 92  
Acid Green 108:1  
Cationic Dye | Anilan Navy Blue RLN  
Acid Navy Blue GGR  
Acid Brown G

---

**IS NIR Spectra – Complete Collection**

6049 NIR spectra

Copyright © NICODOM 2012

NICODOM Ltd., Hlavni 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

[ftirspectra@gmail.com](mailto:ftirspectra@gmail.com)

Acid Brown R  
Acid Blue 352 Boruta  
Acid Red 433  
Direct Catechin GS  
Acid Orange 142 based | Rylan Orange R  
Phthalostane Blue  
Direct Black 9 | Direct Black FF  
Acid Yellow 40 | Rybanyl Yellow NL 5G  
Acid Blue 40 based | Rybanyl Blue EL 2G  
Rybanyl Black GR  
Acid Red 432 | Rylan Red 2B  
Direct Blue 162 | Rybantin Blue RLL  
Direct Blue 120 | Azogen Blue BR  
Vat Red 13 | Ostanthren Rubine R  
Grelit A  
Grelit E  
Pigment Yellow 37 | Multiplas Yellow MP 1036  
Pigment White 6 based | Multiplas White MP 1009  
Solvent Red 8 | Alcovar Red BEC  
Solvent Yellow | Alcovar Yellow RN  
Solvent Yellow 38 | Alcovar Yellow 3G  
Solvent Yellow 83 | Alcovar Yellow PW  
Solvent Orange 19 | Alcovar Orange APD  
Solvent Red | Alcovar Fire Red GN  
Solvent Red 37 | Alcovar Red BN  
Solvent Blue 25 | Alcovar Blue HD Extra  
Solvent Yellow 12 | Oilsol Yellow HA  
Solvent Orange 2 | Oilsol Orange TX  
Solvent Green | Oilsol Green CS  
Food Yellow 3 based | Water Orange 178759  
Water Blue 180694  
Maxilon Flavine BY-4X0  
Solvent Yellow 57 | Lampronol Yellow J  
Solvent Yellow 21 based | Lampronol Yellow 2RN  
Solvent Red 69 | Lampronol Scarlet R  
Solvent Orange 62 | Lampronol Orange R  
Solvent Yellow 30 | Waxoline Yellow A  
Solvent Red 20 | Waxoline Red A  
Solvent Violet 8 | Waxoline Violet R  
Solvent Blue 36 | Waxoline Blue A  
Solvent Green 3 | Waxoline Green G  
Mixture | Waxoline Black 5BP  
Durazol Blue 10GP  
Coomassie Blue SE

---

**IS NIR Spectra – Complete Collection**

6049 NIR spectra

Copyright © NICODOM 2012

NICODOM Ltd., Hlavni 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

[ftirspectra@gmail.com](mailto:ftirspectra@gmail.com)

Acid Black 24 | Coomasie Black B  
Solvent Blue 38 based | Methasol Blue 2B  
Basic Orange 2 based | Chrysoidine Y  
Acid Orange 10 | Naphthalene Fast Orange 2GS  
Diazol Black VP  
Methyl Violet R  
Methic Black GPN  
Basic Orange 1:1  
Acid Yellow 73 based | Fluoresceine LT  
Carbolan Blue 2G  
Carbolan Blue GN  
Acid Blue 139 | Carbolan Blue 2R  
Lissamine Red 4B  
Acid Blue 45 | Lissamine Blue B  
Lissamine Blue 2BR  
Lissamine Blue A  
Lissamine Blue AR  
Acid Blue 40 based | Lissamine Blue 2G  
Lissamine Blue PFN  
Food Green 2 | Lissamine Green SF  
Acid Black 1 based | Lissamine Black 12B  
Lissamine Yellow FF  
Solvent Blue 38 based | Luxol Fast Blue MBSN  
Solvent Orange 25 | Luxol Fast Orange GRL  
Solvent Brown 20 | Luxol Fast Brown K  
Solvent Yellow 47 | Luxol Fast Yellow TN  
Irisol Fast Scarlet GN  
Ceres Brown R  
Ceres Blue R  
Propyl Astra Blue Iodide  
Basic Yellow 13 | Astrazon Yellow 8GL  
Basic Yellow 51 | Astrazon Yellow 5GL  
Basic Orange 27 based  
Solvent Brown 1 based | Ceres Brown RRN  
Solvent Yellow 21 based | Irisol Fast Yellow GRE  
Irisol Fast Scarlet BN  
Green 82 Powder for Inks  
Elvacite 2046, Copolymer Butylmethacrylate/Isobutylmethacrylate  
Umacryl KI, Methylmethacrylate - styrene copolymer  
Methacrylic Acid  
Sulryn random copolymer poly(ethylene-co-methacrylic acid), EMAA  
Fluorescent Brightener | Rylux BSU  
Fluorescent Brightener 49 | Rylux BS  
Plexigum, Isobutyl methacrylate

---

**IS NIR Spectra – Complete Collection**

6049 NIR spectra

Copyright © NICODOM 2012

NICODOM Ltd., Hlavni 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

[ftirspectra@gmail.com](mailto:ftirspectra@gmail.com)

Fluorotelomer Vidax 1000  
Fluorescent Brightener | Blancophor Sample B  
Tungophen B  
Fluorescent Brightener | Blancophor Sample A  
Fluorescent Brightener 140 | Rylux VPA  
Fluorescent Brightener | Uvitex Sample B  
Dammar Aged  
Morez 300, Styrene/Acrylic Resin  
Fluorescence S  
Basic Blue 11 | Victoria Blue R  
Turquoise methylene JSA Extra  
Acid Yellow | Alphanol Fast Yellow GG  
Acid Yellow 42 based | Alphanol Fast Yellow F  
Acid Red | Alphanol Fast Red R  
Acid Red | Alphanol Fast Brilliant Red BL  
Acid Red | Alphanol Fast Red 5B  
Acid Violet | Alphanol Fast Violet FR  
Acid Blue | Alphanol Fast Blue FBL  
Acid Blue 281 | Alphanol Fast Blue 2R  
Acid Blue | Alphanol Fast Blue HFL  
Acid Green | Alphanol Fast Green GWA  
Direct Blue 86 based | Alphanol Fast Turquoise Blue GL  
Acid Brown | Alphanol Fast Brown 3GL  
Acid Brown | Alphanol Fast Brown BC  
Anthralan Yellow G  
Acid Red | Anthralan Red HGG  
Acid Red 42 | Anthralan Red B  
Acid Red 81 | Anthralan Red BBT  
Acid Blue 228 | Alphanol Fast Turquoise Blue SW  
Acid Yellow | Amido Flavine FFP  
Amido Yellow R  
Basic Yellow | Amido Yellow 3GL  
Food Blue 5 based | Patent Blue V  
Remacryl Brilliant Red 4G  
Remacryl Red 3BL  
Remacryl Brilliant Red 7B  
Food Yellow 3 based | Orange RGL  
Acid Orange 7 based | Special Orange H  
Suede Skin Fast Black KHN  
Bayplast Yellow E2G  
Spun Rayon Fast Black B 300  
Acid Brown | Coranil Brown HEMR  
Reflex Blue R 51  
Acid Blue 185 | Turquoise Blue 5G

---

**IS NIR Spectra – Complete Collection**

6049 NIR spectra

Copyright © NICODOM 2012

NICODOM Ltd., Hlavni 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

[ftirspectra@gmail.com](mailto:ftirspectra@gmail.com)

Disperse Red 374  
Acid Yellow 250  
Disperse Red 277 based  
Solvent Brown 1 based | Fat Brown RR  
Solvent Red 1 | Fat Red G  
Vat Red 41 | Hostasol Red 5B  
Hostasol PV Fast Light R01  
Solvent Orange 41 | Savinyl Orange RLS  
Solvent Blue 45 | Savinyl Blue RS  
Solvent Red 91 | Savinyl Red 3BLS  
Acid Violet 47 | Sandolan Fast Violet P-3RL  
Pigment White 6 based | Sandospense White E-TD  
Acid Red | Sandolan Walk Red N-FBL 150  
Acid Blue 182 | Sandolan Blue E-HRL  
Acid Green | Sandolan Green MF-BL  
Solvent Black 27 | Savinyl Black NS  
Pigment Violet 19 based | Permajet Red E5B  
Pigment Yellow 101 | Lumogen Yellow S 0790  
Solvent Orange 3 | Neptun Orange Base 206  
Solvent Red 49 | Neptun Red Base 543  
Solvent Blue 4 | Neptun Blue base 627  
Solvent Yellow 162 | Neptun Yellow 075  
Solvent Blue 64 | Neptun Blue 698  
Neptun Black FF 4247  
Mixture | Basonyl Black X-11  
Solvent Orange 54 | Neopen Orange 252  
Basacryl Navy Blue X2 RN  
Basic Yellow 65 | Basacryl Yellow X2 GL  
Basic Violet 35 | Basacryl Violet XFL  
Basacryl Red YCN  
Basic Yellow 49 | Basacryl Yellow X-RL  
Basacryl Yellow X7 GLS  
Acid Red 114 based | Basacryl Brilliant Red BG  
Basic Blue | Basacryl Blue FR  
Basic Blue 3 based | Basacryl Blue Green X-5G  
Pigment Inorganic | Nacrelon R20  
Pigment Violet 37 based | Microlith Violet  
Pigment Inorganic | Iriodin 300 Gold pearls  
Pigment Inorganic | Iriodin 355 Glitter Gold  
Pigment Inorganic | Iriodin 110 Fine Satin  
Pigment Yellow 83 based | Helio Fast Yellow HRP  
Pigment Yellow 17 based | Helio Fast Yellow GGP  
Solvent Yellow 135 | Macrolex Fluorescence Yellow FP  
Solvent Green 28 | Macrolex Green G

---

**IS NIR Spectra – Complete Collection**

6049 NIR spectra

Copyright © NICODOM 2012

NICODOM Ltd., Hlavni 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

[ftirspectra@gmail.com](mailto:ftirspectra@gmail.com)

Acid Yellow | Special Yellow N-S  
Acid Blue 24 | Acilan Brilliant Blue R  
Acilan Amido Blue GGR  
Chromleder Brilliant Black 2R  
Isolan Black K NGL  
Direct Blue 199 based | Sirius Light Turquoise Blue FBLL  
Nussbaum Brown N  
Acid Yellow 49 | Supracen Yellow GR  
Supranol Cyanine 7BF  
Acid Blue 220 | Supramin Blue FRW  
Blue 63  
Basic Blue 52 | Astra Blue 3R  
Bayscript Special Red T  
Basic Yellow 76 | Astrazon Brilliant Flavine GFL  
Basic Yellow 63 | Astrazon Yellow 8GSL  
Basic Blue 62 | Astrazon Blue 5RL  
Basic Violet 21 | Astrazon Violet F3RL  
Basic Red 115 | Astrazon Red F2BL  
Alchemist Oil of Delft Thick  
Dammar Varnish  
Polishing Oil  
Mastic Varnish  
Cere-Stucco Soap  
Fluorescent Orange  
Marlipal® 1618/25  
Terpineol  
Larch Turpentine Oil  
Perilla Oil  
Fluorescent Yellow  
Mowilith® DM 778  
Latex Milk  
Solvent Red 118 | Red 471  
Orotan 731 K  
Wax-Paint Medium  
Cembran Pine Oil  
Silver Fir Turpentine Balsam  
Orange Turpentine Oil  
Walnut Oil  
Primal® WS 24  
Solvent Yellow 146 | Orasol Yellow 4GN  
Orange Madder  
Rubio Mineral Ground, Varnish  
Solvent Blue 70 based | Blue 807  
Olibanum

---

**IS NIR Spectra – Complete Collection**

6049 NIR spectra

Copyright © NICODOM 2012

NICODOM Ltd., Hlavni 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

[ftirspectra@gmail.com](mailto:ftirspectra@gmail.com)

Poppy Oil  
Mowilith® SDM 5, PVAc  
Tolu Balsam  
Lascaux® Beva 371 Film Thin  
Vulpex Liquid Soap  
Galbanum  
Plastorit® Super  
Mowiol® 4-88, PVA  
Solvent Black 3 | Sudan Black B  
Archäocoll 2000  
Solvent Orange 99 | Orasol Orange RG  
Pigment Green 17 based | Haematite-Chrome Oxide  
Laropal A-81, Ketone based Resin  
Acrylic polymer for emulsification  
Franklin Hyde Glue  
Pigment Brown 33 | Zinc-Iron Chrome Brown  
Manila Copal  
Gum Asafoetida  
Pigment White 1 | Flake White  
Keton Resin N  
Sepiolite  
Gum Benzoin  
Pigment White 20 | Muscovite Mica  
Phenolic Resin  
Colophony, Amber  
Sandarac  
Gum Elemi  
Fu-Nori  
Propolis  
Colophony, Rosin  
Pigment Historic | Atacamite  
Dammar  
Copal Varnish  
Lascaux® Gesso  
Attapulgit  
Pigment Yellow 119 | Zinc-Iron Brown Dark  
Lascaux® Acrylic Glue 498 HV  
Cyclododecan  
Fish Glue  
Canada Balsam  
Burgundy Resin  
Shellac Powder  
Guajac Resin  
MS2A-Resin

---

**IS NIR Spectra – Complete Collection**

6049 NIR spectra

Copyright © NICODOM 2012

NICODOM Ltd., Hlavni 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

[ftirspectra@gmail.com](mailto:ftirspectra@gmail.com)

Laponite® RD  
Paraloid™ B 48 N  
Paraloid™ B 66  
Pigment Glitter | Fish Silver Powder  
Polyethylene(2)-oleyl-ether  
Regalrez® 1126  
Savinyl Red SR 90  
Paraloid™ B 72  
Fluorescent Red  
Solvent Orange 11 | Orasol Orange G  
Solvent Blue 67 | Orasol Blue GN  
Mordant Yellow 1 | Alizarine Yellow GG  
Olive Oil Soap  
Bis(2-hydroxyethyl) cocoalkyl-amine  
Gillot Watercolor Medium  
Fluorescent Violet  
Gneis Green, dried  
Tylose® C 300  
Claw Oil  
Phosphorescent pigment, Blue  
Turpentine Rosin Oil thick, Varnish  
Pigment Glitter | Vermiculite  
Tung Oil Varnish  
Mastic  
Plexigum® PQ 611  
Acronal 500 D  
Oil Paint Medium  
Reactive Orange 91 | Cibacron Yellow F-3R  
Carbopol® EZ 2  
Thickener ASE 60  
Disperse Aid  
Tixogel® VZ  
Copper Blue  
Owatrol Oil  
Pigment Red 259 | Ultramarine Red  
Solvent Yellow 82 | Yellow 157  
Burgundy Red Ochre Deep  
Pigment Green 50 | Cobalt Bottle Green  
Gum Accroides  
Wax V  
Pigment Yellow 42 | Translucent Yellow  
Veronese White, colored marble  
Pigment White 19 | China Clay  
Pigment Green 17 based | Haematite-Chrome Oxide

---

**IS NIR Spectra – Complete Collection**

6049 NIR spectra

Copyright © NICODOM 2012

NICODOM Ltd., Hlavni 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

[ftirspectra@gmail.com](mailto:ftirspectra@gmail.com)

Peru Balsam  
Terra Alba  
Pigment Red 107 | Antimony Red  
Gum Copal  
Fluorescent pigment, White  
Fir Turpentine  
Copaiba Balsam  
Resin Amber  
Dye Natural | Lycopodium  
Pigment Red 83 | Alizarine Crimson dark  
Iwa-Enogu Usukuchi-Murasaki  
Han-Blue Fine  
Pigment Green 19 | Cobalt Green bluish A  
Natural Yellow 24 based | Pipe Gamboge  
Pigment Yellow 184 | Bismuth-Vanadate Yellow Lemon  
Natural Red 22 | Sandalwood  
Pigment Orange 61 | Isoindolor Orange  
Reactive Orange 133 | Cibacron Orange F-BR  
Pigment Violet 15 | Ultramarine Violet Medium  
Pigment Red 7 | Permanent Red  
Blue Bice  
Pigment Mineral | Jade  
Fuchsite  
Turquoise, sky-Blue  
Henna Red  
Pigment Mineral | Jarosite  
Pigment Earth | Jarosite  
Pigment Mineral | Brown Ochre from Elba  
IWA-Enogu® Kibeni  
Tourmaline DENKISEKI-MATSU  
Pigment Japanese | Sodalite SHIUN-MATSU  
Pigment White 25 | Marienglas  
Pigment Black 18 | Pencil Clay  
Dye Natural | Reseda  
Natural Red 25 | Lac Dye  
Pigment Japanese | Garnet SEKURA-NEZUMI  
Dye Natural | Wild Saffron  
Pigment Historic | Han-Purple  
Pigment Violet 16 | Manganese Violet  
Pigment Mineral | Aragonite  
Pigment Yellow 129 | Irgazin Yellow, Greenish  
Pigment Mineral | Basalt Powder Black  
Dye Natural | Safflower  
Pink Coral USUKUTI-SANGO-MATSU

---

**IS NIR Spectra – Complete Collection**

6049 NIR spectra

Copyright © NICODOM 2012

NICODOM Ltd., Hlavni 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

[ftirspectra@gmail.com](mailto:ftirspectra@gmail.com)

Pigment Mineral | Cavansite  
Pigment Historic | Atramentum  
Taunus Ochre, light  
Ploss Blue  
Pigment Historic | Vivianite  
Pigment Orange 48 | Cinquasia Gold, Red-gold  
Pigment Blue 30 based | Azurite Gunjyou  
Blue Verditer  
Pigment Mineral | Vesuvianite  
Pigment Blue 30 based | Azurite Natural  
IWA-Enogu® Iwashikou  
Pigment White 5 | Lithopone  
Natural Violet 1 | Tyrian Purple  
Turmeric Powder  
Dye Natural | Gofun Shirayuki From Oyster Shells  
Pigment Orange 51 | Paliogen Orange  
Pigment Orange 49 | Cinquasia Gold, brown-gold  
Pigment Brown 7 | Pompei Red  
Pigment Green 24 | Ultramarine Green, genuine  
Pigment Blue 71 | Zirconium Cerulean Blue  
Pigment Yellow 173 | Isoindolinon Yellow  
Pigment Black 10 | Asphaltum  
Pigment Blue 35 | Cobalt Blue Light  
Natural Brown 7 | Walnut Hulls  
Ultramarine Ash  
Natural Brown 3 | Catechu  
Natural Red 31 | Dragon's Blood IYDAHA  
Pigment Violet 14 | Cobalt Violet Dark  
Pigment Blue 36 based | Cobalt Blue, Greenish  
Pigment Blue 27 | Prussian Blue  
Pigment Blue 74 | Cobalt Blue Dark  
Pigment Violet 49 | Cobalt Violet Brilliant Light  
Pigment Green 23 | Green Earth  
Natural Black 1 | Logwood, pieces  
Natural Orange 4 | Annatto Seeds  
Pigment Japanese | Tiger-Eye TYOUJICHA  
Pigment Red 254 based | XSL Irgazin Red DPP  
Cibacron Black F-2 B  
Pigment Yellow 40 | Cobalt Yellow  
Pigment Yellow 43 based | Iron Oxide Yellow  
Malachite  
Pigment Blue 36 based | Cobalt Blue Turquoise Dark  
Natural Yellow 2 | Reseda Lemon  
Henna Red

---

**IS NIR Spectra – Complete Collection**

6049 NIR spectra

Copyright © NICODOM 2012

NICODOM Ltd., Hlavni 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

[ftirspectra@gmail.com](mailto:ftirspectra@gmail.com)

Vat Blue 1 based | Indigo, Genuine Indian  
Pigment Yellow 39 | Orpiment, genuine  
Pigment Red 206 | Cinquasia chestnut brown  
Verdigris, Copper(II) Acetate Hydrate  
Verona Green Earth with Virid  
Dye Natural | Aloe  
Natural Red 24 | Redwood, Brazilwood  
Pigment Mineral | Conicalcite  
Pigment Green 18 | Viridian Green  
Dragon's Blood, powder  
Pigment Violet 37 based | Dioxazine Violet  
Pigment Historic | Red Coral Koikuti-Sango-Matsu  
Natural Orange 2 | Resina Kamala  
Pigment Earth | Bohemian Green Earth, imitation  
Bristol Yellow pale  
Pigment Japanese | Malachite MATSUBA-ROKUSYOU  
Rhodochrosite, CI 77748  
Preventol® ON extra  
Pigment Brown 23 | Gubbio Red  
Nicosia Green  
Pigment Green 26 | Cobalt Oxide Green Blue  
Natural Red 20 | Alkanet  
Pigment Blue 32 | Smalt, standard grind  
Natural Red 8 | Madder Lake  
Pigment Yellow 159 | Intensive Yellow  
IWA-Enogu® Shinsia  
Pigment Yellow 43 based | French Ochre JTCLES  
Natural Yellow 24 based | Gamboge  
Pigment Violet 19 based | Hostaperm Red  
Pigment Yellow 31 | Baryte Yellow  
Pigment Red 4:1 | Carmine Naccarat  
Solvent Red 27 | Wax Red 5B02  
Reactive Blue 268 | Cibacron Blue FN-G  
Burgundy Yellow Ochre  
Pigment Earth | Terra Ercolano  
Garnet Powder Red  
Natural Yellow 11 | Yellow Wood  
Ammonium iodide  
Manganese(II) nitrate tetrahydrate  
Ammonium fluoride  
Antimony(IV) oxide  
Diethyl sulfite  
Ethylamine  
Diethyl carbonate

---

**IS NIR Spectra – Complete Collection**

6049 NIR spectra

Copyright © NICODOM 2012

NICODOM Ltd., Hlavni 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

[ftirspectra@gmail.com](mailto:ftirspectra@gmail.com)

Isopentyl nitrite  
Ethyl hydrogen malonate  
Alloxan monohydrate  
Hydroxylamine-o-sulfonic acid  
Trimethylsulfonium iodide  
Sodium cyanoborohydride  
Potassium thiocyanate  
Acetaldehyde  
Poly(ethylene glycol) 6000  
Gadolinium(III) oxide  
Praseodymium(III) oxide  
Holmium(III) oxide  
D-Cellobiose octaacetate  
Sodium tetraoxiodate(VII)  
Lithium borohydride hydrate  
Polyester, adipic Acid  
Poly(methacrylic acid) #2  
Poly(vinyl acetate) #2  
Silipor 600  
Aluminumsilicate Calcium Sodium  
Poly(ethylene glycol) 8000  
Calcium Chloride Powder  
Silicagel  
Carbonate Sodium Decahydrate  
Ethyl Cellulose  
Spheron Oxin 1000 LC  
N-Acetyl-3,5-diiodo-L-tyrosine methyl ester  
D-Glycogen  
Lactose octaacetate  
Ammonium iodide  
Sodium isopropyl xanthate  
1,2-Dibromo-1,1-difluoroethane  
Lithium citrate  
Glyceryl triacetate, triacetin  
Retinol acetate  
Hydrochelidonic acid  
Dextran Blue  
Brown FK Kipper Brown  
Furcelleran  
Polyethylene plasticized #2  
Polyacrylate  
Luviskol K 60  
Luviskol VA 64  
Luviskol K 90

---

**IS NIR Spectra – Complete Collection**

6049 NIR spectra

Copyright © NICODOM 2012

NICODOM Ltd., Hlavni 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

[ftirspectra@gmail.com](mailto:ftirspectra@gmail.com)

Luviskol K 17  
Luviskol K 80  
Lutensit AFK  
Lutonal I 60  
Lutensol FA 12  
Pluronic PE 3100  
Plurafac LF 221  
Pluronic RPE  
Plurafac LF 403  
Disponil O 5  
Disponil NP 307  
Disponil AES 63  
Disponil SUS 87  
Disponil OSS 50  
Aerosil MA 80  
Aerosol 501  
Aerosol A 102  
Hydropalat 759  
Polyphobe 106 E  
Polyolefin CP 164-1  
Synperonic PE/P 103  
BYK 073  
Vinnapas UW-1  
Nopco 8034  
SER/AD FN 225  
Synperonic PE/L 64  
Preventol HS 75  
Piccotac 95 E  
Maranil A  
Eastotack H 100 R  
Staybelite resin E  
Regalite R 7100  
Surfynol PSA-336  
Hydrazine hydrate  
Worléésin MS 265  
Dibasic Ester Mixture  
Optigel SH  
Sulfate Barium  
Busperze 275  
Bentone SD 1  
Exxate 700  
Agitan 731  
Agitan P 813  
Agitan P 823

ir-spectra.com IS NIR Spectra

---

**IS NIR Spectra – Complete Collection**

6049 NIR spectra

Copyright © NICODOM 2012

NICODOM Ltd., Hlavni 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

[ftirspectra@gmail.com](mailto:ftirspectra@gmail.com)

Agitan P 833  
Agitan VP P 804  
Aquacer 535  
Baysilone OL 31  
Baysilone VP AI 3468  
Dow Corning 1-6184  
Dow Corning 84  
BYK 341  
BYK 9076  
CHS Epoxy 200V 55  
Dehydran 1620 1:1BG  
Dapro DF 1161  
Dapro DF 1760  
Dapro DF 4164  
CHS Epoxy 530  
Cardolite NC 541 LV  
Flat-ayd FA-W-34  
Maincote AU - 28  
Bayhydrol UH 135  
Bayhydrol UH 240  
Bayhydrol UH 340  
Luhydran 8975  
Troysan Polyphase EC 17  
Contrapen PR/194  
Contrapen PR/325  
Contrapen PR/388  
Busan 30L  
Slip-Ayd SL 145E  
Aquadlip 671  
Colorol E  
Metolat P 871  
Baragel 3000  
Superchlone CR 40  
Tamol NNO  
Aluminium trihydrate  
Jeffamine M-600  
Jeffamine M-2070  
Disponil FES 77  
Lutensol FSA 10  
Triton CF-32  
Disponil LDBS 55  
Tergitol NP-9  
Acronal LR-8820  
Disponil AES 48 e

---

**IS NIR Spectra – Complete Collection**

6049 NIR spectra

Copyright © NICODOM 2012

NICODOM Ltd., Hlavni 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

[ftirspectra@gmail.com](mailto:ftirspectra@gmail.com)

Disponil FES 61 IS  
Hydropalat 7003  
Disponil AFX 780  
Plioway Ultra 200  
Pliolite VTAC-L  
Plioway Ultra G20  
Pliolite AC3-H  
Luwax AF 32  
Igepal CO 210  
Acresin DS 3500  
Acrysol RM 8  
Telalit 180  
Thixatrol Plus  
Fascat 4350  
Telalit 410  
Disponil AAP-43  
Cymel 303 LF  
Cymel 323  
EFKA 34  
Heucophos - SAPP  
Heucophos - ZPA  
Heucophos - ZPZ  
Rheolate 255  
Bayhydrol VP LS 2235  
Borchigel PW-25  
SER AD FA 179  
Dapro DF 911  
SER AD FA 607  
SER AD FX 504  
SER AD FX 511  
BYK 051  
Heucophos - ZPM  
Epikote 1009  
Atmer 151  
Exxal 18  
Polypropyleneglycolmonoacrylate  
Epikote 1001-X-75  
Epikote DX 874-X-90  
Euredur 43  
Halox BW 111  
Halox CW 2230  
Heucophos - ZPO  
Moly White - 101  
Moly white - 212

---

**IS NIR Spectra – Complete Collection**

6049 NIR spectra

Copyright © NICODOM 2012

NICODOM Ltd., Hlavní 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

[ftirspectra@gmail.com](mailto:ftirspectra@gmail.com)

Moly White - MZAP  
Moly white - ZPN  
Shieldex CP-4  
SER AD DP FS 531  
Epilink DP 700  
Sicor SPO  
Sylodex 04  
Sylodex 07  
Sylodex 08  
Disperbyk 107  
BYK 034  
Disperbyk P 104 S  
Blister Free 3  
Blister Free 45  
Blister Free 75  
Estasol  
BYK 035  
Luviskol K 30  
Agitan 760  
Aquacer 593  
Disponil A 1580  
Lupasol P  
Althosan MB  
Disponil FES 993  
Agitan P 801  
Agitan P 800  
EFKA 4010  
Additol VXW 6210  
SER AD BEZ 75  
SER AD FA 196  
Serdas MOM  
SER AD FN 1566  
EFKA 2022  
EFKA 2720  
EFKA 3232  
SER AD FK 321  
Nuodex Combi HS  
Nuodex Combi QS  
Nuodex Zr 18 %  
Tego Airex 900  
EFKA 3883  
Serdas GBR  
SER AD FX 1070  
SER AD FX 2050

ir-spectra.com IS NIR Spectra

---

**IS NIR Spectra – Complete Collection**

6049 NIR spectra

Copyright © NICODOM 2012

NICODOM Ltd., Hlavni 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

[ftirspectra@gmail.com](mailto:ftirspectra@gmail.com)

Serdas 7540  
SER AD FS 530  
SER AD FA 601  
SER AD FX 9010  
Additol XW 330  
Additol XL 425  
Additol XL 260  
Additol XL 109  
Additol XL 197  
Additol XL 251  
Additol VXL 1105  
Additol VXL 4910  
Additol VXL 4920  
Additol VXL 4929  
Additol VXW 4932  
Additol VXW 4952  
Additol VXW 6203  
EFKA 66  
EFKA 4009  
EFKA 4050  
Bentone 27  
Bentone SD-2  
Bentone SD-3  
EFKA 4520  
Dyhard PC 45  
Wacor ZBP-M  
Zinc Borate  
Vestagon B 68  
Leukonol LBA-2  
Metolat FC 530  
Metolat P 853  
Metolat P 860  
Metolat P 861  
Metolat TH 75  
Nuodex Combi APB  
Nuodex Combi LS  
Versamine I 74  
Bentone 38  
Telalit 160  
Disperbyk 185  
Butrol 23  
Anti-Terra 206  
Ancamine 2143  
Disodium oxalate

ir-spectra.com IS NIR Spectra

---

**IS NIR Spectra – Complete Collection**

6049 NIR spectra

Copyright © NICODOM 2012

NICODOM Ltd., Hlavni 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

[ftirspectra@gmail.com](mailto:ftirspectra@gmail.com)

Nitrate Zinc hexahydrate  
Nitrate Magnesium hexahydrate  
BYK 031  
BYK 033  
EFKA 27  
Butanox M-50  
BYK 366 P  
BYK W 966  
Disperplast 1142  
Disponil STO 120 F1  
Luperox DP10  
Dynasylan BSM 100 W  
BYK 360 P  
Bentone EW  
Ancamide 220 X 70  
Ancamine 1618  
Ancamine 2320  
Ancamine 2379  
BYK 040  
Lanco Disperse 10 UN  
BYK 405  
Desmophen RD 181 X  
Disperbyk 2001  
Setal 169 SS-67  
Setal 90173 SS-50  
Carbonate Barium  
Rhodapex CO-436  
Disponil AES-60  
Feldspar  
Carbonate Magnesium  
Roskydal K68  
Hycar 1300 x 42  
Wollastonite  
Additol XW 329  
BYK 364 P  
BYK 365 P  
BYK 151  
Zinc sulphate  
Lutensol AP-9  
Fatty Acids  
Rhodoline DF 6120  
Polyol TP 200  
Additol VXW 6206  
Addid 900

---

**IS NIR Spectra – Complete Collection**

6049 NIR spectra

Copyright © NICODOM 2012

NICODOM Ltd., Hlavni 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

[ftirspectra@gmail.com](mailto:ftirspectra@gmail.com)

SER AD FX 521  
Disponil TA 430  
Rhodoline DF 311M  
Rhodoline DF 6681  
Rhodoline DF 999  
Blister free 54  
Anti-Terra 203  
Schwego PUR 8050 VP  
Schwego Wett 8075  
Epoxide 8  
Poly{[2-(5-phenyl-1,3,4-oxadiazole-2-yl)-1,4-phenylene]-alt-(9,9-dioctylfluorene-2,7-diyl)}  
Poly[3-(2-heptyl-1,2,4-oxadiazol-5-yl)thiophene-2,5-diyle-2,5-dioctyloxy-1,4-phenylene]  
Poly[9,9-di(4-octyloxyphenyl)fluorene-2,7-diyl-1,4-phenylene]  
Poly[3-(2-heptyl-1,2,4-oxadiazol-5-yl)thiophene-2,5-diyle-9,9-dioctylfluorene-2,7-diyle]  
Dabco DC 193  
Lewatit M610  
DER 321  
BYK 080  
SER AD FA 620  
Hydroxide lithium  
Nitrate Lithium  
Zinc acetate  
Disperbyk 2000  
Zinc Chromate  
Grilonit RV 1814  
Dianol 22  
Laromin A 327  
Barium hydroxide octahydrate  
Ammonium ferric sulfate dodecahydrate  
Polypox R24  
Ammonium sodium phosphate  
Fluorad FC 430  
DER 353  
CoatOSil 2400  
Edaplan VP LA 420  
Cycat 500  
Plastorit M  
SER AD FA 579  
SER AD FN 265  
Surfynol DF 58  
Surfynol 104 S  
SER AD FX 9085  
Sylodex H2  
Nuosperse 657

---

**IS NIR Spectra – Complete Collection**

6049 NIR spectra

Copyright © NICODOM 2012

NICODOM Ltd., Hlavni 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

[ftirspectra@gmail.com](mailto:ftirspectra@gmail.com)

Surfynol SEF  
Carbonate Rubidium  
Ammonium Heptamolybdate Tetrahydrate  
Sodium bismuthate  
Hemoglobine  
Carbonate Zinc  
Thiocyanate barium trihydrate  
Ammonium hexafluorophosphate  
Ammonium oxalate monohydrate  
Lithium acetate  
Lithium Tetrafluoroborate  
Sulfate Lithium  
Carbonate Lithium  
Sodium paraperiodate  
Sodium Azide  
Sodium Thiocyanate Dihydrate  
Sodium Tetrafluoroborate  
Sodium Acetylacetonate  
Sodium potassium tartrate tetrahydrate  
Nitrate Bismuth pentahydrate  
Ammonium metavanadate  
Sodium molybdate dihydrate  
Hexaammonium molybdate  
Sodium tungstate dihydrate  
Tungstic acid  
Carbonate Manganese (II)  
Sulfate Manganese Tetrahydrate  
Nitrate Ferric nonahydrate  
Ammonium Iron(II) Sulfate Hexahydrate  
Sodium hexanitrocobaltate(III)  
Potassium hexacyanocobaltate(III)  
Ferrocenecarboxaldehyde  
(Ferrocenylmethyl)trimethylammonium iodide  
L-Antimony potassium tartrate  
Carbonate Cesium  
Cesium acetate  
Sulfate Cesium  
Nitrate Cesium  
Rubidium Acetate  
Nitrate Rubidium  
Sulfate Rubidium  
2-Methylthiophene  
Ethyl vinyl ether  
Fluorobenzene

---

**IS NIR Spectra – Complete Collection**

6049 NIR spectra

Copyright © NICODOM 2012

NICODOM Ltd., Hlavni 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

[ftirspectra@gmail.com](mailto:ftirspectra@gmail.com)

Food Red 17 | Allura Red AC  
L-meso-Lanthionine  
Calcium carbonate Chalk  
DL-Lanthionine  
Poly(ethylene glycol) 2000 dimethyl ether  
N-Trityl Aminoisopentylphosphonic Acid Diethylester  
Microcrystalline cellulose  
N-Benzyloxycarbonyl Aminobenzylphosphonic Acid Diphenyl Ester  
N-Trityl Aminoethylphosphonic Acid Diethylester  
N-Trityl Aminobutylphosphonic Acid Diisopropylester  
N-Benzyloxycarbonyl Aminomethylphosphonic Acid Diphenyl Ester  
N-Trityl Aminobutylphosphonic Acid Diethylester  
N-Trityl Aminomethylphosphonic Acid Dimethylester  
Aminomethylphosphonic Acid Dimethylester  
(Aminoisobutyl) methylphosphinate  
(Aminoethyl) methylphosphinate  
(Aminomethyl) methylphosphinate  
(Aminomethyl) t-butylphosphinate  
(Morpholin-4-ylmethyl)phenylphosphinic acid  
(Piperazin-1-ylmethyl)phenylphosphinic acid  
Calcium silicate  
Potassium superoxide  
Catalase from bovine liver  
Glucose oxidase  
Cyclopentadienyliron dicarbonyl dimer  
Sephadex G-50  
Sephadex G-15  
Poly(styrene-co-divinylbenzene) aminomethylated  
Polyethyleneimine  
Poly(styrene-co-divinylbenzene) amine functionalized  
Polystyrene-b-poly(4-vinyl pyridine)  
Sephadex G-25  
Ferrocene, [(1E)-3-ethoxy-3-oxo-1-propenyl]- (9CI)  
Sodium Orthovanadate  
Antioxidant SP  
Antioxidant 301  
Nonox SP  
Wingstay T  
Zalba special  
Altofane BB 123  
Altofane M2B  
Sterically hindered bisphenol, butylidenebisxylenol type  
Nonox EXN  
Nonox EXP

---

**IS NIR Spectra – Complete Collection**

6049 NIR spectra

Copyright © NICODOM 2012

NICODOM Ltd., Hlavni 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

[ftirspectra@gmail.com](mailto:ftirspectra@gmail.com)

ASM AFC  
Age Rite Gel  
Agerite Stalite  
Styrenated diphenylamine  
AgeRite Stalite S  
Thermoflex A  
Age-Rite Hipar  
Santoflex 134  
Altofane RA  
Age-Rite Resin  
Agerite Superflex  
Ammonium phosphomolybdate  
Boric Acid-D3  
Aluminium Oxide Hydrated  
Nickel Isothiocyanate  
Ammonium paratungstate  
Ferric ammonium oxalate  
Sulfate Ferric  
Copper Hexafluorosilicate Hexahydrate  
Phosphotungstic acid  
Ammonium Hydrogen Imidodisulphate  
Sodium Dithionate Anhydrous  
Manganese (III) acetate dihydrate  
Calcium hypophosphite hydrate  
Lithium hydroxide with Carbonate  
Nitrate Calcium  
Barium phosphinate  
Potassium pyrosulfate  
Strontium acetate  
Carbonate Strontium  
Phosphinic acid 50% in water  
Tetraethylammonium iodide  
Sulfate Vanadyl  
Sulfate Vanadyl Hydrate  
Amberlite 15  
Titanium Oxysulfate Sulfuric Acid Complex Hydrate  
Ferric ammonium citrate  
Amberlite CG-50  
Amberlite IRA 402  
Dowex 2x8  
Amberlite IRC-50  
Varionap  
Amberlite 410  
Dowex (R) 50 WX8

---

**IS NIR Spectra – Complete Collection**

6049 NIR spectra

Copyright © NICODOM 2012

NICODOM Ltd., Hlavni 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

[ftirspectra@gmail.com](mailto:ftirspectra@gmail.com)

Amberlite IR 45  
Zerolit 225  
Zerolit FF (SRA 70)  
Amberlite IR-48  
Naugard A  
Nonox BL  
Acorox 88  
Advastat E 49  
Styrenated phosphites  
Ferro 541 A  
Phosclere T 210  
Phosclere T 268  
Phosphit P  
Weston 619  
Weston 732  
Alkylated Phenyl Phosphites, n-octylated  
Alkylated Phenyl Phosphites, dilauryl  
Alkylated Phenyl Phosphites, t-octylated  
Mixture of hindered alkylated phenols  
Mark 158  
Mark 522  
Antiozonant AFD  
Antiozonant 70-TBPA  
Wingstay L 177  
Advastab  
Advastab BC-30  
Ca, Zn Carboxylate  
Mark 99  
Mark WS-X  
Beckopox EP 301  
Soybean Oil Epoxy Butylester  
Estabex 2349  
Advastab 15 MS  
Advastab 17 M  
Advastab 17 MOL  
Estabex E  
Meister OZ 2  
Meister Z 25 N  
Stanclere 173  
Citronella Oil From Asia  
Rosemary Oil  
Cobalt(III) hexamine chloride  
Mark 328  
Herax UTS

---

**IS NIR Spectra – Complete Collection**

6049 NIR spectra

Copyright © NICODOM 2012

NICODOM Ltd., Hlavni 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

[ftirspectra@gmail.com](mailto:ftirspectra@gmail.com)

Mark 548  
Mark 550  
Stabaxol PCD  
Mark 1306 A  
Cremophor  
Fatty acid polyglycol ester  
Advastab CZ-59  
Mark 133  
Mark 281 B  
Mark 518  
Mark 553  
Mark 67 D  
Mark 701  
Mark I 68  
Mellite 218 Ba/Cd  
Mellite 220  
Tin Soap Complex  
Ca, Zn Complex  
Meister O-3A  
Meister Z 41  
Stanclere 1250  
Dastib 256  
UV-absorber 318  
Ferro AM 101  
Dowex(R) 50WX4 Hydrogen Form  
Lewatit TP-260  
Dowex 1x2  
Ion exchanger III  
Dowex 50WX4  
Lewatit Monoplus M 500  
Amberlite(R) XAD-1180  
Amberlite XAD-4 (with Ethanol)  
Amberlite IRC-86  
Lewatit TP-207  
Cotton Lot#2  
Polyimide  
Poly(methylphenylsiloxane) fluid  
Tegosoft CT  
Rewoderm LI 420-70  
Epoxy Resin  
Polyurethane #1  
Poly(styrene:ethylacrylate)  
Silicate with cellulose - pigment  
Cellophane

---

**IS NIR Spectra – Complete Collection**

6049 NIR spectra

Copyright © NICODOM 2012

NICODOM Ltd., Hlavni 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

[ftirspectra@gmail.com](mailto:ftirspectra@gmail.com)

Polyisoprene cis #1  
Silicone rubber #1  
Polyester; terephthalate based  
Poly(styrene:acrylonitrile) (70% styrene)  
Nitrocellulose in mixture  
Dimeticone  
Ludipress  
Ethoxylated Lanoline 4000  
Sodium Tetraborate Decahydrate  
PEG 7 Glyceryl Cocoate  
Surfactant, Sodium Alkyl polyglycol ether sulfate  
Maltodextrin KMS X-50  
Brown HT Chocolate brown HT, in mixture  
Hydroxyethyl ether cellulose  
Poloxamer 105 (PEG-PPG copolymer)  
Dimethicone  
Maize Oil  
Sorbitan Oleate, Span 80  
Sodium Chloride  
Poly(ethylene glycol) 14000  
Carboxymethyl cellulose, Na salt  
Casein hydrolysate  
Soda Ash  
Starch Modified  
Opadry OY-B-22932 yellow - coating mixture  
Emetine dihydrochloride heptahydrate  
Silicon dioxide  
Calcium chloride dihydrate  
Iron(II) fumarate  
Opadry OY-S-8720 orange - coating mixture  
Polyethylene Oxide 4000000  
Oxide Aluminium  
Nacrelon S-30 (pearlescent pigment)  
Pareth 15 (ethoxylated alcohols)  
Sodium Metabisulphite  
Silicon dioxide  
Quartz  
Sulfophenylazochromotropic acid  
Dichlorobenzene  
3,5-Diiodo- L-tyrosine  
DL-Phenylserine  
Hydroxylamine sulfate  
Cinchonine sulfate  
N-Phenylglycine

---

**IS NIR Spectra – Complete Collection**

6049 NIR spectra

Copyright © NICODOM 2012

NICODOM Ltd., Hlavni 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

[ftirspectra@gmail.com](mailto:ftirspectra@gmail.com)

FD&C Blue 2, Indigo carmine, Indigotine  
Potassium hexacyanochromate  
Sodium Dihydrogen Phosphate monohydrate  
Sodium Polyphosphate  
Sodium Thiosulfate Pentahydrate  
Potassium Thiosulfate  
Tetramethylammonium hydroxide pentahydrate  
Potassium Tellurite  
Potassium Binoxalate  
Strontium Hydroxide Octahydrate  
Aluminosilicate  
Pumice Stone  
Nitrate Calcium Tetrahydrate  
Ammonium Fluoborate  
Potassium Dicyanoaurate  
Potassium Ferrocyanide Trihydrate  
Sulfate Aluminium Hydrate  
Barium Hydroxide with Carbonate  
Potassium (Poly)Sulfide  
Bismuth Subcarbonate  
Potassium Pyrophosphate  
Phosphorous acid  
Carbonate Magnesium Hydrate  
Disodium Pyrophosphate  
Sodium dichromate dihydrate  
Ammonium Hypophosphite  
Boron Trioxide  
Lithium hydride  
Ammonium dimolybdate  
Potassium Bisulfate  
Potassium Hydrogen Fluoride  
Cesium Hydroxide Monohydrate  
Sodium malonate  
2-Methyl-3-butanol  
Poly(dimethylsiloxane) 60000  
Ethyl-4-bromobutyrate  
Methyl propionate  
Hexamethylphosphorous triamide  
Poly(propylene glycol) 500  
cis-1,4-Cyclohexanedimethanol  
Dibutoxytetraglycol  
Tween 21, POE (4) sorbitan monolaurate  
Diethyl ketone  
Methyl pentyl ketone

---

**IS NIR Spectra – Complete Collection**

6049 NIR spectra

Copyright © NICODOM 2012

NICODOM Ltd., Hlavni 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

[ftirspectra@gmail.com](mailto:ftirspectra@gmail.com)

Butyronitrile  
Tributyl phosphite  
Pigment Rubine, Lithol Rubine BK  
Copper complexes of chlorophyllins and chlorophylls  
Potassium hexafluorophosphate  
Carbonate Potassium  
Hymono 9004  
Tegin 4100  
Cocoa Butter  
Hostacerin WO  
Elfan AT 84 G  
Tego Alcanol 1618  
Cerasynt IP  
MYRJ 52 S  
Polymer JR-400 UCARE  
MYRJ 53  
Alfol  
Montane 83  
Cetiol 1414 E  
Avocado Oil  
Veegum  
Polawax GP 200 (non-animal)  
AF Emulsion  
Cetiol V  
Parsol MCX  
Pecan Nut Oil  
Tego Betain CK KB 5  
Plantacare 2000 UP  
Cetiol OE  
Cetiol MM  
Montanov 68  
Vitazyme C  
Bovexan-P (colagen)  
Crolactil SSL  
Crosilk powder  
Cosmowax D  
Lemon Oil Type 2  
Myrrh Oil Type B  
Sassafras Oil  
Kerosene  
SYNTOPON A 100  
L -Threonine  
Ferric acetate hydrate  
Potassium hypophosphite

---

**IS NIR Spectra – Complete Collection**

6049 NIR spectra

Copyright © NICODOM 2012

NICODOM Ltd., Hlavni 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

[ftirspectra@gmail.com](mailto:ftirspectra@gmail.com)

Cobalt(II) carbonate hydrate  
Magnesium hydrogencitrate  
Ammonium cobalt(II) sulfate hexahydrate  
Ammonium cupric chloride dihydrate  
Chromium nitrate nonahydrate  
Bismuth(III) sulfate  
Calcium Tartrate  
Calcium citrate tetrahydrate  
Aluminium acetate  
Sodium ferrocyanide  
Strontium sulfate  
Strontium chloride hexahydrate  
Rubidium chromate(VI) hydrate  
Phosphonic acid, disodium salt  
2-(Dimethylamino)ethyl 4-butylaminobenzoate .HCl  
Acetophenetidin  
Antipyrine  
Rutine Sodium Sulfate  
Valerianae Dry Extract  
Bismuth Subnitrate  
Sodium Citrate Dihydrate  
Basic Green 1 | Brilliant Green  
Oil of Pine  
Oil of Thyme  
Bergamot Oil  
Cod-Liver Oil  
Oil, Coriander Oil  
Eucalyptus Oil  
Methocel K15  
Gelatin  
Food Red 9 | Amaranth  
Annatto bixin norbixin, liquid in PG  
Brilliant Blue FCF  
Capsanthin, in oil  
Carminic Acid  
Carmoisine B  
Beetroot Red, Betanin, with starch  
Enocyanine  
Erythrosin B  
L-a-Phosphatidylcholine from turkey egg yolk  
Emulsifier E 2149, Stearyl Alcohol and Esters Mixture  
Carboxymethylamylopectine Sodium  
Glycerin  
Tween 20, POE (20) sorbitan monolaurate

---

**IS NIR Spectra – Complete Collection**

6049 NIR spectra

Copyright © NICODOM 2012

NICODOM Ltd., Hlavni 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

[ftirspectra@gmail.com](mailto:ftirspectra@gmail.com)

Agar with maltose extract  
Maize starch extract  
Sucrose  
Marine fish oil  
Sulfate Manganese  
Calcium Chloride Powder  
Sulfate Aluminium Octadecahydrate  
Adonitol  
Neo Dragocid Liquid (mixture)  
Poly(ethylene maleate)  
Poly(1,4-butylene terephthalate)  
Poly(isopropylene endomethylene tetrahydrophthalate)  
Poly(diethylene phthalate)  
Poly(ethylene maleinadipate)  
Epon 828  
Epoxy resin, Bisphenol A + epichlorohydrin  
Acrylic Copolymer with Styrene  
Poly[(acrylic acid):(butyl acrylate):(ethyl acrylate):(2-ethylhexyl acrylate):(methacryl amide)]  
Poly[(acrylic acid):(butyl acrylate):(methyl methacrylate)]  
Poly[(acrylic acid):(butyl acrylate):]  
Poly[(acrylic acid):(butyl acrylate):(methacryl amide):(methyl methacrylate)]  
Poly[(acrylic acid):(2-ethylhexyl acrylate):styrene]  
Poly[(acrylic acid):(butyl acrylate):styrene]  
Poly[acrylonitrile:(butyl acrylate):(methacryl amide)]  
Poly[(acrylic acid):(butyl acrylate):2-(hydroxyethyl methacrylate):styrene]  
Poly[(acrylic acid):(ethyl acrylate):(2-hydroxyethyl methacrylate):styrene]  
Poly[(acrylic acid):(ethyl acrylate):(methacryl amide):(methyl methacrylate)]  
Poly[(ethyl acrylate):(methacrylic acid)]  
Poly[(acrylic acid):(butyl acrylate):(methacryl amide):styrene]  
Poly[(acrylic acid):(butyl acrylate):(ethyl acrylate)]  
Poly[(acrylic acid):(maleic anhydride)]  
Poly[(acrylic acid):(ethyl acrylate)]  
Poly[(acrylic acid):acrylonitrile]  
ARC I-1010  
Irganox 1076  
Irganox 3114  
Irgafos 168  
Antioxidant 2246  
Erucamide  
Color masterbatch Polyethylene + 45% yellow pigment  
Color masterbatch Polystyrene + 50% brown pigment  
Color masterbatch Polyethylene + 60% brown pigment  
Color masterbatch Polyethylene + 60% Fe-oxide

---

**IS NIR Spectra – Complete Collection**

6049 NIR spectra

Copyright © NICODOM 2012

NICODOM Ltd., Hlavni 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

[ftirspectra@gmail.com](mailto:ftirspectra@gmail.com)

Color masterbatch Poly[ethylene-co-(vinyl acetate)] + 45% benzimidazolone  
Color masterbatch Polyethylene + 30% red pigment #3  
Color masterbatch Polystyrene + 30%, Sb, Ni, Ti-oxide  
Color masterbatch Poly[ethylene-co-(vinyl acetate)] + 45% yellow pigment #2  
Color masterbatch Polyethylene + 30% yellow pigment #2  
Tannic acid  
Tween 80, POE (20) sorbitan monooleate  
Caustic sulfite caramel  
Sulfite ammonia caramel  
Sodium potassium and calcium salts of fatty acids  
Ammonia caramel  
Acetylated distarch adipate  
Carnauba wax yellow  
Candelilla wax  
Sulfite ammonia caramel, in mixture  
Methyl Cellulose  
Cellulose Microcrystalline (Vivapur 105)  
Modified Starch  
Gum arabic  
Sorbitan Monolaurate, Span 20  
Starch Tapioca  
Span 40, Sorbitan Monopalmitate  
Polyglycerol Polyricinoleate  
Empiphos Alkylphosphoric acid C16-C18/5632  
Empilan KM11, Alcohol ethoxylate EO11  
Kortamol NNO  
Phthalopal LR8525, Phthalate based Resin  
Laurex CS, Cetyl Stearyl Alcohol  
Tagat CH 60, PEG-60 Hydrogenated Castor Oil  
Preventol CI-8  
Diethylene Glycol Monostearate  
Abimal 125D, Glycerol ester + colophony + maleic  
Abiester 100P, mixed ester of pentaerythritol and resin acids of colophony  
Resin AW-2, Polycyclohexanone based Resin  
Cyclohexanon Formaldehyde resin  
Alpex CK-450, Cyclized rubber  
Alnovol PN-320, Phenol Novolak  
Shellack Aged  
Shellack Angelo 4621 Superblonde  
Modocoll EK-20, Cellulose modified  
Colophony  
Resin MH  
Glyptal Resin  
Resin SK, Modified Acetophenone/Formaldehyde Resin

---

**IS NIR Spectra – Complete Collection**

6049 NIR spectra

Copyright © NICODOM 2012

NICODOM Ltd., Hlavni 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

[ftirspectra@gmail.com](mailto:ftirspectra@gmail.com)

Resin AFS  
Resin CA  
Alresat KM-40, modified maleic resin  
Worleesin MS-240, Maleic resin  
Hostaflex CN-131, Polyvinyl Acetate/Polyvinyl Chloride  
Laropal A-101, Aldehyde based Resin  
Amibal 145  
Acid Blue 93 based | Ink Blue BJTN  
Acid Blue 93 based | Ink Pink Blue RA  
Solvent Yellow 36 | Simpsol Lemon Yellow  
Solvent Blue 70 based  
Solvent Black 37  
Simrect Turquoise Blue  
Solvent Red 45  
Disperse Yellow 5 | Synten Yellow P-5G  
Disperse Yellow 1 | Synten Yellow P 2R  
Disperse Orange 30 based  
Disperse Red 5  
Synten Blue P-RL  
Disperse Blue 73 based | Synten Blue P-BGL  
Direct Yellow 87 | Folan Yellow 7G  
Acid Yellow G  
Acid Blue 83 | Brilliant Blue R  
Verde Brilliant Per Lana GB  
Nero Diretto ANBN  
Arancio Novamina PO  
Acid Dye | Scarlatto Brilliant Per Lana LC  
Ariavit Indigo Carmine  
Pigment Blue 15 based  
Pigment Blue 15:1  
Pigment Blue 15:2  
Pigment Blue 15:3  
Pigment Red 48:5 | Irgalith Red MGP 2344  
Permanent Bordo HFR DH DX 305  
Hostavinyl Yellow GRD DL CT 005  
Pigment Red 208 based | Hostaprint Red HF2B  
Direct Blue 299  
Vat Blue 6 | NoVatic Blue BC  
Vat Green 1 | NoVatic Green XBN  
Vat Green 3 | NoVatic Olive Green B  
Vat Brown 3 based | NoVatic Brown RN  
Vat Blue 6:1 | NoVatic Blue RCL  
Vat Blue 4 | Farbanthren Blue RS  
Reactive Dye | Synozol Yellow HB

---

**IS NIR Spectra – Complete Collection**

6049 NIR spectra

Copyright © NICODOM 2012

NICODOM Ltd., Hlavni 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

[ftirspectra@gmail.com](mailto:ftirspectra@gmail.com)

Reactive Dye | Synozol Red HB  
Reactive Blue 250 | Synozol N/Blue HB  
Reactive Black 5 based | Synozol Deep Black HB  
Acid Orange 67 | SynoAcid Yellow N-3RF  
Acid Red 336 | SynoAcid Red N-2RF  
Acid Blue 350 | SynoAcid Blue N-BF  
Direct Yellow 86 | Celasol Yellow NRL 150%  
Direct Red 83 based | Celasol Rubine 3BLN 200%  
Direct Red 83:1 | Celasol Rubine KC-BL  
Direct Orange 46 | Celasol Orange 7GL 143%  
Direct Red 99 | Celasol Bordo BL  
Direct Red 80 | Celasol Red 3B 200%  
Direct Blue 218 | Celasol Blue 7GL  
Direct Black 112 | Celasol Grau BGL 167%  
Direct Brown 103 | Celasol Braun BL  
Direct Red 227 | CelaDirect Brilliant Pink B  
Direct Red 23 | CelaDirect Scarlet 4BS  
Direct Yellow 142 | Inodir Yellow G  
Direct Orange 39 | Inodir Orange 2GL  
Direct Blue 71 | Inodir Blue RL  
Direct Blue 199 based | Inodir Turquoise Blue JRL  
Direct Green 26 based | Inodir Green BL  
Direct Brown 229 | Inodir Brown BRLL  
Direct Black 22 | Inodir Black VSF  
Direct Blue 293 | Inodir Blue 2G  
Direct Red 224 | Inodir Scarlet F2G  
Direct Red 243 | Inodir Red BWS  
Direct Blue 86 based | Inodir Turquoise Blue FBLL  
Direct Blue 106 | Saturn Blue LB 200  
Acid Yellow 220 | Bestalan Yellow S-2G  
Acid Orange 154 | Bestalan Orange S-2R  
Acid Red 315 | Bestalan Red S-G  
Acid Blue 317 | Bestalan Navy S-G  
Acid Green 73 | Bestalan Oliv S-G  
Acid Brown 282 | Bestalan Brown S-G  
Acid Black 207 | Bestalan Grey S-B  
Acid Black 107 based | Bestalan Black BGL  
Disperse Orange 31 | Ostacet Orange SE-2G 180  
Disperse Blue 56 based | Ostacet Blue E-LR  
Disperse Blue 79 based | Ostacet Blue S-G 200  
Disperse Yellow 23 | Ostacet Yellow E-L5R  
Disperse Yellow 42 | Ostacet Yellow SE-LG  
Disperse Orange 21 | Ostacet Orange E-R  
Disperse Red 72 | Ostacet Red S-LR

---

**IS NIR Spectra – Complete Collection**

6049 NIR spectra

Copyright © NICODOM 2012

NICODOM Ltd., Hlavni 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

[ftirspectra@gmail.com](mailto:ftirspectra@gmail.com)

Disperse Red 60 based | Ostacet Brilliant Red E-LB  
Disperse Red 121 | Ostacet Brilliant Red S-B  
Disperse Red 73 based | Ostacet Rubine SE-LB 150  
Mixture | Ostacet Red E-LG  
Disperse Red 54 | Ostacet Scarlet S-L2G  
Disperse Blue 81 | Ostacet Blue E-LG  
Mixture | Ostacet Black S-B  
Disperse Yellow 60 | Ostacet Yellow E-LR  
Mixture | Ostacet Green PG  
Ostazin Yellow Gold  
Reactive Yellow 3 | Ostazin Yellow H-A 120  
Reactive Orange 72 | Ostazin Brilliant Orange V-3G  
Reactive Violet 1 | Ostazin Brilliant Violet H-3R  
Reactive Yellow 7 | Ostazin Yellow S-GR  
Reactive Yellow 1 | Ostazin Brilliant Yellow S-6G  
Reactive Blue 4 | Ostazin Brilliant Blue S-R  
Reactive Blue 21 based | Ostazin Turquoise V-2G  
Ostazin Green H-3G  
Reactive Orange 13 | Ostazin Brilliant Orange H-2R  
Reactive Brown 2 | Ostazin Brown H-4GR  
Reactive Green 8 | Ostazin Olive H-G  
Reactive Orange 5 | Ostazin Brilliant Orange H-G  
Reactive Orange 86 | Ostazin Yellow S-3R  
Reactive Blue 109 | Ostazin Blue S-2G  
Reactive Orange 111 | Ostazin Brown S-2G  
Mixture | Ostazin Brown HG  
Reactive Black 8 | Ostazin Black H-N  
Reactive Red 2 | Ostazin Brilliant Red S-5B  
Reactive Red 24:1 based | Ostazin Brilliant Red H-B  
Reactive Blue 234 | Ostazin Navy Blue H-5R  
Diadavin UN  
Erkantol NR  
Levalin VKU-F  
Mesitol AS  
Nicca Sunsoft LM-7  
Nicca Sunsoft RF-557  
Nofome AF  
Nofome BLF 01  
Persoftal LU  
Tanawet RC-N  
Tanede RC  
Astragal PAN  
Baypret USV  
Diadavin ANE

---

**IS NIR Spectra – Complete Collection**

6049 NIR spectra

Copyright © NICODOM 2012

NICODOM Ltd., Hlavni 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

[ftirspectra@gmail.com](mailto:ftirspectra@gmail.com)

Diadavin DSP  
Diadavin EPF 01  
Eulan SPA 01, Permethrine based  
Levegal DLP  
Mesito BFN  
Perlit SI/SW  
Persoftal NPF  
Persoftal OE  
Persoftal PEK  
Plexene FED  
Protorez FFO  
Spanscour Spark  
Statexan HA  
Statexan PES  
Tanastat P  
Tanastat PH  
Tanaterge CSU Liquid  
Tanaterge DRC  
Tannex BCA-N  
Methyl-Ethyl Keton Peroxide  
Empilan BQ 100, Fatty acid ethoxylate  
Antimussol WLN  
Slovamix ZF1  
Slovamix TZ  
Lutensol A8, C12-14 fatty alcohol 8EO  
Altaran S-8  
Acid Orange 60 | Ostalan Orange RL-P  
Acid Orange 19 | Egacide Red GG  
Alcopol O-70 PG, Aqueous/PPG contains dioctyl sulfosuccinate sodium  
Empilan KB2/S, Lauryl ethoxylate (2EO)  
Empilan AMT 2, Alkylamine Ethoxylate  
Enorex Print 2030  
Dow Corning 67  
Dow Corning 30, Silicone with Propylene carbonate  
Dow Corning Z6137  
Dow Corning 66  
K-Kat 4205 (pentane 2,4-dione)  
Byk 337  
Disperbyk 193  
Reactive Yellow 167 | Sumifix Supra Brilliant Yellow 3GF 150%  
Reactive Yellow 145 based | Sumifix Supra Yellow 3RF 150%  
Reactive Red 222 | Sumifix Supra Scarlet 2GF 150%  
Reactive Red 223 | Sumifix Supra Brilliant Red GF 150%  
Reactive Red 195 | Sumifix Supra Brilliant Red 3BF 150%

---

**IS NIR Spectra – Complete Collection**

6049 NIR spectra

Copyright © NICODOM 2012

NICODOM Ltd., Hlavni 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

[ftirspectra@gmail.com](mailto:ftirspectra@gmail.com)

Reactive Dye | Sumifix Supra Red 4BNF 150%  
Reactive Blue 221 | Sumifix Supra Blue BRF 150%  
Reactive Blue 231 | Sumifix Supra Turquoise Blue BGF  
Reactive Blue 222 based | Sumifix Supra Navy Blue 3GF 150%  
Reactive Blue 222 based | Sumifix Supra Navy Blue BF  
Reactive Dye | Sumifix Supra Brown RNF  
Reactive Yellow 145 based | Sumifix Supra Yellow E-XF  
Reactive Yellow 145 based | Sumifix Supra Yellow Brown E-XF  
Reactive Dye | Sumifix Supra Rubine E-XF  
Reactive Blue 19 based | Sumifix Supra Blue E-XF  
Reactive Dye | Sumifix Supra Lemon Yellow E-XF  
Reactive Dye | Sumifix Supra Crimson E-XF  
Reactive Dye | Sumifix Supra Scarlet E-XF  
Reactive Blue 222 based | Sumifix Supra Deep Blue E-XF  
Reactive Yellow 37 | Sumifix Yellow 2GL 150%  
Reactive Yellow 15 | Sumifix Yellow GR 150%  
Reactive Orange 56 | Sumifix Brilliant Orange GRS  
Reactive Orange 16 | Sumifix Brilliant Orange 3R 135%  
Reactive Orange 107 | Sumifix Golden Yellow RNL 133%  
Reactive Red 21 | Sumifix Brilliant Red BB 150%  
Reactive Blue 19 based | Sumifix Brilliant Blue R gr.  
Reactive Blue 21 based | Sumifix Turquoise Blue G(N) conc.  
Reactive Blue 220 | Sumifix Brilliant Blue BB 133%  
Reactive Black 5 based | Sumifix Black B 150%  
Reactive Black 5 based | Sumifix Black EX(A)  
Reactive Orange 122 | Inozin Brilliant Orange 2R-FS  
Reactive Red 198 | Itofix Red VSRB 133%  
Reactive Blue 38 | Inozin Brilliant Green V-6B 90%  
Reactive Blue 203 | Itofix Navy Blue VSGG 133%  
Reactive Brown 18 | Inozin Brown V-GR  
Reactive Red 49 | Inozin Bordo V-B  
Reactive Violet 5 | Inozin Violet V-5R 200%  
Reactive Blue 89 | Inozin Dunkel Blau V-HR  
Reactive Dye | Sumifix Supra Sky Blue S-PX  
Mix C.I.Reactive Black 5 | Sumireact Jet Black S-PX  
Reactive Yellow 186 | Sumifix HF Yellow 5GL  
Reactive Dye | Sumifix HF Yellow 3R  
Reactive Dye | Sumifix HF Red G  
Reactive Dye | Sumifix HF Red 2B  
Reactive Dye | Sumifix HF Red 4B  
Reactive Dye | Sumifix HF Blue 2R  
Reactive Dye | Sumifix HF Blue BG  
Reactive Dye | Sumifix HF Dark Blue BR  
Reactive Dye | Sumifix HF Navy 2G

---

**IS NIR Spectra – Complete Collection**

6049 NIR spectra

Copyright © NICODOM 2012

NICODOM Ltd., Hlavni 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

[ftirspectra@gmail.com](mailto:ftirspectra@gmail.com)

Reactive Dye | Tulactiv Red 4B XLE  
Reactive Dye | Synozol Brilliant Yellow CP-F  
Reactive Dye | Synozol Brilliant Red CP-F  
Reactive Dye | Synozol Yellow CP-LP  
Reactive Dye | Synozol Red CP-HL  
Reactive Dye | Synozol Blue CP  
Reactive Dye | Synozol Deep Red HB  
Reactive Dye | Synozol Yellow K-HL  
Reactive Dye | Synozol Red K-HL  
Reactive Dye | Synozol Blue K-HL  
Reactive Dye | Synozol Grey K-HL  
Reactive Dye | Synozol Violet K-HL  
Reactive Dye | Synozol Yellow K-BD conc.  
Reactive Dye | Synozol Red K-BD conc.  
Reactive Dye | Synozol Red K-GD conc.  
Reactive Dye | Synozol Brilliant Orange K-R  
Reactive Dye | Synozol Brilliant Red K-3G  
Reactive Dye | Synozol Brilliant Blue K-RL conc.  
Reactive Dye | Synozol Yellow HF-4GL 150%  
Reactive Red 23 | Synozol Red 3B  
Reactive Yellow 160 | Synozol Yellow 4GL  
Vat Yellow 46 | Farbanthren Brilliant Yellow 2G  
Vat Yellow 33 | Farbanthren Yellow 3GLS  
Vat Orange 11 | Farbanthren Yellow 3RT  
Vat Orange 1 | Farbanthren Golden Yellow RK  
Vat Orange 7 | Farbanthren Brilliant Orange GR  
Vat Dye | Farbanthren Scarlet FB  
Vat Red 14 | Farbanthren Scarlet GG-N  
Vat Red 32 | Farbanthren Brilliant Red GGN  
Vat Red 1 | Farbanthren Brilliant Pink R  
Vat Red 31 | Farbanthren Red F3B  
Vat Violet 3 | Farbanthren Red Violet RRN  
Vat Blue 66 | Farbanthren Blue CLF  
Vat Blue 16 | Farbanthren Marine Blau G  
Vat Dye mixed | Farbanthren Grey CL-N  
Vat Dye mixed | Farbanthren Grey 1277  
Vat Dye mixed | Farbanthren Direct Black 1188  
Vat Orange 15 | NoVatic Orange 3G  
Vat Orange 2 | NoVatic Orange RRT  
Vat Red 10 | NoVatic Red 3B  
Vat Violet 9 | NoVatic Violet 3B  
Vat Violet 13 | NoVatic Violet XBN  
Acetyl tributyl citrate  
Versamine EH 50

---

**IS NIR Spectra – Complete Collection**

6049 NIR spectra

Copyright © NICODOM 2012

NICODOM Ltd., Hlavni 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

[ftirspectra@gmail.com](mailto:ftirspectra@gmail.com)

1,4-Diazabicyclo [2.2.2]octan  
DDBSA Acid catalyst  
Dimethyldodecyl Amine  
Dynasylan HS 2926  
Isononanoic acid  
AMP 95  
Nacure 5076  
O-ethylbenzoin  
Preventol MP 260  
Silquest A 187  
Silquest A 189  
Silquest A 1120  
Silquest A 171  
Polypox 480 (Modified cycloaliphatic polyamine)  
Polypox H 016 (Phenol-free polyamine)  
Polypox H 354 L (modified cycloaliphatic polyamine)  
Cardura E5 (glycidyl neopentanoate)  
Cyasorb UV-416 Light Absorber  
Thiabendazole  
Coleseed Oil  
Long Urethanized Alkyd  
Long Soya Bean Alkyd  
Hydroxyfunctional Alkyd from talow oil  
Tallow Oil  
Linseed Oil  
Disponil SUS 875 IC (Anionic di-2-ethylhexyl sulphosuccinate, sodium salt)  
Halox Xtain L44 (Tannin-Stain Blocker)  
Ostazin Red H-6B  
Midlon Red PRS  
Midlon Red GG  
Acid Blue 129 | Midlon Permanent Blue E  
Acid Blue | Midlon Blue B  
Acid Black 26 based | Midlon Black VL Supra  
Midlon Green P-BS  
Acid Yellow 61 | Midlon Permanent Yellow E  
Midlon Yellow PR  
Egacide Yellow MA 135%  
Acid Yellow 36 | Egacide Yellow M  
Egacide Red  
Acid Blue 41 | Egacide Blue BRL 200  
Acid Blue 78 | Egacide Blue B  
Mixture | Egacide Black A2T  
Acid Black 107 based | Ostalan Black BGL  
Acid Violet 7 based | Egacide Red 6B

---

**IS NIR Spectra – Complete Collection**

6049 NIR spectra

Copyright © NICODOM 2012

NICODOM Ltd., Hlavni 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

[ftirspectra@gmail.com](mailto:ftirspectra@gmail.com)

Food Yellow 5 based | Egacide Yellow 3G  
Acid Yellow | Egacide Yellow 3GL  
Egacide Red B  
Egacide Red 2B  
Navy Blue GGRLL for wool  
Acid Green 27 | Rybacid Green BG  
Acid Black 26 based | Rybacid Black VL  
Acid Yellow 99 | Chromolan Yellow GR 150  
Chromolan Yellow ELN  
Acid Orange 74 | Chromolan Orange G 200  
Acid Violet 58 based | Chromolan Violet R  
Mixture | Chromolan Bordo RM  
Acid Red 195 | Chromolan Pink BE 200  
Mixture | Chromolan Navy Blue R  
Acid Green 12 | Chromolan Green BL  
Acid Blue 158 | Chromolan Blue GG 200  
Mixture | Chromolan Brown G  
Acid Black 52 | Chromolan Black WA 120  
Acid Orange 173 based | Ostalan Orange RLN  
Mixture | Ostalan Brown 2GL Supra  
Acid Yellow 116 | Ostalan Yellow GRL  
Ostalan Scarlet  
Acid Orange 173 based | Ostalan Orange RL  
Acid Red 213 | Ostalan Bordo GRL  
Disperse Red | Ostalan Scarlet S  
Acid Green 43 | Ostalan Green GL  
Acid Brown 19 | Ostalan Brown BL  
Ostalan Brown 5RLN  
Ostalan Brown FS  
Mixture | Ostalan Grey 2GL  
Acid Black 60 | Ostalan Grey BL  
Solvent Black 29 based | Bricosol Black NBC  
Disperse Yellow 184:1 | Itosperse Brilliant Yellow K10GN 200%  
Disperse Orange | Itosperse Brilliant Orange FGR  
Disperse Red 277 based | Papilion Brilliant Red FFG  
Disperse Yellow | Itosperse Yellow RAP (L)  
Disperse Red 60 based | Itosperse Red RAP (L)  
Disperse Blue | Itosperse Blue RAP (L)  
Disperse Yellow 119 | Itosperse Yellow C5G 200%  
Disperse Orange 30 based | Itosperse Yellow Brown SERL 150%  
Disperse Orange 29 | Itosperse Yellow Brown R 240 %  
Disperse Red 153 | Itosperse Scarlet 2RSE 200%  
Disperse Red 152 | Itosperse Red BSE 200%  
Disperse Red 343 | Itosperse Red F3BS 150%

---

**IS NIR Spectra – Complete Collection**

6049 NIR spectra

Copyright © NICODOM 2012

NICODOM Ltd., Hlavni 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

[ftirspectra@gmail.com](mailto:ftirspectra@gmail.com)

Disperse Violet 26 | Itosperse Bordeaux 2R 160%  
Disperse Blue 366 | Itosperse Blue CRE 200%  
Disperse Blue 183 | Itosperse Blue SE2R  
Disperse Blue 73 based | Itosperse Blue SBG  
Disperse Blue 291 | Itosperse Blue C5G 150%  
Disperse Blue 60 | Itosperse Turquoise Blue CGLF 200%  
Disperse Blue 79 based | Itosperse Navy CD2G 200%  
Disperse Dye | Itosperse Navy SEFN 300%  
Disperse Brown 19 | Itosperse Yellow Brown HW  
Disperse Dye | Itosperse Crimson HW  
Disperse Red | Itosperse Red EHW  
Disperse Red 356 | Itosperse Red THW  
Disperse Red | Itosperse Rubine HW  
Disperse Blue 284 | Itosperse Blue HW  
Disperse Blue | Itosperse Navy HW  
Disperse Blue 77 | Itosperse Blue HW  
Disperse Red | Itosperse Auto Red A2G  
Disperse Blue 56 based | Itosperse Auto Blue ABL 150%  
Disperse Yellow | Itosperse Yellow XFE  
Disperse Black | Itosperse Black SERN 300%  
Disperse Orange 30 based | Itosperse Orange ECH conc.  
Disperse Red 167 | Lumacron Rubine 2GFL  
Disperse Red 73 based | Lumacron Rubine GL 200%  
Disperse Blue 165 | Lumacron Blue BBL5 200%  
Disperse Violet 77 | Sumikaron Violet S-4RL conc.  
Disperse Red 135 | Sumikaron Red S-GG 200%  
Mix C.I.Disp..Yellow 54 + Yellow | Sumikaron Yellow E-RPD  
Mix C.I.Disp..Red 60 + Red | Sumikaron Red E-RPD  
Mix C.I.Disp..Blue 56 + Violet | Sumikaron Blue E-RPD  
Disperse Yellow 192 | Sumikaron Yellow SE-RPD  
Mix C.I.Disp..Orange 30 + Yellow | Sumikaron Orange SE-RPD  
Mix C.I.Disp..Red 60 + Red | Sumikaron Red SE-RPD  
Mix C.I.Disp..Red 360 + Red | Sumikaron Rubine SE-RPD  
Mix C.I.Disp..Blue 56 + Blue | Sumikaron Blue SE-RPD  
Mix C.I.Disp..Blue 79 + Blue | Sumikaron Dark Blue SE-RPD  
Disperse Dye Mixed | Sumikaron Yellow S-RPD  
Disperse Dye Mixed | Sumikaron Red S-RPD  
Disperse Dye Mixed | Sumikaron Blue S-RPD  
Acid Black 194 | Basacid Black X 40 fl.  
Acid Blue 7 | Acid Turquoise Blue 100 %  
Acid Blue 9 based | Erioglaurine ABU 9 SF  
Acid Blue 93 based | Duasyn Ink Blue SLK  
Acid Blue 93:1 | Duasyn Ink Blue 2R  
Acid Green 16 | Acid Green W

---

**IS NIR Spectra – Complete Collection**

6049 NIR spectra

Copyright © NICODOM 2012

NICODOM Ltd., Hlavni 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

[ftirspectra@gmail.com](mailto:ftirspectra@gmail.com)

Acid Red 18 based | Basacid Red 400 E  
Acid Red 52 | Acid Rhodamine B-SF  
Acid Red 87 based | Eosine Y  
Solvent Green 7 | Pyranin 120 %  
Acid Red 51 | Erythrosine based  
Acid Violet 17 | Acid Violet 4BH-SF  
Basic Blue 7 | Basonyl Blue 636  
Basic Red 11 | Rhodamine S  
Basic Red 1:1 | Basonyl Red 485  
Doracryl Brilliant Red X-4GS 300%  
Basic Violet 1 based | Basonyl Violet 600  
Basic Yellow 37 | Basonyl Yellow 105  
Fluor Brilliant Flavine BY 4XO  
Basic Blue 26 based | Victoria Blue BB-XX6 non dedus  
Solvent Black 29 based | Neozapon Schwarz X55  
Solvent Red 160 | Neozapon Red 365  
Carmine 20%  
Basic Violet 11:1 | Basonyl Red 560  
Basic Violet 10 based | Basonyl Red 540  
Solvent Blue 44 | Savinyl Blau GLS  
Pigment Red 146 based | Renol Carmine FBB-HW  
Solvent Black 46 | Spectrasol black 46  
Acid Blue 113 | Rybacid Blue 5R 180  
Acid Red 299 | Rybacid Red  
Reactive Red 180 | Brilliant Red R-CI  
Reactive Black 31  
Solvent Red 125 | Spectrasol Red 125  
Pigment Yellow 1:1 | Monolite Yellow 2RE HD (liquid)  
Dispersol Yellow 7G PC (liquid)  
Direct Black 161 | Chlorazol Black 2G (liquid)  
Pigment Yellow 13 | Hostafine Yellow GR (liquid)  
Pigment Yellow 83 based | Hostafine Yellow HR  
Pigment Red 184 | Hostafine Rubine F6B (liquid)  
Direct Black 168 | Duasyn Direct Black HEF SF (liquid)  
Solvent Yellow 79 | Savinyl Yellow 2GLS  
Solvent Red 83:1 | Savinyl Yellow RLS  
Solvent Red 127 | Savinyl Pink 6BLS  
Solvent Red 89 | Savinyl Fire Red GLS  
Solvent Red 122 | Savinyl Red 4GLS  
Acid Yellow 5 | Basacid Yellow 099 (liquid)  
Acid Orange 172 | Levaderm Orange N (liquid)  
Acid Blue | Levaderm Blue FLN (liquid)  
Acid Brown 427 based | Levaderm Light Brown N (liquid)  
Acid Black | Levaderm Black Brown N (liquid)

---

**IS NIR Spectra – Complete Collection**

6049 NIR spectra

Copyright © NICODOM 2012

NICODOM Ltd., Hlavni 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

[ftirspectra@gmail.com](mailto:ftirspectra@gmail.com)

Acid Black | Levaderm Black N (liquid)  
Paper Light Yellow GG (liquid)  
Paper Yellow RRN (liquid)  
Paper Orange WS (liquid)  
Direct Red 255 based | Paper Scarlet 4BF  
Paper Red 4B (liquid)  
Direct Blue 14 | Paper Blue S (liquid)  
Direct Black 167 | Paper Black G (liquid)  
Paper Black PN (liquid)  
Solvent Blue 54 | Astra Blue 4R (liquid)  
Bayscript Yellow CA (liquid)  
Bayscript Cyan BG (liquid)  
Bayscript Cyan VPSP (liquid)  
Basic Violet 7 | Astrazon Red 6B (liquid)  
Basic Blue 69 | Astrazon Blue FRR (liquid)  
Basic Red 22 | Astrazon Red F3BL (liquid)  
Basic Yellow 29 | Astrazon Yellow GRL (liquid)  
Basic Blue 147 | Astrazon Blue F2RL (liquid)  
Bayderm Red A (liquid)  
Bayderm Rubine A  
Bayderm Dark Blue A (liquid)  
Pigment Orange 67 | Levanyl Orange RNLF (liquid)  
Solvent Black 48 | Morfast Black 101 (liquid)  
Solvent Blue 129 | Morfast Blue 100 (liquid)  
Solvent Brown 52 | Morfast Brown 100 (liquid)  
Solvent Red 208 | Morfast Red 102 (liquid)  
Solvent Red 68 | Morfast Red 101 (liquid)  
Solvent Yellow 143 | Morfast Yellow 101 (liquid)  
Solvent Yellow 161 | Morfast Yellow 102 (liquid)  
Solvent Red 124 | Savisol Fire Red 3 GLS (liquid)  
Solvent Yellow 62 | Savisol Yellow 2RLS liquid (liquid)  
Abeson S, dodecylbenzenesulfonic acid based  
Polyquart H81  
Acropol 35  
Pentaerythritol Tetraoleate  
Glycerol Trioleate  
PEG 400 Monooleate  
PPG Dioleate  
Isobutyl Stearate  
Empilan, Ethoxylated Propoxylated fatty alcohol  
Empilan KB-1, Ethoxylated C12-C14 1-EO  
Empicol 0216, Oleoyl etherphosphoric acid  
Empilan BQ4/S, Ethoxylated Propoxylated fatty alcohol  
Etoxon AF-5

---

**IS NIR Spectra – Complete Collection**

6049 NIR spectra

Copyright © NICODOM 2012

NICODOM Ltd., Hlavni 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

[ftirspectra@gmail.com](mailto:ftirspectra@gmail.com)

Disperbyk 194  
Byk Cera Aquacer 526  
Dispercoll U 54  
Pigment Yellow 16 | Flexonyl Yellow NCG-LA (liquid)  
Acid Violet 126 | Sanolin Violet E2R  
Acid Violet 48 based | Sanolin Violet FBL  
Acid Yellow 184 | Sanolin Flavine 8GZ  
Starch, Potato Starch Acetyl Oxidized  
Starch, Corn Starch Pregelatinised  
Sorbitan Stearate, Span 60  
Choline Tartrate  
Polyvidone 30, Polyvinylpolypyrrolidone  
Methocel E15  
Methocel E3  
Methocel E5  
Methocel K100  
Methocel K4  
Chloramphenicol Palmitate  
Carboxycellulose Calcium  
Trimecaine Hydrochloride  
Phenyl Mercury Borate  
Mercury Aminochloride  
Carboxymethyl cellulose, Na salt, ultra low viscosity  
Oil Geranium  
Phenazone with Caffeine and Acid Citric  
Poly(acrylic acid) m.w.4000000  
Cellacefate, Cellulose acetate phthalate  
Hydroxybutyl Methyl Cellulose  
Macrogol 1500, Poly(ethylene glycol) 1500  
Sodium Iodide  
Podophyllum Resin  
Fast Yellow AB  
Hypromellose (4000)  
Methylthionine Chloride  
Paraffin Wax Solid  
Magnesium Chloride Hexahydrate  
Caryophylllic Oil  
Peanut Oil  
Safflower Oil  
Starch, Rice Starch  
Explotab (Sodium Starch Glycolate)  
Povidone 25, Polyvinylpolypyrrolidone  
Povidone 90, Polyvinylpolypyrrolidone  
Polyvidone Insoluble, Polyvinylpyrrolidone

---

**IS NIR Spectra – Complete Collection**

6049 NIR spectra

Copyright © NICODOM 2012

NICODOM Ltd., Hlavni 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

[ftirspectra@gmail.com](mailto:ftirspectra@gmail.com)

Macrogol 4000, Poly(ethylene glycol) 4000  
Macrogol 6000, Poly(ethylene glycol) 6000  
Zinc Lactate  
Opaspray (mixture)  
Hard Fat  
Color masterbatch Polyethylene + 60% yellow pigment #3  
Color masterbatch Polypropylene + 70% white pigment  
Color masterbatch Polyethylene + 30% yellow pigment #4  
Color masterbatch Polyethylene + 75% white pigment  
Color masterbatch Poly(styrene-co-acrylonitrile) + 50% white pigment  
Poly(dipropylene hexachlor endomethylene tetrahydrophthalate)  
Poly(ethylene succinate)  
Poly(neopentylene succinate)  
Poly(4,4-dipropoxy-2,2-diphenyl propane fumarate)  
Poly(1,4-butylene adipate)  
Poly(dianol 22 endomethylene tetrahydrophthalate)  
Poly(trimethylpentylen endomethylene tetrahydrophthalate)  
Poly(dianol 33 endomethylene tetrahydrophthalate)  
Poly(dianol 33 fumarate)  
Poly(dipropylene isophthalate)  
Poly(propylene fumarate)  
Poly(propylene succinate)  
Poly(1,4-butylene succinate)  
Poly(dianol 22 isophthalate)  
Poly(dianol 33 phthalate)  
Poly(dianol 33 maleate)  
Poly(dianol 22 maleate)  
Poly(neopentylene fumarate)  
Poly(neopentylene adipate)  
Poly(neopentylene hexachlor endomethylene tetrahydrophthalate)  
Poly(trimethylpentylen maleate)  
Poly(dianol 22 succinate)  
Poly(diethylene tetrahydrophthalate)  
Poly(dianol 22 adipate)  
Poly(propylene endomethylene tetrahydrophthalate)  
Poly(dianol 33 isophthalate)  
Poly(1,3-butylene phthalate)  
Poly(ethylene phthalate)  
Poly(propylene hexachlor endomethylene tetrahydrophthalate)  
Poly(trimethylpentylen adipate)  
Poly(neopentylene isophthalate)  
Poly(1,3-butylene endomethylene tetrahydrophthalate)  
Poly(neopentylene terephthalate)  
Poly(ethylene endomethylene tetrahydrophthalate)

---

**IS NIR Spectra – Complete Collection**

6049 NIR spectra

Copyright © NICODOM 2012

NICODOM Ltd., Hlavni 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

[ftirspectra@gmail.com](mailto:ftirspectra@gmail.com)

Poly(1,4-butylene hexachlor endomethylene tetrahydrophthalate)  
Dicapryl Adipate  
Poly(dianol 33 adipate)  
Poly(ethylene hexachlor endomethylene tetrahydrophthalate)  
Poly(1,3-butylene adipate)  
Poly(trimethylpentylene sebacate)  
Poly(propylene isophthalate)  
Poly(neopentylene phthalate)  
Poly(dianol 33 tetrahydrophthalate)  
Poly(1,4-butylene tetrahydrophthalate)  
Poly(propylene maleate)  
Poly(1,3-butylene isophthalate)  
Color masterbatch Polystyrene + 50% white pigment  
Poly(propylene adipate)  
Poly(dianol 22 tetrahydrophthalate)  
Poly(1,3-butylene hexachlor endomethylene tetrahydrophthalate)  
Poly(diethylene hexachlor endomethylene tetrahydrophthalate)  
Poly(dianol 33 succinate)  
Poly(ethylene isophthalate)  
Streptomycin Sulfate  
Dextran 60000-80000  
Sulphur Colloidal  
Sodium Hydrogencitrate  
Cellulose Beaded (Perloza)  
Glycyl-L-Tyrosine Dihydrate  
Procaine Hydrochloride  
Neomycine Sulfate  
Vinylpyrrolidone-Vinylacetate Copolymer  
Dextran 40000  
Starch, Maize Starch Pregelatinised  
Poly(1,3-butylene fumarate)  
Poly(1,4-butylene fumarate)  
Poly(hexandiol 1,6 maleinadipate)  
Poly(dipropylene endomethylene tetrahydrophthalate)  
Poly(1,3-butylene succinate)  
Poly(diallylphthalate)

## List of Spectra /CAS Numbers

100-00-5  
100-01-6  
100-02-7  
100-06-1

---

**IS NIR Spectra – Complete Collection**

6049 NIR spectra

Copyright © NICODOM 2012

NICODOM Ltd., Hlavni 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

[ftirspectra@gmail.com](mailto:ftirspectra@gmail.com)

100-09-4  
100-10-7  
100-11-8  
10016-20-3  
100-16-3  
100-19-6  
100208-56-8, 6369-51-3, 8001-76-1  
100208-57-9  
100208-59-1, 1397-68-8, 1397-69-9, 1397-70-2, 6771-96-6  
100208-62-6  
100-21-0  
10022-31-8  
10022-68-1  
100239-44-9  
100-25-4  
10025-64-6  
10025-70-4  
10026-22-9  
10026-23-0  
10026-24-1  
10028-22-5  
10028-24-7  
1002-84-2  
10031-49-9  
10031-51-3  
1003-29-8  
1003-31-2  
10034-81-8  
10034-88-5  
10034-93-2  
10034-99-8  
10035-04-8  
10035-06-0  
100-36-7  
1003-67-4  
100-37-8  
100-39-0  
10039-54-0  
10040-45-6  
100-41-4  
100-42-5  
10042-76-9  
10043-01-3  
10043-52-4

ir-spectra.com IS NIR Spectra

---

**IS NIR Spectra – Complete Collection**

6049 NIR spectra

Copyright © NICODOM 2012

NICODOM Ltd., Hlavní 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

[ftirspectra@gmail.com](mailto:ftirspectra@gmail.com)

10043-52-4  
1004-38-2  
100-44-7  
10045-25-7  
100-46-9  
100-47-0  
100-48-1  
10048-95-0  
10048-98-3  
10049-01-1  
10049-21-5  
100-51-6  
100-52-7  
100-54-9  
100-60-7  
1006102-75-5  
100-63-0  
100641-16-5  
100-68-5  
10075-85-1  
100-79-8  
10081-67-1  
100-83-4  
100-93-6  
100-97-0  
100986-85-4  
10098-89-2  
10099-74-8  
10101-41-4  
10101-53-8  
10101-66-3  
10101-68-5  
10101-89-0  
10101-97-0  
10101-98-1  
101-02-0  
101020-29-5  
10102-06-4  
10102-17-7  
10102-18-8  
10102-40-6  
10102-45-1  
101033-44-7  
10103-46-5

ir-spectra.com IS NIR Spectra

---

**IS NIR Spectra – Complete Collection**

6049 NIR spectra

Copyright © NICODOM 2012

NICODOM Ltd., Hlavni 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

[ftirspectra@gmail.com](mailto:ftirspectra@gmail.com)

101098-23-1  
10114-58-6  
10124-27-3  
10124-36-4  
10124-37-5  
10124-53-5  
10127-27-2  
101348-81-6  
101-38-2  
10138-94-0  
10147-11-2  
1015236-20-0  
101-61-1  
101-72-4  
101-77-9  
101-81-5  
101-84-8  
101-87-1  
1018787-34-2  
10196-18-6  
101-96-2  
10199-89-0  
102-01-2  
102-02-3  
1020-31-1  
102-07-8  
102-08-9  
10213-10-2  
10217-52-4  
1021910-89-3  
10238-21-8  
10241-21-1  
102-54-5  
10262-69-8  
102-69-2  
1027-14-1  
102-71-6  
102-76-1  
102767-28-2  
10277-43-7  
102-81-8  
102-85-2  
102-86-3  
10293-06-8

ir-spectra.com IS NIR Spectra

---

**IS NIR Spectra – Complete Collection**

6049 NIR spectra

Copyright © NICODOM 2012

NICODOM Ltd., Hlavní 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

[ftirspectra@gmail.com](mailto:ftirspectra@gmail.com)

10293-10-4  
10294-26-5  
10294-39-0  
10294-40-3  
10294-41-4  
10294-42-5  
10294-46-9  
10294-54-9  
10294-56-1  
10294-66-3  
102-96-5  
10297-06-0  
102-98-7  
103-01-5  
10304-39-9  
103-11-7  
103-23-1  
10323-20-3  
10325-94-7  
10326-27-9  
10326-28-0  
103-26-4  
10326-41-7  
10328-92-4  
103-30-0  
103-33-3  
103-36-6  
10343-58-5  
103-45-7  
103-47-9  
103-49-1  
103577-45-3  
103-60-6  
103628-48-4  
103-64-0  
103-67-3  
10368-90-8  
103-69-5  
103-71-9  
103-72-0  
103-76-4  
10377-48-7  
10377-60-3  
103-81-1

ir-spectra.com IS NIR Spectra

---

**IS NIR Spectra – Complete Collection**

6049 NIR spectra

Copyright © NICODOM 2012

NICODOM Ltd., Hlavni 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

[ftirspectra@gmail.com](mailto:ftirspectra@gmail.com)

103-82-2  
103-83-3  
103-84-4  
103-85-5  
10387-40-3  
103890-78-4  
103-89-9  
103-90-2  
103-92-4  
103939-31-7  
103939-32-8  
104074-25-1  
10415-75-5  
10419-77-9  
104244-10-2  
10424-65-4  
104344-23-2  
10441-27-7  
10442-39-4  
1044571-64-3  
104-46-1  
10450-55-2  
104-54-1  
1047-16-1  
1047-16-1  
1047-16-1, 1503-48-6, 71819-74-4  
1047-25-2  
104-75-6  
104810-47-1  
104-85-8  
10486-00-7  
104-88-1  
104-90-5  
104-92-7  
104934-52-3  
104-94-9  
105-07-7  
105-08-8  
10510-54-0  
105155-66-6  
105157-03-7  
105-34-0  
10534-89-1  
105-36-2

ir-spectra.com IS NIR Spectra

---

**IS NIR Spectra – Complete Collection**

6049 NIR spectra

Copyright © NICODOM 2012

NICODOM Ltd., Hlavni 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

[ftirspectra@gmail.com](mailto:ftirspectra@gmail.com)

105-39-5  
10544-63-5  
105-50-0  
105-53-3  
105-54-4  
105-56-6  
105-58-8  
105-59-9  
105-60-2  
10563-26-5  
105-76-0  
105832-38-0  
1058-92-0  
10605-21-7  
106-11-6  
106-18-3  
106-22-9  
106266-06-2  
1062-96-0  
106-32-1  
106-33-2  
106-34-3  
106-37-6  
106-39-8  
106-40-1  
106-41-2  
106-43-4  
106-44-5  
1064-48-8  
106-46-7  
106-46-7  
106-47-8  
106-48-9  
106-49-0  
106-50-3  
106-51-4  
1066-30-4  
106-65-0  
1066-51-9  
1067-33-0  
1067-71-6  
1068-22-0  
106-87-6  
106-89-8

ir-spectra.com IS NIR Spectra

---

**IS NIR Spectra – Complete Collection**

6049 NIR spectra

Copyright © NICODOM 2012

NICODOM Ltd., Hlavni 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

[ftirspectra@gmail.com](mailto:ftirspectra@gmail.com)

106-91-2  
106-93-4  
106-94-5  
106990-43-6  
107-05-1  
107-06-2  
107-07-3  
107-11-9  
107-13-1  
107-14-2  
1071-46-1  
107-15-3  
107-18-6  
107-19-7  
107-20-0  
107-21-1  
1072-35-1  
1072-71-5  
1072-83-9  
1072845-93-2  
107-31-3  
107-35-7  
1073653-89-0  
1073-70-7  
107-41-5  
107-43-7  
1075258-11-5  
107-54-0  
107720-81-0  
107724-20-9  
1077-28-7  
107-81-3  
1078-19-9  
107868-30-4  
107-87-9  
107-88-0  
107-91-5  
107-93-7  
107-95-9  
107-96-0  
1079-66-9  
107-98-2  
108-00-9  
108-01-0

ir-spectra.com IS NIR Spectra

---

**IS NIR Spectra – Complete Collection**

6049 NIR spectra

Copyright © NICODOM 2012

NICODOM Ltd., Hlavni 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

[ftirspectra@gmail.com](mailto:ftirspectra@gmail.com)

1080-32-6  
108-05-4  
108-10-1  
108-18-9  
108-19-0  
108-20-3  
108-24-7  
108-29-2  
108-30-5  
108-31-6  
108-32-7  
1083-48-3  
108-39-4  
108-42-9  
108-43-0  
108-44-1  
108-46-3  
108-47-4  
108-48-5  
1084983-90-3  
108-59-8  
108-65-6  
108-68-9  
108-69-0  
108-73-6  
108-74-7  
108-75-8  
108-77-0  
108-78-1  
108-80-5  
108-83-8  
108-88-3  
108-89-4  
108-91-8  
108-93-0  
108-94-1  
108-95-2  
108-99-6  
109-00-2  
109-02-4  
109-04-6  
109-05-7  
109-06-8  
109-07-9

ir-spectra.com IS NIR Spectra

---

**IS NIR Spectra – Complete Collection**

6049 NIR spectra

Copyright © NICODOM 2012

NICODOM Ltd., Hlavni 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

[ftirspectra@gmail.com](mailto:ftirspectra@gmail.com)

109-12-6  
109-16-0  
1092543-29-7  
1092545-73-7  
109-30-8  
109488-50-8  
109-49-9  
109-52-4  
109-55-7  
109-57-9  
109-60-4  
109-63-7  
109-64-8  
109-65-9  
109-69-3  
109-70-6  
109-74-0  
109-76-2  
109-77-3  
109-78-4  
109-83-1  
109-85-3  
109-86-4  
109-89-7  
109-92-2  
109-94-4  
109945-04-2  
109-97-7  
109-99-9  
110-00-9  
110-02-1  
11006-55-6  
110119-84-1, 120011-70-3  
110-13-4  
110-14-5  
110-15-6  
110-16-7  
110-17-8  
110-18-9  
110-26-9  
110-27-0  
110-30-5  
110342-29-5  
110-43-0

ir-spectra.com IS NIR Spectra

---

**IS NIR Spectra – Complete Collection**

6049 NIR spectra

Copyright © NICODOM 2012

NICODOM Ltd., Hlavni 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

[ftirspectra@gmail.com](mailto:ftirspectra@gmail.com)

110-44-1  
110-45-2  
110-46-3  
110-52-1  
110553-27-0  
110-58-7  
110-64-5  
110-65-6  
11071-15-1  
110-71-4  
110-73-6  
11075-16-4  
110-80-5  
11081-40-6  
110-82-7  
110-83-8  
11084-87-0  
11084-87-0,10060-13-6  
110-85-0  
110-86-1  
110-89-4  
110-91-8  
11092-32-3  
110-93-0  
110-94-1  
11099-07-3  
11099-79-9  
11113-50-1  
11113-61-4  
11114-17-3  
111-15-9  
111-16-0  
11119-67-8  
11120-25-5  
111-20-6  
111-24-0  
111-25-1  
1112-67-0  
111-29-5  
111-34-2  
11136-36-0  
11137-99-8  
11138-66-2  
111-40-0

ir-spectra.com IS NIR Spectra

---

**IS NIR Spectra – Complete Collection**

6049 NIR spectra

Copyright © NICODOM 2012

NICODOM Ltd., Hlavní 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

[ftirspectra@gmail.com](mailto:ftirspectra@gmail.com)

111-41-1  
111-42-2  
111-44-4  
111-46-6  
111470-99-6  
111-48-8  
111-50-2  
1115-20-4  
111-65-9  
1116-76-3  
111-70-6  
111-71-7  
1117-31-3  
111-76-2  
111-78-4  
111-83-1  
111-86-4  
1118-71-4  
111-87-5  
111-90-0  
1119-46-6  
111-96-6  
111974-69-7  
111974-72-2  
1120-01-0  
1120-36-1  
112-05-0  
112068-01-6  
112-07-2  
112099-31-7  
112099-32-8  
1121-25-1  
1121-60-4  
112-18-5  
112193-83-6  
112-20-9  
112-24-3  
1122-54-9  
1122-58-3  
1122-60-7  
1122-62-9  
112-27-6  
1122-91-4  
112-30-1

ir-spectra.com IS NIR Spectra

---

**IS NIR Spectra – Complete Collection**

6049 NIR spectra

Copyright © NICODOM 2012

NICODOM Ltd., Hlavní 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

[ftirspectra@gmail.com](mailto:ftirspectra@gmail.com)

112-30-1  
112-31-2  
112-34-5  
112-36-7  
112-38-9  
112-39-0  
112-40-3  
112-41-4  
112-42-5  
1124-33-0  
112-43-6  
112-47-0  
112529-15-4  
112-53-8  
112-54-9  
1125-78-6  
112-59-4  
112-60-7  
112-61-8  
112633-26-8  
112-66-3  
112-69-6  
112-71-0  
112-73-2  
112-77-6  
112-80-1  
112809-51-5  
112-82-3  
112-84-5  
112-88-9  
112-89-0  
112-92-5  
1129-30-2  
1129-65-3  
112-96-9  
112-98-1  
113-00-8  
1131-62-0  
113-24-6  
1132-61-2  
1132819-30-7  
113299-38-0  
113299-40-4  
1134-35-6

ir-spectra.com IS NIR Spectra

---

**IS NIR Spectra – Complete Collection**

6049 NIR spectra

Copyright © NICODOM 2012

NICODOM Ltd., Hlavní 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

[ftirspectra@gmail.com](mailto:ftirspectra@gmail.com)

113-53-1  
1135-40-6  
113-59-7  
1139-30-6  
114060-55-8, 532413-17-5  
1141-59-9  
1143-38-0  
1144100-80-0  
114-49-8  
114-83-0  
114873-52-8  
1149-16-2  
114921-13-0  
114921-13-0  
1149-23-1  
1149-24-2  
114971-52-7  
115-19-5  
1153-06-6  
115340-31-3  
115-39-9  
115-40-2  
1154-12-7  
115436-72-1  
115682-09-2  
115682-10-5  
115-77-5  
115-84-4  
115857-31-3  
115-86-6  
115-95-7  
115959-99-4  
115965-77-0  
1160183-23-2  
1160361-58-9  
116-17-6  
116205-27-7  
116362-26-6  
116-54-1  
1165-91-9  
116-63-2  
1166-52-5  
116711-93-4  
116971-11-0

ir-spectra.com IS NIR Spectra

---

**IS NIR Spectra – Complete Collection**

6049 NIR spectra

Copyright © NICODOM 2012

NICODOM Ltd., Hlavni 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

[ftirspectra@gmail.com](mailto:ftirspectra@gmail.com)

117-08-8  
1173710-58-1  
1174-72-7  
117-59-9  
1177-87-3  
117-80-6  
117-81-7  
118-00-3  
118-02-5  
118-10-5  
1182-07-6  
1182-66-7  
118-29-6  
1183-04-6  
118-32-1  
118337-09-0  
1185-53-1  
1185-57-5  
118-55-8  
118-57-0  
1185-81-5  
1186-10-3  
118-61-6  
118-71-8  
118-75-2  
118-79-6  
1188-02-9  
118-82-1  
118-90-1  
118-91-2  
118-92-3  
118-93-4  
118-96-7  
1189763-33-4  
119-15-3  
119-15-3  
1192-62-7  
119-26-6  
119-28-8  
1193-62-0  
119-36-8  
119392-95-9  
119445-90-8  
119-47-1

ir-spectra.com IS NIR Spectra

---

**IS NIR Spectra – Complete Collection**

6049 NIR spectra

Copyright © NICODOM 2012

NICODOM Ltd., Hlavní 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

[ftirspectra@gmail.com](mailto:ftirspectra@gmail.com)

119-52-8  
119-53-9  
1195563-23-5  
1195-59-1  
119-61-9  
119-64-2  
119-65-3  
1197-18-8  
119-79-9  
1198-77-2  
119-90-4  
119-91-5  
119930-52-8  
119-93-7  
1199-95-7  
120011-70-3  
12001-26-2  
120014-06-4  
12001-98-8  
12001-99-9  
120060-80-2  
120-07-0  
12011-78-8  
120-12-7  
120-14-9  
120-23-0  
12026-66-3  
12027-06-4  
12027-06-4  
12027-67-7  
12030-88-5  
120-32-1  
12035-88-0  
12037-29-5  
120-46-7  
120-47-8  
120-51-4  
12054-48-7  
12054-85-2  
12055-62-8  
120-55-8  
1205-62-5  
120-57-0  
120-58-1

ir-spectra.com IS NIR Spectra

---

**IS NIR Spectra – Complete Collection**

6049 NIR spectra

Copyright © NICODOM 2012

NICODOM Ltd., Hlavni 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

[ftirspectra@gmail.com](mailto:ftirspectra@gmail.com)

120-61-6  
12062-24-7  
12064-62-9  
120-66-1  
12067-99-1  
12068-03-0  
12069-69-1  
120-72-9  
120-78-5  
12079-65-1  
12080-32-9  
120-80-9  
12081-18-4  
120-83-2  
12086-40-7  
12088-57-2  
120919-08-6  
120-92-3  
12092-47-6  
12093-10-6  
120-93-4  
12098-17-8  
121-00-6  
1210-12-4  
1210-35-1  
12107-56-1  
121-14-2  
12124-97-9  
12125-01-8  
12125-02-9  
12125-06-3  
121-32-4  
121-34-6  
1214-39-7  
121-44-8  
12145-47-0  
121-45-9  
121-47-1  
12150-46-8  
121-54-0  
121-57-3  
121-69-7  
12173-47-6  
1217358-08-1

ir-spectra.com IS NIR Spectra

---

**IS NIR Spectra – Complete Collection**

6049 NIR spectra

Copyright © NICODOM 2012

NICODOM Ltd., Hlavní 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

[ftirspectra@gmail.com](mailto:ftirspectra@gmail.com)

1217-45-4  
121-79-9  
121-82-4  
121-87-9  
121888-67-3  
121888-68-4  
121-89-1  
121-92-6  
12209-98-2  
12210-47-8  
122-11-2  
12217-35-5, 83027-46-7  
12217-36-6  
12217-38-8  
12217-48-0, 65122-06-7  
12217-50-4  
12217-52-6  
12217-79-7  
12217-80-0  
12217-83-3, 10000-42-7  
12217-95-7  
12218-95-0  
12218-96-1  
12218-96-1  
12219-20-4  
12219-26-0  
12219-31-7  
12219-35-1  
12219-65-7  
12219-88-4  
12219-93-1  
12220-24-5  
12220-29-0  
12220-51-8  
12220-51-8  
12220-6-3  
12221-36-2  
12221-61-3  
12222-04-7  
12222-48-9  
12222-51-4  
12222-78-5  
12222-78-5, 15114-15-5, 13716-91-1, 13698-89-0  
12222-79-6

---

**IS NIR Spectra – Complete Collection**

6049 NIR spectra

Copyright © NICODOM 2012

NICODOM Ltd., Hlavni 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

[ftirspectra@gmail.com](mailto:ftirspectra@gmail.com)

12223-23-3  
12223-23-3  
12223-39-1  
12223-62-0  
12225-06-8  
12225-26-2  
12225-68-2  
12225-73-9  
12225-77-3, 70210-47-8  
12225-88-6  
12226-38-9  
12226-46-9  
12226-47-0  
12226-74-3  
12226-77-6, 58104-34-0  
12226-78-7  
12226-82-3, 28198-05-2  
12227-04-2, 97392-94-4, 71077-14-0  
12227-47-3  
12227-50-8  
12227-55-3  
12227-60-0  
12227-67-7  
12227-89-3  
122-28-1  
12230-71-6  
122-31-6  
122-32-7  
12232-99-4  
12234-64-9  
12234-73-0  
12235-16-4  
12235-47-1  
12236-62-3  
12236-86-1  
12236-86-1  
12236-90-7  
12236-92-9  
12236-93-0, 70210-17-2  
12237-02-4  
12237-16-0  
12237-22-8  
12237-24-0  
12237-24-0

---

**IS NIR Spectra – Complete Collection**

6049 NIR spectra

Copyright © NICODOM 2012

NICODOM Ltd., Hlavní 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

[ftirspectra@gmail.com](mailto:ftirspectra@gmail.com)

12237-26-2  
12237-30-8  
12237-31-9  
12237-50-2  
12239-11-1  
12239-15-5, 69762-08-9  
12239-18-8  
12239-29-1, 70306-36-8  
12239-34-8  
122-39-4  
12239-45-1  
12239-58-6, 27425-55-4  
12239-74-6  
12239-75-7  
12240-03-8  
12240-15-2, 14038-43-8, 25869-00-5  
12244-51-8  
122-51-0  
122-52-1  
122555-91-3  
122560-63-8, 9002-24-8  
122-57-6  
122-62-3  
12269-90-8  
12270-13-2  
12270-30-3  
12271-00-0  
122-75-8  
122-80-5  
122-85-0  
12286-66-7  
1229-55-6  
1229-55-6  
122-96-3  
122-98-5  
122-99-6  
123-03-5  
123-04-6  
123-08-0  
123-11-5  
123-19-3  
123-25-1  
123-30-8  
123-31-9

ir-spectra.com IS NIR Spectra

---

**IS NIR Spectra – Complete Collection**

6049 NIR spectra

Copyright © NICODOM 2012

NICODOM Ltd., Hlavni 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

[ftirspectra@gmail.com](mailto:ftirspectra@gmail.com)

123-33-1  
123333-66-4  
123-38-6  
123-41-1  
123-42-2  
123-51-3  
123-54-6  
12354-85-7  
123-56-8  
123-62-6  
123-63-7  
123-66-0  
123-72-8  
123-76-2  
123-79-5  
123-86-4  
123-91-1  
123-92-2  
12392-64-2  
123-94-4  
1239-45-8  
123-96-6  
123968-25-2  
123-99-9  
124-04-9  
124-09-4  
124-13-0  
124172-53-8  
124-17-4  
124-18-5  
1242003-07-1  
124-22-1  
124-26-5  
124-30-1  
12439-96-2  
1244956-73-7  
124529-02-8  
1246271-34-0  
124-63-0  
124-68-5  
124750-99-8  
124755-24-4  
1248-42-6  
124944-40-7

ir-spectra.com IS NIR Spectra

---

**IS NIR Spectra – Complete Collection**

6049 NIR spectra

Copyright © NICODOM 2012

NICODOM Ltd., Hlavní 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

[ftirspectra@gmail.com](mailto:ftirspectra@gmail.com)

125-20-2  
1252679-07-4  
1253-84-5  
125-60-0  
125700-67-6  
125878-86-6  
1260-17-9  
12602-23-2  
12609-80-2  
126-11-4  
126-14-7  
126161-18-0  
126235-76-5  
12626-73-2  
126-30-7  
126-33-0  
12671-74-8, 27870-92-4  
126-73-8  
12674-05-4  
126-81-8  
126-86-3  
127-06-0  
127-08-2  
127-09-3  
1271-19-8  
1271-42-7  
1271-47-2  
1271-55-2  
12715-60-5  
1271-56-3  
127-19-5  
12731-63-4  
1273-86-5  
1273-94-5  
127-47-9  
127-52-6  
127-63-9  
127-65-1  
127-68-4  
12768-78-4  
12769-07-2  
12769-96-9  
127-91-3  
127-95-7

ir-spectra.com IS NIR Spectra

---

**IS NIR Spectra – Complete Collection**

6049 NIR spectra

Copyright © NICODOM 2012

NICODOM Ltd., Hlavní 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

[ftirspectra@gmail.com](mailto:ftirspectra@gmail.com)

128-08-5  
128-09-6  
128-13-2  
128-37-0  
128-39-2  
128416-19-3  
128-42-7  
128-44-9  
128513-84-8  
128-58-5  
128-69-8  
1287-17-8  
128-80-3  
128-95-0  
128-97-2  
129-09-9  
129101-54-8  
129113-00-4  
1291-32-3  
129-17-9  
129-17-9, 20262-76-4  
129318-43-0  
1293-73-8  
1293-87-4  
129423-54-7  
129-44-2  
129-66-8  
129722-12-9  
129722-25-4  
129893-27-2  
129898-77-7  
130018-87-0  
13007-92-6  
130-13-2  
130-20-1  
130-22-3  
1302-78-9  
13031-60-2  
1303-33-9  
1303-86-2  
130390-90-8  
1303-96-4  
1304-85-4  
130500-95-7

ir-spectra.com IS NIR Spectra

---

**IS NIR Spectra – Complete Collection**

6049 NIR spectra

Copyright © NICODOM 2012

NICODOM Ltd., Hlavni 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

[ftirspectra@gmail.com](mailto:ftirspectra@gmail.com)

1305-62-0  
13057-17-5  
1305-78-8  
1306-23-6  
1308-14-1  
1308-38-9, 68909-79-5  
1308-38-9, 68909-79-5  
130-89-2  
130-95-0  
1309-64-4  
130995-13-0  
1310-58-3  
1310-65-2  
1310-66-3  
1310-73-2  
131-08-8  
1311-10-0  
131-11-3  
131-22-6  
13126-12-0  
131-27-1  
1312-81-8  
1313-60-6  
13138-36-8  
13139-14-5  
13139-15-6  
1315-04-4  
131-53-3  
131-55-5  
131-56-6  
131-57-7  
131594-92-8  
131640-49-8, 58128-18-0, 28316-43-0  
1317-36-8  
131-73-7  
1317-38-0  
131-74-8  
1317-60-8  
1318-02-1  
13185-73-4  
131-91-9  
131-92-0  
1319-46-6  
1320-06-5

ir-spectra.com IS NIR Spectra

---

**IS NIR Spectra – Complete Collection**

6049 NIR spectra

Copyright © NICODOM 2012

NICODOM Ltd., Hlavní 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

[ftirspectra@gmail.com](mailto:ftirspectra@gmail.com)

13204-98-3  
1321-14-8  
1321-73-9  
13236-02-7  
132396-99-7  
1324-04-5  
1324-11-4  
1324-17-0  
1324-35-2  
1324-53-4  
132-53-6  
132539-06-1  
1325-54-8  
1325-62-8  
1325-63-9  
13257-51-7  
132679-61-9  
13268-42-3  
132723-93-4  
1327-41-9  
1328-51-4  
1328-51-4  
1328-53-6  
1328-59-2, 8005-33-2, 90294-88-5, 475-25-2  
132-86-5  
132893-47-1  
13291-61-7  
13294-71-8  
1330-20-7  
1330-38-7  
1330-38-7  
13303-98-5  
1330-78-5  
133-10-8  
133-13-1  
1332-09-8  
13324-20-4, 4499-01-8  
1332-58-7  
1332-81-6  
133-32-4  
1333-82-0  
133404-94-1, 11140-65-1  
1334-76-5  
1335-12-2

---

**IS NIR Spectra – Complete Collection**

6049 NIR spectra

Copyright © NICODOM 2012

NICODOM Ltd., Hlavni 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

[ftirspectra@gmail.com](mailto:ftirspectra@gmail.com)

133556-24-8  
1338-23-4  
1338-39-2  
1338-41-6  
1338-43-8  
13395-16-9  
1339-83-9  
1339-83-9  
134-03-2  
13403-73-1  
134050-26-3  
1341-49-7  
13419-61-9  
134-31-6  
13438-45-4  
1343-98-2  
1343-98-2  
13444-71-8  
13445-51-7  
13446-18-9  
13446-19-0  
13446-28-1  
13446-72-5  
1344-95-2  
1344-98-5  
1345-00-2, 68186-90-3, 57455-37-5  
1345-05-7  
134-50-9  
1345-16-0  
1345-19-3, 68187-05-3, 83712-59-8  
13453-78-6  
13454-73-4  
13455-36-2  
134595-59-8  
13462-86-7  
134633-09-3  
13463-40-6  
13463-67-7  
13463-67-7  
13463-67-7  
13465-95-7  
134687-49-3  
134695-74-2  
13477-34-4

ir-spectra.com IS NIR Spectra

---

**IS NIR Spectra – Complete Collection**

6049 NIR spectra

Copyright © NICODOM 2012

NICODOM Ltd., Hlavní 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

[ftirspectra@gmail.com](mailto:ftirspectra@gmail.com)

13478-33-6  
134-81-6  
134910-54-6  
13494-90-1  
13494-98-9  
135-02-4  
13510-49-1  
135-19-3  
13520-56-4  
13520-61-1  
13520-83-7  
13521-83-0  
13528-93-3  
13530-65-9  
13536-84-0  
13537-15-0  
13548-38-4  
13550-47-5  
135-57-9  
13560-49-1  
13573-16-5  
13586-38-4  
135-88-6  
13590-82-4  
13596-22-0  
13597-54-1  
13598-36-2  
13600-88-9  
13600-98-1  
13601-11-1  
13613-65-5  
136-35-6  
136434-34-9  
136-47-0  
136496-87-2  
13668-61-6  
13681-82-8  
13689-92-4  
136-94-7  
136-95-8  
13701-59-2  
13701-64-9  
137-07-5  
137-08-6

ir-spectra.com IS NIR Spectra

---

**IS NIR Spectra – Complete Collection**

6049 NIR spectra

Copyright © NICODOM 2012

NICODOM Ltd., Hlavni 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

[ftirspectra@gmail.com](mailto:ftirspectra@gmail.com)

13708-85-5  
137-09-7  
13710-19-5  
137-16-6  
13717-00-5  
13718-59-7  
13721-39-6  
137-29-1  
137-30-4  
13731-82-3  
137-32-6  
13732-62-2  
13734-41-3  
137-40-6  
137-43-9  
13746-66-2  
13746-89-9  
137506-60-6  
13755-29-8  
13755-38-9  
13759-83-6  
13765-49-6  
137-66-6  
13768-11-1  
13770-16-6  
13773-69-8  
13776-84-6  
13782-01-9  
138067-74-0  
138-09-0  
13814-33-0  
13814-62-5  
13814-75-0  
138149-64-1  
138-15-8  
138184-36-8  
138-24-9  
13825-74-6  
138265-88-0  
13826-66-9  
13826-83-0  
138-32-9  
13838-08-9  
138402-11-6

ir-spectra.com IS NIR Spectra

---

**IS NIR Spectra – Complete Collection**

6049 NIR spectra

Copyright © NICODOM 2012

NICODOM Ltd., Hlavní 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

[ftirspectra@gmail.com](mailto:ftirspectra@gmail.com)

138-52-3  
13859-51-3  
138-60-3  
138729-47-2  
13881-91-9  
138-86-3  
138-89-6  
13894-21-8  
138988-65-5  
139-05-9  
1390-65-4  
1390-65-4, 1260-17-9  
139-12-8  
139-12-8  
139-13-9  
139262-23-0  
13927-77-0  
139-33-3  
1393-63-1  
1393-63-1, 6983-79-5  
13939-06-5  
13940-38-0  
13943-58-3  
1395-35-6  
139-59-3  
139-60-6  
13963-57-0  
13963-58-1  
13965-03-2  
13967-50-5  
139691-02-4  
13983-17-0  
139-85-5  
140-10-3  
140-11-4  
1401-55-4  
1401-55-4  
140-22-7  
14024-18-1  
14024-48-7  
14024-61-4  
14024-63-6  
140-29-4  
140-31-8

ir-spectra.com IS NIR Spectra

---

**IS NIR Spectra – Complete Collection**

6049 NIR spectra

Copyright © NICODOM 2012

NICODOM Ltd., Hlavni 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

[ftirspectra@gmail.com](mailto:ftirspectra@gmail.com)

14040-11-0  
14047-23-5  
1405-10-3  
1405-87-4  
1405-89-6  
14059-33-7  
140-66-9  
140679-55-6  
140-72-7  
14074-80-7  
140-87-4  
140-89-6  
140-93-2  
14097-03-1  
141-01-5  
141093-26-7  
14114-05-7  
141-28-6  
141-32-2  
141-43-5  
14149-58-7  
14152-97-7  
141-53-7  
141-78-6  
141-79-7  
141-82-2  
141-84-4  
141-86-6  
14187-32-7  
141-90-2  
141-95-7  
141-97-9  
142-04-1  
14205-39-1  
142106-29-4  
14217-01-7  
14220-21-4  
14221-01-3  
142227-49-4  
142227-51-8  
142-29-0  
14233-37-5  
142-47-2  
142-62-1

ir-spectra.com IS NIR Spectra

---

**IS NIR Spectra – Complete Collection**

6049 NIR spectra

Copyright © NICODOM 2012

NICODOM Ltd., Hlavní 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

[ftirspectra@gmail.com](mailto:ftirspectra@gmail.com)

14264-16-5  
142-64-3  
14267-17-5  
142-71-2  
142-73-4  
142-82-5  
14283-07-9  
14285-68-8  
142902-13-4  
14295-43-3  
142-96-1  
142986-44-5  
14302-13-7  
143-07-7  
143-08-8  
143-15-7  
143179-86-6  
143-19-1  
14320-04-8  
143-22-6  
14324-55-1  
143-24-8  
143-28-2  
14348-75-5  
143-66-8  
143687-46-1  
14371-10-9  
143-74-8  
143860-04-2  
1439-41-4  
143-98-6  
144-02-5  
14402-75-6  
144-19-4  
144-23-0  
144-23-0  
144246-02-6  
14426-21-2  
14431-43-7  
144-33-2  
14434-22-1  
144-48-9  
144-55-8  
1445-73-4

ir-spectra.com IS NIR Spectra

---

**IS NIR Spectra – Complete Collection**

6049 NIR spectra

Copyright © NICODOM 2012

NICODOM Ltd., Hlavní 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

[ftirspectra@gmail.com](mailto:ftirspectra@gmail.com)

144598-75-4  
14459-95-1  
144701-48-4  
144-82-1  
14483-18-2  
145017-98-7  
145-13-1  
14516-47-3  
14516-54-2  
145483-63-2  
1454-85-9  
1455-77-2  
14590-13-7  
1460-02-2  
1460-05-5  
1460-57-7  
146062-45-5  
1461-15-0  
1461-22-9  
1461-23-0  
146-17-8  
14634-93-6  
1465-25-4  
14660-52-7  
146-68-9  
1468-83-3  
14694-95-2  
147098-20-2  
147-14-8  
147-14-8  
147-14-8  
147-14-8  
147-14-8  
14726-36-4  
147315-50-2  
147-47-7  
147-58-0  
147-71-7  
1477-55-0  
147-82-0  
147826-71-9  
147-85-3  
147-93-3  
14807-96-6

ir-spectra.com IS NIR Spectra

---

**IS NIR Spectra – Complete Collection**

6049 NIR spectra

Copyright © NICODOM 2012

NICODOM Ltd., Hlavni 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

[ftirspectra@gmail.com](mailto:ftirspectra@gmail.com)

14808-60-7  
14808-60-7  
14808-60-7  
148-18-5  
148-24-3  
148-25-4  
14836-60-3  
1484-80-6  
148618-29-5  
148619-01-6  
14866-33-2  
14871-79-5  
14882-18-9  
14918-21-9  
14918-39-9  
149-30-4  
1493-13-6  
149-32-6  
1493-27-2  
14937-45-2  
149-44-0  
149-45-1  
1499-10-1  
149-91-7  
15000-59-6  
150-13-0  
1501-82-2  
150-30-1  
1503-48-6, 71819-76-6  
15060-62-5  
150-76-5  
150-86-7  
150-90-3  
15120-17-9  
151-21-3  
15131-43-8  
151-56-4  
151616-42-1  
15168-20-4  
1518-58-7  
1518-86-1  
15205-66-0  
152-11-4  
15214-89-8

ir-spectra.com IS NIR Spectra

---

**IS NIR Spectra – Complete Collection**

6049 NIR spectra

Copyright © NICODOM 2012

NICODOM Ltd., Hlavni 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

[ftirspectra@gmail.com](mailto:ftirspectra@gmail.com)

152165-67-8  
1523-11-1  
15233-47-3  
15243-33-1  
15244-37-8  
152459-95-5  
15280-53-2  
1529-68-6  
1530-32-1  
15307-79-6  
15307-81-0  
153-18-4  
15321-51-4  
15356-70-4  
15364-94-0  
153-78-6  
153-94-6  
15418-16-3  
154-23-4  
154336-20-6  
15435-71-9  
154361-50-9  
15475-27-1  
155-04-4  
15507-52-5  
155141-29-0  
1552-12-1  
15521-65-0  
15530-47-9  
15546-12-0  
155684-37-0  
15571-58-1  
155773-72-1  
155775-31-8  
15593-52-9  
15593-61-0  
15598-34-2  
156028-40-9  
1560-28-7  
156260-01-4  
1562-85-2  
15629-92-2  
156-43-4  
1564-64-3

ir-spectra.com IS NIR Spectra

---

**IS NIR Spectra – Complete Collection**

6049 NIR spectra

Copyright © NICODOM 2012

NICODOM Ltd., Hlavní 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

[ftirspectra@gmail.com](mailto:ftirspectra@gmail.com)

156-54-7  
156-57-0  
15658-60-3  
156-59-2  
15676-16-1  
15680-42-9  
15687-27-1  
156-87-6  
1569-69-3  
156970-79-5  
15702-36-0  
1570-45-2  
157-06-2  
15708-48-2  
15709-76-9  
1571-08-0  
1571-33-1  
157185-56-3  
15727-65-8  
15761-38-3  
15761-39-4  
1576-35-8  
15782-05-5  
157856-80-9  
15790-07-5  
15790-15-5  
158129-68-1  
158226-29-0  
15827-60-8  
15878-77-0  
15879-93-3  
15901-08-3  
15901-11-8  
1592-23-0  
15988-96-2  
1603-40-3  
1603-79-8  
16039-53-5  
16039-53-5  
16039-64-8  
16045-72-0  
16051-77-7  
16069-36-6  
1608-26-0

ir-spectra.com IS NIR Spectra

---

**IS NIR Spectra – Complete Collection**

6049 NIR spectra

Copyright © NICODOM 2012

NICODOM Ltd., Hlavni 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

[ftirspectra@gmail.com](mailto:ftirspectra@gmail.com)

16091-04-6  
16096-31-4  
1611-35-4  
16115-80-3  
16116-78-2  
16182-70-0  
162148-48-3  
162157-03-1  
1622-62-4  
162359-56-0  
16245-77-5  
16283-36-6  
16294-75-0  
1631-83-0  
1633-05-2  
163442-67-9  
1634-82-8  
163515-54-6  
16357-59-8  
1637-73-6  
1641-09-4  
16419-60-6  
16423-68-0  
16423-68-0  
1643-19-2  
16433-88-8  
16433-96-8  
16454-60-7  
164578-37-4  
164858-92-8  
16545-54-3  
1657-16-5  
16577-73-4  
1658-56-6, 18268-54-7freeacid  
16588-74-2  
16593-81-0  
1660-94-2  
1662-01-7  
1663-45-2  
1667-01-2  
16674-78-5  
166798-98-7  
1668-10-6  
16695-22-0

ir-spectra.com IS NIR Spectra

---

**IS NIR Spectra – Complete Collection**

6049 NIR spectra

Copyright © NICODOM 2012

NICODOM Ltd., Hlavni 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

[ftirspectra@gmail.com](mailto:ftirspectra@gmail.com)

1670-14-0  
16731-55-8  
1675-54-3  
1675-54-3  
1675-54-3  
1675-54-3  
1675-54-3  
16773-42-5  
16774-21-3  
1679-64-7  
16836-95-6  
16858-02-9  
16867-03-1  
16872-11-0  
1687-30-5  
16889-10-4  
16889-10-4  
1689-82-3  
16909-78-7  
16940-66-2  
16940-66-2  
1694-09-3  
16941-11-0  
16949-15-8  
169557-48-6  
169830-46-0  
170161-27-0  
17026-42-5  
17032-40-5  
1707-03-5  
17084-13-8  
17095-24-8  
17095-24-8  
17095-24-8  
1709-70-2  
171290-52-1  
1714-29-0  
17146-08-6  
171596-29-5  
171669-30-0  
17176-77-1  
1719-71-7  
17199-29-0  
172201-26-2

ir-spectra.com IS NIR Spectra

---

**IS NIR Spectra – Complete Collection**

6049 NIR spectra

Copyright © NICODOM 2012

NICODOM Ltd., Hlavni 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

[ftirspectra@gmail.com](mailto:ftirspectra@gmail.com)

172399-10-9  
17261-34-6  
172672-69-4  
17282-04-1  
17309-53-4  
1733-12-6  
173346-69-5  
1736-70-5  
17372-87-1  
1738-02-9  
1738-25-6  
174137-97-4  
17418-58-5  
17418-58-5, 54243-60-6  
174264-50-7  
1742-95-6  
17455-13-9  
174600-19-2  
17464-88-9  
17471-45-3  
174820-49-6  
1750-12-5  
17501-44-9  
17557-23-2  
17560-51-9  
175705-31-4  
1758-68-5  
17593-73-6  
1759-53-1  
1760-24-3  
1761-56-4  
1761-61-1  
17629-30-0  
1762-95-4  
17636-10-1  
17637-11-5  
17696-69-4  
177256-40-5  
17741-63-8  
17741-63-8  
1774-47-6  
177937-60-9  
1779-48-2  
1779-51-7

ir-spectra.com IS NIR Spectra

---

**IS NIR Spectra – Complete Collection**

6049 NIR spectra

Copyright © NICODOM 2012

NICODOM Ltd., Hlavni 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

[ftirspectra@gmail.com](mailto:ftirspectra@gmail.com)

178037-57-5  
17804-49-8  
1783-96-6  
1787-61-7  
17927-65-0  
18014-00-1  
180200-66-2  
1806-34-4  
1809-20-7  
18108-24-2  
18109-81-4  
18156-74-6  
18162-48-6  
18163-47-8  
18173-64-3  
1820-81-1  
182507-22-8  
1829-00-1  
1830-54-2  
1837-57-6  
18419-53-9  
1843-03-4  
1843-05-6  
1844-09-3  
18472-51-0  
1847-55-8  
18496-40-7  
18497-13-7  
18531-99-2  
18624-44-7  
186286-86-8  
1863-63-4  
186397-64-8  
1871-22-3  
1878-68-8  
1879-09-0  
18851-33-7  
1885-14-9  
1885-29-6  
18917-93-6  
1892-57-5  
189441-17-6  
189446-58-0  
1896-62-4

ir-spectra.com IS NIR Spectra

---

**IS NIR Spectra – Complete Collection**

6049 NIR spectra

Copyright © NICODOM 2012

NICODOM Ltd., Hlavní 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

[ftirspectra@gmail.com](mailto:ftirspectra@gmail.com)

18987-59-2  
18996-35-5  
19034-13-0  
191217-81-9  
1912-28-3  
1913-93-5  
1914-99-4  
19201-53-7  
19222-41-4  
19230-28-5  
1923-70-2  
19237-84-4  
1928-81-0  
19321-40-5  
1934-20-9  
1934-21-0  
1937-34-4  
19381-50-1  
1941-27-1  
1944-12-3  
1945-84-2  
194800-56-1  
194800-57-2  
1948-33-0  
19485-76-8  
19513-05-4  
19598-90-4  
197252-55-4  
19780-11-1  
19794-93-5  
19816-88-7  
1984-06-1  
1984-15-2  
1989-52-2  
198964-46-4  
1990-11-0  
199015-36-6  
199620-18-3  
19992-12-2  
2001-45-8  
20017-83-8  
200630-46-2  
20123-80-2  
201363-83-9

ir-spectra.com IS NIR Spectra

---

**IS NIR Spectra – Complete Collection**

6049 NIR spectra

Copyright © NICODOM 2012

NICODOM Ltd., Hlavní 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

[ftirspectra@gmail.com](mailto:ftirspectra@gmail.com)

2016-42-4  
20213-65-4  
20227-30-9  
20241-76-3  
202483-55-4  
202875-89-6  
2032-35-1  
2033-24-1  
20344-49-4  
20349-39-7  
204143-01-1  
2042-37-7  
2043-61-0  
2050-76-2  
20526-58-3  
20624-25-3  
206-44-0  
20664-98-6  
20694-39-7  
20702-77-6  
2071-20-7  
207386-92-3  
207683-19-0  
20768-67-6  
2078-12-8#  
2078-54-8  
208252-54-4  
2082-79-3  
2082-81-7  
20830-75-5  
208350-19-0  
20837-86-9  
2092-17-3  
21041-95-2  
21083-47-6  
211629-88-8  
2117-69-3  
211811-37-9  
21187-98-4  
212516-19-3  
212652-59-0  
21286-54-4  
214078-95-2  
214360-73-3

ir-spectra.com IS NIR Spectra

---

**IS NIR Spectra – Complete Collection**

6049 NIR spectra

Copyright © NICODOM 2012

NICODOM Ltd., Hlavni 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

[ftirspectra@gmail.com](mailto:ftirspectra@gmail.com)

2150-60-9  
215462-83-2  
21564-17-0  
215929-30-9  
2163-48-6  
21679-31-2  
2172-33-0  
2175-86-2  
2175-90-8  
217-59-4  
217797-14-3  
217819-09-5  
2181-42-2  
21829-25-4  
2185-86-6  
2192-20-3  
21959-36-4  
2197-63-9  
219861-08-2  
2199-58-8  
220182-19-4  
220617-96-9  
22071-15-4  
2210-79-9  
22144-60-1  
2215-20-5  
221628-41-7  
2216-51-5  
2219-31-0  
22204-53-1  
22204-88-2  
22206-57-1  
2223-93-0  
2223-95-2  
2226-96-2  
2227-79-4  
223251-74-9  
223251-86-3  
22326-55-2  
223446-13-7  
223745-37-7  
22396-48-1  
2243-76-7  
2244-11-3#

ir-spectra.com IS NIR Spectra

---

**IS NIR Spectra – Complete Collection**

6049 NIR spectra

Copyright © NICODOM 2012

NICODOM Ltd., Hlavni 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

[ftirspectra@gmail.com](mailto:ftirspectra@gmail.com)

2246-46-0  
22483-09-6  
22664-55-7  
226882-93-5  
226923-51-9  
22722-98-1  
2273-43-0  
22743-80-2  
22767-50-6  
2280-01-5#  
229485-33-0  
2295-31-0  
229-87-8  
229959-20-0  
230-46-6  
2309-94-6  
23141-83-5  
23165-29-9  
2321-07-5  
2321-07-5  
2328-96-3  
2338-05-8  
23389-75-5  
23539-10-8  
2362-14-3  
23624-49-9  
23647-14-5  
23647-14-5  
23767-99-9  
2379-74-0  
2379-75-1  
2379-77-3  
2379-78-4  
2379-79-5  
23828-92-4  
23836-43-3  
23850-94-4  
2387-23-7  
238756-91-7  
2390-60-5  
23906-97-0  
23911-80-0  
2391-30-2  
2396-63-6

ir-spectra.com IS NIR Spectra

---

**IS NIR Spectra – Complete Collection**

6049 NIR spectra

Copyright © NICODOM 2012

NICODOM Ltd., Hlavni 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

[ftirspectra@gmail.com](mailto:ftirspectra@gmail.com)

2396-65-8  
2403-88-5  
2403-89-6  
24057-28-1  
2409-55-4  
2411-89-4  
2425-79-8  
2425-85-6  
24280-93-1  
24287-95-4  
2438-80-4  
2439-99-8  
2440-22-4  
24424-99-5  
2444-46-4  
2444-90-8  
24463-19-2  
2451-01-6  
24516-39-0  
24544-04-5  
2456-81-7  
2461-15-6  
24623-77-6  
2465-27-2  
2465-93-2  
2469-45-6  
2475-31-2  
2476-37-1  
2481-94-9  
248606-45-3  
24886-73-5  
2491-32-9  
24937-16-4  
24937-78-8  
24937-78-8  
24937-78-8  
24937-78-8  
24937-79-9  
24937-79-9  
24968-12-5  
24991-55-7  
2500-88-1  
25013-16-5  
25036-13-9

ir-spectra.com IS NIR Spectra

---

**IS NIR Spectra – Complete Collection**

6049 NIR spectra

Copyright © NICODOM 2012

NICODOM Ltd., Hlavni 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

[ftirspectra@gmail.com](mailto:ftirspectra@gmail.com)

25036-13-9  
25037-45-0  
25037-45-0  
25038-36-2  
25038-36-2  
25038-36-2  
25038-54-4  
25038-54-4  
25038-54-4  
25038-59-9  
25038-59-9  
25038-59-9  
25053-15-0  
25053-15-0  
25067-59-8  
25086-89-9  
25086-89-9  
25087-26-7  
25101-03-5  
2512-29-0  
25135-51-7, 25154-01-2  
251474-95-0  
25151-00-2  
25154-52-3  
25168-24-5  
2517-04-6  
2517-43-3  
25189-55-3  
25190-06-1  
2519-30-4  
25213-24-5  
25213-39-2  
25249-54-1  
25249-54-1  
25249-54-1  
25265-77-4  
2530-83-8  
2530-85-0  
25322-68-3  
25322-68-3  
25322-68-3  
25322-68-3  
25322-68-3  
25322-68-3

ir-spectra.com IS NIR Spectra

---

**IS NIR Spectra – Complete Collection**

6049 NIR spectra

Copyright © NICODOM 2012

NICODOM Ltd., Hlavni 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

[ftirspectra@gmail.com](mailto:ftirspectra@gmail.com)

25322-68-3  
25322-68-3  
25322-68-3  
25322-68-3  
25322-69-4  
25322-69-4  
25332-39-2  
25360-72-9  
25383-99-7  
2538-85-4  
25425-97-2  
254436-04-9  
25448-25-3  
25470-96-6  
25498-49-1  
2550-06-3#  
2553-08-4#  
25550-98-5  
25587-80-8  
25637-99-4  
25667-11-2  
2568-33-4  
25704-18-1  
25713-23-9  
2576-47-8  
25805-17-8  
2580-56-5  
2580-78-1  
2580-78-1  
2581-34-2  
2582-30-1  
25826-34-0, 12237-86-4  
25852-47-5  
2586-62-1  
25895-60-7  
25910-85-4  
2592-95-2  
25952-53-8  
25955-51-5  
25956-17-6  
2595-98-4  
259731-40-3  
25973-55-1  
25987-06-8

ir-spectra.com IS NIR Spectra

---

**IS NIR Spectra – Complete Collection**

6049 NIR spectra

Copyright © NICODOM 2012

NICODOM Ltd., Hlavni 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

[ftirspectra@gmail.com](mailto:ftirspectra@gmail.com)

2602-36-0  
26027-38-3  
260-94-6  
26098-55-5  
2610-05-1  
2610-10-8  
2610-11-9, 25188-42-5  
2611-82-7  
2611-82-7  
26156-84-3  
2622-14-2  
26222-40-2  
26266-57-9  
26279-24-3  
26401-97-8  
2644-70-4  
2646-17-5  
26472-00-4  
264917-73-9  
264917-87-5  
2650-17-1  
26523-78-4  
26544-38-7  
26628-22-8  
26635-93-8  
2664-63-3  
26658-19-5  
2666-17-3  
26677-99-6  
267241-27-0  
26761-45-5  
2679-01-8  
26807-65-8  
26826-40-4  
2687-12-9#  
2687-91-4  
2687-94-7  
2687-96-9  
26896-18-4  
2691-41-0  
2706-28-7  
27107-88-6  
271261-26-8  
27176-87-0

ir-spectra.com IS NIR Spectra

---

**IS NIR Spectra – Complete Collection**

6049 NIR spectra

Copyright © NICODOM 2012

NICODOM Ltd., Hlavni 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

[ftirspectra@gmail.com](mailto:ftirspectra@gmail.com)

27177-63-5  
27206-35-5  
27208-38-4  
27312-17-0, 6408-72-6  
273-13-2  
27318-90-7  
27329-60-8  
2743-38-6  
2751-90-8  
27532-96-3  
27546-07-2  
27565-41-9  
2758-06-7#  
2758-18-1  
2766-77-0  
27676-62-6  
2768-02-7#  
2782-07-2  
2783-17-7  
2783-94-0  
2783-94-0  
27848-81-3  
2786-76-7  
27896-84-0  
28020-73-7  
280-57-9  
2809-21-4  
281-23-2  
2814-77-9  
28179-44-4  
28258-59-5  
2829-43-8, 25255-10-1  
2832-45-3  
2834-05-1#  
2835-06-5#  
2835-77-0  
2835-81-6  
28406-56-6  
28407-37-6  
28539-02-8  
28553-12-0  
2861-02-1  
28631-66-5  
28697-53-2

ir-spectra.com IS NIR Spectra

---

**IS NIR Spectra – Complete Collection**

6049 NIR spectra

Copyright © NICODOM 2012

NICODOM Ltd., Hlavni 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

[ftirspectra@gmail.com](mailto:ftirspectra@gmail.com)

287-92-3  
28804-96-8  
288-13-1  
288-32-4  
288-88-0  
2893-78-9  
28940-11-6  
2897-60-1  
28981-97-7  
28983-56-4  
28983-56-4  
28983-56-4  
2899-60-7  
290-37-9  
29079-00-3  
29094-61-9  
29117-01-9  
29122-68-7  
2917-26-2  
2920-38-9  
2923-16-2  
29364-29-2  
29387-86-8  
29412-85-9  
2944-58-3  
2944-59-4  
2945-03-1#  
2945-96-2, 25180-12-5  
294-62-2  
294-90-6  
294-93-9  
2950-43-8  
29556-33-0  
29598-76-3  
2969-81-5  
29712-28-5  
29712-30-9  
29712-32-1  
2971-90-6  
298-12-4  
2985-59-3  
298-57-7  
29858-07-9  
29860-47-7, 50586-59-9, 99488-65-0

ir-spectra.com IS NIR Spectra

---

**IS NIR Spectra – Complete Collection**

6049 NIR spectra

Copyright © NICODOM 2012

NICODOM Ltd., Hlavní 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

[ftirspectra@gmail.com](mailto:ftirspectra@gmail.com)

298-83-9  
298-95-3  
298-96-4  
29898-32-6  
299199-23-8  
299-39-8  
300-39-0  
3005-32-1  
3006-15-3  
300-87-8  
301-04-2  
30112-70-0  
30112-70-0  
30112-70-0, 12731-53-2  
3011-34-5  
3012-65-5  
30211-73-5  
302-17-0  
302-72-7  
30282-54-3  
302-84-1  
302-95-4  
303-49-1  
304-55-2  
304-59-6  
304-88-1  
30499-70-8  
3051-09-0  
3051-11-4  
30525-89-4  
3056-93-7  
305-84-0  
306-08-1  
3061-42-5  
30653-05-5  
3068-39-1  
3071-73-6  
30737-24-7  
307531-90-4  
307531-92-6  
308068-88-4  
308069-26-3  
308075-91-4  
308083-59-2

ir-spectra.com IS NIR Spectra

---

**IS NIR Spectra – Complete Collection**

6049 NIR spectra

Copyright © NICODOM 2012

NICODOM Ltd., Hlavni 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

[ftirspectra@gmail.com](mailto:ftirspectra@gmail.com)

3084-40-0  
3087-16-9  
3087-36-3  
30899-19-5  
309-00-2  
309252-26-4  
309252-29-7  
3094-09-5#  
30947-30-9  
30948-01-7  
3095-65-6  
3109-63-5  
31121-93-4  
311-28-4  
3115-68-2  
31183-11-6  
3118-97-6  
3121-71-9  
312-30-1  
312693-80-4  
312-84-5  
3130-75-4  
3130-87-8  
313-72-4  
3141-26-2  
3142-42-5  
3143-02-0  
3144-16-9  
31566-31-1  
31570-04-4  
31570-04-4  
3164-29-2  
3164-34-9  
31694-16-3  
317-34-0  
3173-53-3  
31772-06-2  
31778-10-6  
31778-10-6  
317802-08-7  
3179-63-3  
31799-71-0  
31810-89-6, 12217-79-7  
3183-08-2#

ir-spectra.com IS NIR Spectra

---

**IS NIR Spectra – Complete Collection**

6049 NIR spectra

Copyright © NICODOM 2012

NICODOM Ltd., Hlavni 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

[ftirspectra@gmail.com](mailto:ftirspectra@gmail.com)

3184-13-2  
31879-05-7  
318-98-9  
3196-73-4  
319-78-8  
3200-06-4  
32005-36-0  
321-14-2  
321-30-2  
32131-17-2  
3214-47-9, 25738-24-3  
32220-82-9  
32221-81-1  
3234-85-3  
3244-88-0  
32503-27-8  
32509-66-3  
3251-23-8  
3264-82-2  
32660-24-5  
32687-78-8  
3268-79-9  
326-91-0  
32708-95-5  
3271-76-9  
327-56-0  
327-57-1  
3279-54-7  
3282-30-2  
3283-12-3#  
328-38-1  
328-50-7  
32870-98-7  
3290-92-4  
329-71-5  
33005-95-7  
33008-06-9  
33085-32-4  
331-39-5  
33203-82-6  
33227-76-8  
33270-70-1  
332-77-4  
33286-22-5

ir-spectra.com IS NIR Spectra

---

**IS NIR Spectra – Complete Collection**

6049 NIR spectra

Copyright © NICODOM 2012

NICODOM Ltd., Hlavni 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

[ftirspectra@gmail.com](mailto:ftirspectra@gmail.com)

332863-88-4  
333-20-0  
3333-62-8  
3333-67-3  
33342-05-1  
334-48-5  
3349-06-2#  
33504-03-9  
3351-05-1  
3353-05-7  
3354-97-0, 25180-30-7  
33577-16-1  
3375-31-3  
33787-28-9  
3380-34-5  
338-69-2  
33893-89-9  
3392-07-2#  
3396-11-0  
33973-59-0  
3398-22-9  
34002-94-3  
34381-68-5  
3441-14-3  
34433-71-1  
3445-11-2#  
34523-28-9  
3453-00-7  
34552-84-6  
3458-28-4  
3458-72-8  
34590-94-8  
345989-86-8  
3466-75-9  
3468-63-1  
3469-26-9  
34699-71-3  
3470-34-9  
34722-90-2  
34807-26-6  
348083-27-2  
348-52-7  
348-54-9  
3485-80-1

ir-spectra.com IS NIR Spectra

---

**IS NIR Spectra – Complete Collection**

6049 NIR spectra

Copyright © NICODOM 2012

NICODOM Ltd., Hlavni 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

[ftirspectra@gmail.com](mailto:ftirspectra@gmail.com)

3486-30-4  
3486-35-9  
34946-82-2  
350-03-8  
35099-99-1  
351002-92-1  
35103-79-8  
3520-42-1  
3520-72-7  
352680-53-6  
352-97-6  
3536-49-0  
354526-75-3  
3546-41-6  
3547-38-4  
3550-43-4  
35602-69-8  
3564-18-9  
3567-66-6  
3567-69-9  
3567-69-9, 53026-69-9  
3568-26-1  
3570-80-7  
35725-34-9  
357413-27-5  
357-57-3  
358-23-6  
35869-64-8  
35948-25-5  
35963-20-3  
3612-20-2  
3615-37-0  
3615-41-6  
3615-44-9  
36177-92-1  
3618-12-0  
3625-57-8  
3626-36-6, 25188-23-2  
3626-40-2, 25180-33-0  
36322-90-4  
36330-85-5  
36339-47-6  
3634-89-7  
363-52-0

ir-spectra.com IS NIR Spectra

---

**IS NIR Spectra – Complete Collection**

6049 NIR spectra

Copyright © NICODOM 2012

NICODOM Ltd., Hlavni 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

[ftirspectra@gmail.com](mailto:ftirspectra@gmail.com)

36393-42-7  
36443-68-2  
3648-36-0  
36548-87-5  
36653-82-4  
367-12-4  
3671-99-6  
3674-13-3  
36768-62-4  
36791-04-5  
367-93-1  
3682-35-7  
368-39-8  
3687-18-1  
36876-02-5  
36888-99-0  
3695-00-9  
36951-72-1  
369632-52-0  
37052-78-1  
370-81-0  
37191-69-8  
37206-32-9  
37211-54-4  
37229-23-5  
37279-47-3  
37294-49-8  
3731-51-9  
3731-52-0  
3731-53-1  
37338-62-8  
373-44-4  
3734-67-6  
3734-67-6  
37353-59-6  
37372-87-5, 81064-47-3  
375-22-4  
37625-56-2  
3769-57-1  
3770-50-1  
3772-44-9  
3775-90-4  
379-79-3  
38002-45-8

ir-spectra.com IS NIR Spectra

---

**IS NIR Spectra – Complete Collection**

6049 NIR spectra

Copyright © NICODOM 2012

NICODOM Ltd., Hlavni 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

[ftirspectra@gmail.com](mailto:ftirspectra@gmail.com)

3810-63-7  
3810-74-0  
38117-54-3  
38222-83-2  
38416-67-0  
3844-45-9  
3844-45-9  
3844-45-9  
3844-94-8  
3846-71-7  
38613-77-3  
3861-73-2  
3864-99-1  
38662-28-1  
38668-48-3  
3891-59-6  
38916-42-6  
3896-11-5#  
3905-19-9  
392230-83-0  
3923-52-2  
39236-46-9  
3926-62-3  
39277-77-5  
39279-59-9, 38151-74-5  
39279-69-1  
39336-74-8  
39365-88-3  
3937-56-2  
39389-20-3, 52038-26-3  
39390-62-0  
39393-39-0  
39393-70-9  
39409-82-0  
39423-51-3  
39429-60-2  
39443-66-8  
39456-78-5  
3946-91-6  
3956-55-6, 12239-34-8  
3964-18-9  
3992-42-5  
39968-33-7  
40018-26-6

---

**IS NIR Spectra – Complete Collection**

6049 NIR spectra

Copyright © NICODOM 2012

NICODOM Ltd., Hlavni 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

[ftirspectra@gmail.com](mailto:ftirspectra@gmail.com)

40070-59-5  
4009-98-7  
40149-83-5  
4016-14-2  
40231-03-6  
40243-08-1  
4051-63-2  
4053-08-1  
40618-31-3  
4065-45-6  
4066-02-8#  
4070-80-8  
40716-47-0  
40850-57-5  
4101-68-2  
4106-67-6  
41100-52-1  
41137-60-4  
4115-52-0  
4129-84-4  
41337-81-9  
41444-62-6  
41484-35-9  
4152-09-4#  
41556-26-7  
4160-72-9  
41642-51-7  
41669-30-1  
4170-30-3  
417708-80-6  
4180-23-8  
4187-87-5  
4196-86-5  
4197-25-5  
4199-88-6  
4200-06-0  
4202-14-6  
4203-77-4  
4208-80-4, 67815-70-7  
4210-32-6  
4211-67-0  
4215-80-9  
4216-01-7  
4216-02-8

ir-spectra.com IS NIR Spectra

---

**IS NIR Spectra – Complete Collection**

6049 NIR spectra

Copyright © NICODOM 2012

NICODOM Ltd., Hlavní 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

[ftirspectra@gmail.com](mailto:ftirspectra@gmail.com)

4221-80-1  
4224-87-7  
4224-96-8  
4252-78-2  
4254-17-5  
425659-58-1  
42612-12-4  
42616-96-6  
4271-30-1#  
42765-77-5  
42783-06-2, 70729-65-6, 59709-38-5  
42844-93-9  
42882-31-5  
4292-19-7  
431-03-8  
4310-80-9  
4314-14-1  
4316-74-9  
4317-14-0  
43200-80-2  
4321-69-1  
4342-61-4  
4360-12-7  
4360-63-8  
4368-56-3  
436-910-3  
437713-74-1  
437-74-1  
4378-61-4  
4394-85-8  
4399-55-7  
4403-90-1  
4405-13-4  
4409-98-7  
441-38-3  
4424-6-0  
4424-87-7  
4430-20-0  
4432-31-9  
443-48-1  
443-72-1  
4450-94-6  
44813-14-1  
4514-87-8

ir-spectra.com IS NIR Spectra

---

**IS NIR Spectra – Complete Collection**

6049 NIR spectra

Copyright © NICODOM 2012

NICODOM Ltd., Hlavni 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

[ftirspectra@gmail.com](mailto:ftirspectra@gmail.com)

4526-07-2#  
4531-49-1  
4531-49-1  
4548-53-2  
455-18-5  
4553-89-3  
4553-89-3  
456-22-4  
457931-26-9  
459-57-4  
460740-96-9  
460741-05-3  
461-58-5  
461-72-3  
461-84-7  
461-89-2  
462-06-6  
462-08-8  
4637-24-5  
465-42-9  
465-65-6  
4657-00-5  
4672-35-9  
4680-78-8  
4702-13-0  
470-82-6  
471-25-0  
4712-55-4  
471-34-1  
471-34-1  
471-46-5  
4720-82-5  
47230-38-6  
4727-50-8  
4730-75-0  
4731-84-4  
4733-39-5  
473711-70-5  
474-25-9  
475-31-0  
4755-77-5  
4767-03-7#  
477572-63-7  
477-73-6

ir-spectra.com IS NIR Spectra

---

**IS NIR Spectra – Complete Collection**

6049 NIR spectra

Copyright © NICODOM 2012

NICODOM Ltd., Hlavni 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

[ftirspectra@gmail.com](mailto:ftirspectra@gmail.com)

477740-25-3  
479-27-6  
479-33-4  
479-45-8  
479-92-5  
480-16-0  
4812-40-2  
482-89-3  
483-55-6  
483-84-1  
484-11-7  
4851-50-7  
485-47-2  
485-49-4  
486-25-9  
487-89-8  
488-38-0  
488791-77-1  
488-81-3  
488-81-3  
4891-38-7  
4913-13-7  
491-35-0  
492462-76-7  
492-94-4  
492-99-9  
493-52-7  
4940-11-8  
494-38-2  
4948-15-6  
49540-00-3  
495-54-5  
49562-28-9  
495-69-2  
496-15-1  
496-67-3  
497-19-8  
497-76-7  
498-21-5  
498-23-7  
498-59-9  
499-83-2  
50-01-1  
5001-51-4

ir-spectra.com IS NIR Spectra

---

**IS NIR Spectra – Complete Collection**

6049 NIR spectra

Copyright © NICODOM 2012

NICODOM Ltd., Hlavni 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

[ftirspectra@gmail.com](mailto:ftirspectra@gmail.com)

50-03-3  
50-06-6  
501-53-1  
501-65-5  
501-97-3  
50-21-5  
502-41-0  
502-49-8  
50-29-3  
50370-12-2  
503-87-7  
504-02-9  
504-17-6  
504-24-5  
504-29-0  
5045-40-9  
504-63-2  
50-48-6  
50-52-2  
5053-08-7#  
50-54-4  
505-48-6  
50-55-5  
505-66-8  
505-95-3  
50-59-9  
506-12-7  
50-63-5  
506-46-7  
506-59-2  
506-61-6  
506-68-3  
50679-08-8  
506-87-6  
50-69-1  
506-93-4  
50-70-4  
507-09-5  
50-71-5  
507-28-8  
507-70-0  
50-78-2  
50814-31-8  
50-81-7

ir-spectra.com IS NIR Spectra

---

**IS NIR Spectra – Complete Collection**

6049 NIR spectra

Copyright © NICODOM 2012

NICODOM Ltd., Hlavni 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

[ftirspectra@gmail.com](mailto:ftirspectra@gmail.com)

50858-51-0  
50890-67-0  
5089-16-7, 25489-31-0  
50-89-5  
50918-71-3  
50925-42-3  
509-34-2  
50-98-6  
51012-33-0  
5102-83-0  
5104-49-4  
51-05-8  
5108-96-3  
51-17-2  
5124-25-4  
512-56-1  
51274-00-1  
51274-00-1  
51-28-5  
5131-60-2  
5131-66-8  
51333-22-3  
51-35-4  
51373-56-9  
5137-55-3  
513-77-9  
513-78-0  
513-85-9  
514-10-3  
5141-20-8  
51-42-3  
5142-76-7  
51-43-4  
51-45-6  
5153-67-3  
51-55-8  
51-56-9  
515-69-5  
515-74-2  
51595-55-2  
516-02-9  
51-60-5  
516-06-3  
51-61-6

ir-spectra.com IS NIR Spectra

---

**IS NIR Spectra – Complete Collection**

6049 NIR spectra

Copyright © NICODOM 2012

NICODOM Ltd., Hlavni 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

[ftirspectra@gmail.com](mailto:ftirspectra@gmail.com)

517-28-2  
517-51-1  
5175-83-7  
517-88-4, 517-90-8  
517906-55-7  
51-79-6  
518-47-8  
518-47-8  
51-85-4  
51868-46-3, 52697-38-8, 13301-61-6  
51877-53-3  
51934-41-9  
519-34-6, 70892-58-9, 480-16-0, 99328-49-1, 85595-40-0  
519-73-3  
52001-91-9  
520-33-2  
520-34-3  
520-45-6  
52080-58-7  
521-13-1  
521-24-4  
521-31-3  
521-74-4  
52218-73-2  
52256-37-8  
522608-60-2  
522-75-8  
52315-92-1  
52320-66-8  
523-27-3  
5232-99-5  
523-44-4  
52439-77-7  
52462-29-0  
52-49-3  
525-05-3  
52-51-7  
525573-16-4  
525-79-1  
52596-06-2  
52598-24-0  
5261-31-4  
52627-00-6  
52688-46-7

ir-spectra.com IS NIR Spectra

---

**IS NIR Spectra – Complete Collection**

6049 NIR spectra

Copyright © NICODOM 2012

NICODOM Ltd., Hlavní 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

[ftirspectra@gmail.com](mailto:ftirspectra@gmail.com)

526-94-3  
526-99-8  
527-07-1  
5271-38-5  
52722-53-9  
5272-36-6  
52746-49-3  
527-69-5  
527697-23-0  
5280-68-2  
5280-68-2  
5280-78-4  
5280-80-8  
5281-04-9  
52829-07-9  
528-50-7  
52-89-1  
52-90-4  
5290-62-0  
529-12-4  
529-19-1  
5292-43-3  
529-28-2  
5293-84-5  
5296-78-6  
53028-40-3  
530-43-8  
53118-99-3  
53124-00-8  
53124-74-6  
531-75-9  
531-85-1  
531-86-2  
531-91-9  
53193-09-2  
532-02-5  
532-32-1  
532-82-1  
5329-14-6  
53293-00-8  
532-94-5  
5332-26-3  
5332-73-0  
5337-17-7

ir-spectra.com IS NIR Spectra

---

**IS NIR Spectra – Complete Collection**

6049 NIR spectra

Copyright © NICODOM 2012

NICODOM Ltd., Hlavni 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

[ftirspectra@gmail.com](mailto:ftirspectra@gmail.com)

5337-19-9  
533-87-9  
534-03-2  
5341-61-7  
534-17-8  
53-41-8  
53-43-0  
5344-90-1  
5345-01-7  
534-52-1  
5346-90-7  
5349-80-4  
5350-57-2  
535-34-2  
5355-16-8  
536-17-4  
53626-51-0  
53637-25-5  
53639-82-0, 6856-31-1  
536-46-9  
536-90-3  
53716-49-7  
537-45-1  
537-55-3  
538-51-2  
538-58-9  
53-86-1  
538-65-8  
538-74-9  
538-75-0  
538-93-2  
5390-04-5#  
5392-12-1  
539-47-9  
539-48-0  
5394-83-2  
539-74-2  
539-88-8  
53994-73-3  
540-10-3  
54018-40-5  
540-23-8  
5402-73-3  
540-37-4

ir-spectra.com IS NIR Spectra

---

**IS NIR Spectra – Complete Collection**

6049 NIR spectra

Copyright © NICODOM 2012

NICODOM Ltd., Hlavni 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

[ftirspectra@gmail.com](mailto:ftirspectra@gmail.com)

540-38-5  
54039-38-2  
540-51-2  
54060-92-3  
54063-53-5  
540-69-2  
540-84-1  
540-88-5  
541-09-3  
54-11-5  
541-16-2  
54-12-6  
5413-75-2  
541-41-3  
541-47-9  
54-16-0  
541-69-5  
54187-96-1  
542-16-5  
54-21-7  
542-32-5  
542-59-6  
542-69-8  
5429-56-1  
5431-44-7  
543-24-8  
5437-45-6  
543-80-6  
543-81-7  
543-90-8  
543-94-2  
544-01-4  
544-10-5  
544-17-2  
544-40-1  
544-63-8  
544-85-4  
544-92-3  
54527-84-3  
54553-90-1  
5459-58-5  
5460-09-3  
54-64-8  
54660-00-3

ir-spectra.com IS NIR Spectra

---

**IS NIR Spectra – Complete Collection**

6049 NIR spectra

Copyright © NICODOM 2012

NICODOM Ltd., Hlavni 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

[ftirspectra@gmail.com](mailto:ftirspectra@gmail.com)

546-89-4  
546-93-0  
5469-45-4  
5470-11-1  
5471-63-6  
54-71-7  
547-58-0  
547-91-1  
548-00-5  
54827-17-7  
54-85-3  
548-62-9  
548-80-1  
549-18-8  
550-74-3  
550-97-0  
55142-85-3  
5521-31-3  
55216-11-0  
55-22-1  
552-89-6  
553-12-8  
553-24-2  
553-26-4  
55349-16-1  
5536-61-8  
553-90-2  
554-00-7  
55406-53-6  
554-12-1  
554-13-2  
554-14-3  
55442-07-4  
55-48-1  
554-84-7  
55488-87-4  
55488-87-4  
554-95-0  
5550-12-9  
555-30-6  
555-44-2  
555-45-3  
55-55-0  
55-56-1

ir-spectra.com IS NIR Spectra

---

**IS NIR Spectra – Complete Collection**

6049 NIR spectra

Copyright © NICODOM 2012

NICODOM Ltd., Hlavni 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

[ftirspectra@gmail.com](mailto:ftirspectra@gmail.com)

555-75-9  
55589-62-3  
556-03-6  
556-48-9  
556-50-3  
5567-15-7  
5567-15-7  
557-04-0  
557-05-1  
557-08-4  
557-09-5  
557-19-7  
557-34-6  
557-39-1  
5574-97-0  
557-66-4  
55777-80-5  
5580-57-4  
5580-58-5  
55840-60-3  
55904-02-4  
55913-03-6  
55913-05-8  
55913-06-9  
55913-07-0  
55934-93-5  
5600-62-4  
5601-29-6  
5601-29-6  
56090-78-9  
56090-80-3  
5610-64-0  
561-20-6  
56-12-2  
5613-46-7  
56-17-7  
56187-04-3  
56-18-8  
56189-09-4  
56238-63-2  
56257-90-0  
563-03-1  
5630-53-5  
56317-51-2

ir-spectra.com IS NIR Spectra

---

**IS NIR Spectra – Complete Collection**

6049 NIR spectra

Copyright © NICODOM 2012

NICODOM Ltd., Hlavni 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

[ftirspectra@gmail.com](mailto:ftirspectra@gmail.com)

563-41-7  
56-35-9  
563-63-3  
563-67-7  
563-68-8  
56-37-1  
563-76-8  
563-80-4  
56392-17-7  
56-40-6  
56-41-7  
564-25-0  
56-45-1  
56491-86-2  
56512-49-3  
56548-64-2  
56731-19-2  
56-75-7  
56773-42-3  
5680-79-5  
56-81-5  
56-82-6  
5683-31-8  
56-84-8  
56-85-9  
56-86-0  
56-89-3  
56-92-8  
56-93-9  
569-58-4  
569-61-9  
56962-08-4  
56996-68-0  
570-08-1  
570-24-1  
5704-04-1#  
57-08-9  
57-09-0  
57-10-3  
57-11-4  
57-13-6  
57-15-8  
571-58-4  
571-61-9

ir-spectra.com IS NIR Spectra

---

**IS NIR Spectra – Complete Collection**

6049 NIR spectra

Copyright © NICODOM 2012

NICODOM Ltd., Hlavni 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

[ftirspectra@gmail.com](mailto:ftirspectra@gmail.com)

57248-88-1  
57285-14-0  
57-30-7  
57308-41-5  
573-56-8  
573-58-0  
57-41-0  
5743-26-0  
5743-49-7  
57-44-3  
57455-37-5  
57-48-7  
57-50-1  
57-50-1  
57516-68-4  
57-55-6  
57-56-7  
57584-27-7  
57637-97-5  
57-67-0  
57-68-1  
57-71-6  
5785-44-4  
578-57-4  
57-87-4  
57-88-5  
5794-13-8  
579-44-2  
57971-97-8  
579-75-9  
58051-96-0  
58-08-2  
581-47-5  
58-15-1  
582-17-2  
582-22-9  
5824-40-8  
58-27-5  
58-33-3  
583-39-1  
583-55-1  
583-78-8  
58392-07-7  
58392-11-3

ir-spectra.com IS NIR Spectra

---

**IS NIR Spectra – Complete Collection**

6049 NIR spectra

Copyright © NICODOM 2012

NICODOM Ltd., Hlavni 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

[ftirspectra@gmail.com](mailto:ftirspectra@gmail.com)

584-08-7  
584-09-8  
5844-07-5  
584-42-9  
584-48-5  
58-46-8  
584-84-9  
5850-39-5  
5850-44-2  
58-55-9  
5858-39-9  
585-86-4  
58613-37-9  
58-61-7  
5863-46-7, 25305-80-0  
58-63-9  
586-76-5  
586-78-7  
586-98-1  
587-04-2  
5872-08-2  
58-74-2  
587-85-9  
587-98-4  
58816-01-6  
58-86-6  
589-15-1  
589-17-3  
5892-10-4#  
58-95-7  
589-87-7  
589-92-4  
590-00-1  
59-02-9  
590-29-4  
590-36-3  
590-46-5  
590-67-0  
590-86-3  
5908-80-5  
591-01-5  
591-08-2  
591-19-5  
591-27-5

ir-spectra.com IS NIR Spectra

---

**IS NIR Spectra – Complete Collection**

6049 NIR spectra

Copyright © NICODOM 2012

NICODOM Ltd., Hlavni 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

[ftirspectra@gmail.com](mailto:ftirspectra@gmail.com)

59130-69-7  
591-31-1  
5913-70-2  
591-50-4  
59160-29-1  
591-80-0  
59-23-4  
592-76-7  
59277-89-3  
592-84-7  
592-85-8  
592-87-0  
592-88-1  
5934-29-2  
593-51-1  
593-85-1  
59-43-8  
594-60-5  
59-46-1  
594-65-0  
5949-16-6  
5949-29-1  
59-51-8  
59568-28-4  
595-89-1  
595-90-4  
596-27-0  
59-66-5  
59-67-6  
59686-68-9  
5970-45-6  
5970-47-8  
5972-73-6  
59729-32-7  
59786-91-3  
5979-28-2  
59804-37-4  
598-04-9  
598-62-9  
598-63-0  
598-72-1  
598-75-4  
598-78-7  
59-88-1

ir-spectra.com IS NIR Spectra

---

**IS NIR Spectra – Complete Collection**

6049 NIR spectra

Copyright © NICODOM 2012

NICODOM Ltd., Hlavni 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

[ftirspectra@gmail.com](mailto:ftirspectra@gmail.com)

598-82-3  
5989-81-1  
59918-86-4  
59-92-7  
5995-25-5  
599-64-4  
59-97-2  
60001-11-8  
60-00-4  
60-01-5  
60-09-3  
6009-70-7  
60-10-6  
60-11-7  
60-13-9  
60147-22-0  
60177-39-1  
60181-77-3  
60-18-4  
6018-89-9  
60-19-5  
60223-95-2  
602302-04-5  
60239-18-1  
60-24-2  
60-29-7  
603-00-9  
603-11-2  
60-32-2  
603-32-7  
60-33-3  
603-34-9  
603-35-0  
603-36-1  
603-45-2  
603-52-1  
60-35-5  
603-76-9  
603-79-2  
6038-19-3  
603-83-8  
6041-94-7  
604-32-0  
604-35-3

ir-spectra.com IS NIR Spectra

---

**IS NIR Spectra – Complete Collection**

6049 NIR spectra

Copyright © NICODOM 2012

NICODOM Ltd., Hlavni 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

[ftirspectra@gmail.com](mailto:ftirspectra@gmail.com)

604-59-1  
6046-93-1  
6046-93-1  
604766-57-6  
605-65-2  
60668-06-6  
60668-11-3  
60668-24-8  
6066-82-6  
60668-50-0  
60687-93-6  
60703-84-6  
607-88-5  
60-80-0  
6080-56-4  
608-08-2  
60-81-1  
608-50-4  
608-66-2  
608-68-4  
6087-13-4  
609-06-3  
609-08-5  
609-14-3  
609-15-4  
609844-11-3  
6099-90-7  
609-99-4  
6100-05-6  
6100-16-9  
610-02-6  
610-27-5  
61036-28-0  
6104-58-1  
6104-59-2  
6106-21-4  
6106-24-7  
610-97-9  
611-10-9  
611-13-2  
611-14-3  
61116-28-7  
61-12-1  
61128-46-9

ir-spectra.com IS NIR Spectra

---

**IS NIR Spectra – Complete Collection**

6049 NIR spectra

Copyright © NICODOM 2012

NICODOM Ltd., Hlavni 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

[ftirspectra@gmail.com](mailto:ftirspectra@gmail.com)

61167-58-6  
61177-45-5  
6119-47-7  
6119-70-6  
612-22-6  
61-25-6  
612-82-8  
613-03-6  
6130-98-9  
6131-90-4  
6132-02-1  
6132-04-3  
6132-04-3  
6132-05-4#  
61336-70-7  
613-45-6  
6134-79-8  
613-78-5  
6138-23-4  
613-94-5  
614-16-4  
614-18-6  
6147-53-1  
61512-61-6, 73309-48-5  
6151-40-2  
6153-39-5  
6153-56-6  
615-37-2  
615-42-9  
615-43-0  
615-50-9  
6155-57-3  
615-67-8  
61583-60-6  
6159-44-0  
616-05-7  
616-34-2  
6163-58-2  
616-45-5  
616-47-7  
616-68-2  
616-91-1  
61703-04-6, 6441-84-5  
61703-12-6

ir-spectra.com IS NIR Spectra

---

**IS NIR Spectra – Complete Collection**

6049 NIR spectra

Copyright © NICODOM 2012

NICODOM Ltd., Hlavni 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

[ftirspectra@gmail.com](mailto:ftirspectra@gmail.com)

61-72-3  
61724-72-9, 70210-31-0  
61724-78-5, 70304-38-0  
61725-71-1  
61725-76-6  
61725-81-3  
61-73-4  
617-35-6  
617-48-1  
61-76-7  
61789-32-0  
61789-91-1  
61790-81-6  
61791-14-8  
61791-31-9  
61813-61-4, 73507-36-5  
61813-62-5  
61813-86-3, 72208-28-7  
61813-90-9  
61814-52-6  
61815-09-6  
61-82-5  
61827-42-7  
61869-08-7  
618-98-4  
61901-39-1  
61901-87-9  
61901-87-9  
61901-91-5  
61901-92-6  
61901-95-9  
61902-11-2  
619-04-5  
61-90-5  
619-08-9  
61916-41-4  
6192-52-5  
61931-02-0  
61931-55-3  
61931-83-7  
61932-63-6  
619-42-1  
61951-47-1  
61951-76-6, 70865-31-5

ir-spectra.com IS NIR Spectra

---

**IS NIR Spectra – Complete Collection**

6049 NIR spectra

Copyright © NICODOM 2012

NICODOM Ltd., Hlavni 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

[ftirspectra@gmail.com](mailto:ftirspectra@gmail.com)

619-58-9  
61968-07-8  
61968-12-5  
61968-26-1  
61968-38-5, 68391-42-4  
61968-52-3  
61968-71-6  
61968-98-7  
61969-48-0  
61969-51-5  
6197-30-4  
62005-64-5  
620-13-3  
620-40-6  
62068-38-6  
620-89-3  
621-82-9  
621-84-1  
6219-63-2, 643-56-1  
6219-73-4  
622-08-2  
622-24-2  
62-23-7  
6223-83-2  
622-40-2  
622-64-0  
6226-79-5  
6226-87-5  
622-88-8  
6229-00-1  
622-93-5  
623-00-7  
623-03-0  
623-24-5  
623-25-6  
623-26-7  
623-27-8  
623-37-0  
623-43-8  
623-47-2  
623-50-7  
623-81-4  
62-38-4  
623-84-7

ir-spectra.com IS NIR Spectra

---

**IS NIR Spectra – Complete Collection**

6049 NIR spectra

Copyright © NICODOM 2012

NICODOM Ltd., Hlavni 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

[ftirspectra@gmail.com](mailto:ftirspectra@gmail.com)

623-91-6  
624-18-0  
624-38-4  
62-44-2  
624-48-6  
624-49-7  
624-54-4  
6245-60-9, 25317-07-1  
62-50-0  
6250-23-3  
625-23-0  
62-53-3  
62-54-4  
625-52-5  
625534-03-4  
62-55-5  
62-56-6  
62571-86-2  
62-57-7  
625-84-3  
625853-76-1  
625-86-5  
626-17-5  
626-19-7  
6262-07-3  
6262-07-3  
62625-32-5  
62637-93-8  
626-39-1  
626-43-7  
626-51-7  
62668-02-4  
626-93-7  
626-95-9  
627-05-4  
62708-56-9  
627-30-5  
627-41-8  
62-76-0  
627-75-8  
628-13-7  
6282-02-6#  
628-20-6  
6284-40-8

ir-spectra.com IS NIR Spectra

---

**IS NIR Spectra – Complete Collection**

6049 NIR spectra

Copyright © NICODOM 2012

NICODOM Ltd., Hlavni 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

[ftirspectra@gmail.com](mailto:ftirspectra@gmail.com)

628-61-5  
628-63-7  
628-97-7  
629-03-8  
629-04-9  
629-11-8  
6291-42-5  
629-15-2  
629-30-1  
629-73-2  
629-78-7  
629-82-3  
62992-46-5  
630-02-4  
630-06-8  
630-18-2  
63026-01-7  
6303-21-5  
6305-43-7  
6309-30-4  
6309-51-9  
63126-87-4  
63148-57-2  
63148-62-9  
63148-62-9  
63148-62-9  
63148-65-2  
63148-69-6  
6315-52-2  
631-61-8  
6320-14-5  
632-20-2  
632-22-4  
6324-00-1  
632-69-9  
632-79-1  
632-99-5  
633-03-4  
633-31-8  
6334-25-4  
63394-02-5  
633-96-5  
63409-48-3  
63-42-3

ir-spectra.com IS NIR Spectra

---

**IS NIR Spectra – Complete Collection**

6049 NIR spectra

Copyright © NICODOM 2012

NICODOM Ltd., Hlavni 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

[ftirspectra@gmail.com](mailto:ftirspectra@gmail.com)

63449-39-8,8029-39-8  
63451-28-5  
634558-40-0  
634909-88-9  
634-97-9  
635-65-4  
6358-30-1  
6358-31-2  
6358-36-7  
6358-69-6  
6359-82-6, 25739-63-3  
6359-85-9, 25739-67-7  
6359-98-4  
6359-98-4  
6360-07-2  
636-30-6  
63641-88-3  
636-54-4  
636-61-3  
636-72-6  
6367-75-5  
63-68-3  
6368-72-5  
6370-08-7  
6370-43-0  
637-12-7  
6372-96-9  
6373-74-6  
63-74-1  
6375-55-9  
637-87-6  
637-88-7  
63800-37-3  
638-02-8  
6381-59-5  
6381-79-9  
6381-92-6  
638-23-3  
63-84-3  
63843-89-0  
638-45-9  
6386-38-5  
6388-26-7  
6388-26-7, 25180-48-7

ir-spectra.com IS NIR Spectra

---

**IS NIR Spectra – Complete Collection**

6049 NIR spectra

Copyright © NICODOM 2012

NICODOM Ltd., Hlavni 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

[ftirspectra@gmail.com](mailto:ftirspectra@gmail.com)

63-89-8  
63-91-2  
639-58-7  
6397-02-0  
64-02-8  
64044-51-5  
640-68-6  
6407-99-4  
6408-50-0  
6408-57-7  
6408-78-2  
6410-30-6  
6410-32-8  
6410-38-4  
64-10-8  
6416-57-5  
6416-57-5  
64-17-5  
64-18-6  
64-20-0  
6420-33-3  
6420-65-1  
64211-46-7  
642-31-9  
6424-75-5  
6424-76-6  
6424-85-7  
6424-85-7, 16247-34-3, 6247-34-3  
6428-31-5, 7518-68-5  
64294-91-3  
64294-91-3  
64338-16-5  
64365-13-5  
643-79-8  
6439-53-8  
644-08-6  
6441-82-3  
6441-91-4  
6441-93-6  
644-62-2  
64485-93-4  
645-49-8  
6459-94-5  
646-06-0

ir-spectra.com IS NIR Spectra

---

**IS NIR Spectra – Complete Collection**

6049 NIR spectra

Copyright © NICODOM 2012

NICODOM Ltd., Hlavní 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

[ftirspectra@gmail.com](mailto:ftirspectra@gmail.com)

646-13-9  
646-31-1  
64-69-7  
6471-09-6  
6471-49-4  
6471-51-8  
6472-38-4  
6473-13-8  
64-75-5  
64-77-7  
6483-64-3  
6483-73-4  
6486-23-3  
6487-48-5  
6493-05-6  
6493-06-7#  
6493-07-8#  
650-51-1  
6505-45-9  
65140-91-2  
65-19-0  
65212-77-3  
652-29-9  
65-23-6  
6527-70-4  
65277-42-1  
65294-17-9  
65-31-6  
6533-73-9  
6535-42-8  
6535-46-2  
6539-67-9, 4988-30-1  
65423-44-1  
65447-77-0  
65-46-3  
6553-96-4  
6556-12-3  
65-61-2  
657-27-2  
65753-47-1  
657-84-1  
65800-06-8  
65-85-0  
65853-65-8

ir-spectra.com IS NIR Spectra

---

**IS NIR Spectra – Complete Collection**

6049 NIR spectra

Copyright © NICODOM 2012

NICODOM Ltd., Hlavni 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

[ftirspectra@gmail.com](mailto:ftirspectra@gmail.com)

65997-06-0  
6600-40-4  
66023-94-7  
6602-54-6  
660392-09-6  
660392-10-9  
660392-11-0  
661-19-8  
66-22-8  
66-23-9  
66-27-3  
66335-38-4  
66376-36-1  
6642-29-1  
66449-66-9  
66455-14-9  
6651-34-9  
66516-09-4  
665-66-7  
6657-37-0, 6021-61-0  
66-71-7  
6674-22-2  
66-79-5  
6683-19-8  
66-84-2  
66851-10-3  
66904-06-1  
66-99-9  
67-03-8  
670-54-2  
6706-59-8  
6706-96-3  
67-20-9  
67209-29-4  
673-32-5  
6737-42-4  
67398-15-6  
67421-02-7  
67-43-6  
67462-91-3  
67-52-7  
67-56-1  
67577-84-8  
67-63-0

ir-spectra.com IS NIR Spectra

---

**IS NIR Spectra – Complete Collection**

6049 NIR spectra

Copyright © NICODOM 2012

NICODOM Ltd., Hlavní 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

[ftirspectra@gmail.com](mailto:ftirspectra@gmail.com)

67-64-1  
676448-30-9  
67-66-3  
67674-67-3  
67-68-5  
677335-18-1  
677335-21-6  
67762-27-0  
67845-93-6  
67861-96-5  
6786-83-0  
67876-63-5  
67906-42-7  
67938-21-0  
68016-36-4  
680-31-9  
68-04-2  
68-05-3  
68071-75-0  
68100-20-9  
68-11-1  
681132-21-8  
68-12-2  
68130-14-3  
68155-09-9  
68155-92-0  
681-84-5  
68186-85-6  
68186-88-9  
68186-95-8  
68187-11-1  
68187-11-1  
68187-12-2  
68187-15-5  
68187-49-5  
68187-51-9  
68-19-9  
68201-46-7  
68213-23-0  
682-30-4  
68247-76-7  
6834-92-0  
6837-87-2  
68399-81-5

ir-spectra.com IS NIR Spectra

---

**IS NIR Spectra – Complete Collection**

6049 NIR spectra

Copyright © NICODOM 2012

NICODOM Ltd., Hlavní 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

[ftirspectra@gmail.com](mailto:ftirspectra@gmail.com)

68399-99-5  
68412-54-4  
68412-74-8  
68440-66-4  
6846-33-9  
68476-25-5  
68476-50-6  
68511-62-6  
68516-73-4  
68541-71-9  
6856-08-2  
6856-31-1  
685-73-4  
68584-77-0  
685-87-0  
685-91-6  
68610-51-5  
687-47-8  
68784-69-0?  
68859-25-6, 1306-23-6  
688-73-3  
68-89-3  
68912-12-9  
68915-31-1  
68916-42-7  
68917-82-8  
68937-54-2  
68937-55-3  
68938-54-5  
68-94-0  
68954-91-6, 115965-78-1  
68954-91-6, 58450-52-5  
68957-94-8  
68988-92-1  
68990-67-0  
69011-20-7  
69071-34-7  
69-09-0  
691397-13-4  
6921-29-5  
69219-70-1  
69249-61-2  
69257-04-1  
693-23-2

ir-spectra.com IS NIR Spectra

---

**IS NIR Spectra – Complete Collection**

6049 NIR spectra

Copyright © NICODOM 2012

NICODOM Ltd., Hlavni 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

[ftirspectra@gmail.com](mailto:ftirspectra@gmail.com)

693-36-7  
693-36-7  
693-58-3  
693-98-1  
69458-20-4  
69458-70-4  
695-06-7  
69-52-3  
69-53-4  
6959-47-3  
696-07-1  
696-54-8  
69-65-8  
696-62-8  
69671-47-2  
69-72-7  
69-79-4  
6985-92-8  
698-63-5  
69-89-6  
69899-68-9  
6992-11-6  
69-96-5  
699-98-9  
700-58-3  
70-07-5  
700842-60-0  
70-11-1  
70-18-8  
70209-87-9  
70210-21-8, 12237-97-7  
70-26-8  
70294-19-8  
70294-19-8  
70321-86-7  
70-34-8  
70356-09-1  
70-47-3  
7050-19-3  
7051-16-3  
705-31-7  
70-55-3  
70563-36-9  
70616-89-6

ir-spectra.com IS NIR Spectra

---

**IS NIR Spectra – Complete Collection**

6049 NIR spectra

Copyright © NICODOM 2012

NICODOM Ltd., Hlavni 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

[ftirspectra@gmail.com](mailto:ftirspectra@gmail.com)

70694-55-2  
70-70-2  
70729-65-6, 61355-92-8  
70782-12-6  
70851-08-0  
70851-34-2  
7085-55-4  
7087-68-5  
71-00-1  
71011-26-2  
71011-27-3  
71125-38-7  
71130-06-8  
712304-21-7  
71-23-8  
71243-80-6  
71-30-7  
71320-77-9  
71396-29-7  
7142-69-0  
71-43-2  
7144-08-3#  
7144-37-8  
71-48-7  
71538-26-6, 476-41-5  
71-55-6  
71637-34-8  
71759-92-7  
7177-48-2  
71775-27-4  
7179-49-9  
71819-75-5  
7182-77-6  
71838-82-9  
71872-43-0  
71872-49-6  
71872-67-8  
71878-19-8  
71902-08-4  
71963-77-4  
72-14-0  
72152-49-9  
72-17-3  
72-18-4

ir-spectra.com IS NIR Spectra

---

**IS NIR Spectra – Complete Collection**

6049 NIR spectra

Copyright © NICODOM 2012

NICODOM Ltd., Hlavni 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

[ftirspectra@gmail.com](mailto:ftirspectra@gmail.com)

72-19-5  
72-19-5  
7220-79-3  
72287-26-4  
7232-21-5  
7240-90-6  
72428-99-0  
72-48-0  
72-57-1  
7275-43-6  
72828-03-6  
72829-25-5  
72906-25-3  
72914-19-3  
7292-14-0  
72956-09-3  
72968-14-0  
72968-82-2  
730-11-0  
73049-39-5  
7311-70-8  
7320-34-5  
73-22-3  
73-24-5  
73270-46-9  
73-32-5  
73398-77-3  
73-40-5  
736-30-1  
7365-45-9  
7369-51-9  
73713-75-4  
7379-35-3  
738-70-5  
73987-94-7  
7400-08-0  
74-11-3  
7411-49-6  
741240-85-7  
741258-58-2  
7425-14-1  
7429-37-0  
74311-50-5  
74-31-7

ir-spectra.com IS NIR Spectra

---

**IS NIR Spectra – Complete Collection**

6049 NIR spectra

Copyright © NICODOM 2012

NICODOM Ltd., Hlavni 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

[ftirspectra@gmail.com](mailto:ftirspectra@gmail.com)

74336-59-7  
7434-89-1  
74-39-5  
7440-41-7  
7446-07-3#  
7446-18-6  
7446-20-0  
7446-27-7  
7446-31-3  
7447-40-7  
74568-07-3  
7467-90-5  
74-79-3  
74842-55-0  
7487-88-9  
74-88-4  
7488-54-2  
7488-55-3  
74-89-5  
7491-74-9  
74-94-2  
74-95-3  
749836-65-5  
75-04-7  
75-05-8  
75-07-0  
75-09-2  
75-11-6  
75-12-7  
75-13-8  
75199-01-8  
75-26-3  
7529-22-8  
75-30-9  
75330-75-5  
7534-94-3  
75-36-5  
754-10-9  
75-47-8  
7550-35-8  
75-52-5  
75-57-0  
75601-68-2  
75-64-9

ir-spectra.com IS NIR Spectra

---

**IS NIR Spectra – Complete Collection**

6049 NIR spectra

Copyright © NICODOM 2012

NICODOM Ltd., Hlavní 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

[ftirspectra@gmail.com](mailto:ftirspectra@gmail.com)

75-65-0  
756-79-6  
75708-92-8  
75-75-2  
7576-65-0, 10319-14-9  
7580-46-3  
7580-67-8  
75-82-1  
75-84-3  
7585-39-9  
75-85-4  
75-86-5  
75-89-8  
75920-10-4  
759-24-0  
7598-70-1  
760157-58-2  
7601-89-0  
76-03-9  
76-09-5  
76095-16-4  
76168-75-7  
76199-85-4  
762-04-9  
76-22-2  
762-42-5  
76-24-4  
76-25-5  
7631-90-5  
7631-94-9  
763-69-9  
76-38-0  
76420-72-9  
764-35-2  
7643-75-6  
764-42-1  
7646-93-7  
7647-14-5  
76486-10-7  
765-13-9  
7651-99-2  
765-30-0  
76-54-0  
765-70-8

ir-spectra.com IS NIR Spectra

---

**IS NIR Spectra – Complete Collection**

6049 NIR spectra

Copyright © NICODOM 2012

NICODOM Ltd., Hlavni 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

[ftirspectra@gmail.com](mailto:ftirspectra@gmail.com)

76-59-5  
7659-95-2  
76-60-8  
76-61-9  
7665-72-7  
7677-24-9  
767-91-9  
7681-38-1  
7681-82-5  
76824-35-6  
768-33-2  
76-83-5  
76-86-8  
76-87-9  
768-90-1  
768-94-5  
768-95-6  
769174-08-5  
76-93-7  
769-42-6  
76947-02-9  
76947-39-2  
770-35-4  
7704-34-9  
77-09-8  
771-50-6  
771-60-8  
7720-78-7  
7722-76-1  
7722-88-5  
772-38-3  
7726-03-6  
7727-21-1  
7727-54-0  
77-28-1  
7732-18-5  
7733-02-0  
77337-73-6  
77-36-1  
773-76-2  
77393-49-8  
77-40-7  
77-48-5  
77538-14-8

ir-spectra.com IS NIR Spectra

---

**IS NIR Spectra – Complete Collection**

6049 NIR spectra

Copyright © NICODOM 2012

NICODOM Ltd., Hlavni 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

[ftirspectra@gmail.com](mailto:ftirspectra@gmail.com)

7757-69-9  
7757-79-1  
7757-82-6  
7757-93-9  
7758-09-0  
7758-11-4  
7758-16-9  
7758-23-8  
7758-29-4  
77-58-7  
7758-87-4  
7758-99-8  
7759-02-6#  
77-61-2  
7761-88-8  
77-62-3  
776-74-9?  
77-67-8  
77-71-4  
7772-98-7  
7772-99-8  
77-73-6  
7775-14-6  
7775-27-1  
7778-50-9  
7778-53-2  
7778-74-7  
7778-77-0  
7778-80-5  
77804-81-0  
7782-61-8  
7782-68-5  
7782-87-8  
7783-00-8  
7783-03-1#  
7783-13-3  
7783-20-2  
7783-28-0  
7783-34-8  
7783-35-9  
7783-36-0  
7783-83-7  
7783-85-9  
7784-24-9

ir-spectra.com IS NIR Spectra

---

**IS NIR Spectra – Complete Collection**

6049 NIR spectra

Copyright © NICODOM 2012

NICODOM Ltd., Hlavní 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

[ftirspectra@gmail.com](mailto:ftirspectra@gmail.com)

7784-26-1  
7784-31-8  
7784-46-5  
77-85-0  
7785-20-8  
7785-26-4  
7785-87-7  
7786-17-6  
7786-30-3  
7787-56-6  
7787-68-0  
7788-98-9  
7788-99-0  
7789-00-6  
7789-02-8#,13548-38-4  
7789-09-5#  
7789-12-0  
7789-18-6  
7789-29-9  
7789-77-7  
7789-79-9  
7790-28-5  
7790-62-7  
7790-69-4  
77-90-7  
7790-78-5  
7790-98-9  
7791-07-3#  
7791-13-1  
7791-18-6  
7791-20-0  
7791-25-5  
77-92-9  
77-93-0  
77-95-2  
7798-23-4  
779-84-0  
77-98-5  
77987-49-6  
779-90-8  
77-99-6  
780015-20-5  
7803-55-6  
7803-58-9

ir-spectra.com IS NIR Spectra

---

**IS NIR Spectra – Complete Collection**

6049 NIR spectra

Copyright © NICODOM 2012

NICODOM Ltd., Hlavni 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

[ftirspectra@gmail.com](mailto:ftirspectra@gmail.com)

7803-65-8  
7803-68-1  
78-04-6  
78-11-5  
781560-61-0  
78246-49-8  
78-27-3  
78-38-6  
78-39-7  
784-04-3  
78564-86-0  
78564-87-1  
78-67-1  
78-75-1  
78769-04-7  
78-77-3  
78-79-5  
78-82-0  
78-83-1  
78-90-0  
789-25-3  
78-93-3  
78-94-4  
78-95-5  
78-96-6  
78-98-8  
79-01-6  
79-05-0  
79-06-1  
79-07-2  
79-09-4  
791-28-6  
791-31-1  
79-14-1  
791573-87-0  
79-16-3  
79-19-6  
79-24-3  
79-31-2  
793-24-8  
79-33-4  
79-34-5  
79350-37-1  
79-40-3

ir-spectra.com IS NIR Spectra

---

**IS NIR Spectra – Complete Collection**

6049 NIR spectra

Copyright © NICODOM 2012

NICODOM Ltd., Hlavní 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

[ftirspectra@gmail.com](mailto:ftirspectra@gmail.com)

79-41-4  
79-43-6  
79-46-9  
79-74-3  
79771-28-1, 7757-82-6  
79785-97-0  
79785-99-2  
79817-88-2  
79879-38-2  
79887-10-8  
79902-63-9  
79-97-0  
8000-34-8  
8000-48-4  
8000-75-7  
8000-95-1  
8001-22-7  
8001-23-8  
8001-26-1  
8001-30-7  
8001-61-4  
8001-69-2  
8001-75-0  
8001-79-4  
8002-03-7  
8002-09-3  
8002-11-7  
8002-43-5  
8002-74-2  
8003-22-3, 6365-50-0  
8003-58-5  
8003-69-8  
8004-87-3  
8004-92-0, 100208-62-6  
8005-32-1  
8005-52-5, 7248-45-5  
8005-56-9  
80-05-7  
8005-78-5, 104744-50-5, 68425-18-3, 4482-25-1  
8006-40-4  
8006-44-8  
8006-54-0  
8006-64-2  
8006-80-2

---

**IS NIR Spectra – Complete Collection**

6049 NIR spectra

Copyright © NICODOM 2012

NICODOM Ltd., Hlavní 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

[ftirspectra@gmail.com](mailto:ftirspectra@gmail.com)

8006-90-4  
8007-18-9  
8007-46-3  
8007-47-4  
8007-75-8  
80-07-9  
80-08-0  
8008-52-4  
8008-56-8  
80-11-5  
8011-86-7, 6637-87-2, 68155-63-5  
8011-87-8  
8012-89-3  
8012-95-1  
8013-07-8#  
80157-00-2, 93050-80-7  
80157-00-2, 93050-80-7  
80157-00-2, 93050-80-7  
8015-86-9  
8015-86-9  
8015-91-6  
8016-13-5  
8016-37-3  
8020-83-5  
80210-62-4  
80214-83-1  
80223-99-0, 106463-17-6  
8022-87-5  
8024-09-7  
8024-32-6  
8028-48-6  
8028-89-5  
8028-89-5  
8028-89-5  
8028-89-5  
8028-89-5  
8042-47-5  
8047-15-2  
80474-97-1  
8047-76-5  
80-48-8  
8050-09-7#  
8051-02-3  
80-52-4

ir-spectra.com IS NIR Spectra

---

**IS NIR Spectra – Complete Collection**

6049 NIR spectra

Copyright © NICODOM 2012

NICODOM Ltd., Hlavni 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

[ftirspectra@gmail.com](mailto:ftirspectra@gmail.com)

80-56-8  
80-58-0  
80-59-1  
80-62-6  
8063-82-9  
8064-00-4  
80-97-7  
81012-89-7  
81-13-0  
81-16-3  
81-23-2  
81235-32-7  
81-25-4  
813-80-2  
814-49-3  
814-80-2  
81-48-1  
814-89-1  
814-95-9  
815-17-8  
81-54-9  
81-64-1  
816420-72-1  
81-77-6  
818-08-6  
81-81-2  
818-38-2  
81-84-5  
81872-10-8  
81-88-9  
81-88-9  
819045-92-6  
81-93-6  
81988-97-8  
81990-04-7  
82029-77-4  
82-07-5  
82192-57-2  
82199-12-0  
82248-59-7  
82419-36-1  
82451-48-7  
82452-93-5  
82-46-2

ir-spectra.com IS NIR Spectra

---

**IS NIR Spectra – Complete Collection**

6049 NIR spectra

Copyright © NICODOM 2012

NICODOM Ltd., Hlavní 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

[ftirspectra@gmail.com](mailto:ftirspectra@gmail.com)

824-72-6  
824-79-3  
82586-55-8  
82640-04-8, 84449-90-1  
82-71-3  
827-54-3  
82-75-7  
82-76-8  
828-27-3  
828-51-3  
830-03-5  
830-13-7  
83-07-8  
830-96-6  
83-32-9  
83-33-0  
83382-48-3  
83-46-5  
83-56-7  
83567-04-8  
83567-04-8  
83-67-0  
83713-01-3  
83-72-7  
83-73-8  
83834-59-7  
83881-52-1  
83-88-5  
83929-91-3  
83929-92-4  
83949-75-1  
84000-63-5  
84057-84-1  
84-15-1  
84485-00-7  
84-61-7  
84625-61-6  
84632-65-5  
84632-65-5  
84-65-1  
84-66-2  
846-70-8  
84-74-2  
847831-15-6

ir-spectra.com IS NIR Spectra

---

**IS NIR Spectra – Complete Collection**

6049 NIR spectra

Copyright © NICODOM 2012

NICODOM Ltd., Hlavni 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

[ftirspectra@gmail.com](mailto:ftirspectra@gmail.com)

84-80-0  
84843-33-4, 94246-75-0  
84-88-8  
85-01-8  
85-18-7  
85190-63-2  
853-34-9  
853-39-4  
85-41-6  
85-44-9  
85-48-3  
85496-36-2  
855127-25-2  
85536-87-4, 71902-15-3  
85-60-9  
85631-88-5, 83046-38-2  
85666-92-8  
85-68-7  
85721-33-1  
85-73-4  
85828-27-9  
85-83-6  
85-85-8  
85959-60-0  
85959-83-7  
860-22-0  
860-22-0  
86-34-0  
863-57-0  
86386-73-4  
86-48-6  
86-53-3  
86543-85-3  
865451-15-6  
865758-15-2  
86-57-7  
86-65-7  
866-82-0  
866-84-2  
867-13-0  
86-73-7  
867-44-7  
86-74-8  
867-56-1

ir-spectra.com IS NIR Spectra

---

**IS NIR Spectra – Complete Collection**

6049 NIR spectra

Copyright © NICODOM 2012

NICODOM Ltd., Hlavni 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

[ftirspectra@gmail.com](mailto:ftirspectra@gmail.com)

868-14-4  
868-18-8  
868-59-7  
868-85-9  
86937-80-6  
86-98-6  
870-24-6  
87-02-5  
870-93-9  
87-13-8  
871-84-1  
872-05-9  
87-20-7  
87210-32-0,9011-14-7  
87210-32-0,9011-14-7  
87239-81-4  
872469-60-8  
872-50-4  
87269-97-4  
87-32-1  
873-32-5  
87333-19-5  
873-55-2  
873-75-6  
873-83-6  
87-41-2  
874-42-0  
87-48-9  
874-90-8  
87-51-4  
87-62-7  
876301-33-6  
87630-35-1  
87-65-0  
87-66-1  
87-67-2  
87-67-2  
87-69-4  
877-06-5  
877-24-7  
87749-50-6  
87-78-5  
87-79-6  
87-81-0

ir-spectra.com IS NIR Spectra

---

**IS NIR Spectra – Complete Collection**

6049 NIR spectra

Copyright © NICODOM 2012

NICODOM Ltd., Hlavní 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

[ftirspectra@gmail.com](mailto:ftirspectra@gmail.com)

87-85-4  
87-86-5  
87-89-8  
879549-45-8  
87-99-0  
88025-88-1  
88-06-2  
88-12-0  
88-13-1  
88-14-2  
88144-23-4  
88150-47-4  
88-15-3  
88-18-6  
881-90-3  
88-24-4  
88-26-6  
88-27-7  
88385-22-2  
883-93-2  
88544-01-8  
88554-07-8  
88-58-4  
886-65-7  
88-67-5  
886986-20-5  
886986-21-6  
88-72-2  
88-74-4  
88-75-5  
88-88-0  
88-89-1  
88935-43-7  
889462-54-8  
88985-18-6  
88-99-3  
89-05-4  
89073-36-9  
89-25-8  
89382-86-5  
894810-12-9  
89551-24-6  
89-55-4  
89-57-6

ir-spectra.com IS NIR Spectra

---

**IS NIR Spectra – Complete Collection**

6049 NIR spectra

Copyright © NICODOM 2012

NICODOM Ltd., Hlavni 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

[ftirspectra@gmail.com](mailto:ftirspectra@gmail.com)

89599-47-3  
89-62-3  
896-33-3  
89-65-6  
89-68-9  
897027-45-1  
89749-77-9  
89759-80-8  
89-82-7  
89830-32-0  
89-83-8  
89-84-9  
89959-98-8, 12221-69-1  
89-98-5  
9000-01-5  
9000-01-5  
9000-07-1  
9000-11-7  
9000-16-2  
9000-16-2  
9000-55-9  
9000-64-0  
9000-65-1  
9000-69-5  
9000-70-8  
9000-71-9  
9000-71-9  
9001-05-2#  
9001-37-0  
90-01-7  
9002-18-0  
9002-22-6  
9002-26-0  
9002-29-3  
90-02-8  
9002-81-7  
9002-84-0  
9002-86-2  
9002-88-4  
9002-89-5  
9002-89-5  
9002-92-0  
9002-93-1  
9002-93-1

ir-spectra.com IS NIR Spectra

---

**IS NIR Spectra – Complete Collection**

6049 NIR spectra

Copyright © NICODOM 2012

NICODOM Ltd., Hlavni 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

[ftirspectra@gmail.com](mailto:ftirspectra@gmail.com)

9002-93-1, 9016-45-9  
9002-98-6  
9003-01-4  
9003-03-6#  
9003-05-8  
9003-07-0  
9003-07-0  
9003-07-0  
9003-07-0  
9003-09-2#  
9003-11-6  
9003-11-6#  
9003-17-2  
9003-17-2  
9003-17-2  
9003-20-7  
9003-31-0  
9003-39-8  
9003-44-5  
9003-49-0  
9003-54-7  
9003-54-7  
9003-54-7  
9003-55-8  
9003-55-8  
9003-55-8  
9003-55-8  
9003-55-8  
9003-55-8  
9003-56-9  
9003-70-7  
9003-70-7,9003-53-6  
9003-70-7,9003-53-6  
90-04-0  
9004-32-4  
9004-32-4  
9004-34-6  
9004-34-6  
9004-34-6  
9004-35-7  
9004-35-7  
9004-38-0  
9004-53-9  
9004-54-0

ir-spectra.com IS NIR Spectra

---

**IS NIR Spectra – Complete Collection**

6049 NIR spectra

Copyright © NICODOM 2012

NICODOM Ltd., Hlavni 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

[ftirspectra@gmail.com](mailto:ftirspectra@gmail.com)

9004-54-0  
9004-57-3  
9004-58-4  
9004-62-0  
9004-65-3  
9004-65-3  
9004-65-3  
9004-65-3  
9004-65-3  
9004-65-3  
9004-65-3  
9004-67-5  
9004-70-0  
9004-70-0  
9004-70-0  
9004-96-0  
9004-98-2  
9004-99-3  
9005-00-9  
90-05-1  
9005-12-3  
9005-25-8  
9005-25-8  
9005-25-8  
9005-25-8  
9005-25-8  
9005-25-8  
9005-25-8  
9005-25-8  
9005-25-8  
9005-38-3  
9005-46-3  
9005-49-6  
9005-53-2  
9005-64-5  
9005-64-5  
9005-65-6  
9005-65-6  
9005-67-8  
9005-70-3  
9005-71-4  
9005-79-2  
9005-81-6  
9006-03-5#  
9006-52-4

ir-spectra.com IS NIR Spectra

---

**IS NIR Spectra – Complete Collection**

6049 NIR spectra

Copyright © NICODOM 2012

NICODOM Ltd., Hlavni 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

[ftirspectra@gmail.com](mailto:ftirspectra@gmail.com)

9006-65-9  
9007-48-1  
9008-66-6  
9009-54-5  
900-95-8  
9009-90-9  
9010-66-6  
90106-86-8  
9010-79-1  
9010-79-1  
9010-79-1  
90-11-9  
90-12-0  
9012-36-6  
9012-76-4  
90-13-1  
9014-90-8, 37226-44-1, 54018-84-7  
90-15-3  
9016-00-6  
9016-45-9  
90170-93-7  
90-20-0  
90268-23-8  
90-27-7  
90281-17-7  
90295-11-7, 12222-04-7, 63950-02-7  
90-30-2  
90-33-5  
90357-06-5  
90367-28-5  
9037-22-3  
9041-35-4  
9041-56-9  
9042-11-9#  
90-43-7  
90438-79-2  
90-44-8  
9047-09-0  
9048-46-8  
9048-71-9  
9049-11-0  
9049-32-5  
90498-90-1  
9050-36-6

ir-spectra.com IS NIR Spectra

---

**IS NIR Spectra – Complete Collection**

6049 NIR spectra

Copyright © NICODOM 2012

NICODOM Ltd., Hlavni 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

[ftirspectra@gmail.com](mailto:ftirspectra@gmail.com)

9051-57-4  
90597-79-8  
90604-31-2  
90604-90-3  
9060-74-6  
9063-38-1  
906371-67-3  
90-64-2  
90-66-4  
907215-89-8  
90-72-2  
90-80-2  
9083-41-4  
9085-42-1  
90880-77-6  
90-94-8  
91-00-9  
91-01-0  
910-31-6  
91-15-6  
91-16-7  
91-17-8  
91-20-3  
91254-18-1  
912-60-7  
91-30-5  
91374-20-8  
91-40-7  
91435-84-9  
91-56-5  
915-67-3  
915-67-3  
91-57-6  
91585-49-8  
91-63-4  
91-64-5  
91-66-7  
91722-66-6  
917-23-7  
91-75-8  
917-61-3  
919-16-4,6080-58-6  
919-30-2  
92-13-7

ir-spectra.com IS NIR Spectra

---

**IS NIR Spectra – Complete Collection**

6049 NIR spectra

Copyright © NICODOM 2012

NICODOM Ltd., Hlavni 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

[ftirspectra@gmail.com](mailto:ftirspectra@gmail.com)

921-53-9  
921-60-8  
92228-81-4  
92308-56-0  
92-43-3  
92-44-4  
92462-68-5  
92-52-4  
92-66-0  
92-69-3  
92-70-6  
92-71-7  
92-84-2  
928-49-4  
92-87-5  
928-90-5  
92-91-1  
92-94-4  
929-59-9  
93-00-5  
93-02-7  
93-03-8  
93040-61-6  
93-04-9  
93-05-0  
93050-79-4  
93051-41-3  
93051-43-5  
93051-44-6  
93051-44-6  
93051-44-6  
93051-44-6  
93051-45-7  
93-07-2  
93-08-3  
93-09-4  
93-10-7  
931-35-1  
93-14-1  
931-48-6  
93-18-5  
931-88-4  
932-53-6  
93-35-6

ir-spectra.com IS NIR Spectra

---

**IS NIR Spectra – Complete Collection**

6049 NIR spectra

Copyright © NICODOM 2012

NICODOM Ltd., Hlavni 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

[ftirspectra@gmail.com](mailto:ftirspectra@gmail.com)

93390-81-9  
93-40-3  
93407-96-6  
93413-69-5  
934-32-7  
93-46-9  
935-30-8  
93-56-1  
93-58-3  
936101-58-5  
93-61-8  
93623-34-8  
93-67-4  
937-20-2  
93-72-1  
93-89-0  
93951-21-4  
93-98-1  
94035-02-6  
940-71-6  
94-09-7  
94-13-3  
94158-79-9  
94158-79-9  
941-98-0  
94-25-7  
94313-91-4  
943319-37-7  
94-36-0  
94-41-7  
94-44-0  
946616-14-4  
94-67-7  
947-42-2  
94-75-7  
94765-62-5  
94-91-7  
94-93-9  
949-99-5  
95-01-2  
95040-74-8  
95-13-6  
95-14-7  
95-33-0

ir-spectra.com IS NIR Spectra

---

**IS NIR Spectra – Complete Collection**

6049 NIR spectra

Copyright © NICODOM 2012

NICODOM Ltd., Hlavni 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

[ftirspectra@gmail.com](mailto:ftirspectra@gmail.com)

95-45-4  
95464-05-4  
95-47-6  
95481-62-2  
95-48-7  
95-50-1  
95507-31-6  
95-51-2  
95-53-4  
95-54-5  
95-55-6  
95-57-8  
956-02-5  
956-04-7  
95-64-7  
956-48-9  
95-68-1  
95-69-2  
956931-56-9  
95713-52-3  
95716-98-6  
95-74-9  
95-76-1  
95-78-3  
95-80-7  
958-09-8  
95-82-9  
95-84-1  
95-85-2  
95-92-1  
959-23-9  
95-93-2  
959-33-1  
95-95-4  
96-02-6  
961-07-9  
96107-95-8  
96-22-0  
96-23-1  
96-24-2  
96-29-7  
96-33-3  
96352-23-7  
96-41-3

ir-spectra.com IS NIR Spectra

---

**IS NIR Spectra – Complete Collection**

6049 NIR spectra

Copyright © NICODOM 2012

NICODOM Ltd., Hlavni 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

[ftirspectra@gmail.com](mailto:ftirspectra@gmail.com)

96433-09-9  
96478-09-0  
96-48-0  
96-50-4  
96-53-7  
96549-95-0  
96-69-5  
96-76-4  
96-81-1  
97-00-7  
97-02-9  
97-05-2  
97-23-4  
97396-56-0  
975-17-7  
97-53-0  
97553-90-7  
97-59-6  
97-65-4  
976-56-7  
97-67-6  
97753-82-7  
97-88-1  
97-90-5  
97-97-2  
97-99-4  
98-01-1  
980-26-7  
98-03-3  
98-04-4  
98-05-5  
98-11-3  
98185-21-8  
98-29-3  
98-36-2  
98418-47-4  
98432-80-5  
98-47-5  
98-54-4  
98-59-9  
98725-14-5  
98-74-8  
98-80-6  
98-82-8

ir-spectra.com IS NIR Spectra

---

**IS NIR Spectra – Complete Collection**

6049 NIR spectra

Copyright © NICODOM 2012

NICODOM Ltd., Hlavni 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

[ftirspectra@gmail.com](mailto:ftirspectra@gmail.com)

98-86-2  
98-87-3  
98-92-0  
989-38-8  
98-95-3  
98-96-4  
98-97-5  
98-98-6  
99031-78-6  
99035-77-5  
99-04-7  
99-06-9  
99085-43-5  
99-09-2  
99108-56-2  
991-84-4  
99241-24-4  
992-53-0  
992-59-6  
99-26-3  
99270-16-3  
99280-30-5  
99294-93-6  
99300-78-4  
99-33-2  
99-34-3  
99-35-4  
99400-01-8  
99402-80-9  
99-55-8  
99614-02-5  
996-19-0  
996-34-9  
99-65-0  
996-50-9  
99727-20-5  
99-73-0  
997-55-7  
99-76-3  
99-77-4  
99-86-5  
99-87-6  
99-90-1  
99-91-2

ir-spectra.com IS NIR Spectra

---

**IS NIR Spectra – Complete Collection**

6049 NIR spectra

Copyright © NICODOM 2012

NICODOM Ltd., Hlavni 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

[ftirspectra@gmail.com](mailto:ftirspectra@gmail.com)

99922-43-7  
99-92-3  
99-93-4  
99-94-5  
99-96-7  
99-97-8  
99-99-0  
999-97-3

ir-spectra.com IS NIR Spectra

---

**IS NIR Spectra – Complete Collection**  
6049 NIR spectra  
Copyright © NICODOM 2012  
NICODOM Ltd., Hlavni 2727, CZ-14100 Praha 4, Czech Republic, EU  
Tel: +420-281914970, Fax: +420-281914971  
<http://www.ir-spectra.com/>  
[ftir@ftir.cz](mailto:ftir@ftir.cz)  
[ftirspectra@gmail.com](mailto:ftirspectra@gmail.com)