



## NICODOM IR Coatings – List of Spectra

The infrared spectral library „NICODOM IR Coatings“ is a unique collection of 1816 FTIR (ATR) digital spectra of Coatings and related materials.

„NICODOM IR Coatings“ can be purchased separately (1816 spectra) or as a cost effective package of 3216 spectra together with „NICODOM IR Dyes and Pigments“ (1400 spectra)  
This list contents sample information of „NICODOM IR Coatings“ only.

The spectra were collected using the FTIR Spectrometer Nexus 670™ (Thermo) and the Miracle™ single bounce ATR accessory (Pike Technologies) equipped by Si crystal.  
The spectral library features resolution of 4cm<sup>-1</sup> in the spectral range 600-4000 cm<sup>-1</sup>, measurement time 1 minute, the instrument was purged by dried air.

Typical information about samples like Name, Manufacturer, Category (sample type), chemical composition, CAS number is also included (when available).

*Example of included information:*

*Byk 011, Additive, Silicone-Free Defoamer, Byk Chemie*

*Pentaerythritol tetra(3-mercaptopropionate), PTMP, PETMP, Additive, Polymer crosslinking, Bruno Bock, CAS: 7575-23-7*

This library is available as digital searchable library compatible with your spectroscopic software. Also it can be delivered with SEARCH software or in \*.pdf format.

To download the list of spectra from this library, to download the free Nicodomo IR Demo Library (50 spectra) compatible with your software, to check the format compatibility as well as for pricelist, ordering info and information about all NICODOM IR/NIR Libraries and other NICODOM products please visit our webpage.

---

NICODOM IR Coatings, 1816 FTIR spectra

Copyright © NICODOM 2008

NICODOM Ltd., Hlavni 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

## List of Categories

The infrared spectral library „NICODOM IR Coatings“ includes following types of samples.

Acidity regulator  
Acrylate  
Acrylic Copolymer  
Acrylic/alkyd emulsion  
Acrylic/styrene emulsion  
Acrylic/VA dispersion  
Additive for Adhesives  
Additive in alkyd coatings  
Additive Inorganic  
Adhesion Promotor  
Aftertreatment agent for dyeings  
Aliphatic fatty acid modified PU  
Aliphatic PU dispersion  
Aliphatic PU/Acrylic dispersion  
Alkyd  
Alkyd (solvent)  
Alkyd emulsion  
Aminofunctional reactive thinner  
Aminopolysiloxane based softening agent  
Anionic Wetting and dispersing agent  
Anti bleeding agent  
Anti gelling agent  
Anti separation Additive  
Anti settling Additive  
Anti skinning agent  
Antistatic Agent  
Aqueous polyurethane dispersion  
Aromatic fatty acid modified PU  
Aromatic PU/acrylic dispersion  
Binder  
Biodegradable reducing agent  
Butadiene-styrene based dispersion binder  
Carbohydrates  
Castor oil based aliphatic PU dispersion  
Catalyst  
Catalyst for EP  
Catalyst for NCO  
Coalescent  
Combination drier  
Copolymer  
Copolymer/aliphatic PU dispersion  
Crease-prevention agent

---

NICODOM IR Coatings, 1816 FTIR spectra

Copyright © NICODOM 2008

NICODOM Ltd., Hlavní 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

Crosslinker  
Cycloaliphatic EP resin  
Deaerator  
Defoamer  
Dispersing Additive  
Dyeing Auxiliaries  
Dyes/Pigments (see details later in this list)  
Emulsifier  
EP reactive diluent  
Epoxyester emulsion  
Epoxy-silicone resin  
Filler  
Fixing agent  
Flow agent and defoamer  
Flow and Leveling Additive  
Foam Control Agents  
Glue  
Graphitti Spray  
Hair Dye  
Hair Spray  
Hardener  
High Mw Dispersant  
High temperature grease  
Hydrophilic polymer  
Hydrophobizing agent  
Chelating agent  
Improves colorants compatibility  
Improves conductivity  
Improves Scratch Resistance  
Improves wetfastness of dyeings  
Inhibitor and chelating agent  
Interfacial Tension Modifier  
Isocyanate  
Leveling and surface wetting Additive -silicone  
Levelling Additive  
Low foaming dispersing agent, surfactants  
Low MW Epoxy resin  
Melamine Formaldehyde Resin  
Modified acrylic emulsion  
Modified Amine  
Modified polyester resin  
Moisture Scavenger  
Moth and Beetle-proofing Agents for Wool  
Nail Lacquere  
Natural based  
Neutralization Agent  
Offset Print

---

NICODOM IR Coatings, 1816 FTIR spectra

Copyright © NICODOM 2008

NICODOM Ltd., Hlavni 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

Oil  
Pad-stream dyeing of cellulose, Aryl Polyglycol Ether and Alkanesulfonates  
Paint  
Paint Industrial/Commercial  
Pigment wetting Additive  
Plasticizer  
Polyamine adduct  
Polycarbonate  
Polyester  
Polyether Siloxane Flow and Levelling Additive  
Polymer  
Polymer crosslinking  
Polymer dispersant  
Polymeric thickener  
Polyurethane dispersion  
Pot-life retarder  
Powder Coating  
Preservative  
Preservative  
PU dispersion  
PU High Mw Dispersant  
PU thickener  
PU water dispersion  
PU Wetting and dispersing Agent  
Raw Material Ceramics  
Reactive Diluent  
Reactive Flexibizer/Diluent  
Reactive Silicone Levelling Additive  
Resin  
Retarder for dyeing acrylic fibers  
Rheology Additive  
Saturated polyester  
Scouring/wetting agent  
Scratch and slip Additive  
Short Alkyd  
Silane  
Silicone Additive  
Silicone modified alkyd  
Silicone Resin  
Slip Additive  
Slip and anti blocking agent  
Slip and Flow Additive  
Slip, levelling, soft touch  
Softener  
Solvent  
Stabiliser  
Stoving Varnish

---

NICODOM IR Coatings, 1816 FTIR spectra

Copyright © NICODOM 2008

NICODOM Ltd., Hlavní 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

Styren/acrylic  
Substrate Wetting Additive  
Surfactant  
Thickener, Retarder  
Tin catalyst  
Tinting Color  
Toner  
Urethane acrylic copolymer  
UV stabilizer  
Vinylacetate copolymer dispersion  
Water Acrylic Dispersion  
Water Dispersion VAc/Ethylen  
Water Dispersion VAc/VC/Ethylen  
Wax  
Wetting Agent  
Wetting and anti settling Additive  
Wetting and Dispersing Agent  
Wetting and Flow Agent

## **List of Manufacturers**

3M  
Air Products  
AKZO  
Alberdingk Boley  
Albrecht  
Albright&Wilson  
Aldrich  
Alfred Gasterstaedt  
Allied Colors  
Angus  
Arets Graphics  
Arsonsisi  
Asuso  
Austis  
BAL  
Balak  
Balakom  
BASF  
Bauhaus  
Bayer  
Belton  
Berol  
Biopol

---

NICODOM IR Coatings, 1816 FTIR spectra

Copyright © NICODOM 2008

NICODOM Ltd., Hlavní 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

BK-Guilini  
Bondex  
BPB Furmuia  
Brillux  
Bruno Bock  
Buechner  
Büfa  
Byk  
C.H. Erbslöh  
Cardolite  
CDF Chemie  
Cedro  
Celanese  
Ciba  
Clariant  
Clarins  
Cognis  
Collano  
Colorlak  
Cray Valley  
Croda  
CVC  
CWS  
Cytec  
De Rem  
Degronit  
Degussa  
Desowag  
Detecha  
Dinova  
Dior  
DM  
Doerken  
Dow  
DSM  
Du Pont  
Dufa  
Dupont  
DWA  
Dyrup  
Eastman  
Eckart  
Elementis  
Epifanes  
Essence  
Estee Lauder  
Europolveri

---

NICODOM IR Coatings, 1816 FTIR spectra

Copyright © NICODOM 2008

NICODOM Ltd., Hlavní 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

Fingrs  
Flint Group  
Gabiella Salvete  
Garnier  
GE  
Geyeral  
Goldschmidt  
Guerlain  
Hagesan HG  
HB Lak  
Helena Rubinstein  
Henkel  
HET  
Hexion  
Hoechst  
Hoeveling  
Hornbach  
Huber  
Hüls  
Hyundai  
Chanel  
Chemolak  
Chemotex  
ICI  
IGP  
Industrial Copolymers  
Jaeger  
Janecke-Schneemann  
Jansen  
Johnstone  
Jordana  
Jotun Norway  
King Industries  
Kluthe  
Kremer  
Kreul  
Lascaux  
Lechler  
Linden  
Loewe  
Lubrizol  
Marchnon  
Matador Fix  
Maybelline  
Mazda  
Mechelen  
Merck

---

NICODOM IR Coatings, 1816 FTIR spectra

Copyright © NICODOM 2008

NICODOM Ltd., Hlavni 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

Minolta  
Miss Sporty  
Montana Cans  
Morton  
Motip Dupli  
Münzing  
Murexin  
Nanoresins  
Naomi  
Nicca  
Nissan  
Noroxo  
Nuplex  
Oleofina and Synfina  
Osi specialties  
Pentel  
Polymerlatex  
Primalex  
Reinol  
RHG  
Rica  
Ripol  
Rohm&Haas  
Rothirsch  
Rühl  
Sadolin  
Sandoz  
Sephora  
Schoenfeld  
Schülke&Mayr  
Sigma  
Sigma Hungary  
Sika  
Slovchema  
Sokrates  
Solange  
Spangenberg  
SPI-Chem  
Spolchemie  
Subaru  
Suzuki  
Synpo  
Synthesia  
Synthopol  
Tanatex  
Tebas  
Teneco

---

NICODOM IR Coatings, 1816 FTIR spectra

Copyright © NICODOM 2008

NICODOM Ltd., Hlavní 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)



Thor  
Thywissen  
Triga  
Umton  
Univar  
UPPC  
Van Son  
Vettewinkel  
Vianova  
VW  
Wacker  
Wagemakers  
Wella  
Welte  
Weyl  
Wiederhold  
Wilckens  
Worlee  
Yves Saint Laurent

## **List of spectra**

The following list contains only common names. The original list delivered with the library will contain the complete information containing CI name, Trade name, Chemical composition, and Manufacturer (when available).

## **Additives**

Anti Terra - 204  
Anti Terra - 204 (dried)  
Anti Terra - U  
Anti Terra - U (dried)  
Byk - Dynwet 800  
Byk 011  
Byk 011 (dried)  
Byk 012  
Byk 017  
Byk 018  
Byk 019  
Byk 020  
Byk 020 (dried)  
Byk 022

---

NICODOM IR Coatings, 1816 FTIR spectra

Copyright © NICODOM 2008

NICODOM Ltd., Hlavní 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

Byk 024  
Byk 025  
Byk 025 (dried)  
Byk 028  
Byk 037  
Byk 037 (dried)  
Byk 044  
Byk 044 (dried)  
Byk 052  
Byk 052 (dried)  
Byk 055  
Byk 055 (dried)  
Byk 066  
Byk 066 (dried)  
Byk 070  
Byk 070 (dried)  
Byk 077  
Byk 077 (dried)  
Byk 088  
Byk 088 (dried)  
Byk 140  
Byk 140 (dried)  
Byk 141  
Byk 142  
Byk 142 (dried)  
Byk 154  
Byk 154 (dried)  
Byk 220S  
Byk 220S (dried)  
Byk 300 (dried)  
Byk 301  
Byk 301 (dried)  
Byk 302  
Byk 306  
Byk 306 (dried)  
Byk 307  
Byk 310  
Byk 310 (dried)  
Byk 320  
Byk 320 (dried)  
Byk 322  
Byk 323  
Byk 325  
Byk 325 (dried)  
Byk 330  
Byk 330 (dried)  
Byk 331

---

NICODOM IR Coatings, 1816 FTIR spectra

Copyright © NICODOM 2008

NICODOM Ltd., Hlavní 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

Byk 333  
Byk 335  
Byk 335 (dried)  
Byk 337  
Byk 337 (dried)  
Byk 344  
Byk 345  
Byk 345 (dried)  
Byk 346  
Byk 346 (dried)  
Byk 348  
Byk 348 (dried)  
Byk 354  
Byk 354 (dried)  
Byk 355  
Byk 358  
Byk 358 (dried)  
Byk 361  
Byk 373  
Byk 373 (dried)  
Byk 380  
Byk 380 (dried)  
Byk 381  
Byk 381 (dried)  
Byk 390  
Byk 390 (dried)  
Byk 392  
Byk 392 (dried)  
Byk 410  
Byk 410 (dried)  
Byk 420  
Byk 420 (dried)  
Byk 425  
Byk 4510  
Byk 4510 (dried)  
Byk A530  
Byk A530 (dried)  
Byk P 104  
Byk P 104 (dried)  
Byk P-105  
Byk S 706  
Byk S 706 (dried)  
Byketol OK  
Byketol OK (dried)  
Byketol WS  
Byketol WS (dried)  
Disperbyk 108

---

NICODOM IR Coatings, 1816 FTIR spectra

Copyright © NICODOM 2008

NICODOM Ltd., Hlavní 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

Disperbyk 110  
Disperbyk 110 (dried)  
Disperbyk 111  
Disperbyk 111 (dried)  
Disperbyk 116  
Disperbyk 161  
Disperbyk 161 (dried)  
Disperbyk 162  
Disperbyk 162 (dried)  
Disperbyk 163  
Disperbyk 163 (dried)  
Disperbyk 164  
Disperbyk 164 (dried)  
Disperbyk 166  
Disperbyk 166 (dried)  
Disperbyk 167  
Disperbyk 167 (dried)  
Disperbyk 170  
Disperbyk 170 (dried)  
Disperbyk 180  
Disperbyk 180 (dried)  
Disperbyk 181  
Disperbyk 181 (dried)  
Disperbyk 182  
Disperbyk 184  
Disperbyk 184 (dried)  
Disperbyk 190  
Disperbyk 190 (dried)  
Disperbyk 192  
Disperbyk 193  
Disperbyk 193 (dried)  
Disperbyk 194  
Disperbyk 194 (dried)  
Lactimon - WS  
Lactimon - WS (dried)  
Nanobyk 3600  
Nanobyk 3600 (dried)  
Efka 2018 (dried)  
Efka 2021  
Efka 2021 (dried)  
Efka 2023  
Efka 2023 (dried)  
Efka 2025  
Efka 2025 (dried)  
Efka 2028  
Efka 2028 (dried)  
Efka 2038

---

NICODOM IR Coatings, 1816 FTIR spectra

Copyright © NICODOM 2008

NICODOM Ltd., Hlavní 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

Efka 2040  
Efka 2040 (dried)  
Efka 2526  
Efka 2526 (dried)  
Efka 2532  
Efka 2532 (dried)  
Efka 2550  
Efka 2722  
Efka 2722 (dried)  
Efka 3030  
Efka 3031  
Efka 3031  
Efka 3031 (dried)  
Efka 3032  
Efka 3033  
Efka 3033 (dried)  
Efka 3034  
Efka 3034 (dried)  
Efka 3236  
Efka 3239  
Efka 3288  
Efka 3500  
Efka 3500 (dried)  
Efka 3570  
Efka 3570 (dried)  
Efka 3580  
Efka 3580 (dried)  
Efka 3600  
Efka 3772  
Efka 3772 (dried)  
Efka 3778  
Efka 3778 (dried)  
Efka 4046  
Efka 4046 (dried)  
Efka 4047  
Efka 4047 (dried)  
Efka 4060  
Efka 4060 (dried)  
Efka 4400  
Efka 4400 (dried)  
Efka 4401  
Efka 4401 (dried)  
Efka 4402  
Efka 4402 (dried)  
Efka 4403  
Efka 4500  
Efka 4500 (dried)

---

NICODOM IR Coatings, 1816 FTIR spectra

Copyright © NICODOM 2008

NICODOM Ltd., Hlavní 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

Efka 4530  
Efka 4530 (dried)  
Efka 4550  
Efka 4550 (dried)  
Efka 4560  
Efka 4560 (dried)  
Efka 4580  
Efka 4580 (dried)  
Efka 5044  
Efka 5044 (dried)  
Efka 5054  
Efka 5054 (dried)  
Efka 5065  
Efka 5066  
Efka 5066 (dried)  
Efka 5071  
Efka 5071 (dried)  
Efka 5207  
Efka 5244  
Efka 6220  
Efka 6230  
Efka 6700  
Efka 6700 (dried)  
Efka 6780  
Efka 6780 (dried)  
Efka 7022  
Efka 7022 (dried)  
Efka 8400  
Efka 8400 (dried)  
Efka 8503  
Efka 8510  
Efka 8510 (dried)  
Efka 8530  
Efka 8832  
Efka 8832 (dried)  
Efka 9009  
Efka 9009 (dried)  
Additol VXL 4930  
Additol VXL 4930 (dried)  
Additol VXL 6237  
Additol VXL 6237 (dried)  
Additol VXW 4909  
Additol VXW 4909 (dried)  
Additol VXW 4926  
Additol VXW 4973  
Additol VXW 4973 (dried)  
Additol VXW 6200

---

NICODOM IR Coatings, 1816 FTIR spectra

Copyright © NICODOM 2008

NICODOM Ltd., Hlavní 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

Additol VXW 6200 (dried)  
Additol VXW 6208  
Additol VXW 6208 (dried)  
Additol VXW 6209  
Additol VXW 6209 (dried)  
Additol VXW 6211  
Additol VXW 6213  
Additol VXW 6213 (dried)  
Additol VXW 6360  
Additol VXW 6387  
Additol VXW 6396  
Additol VXW 6396 (dried)  
Additol WL 490/50  
Additol WL 490/50 (dried)  
Additol XL 018  
Additol XL 018 (dried)  
Additol XL 123  
Additol XL 123 (dried)  
Additol XL 125  
Additol XL 252  
Additol XL 252 (dried)  
Additol XL 255  
Additol XL 255 (dried)  
Additol XL 270  
Additol XL 297, Oxime Based  
Additol XL 6526  
Additol XW 314  
Additol XW 314 (dried)  
Additol XW 390  
Additol XW 390 (dried)  
Additol XW 395  
Additol XW 395 (dried)  
Abiester 100P  
Acropol 35  
Acropol 35 (dried)  
Amibal 145  
AMP 95  
Color Naturals 4.3 Goldbrown  
Color Naturals 4.3 Goldbrown (after oxidation)  
Color Naturals 4.3 Goldbrown (after oxidation) (dried)  
Color Naturals 4.3 Goldbrown (dried)  
Desmophen VP LS 2973  
Desmophen VP LS 2973 (dried)  
Dibenzoyl peroxide  
Diocetyl Sodium Sulphosuccinate  
Dow Corning 30, Silicone, Propylene carbonate  
Dow Corning 66

---

NICODOM IR Coatings, 1816 FTIR spectra

Copyright © NICODOM 2008

NICODOM Ltd., Hlavni 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

Dow Corning 66  
Edolan XCIC  
Edolan XCIC (dried)  
Ethyl palmitate  
Halox Xtain L44 (Tannin-Stain Blocker)  
Halox Xtain L44 (Tannin-Stain Blocker) (dried)  
Isobutyl Stearate  
Laurex CS, Cetyl Stearyl Alcohol  
Lignine  
Lukas Tinting Color Chrom Yellow Light 9507, Dried  
Lukas Tinting Color Chrom Yellow Light 9507, Liquid  
Lukas Tinting Color Light Ochre Yellow 9510, Dried  
Lukas Tinting Color Umbra Greenish 9522, Dried  
Lukas Tinting Color Umbra Greenish 9522, Liquid  
Maize Starch  
Methyl-Ethyl Keton Peroxide  
Metolat LA 571  
Metolat LA 571 (dried)  
Microcrystalline cellulose  
Monoethanolamine  
Neofix R-800  
Nitrocellulose  
Pentaerytritol tetra mercapto propionate  
Perlit SI/SW  
Perlit SI/SW  
Persoftal LU  
Polyoxyethylene Polymethylene Siloxane DC470A (dried)  
Polyquart H81  
Polyquart H81 (dried)  
Protorez FFO  
Protorez FFO  
Silquest A 1100  
Silquest A 1100 (dried)  
Silquest A 1120  
Silquest A 171  
Silquest A 171 (dried)  
Silquest A 187  
Silquest A 187 (dried)  
Silquest A 189  
Silquest A 189 (dried)  
Syntegal V-7  
Tanede Nova  
Tanede RC  
Tinting Oil Color Brilliant Ochre Light Nr. 1507, Liquid  
Tinting Oil Color Burned Sienna Nr. 1510, Dried  
Tinting Oil Color Burned Sienna Nr. 1510, Liquid  
Tinting Oil Color Chrom Yellow Light Nr. 1504, Dried

---

NICODOM IR Coatings, 1816 FTIR spectra

Copyright © NICODOM 2008

NICODOM Ltd., Hlavní 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)



Tinting Oil Color Chrom Yellow Light Nr. 1504, Liquid  
Tinting Oil Color Ivory Extra Burning, Dried  
Tinting Oil Color Ivory Extra Burning, Liquid  
Tinting Oil Color Real Orange Nr. 1600, Dried  
Tinting Oil Color Real Yellow Light Nr. 1598, Dried  
Tinting Oil Color Real Yellow Light Nr. 1598, Liquid  
Tinting Oil Color Umbra Brown Nr. 1516, Dried  
Tinting Oil Color Umbra Brown Nr. 1516, Liquid  
Tinting Oil Color Umbra Greenish Nr. 1648, Dried  
Tinting Oil Color Umbra Greenish Nr. 1648, Liquid  
Tinting Oil Color Zinc Green Light Nr.1533, Dried  
Tinting Oil Color Zinc Green Light Nr.1533, Liquid  
Tinting Oil Color Zinc Green Middle Nr. 1532, Dried  
Triethanolamine  
Tungophen B  
Zoldine MS Plus  
Dow Corning Z6137  
Dow Corning Z6137  
Lubrizol 2061  
Lubrizol 2061 (dried)  
Lubrizol 2062  
Lubrizol 2062 (dried)  
Lubrizol 2063  
Lubrizol 2063 (dried)  
Solthix 250  
Solthix 250 (dried)  
Ancamide 2353  
Ancamide 2353 (dried)  
Ancamine 2280  
Ancamine 2280 (dried)  
Ancamine 2432  
Ancamine 2432 (dried)  
Vestamin TMD (branched aliphatic diamine)  
Etoxon AF-5  
Etoxon AF-5 (dried)  
Statexan HA  
Statexan HA (dried)  
Statexan PES  
Statexan PES  
Tanastat P  
Tanastat P  
Tanastat PH  
Tanastat PH  
1,4-Diazabicyclo [2.2.2]octan  
Dabco T 12 CL (Dibutyltin Dilaurate)  
DDBSA Acid catalyst (Dodecylbenzene Sulfonic Acid)  
DDBSA Acid catalyst (Dodecylbenzene Sulfonic Acid) (dried)

---

NICODOM IR Coatings, 1816 FTIR spectra

Copyright © NICODOM 2008

NICODOM Ltd., Hlavni 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

DEEA (Diethyl Ethanolamine)  
DMDA (Dimethyldodecyl Amine)  
DMP 30, 2,4,6-tris (dimethylaminomethyl) phenol  
DMP 30, 2,4,6-tris (dimethylaminomethyl) phenol (dried)  
Isononanoic acid  
K Kat XC 5218  
K Kat XC 5218 (dried)  
K-Kat 4205 (pentane 2,4-dione)  
Nacure 5076  
Nacure 5076 (dried)  
Tegokat 226  
Tegokat 226 (dried)  
Versamine EH 50  
Versamine EH 50 (dried)  
Lusolvan FBH (Di-isobutyl ester of a mixture of dicarboxylic acids) (dried)  
Optifilm Enhancer 300 (dried)  
Plastilit 3060 (Poly(propylene glycol alkylphenyl ether) (dried)  
Nofome AF  
Nofome AF  
Tego Airex 1435  
Tego Airex 1435 (dried)  
Tego Airex 902 W  
Tego Airex 902 W (dried)  
Tego Airex 907 W  
Tego Airex 907 W (dried)  
Tego Airex 935  
Tego Airex 935 (dried)  
Tego Airex 980  
Tego Airex 985  
Tego Airex 985 (dried)  
Tego Airex 986  
Tego Airex 986 (dried)  
Agitan 100  
Agitan 100 (dried)  
Agitan 120  
Agitan 120 (dried)  
Agitan 150 (dried)  
Agitan 151  
Agitan 151 (dried)  
Agitan 217  
Agitan 218  
Agitan 232  
Agitan 255  
Agitan 255 (dried)  
Agitan 260  
Agitan 260 (dried)  
Agitan 280

---

NICODOM IR Coatings, 1816 FTIR spectra

Copyright © NICODOM 2008

NICODOM Ltd., Hlavní 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

Agitan 281  
Agitan 281 (dried)  
Agitan 285  
Agitan 285 (dried)  
Agitan 288  
Agitan 288 (dried)  
Agitan 290  
Agitan 295  
Agitan 295 (dried)  
Agitan 296  
Agitan 296 (dried)  
Agitan 305  
Agitan 315  
Agitan 315 (dried)  
Agitan 380  
Agitan 380 (dried)  
Agitan 633  
Agitan 650  
Agitan 700  
Agitan 702  
Agitan 703 N  
Agitan 770  
Agitan E 257 (dried)  
Antimussol WLN  
Antimussol WLN  
Dapro DF 1181  
Dapro DF 1492  
Dapro DF 2162  
Dapro DF 880  
Dapro DF 881 S  
Dapro DF 881 S (dried)  
Dapro DF 900  
Dapro DF 944  
Dapro DF 975  
Dapro NA 1621  
Dapro NA 1621 (dried)  
Dehydran 1208  
Dehydran 1208 (dried)  
Perenol E 1  
Tego Foamex 1495  
Tego Foamex 1495 (dried)  
Tego Foamex 7447 (dried)  
Tego Foamex 8030  
Tego Foamex 8030 (dried)  
Tego Foamex 805  
Tego Foamex 805 (dried)  
Tego Foamex 815 (dried)

---

NICODOM IR Coatings, 1816 FTIR spectra

Copyright © NICODOM 2008

NICODOM Ltd., Hlavní 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

Tego Foamex 825  
Tego Foamex 825 (dried)  
Tego Foamex 830  
Dapro W 95 HS  
Disperse AYD W 22  
Disperse AYD W 22 (dried)  
Dispex N 40  
Dispex N 40 (dried)  
Kortamol NNO  
Lopon 890 (dried)  
Lopon PO  
Lopon PO (dried)  
Metolat FC 514  
Metolat FC 514 (dried)  
Metolat FC 515  
Metolat FC 515 (dried)  
Orotan 681  
Orotan 681 (dried)  
Orotan 850  
Sol Plus D 540  
Sol Plus D 540 (dried)  
Solsperse 27000  
Solsperse 41090  
Zetasperse 1200  
Zetasperse 1200 (dried)  
Zetasperse 1600  
Zetasperse 1600 (dried)  
Zetasperse 2100  
Zetasperse 2300  
Levegal DLP  
Levegal DLP  
Tanasperse OLG  
Tanasperse OLG (dried)  
Additol VXW 4940  
Additol VXW 4940 (dried)  
Additol VXW 6240  
Additol VXW 6240 (dried)  
Astragal PAN  
Astragal PAN (dried)  
Levalin VKU-F  
Levalin VKU-F  
Levogen BF 01  
Levogen BF 01 (dried)  
Mesitol AS  
Mesitol AS  
Mesitol BFN  
Mesitol BFN

---

NICODOM IR Coatings, 1816 FTIR spectra

Copyright © NICODOM 2008

NICODOM Ltd., Hlavní 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

Persoftal L  
Persoftal L  
Aluminium stearate  
Brij 35, Polyoxyethylene (35) Lauryl Ether  
Diethyl adipate  
Diethylene Glycol Monostearate  
Disponil AAP 436  
Disponil AAP 436 (dried)  
Disponil SUS 875 IC (Anionic di-2-ethylhexyl sulphosuccinate, sodium salt)  
Disponil SUS 875 IC (Anionic di-2-ethylhexyl sulphosuccinate, sodium salt) (dried)  
Glycerol  
Glyceryl monostearate  
Linoleic Acid  
Magnesium stearate  
Methyl esters of fatty acid  
Methyl Linoleate  
Methyl Palmitate  
Oleic acid  
Oleyl alcohol  
Palmitic acid  
PEG 400 Monooleate  
PEG 400 Monostearate  
Pentaerythritol Tetraoleate  
Polyglycerol esters of fatty acids  
Polyoxyethene (40) stearate  
Polysorbate 20  
PPG Dioleate  
Silwet L-77  
Slovamine SK-7  
Slovamix PC (dried)  
Slovamix TZ  
Slovamix TZ (dried)  
Slovamix ZF1  
Sodium oleate  
Sodium potassium and calcium salts of fatty acids  
Sorbitan monolaurate  
Sorbitan monooleate  
Sorbitan monopalmitate  
Sorbitan monostearate  
Stearic acid  
Tween 60  
Plexene CA  
Plexene FED  
Plexene FED  
Aluminium potassium sulfate  
Aluminium Sulfate  
Ammonium carbonate

---

NICODOM IR Coatings, 1816 FTIR spectra

Copyright © NICODOM 2008

NICODOM Ltd., Hlavní 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

Ammonium chloride  
Ammonium Phosphate  
Bentonite  
Bentonite  
Bicarbonate of soda  
Calcium aluminium silicate  
Calcium chloride  
Calcium oxide  
Calcium silicate  
Calcium sulfate  
Concrete Screed  
Copper Sulfate  
Glass wool  
Granite Grey, dried  
Gypsum white  
Chamotte  
Kaolin  
Kaolin, Aluminium silicate  
Lime Extra White  
Magnesium carbonate  
Magnesium oxide  
Monopotassium Phosphate  
Monosodium phosphate  
Oxide Aluminium  
Plaster of Paris  
Potassium carbonate  
Potassium hydrogen carbonate  
Potassium silicate  
Sand, Silicon dioxide  
Silica Gel  
Sodium aluminium silicate  
Sodium carbonate  
Sodium metasilicate  
Sodium nitrate  
Sodium orthosilicate  
Sodium tetraborate (Borax)  
Talc  
Tricalcium phosphate  
Bayhydur 2336 (Hydrophylic Isocyanate)  
Bayhydur 3100, based on hexamethylene diisocyanate (dried)  
Bayhydur VP LS 2150, based on Isophorone Diisocyanate (IPDI)  
Bayhydur VP LS 2150, based on Isophorone Diisocyanate (IPDI) (dried)  
Cyclohexyl isocyanate  
Desmodur VPLS 2098  
Desmodur VPLS 2098 (dried)  
Desmodur W (cycloaliphatic diisocyanate) (dried)  
Desmodur XP 2410 (aliphatic polyisocyanate)

---

NICODOM IR Coatings, 1816 FTIR spectra

Copyright © NICODOM 2008

NICODOM Ltd., Hlavní 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

Desmodur XP 2410 (aliphatic polyisocyanate) (dried)  
HDI - monomer Hexamethylene diisocyanate (aged) (dried)  
HDI - monomer Hexamethylene diisocyanate (dried)  
Tego Flow 300  
Tego Flow 300 (dried)  
Tego Flow 425  
Acetyl tributyl citrate  
Dibutylphthalate  
Diisodecyl phthalate  
Dioctyl adipate  
Kodaflex TXIB  
Persoftal NOBLE  
Persoftal NOBLE (dried)  
Persoftal OE  
Persoftal OE (dried)  
Persoftal PEK  
Persoftal PEK (dried)  
Persoftal PLUS  
Acetic acid  
Acticide MKC  
Acticide MKC (dried)  
Acticide SR 2278  
Acticide SR 2278 (dried)  
Benzoic acid  
Eulan SPA 01  
Eulan SPA 01 (dried)  
Formaldehyde  
Naphthalene  
Parmetol DF 18  
Parmetol DF 18 (dried)  
Parmetol DF 35  
Parmetol DF 35 (dried)  
Parmetol HF 25  
Parmetol HF 25 (dried)  
Potassium acetate  
Preventol CI-7  
Preventol CI-8  
Preventol D7 (Isothiazolinones) (dried)  
Preventol MP 260  
Preventol MP 260 (dried)  
Preventol® ON extra  
Triclosan  
Cardolite NC 513  
Cardura E5 (glycidyl neopentanoate)  
Erisys GE-22S (dried)  
ChS Epoxy RR 700 (dried)  
K Flex M B 301

---

NICODOM IR Coatings, 1816 FTIR spectra

Copyright © NICODOM 2008

NICODOM Ltd., Hlavní 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

K Flex M B 301 (dried)  
Polypox 480 (Modified cycloaliphatic polyamine)  
Polypox 480 (Modified cycloaliphatic polyamine) (dried)  
Polypox H 016 (Phenol-free polyamine)  
Polypox H 016 (Phenol-free polyamine) (dried)  
Polypox H 354 L (modified cycloaliphatic polyamine)  
Polypox H 354 L (modified cycloaliphatic polyamine) (dried)  
Polypox R 6 (o-cresyl glycidyl ether)  
Polypox R 6 (o-cresyl glycidyl ether) (dried)  
t-Butyl Glycidyl Ether  
Zoldine RD 20  
Acrysol RM 1020  
Acrysol RM 1020 (dried)  
Acrysol RM 12W  
Acrysol RM 12W (dried)  
Acrysol RM 2020  
Acrysol RM 2020 (dried)  
Acrysol RM 825  
Acrysol RM 825 (dried)  
Collacral PU 85  
Collacral PU 85 (dried)  
Rheolate 266 (dried)  
Rheolate 278  
Rheolate 310  
Rheolate 310 (dried)  
Rohagit® SD 15  
Tafigel AP 1/1  
Tafigel AP 1/1 (dried)  
Tafigel AP10  
Tafigel AP10 (dried)  
Tafigel PUR 40  
Tafigel PUR 45 (dried)  
Tafigel PUR 50  
Tafigel PUR 50 (dried)  
Tafigel PUR 60 (dried)  
Tego Visco Plus 3000  
Tego Visco Plus 3000 (dried)  
Tego Visco Plus 3010  
Tego Visco Plus 3010 (dried)  
Tego Visco Plus 3060 (dried)  
Diadavin ANE  
Diadavin ANE (dried)  
Diadavin DSP  
Diadavin DSP (dried)  
Diadavin EPF 01  
Diadavin EPF 01 (dried)  
Diadavin UN

---

NICODOM IR Coatings, 1816 FTIR spectra

Copyright © NICODOM 2008

NICODOM Ltd., Hlavní 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)



Diadavin UN (dried)  
Tanaterge SD  
Tanaterge SD (dried)  
Edaplan LA 402  
Edaplan LA 402 (dried)  
Edaplan LA 410  
Edaplan LA 411  
Edaplan LA 412  
Edaplan LA 431  
Edaplan LA 431 (dried)  
Tego Glide 406  
Tego Glide 406 (dried)  
Tego Glide 410  
Tego Glide 420  
Tego SMO  
2,6-Di-tert-Butyl-p-Cresol  
4,4-thiobis(t-butyl-m-cresol)  
Antioxidant 2246  
Cyasorb UV 3346  
Didodecyl 3,3'-thiodipropionate, DLTPD  
Dioctadecyl 3,3'-thiodipropionate  
Hostanox 03  
Hostavin N30 (HALS)  
Chimasorb 81, 2-hydroxy-4-n-octyloxybenzophenone (UV-531)  
Chimasorb 944  
Irgafos 168  
Irganox 1010  
Irganox 1010 + Irganox 168  
Irganox 1076  
Irganox 1330  
Irganox 3114  
Leucopur EGM  
Millad 3988  
Polyvinylpyrrolidone  
Tinuvin 622  
Uvinul 4050 H  
Byk - LD N 20048  
Byk - LD N 20048 (dried)  
Dapro S 65  
Dapro S 65 (dried)  
Dapro U 99  
Dapro U 99 (dried)  
2-(2-Hydroxy-5-methyl phenyl) Benzotriazole 97%  
2,2,4,4 -Tetrahydroxy-Benzophone 97%  
2-[3-(2H-Benzotriazol-2-yl)-4-hydroxyphenyl]ethyl methacrylate  
3-(2H-Benzotriazol-2-yl)-4-Hydroxy Phenetyl Alcohol  
Cyasorb UV-416 Light Absorber

---

NICODOM IR Coatings, 1816 FTIR spectra

Copyright © NICODOM 2008

NICODOM Ltd., Hlavní 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

Tinuvin 1130 (dried)  
Tinuvin 123-DW, HALS (dried)  
Tinuvin 123-DW, HALS (water dispersion)  
Tinuvin 292 (dried)  
Tinuvin 477-DW (dried)  
Tinuvin 477-DW (water dispersion)  
Tinuvin 5155  
Abeson S, dodecylbenzenesulfonic acid based  
Alresat KM-40, modified maleic resin  
Dow Corning 67  
Edaplan 472  
Edaplan 472 (dried)  
Edaplan 480  
Edaplan 480 (dried)  
Edaplan 482  
Edaplan VP 470 (dried)  
Erkantol NR  
Erkantol NR (dried)  
Lutensol A8, C12-14 fatty alcohol 8EO  
Metolat FC 355  
Metolat FC 388  
Metolat LA 524  
Metolat LA 524 (dried)  
Nofome BLF 01  
Nofome BLF 01 (dried)  
Tanaterge CSU (dried)  
Tanaterge CSU Liquid  
Tanawet RC-N  
Tanawet RC-N (dried)  
Tego Dispers 710  
Tego Dispers 710 (dried)  
Tego Dispers 725 W  
Tego Dispers 725 W (dried)  
Tego Dispers 735 W  
Tego Dispers 735 W (dried)  
Tego Dispers 745 W  
Tego Dispers 745 W (dried)  
Tego Dispers 750 W  
Tego Dispers 750 W (dried)  
Tego Wet 265  
Tego Wet 265 (dried)  
Tego Wet 270  
Tego Wet 500  
Tego Wet 505  
Tego Wet 510

---

NICODOM IR Coatings, 1816 FTIR spectra

Copyright © NICODOM 2008

NICODOM Ltd., Hlavní 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

## Resin/Polymer

Acramin BS-N  
Acramin BS-N (dried)  
Acronal 290D  
Acronal 500 D  
Acronal A 603 (dried)  
Acronal DS 6219  
Acronal DS 6219 (dried)  
Acronal LR 8977 (dried)  
Acronal LR 8988  
Acronal LR 8998  
Acronal S 533  
Acronal S 533 (dried)  
Acronal S 726  
Acronal S 726 (dried)  
Alberdingk AC 2514 (dried)  
Alberdingk AC 2518 (dried)  
Alberdingk AC 2522 (dried)  
Alberdingk AC 2523  
Alberdingk AC 2523 (dried)  
Alberdingk AC 31  
Alberdingk AC 31 (dried)  
Alberdingk AC 5530  
Alberdingk AC 5530 (dried)  
Alberdingk ACU 2000  
Alberdingk ACU 2000 (dried)  
Alberdingk AS 2688  
Alberdingk AS 2688 (dried)  
Alberdingk CUR 69  
Alberdingk CUR 69 (dried)  
Alberdingk CUR 991 VP  
Alberdingk CUR 991 VP (dried)  
Alberdingk MAC 34 (dried)  
Alberdingk U 520  
Alberdingk U 520 (dried)  
Alberdingk U 800  
Alberdingk U 800 (dried)  
Alberdingk U 9800 VP  
Alberdingk U 9800 VP (dried)  
Alberdingk UC 90  
Alberdingk UC 90 (dried)  
Albidur EP 2240 A  
Albidur EP 2240 A (dried)  
Albiflex 348

---

NICODOM IR Coatings, 1816 FTIR spectra

Copyright © NICODOM 2008

NICODOM Ltd., Hlavní 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

Alkydal F 26  
Alkydal F 26 (dried)  
Alnovol PN-320, Phenol Novolak  
Alpex CK-450, Cyclized rubber  
Ancamide 2445  
Ancamide 2445 (dried)  
Ancamide 2573  
Ancamide 2573 (dried)  
Asphalt  
Bakelite  
Bay 2427  
Bay 2427 (dried)  
Bay 2470  
Bay 2470 (dried)  
Bay 2500  
Bay 2500 (dried)  
Bay XP 2593/1  
Bayderm Finish 90 UD  
Bayderm Grund  
Bayderm Grund (dried)  
Baypret USV  
Baypret USV  
Carboset CR 717  
Carboset CR 717 (dried)  
Cardolite LX 5206  
Cardolite LX 5206 (dried)  
Cardolite NC 556X80  
Cardolite NC 556X80 (dried)  
Cardolite NC 562  
Cardolite NC 562 (dried)  
Cellulose Acetobutyrate, Celidor  
Colophony  
Colophony, Amber  
Colophony, Rosin  
Copolymer ABS type  
Copolymer EPDM Type  
Copolymer EVA type  
Copolymer SAN type  
Crodasil SC 650W  
Crodasil SC 650W (dried)  
Crosslinker CX 100 (polyfunctional aziridine) (dried)  
Cyclohexanon Formaldehyde resin  
Cymel 301  
Cymel 301 (dried)  
Cymel 327 (dried)  
Dammar  
Dammar Aged

---

NICODOM IR Coatings, 1816 FTIR spectra

Copyright © NICODOM 2008

NICODOM Ltd., Hlavní 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

Daotan VTW 1265/36WA  
Daotan VTW 1265/36WA (dried)  
DER 660 B-80 (dried)  
Desmodur N 3200, hexamethylene diisocyanate based resin  
Desmodur N 3200, hexamethylene diisocyanate based resin (dried)  
Desmodur N 3600, hexamethylene diisocyanate based resin  
Desmodur N 3600, hexamethylene diisocyanate based resin (dried)  
Desmophen 1150  
Desmophen 2009/1  
Desmophen 2249/1  
Desmophen 2249/1 (dried)  
Desmophen 670 (dried)  
Desmophen 680/1  
Desmophen 680/1 (dried)  
Desmophen 800  
Desmophen 870 (dried)  
Desmophen A 160 SN  
Desmophen A 160 SN (dried)  
Desmophen A 365  
Desmophen A 365 (dried)  
Desmophen A VP LS 2350  
Desmophen A VP LS 2350 (dried)  
Desmophen A VPLS 2158  
Desmophen A VPLS 2158 (dried)  
Desmophen A VPLS 2177  
Desmophen A VPLS 2177 (dried)  
Desmophen C1200  
Desmophen VP LS 2008/1  
Desmophen VP LS 2008/1 (dried)  
Dinomin MI-11-1  
Dinomin MI-11-1 (dried)  
Dispercoll U 54  
Dispercoll U 54 (dried)  
Duroftal VPI 2801/78 BAC  
Duroftal VPI 2801/78 BAC (dried)  
Duroftal VPI 2803/78 BAC  
Duroftal VPI 2803/78 BAC (dried)  
Duroxyn VAX 6127W  
Duroxyn VAX 6127W (dried)  
Efka 1500  
Elvacite 2046, Copolymer Butylmethacrylate/Isobutylmethacrylate  
Enorex M 1630AV  
Enorex M 1630AV (dried)  
Enorex Print 2030  
Enorex Print 2030  
Enorex TE 221  
Enorex TE 221 (dried)

---

NICODOM IR Coatings, 1816 FTIR spectra

Copyright © NICODOM 2008

NICODOM Ltd., Hlavní 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

Epalloy 5200 (dried)  
Epon 828  
Epon 828  
Fluorosilicone Rubber  
Fluorotelomer Vidax 1000  
Glyptal Resin  
Halweftal B 32 75 BAC  
Halweftal B 32 75 BAC (dried)  
Hostaflex CN-131, Polyvinyl Acetate/Polyvinyl Chloride  
Hydrocarbon resin  
Hydroxyfunctional Alkyd from talow oil  
Hydroxyfunctional Alkyd from talow oil (dried)  
Chs Alkyd L56 (dried)  
CHS Epoxy 210X75  
ChS Epoxy Resin 590  
Incorez W 2200  
Incorez W 2200 (dried)  
Incorez W 2600  
Incorez W 2600 (dried)  
Incorez W 830/140  
Incorez W 830/140 (dried)  
Incorez W 830/256  
Incorez W 830/256 (dried)  
Joncryl 1915  
Joncryl 1915 (dried)  
Joncryl 8224  
Joncryl 8224 (dried)  
Joncryl 8226  
Joncryl 8226 (dried)  
Joncryl 8284  
Joncryl 8284 (dried)  
Joncryl 8330  
Joncryl 8330 (dried)  
Joncryl HPD 71 E  
Joncryl HPD 71 E (dried)  
K Flex 188  
Keton Resin N  
Laropal A-101, Aldehyde based Resin  
Laropal A-81, Ketone based Resin  
Levaflex EP-390, polypropylene based  
Liocryl AS 680  
Liocryl AS 680 (dried)  
Liopur 2004-119  
Liopur 2004-119 (dried)  
Long Soya Bean Alkyd  
Long Soya Bean Alkyd (dried)  
Long Urethanized Alkyd

---

NICODOM IR Coatings, 1816 FTIR spectra

Copyright © NICODOM 2008

NICODOM Ltd., Hlavní 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

Long Urethanized Alkyd (dried)  
Luhydran A 848S  
Luhydran A 848S (dried)  
Macrynal SM 51ON 60 LG  
Macrynal SM 51ON 60 LG (dried)  
Maincote HG-86ER  
Maincote HG-86ER (dried)  
Maincote PR 71 (dried)  
Manila Copal  
Melamine Resin  
Melamine-formaldehyde Resin  
Melamine-urea-formaldehyde resin  
Melform 200 I  
Methacrylic Acid  
Modocoll EK-20, Cellulose modified  
Morez 300, Styrene/Acrylic Resin  
Mowilith 611  
Mowilith 611 (dried)  
Mowilith 7410 (dried)  
Mowilith 7416  
Mowilith 7416 (dried)  
Mowilith 7460  
Mowilith 7460 (dried)  
Mowilith 760  
Mowilith 760 (dried)  
Mowilith 7714 (dried)  
Mowilith 7717  
Mowilith 7729  
Mowilith 7729 (dried)  
Mowilith DM 2452  
Mowilith® 20, PVAc  
Mowilith® DM 778  
Mowilith® SDM 5  
Mowilith® SDM 5 (dried)  
Mowiol® 4/98  
Mowiol® 4-88, PVA  
MS2A-Resin  
Natural Rubber  
NeoCryl A-1049  
NeoCryl A-1049 (dried)  
NeoCryl A-1091  
NeoCryl A-1091 (dried)  
NeoCryl A-1120  
NeoCryl A-1120 (dried)  
NeoCryl BT-24  
NeoCryl BT-26 (dried)  
NeoCryl SC-160

---

NICODOM IR Coatings, 1816 FTIR spectra

Copyright © NICODOM 2008

NICODOM Ltd., Hlavní 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

NeoCryl SC-160 (dried)  
NeoCryl XK-12  
NeoCryl XK-15  
NeoCryl XK-205  
NeoCryl XK-205 (dried)  
NeoCryl XK-37  
NeoCryl XK-52  
NeoCryl XK-52 (dried)  
NeoCryl XK-62  
NeoCryl XK-62 (dried)  
NeoCryl XK-75  
NeoCryl XK-75 (dried)  
NeoCryl XK-78  
NeoCryl XK-87  
NeoCryl XK-87 (dried)  
NeoCryl XK-90  
NeoCryl XK-95 (dried)  
NeoCryl XK-98  
NeoPac E-106  
NeoPac E-106 (dried)  
NeoPac E-111  
NeoPac E-111 (dried)  
NeoPac E-114  
NeoPac E-114 (dried)  
NeoPac E-121  
NeoPac E-121 (dried)  
NeoPac E-125  
NeoPac E-125 (dried)  
NeoPac E-130  
NeoPac E-130 (dried)  
NeoPac E-180  
NeoRez R-2100  
NeoRez R-2100 (dried)  
NeoRez R-9403  
NeoRez R-9407  
NeoRez R-9407 (dried)  
Paraformaldehyde  
Paraldehyde  
Paraloid B66  
Paraloid B66 (dried)  
Paraloid™ B 48 N  
Paraloid™ B 66  
Paraloid™ B 72  
PEG 7 Glyceryl Cocoate  
PEG-200 Glyceryl Tallowate  
Phenol resin  
Phenolic Resin

---

NICODOM IR Coatings, 1816 FTIR spectra

Copyright © NICODOM 2008

NICODOM Ltd., Hlavní 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)



Phthalopal LR8525, Phthalate based Resin  
Plexigum® PQ 611  
Plexisol® P 550-40  
Poloxamer 105 (PEG-PPG copolymer)  
Poly(acrylonitrile:butadiene)  
Poly(butadiene::styrene)  
Poly(butyl acrylate)  
Poly(butyl methacrylate)  
Poly(ethyl acrylate)  
Poly(ethyl methacrylate)  
Poly(ethylene adipate)  
Poly(ethylene dimethacrylate)  
Poly(ethylene glycol)  
Poly(ethylene glycol) 1000  
Poly(ethylene terephthalate)  
Poly(ethylene:propylene)  
Poly(isodecyl methacrylate)  
Poly(methacrylic acid)  
Poly(methyl methacrylate)  
Poly(methyl methacrylate)::(butyl methacrylate)  
Poly(propylene glycol) 2000  
Poly(styrene:ethylacrylate)  
Poly(styrene-ethylene-butylene)  
Poly(vinyl butyral)  
Poly(vinyl isobutyl ether)  
Poly(vinyl methyl ketone)  
Polyacrylamide  
Polyacrylonitrile  
Polyamide resin  
Polybutadiene  
Polydimethylsiloxane  
Polyethylene  
Polychloroprene  
Polyisobutylene  
Polyphenylene Oxide, Noryl  
Polypropylene  
Polystyrene  
Polyurethane  
Polyvinylbutyral  
Polyvinylchloride - hard  
Polyvinylchloride - soft  
Primal CL-3371 (dried)  
Primal E 2573  
Primal E 2573 (dried)  
Primal E-3188  
Primal E-3188 (dried)  
Primal EP 6060

---

NICODOM IR Coatings, 1816 FTIR spectra

Copyright © NICODOM 2008

NICODOM Ltd., Hlavní 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

Primal EP 6060 (dried)  
Primal HG 800  
Primal WS 24E  
Primal® AC 35 (dried)  
Primal® WS 24  
Primal® WS 24 (dried)  
PUR Desmophan 460  
Regalrez® 1126, Hydrocarbon based Resin  
Resin AFS  
Resin AW-2, Polycyclohexanone based Resin  
Resin CA  
Resin MH  
Resin MH (dried)  
Resin SK, Modified Acetophenone/Formaldehyde Resin  
Resydrol AX 237W  
Resydrol AX 237W (dried)  
Resydrol AY 241W/40WA  
Resydrol AY 241W/40WA (dried)  
Resydrol AY 430W  
Resydrol AY 430W (dried)  
Resydrol AY 586W/38WA  
Resydrol AY 586W/38WA (dried)  
Resydrol UEA 209  
Resydrol UEA 209 (dried)  
Rheovis CR2, Polyalkyl acrylic acid based Resin  
Rhodopas DS 850  
Rhodopas DS 850 (dried)  
Sancure 20041  
Sancure 20041 (dried)  
SCX 910  
SCX 910 (dried)  
SCX 920  
SCX 920 (dried)  
Setal 168 SS 80  
Setal 168 SS 80 (dried)  
Setalux 6741 AQ 43  
Setalux 6741 AQ 43 (dried)  
Setalux 6742 AQ 40  
Setalux 6742 AQ 40 (dried)  
Setalux C 1151 XX 51  
Setalux C 1151 XX 51 (dried)  
Setalux C 1152 XX 51  
Setalux C 1156 XS 60  
Setalux C 1156 XS 60 (dried)  
Setalux C 1182 SS 55  
Setalux C 1182 SS 55 (dried)  
Setaqua 6520

---

NICODOM IR Coatings, 1816 FTIR spectra

Copyright © NICODOM 2008

NICODOM Ltd., Hlavní 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

Setaqua 6520 (dried)  
Setaqua 6755 (dried)  
Setaqua 6756  
Setaqua 6763  
Setaqua 6763 (dried)  
Setaqua 6771  
Setaqua 6771 (dried)  
Setaqua 6774  
Setaqua 6774 (dried)  
Setaqua 6776 (dried)  
Setyrene 78 XS 55  
Setyrene 78 XS 55 (dried)  
Shellack Aged  
Shellack Angelo 4621 Superblonde  
Short Soya Bean Alkyd  
Short Soya Bean Alkyd (dried)  
Silikophen N P 50/300  
Silikopon EF, Silicone-epoxy Resin  
Silikopur 8001 F  
Silikopur 8001 F (dried)  
Sokrat 412  
Styrofan D 422  
Styrofan D 422 (dried)  
Sulryn random copolymer poly(ethylene-co-methacrylic acid), EMAA  
Synocure 220 HS  
Synocure 220 HS (dried)  
Synolac 1333  
Synolac 1333 (dried)  
Synthalan LS 768  
Synthalan LS 768 (dried)  
Teflon  
Umacryl KI, Methylmethacrylate - styrene copolymer  
Uracr ON CY 474 E (dried)  
Uracron CY 474 E  
Uralac SN 805 X-70  
Uralac SN 805 X-70 (dried)  
Uramex MF 821 B (dried)  
Viacryl SC 166  
Viacryl VSC 6279W  
Viacryl VSC 6279W (dried)  
Viacryl VSC 6295W  
Viacryl VSC 6295W (dried)  
Vialkyd AC 406/75X  
Vialkyd AC 406/75X (dried)  
Vialkyd VAF 7200/70 BAC (dried)  
Vinnapas CEP 15  
Vinnapas CEP 15 (dried)

---

NICODOM IR Coatings, 1816 FTIR spectra

Copyright © NICODOM 2008

NICODOM Ltd., Hlavní 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

Vinnapas LL 6484  
Vinyl Short Alkyd  
Vinyl Short Alkyd (dried)  
WorléeCryl 7930  
WorléeCryl 7930 (dried)  
Worleesin MS-240, Maleic resin  
Zytel E103 HSLC NC-10  
Abimal 125D, Glycerol ester + colophony + maleic  
Alkaprene, Chloroprene based Adhesive  
Refill Glue ER-S

## Oil

Castor Oil  
Coleseed Oil  
Dehydrogenated Castor Oil (dried)  
Dow Corning 41  
Fluorolube  
Linseed Oil  
Mineral Oil, light and heavy  
Oil of Pine  
Oil of Thyme  
Bergamot Oil  
Caraway oil  
Cedar Wood Oil  
Cod-Liver Oil  
Eucalyptus Oil  
Lavender Oil  
Marine fish oil  
Sesame Oil  
Yoyoba Oil  
Owatrol Oil  
Paraffin Oil  
Perilla Oil  
Polishing Oil  
Poppy Oil  
Soya Bean Oil  
Tagat CH 60, PEG-60 Hydrogenated Castor Oil  
Tallow Oil  
Tung Oil  
Walnut Oil

## Solvent

1,1,2,2-Tetrachloroethane

---

NICODOM IR Coatings, 1816 FTIR spectra  
Copyright © NICODOM 2008  
NICODOM Ltd., Hlavní 2727, CZ-14100 Praha 4, Czech Republic, EU  
Tel: +420-281914970, Fax: +420-281914971  
<http://www.ir-spectra.com/>  
[ftir@ftir.cz](mailto:ftir@ftir.cz)

1,1,2-Trichlorotrifluoroethane  
1,2-Dichloroethane  
1,6-Hexandiol  
1-Propanol  
Acetone  
Acetonitrile  
Benzene  
Benzyl acetate  
Benzyl alcohol  
Benzyl benzoate  
Butyl acetate  
Carbon Tetrachloride  
Cetyl alcohol  
Cyclododecan  
Cyclohexane  
Diethyl ether  
Diethylene Glycol Dibutyl Ether  
Dichloromethane  
Dimethylsulfoxide  
Dodecane  
Dowanol DPM, Dipropylene Glycol Methyl Ether  
Dowanol EPH, 2-phenoxyethanol  
Dowanol PNP, Propylene glycol n-propyl ether  
Dowanol PPH, Propylene Glycol Phenyl Ether  
Ethanol  
Ethyl acetate  
Ethylbenzene  
Ethylene glycol  
Ethylene glycol dimethyl ether  
Ethylene glycol monoethyl ether  
Heptane  
Hexane  
Hexanol  
Chloroform  
Isobutyl alcohol  
m-Cresol  
N-methylpyrrolidone  
o-Cresol  
o-Xylene  
p-Cresol  
Pentane  
Petroleum ether  
Phenol  
Propylene glycol  
Tetrahydrofuran  
Tetrachloroethylene  
Toluene

---

NICODOM IR Coatings, 1816 FTIR spectra

Copyright © NICODOM 2008

NICODOM Ltd., Hlavní 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

TPM, Tripropylene Glycol (Mono) Methyl Ether  
TPNB, Tripropylene Glycol n-Butyl Ether  
Triethylene Glycol  
Triethylene Glycol Monobutylether  
Water  
Xylenes Mixed

## **Surfactant**

2,4,7,9-Tetramethyl-5-Decyne-4,7-Diol  
Alcopol O-70 PG, Aqueous/PPG contains dioctyl sulfosuccinate sodium  
Altaran S-8  
Altaran S-8 (dried)  
Disponil A 3065 (modified ethoxylated fatty alcohols POE 30)  
Disponil A 3065 (modified ethoxylated fatty alcohols POE 30) (dried)  
Disponil A 4065 EXP (modified ethoxylated fatty alcohols POE 40)  
Disponil A 4065 EXP (modified ethoxylated fatty alcohols POE 40) (dried)  
Disponil A 5560 EXP (modified ethoxylated fatty alcohols) (dried)  
Dow Corning 57  
Dynasylan HS 2926  
Dynasylan HS 2926 (dried)  
Empicol 0216, Oleoyl etherphosphoric acid  
Empilan AMT 2, Alkylamine Ethoxylate  
Empilan BQ 100, Fatty acid ethoxylate  
Empilan BQ4/S, Ethoxylated Propoxylated fatty alcohol  
Empilan KB-1, Ethoxylated C12-C14 1-EO  
Empilan KB2/S, Lauryl ethoxylate (2EO)  
Empilan KM11, Alcohol ethoxylate EO11  
Empilan, Ethoxylated Propoxylated fatty alcohol  
Empiphos Alkylphosphoric acid C16-C18/5632  
Emulgator WNS, Arylpolyglykolether (dried)  
Fluorad FC-129  
Fluowet OTN  
Fluowet OTN (dried)  
Glycerol Trioleate  
Marlophen 87 P  
Nicca Sunsolt LM-7  
Nicca Sunsolt RF-557  
Nicca Sunsolt RF-557 (dried)  
Sodium Dodecyl Sulfate  
Spanscour Spark  
Spanscour Spark (dried)  
Surfynol 440 (dried)  
Tanaterge DRC  
Tanaterge DRC

---

NICODOM IR Coatings, 1816 FTIR spectra

Copyright © NICODOM 2008

NICODOM Ltd., Hlavní 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

TANEDE CL  
Tannex BCA-N  
Tannex BCA-N (dried)  
Tego Twin 4000  
Triton X-100 (dried)  
Zonyl 9361  
Zonyl 9361 (dried)  
Zonyl FS 62  
Zonyl FS 62 (dried)  
Zonyl FSA  
Zonyl FSA (dried)  
Zonyl FSH  
Zonyl FSH (dried)  
Zonyl FSJ  
Zonyl FSJ (dried)  
Zonyl FSN  
Zonyl FSN (dried)  
Zonyl FSO 100  
Zonyl FSO 100 (dried)  
Zonyl FSP  
Zonyl FSP (dried)

## **Wax**

Byk Cera Aquacer 526  
Byk Cera Aquacer 526 (dried)  
Fluowax Pink  
Fluowax Pink (dried)  
Glycerol Tristearate  
Chlorinated paraffin  
Microcrystalline wax  
Persoftal NPF  
Persoftal NPF (dried)  
Triacetin  
Trilaurin  
Trimyristin  
Wax V

## **Paint Industrial/Commercial**

1K Clearcoat  
1K Clearcoat (dried)

---

NICODOM IR Coatings, 1816 FTIR spectra  
Copyright © NICODOM 2008

NICODOM Ltd., Hlavní 2727, CZ-14100 Praha 4, Czech Republic, EU  
Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>  
[ftir@ftir.cz](mailto:ftir@ftir.cz)

Acryl White Universal Paint  
Activ Wood Stain Colorless  
Alkyd Lacquere Transparent  
Amarol Triol Paint for Wood Aldehyde Resin based  
Anticorrosion Primer Paint, Binder part, Dried  
Anticorrosion Primer Paint, Mixed, Dried  
Anticorrosive Priming Paint  
Antikor 200 Paint  
Antirezin Black Paint  
Asphaltum Laquer  
Asphaltum Laquer (dried)  
Balakryl Uni Mat Paint 0750 Orange  
Balakryl Uni Mat Paint 0100 White  
Balakryl Uni Mat Paint 0815 Red  
Balakryl Uni Mat Paint 0510 Light Green  
Balakryl Uni Mat Paint 0660 Ochre  
Balakryl Uni Mat Paint 0220 Brown Light  
Balakryl Uni Mat Paint 0250 Palisander  
Balatix Paint, Tixotropic Enamel  
Base Lissante Sheer Ivory Base Coat (for nails)  
Base Lumere Protective Basecoat (for nails)  
Base Paint for Metall, Glossy part, dried, from 1960  
Base Paint for Metall, Matt Part, Dried, from 1960  
Base Paint, from 1960  
Bisil Light Blue Paint  
Black Matt Paint for Tennis Tables, from 1960  
Black Protective and Isolation Paint, Dried  
Black Protective and Isolation Paint, Liquid  
Black Toner MT 105B  
Bondex Profi Wood Stain  
Bondex Wood Care Product Natural Mahagony  
Bondex Wood Care Product, Dried  
Bondex Wood Care Product, Liquid  
Bootlack, Paint for Boats  
Brass Clear Varnish for Metal, Dried  
Brass Clear Varnish for Metal, Liquid  
Brillux Lacryl Glossy Finish Paint 275 RAL 1015  
Bronze Goddess AJ (for nails)  
C1036 Paint 0000, Nitrocellulose based  
C1037 Celolesk 0000, Paint for wood, nitrocellulose based  
C1038 Celomat Paint 0000, Wood Stain, nitrocellulose based  
C2000 Paint  
C2000 Paint (dried)  
C5000 Sealant 110, Nitrocellulose based  
Celox Paint  
Cewesol Enamel Binder part, Liquid  
Cewesol Enamel Settlement part, Dried

---

NICODOM IR Coatings, 1816 FTIR spectra

Copyright © NICODOM 2008

NICODOM Ltd., Hlavni 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)



Cewesol Enamel Settlement part, Liquid  
Cleaner 205 Bonding Agent  
Clipper I Penguin Yachting Paint  
Color-Matic Spray Paint Red  
Consolan Wood Stain  
Copper Rose Nail Lacquere  
Cova Deep, Wall Paint Interior  
Cova Medium, Wall Paint Interior  
Cova Pastel, Wall Paint Interior  
Cova White, Wall Paint Interior  
Cova Yellow, Wall Paint Interior  
Cowacryl Acrylic Paint Waterthinnable  
Cowidur Glossy Enamel, PU based, Component A+B, Aged A  
Cowidur Glossy Enamel, PU based, Component A+B, Aged B  
Cowidur Glossy Enamel, PU based, Component A+B, Liquid  
Cowidur Hardener  
DD Primer Paint for Yachts Grey  
DD-Hartlack, PU based Lacquer for Yachts, Component A, Dried  
DD-Hartlack, PU based Lacquer for Yachts, Component A, Liquid  
DD-Hartlack, PUR based Lacquer for Yachts, Component B, Dried  
DD-Hartlack, PUR based Lacquer for Yachts, Component B, Liquid  
DD-Hartlack, PUR based Lacquer for Yachts, Components A+B mixed, Dried  
DD-Hartlack, PUR based Lacquer for Yachts, Components A+B mixed, Liquid  
Degrolux Synthetic Resin Based Enamel Paint, from 1960  
Degronit 767 Floor Paint, from 1960  
Deltalac Paint with Epicote, Binder part, Dried  
Deremaqua Paint (dried)  
Dermaline Nail Lacquere  
Deyellowing Serum (for nails)  
Dibenol (binder), from 1960  
Dinocryl Dispersion Paint  
Ducolux Enamel  
Ducolux Metal Paint, Binder part  
Ducolux Metal Paint, Settlement  
Dulux Enamel Plus Red, Paint Silicone based  
Duralit 1516 - 029, Industrial Paint for Metal  
Duralit 1516 - 029, Industrial Paint for Metal (dried)  
Durodur 3057-025, Industrial Paint for Metal  
Durodur 3057-025, Industrial Paint for Metal (dried)  
Eco Super White Enamel Inside, from 1960  
Eggshell Paint  
Ekohost Colorless Matt 0001, Paint urethane-acrylate based, for wood inside  
Ekohost Colorless Semigloss 0000, Paint urethane-acrylate based, for wood inside  
Enamel Spray for Cars Acrylate Combination Blue  
Enamel Spray for Cars Acrylate Combination Blue (dried)  
Enamel Spray for Cars Acrylate Combination Yellow (dried)  
Enamel Spray for Cars Acrylate Green

---

NICODOM IR Coatings, 1816 FTIR spectra

Copyright © NICODOM 2008

NICODOM Ltd., Hlavni 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

Enamel Spray for Cars Acrylate Red  
Enamel Spray for Cars Acrylate Red (dried)  
Epolex Colorless Glossy Paint 0000  
Epoxide Based Lacquere Two componets  
Erbedol Semigloss Paint A + I  
Erbedol Synthetic Resin Based, Impact-Proof Paint  
Eternal Mat Revital Paint 217/201/218 Orange  
Eternal Mat Revital Paint 217/216 Green  
Eternal Mat Revital Paint 206/216 Dark Green  
Eternal Mat Revital Paint 205/216 Blue-Green  
Eternal Mat Revital Paint 205/202/217 Brown-Yellow  
Eternal Mat Revital Paint 201/216/218 Purple  
Eternal Mat Revital Paint 201/213 Grey  
Eternal Paint 01 White  
Eternal Paint 05 Yellow  
Eternal Paint 07 Red-Brown  
Eternal Paint 08 Brick Red  
Eternal Paint 14 Ivory  
Eternal Paint 17 Light yellow  
Eternal Paint 18 Stawberry Red  
Eternal Paint 22 Dark Green  
Euro Alkyd Based Lacquere Glossy  
Express Finish 260 Exotic nail Lacquere  
Exterior Car Car Paint Chip  
Facade Paint TD  
Fasadin, Facade Paint  
Finadur 781 Stoving Varnish (dried)  
Finadur 781, Stoving Varnish  
Finishing Paint Alkyd Based, Matt, P199M Black  
Finishing Paint Alkyd Based, Gloss, P6700 Ochre  
Finishing Paint Alkyd Based, Matt, P100M White  
Finishing Paint Alkyd Based, Gloss, P5321 Green  
Finishing Paint Alkyd Based, Gloss, P6006 Ivory  
Finishing Paint Alkyd based, wood and metal  
Formex Extra Palisander 0280, Anticorrosion Single Layer Paint, for zinc  
Formex Extra Sandy 0670, Anticorrosion Single Layer Paint, for zinc  
Formex Extra White-Grey 0101, Anticorrosion Single Layer Paint, for zinc  
Fortissimo, Wall Paint Interior  
Glasurit Oil Color with Synthetic Resin  
Glasurit Oil Color with Synthetic Resin  
Gloss Paint  
Glossy Colour Enamel Acrylic based, Dried  
Glossy Colour Enamel Acrylic based, Liquid  
Glossy Colour Enamel Alkyd based, Binder part, Dried  
Glossy Colour Enamel Alkyd based, Mixed, Dried  
Goldlack Synthetic Resin based Lacquere for effect, Binder part, Dried  
Goldlack Synthetic Resin based Lacquere for effect, Binder part, Liquid

---

NICODOM IR Coatings, 1816 FTIR spectra

Copyright © NICODOM 2008

NICODOM Ltd., Hlavni 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

Goldlack Synthetic Resin based Lacquere for effect, Mixed, Dried  
Goldlack Synthetic Resin based Lacquere for effect, Mixed, Liquid  
Graphitti Spray Black - White Mate  
Graphitti Spray Gold - White Creme  
Graphitti Spray Molotow Premium - White  
Graphitti Spray Pocket - White  
GS Extra Mat Nail Lacquere  
Gussofarb 66 Boat Paint, Dried  
Gussofarb 66 Boat Paint, Liquid  
H2203 Paint  
Hairstyling Spray Ultralong Lasting  
Hammerite Paint for Metal Green  
Hammerite Paint for Metal Grey  
Hammerite Paint for Metal Silver  
Hammerite Paint for Metal Red  
Hammerite Paint for Metal Blue  
Hammerite Paint for Metal Dark Green  
Herbol - Offenporig, Glazing Paint  
Hetline Paint  
High Build Stain  
High Build Wood Stain  
High Glossy Lacquere White  
Hostagrund Aluminium 0910, Anticorrosion Single Layer Paint  
Hostagrund Black 0199, Anticorrosion Single Layer Paint  
Hostagrund Brown 0243, Anticorrosion Single Layer Paint  
Hostagrund Dark Blue 0455, Anticorrosion Single Layer Paint  
Hostagrund Dark Green 0540, Anticorrosion Single Layer Paint  
Hostagrund Light Blue 0440, Anticorrosion Single Layer Paint  
Hostagrund Red-Brown 0840, Anticorrosion Single Layer Paint  
Hostagrund White 0100, Anticorrosion Single Layer Paint  
Hostagrund Yellow 6200, Anticorrosion Single Layer Paint  
Hyundai Bumper ZW Dark Grey Effect Car Paint Chip  
Icosit Dickschicht 5530 RAL 1014, Anticorrosion Paint  
Icosit Dickschicht 5530 RAL 1014, Anticorrosion Paint  
Impra Color Colorless Wood Stain  
Impredur Glossy Lacquer 840 Alkyd Resin based, Dried  
Impredur Glossy Lacquer 840 Alkyd Resin based, Liquid  
Impredur Silk Glossy Lacquer 880 Alkyd Resin based, Dried  
Impredur Silk Glossy Lacquer 880 Alkyd Resin based, Liquid  
Industral Synthetic Paint for Outside - Brown  
Industrol Black 1999, Paint Alkyd based  
Industrol Blue 4265, Paint Alkyd based  
Industrol Dark Red 8191, Paint Alkyd based  
Industrol Green 5300, Paint Alkyd based  
Industrol Ivory 6003, Paint Alkyd based  
Industrol Light Red 8140, Paint Alkyd based  
Industrol Orange 7551, Paint Alkyd based

---

NICODOM IR Coatings, 1816 FTIR spectra

Copyright © NICODOM 2008

NICODOM Ltd., Hlavni 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

Industrol Red-Brown 8440, Paint Alkyd based  
Industrol Synthetic Paint for Outside - White  
Industrol Yellow 6201, Paint Alkyd based  
Inertol Paint for steel and concrete  
Interior LK Rich Black Car Paint Chip  
Interior Primer Paint, Dried, from 1960  
Interior Primer Paint, Liquid, from 1960  
Interior TP Paint  
Interior ZK Beige Car Paint Chip  
Izoban, Paint for swimming pools  
Izoprim, Wall Paint Interior  
Kronalux Street Marking Paint, Liquid PVC  
Kunstschmiede Mattlack, PU-Alkyd Resin based Lacquer, Settlement part, Dried  
Laksil Semiglossy Paint  
Lazurol Colorless 0000, Alkyd Based Wood Stain with Oil and Wax  
Lazurol Ebony 0099, Alkyd Based Wood Stain with Oil and Wax  
Lazurol Green 0051, Alkyd Based Wood Stain with Oil and Wax  
Lazurol Chestnut 0020, Alkyd Based Wood Stain with Oil and Wax  
Lazurol Mahagony 0080, Alkyd Based Wood Stain with Oil and Wax  
Lazurol Nut 0021, Alkyd Based Wood Stain with Oil and Wax  
Lazurol Palisander 0022, Alkyd Based Wood Stain with Oil and Wax  
Lazurol Pine Dark 0060, Alkyd Based Wood Stain with Oil and Wax  
Lazurol Pine Light 0062, Alkyd Based Wood Stain with Oil and Wax  
Lazurol Sipo 0025, Alkyd Based Wood Stain with Oil and Wax  
Lazurol White 0010, Alkyd Based Wood Stain with Oil and Wax  
Lazurol, Wood Stain, Light Oak  
Le Vernis Multi Eclat 222 (for nails)  
Le Vernis Nail Colour 461 Blue Satin  
Le Vernis Nail Colour 465 Azur  
Long Lasting Nail Lacquere 30  
Luxol Extra Paint S1023, Wood Stain  
Malak Primer Paint for Windows, from 1960  
Malvena, Facade Paint  
Matt, Synthetic Resin Based Paint, from 1960  
Meisterpreis Clear Varnish Seidenmatt with Polyurethane, Gel part, Dried  
Meisterpreis Euro Ventilations, White Paint for Wood, Alkyd Resin Based  
Meisterpreis Super Silky Polyurethane based Paint, Dried  
Meisterpreis Super Silky Polyurethane based Paint, Liquid  
Meisterpreis Tixoton High Build Wood Stain, Dried  
Meisterpreis Tixoton High Build Wood Stain, Liquid  
Meisterpreis Tixoton, High Build Stain  
Meisterpreis White Enamel  
ML AC 300 Stoving Varnish (dried)  
ML AC 300, Stoving Varnish  
Nail Art Multiglitter Tip Painter  
Nail Enamel White  
Nail Hardener

---

NICODOM IR Coatings, 1816 FTIR spectra

Copyright © NICODOM 2008

NICODOM Ltd., Hlavni 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

Nail Lacquere for Childern Strawberry  
Nail Polish No. 48  
Nail Stimulating Serum  
Nissan Infinity Bumper AG05 Medium Grey Effect Car Paint Chip  
O2004 Paint Undercoat for metals, industrial, outside, oil based  
Oil Ground - Preparation for Oil Painting  
Oknobal Palisander 2800, Paint Urethanised Alkyd Glossy  
Oknobal White 1000, Paint Urethanised Alkyd Glossy  
Oldorit Anticorrosive Paint Redbrown, Binder part, Dried  
Oldorit Anticorrosive Paint Redbrown, Mixed, Dried  
Oldorit Anticorrosive Paint Redbrown, Mixed, Liquid  
Original French manicure for Nails  
P05 Paint Remover  
Pinotex Mahagony, Wood Care Product, Dried  
Pinotex Mahagony, Wood Care Product, Liquid  
Plus, Wall Paint Interior  
Polar, Wall Paint Interior  
Pragoprimer, Anticorrosion Paint Primer  
Praktik Gloss Paint 9110  
Praktik Primer Paint 100  
Primalex Procolor S-PL, Paint for wood, base without pigments  
Primer Paint 609 White, Dried, from 1960  
Primer Paint 609 White, Liquid, from 1960  
Primer Paint TP  
Priming Paint For Wood Alkyd Based  
Protective Concrete Paint, Dried  
Protective Concrete Paint, Liquid  
Quick Dry Top Coat For Nails  
Quick Shine For Nails 510 Alin  
Radbal White 1000, Paint Acrylate Based for Radiators  
Rampolux Synthetic Resin Based Paint  
Rampolux Synthetic Resin Based Paint, Dried  
Rampolux Synthetic Resin Based Paint, Liquid  
Rampotex Chlorinated Rubber based Paint  
Rampozet Liquid Plastics  
Ridgefiller (for nails)  
Ripolin Ripoderm Primer Paint for Wood, Dried  
Ripolin Ripoderm Primer Paint for Wood, Settlement part, Liquid  
Ritual Color 13 Fatale (for nails)  
Rostosam Anticorrosion Primer Paint, Dried  
Rostosam Anticorrosion Primer Paint, Liquid  
Rustic Satin High Build Stain, Binder part, Liquid  
Rustic Wood Stain, Liquid, from 1960  
S2028 Paint  
S2070 - Syntecol Primer 0100, Paint for wood, Alkyd based  
S2222 Profi Radiator 1000, Alkyd/urethane based paint, for radiators  
S2819 Paint

---

NICODOM IR Coatings, 1816 FTIR spectra

Copyright © NICODOM 2008

NICODOM Ltd., Hlavni 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

S2852 Krastenol 0100, Paint Polystyrene based for swimming pools  
Sadosol, Wood Care Product  
Sadotec, Wood Care Product  
Seidenmattlack Matte Lacquere for Outside Black  
Semigloss Enamel, Binder part, Liquid  
Semigloss Enamel, Settlement part, Liquid  
Semigloss Enamel, Settlement part, Liquid  
SemiGloss Paint 530  
S-he Stylezone Nail Polish  
Sheer Shimmer Nail Colour 105  
Sheer Strength Nail Foundation  
Sigmastar Alkyde Resin based Paint  
Specolor Priming Undercoat Interior  
Specolux Titan Glossy Enamel  
Specosal Glossy Color Enamel RAL 7035  
Stove Pipe Paint Black Glossy for high Temperatures, Dried  
Stove Pipe Paint Black Glossy for high Temperatures, Liquid  
Stove Pipe Paint Black Glossy, Dried  
Stove Pipe Paint Black Glossy, Liquid  
Strahlenweiss Alkyd Resin Anticorrosion Paint  
Supra Lux Synthetic Resin based Lacquere, Dried  
Supra Lux Synthetic Resin based Lacquere, Liquid  
Supra Oil Colors Glossy Enamel, Dried  
Supra Oil Colors Glossy Enamel, Liquid  
Swing Color, Anticorrosion Paint  
Sylabal Colorless Matt 0001, Paint Alkyd based, transparent, for metal and wood inside  
Teolin High Glossy Lacquere  
Thick Layer Glazing Paint for Wood  
Thin Wood Stain  
Top Coat for Nails  
Topcoat Paint 100 Synthetic  
Topcoat Paint 100 Water Based  
Transparent Sealer for Wood, from 1960  
Trinity Nail Lacquere  
Twinkle Shine Nail Lacquere  
Undercoat Paint  
Underhood Car Paint Chip 6517 Vermillon  
Underhood Car Paint Chip 6519 Black  
Underhood Car Paint Chip 6534 Titanium Frost  
Underhood Car Paint Chip 6579 Grey  
Underhood Car Paint Chip 6726 Red  
Underhood Car Paint Chip 6901 Spruce Green  
Underhood Car Paint Chip 6980 Chestnut  
Underhood Car Paint Chip 7058 Dark SVT Blue  
Underhood Car Paint Chip CS30 Silky White  
Underhood Car Paint Chip Z26 Red  
Underhood Car Paint Chip ZOC Light Grey

---

NICODOM IR Coatings, 1816 FTIR spectra

Copyright © NICODOM 2008

NICODOM Ltd., Hlavni 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

Unibunt Synthetic Resin Based Paint RAL 9001, from 1960  
Unimat Black Paint 1999, Topcoat Matt  
Unimat White Paint 1000, Topcoat Matt  
Universal Brown Paint Urethane/Alkyd based  
Universal Dispersion for Paint to enhance Elasticity  
V 1314 Stoving Varnish  
V 1314 Stoving Varnish (dried)  
V1600 Paint  
V1605 Polyurex Mat 0000, Paint for wood, PU/acrylate based, flat  
V1640 Paint  
V2018 Latex Universal 1000, Wall Paint Interior  
V2020 Paint  
V2025 Okapol 842, Paint for Eaves  
V2026 Paint  
V2045 Balakryl Uni Mat 220, Paint for wood and walls, acrylate based, flat  
V2068 Balakryl Gloss 1000, Paint for wood and walls, acrylate based  
Velvet Orchid Nail Lacquere  
Venticolor Window Enamel  
Ventilux Paint for Windows, Dried  
Ventilux Paint for Windows, Liquid  
Vernis 001 (for nails)  
Vernis 019 (for nails)  
Viscolor Paint, from 1960  
Viscolor Synthetic Resin based Enamel, Binder part, Dried  
VP 169 Stoving Varnish  
VP 169 Stoving Varnish (dried)  
Waterproof Paint for Tiles, Isobutyl methacrylate based  
Weissglut Paint 650  
Wellaflex Hair Spray  
Wilkenit Synthetic Resin Lacquere, Binder part, Dried  
Wilkenit Synthetic Resin Lacquere, Mixed, Dried  
Window Enamel White, from 1960  
Wood Care Product Colorless  
Wood Primer Paint  
Wood Stain TC  
Xylabrilant Wood Stain  
Xyladecor 200 Wood Stain  
Xyladecor Blue Wood Stain  
Xylamon Colorless Wood Worm Killer  
Yacht Systeme, Paint for Yachts

## **Powder Coatings**

Powder Coating 19L, Epoxy polyester, glossy, RAL 1021 (powder)  
Powder Coating 19L, Epoxy polyester, glossy, RAL 3020 (powder)  
Powder Coating 3303, Epoxy polyester, matt, RAL 9005 (powder)

---

NICODOM IR Coatings, 1816 FTIR spectra

Copyright © NICODOM 2008

NICODOM Ltd., Hlavni 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)

Powder Coating 4007, Acrylate, silk glossy, RAL 6005 (cured)  
Powder Coating 4007, Acrylate, silk glossy, RAL 6005 (powder)  
Powder Coating 500, Polyester, silk glossy, transparent (powder)  
Powder Coating 504, Polyester, matt, RAL 9005 (powder)  
Powder Coating 58L, Polyester, silk glossy, RAL 1023 (powder)  
Powder Coating 58L, Polyester, silk glossy, RAL 3002 (powder)  
Powder Coating 58L, Polyester, silk glossy, RAL 6005 (powder)  
Powder Coating 58L, Polyester, silk glossy, RAL 9006 (powder)  
Powder Coating 8009, Polyurethane, glossy, transparent (powder)  
Powder Coating B23, Epoxide, matt, RAL 9005 (cured)  
Powder Coating B23, Epoxide, matt, RAL 9005 (powder)  
Powder Coating B44, Epoxide, matt, RAL 9005 (cured)  
Powder Coating B44, Epoxide, matt, RAL 9005 (powder)  
Powder Coating H50, Epoxypolyester, glossy, transparent (powder)

## **Printing**

Matte Overprint Varnish V5080 (dried), Offset Print Color  
Matte Overprint Varnish V5080, Offset Print Color  
Novaspot Pantone Black (dried), Offset Print Color  
Novavit Bio HKS44K (dried), Offset Print Color  
Novavit Bio HKS44K, Offset Print Color  
Pantone 186GA Red (dried), Offset Print Color  
Pantone 186GA Red, Offset Print Color  
Pantone 8383 Gold (dried), Offset Print Color  
Pantone Black Natural 10009 (dried), Offset Print Color  
Pantone Green Natural 10008 (dried), Offset Print Color  
Pantone Green Natural 10008, Offset Print Color  
Pantone Orange 10021 (dried), Offset Print Color  
Pantone Orange 10021, Offset Print Color  
Pantone Process Blue 10007 (dried), Offset Print Color  
Pantone Process Blue 10007, Offset Print Color  
Pantone Purple Oxy (dried), Offset Print Color  
Pantone Purple Oxy, Offset Print Color  
Pantone Red 10032 (dried), Offset Print Color  
Pantone Red 10032, Offset Print Color  
Pantone Reflex Blue 10006 (dried), Offset Print Color  
Pantone Reflex Blue 10006, Offset Print Color  
Pantone Transparent White (dried), Offset Print Color  
Pantone Transparent White, Offset Print Color  
Pantone Violet Natural 10005 (dried), Offset Print Color  
Pantone Violet Natural 10005, Offset Print Color  
Pantone Warm Red 10001 (dried), Offset Print Color  
Pantone Warm Red 10001, Offset Print Color  
Pantone Yellow Natural 10000 (dried), Offset Print Color  
Pantone Yellow Natural 10000, Offset Print Color

---

NICODOM IR Coatings, 1816 FTIR spectra

Copyright © NICODOM 2008

NICODOM Ltd., Hlavni 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

[ftir@ftir.cz](mailto:ftir@ftir.cz)



Print Varnish Glossy 10LW9500 (dried), Offset Print Color  
Print Varnish Glossy 10LW9500, Offset Print Color  
Turbo Board Plus Cyan 80B32 (dried), Offset Print Color  
Turbo Board Plus Magenta 80B31 (dried), Offset Print Color  
Turbo Board Plus Magenta 80B31, Offset Print Color  
Turbo Board Plus Yellow 80B30 (dried), Offset Print Color  
Turbo Board Plus Yellow 80B30, Offset Print Color

Nicodrom IR Coatings. [www.ir-spectra.com](http://www.ir-spectra.com)

---

NICODOM IR Coatings, 1816 FTIR spectra  
Copyright © NICODOM 2008  
NICODOM Ltd., Hlavní 2727, CZ-14100 Praha 4, Czech Republic, EU  
Tel: +420-281914970, Fax: +420-281914971  
<http://www.ir-spectra.com/>  
[ftir@ftir.cz](mailto:ftir@ftir.cz)