



## NICODOM IR Polymers Plus (3431 Spectra)

The infrared spectral library „NICODOM IR Polymers Plus“ is a unique collection of 3431 FTIR (ATR) digital spectra of polymers, polymers additives and other related compounds used in polymer industry.

Examples of compound types:

- Polymers Additives (Acidity regulator, Inorganics, Adhesion Promotor, Aftertreatment agent for dyeings, Thinner, Wetting and dispersing agent, Anti bleeding agent, Anti gelling agent, Anti separation Additive, Anti settling Additive, Anti skinning agent, Antistatic Agent, Binder, Catalyst, Coalescent, Combination drier, Crease-prevention agent, Crosslinker, Deaerator, Defoamer, Dispersing Additive, Dyeing Auxiliaries, Emulsifier, Filler, Fixing agent, Flow and Leveling Additive, Foam Control Agents, Glue, Hardener, High Mw Dispersant, High temperature grease, Hydrophobizing agent, Chelating agent, Inhibitor and chelating agent, Interfacial Tension Modifier, Levelling Additive, Low foaming dispersing agent, surfactants, Moisture Scavenger, Moth and Beetle-proofing Agents for Wool, Oil, Pigment wetting Additive, Plasticizer, dispersant, Pot-life retarder, Preservative, Reactive Flexibizer/Diluent, Retarder for dyeing acrylic fibers, Rheology Additive, Scouring/wetting agent, Scratch and slip Additive, Slip and anti blocking agent, Slip and Flow Additive, Softener, Stabiliser, Surfactant, Thickener, Retarder, UV stabilizer, Wax, Wetting Agent
- Natural and synthetic polymers (Cellulose derivatives, PES, PEG/PPG, Acrylates, PA, PB, PE, PP, PS, PU, Siloxanes, Halogenated polymers, Vinyl polymers, Resins synthetic and natural, Rubbers, Alkyds...)
- Copolymers (ABS, Acrylates, BS, EPDM, EVA, PEPP, SAN...)
- Monomers

The spectra were collected using the FTIR Spectrometer Nexus 670™ (Thermo) and the Miracle™ single bounce ATR accessory (Pike Technologies) equipped by ZnSe or Si crystal.

The spectral library features resolution of 4cm<sup>-1</sup> in the spectral range 600-4000 cm<sup>-1</sup>, measurement time 1 minute, the instrument was purged by dried air.

Typical information about the sample contains:

Name  
CAS Number  
Chemical Formula/Composition  
Sample Type  
Manufacturer  
Brand Name

This library is available as digital searchable library compatible with your spectroscopic software. Also it can be delivered with SEARCH software or in \*.pdf format.

To download the list of spectra from this library, to download a free demo library compatible with your software, to check the format compatibility as well as for pricelist, ordering info and information about all NICODOM IR/NIR Libraries and other NICODOM products please visit our webpage.

The infrared spectral library package „NICODOM IR Polymers Plus“ is delivered as two separate spectral libraries. This allows you to search in separate libraries, in group of libraries or in all libraries.

“NICODOM IR Polymers” (1615 spectra)

“NICODOM IR Coatings” (1816 spectra)

## Part I Polymers

The Following list contains only the compound names. The original list delivered with the library will contain the complete information containing Compound Name, CAS Number, Chemical Formula, Sample Type (when available).

MP0001 Polyethylene high density, Eraclene MP 90  
MP0002 Polyethylene linear low density, Escorene LLN100 4YB  
MP0003 Polypropylene  
MP0004 Polystyrene  
MP0005 Polystyrene general purpose, Krasten 127/9001  
MP0006 Polyvinylchloride  
MP0007 Teflon  
MP0008 Poly(vinyl acetate)  
MP0009 Poly(vinyl alcohol)  
MP0010 Poly(vinyl butyral)  
MP0011 Cellulose Microcrystalline (Vivapur 105)  
MP0012 Cellophane  
MP0013 Cellulose Acetate  
MP0014 Ethyl Cellulose  
MP0015 Hydroxybutyl Methyl Cellulose  
MP0016 Hydroxypropylmethylcellulose  
MP0017 Methyl Cellulose  
MP0018 Nitrocellulose  
MP0019 Starch  
MP0020 Starch, Maize starch  
MP0021 Starch, Potato Starch  
MP0022 Zein  
MP0023 Leather  
MP0024 Poly(1,4-butylene terephthalate)  
MP0025 Poly(diallylphthalate)  
MP0026 Poly(ethylene terephthalate)  
MP0027 Nylon 6  
MP0028 Nylon 11  
MP0029 Bisphenol-A-polycarbonate  
MP0030 Polyurethane  
MP0031 Phenol-formaldehyde resin  
MP0032 Methacrylate oligomer  
MP0033 Paragum 501 (Polyacrylate)  
MP0034 Poly(2-octyl acrylate)  
MP0035 Poly(acrylamide) #1  
MP0036 Poly(acrylamide) #3  
MP0037 Poly(acrylamide) #2  
MP0038 Poly(acrylic acid) m.w.4000000  
MP0039 Poly(butyl acrylate)  
MP0040 Poly(butyl methacrylate)  
MP0041 Poly(ethyl acrylate) #1  
MP0042 Poly(ethyl methacrylate)  
MP0043 Poly(ethylene dimethacrylate)  
MP0044 Poly(hexadecyl methacrylate)  
MP0045 Poly(isodecyl methacrylate)  
MP0046 Poly(methacrylic acid) #1  
MP0047 Poly(methacrylic acid) #2  
MP0048 Poly(methyl methacrylate)  
MP0049 Poly(methyl methacrylate) m.w. 30000  
MP0050 Polyacrylate

MP0051 Polymethacrylate  
MP0052 Chemlon  
MP0053 Nylon 12 #1  
MP0054 Nylon 12 #2  
MP0055 Nylon 6 #1  
MP0056 Nylon 6 #2  
MP0057 Nylon 610  
MP0058 Nylon 612  
MP0059 Nylon 66  
MP0060 Poly(amide-6,6:amide-6)  
MP0061 Polyamide  
MP0062 Polyamide 6 #1  
MP0063 Polyamide 6 #2  
MP0064 Polyamide 6 #3  
MP0065 Ultramid C-35  
MP0066 Buna CB HC 502 TC, Polymer PB type  
MP0067 Buna CB HX 529 SIC, Polymer PB type  
MP0068 Buna CB HX 529 TC, Polymer PB type  
MP0069 Diene 55 AC, Polymer PB type  
MP0070 Intene 40 A, Polymer PB type  
MP0071 Intene 50 AM, Polymer PB type  
MP0072 Poly(butadiene) (45% cis & trans, 55% vinyl)  
MP0073 Poly(butadiene) (55% cis & trans, 45% vinyl)  
MP0074 Poly(butadiene) (36% cis, 55% trans, 9% vinyl)  
MP0075 Poly(butadiene), phenyl terminated  
MP0076 Polybutadiene CIS  
MP0077 Polybutadiene, Asadene 70 X  
MP0078 Polybutadiene, Asaprene 700 A  
MP0079 Polyisoprene cis #1  
MP0080 Techlen BR 1209, Polymer PB type  
MP0081 Poly(1,3-butylene adipate)  
MP0082 Poly(1,3-butylene endomethylene tetrahydrophthalate)  
MP0083 Poly(1,3-butylene fumarate)  
MP0084 Poly(1,3-butylene hexachlor endomethylene tetrahydrophthalate)  
MP0085 Poly(1,3-butylene isophthalate)  
MP0086 Poly(1,3-butylene phthalate)  
MP0087 Poly(1,3-butylene succinate)  
MP0088 Poly(1,4-butylene adipate)  
MP0089 Poly(1,4-butylene fumarate)  
MP0090 Poly(1,4-butylene hexachlor endomethylene tetrahydrophthalate)  
MP0091 Poly(1,4-butylene succinate)  
MP0092 Poly(1,4-butylene terephthalate) #1  
MP0093 Poly(1,4-butylene tetrahydrophthalate)  
MP0094 Poly(cyclohexylene terephthalate)  
MP0095 Poly(diethylene hexachlor endomethylene tetrahydrophthalate)  
MP0096 Poly(diethylene phthalate)  
MP0097 Poly(diethylene tetrahydrophthalate)  
MP0098 Poly(dipropylene endomethylene tetrahydrophthalate)  
MP0099 Poly(dipropylene hexachlor endomethylene tetrahydrophthalate)  
MP0100 Poly(dipropylene isophthalate)  
MP0101 Poly(ethylene adipate)  
MP0102 Poly(ethylene endomethylene tetrahydrophthalate)  
MP0103 Poly(ethylene hexachlor endomethylene tetrahydrophthalate)  
MP0104 Poly(ethylene isophthalate)  
MP0105 Poly(ethylene maleate)  
MP0106 Poly(ethylene maleinadipate)  
MP0107 Poly(ethylene phthalate)

MP0108 Poly(ethylene succinate)  
MP0109 Poly(ethylene terephthalate) #1  
MP0110 Poly(hexandiol 1,6 maleinadipate)  
MP0111 Poly(isopropylene endomethylene tetrahydrophthalate)  
MP0112 Poly(neopentylene adipate)  
MP0113 Poly(neopentylene fumarate)  
MP0114 Poly(neopentylene hexachlor endomethylene tetrahydrophthalate)  
MP0115 Poly(neopentylene isophthalate)  
MP0116 Poly(neopentylene phthalate)  
MP0117 Poly(neopentylene succinate)  
MP0118 Poly(neopentylene terephthalate)  
MP0119 Poly(neopentylene tetrahydrophthalate)  
MP0120 Poly(propylene adipate)  
MP0121 Poly(propylene endomethylene tetrahydrophthalate)  
MP0122 Poly(propylene fumarate)  
MP0123 Poly(propylene hexachlor endomethylene tetrahydrophthalate)  
MP0124 Poly(propylene isophthalate)  
MP0125 Poly(propylene maleate)  
MP0126 Poly(propylene succinate)  
MP0127 Poly(propylene tertahydrophthalate))  
MP0128 Poly(trimethylpentylen adipate)  
MP0129 Poly(trimethylpentylen endomethylene tetrahydrophthalate)  
MP0130 Poly(trimethylpentylen maleate)  
MP0131 Poly(trimethylpentylen phthalate)  
MP0132 Poly(trimethylpentylen sebacate)  
MP0133 Poly(trimethylpentylen succinate)  
MP0134 Poly(trimethylpentylen tetrahydrophthalate)  
MP0135 Polyester based mixture  
MP0136 Polyester, adipic Acid  
MP0137 Polyester; terephthalate based  
MP0138 Perdeuteropolyethylene  
MP0139 Poly(ethylene:propylene) (60% ethylene)  
MP0140 Polyethylene chlorosulphonated  
MP0141 Polyethylene low density  
MP0142 Polyethylene medium density  
MP0143 Polyethylene oxidized  
MP0144 Polyethylene plasticized #1  
MP0145 Polyethylene plasticized #2  
MP0146 Polyethylene, Daplex BF 4572  
MP0147 Polyethylene, Eltex K 52-05-159  
MP0148 Polyethylene, Engage 8180  
MP0149 Polyethylene, Eraclene 80  
MP0150 Polyethylene, Finathene 3802 Y  
MP0151 Polyethylene, Hostalen GM 7040  
MP0152 Polyethylene, Liten MB 62  
MP0153 Polyethylene, Lupolen 6021 D  
MP0154 Polyethylene, Moplen ZB 5015 F  
MP0155 Polyethylene, Scolefin PE AG 62 BA  
MP0156 Polyethylene, Stavrolen PE 4FE69  
MP0157 Acacia gum  
MP0158 Alkylaryl poly(ethylene glycol)  
MP0159 Carrageenan  
MP0160 Gum arabic  
MP0161 Paraform  
MP0162 Paraformaldehyde  
MP0163 Paraldehyde  
MP0164 PEG 7 Glyceryl Cocoate

MP0165 PEG-200 Glyceryl Tallowate  
MP0166 Poly(1-butene)  
MP0167 Poly(acrylonitrile)  
MP0168 Poly(diallylphthalate) #1  
MP0169 Poly(ethyl acrylate) #2  
MP0170 Poly(ethyl cyanoacrylate)  
MP0171 Poly(m-phenylene isophthalamide)  
MP0172 Poly(N-vinyl carbazole), Luvican  
MP0173 Poly(phenylmethylsiloxane)  
MP0174 Poly(vinyl acetate) #2  
MP0175 Poly(vinyl formal)  
MP0176 Poly(vinyl isobutyl ether)  
MP0177 Poly(vinyl methyl ketone)  
MP0178 Polyether polyol  
MP0179 Polyimide  
MP0180 Polyisobutylene  
MP0181 Polyolefinamid Alkeneamine  
MP0182 Polyoxymethylene #1  
MP0183 Polyoxymethylene #2  
MP0184 Polyvidone 17, Polyvinylpolypyrrolidone  
MP0185 Polyvidone 25, Polyvinylpolypyrrolidone  
MP0186 Polyvidone 60, Polyvinylpolypyrrolidone  
MP0187 Polyvidone 90, Polyvinylpolypyrrolidone  
MP0188 Polyvidone Insoluble, Polyvinylpyrrolidone  
MP0189 HIPS/HDPE alloy #1  
MP0190 HIPS/HDPE alloy #2  
MP0191 HIPS/PP alloy  
MP0192 HIPS/PP/PA6 alloy  
MP0193 Polypropylene with 20% Talc  
MP0194 Polypropylene with rubber content  
MP0195 Polypropylene, atactic #1  
MP0196 Polypropylene, atactic #2  
MP0197 Polypropylene, isotactic #1  
MP0198 Polypropylene, isotactic #2  
MP0199 Polypropylene, Moplen EPQ 30 RF  
MP0200 Polypropylene, Stamyran P 512MN10  
MP0201 Polypropylene+Vistalon  
MP0202 Polypropylene+Vistalon 404, 1:1  
MP0203 Polypropylene+Vistalon 719, 1:1  
MP0204 Polystyrene high impact #1  
MP0205 Polystyrene high impact #2  
MP0206 Polystyrene high impact #3  
MP0207 Polystyrene, dicarboxy terminated  
MP0208 Polystyrene, monocarboxy terminated  
MP0209 Polystyrene, Nope PS-Netsark 336M  
MP0210 Polyurethane #1  
MP0211 Polyurethane #3  
MP0212 Polyurethane foam  
MP0213 Polyurethane resin based polymer #2  
MP0214 Poly(propylene glycol) 500  
MP0215 Bakelite  
MP0216 Bisphenol-A-polycarbonate, 20% glass reinforced  
MP0217 Diglycidyl ether of Bisphenol A  
MP0218 Epon 828  
MP0219 Epoxy Resin  
MP0220 Epoxy Resin Novolac  
MP0221 Epoxy resin, Bisphenol A + epichlorohydrin

MP0222 Hydrocarbon resin  
MP0223 Melamine Resin  
MP0224 Melamine-formaldehyde condensate  
MP0225 Melamine-formaldehyde Resin  
MP0226 Melamine-urea-formaldehyde resin #1  
MP0227 Melamine-urea-formaldehyde resin #2  
MP0228 Phenol resin #1  
MP0229 Phenol resin #2  
MP0230 Phenolic Resin  
MP0231 Phenolic Resol  
MP0232 Phenoxy resin  
MP0233 Phenoxy resin Dow Chemical  
MP0234 Phenoxy resin UCAR  
MP0235 Polyamide resin  
MP0236 Polycarbonate  
MP0237 Polyurethane resin  
MP0238 Polyurethane resin based polymer #1  
MP0239 Versamid 100  
MP0240 Fluorosilicone Rubber  
MP0241 Natural Rubber  
MP0242 Rubber Bike break  
MP0243 Rubber Bike Tyre  
MP0244 Rubber Car Tyre 1  
MP0245 Rubber Car Tyre winter  
MP0246 Rubber Carbon Filled #1  
MP0247 Rubber Carbon Filled #2  
MP0248 Synthetic Rubber  
MP0249 White Rubber  
MP0250 Dimethicone  
MP0251 Dimethyl Polysiloxane  
MP0252 Dimeticone  
MP0253 Poly(dimethylsiloxane) #1  
MP0254 Poly(dimethylsiloxane) #2  
MP0255 Poly(dimethylsiloxane) (n=1,000 ctk) s  
MP0256 Poly(dimethylsiloxane) 10000  
MP0257 Poly(dimethylsiloxane) 60000  
MP0258 Poly(methylphenylsiloxane) fluid  
MP0259 Poly(methylphenylsiloxane) viscose fluid  
MP0260 Polysiloxane  
MP0261 Silicone rubber #1  
MP0262 Silicone rubber #2  
MP0263 Silicone rubber #3  
MP0264 Acrylic Copolymer with Styrene  
MP0265 Lotryl 20MA08  
MP0266 Lotryl 35MA08  
MP0267 Ploxadmer 105 (PEG-PPG copolymer)  
MP0268 Poly(2,6-dimethyl-1,4-phenylene ether) + polystyrene  
MP0269 Poly(acrylonitrile:butadiene)  
MP0270 Poly(acrylonitrile:divinylbenzene)  
MP0271 Poly(acrylonitrile:ethylene:propylene:styrene) #1  
MP0272 Poly(acrylonitrile:ethylene:propylene:styrene) #2  
MP0273 Poly(alpha-methylstyrene-block-butadiene)  
MP0274 Poly(chlorotrifluoroethylene:propene)  
MP0275 Poly(ethylene:tetrafluoroethylene)  
MP0276 Poly(ethylene-chlorotrifluoroethylene)  
MP0277 Poly(methyl methacrylate:butadiene)  
MP0278 Poly(m-phenylene isophthalamide)/(viscose)

MP0279 Poly(styrene:ethylacrylate)  
 MP0280 Poly(styrene:n-butyl methacrylate)  
 MP0281 Poly(styrene-ethylene-butylene)  
 MP0282 Poly(vinylidene chloride:acrylonitrile)  
 MP0283 Poly(vinylidene chloride:vinyl chloride)  
 MP0284 Poly[(2,3-epoxypropyl methacrylate):styrene:(ethylene dimethacrylate)]  
 MP0285 Poly[(acrylic acid):(2-ethylhexyl acrylate):styrene]  
 MP0286 Poly[(acrylic acid):(acrylate):styrene]  
 MP0287 Poly[(acrylic acid):(butyl acrylate):(methacryl amide):styrene]  
 MP0288 Poly[(acrylic acid):(butyl acrylate):2-(hydroxyethyl methacrylate):styrene]  
 MP0289 Poly[(acrylic acid):(butyl acrylate):styrene]  
 MP0290 Poly[(acrylic acid):(ethyl acrylate):(2-hydroxyethyl methacrylate):styrene]  
 MP0291 Poly[(acrylic acid):(maleic anhydride)]  
 MP0292 Poly[(acrylic acid):ethylene], Surlyn  
 MP0293 Poly[(divinyl benzene):styrene]  
 MP0294 Poly[ethylene:(vinyl acetate)]  
 MP0295 Poly[styrene-co-acrylonitrile-co-(divinyl benzene)]  
 MP0296 Polypropylene-graft-maleic anhydride  
 MP0297 Ultramid 6A, Poly[caprolactam:(hexamethylene adipamide)]  
 MP0298 Vinylpyrrolidone-Vinylacetate Copolymer  
 MP0299 Cycolac AM-AS 48879, Copolymer ABS type  
 MP0300 Cycolac G 101, Copolymer ABS type  
 MP0301 Cycolac TIAS, Copolymer ABS type  
 MP0302 Forsan 449/9002, Copolymer ABS type  
 MP0303 Forsan 548/9002, Copolymer ABS type  
 MP0304 Forsan 573, Copolymer ABS type  
 MP0305 Lustran HH 1655, Copolymer ABS type  
 MP0306 Magnum 3404, Copolymer ABS type  
 MP0307 Magnum 3416 SC, Copolymer ABS type  
 MP0308 Magnum 3453, Copolymer ABS type  
 MP0309 Novodur KU 2, Copolymer ABS type  
 MP0310 Novodur P2H AT, Copolymer ABS type  
 MP0311 Novodur P2T, Copolymer ABS type  
 MP0312 Poly(acrylonitrile:butadiene:styrene) #1  
 MP0313 Poly(acrylonitrile:butadiene:styrene) #2  
 MP0314 Poly(acrylonitrile-co-butadiene-co-styrene)  
 MP0315 Polylac 737 PA, Copolymer ABS type  
 MP0316 Sinkral CA24Y, Copolymer ABS type  
 MP0317 Starex HR 0390, Copolymer ABS type  
 MP0318 Starex SM 0150, Copolymer ABS type  
 MP0319 Starex VH 0800, Copolymer ABS type  
 MP0320 Terluran 967 K, Copolymer ABS type  
 MP0321 Terluran KR 2893, Copolymer ABS type  
 MP0322 Terluran KR P 31147, Copolymer ABS type  
 MP0323 Urtal B 23, Copolymer ABS type  
 MP0324 Poly(methyl methacrylate):(butyl methacrylate)  
 MP0325 Poly[(2,3-epoxypropyl methacrylate):(ethylene dimethacrylate)]  
 MP0326 Poly[(acrylic acid):(2-hydroxyethyl methacrylate)]  
 MP0327 Poly[(acrylic acid):(butyl acrylate):(ethyl acrylate):(2-ethylhexyl acrylate):(methacryl amide)]  
 MP0328 Poly[(acrylic acid):(butyl acrylate):(ethyl acrylate)]  
 MP0329 Poly[(acrylic acid):(butyl acrylate):(methacryl amide):(methyl methacrylate)]  
 MP0330 Poly[(acrylic acid):(butyl acrylate):(methyl methacrylate)]  
 MP0331 Poly[(acrylic acid):(butyl acrylate):]  
 MP0332 Poly[(acrylic acid):(ethyl acrylate):(methacryl amide):(methyl methacrylate)]  
 MP0333 Poly[(acrylic acid):(ethyl acrylate)]  
 MP0334 Poly[(acrylic acid):acrylonitrile]

MP0335 Poly[(ethyl acrylate):(methacrylic acid)]  
MP0336 Poly[acrylonitrile:(butyl acrylate):(methacryl amide)]  
MP0337 Poly[acrylonitrile:(ethyl acrylate):(methacrylic acid):(methyl methacrylate)]  
MP0338 Poly(butadiene::styrene)  
MP0339 Poly(butadiene:styrene), Buna BL 6533  
MP0340 Poly(butadiene:styrene), Cariflex 1500  
MP0341 Poly(butadiene:styrene), Cariflex 1712  
MP0342 Poly(butadiene:styrene), Solprene 1205  
MP0343 Poly(butadiene:styrene), Solprene 1322  
MP0344 Poly(styrene::butadiene) ABA block (28% styrene), Kraton D 1101  
MP0345 Poly(styrene::butadiene) ABA block (30% styrene), Solprene 411  
MP0346 Poly(styrene::butadiene) block  
MP0347 Poly(styrene::butadiene) graft  
MP0348 Poly(styrene::butadiene::divinylbenzene)  
MP0349 Poly(styrene::butadiene::styrene)  
MP0350 Poly(styrene::butadiene::styrene::butadiene::styrene)  
MP0351 Poly(styrene::isoprene) ABA block (14% styrene)  
MP0352 Poly(styrene:butadiene)  
MP0353 Buna AP 251, Copolymer EPDM type  
MP0354 Buna AP 447, Copolymer EPDM type  
MP0355 Buna AP 541, Copolymer EPDM type  
MP0356 Dutral PMX 9705, Copolymer EPDM type  
MP0357 Dutral PMX 9710, Copolymer EPDM type  
MP0358 Dutral TER 4038 EP, Copolymer EPDM type  
MP0359 Dutral TER 9046, Copolymer EPDM type  
MP0360 Keltan 1446 A, Copolymer EPDM type  
MP0361 Keltan 314, Copolymer EPDM type  
MP0362 Keltan 480X100, Copolymer EPDM type  
MP0363 Keltan 509 X 100, Copolymer EPDM type  
MP0364 Keltan 512, Copolymer EPDM type  
MP0365 Poly(ethylene:propylene:diene) (50% ethylene, 4% diene)  
MP0366 Poly(ethylene:propylene:diene) (50% ethylene, 8% diene)  
MP0367 Poly(ethylene:propylene:diene) (70% ethylene, 4% diene)  
MP0368 Polysar EPDM 345, Copolymer EPDM type  
MP0369 Royalene 501, Copolymer EPDM type  
MP0370 Vistalon 5630, Copolymer EPDM type  
MP0371 Vistalon 6505, Copolymer EPDM type  
MP0372 Vistalon 6630, Copolymer EPDM type  
MP0373 Vistalon 8510, Copolymer EPDM type  
MP0374 Vistalon 9500, Copolymer EPDM type  
MP0375 Elvax 265, Copolymer EVA type  
MP0376 Elvax 360, Copolymer EVA type  
MP0377 Elvax 460, Copolymer EVA type  
MP0378 Escorene ultra 02020, Copolymer EVA type  
MP0379 Escorene ultra FL00728, Copolymer EVA type  
MP0380 Escorene ultra UL 00723, Copolymer EVA type  
MP0381 Evatane 1005 VN 3, Copolymer EVA type  
MP0382 Evatane 1005 VN 4, Copolymer EVA type  
MP0383 Evatane 1020 VB 2, Copolymer EVA type  
MP0384 Levapren 8338, Copolymer EVA type  
MP0385 Miravithen D23EA, Copolymer EVA type  
MP0386 Poly(ethylene:propylene), Daplen KSR 4525  
MP0387 Poly(ethylene:propylene), Keltan 740  
MP0388 Poly(propylene:ethylene), Novolen 2900NC  
MP0389 Polypropylene+poly(ethylene:propylene)  
MP0390 JSR SAN 230 NP, Copolymer SAN type  
MP0391 Kibisan 127, Copolymer SAN type

MP0392 Kibisan, Copolymer SAN type  
MP0393 Kostil AP, Copolymer SAN type  
MP0394 Kostil B555, Copolymer SAN type  
MP0395 Poly(acrylonitrile:styrene), PD Alfa 468  
MP0396 Poly(acrylonitrile:styrene), Polykemi  
MP0397 Poly(acrylonitrile:styrene), SAN TOYO 100 PCF  
MP0398 Poly(acrylonitrile:styrene), SAN TOYO 290  
MP0399 Poly(acrylonitrile:styrene), SAN TOYO 330  
MP0400 Poly(styrene:acrylonitrile) (70% styrene)  
MP0401 Kralex 1502, Copolymer SBR type  
MP0402 Kralex 1507, Copolymer SBR type  
MP0403 Poly(vinyl chloride:vinyl acetate:vinyl alcohol)  
MP0404 Poly[(vinyl acetate):(vinyl chloride)]  
MP0405 Reilex 425, Poly[(divinyl benzene):(vinyl pyridine)]  
MP0406 Poly(chloroprene)  
MP0407 Poly(hexafluoropropylene oxide) (mw: 7,000)  
MP0408 Poly(tetrafluoroethylene)  
MP0409 Poly(vinylidene fluoride) #1  
MP0410 Poly(vinylidene fluoride) #2  
MP0411 Polychlorotrifluorethylen, Fluorolube  
MP0412 Polyvinyl Fluoride  
MP0413 Polyvinylchloride #2  
MP0414 Polyvinylchloride, plasticized  
MP0415 Rubber, chlorinated  
MP0416 Elastin (protein structure)  
MP0417 Xanthan gum  
MP0418 Acetylated distarch adipate  
MP0419 Agar  
MP0420 Agar powder  
MP0421 Agar with maltose extract  
MP0422 Amylopectine  
MP0423 Amylopectine Ultra N  
MP0424 Asphalt  
MP0425 Caramel  
MP0426 Carboxymethyl cellulose, Na salt  
MP0427 Carboxymethyl cellulose, Na salt, ultra low viscosity  
MP0428 Carboxymethylamylopectine Sodium  
MP0429 Carboxymethylcellulose (Carmellose)  
MP0430 Casein  
MP0431 Casein hydrolysate  
MP0432 Casein sodium salt  
MP0433 Cellacefate, Cellulose acetate phthalate  
MP0434 Cellulose #1  
MP0435 Cellulose #2  
MP0436 Cellulose (MDOC), Ca salt  
MP0437 Cellulose Acetate (39.8% acetyl)  
MP0438 Cellulose Acetate Butyrate  
MP0439 Cellulose Acetate cine film  
MP0440 Cellulose Beaded (Perloza)  
MP0441 Cellulose Microcrystalline (Avicel)  
MP0442 Cellulose Modified  
MP0443 Cellulose Modified  
MP0444 Cellulose swabs  
MP0445 Cork  
MP0446 Deerskin  
MP0447 Dextran 40000  
MP0448 Dextran 60000-80000

MP0449 Explotab (Sodium Starch Glycolate)  
MP0450 Gelatin  
MP0451 Gelatin foam  
MP0452 Human Skin  
MP0453 Hydroxy propyl distarch phosphate  
MP0454 Hydroxybutyl Methyl Cellulose  
MP0455 Hydroxyethyl ether cellulose  
MP0456 Hydroxyethylcellulose  
MP0457 Hydroxymethylcellulose  
MP0458 Hypromellose (4000)  
MP0459 Kollidon VA 64  
MP0460 Leather  
MP0461 Lignine  
MP0462 Maltodextrin KMS X-50  
MP0463 Methocel E15  
MP0464 Methocel E3  
MP0465 Methocel E5  
MP0466 Methocel K100  
MP0467 Methocel K15  
MP0468 Methocel K4  
MP0469 Modified Starch  
MP0470 Nitrocellulose in mixture  
MP0471 Paper coated  
MP0472 Paper coated office  
MP0473 Paper coated, white card  
MP0474 Paper coated, white card matt  
MP0475 Paper filled with kaoline pigment  
MP0476 Paper glossy plotter  
MP0477 Paper napkin  
MP0478 Paper with polymer coating  
MP0479 Paper with polymer coating #1  
MP0480 Paper with polymer coating #2  
MP0481 Paper with polymer coating #3  
MP0482 Paper with polymer coating #4  
MP0483 Paper, beige recycled  
MP0484 Paper, brown box unbleached  
MP0485 Paper, brown unbleached for packing  
MP0486 Paper, colored, purple  
MP0487 Paper, extra thin yellow  
MP0488 Paper, matt  
MP0489 Paper, white  
MP0490 Paper, white box  
MP0491 Paper, white office  
MP0492 Papyrus Egypt  
MP0493 Pectine  
MP0494 Pectine Citrus  
MP0495 Pine Resin  
MP0496 Polyvidone 30, Polyvinylpolypyrrolidone  
MP0497 Scallop  
MP0498 Starch  
MP0499 Starch Modified  
MP0500 Starch sodium octenyl succinate  
MP0501 Starch Tapioca  
MP0502 Starch, Corn Starch Pregelatinised  
MP0503 Starch, Maize Starch Pregelatinised  
MP0504 Starch, Potato Starch Acetyl Oxidized  
MP0505 Starch, Rice Starch

MP0506 Macrogol 1500, Poly(ethylene glycol) 1500  
MP0507 Macrogol 200, Poly(ethylene glycol) 200  
MP0508 Macrogol 4000, Poly(ethylene glycol) 4000  
MP0509 Macrogol 6000, Poly(ethylene glycol) 6000  
MP0510 Poly(ethylene glycol) 1000  
MP0511 Poly(ethylene glycol) 14000  
MP0512 Poly(ethylene glycol) 2000  
MP0513 Poly(ethylene glycol) 300  
MP0514 Poly(ethylene glycol) 400  
MP0515 Poly(ethylene glycol) 4000  
MP0516 Poly(ethylene glycol) 6000  
MP0517 Poly(ethylene glycol) 8000  
MP0518 Poly(ethylene oxide)  
MP0519 Poly(propylene glycol) 2000  
MP0520 Polyethylene Oxide 4000000  
MP0521 Polyethylene Oxide 5000000  
MP0522 Polyethylene Oxide 900000  
MP0523 Poly(4,4-dipropoxy-2,2-diphenyl propane fumarate)  
MP0524 Poly(dianol 22 adipate)  
MP0525 Poly(dianol 22 endomethylene tetrahydrophthalate)  
MP0526 Poly(dianol 22 isophthalate)  
MP0527 Poly(dianol 22 maleate)  
MP0528 Poly(dianol 22 succinate)  
MP0529 Poly(dianol 22 tetrahydrophthalate)  
MP0530 Poly(dianol 33 adipate)  
MP0531 Poly(dianol 33 endomethylene tetrahydrophthalate)  
MP0532 Poly(dianol 33 fumarate)  
MP0533 Poly(dianol 33 isophthalate)  
MP0534 Poly(dianol 33 maleate)  
MP0535 Poly(dianol 33 phthalate)  
MP0536 Poly(dianol 33 succinate)  
MP0537 Poly(dianol 33 tetrahydrophthalate)  
MP0538 2,6-Di-tert-Butyl-p-Cresol  
MP0539 4,4-thiobis(t-butyl-m-cresol)  
MP0540 Antioxidant 2246  
MP0541 ARC I-1010  
MP0542 BHT  
MP0543 Chimisorb 81, 2-hydroxy-4-n-octyloxybenzophenone (UV-531)  
MP0544 Chimisorb 944  
MP0545 Cyasorb UV 3346  
MP0546 Didodecyl 3,3'-thiodipropionate, DLTDP  
MP0547 Dioctadecyl 3,3'-thiodipropionate  
MP0548 DSTDP  
MP0549 Hostanox 03  
MP0550 Hostavin N30 (HALS)  
MP0551 Irgafos 168  
MP0552 Irganox 1010  
MP0553 Irganox 1010 + Irganox 168  
MP0554 Irganox 1076  
MP0555 Irganox 1330  
MP0556 Irganox 168, Tris(2,4-di-terc.butylphenyl)phosphite  
MP0557 Irganox 3114  
MP0558 Leucopur EGM  
MP0559 Millad 3988  
MP0560 Santonox R  
MP0561 Tinuvin 622  
MP0562 Uvinul 4050 H

MP0563 Candelilla wax  
MP0564 Carnauba wax  
MP0565 Carnauba wax yellow  
MP0566 Ceresin  
MP0567 Cetaceum, Cetyl Palmitate 15  
MP0568 Chlorinated paraffin  
MP0569 Chlorowax  
MP0570 Ethoxylated Lanoline 4000  
MP0571 Glycerol Tristearate  
MP0572 Lanolin Anhydrous  
MP0573 Lanolin Wax  
MP0574 Paraffin  
MP0575 Paraffin wax  
MP0576 Paraffin, Flakes  
MP0577 Trimyristin  
MP0578 Tripalmitin  
MP0579 Wax  
MP0580 Wax, Emulsifying  
MP0581 Wax, White Bees Wax  
MP0582 Wax, Yellow Bees Wax  
MP0583 Sorbitan Oleate, Span 80  
MP0584 Sorbitan Stearate, Span 60  
MP0585 Tween 20, POE (20) sorbitan monolaurate  
MP0586 Tween 21, POE (4) sorbitan monolaurate  
MP0587 Tween 40, POE (20) sorbitan monostearate  
MP0588 Tween 60, POE (20) Sorbitan Monopalmitate  
MP0589 Tween 61, POE (4) sorbitan monopalmitate  
MP0590 Tween 80, POE (20) sorbitan monooleate  
MP0591 Tween 85, POE (20) sorbitan trioleate  
MP0592 1,3-Diphenylguanidine  
MP0593 1,4-Phenylenediamine  
MP0594 1-Naphthol  
MP0595 2-Mercaptobenzothiazole MBT  
MP0596 Aerosil 200  
MP0597 Aerosil R972  
MP0598 Aliquat 336  
MP0599 Alkanol XC  
MP0600 Alkylated Diphenylamine, Irganox L57  
MP0601 Aluminium acetate basic  
MP0602 Aluminium Distearate  
MP0603 Aluminium Hydroxide  
MP0604 Aluminium Oxide - pigment  
MP0605 Aluminium Oxide 90  
MP0606 Aluminium Stearate  
MP0607 Aluminosilicate  
MP0608 Aluminosilicate Sodium  
MP0609 Ammonium Bromide  
MP0610 Ammonium Chloride  
MP0611 Bicarbonate Ammonium  
MP0612 Bicarbonate Potassium  
MP0613 Bis (2-hydroxyethyl) tallow amine  
MP0614 Brij 35, Polyoxyethylene (35) Lauryl Ether  
MP0615 Bromochlorinated paraffin  
MP0616 Cabosil HS-5, SiO<sub>2</sub> based  
MP0617 Cabosil LM-150, SiO<sub>2</sub> based  
MP0618 Calcium Carbonate  
MP0619 Calcium Chloride Crystalline

MP0620 Calcium Chloride Powder  
MP0621 Calcium Oxide  
MP0622 Calcium Stearate  
MP0623 Aluminiumsilicate Calcium Sodium  
MP0624 Carbonate Ammonium  
MP0625 Carbonate Magnesium  
MP0626 Carbonate Potassium  
MP0627 Carbonate Sodium anh.  
MP0628 Castor oil, ethoxylated  
MP0629 Chloride Magnesium  
MP0630 Coal Tar  
MP0631 Colloid Defoamer  
MP0632 Dapro DF-911, Defoamer  
MP0633 Decabromdiphenyloxide  
MP0634 Di-2-ethylhexyl phthalate  
MP0635 Dibenzoyl peroxide  
MP0636 Dibutyl phthalate  
MP0637 Diethylenetriamine  
MP0638 Diisodecyl phthalate  
MP0639 Di-isononyl phthalate  
MP0640 Disodium Calcium Edetate  
MP0641 Disodium EDTA  
MP0642 Ditridecyl Phthalate  
MP0643 Docosanoic Acid, Behenic Acid  
MP0644 Durafur developer D  
MP0645 Eastobrite OB-1  
MP0646 Edetate Dipotassium Magnesium  
MP0647 Emulsifier E 2149, Stearyl Alcohol and Esters Mixture  
MP0648 Erucamide  
MP0649 Ethoxylated tallow amine  
MP0650 Ethyl palmitate  
MP0651 Ethylenediaminetetraacetic acid = EDTA  
MP0652 Fatty Amines Modified  
MP0653 Ferric Hydroxide  
MP0654 Glass fiber  
MP0655 Glass wool  
MP0656 Glycerides of fatty acids C12-C22  
MP0657 Glycerides of natural fatty acids  
MP0658 Glycerol Stearate  
MP0659 Glyceryl Stearate  
MP0660 Gum Tragacanth  
MP0661 Hexadecanoic Acid, Palmitic acid  
MP0662 Hexamethylenetetramine  
MP0663 Hydroxide Aluminium  
MP0664 Iron (III) Hydroxide  
MP0665 Iron (III) Oxide Hydroxide  
MP0666 Iron (III) Oxide monohydrate, yellow  
MP0667 Iron (III) Oxide, Hematite  
MP0668 Iron (III) Oxide, Red  
MP0669 Iron Oxide  
MP0670 Iron Oxide Black  
MP0671 Iron(II) oxalate dihydrate  
MP0672 Kaolin, Aluminum silicate dihydrate  
MP0673 Loxamid E, Stearoamides mixture  
MP0674 Loxamid OHA, Oleo and stearoamides  
MP0675 Magnesium acetate tetrahydrate  
MP0676 Magnesium Carbonate

MP0677 Magnesium Chloride  
MP0678 Magnesium Chloride Hexahydrate  
MP0679 Magnesium Distearate  
MP0680 Magnesium hydroxide  
MP0681 Magnesium Oxide  
MP0682 Magnesium Stearate  
MP0683 Manganese stearate  
MP0684 Methyl Linoleate  
MP0685 Methyl stearate  
MP0686 Molybdenum Sulfide Lubricant  
MP0687 Mono and Diglycerides of fatty acids  
MP0688 Nacrelon S-30 (pearlescent pigment)  
MP0689 Nacure 1040  
MP0690 Nonadecanoic Acid  
MP0691 Octadecanoic Acid, Stearic Acid  
MP0692 Oil, Caraway oil  
MP0693 Opadry OY-B-22932 yellow - coating mixture  
MP0694 Opadry OY-S-8720 orange - coating mixture  
MP0695 Opaspray (mixture)  
MP0696 Opex, amine based polymer additive  
MP0697 Oxide Aluminium  
MP0698 Pareth 15 (ethoxylated alcohols)  
MP0699 Plasticizer, Tegosoft EE  
MP0700 Polynol A, mono and diglycerides  
MP0701 Polysorbate 20  
MP0702 Quartz  
MP0703 Silica Gel, 70-230 mesh  
MP0704 Silicagel  
MP0705 Silicate - pigment  
MP0706 Silicate with cellulose - pigment  
MP0707 Silicon dioxide  
MP0708 Silicon dioxide, highly dispersed  
MP0709 Silipor 600  
MP0710 Silwet L-77  
MP0711 Soda Ash  
MP0712 Sodium Alkyl Aryl Sulfonate  
MP0713 Sodium Dodecylbenzene Sulfonate  
MP0714 Sodium Metasilicate  
MP0715 Sodium Silicate  
MP0716 Sorbitan Monolaurate, Span 20  
MP0717 Sorbitan tristearate  
MP0718 Span 40, Sorbitan Monopalmitate  
MP0719 Stearamide  
MP0720 Stearate Calcium  
MP0721 Stearate Zinc  
MP0722 Succinimide  
MP0723 Sulfate Calcium  
MP0724 Surfactant, Sodium Alkyl polyglycol ether sulfate  
MP0725 Surfactant, Syntopon A 100  
MP0726 Talc  
MP0727 Tetrasodium Edetate  
MP0728 Titanium Dioxide  
MP0729 Titanium Ethoxide  
MP0730 Tri(n-octyl,n-decyl)trimellitate  
MP0731 Trilaurin  
MP0732 Triton, X-100  
MP0733 Tween 65, POE (20) sorbitan tristearate

MP0734 Tween 81, POE (5) sorbitan monooleate  
 MP0735 Vanfre IL-1  
 MP0736 Zinc Oxide  
 MP0737 Zinc Stearate  
 MP0738 2,6-di-t-butyl-4-methylphenol, BHT  
 MP0739 Alkyd with Castor Oil  
 MP0740 Alkyd, 40% soybean oil  
 MP0741 Alkyd, 45% linseed oil  
 MP0742 Alkyd, 65% soybean oil  
 MP0743 Phthalic Alkyd  
 MP0744 Color masterbatch Poly[ethylene-co-(vinyl acetate)] + 45% benzimidazolone  
 MP0745 Color masterbatch Poly[ethylene-co-(vinyl acetate)] + 45% yellow pigment #3  
 MP0746 Color masterbatch Poly[ethylene-co-(vinyl acetate)] + 45% yellow pigment #4  
 MP0747 Color masterbatch Poly[ethylene-co-(vinyl acetate)] + 45% yellow pigment #1  
 MP0748 Color masterbatch Poly[ethylene-co-(vinyl acetate)] + 45% yellow pigment #2  
 MP0749 Color masterbatch Polyethylene + 30% red pigment #1  
 MP0750 Color masterbatch Polyethylene + 30% red pigment #2  
 MP0751 Color masterbatch Polyethylene + 30% red pigment #3  
 MP0752 Color masterbatch Polyethylene + 30% red pigment #4  
 MP0753 Color masterbatch Polyethylene + 30% yellow pigment #1  
 MP0754 Color masterbatch Polyethylene + 30% yellow pigment #2  
 MP0755 Color masterbatch Polyethylene + 30% yellow pigment #3  
 MP0756 Color masterbatch Polyethylene + 30% yellow pigment #4  
 MP0757 Color masterbatch Polyethylene + 30% yellow pigment #5  
 MP0758 Color masterbatch Polyethylene + 40% yellow pigment  
 MP0759 Color masterbatch Polyethylene + 45% yellow pigment  
 MP0760 Color masterbatch Polyethylene + 60% brown pigment  
 MP0761 Color masterbatch Polyethylene + 60% Fe-oxide  
 MP0762 Color masterbatch Polyethylene + 60% yellow pigment #1  
 MP0763 Color masterbatch Polyethylene + 60% yellow pigment #2  
 MP0764 Color masterbatch Polyethylene + 60% yellow pigment #3  
 MP0765 Color masterbatch Polyethylene + 75% white pigment  
 MP0766 Color masterbatch Polyethylene + Zn, Fe oxide  
 MP0767 Color masterbatch Polypropylene + 70% white pigment  
 MP0768 Color masterbatch Polypropylene + brown pigment  
 MP0769 Color masterbatch Polypropylene + white pigment  
 MP0770 Color masterbatch Polypropylene + yellow pigment #1  
 MP0771 Color masterbatch Polystyrene + 30%, Sb, Ni, Ti-oxide  
 MP0772 Color masterbatch Polystyrene + 50% brown pigment  
 MP0773 Color masterbatch Polystyrene + 50% white pigment  
 MP0774 Polypropylene + yellow pigment #2  
 MP0775 Color masterbatch Poly(styrene-co-acrylonitrile) + 50% white pigment  
 MP0776 (+/-)-Camphor  
 MP0777 (E)-Anethole  
 MP0778 1,12-Diaminododecane  
 MP0779 1,2,4-Triacetoxymethylbenzene  
 MP0780 1,2-Bis(diphenylphosphino)ethane  
 MP0781 1,2-Dicyanobenzene  
 MP0782 1,3,5-Triethylhexahydro-1,3,5-triazine  
 MP0783 1,4-Phenylenediacetic Acid  
 MP0784 1,5-Dimethylnaphthalene  
 MP0785 12-Hydroxydodecanoic acid  
 MP0786 1-Adamantanecarboxylic acid  
 MP0787 1-Bromohexadecane  
 MP0788 1-Bromotetradecane  
 MP0789 1-Chloro-2,4-dinitrobenzene  
 MP0790 1-Docosanol

MP0791 1-Heptadecanol  
MP0792 1-Heptene  
MP0793 1-Hexadecanol  
MP0794 1-Iodobutane  
MP0795 1-Methyl-4-piperidone  
MP0796 1-Methylnaphthalene  
MP0797 1-Nonadecanol  
MP0798 1-Octacosane  
MP0799 1-Octadecanol  
MP0800 1-Phenyl-1,2-ethanediol  
MP0801 1-Phenyl-2-thiourea  
MP0802 2-(Dibutylamino)ethanol  
MP0803 2,2,4-Trimethylpentane  
MP0804 2,2'-Azobis(isobutyronitrile)  
MP0805 2,4,6-Trimethylpyridine  
MP0806 2,4-Dihydroxybenzaldehyde  
MP0807 2,4-Hexadienoic acid, potassium salt  
MP0808 2'-Acetonaphthone  
MP0809 2-Amino-2-(hydroxymethyl)-1,3-propanediol; TRIS  
MP0810 2-Azacyclotridecanone  
MP0811 2-Butenoic acid, ethyl ester  
MP0812 2-Chlorobenzoic acid  
MP0813 2-Chlorobenzonitrile  
MP0814 2-Ethoxynaphthalene  
MP0815 2-Ethylbutyric acid  
MP0816 2-Ethylcaproic Acid  
MP0817 2-Ethylhexyl nitrate  
MP0818 2-Fluoropropionic acid  
MP0819 2-Hydroxyisobutyric acid  
MP0820 2-Methoxybenzaldehyde  
MP0821 2-Methoxynaphthalene  
MP0822 2-Methylnaphthalene  
MP0823 2-Naphthoic acid  
MP0824 2-Naphthol  
MP0825 2-Naphthyl acetate  
MP0826 2-Phenoxyethanol  
MP0827 2-Thiophenecarboxaldehyde  
MP0828 3,5-Dimethylphenol = Xylenol  
MP0829 3-Aminophenol  
MP0830 3-Chlorobenzaldehyde  
MP0831 3-Hexyne  
MP0832 3-Methyl-1-butanol  
MP0833 3-Nitrophenol  
MP0834 3-Picoline = 3-Methylpyridine  
MP0835 4,4'-Methylenebis(N,N-dimethylaniline)  
MP0836 4-Aminobenzoic acid  
MP0837 4-Aminosalicylic acid sodium salt  
MP0838 4-Bromoaniline  
MP0839 4-Bromophenol  
MP0840 4-Chlorobenzoic acid  
MP0841 4-Chlorobutyronitrile  
MP0842 4-Fluorobenzoic acid  
MP0843 4-Iodophenol  
MP0844 4-Methoxybenzaldehyde  
MP0845 4-Methoxybenzonitrile  
MP0846 4-Methyl-2-pentanone  
MP0847 4-Vinylcyclohexene dioxide

MP0848 5,5-Dimethyl-1,3-cyclohexanedione  
MP0849 6-Aminohexanoic acid  
MP0850 9,10-Dicyanoanthracene  
MP0851 Acenaphthene  
MP0852 Acetamide  
MP0853 Acetanilide  
MP0854 Acetoacetanilide  
MP0855 Acetylsalicylic acid  
MP0856 Allyl sulfide  
MP0857 Aluminium Chlorohydrate  
MP0858 Ammonium Acetate  
MP0859 Ammonium Citrate Monobasic  
MP0860 Ammonium iodide  
MP0861 Anthracene  
MP0862 Anthraquinone  
MP0863 Azelaic acid  
MP0864 Bentonite  
MP0865 Benzamide  
MP0866 Benzenesulfonic acid  
MP0867 Benzil  
MP0868 Benzimidazole  
MP0869 Benzoin  
MP0870 Benzophenone  
MP0871 Benzotriazole  
MP0872 Benzyl bromide  
MP0873 Benzyltrimethylammonium chloride  
MP0874 Beta-cyclodextrine  
MP0875 Bicarbonate Sodium  
MP0876 Biphenyl  
MP0877 Biuret  
MP0878 Boric Acid  
MP0879 Bromoacetic acid  
MP0880 Bromobutyrolactone  
MP0881 Butylated hydroxyanisole, BHA  
MP0882 Butyltriphenylphosphonium bromide  
MP0883 Calcium Acetate  
MP0884 Calcium chloride dihydrate  
MP0885 Calcium Diacetate  
MP0886 Calcium Hydroxide  
MP0887 Calcium Lactate  
MP0888 Calcium Petroleum Sulfonate  
MP0889 Camphene  
MP0890 Carbonate Sodium Decahydrate  
MP0891 Carbonate Zinc Basic  
MP0892 Cetyldimethylamine  
MP0893 Chloral hydrate  
MP0894 Chloroacetamide  
MP0895 Chloroacetic acid  
MP0896 Chlorodiphenylphosphine  
MP0897 Cinnamic Acid  
MP0898 Citric Acid  
MP0899 Citric Acid Monohydrate  
MP0900 Cocoamidopropyl Hydroxysultaine  
MP0901 Cupric acetate .H2O  
MP0902 Cyclododecanone  
MP0903 D-(+)-Malic acid  
MP0904 Dextrin

MP0905 D-Galactonic acid g-lactone  
MP0906 Di(ethylene glycol) distearate  
MP0907 Dicyandiamide sulfate  
MP0908 Diethyl dithiophosphate, ammonium salt  
MP0909 Diethyl glutarate  
MP0910 Diethyl isoamylmalonate  
MP0911 Diethyl Sulfate  
MP0912 Diethyl t-butylmalonate  
MP0913 Dimethyl Sulfate  
MP0914 Diphenyl  
MP0915 Diphenylacetylene  
MP0916 Diphosphate Tetrasodium  
MP0917 Diphosphate Tetrasodium Decahydrate  
MP0918 Dipropylene Glycol  
MP0919 Disodium Orthophosphate  
MP0920 Di-t-butyl dicarbonate  
MP0921 Dithionite Sodium  
MP0922 DL-Camphor  
MP0923 DL-Menthol  
MP0924 DL-Tartaric Acid  
MP0925 Dodecanoic acid  
MP0926 Dodecyl aldehyde  
MP0927 D-Tartaric Acid  
MP0928 Ethyl caprylate  
MP0929 Ethyl maltol  
MP0930 Ethyl nicotinate  
MP0931 Ethylene glycol diformate  
MP0932 Gallic acid  
MP0933 Gallotannin  
MP0934 Glycerol triacetate  
MP0935 Glyceryl monooleate  
MP0936 Glyceryl monostearate  
MP0937 Heneicosane  
MP0938 Hexachlorobenzene  
MP0939 Hexamethylbenzene  
MP0940 Hexatriacontane  
MP0941 Hydroquinone  
MP0942 Hypophosphite Sodium  
MP0943 Imidazole  
MP0944 Indole  
MP0945 Inositol  
MP0946 Iodoform  
MP0947 Iron Citrate  
MP0948 Isopropyl Glycidyl Ether  
MP0949 Isostearyl isostearate  
MP0950 Itaconic acid  
MP0951 Lead(II) nitrate  
MP0952 L-Lactic acid  
MP0953 L-Lactic acid, Sodium salt  
MP0954 L-Malic acid  
MP0955 L-Menthol  
MP0956 Ludipress (mixture)  
MP0957 Macrogol Sodium Sulphosuccinate 30%  
MP0958 Magnesium Carbonate Basic Hydrate  
MP0959 Maleic acid  
MP0960 Menthol  
MP0961 Metabisulfite Potassium

MP0962 Metabisulfite Sodium  
MP0963 Methanesulfonic acid  
MP0964 Methyl 2-furoate  
MP0965 Methylenedilaurate  
MP0966 Methylhydroquinone  
MP0967 Methyltin Mercaptoester Sulfides mixed  
MP0968 Methyltriphenylphosphonium bromide  
MP0969 Michler's Ketone  
MP0970 Morpholine hydrobromide  
MP0971 m-Phthalodinitrile  
MP0972 Myo-Inositol  
MP0973 Myristic acid  
MP0974 N,N-Dimethyl-1-naphthylamine  
MP0975 N,N-Dimethyl-p-toluidine  
MP0976 N,N'-Methylenebisacrylamide  
MP0977 Naphthalene  
MP0978 N-Benzyl-N,N,N-triethylammonium chloride  
MP0979 N-Bromosuccinimide  
MP0980 n-Dotriacontane  
MP0981 N-Hexadecyltrimethylammonium chloride  
MP0982 Nicotinic acid = Niacin  
MP0983 Nonacosane  
MP0984 Nonanedioic acid  
MP0985 Nonanoic acid  
MP0986 n-Pentyl Formate  
MP0987 Octadecanoic Acid, Calcium Salt  
MP0988 Octylbenzene  
MP0989 Oil of Peppermint, Menthol  
MP0990 Orthoboric Acid  
MP0991 Ostion (HPLC sorbent)  
MP0992 o-Tolidine  
MP0993 o-Tolunitrile  
MP0994 Oxalic Acid  
MP0995 Pentaerythritol  
MP0996 Pentasodium triphosphate  
MP0997 Pentyl propionate  
MP0998 Persulfate Ammonium  
MP0999 Persulfate Potassium  
MP1000 Persulfate Sodium  
MP1001 Petroleum Sulfonate  
MP1002 Phenyl Salicylate  
MP1003 Phenylboronic acid  
MP1004 Phenylmercuric acetate  
MP1005 Phenylurea  
MP1006 Phosphate Ammonium Dibasic  
MP1007 Phosphate Ammonium Monobasic  
MP1008 Phosphate Calcium Dibasic  
MP1009 Phosphate Calcium Monobasic  
MP1010 Phosphate Calcium Tribasic  
MP1011 Phosphate Calcium, Hydroxylapatite  
MP1012 Phosphate Potassium Dibasic  
MP1013 Phosphate Potassium Dibasic Trihydrate  
MP1014 Phosphate Potassium Monobasic  
MP1015 Phosphate Potassium Tribasic n-hydrate  
MP1016 Phosphate Sodium Dibasic  
MP1017 Phosphate Sodium Dibasic Dihydrate  
MP1018 Phosphate Sodium Dibasic Heptahydrate

MP1019 Phosphate Sodium Monobasic  
MP1020 Phosphate Sodium Monobasic Dihydrate  
MP1021 Phosphate Sodium Tribasic  
MP1022 Phthalide  
MP1023 Piperazine dihydrochloride  
MP1024 Pivalic acid t-butylamide  
MP1025 Poly(sodium 4-styrenesulfonate)  
MP1026 Polyglycerol Polyricinoleate  
MP1027 Polyglyceryl Oleate  
MP1028 Potassium Acetate  
MP1029 Potassium bromate  
MP1030 Potassium carbonate, sesquihydrate  
MP1031 Potassium chloride  
MP1032 Potassium Ferricyanide  
MP1033 Potassium Ferrocyanide  
MP1034 Potassium Nitrate  
MP1035 Potassium permanganate  
MP1036 Propionic Acid  
MP1037 Propyl Gallate  
MP1038 p-Tolunitrile  
MP1039 Pyridine hydrochloride  
MP1040 Quinhydrone  
MP1041 Resorcinol  
MP1042 Salicylic acid  
MP1043 Silver Nitrate  
MP1044 Sodium Acetate  
MP1045 Sodium Acetate Trihydrate  
MP1046 Sodium Cetylsulfate  
MP1047 Sodium Dodecyl Sulfate  
MP1048 Sodium Nitrate  
MP1049 Sodium Octanoate  
MP1050 Sodium Perborate\* 4H<sub>2</sub>O  
MP1051 Sodium Stearyl Fumarate  
MP1052 Sodium Tetraborate  
MP1053 Sodium Tetraborate Decahydrate  
MP1054 Sodium Toluenesulfonate  
MP1055 Sodium trichloroacetate  
MP1056 Spheron Oxin 1000 LC  
MP1057 Stannous Chloride  
MP1058 Strontium oxalate  
MP1059 Sucrose (crystals)  
MP1060 Sulfamic acid = Amidosulfonic acid  
MP1061 Sulfate Aluminium  
MP1062 Sulfate Aluminium Potassium  
MP1063 Sulfate Ammonium  
MP1064 Sulfate Copper  
MP1065 Sulfate Ferrous  
MP1066 Sulfate Magnesium  
MP1067 Sulfate Sodium Monobasic  
MP1068 Tannic acid  
MP1069 Tetrabutylammonium bromide  
MP1070 Tetrabutylammonium chloride  
MP1071 Tetrabutylammonium cyanide  
MP1072 Tetrabutylammonium fluoride trihydrate  
MP1073 Tetrabutylammonium hexafluorophosphate  
MP1074 Tetrabutylammonium hydrogen sulfate  
MP1075 Tetrabutylammonium iodide

MP1076 Tetrabutylammonium perchlorate  
MP1077 Tetrabutylammonium sulfamate  
MP1078 Tetracosane  
MP1079 Tetraethylammonium bromide  
MP1080 Tetraethylammonium chloride  
MP1081 Tetraethylammonium tetrafluoroborate  
MP1082 Tetramethylammonium bromide  
MP1083 Tetramethylammonium fluoride  
MP1084 Tetramethylammonium hydroxide  
MP1085 Tetramethylammonium iodide  
MP1086 Tetraoctylammonium bromide  
MP1087 Tetraphenylboron Sodium  
MP1088 Thioacetic acid  
MP1089 Thiosulfate Sodium  
MP1090 Thiourea  
MP1091 trans-Cinnamic acid  
MP1092 Triaryl Phosphate Ester  
MP1093 Tribenzylamine  
MP1094 Tributyltin bromide  
MP1095 Tributyltin hydride  
MP1096 Tricalcium Orthophosphate  
MP1097 Triethanolamine hydrochloride  
MP1098 Triethyl phosphonoacetate  
MP1099 Triethylene Glycol  
MP1100 Trimethylamine hydrochloride  
MP1101 Trimethylolpropane Trinonanoate  
MP1102 Tripelargonin  
MP1103 Triphenyl phosphate  
MP1104 Triphenyl Phosphine  
MP1105 Triphenylantimony  
MP1106 Tritylalcohol  
MP1107 Vestamid E40  
MP1108 Watermelon Ketone  
MP1109 Zinc acetate dihydrate  
MP1110 Zinc salicylate  
MP1111 (+/-)-2-Butanol  
MP1112 (2R,3R)-(-)-2,3-Butanediol  
MP1113 (2S,3S)-(+)-2,3-Butanediol  
MP1114 (R)-(-)-2-Butanol  
MP1115 (R)-(-)-2-Octanol  
MP1116 (R)-(-)-Ethyl mandelate  
MP1117 (S)-(-)-2-Methyl-1-butanol  
MP1118 (S)-(+)-1,2-Propanediol  
MP1119 (S)-(+)-1,3-Butanediol  
MP1120 (S)-(+)-2-Butanol  
MP1121 (S)-(+)-2-Octanol  
MP1122 (S)-(+)-Ethyl mandelate  
MP1123 1,1,1,3,3,3-Hexafluoro-2-propanol  
MP1124 1,1,2-Trichlorotrifluoroethane  
MP1125 1,1,2-Trimethoxyethane  
MP1126 1,10-Dibromodecane  
MP1127 1,2,3,4-Tetrahydronaphthalene  
MP1128 1,2-Dibromoethane  
MP1129 1,2-Dibromopropane  
MP1130 1,2-Dichlorobenzene  
MP1131 1,2-Dichloroethane  
MP1132 1,3-Butanediol

MP1133 1,3-Butanediol diacetate  
MP1134 1,3-Dibromopropane  
MP1135 1,3-Propanediol  
MP1136 1,4-Butanediol  
MP1137 1,4-Butanediol diglycidyl ether  
MP1138 1,4-Dibromobutane  
MP1139 1,4-Dioxane  
MP1140 1,4-Pentanediol  
MP1141 1,5-Cyclooctadiene  
MP1142 1,5-Pentanediol  
MP1143 1,6-Dibromohexane  
MP1144 1-Bromo-3-chloropropane  
MP1145 1-Bromobutane  
MP1146 1-Bromododecane  
MP1147 1-Bromoheptane  
MP1148 1-Bromohexane  
MP1149 1-Bromonaphthalene  
MP1150 1-Bromononane  
MP1151 1-Bromooctane  
MP1152 1-Butanol  
MP1153 1-Chlorohexane  
MP1154 1-Chloronaphthalene  
MP1155 1-Decanol  
MP1156 1-Decene  
MP1157 1-Dodecene  
MP1158 1-Fluoro-2-nitrobenzene  
MP1159 1-Heptanol  
MP1160 1-Hexadecene  
MP1161 1-Hexanol  
MP1162 1-Iodohexane  
MP1163 1-Nonanol  
MP1164 1-Octadecene  
MP1165 1-Octanol  
MP1166 1-Pentanol = Amyl alcohol  
MP1167 1-Phenyl-1-propanol  
MP1168 1-Propanol  
MP1169 1-Tetradecene  
MP1170 2-(2-Butoxyethoxy)ethanol  
MP1171 2-(Ethylamino)ethanol  
MP1172 2,2,2-Trifluoroethanol  
MP1173 2,2'-Dichloroethyl ether  
MP1174 2,2'-Thiodiethanol  
MP1175 2,3-Butanediol  
MP1176 2,3-Dimethyl-2-butanol  
MP1177 2,6-Lutidine  
MP1178 2-Amino-2-methyl-1-propanol  
MP1179 2-Benzyloxyethanol  
MP1180 2-Bromohexanoic acid  
MP1181 2-Bromopentane  
MP1182 2-Butanol  
MP1183 2-Butanone  
MP1184 2-Butene-1,4-diol (cis and trans)  
MP1185 2-Chlorobenzaldehyde  
MP1186 2-Chloroethanol  
MP1187 2-Ethoxyethyl acetate  
MP1188 2-Ethyl-1-hexanol  
MP1189 2-Ethylhexylamine

MP1190 2-Ethylpiperidine  
MP1191 2-Hexanol  
MP1192 2-Mercaptoethanol  
MP1193 2-Methyl-1-butanol  
MP1194 2-Methyl-1-propanol  
MP1195 2-Methyl-2-butanol  
MP1196 2-Methyl-2-hexanol  
MP1197 2-Methyl-2-pentanol  
MP1198 2-Methylcyclohexanone  
MP1199 2-Nitropropane  
MP1200 2-Octanol  
MP1201 2-Picoline  
MP1202 3-Acetylpyridine  
MP1203 3-Amino-1-propanol  
MP1204 3-Chloro-1,2-propanediol  
MP1205 3-Hydroxypropionitrile  
MP1206 4-Chlorotoluene  
MP1207 4-Isopropyltoluene, p-Cymene  
MP1208 Acetic Acid  
MP1209 Acetic anhydride  
MP1210 Acetone  
MP1211 Acetonitrile  
MP1212 Acetophenone  
MP1213 Acetylacetone  
MP1214 Acrylonitrile  
MP1215 Amyl Acetate  
MP1216 Amylamine  
MP1217 Anethole  
MP1218 Aniline  
MP1219 *n*-Bromostyrene  
MP1220 Benzaldehyde  
MP1221 Benzene  
MP1222 Benzonitrile  
MP1223 Benzyl acetate  
MP1224 Benzyl alcohol  
MP1225 Benzyl benzoate  
MP1226 Benzyl chloride  
MP1227 Benzyl cyanide  
MP1228 Benzylamine  
MP1229 Butyl acetate  
MP1230 Butyl ether  
MP1231 Butyl formate  
MP1232 Butyl laurate  
MP1233 Butyl methacrylate  
MP1234 Butyl vinyl ether  
MP1235 Butyraldehyde  
MP1236 Carbon disulfide  
MP1237 Carbon Tetrachloride  
MP1238 Chloroform  
MP1239 *cis*-2-Butene-1,4-diol  
MP1240 Crotonic acid  
MP1241 Cumene  
MP1242 Cyclododecene  
MP1243 Cycloheptanol  
MP1244 Cyclohexane  
MP1245 Cyclohexanol  
MP1246 Cyclohexanone

MP1247 Cyclohexene  
MP1248 Cyclohexyl acetate  
MP1249 Cyclohexyl isocyanate  
MP1250 Cyclohexylamine  
MP1251 Cyclohexylmercaptan  
MP1252 Cyclooctanone  
MP1253 Cyclooctene  
MP1254 Cyclopentane  
MP1255 Cyclopentanol  
MP1256 Cyclopentanone  
MP1257 D-(-)-Lactic acid  
MP1258 Decane  
MP1259 Diacetone Alcohol  
MP1260 Diallyl Phthalate  
MP1261 Dibutyl ketone  
MP1262 Dibutyl sulfide  
MP1263 Dibutylamine  
MP1264 Dichloromethane  
MP1265 Dicyclohexylamine  
MP1266 Dicyclopentadiene  
MP1267 Diethanolamine  
MP1268 Diethyl acetylmalonate  
MP1269 Diethyl adipate  
MP1270 Diethyl carbonate  
MP1271 Diethyl ether  
MP1272 Diethyl ethylmalonate  
MP1273 Diethyl malonate  
MP1274 Diethyl methylmalonate  
MP1275 Diethyl oxalate  
MP1276 Diethyl phthalate  
MP1277 Diethyl pimelate  
MP1278 Diethyl propylmalonate  
MP1279 Diethyl suberate  
MP1280 Diethyl succinate  
MP1281 Diethylamine  
MP1282 Diethylene glycol  
MP1283 Diethylene Glycol Dibutyl Ether  
MP1284 Diethylene glycol dimethyl ether  
MP1285 Dihydropyran  
MP1286 Diisopropyl phosphite  
MP1287 Diisobutyl carbonate  
MP1288 Diisobutylketone  
MP1289 Diisopropyl ether  
MP1290 Diisopropylamine  
MP1291 Dimethyl phthalate  
MP1292 Dimethyl succinate  
MP1293 Dimethylaminoethanol  
MP1294 Dimethylsulfoxide  
MP1295 Dioctyl phthalate  
MP1296 Diundecyl phthalate  
MP1297 DL-Lactic acid  
MP1298 Dodecane  
MP1299 Dodecenylsuccinic anhydride  
MP1300 Dodecyl acetate  
MP1301 e-Caprolactam  
MP1302 Epichlorhydrin  
MP1303 Ethanol

MP1304 Ethoxyacetic acid  
MP1305 Ethyl acetate  
MP1306 Ethyl acetoacetate  
MP1307 Ethyl Acrylate  
MP1308 Ethyl benzoate  
MP1309 Ethyl bromoacetate  
MP1310 Ethyl butyrate  
MP1311 Ethyl caproate  
MP1312 Ethyl chloroacetate  
MP1313 Ethyl cinnamate  
MP1314 Ethyl crotonate  
MP1315 Ethyl D-lactate  
MP1316 Ethyl DL-mandelate  
MP1317 Ethyl dodecanoate  
MP1318 Ethyl formate  
MP1319 Ethyl Lactate  
MP1320 Ethyl salicylate  
MP1321 Ethyl-2-methylacetoacetate  
MP1322 Ethylbenzene  
MP1323 Ethylene carbonate  
MP1324 Ethylene Diglycol Monoethyl Ether  
MP1325 Ethylene glycol  
MP1326 Ethylene glycol dimethacrylate  
MP1327 Ethylene glycol dimethyl ether  
MP1328 Ethylene glycol monobutyl ether  
MP1329 Ethylene glycol monoethyl ether  
MP1330 Ethylene glycol monomethyl ether  
MP1331 Formamide  
MP1332 Formic acid  
MP1333 Furfuryl alcohol  
MP1334  $\gamma$ -Butyrolactone  
MP1335  $\gamma$ -Caprolactone  
MP1336 Glycerol  
MP1337 Heptadecane  
MP1338 Heptane  
MP1339 Hexachloro-1,3-butadiene  
MP1340 Hexamethylphosphoramide  
MP1341 Hexane  
MP1342 Hexanoic acid  
MP1343 Hexanol  
MP1344 Indene  
MP1345 Isoamyl acetate  
MP1346 Isoamyl alcohol  
MP1347 Isoamyl ether  
MP1348 Isoamyl formate  
MP1349 Isobutyl alcohol  
MP1350 Isobutyronitrile  
MP1351 Isopropenyl acetate  
MP1352 Isopropylamine  
MP1353 Isovaleric acid  
MP1354 L-(+)-Lactic acid, ethyl ester  
MP1355 m-Cresol  
MP1356 Mesityl oxide  
MP1357 Methyl alcohol  
MP1358 Methyl benzoate  
MP1359 Methyl chloroacetate  
MP1360 Methyl chloroformate

MP1361 Methyl Crotonate  
MP1362 Methyl dichloroacetate  
MP1363 Methyl formate  
MP1364 Methyl isobutyl ketone  
MP1365 Methyl Oleate  
MP1366 Methyl Palmitate  
MP1367 Methyl p-toluenesulfonate  
MP1368 Methyl Salicylate  
MP1369 Methylcyclohexanol  
MP1370 Methylcyclohexanols  
MP1371 Monoethanolamine  
MP1372 Morpholine  
MP1373 m-Toluidine  
MP1374 N,N,N',N'-Tetramethylethylenediamine  
MP1375 N,N-Diethylethanolamine  
MP1376 N,N-Dimethylacetamide  
MP1377 N,N-Dimethylformamide  
MP1378 Nitrobenzene  
MP1379 Nitrocyclohexane  
MP1380 Nitroethane  
MP1381 Nitromethane  
MP1382 N-Methylacetamide  
MP1383 N-Methylcyclohexylamine  
MP1384 N-Methylformamide  
MP1385 Nonylamine  
MP1386 Nonylphenol  
MP1387 o-Anisidine  
MP1388 o-Bromaniline  
MP1389 o-Cresol  
MP1390 Octanoic acid = Caprylic acid  
MP1391 Oleic acid  
MP1392 o-Nitrotoluene  
MP1393 o-Toluidine  
MP1394 o-Xylene  
MP1395 p-Cresol  
MP1396 p-Dioxane  
MP1397 Pentane  
MP1398 Petroleum ether  
MP1399 Petroleum naphtha  
MP1400 Phenyl Ether  
MP1401 Phenyl Isocyanate  
MP1402 Piperidine  
MP1403 Propargyl alcohol  
MP1404 Propionaldehyde  
MP1405 Propyl acetate  
MP1406 Propylene carbonate  
MP1407 Propylene Glycol  
MP1408 Propylene Oxide  
MP1409 p-Toluidine  
MP1410 Pyridine  
MP1411 R-(-)-1,3-Butanediol  
MP1412 Styrene  
MP1413 t-Butyl Acetate  
MP1414 t-Butyl alcohol  
MP1415 t-Butylamine  
MP1416 t-Butylbenzene  
MP1417 Tetrachloroethylene

MP1418 Tetrahydrofuran  
MP1419 Toluene  
MP1420 Triacetin  
MP1421 Tributyl phosphate  
MP1422 Tributylamine  
MP1423 Tributyrin  
MP1424 Trichloroethylene  
MP1425 Tridecane  
MP1426 Triethanolamine  
MP1427 Triethyl citrate  
MP1428 Triethylamine  
MP1429 Trioctylamine  
MP1430 Water  
MP1431 Xylenes Mixed  
MP1432 Castor Oil  
MP1433 Castor oil hydrogenated  
MP1434 Cedar Wood Oil  
MP1435 Clove oil  
MP1436 Coconut oil  
MP1437 Cottonseed Oil  
MP1438 Linseed Oil  
MP1439 Maize Oil  
MP1440 Mineral Oil  
MP1441 Mineral Oil, light and heavy  
MP1442 Oil of Pine  
MP1443 Oil of Thyme  
MP1444 Olive Oil  
MP1445 Paraffin Oil  
MP1446 Peanut Oil  
MP1447 Safflower Oil  
MP1448 Soybean Oil  
MP1449 Sunflower oil  
MP1450 (Methylphenyl)diphenylsilanol  
MP1451 1,10-Decanediol = Decamethylene glycol  
MP1452 1,2-Cycloheptanediol  
MP1453 1,2-Phenylenediamine  
MP1454 1,4-butanediol dimethacrylate  
MP1455 1,4-Cyclohexadiene  
MP1456 1,4-Cyclohexanedimethanol  
MP1457 1,7-Heptanediol = Heptamethylene glycol  
MP1458 1,9-Nonanediol  
MP1459 1-Butyldiphenylsilanol  
MP1460 1-Naphthylamine  
MP1461 2-Butyldiphenylsilanol  
MP1462 2-Butyne-1,4-Diol  
MP1463 2-Vinylnaphthalene  
MP1464 3,3-Dimethylacrylic acid  
MP1465 Acrylamide  
MP1466 Acrylic acid  
MP1467 Adipic acid  
MP1468 Allyl alcohol  
MP1469 Aluminum Ethoxide  
MP1470 Benzoic acid  
MP1471 Benzoyl Chloride  
MP1472 Bisphenol A  
MP1473 Calcium formate  
MP1474 Cetyl alcohol

MP1475 Chloromaleic anhydride  
MP1476 cis-1,4-Cyclohexanedimethanol  
MP1477 cis-1,5-Cyclooctanediol  
MP1478 Cyclohexyldiphenylsilanol  
MP1479 D-Glucose  
MP1480 Dicapryl Adipate  
MP1481 Dicyandiamide = Cyanoguanidine  
MP1482 Diethylenediamine  
MP1483 Diphenylamine  
MP1484 Diurethane Dimethacrylate  
MP1485 Dodecanedioic acid  
MP1486 Dodecylamine  
MP1487 D-Tyrosine  
MP1488 Ethyl Carbamate  
MP1489 Ethylenediamine  
MP1490 Formaldehyde 35%  
MP1491 Fumaric acid  
MP1492 Glucose  
MP1493 Glutaric acid  
MP1494 Glycerin  
MP1495 g-Valerolactone  
MP1496 Isocyanic acid  
MP1497 Isophthalic acid  
MP1498 Linoleic Acid  
MP1499 Maleic anhydride  
MP1500 Malonic Acid  
MP1501 Melamine = sym-Triaminotriazine  
MP1502 Methyl methacrylate  
MP1503 Methylamine Hydrochloride  
MP1504 mono-Methyl terephthalate  
MP1505 Neopentyl glycol  
MP1506 Octadecylamine  
MP1507 Octylamine  
MP1508 Oleoyl chloride  
MP1509 Oleyl alcohol, Cis-9-octadecen-1-ol  
MP1510 o-Phthalic acid  
MP1511 Phthalic acid  
MP1512 Phthalic acid, K salt  
MP1513 Phthalimide  
MP1514 Sebacic Acid  
MP1515 Sodium benzoate  
MP1516 Sodium butyrate  
MP1517 Sodium formate  
MP1518 Sodium oleate  
MP1519 Styrene oxide  
MP1520 Suberic acid  
MP1521 Succinic acid  
MP1522 Succinic anhydride  
MP1523 t-Butyldimethylsilyl chloride  
MP1524 Terephthalaldehyde  
MP1525 Terephthalic acid  
MP1526 Terephthalonitrile  
MP1527 Tetrahydrofurfurylalcohol  
MP1528 Titanium(IV) butoxide  
MP1529 trans-1,2-Cyclohexanediol  
MP1530 Triethylene glycol dimethacrylate  
MP1531 Triethylenetetramine

MP1532 Trimethylsilyl cyanide  
MP1533 Triphenylsilane  
MP1534 Triphenylsilanol  
MP1535 Urea  
MP1536 Urethane  
MP1537 Urotropine  
MP1538 Glue - Epoxy  
MP1539 Glue - Office White  
MP1540 Glue - Patex  
MP1541 Glue - Polychloroprene based  
MP1542 Glue - PVA based for wood  
MP1543 Glue - PVAC based for paper  
MP1544 Glue - PVAC based for wall tile  
MP1545 Cotton #1  
MP1546 Cotton #2  
MP1547 Cotton+Acrylic  
MP1548 Cotton+Elastane (92:8)  
MP1549 Cotton+Flax (60:40)  
MP1550 Cotton+Lycra (95:5)  
MP1551 Cotton+Polyester (1:1)  
MP1552 Cotton+Polyester (60:40)  
MP1553 Cotton+Polyester (65:35)  
MP1554 Cotton+Viscose (82:18)  
MP1555 Fabric Acrylic+Cotton (1:1)  
MP1556 Fabric Acrylic+Cotton+Polyester (44:42:14)  
MP1557 Fabric Angora+Wool+Nylon (20:60:20)  
MP1558 Fabric Cellophane based  
MP1559 Fabric Cellulose Acetate based  
MP1560 Fabric Poly(ethylene terephthalate)  
MP1561 Fabric Polyamide  
MP1562 Fabric Polyamide+Elastane (83:17)  
MP1563 Fabric Polyester  
MP1564 Fabric Polyester+Cotton (65:35)  
MP1565 Fabric Polyester+Lycra (9:1)  
MP1566 Fabric Polyester+Viscose (80:20)  
MP1567 Fabric Polymide  
MP1568 Fabric Viscon+Polyester+Lycra (5:4:1)  
MP1569 Moira  
MP1570 Silk I  
MP1571 Silk II  
MP1572 Silk III  
MP1573 Wool  
MP1574 Wool+Cotton+Acrylic+Nylon (32:23:25:20)  
MP1575 Adhesive strip napkin  
MP1576 Adhesive tape yellow  
MP1577 Artificial Fur  
MP1578 Artificial Leather #1  
MP1579 Artificial Leather #2  
MP1580 Artificial Leather #3  
MP1581 Artificial Leather #4  
MP1582 Artificial Leather #5  
MP1583 Artificial Leather #6  
MP1584 Artificial Leather #7  
MP1585 Beach Shoe China #1  
MP1586 Beach Shoe China #2  
MP1587 Beach Shoe China #3  
MP1588 Electronic board

MP1589 Furniture veneer  
MP1590 Insulation material car engine  
MP1591 Latex glove CaCO<sub>3</sub> coated  
MP1592 Lip salve polyvidone based  
MP1593 Medicinal capsule  
MP1594 Nail lack  
MP1595 Packing Material Ultrasoft  
MP1596 Plasticine for children  
MP1597 Polyurethane roller skate wheel  
MP1598 Rubber roller skate break  
MP1599 Rubber seal  
MP1600 Rubber white CaCO<sub>3</sub> filled  
MP1601 Shell button dark  
MP1602 Shoe cream  
MP1603 Silichrom  
MP1604 Ski running surface  
MP1605 Ski surface  
MP1606 Spaghetti insulation CaCO<sub>3</sub> filled #1  
MP1607 Spaghetti insulation CaCO<sub>3</sub> filled #2  
MP1608 Sympatex Dark  
MP1609 Sympatex White  
MP1610 Wadding  
MP1611 Waterpipe seal black  
MP1612 White Lack CD surface  
MP1613 White painting for walls  
MP1614 Window seal black  
MP1615 Yellow painting for walls

## **Part II – Additives, Coatings**

CT0001 Anti Terra - 204  
CT0002 Anti Terra - 204 (dried)  
CT0003 Anti Terra - U  
CT0004 Anti Terra - U (dried)  
CT0005 Byk - Dynwet 800  
CT0006 Byk 011  
CT0007 Byk 011 (dried)  
CT0008 Byk 012  
CT0009 Byk 017  
CT0010 Byk 018  
CT0011 Byk 019  
CT0012 Byk 020  
CT0013 Byk 020 (dried)  
CT0014 Byk 022  
CT0015 Byk 024  
CT0016 Byk 025  
CT0017 Byk 025 (dried)  
CT0018 Byk 028  
CT0019 Byk 037  
CT0020 Byk 037 (dried)  
CT0021 Byk 044  
CT0022 Byk 044 (dried)  
CT0023 Byk 052  
CT0024 Byk 052 (dried)

CT0025 Byk 055  
CT0026 Byk 055 (dried)  
CT0027 Byk 066  
CT0028 Byk 066 (dried)  
CT0029 Byk 070  
CT0030 Byk 070 (dried)  
CT0031 Byk 077  
CT0032 Byk 077 (dried)  
CT0033 Byk 088  
CT0034 Byk 088 (dried)  
CT0035 Byk 140  
CT0036 Byk 140 (dried)  
CT0037 Byk 141  
CT0038 Byk 142  
CT0039 Byk 142 (dried)  
CT0040 Byk 154  
CT0041 Byk 154 (dried)  
CT0042 Byk 220S  
CT0043 Byk 220S (dried)  
CT0044 Byk 300 (dried)  
CT0045 Byk 301  
CT0046 Byk 301 (dried)  
CT0047 Byk 302  
CT0048 Byk 306  
CT0049 Byk 306 (dried)  
CT0050 Byk 307  
CT0051 Byk 310  
CT0052 Byk 310 (dried)  
CT0053 Byk 320  
CT0054 Byk 320 (dried)  
CT0055 Byk 322  
CT0056 Byk 323  
CT0057 Byk 325  
CT0058 Byk 325 (dried)  
CT0059 Byk 330  
CT0060 Byk 330 (dried)  
CT0061 Byk 331  
CT0062 Byk 333  
CT0063 Byk 335  
CT0064 Byk 335 (dried)  
CT0065 Byk 337  
CT0066 Byk 337 (dried)  
CT0067 Byk 344  
CT0068 Byk 345  
CT0069 Byk 345 (dried)  
CT0070 Byk 346  
CT0071 Byk 346 (dried)  
CT0072 Byk 348  
CT0073 Byk 348 (dried)  
CT0074 Byk 354  
CT0075 Byk 354 (dried)  
CT0076 Byk 355  
CT0077 Byk 358  
CT0078 Byk 358 (dried)  
CT0079 Byk 361  
CT0080 Byk 373  
CT0081 Byk 373 (dried)

CT0082 Byk 380  
CT0083 Byk 380 (dried)  
CT0084 Byk 381  
CT0085 Byk 381 (dried)  
CT0086 Byk 390  
CT0087 Byk 390 (dried)  
CT0088 Byk 392  
CT0089 Byk 392 (dried)  
CT0090 Byk 410  
CT0091 Byk 410 (dried)  
CT0092 Byk 420  
CT0093 Byk 420 (dried)  
CT0094 Byk 425  
CT0095 Byk 4510  
CT0096 Byk 4510 (dried)  
CT0097 Byk A530  
CT0098 Byk A530 (dried)  
CT0099 Byk P 104  
CT0100 Byk P 104 (dried)  
CT0101 Byk P-105  
CT0102 Byk S 706  
CT0103 Byk S 706 (dried)  
CT0104 Byketol OK  
CT0105 Byketol OK (dried)  
CT0106 Byketol WS  
CT0107 Byketol WS (dried)  
CT0108 Disperbyk 108  
CT0109 Disperbyk 110  
CT0110 Disperbyk 110 (dried)  
CT0111 Disperbyk 111  
CT0112 Disperbyk 111 (dried)  
CT0113 Disperbyk 116  
CT0114 Disperbyk 161  
CT0115 Disperbyk 161 (dried)  
CT0116 Disperbyk 162  
CT0117 Disperbyk 162 (dried)  
CT0118 Disperbyk 163  
CT0119 Disperbyk 163 (dried)  
CT0120 Disperbyk 164  
CT0121 Disperbyk 164 (dried)  
CT0122 Disperbyk 166  
CT0123 Disperbyk 166 (dried)  
CT0124 Disperbyk 167  
CT0125 Disperbyk 167 (dried)  
CT0126 Disperbyk 170  
CT0127 Disperbyk 170 (dried)  
CT0128 Disperbyk 180  
CT0129 Disperbyk 180 (dried)  
CT0130 Disperbyk 181  
CT0131 Disperbyk 181 (dried)  
CT0132 Disperbyk 182  
CT0133 Disperbyk 184  
CT0134 Disperbyk 184 (dried)  
CT0135 Disperbyk 190  
CT0136 Disperbyk 190 (dried)  
CT0137 Disperbyk 192  
CT0138 Disperbyk 193

CT0139 Disperbyk 193 (dried)  
CT0140 Disperbyk 194  
CT0141 Disperbyk 194 (dried)  
CT0142 Lactimon - WS  
CT0143 Lactimon - WS (dried)  
CT0144 Nanobyk 3600  
CT0145 Nanobyk 3600 (dried)  
CT0146 Efka 2018 (dried)  
CT0147 Efka 2021  
CT0148 Efka 2021 (dried)  
CT0149 Efka 2023  
CT0150 Efka 2023 (dried)  
CT0151 Efka 2025  
CT0152 Efka 2025 (dried)  
CT0153 Efka 2028  
CT0154 Efka 2028 (dried)  
CT0155 Efka 2038  
CT0156 Efka 2040  
CT0157 Efka 2040 (dried)  
CT0158 Efka 2526  
CT0159 Efka 2526 (dried)  
CT0160 Efka 2532  
CT0161 Efka 2532 (dried)  
CT0162 Efka 2550  
CT0163 Efka 2722  
CT0164 Efka 2722 (dried)  
CT0165 Efka 3030  
CT0166 Efka 3031  
CT0167 Efka 3031  
CT0168 Efka 3031 (dried)  
CT0169 Efka 3032  
CT0170 Efka 3033  
CT0171 Efka 3033 (dried)  
CT0172 Efka 3034  
CT0173 Efka 3034 (dried)  
CT0174 Efka 3236  
CT0175 Efka 3239  
CT0176 Efka 3288  
CT0177 Efka 3500  
CT0178 Efka 3500 (dried)  
CT0179 Efka 3570  
CT0180 Efka 3570 (dried)  
CT0181 Efka 3580  
CT0182 Efka 3580 (dried)  
CT0183 Efka 3600  
CT0184 Efka 3772  
CT0185 Efka 3772 (dried)  
CT0186 Efka 3778  
CT0187 Efka 3778 (dried)  
CT0188 Efka 4046  
CT0189 Efka 4046 (dried)  
CT0190 Efka 4047  
CT0191 Efka 4047 (dried)  
CT0192 Efka 4060  
CT0193 Efka 4060 (dried)  
CT0194 Efka 4400  
CT0195 Efka 4400 (dried)

CT0196 Efka 4401  
CT0197 Efka 4401 (dried)  
CT0198 Efka 4402  
CT0199 Efka 4402 (dried)  
CT0200 Efka 4403  
CT0201 Efka 4500  
CT0202 Efka 4500 (dried)  
CT0203 Efka 4530  
CT0204 Efka 4530 (dried)  
CT0205 Efka 4550  
CT0206 Efka 4550 (dried)  
CT0207 Efka 4560  
CT0208 Efka 4560 (dried)  
CT0209 Efka 4580  
CT0210 Efka 4580 (dried)  
CT0211 Efka 5044  
CT0212 Efka 5044 (dried)  
CT0213 Efka 5054  
CT0214 Efka 5054 (dried)  
CT0215 Efka 5065  
CT0216 Efka 5066  
CT0217 Efka 5066 (dried)  
CT0218 Efka 5071  
CT0219 Efka 5071 (dried)  
CT0220 Efka 5207  
CT0221 Efka 5244  
CT0222 Efka 6220  
CT0223 Efka 6230  
CT0224 Efka 6700  
CT0225 Efka 6700 (dried)  
CT0226 Efka 6780  
CT0227 Efka 6780 (dried)  
CT0228 Efka 7022  
CT0229 Efka 7022 (dried)  
CT0230 Efka 8400  
CT0231 Efka 8400 (dried)  
CT0232 Efka 8503  
CT0233 Efka 8510  
CT0234 Efka 8510 (dried)  
CT0235 Efka 8530  
CT0236 Efka 8832  
CT0237 Efka 8832 (dried)  
CT0238 Efka 9009  
CT0239 Efka 9009 (dried)  
CT0240 Additol VXL 4930  
CT0241 Additol VXL 4930 (dried)  
CT0242 Additol VXL 6237  
CT0243 Additol VXL 6237 (dried)  
CT0244 Additol VXW 4909  
CT0245 Additol VXW 4909 (dried)  
CT0246 Additol VXW 4926  
CT0247 Additol VXW 4973  
CT0248 Additol VXW 4973 (dried)  
CT0249 Additol VXW 6200  
CT0250 Additol VXW 6200 (dried)  
CT0251 Additol VXW 6208  
CT0252 Additol VXW 6208 (dried)

CT0253 Additol VXW 6209  
CT0254 Additol VXW 6209 (dried)  
CT0255 Additol VXW 6211  
CT0256 Additol VXW 6213  
CT0257 Additol VXW 6213 (dried)  
CT0258 Additol VXW 6360  
CT0259 Additol VXW 6387  
CT0260 Additol VXW 6396  
CT0261 Additol VXW 6396 (dried)  
CT0262 Additol WL 490/50  
CT0263 Additol WL 490/50 (dried)  
CT0264 Additol XL 018  
CT0265 Additol XL 018 (dried)  
CT0266 Additol XL 123  
CT0267 Additol XL 123 (dried)  
CT0268 Additol XL 125  
CT0269 Additol XL 252  
CT0270 Additol XL 252 (dried)  
CT0271 Additol XL 255  
CT0272 Additol XL 255 (dried)  
CT0273 Additol XL 270  
CT0274 Additol XL 297, Oxime Based  
CT0275 Additol XL 6526  
CT0276 Additol XW 314  
CT0277 Additol XW 314 (dried)  
CT0278 Additol XW 390  
CT0279 Additol XW 390 (dried)  
CT0280 Additol XW 395  
CT0281 Additol XW 395 (dried)  
CT0282 Abiester 100P  
CT0283 Acropol 35  
CT0284 Acropol 35 (dried)  
CT0285 Amibal 145  
CT0286 AMP 95  
CT0287 Color Naturals 4.3 Goldbrown  
CT0288 Color Naturals 4.3 Goldbrown (after oxidation)  
CT0289 Color Naturals 4.3 Goldbrown (after oxidation) (dried)  
CT0290 Color Naturals 4.3 Goldbrown (dried)  
CT0291 Desmophen VP LS 2973  
CT0292 Desmophen VP LS 2973 (dried)  
CT0293 Dibenzoyl peroxide  
CT0294 Dioctyl Sodium Sulphosuccinate  
CT0295 Dow Corning 30, Silicone, Propylene carbonate  
CT0296 Dow Corning 66  
CT0297 Dow Corning 66  
CT0298 Edolan XCIC  
CT0299 Edolan XCIC (dried)  
CT0300 Ethyl palmitate  
CT0301 Halox Xtain L44 (Tannin-Stain Blocker)  
CT0302 Halox Xtain L44 (Tannin-Stain Blocker) (dried)  
CT0303 Isobutyl Stearate  
CT0304 Laurex CS, Cetyl Stearyl Alcohol  
CT0305 Lignine  
CT0306 Lukas Tinting Color Chrom Yellow Light 9507, Dried  
CT0307 Lukas Tinting Color Chrom Yellow Light 9507, Liquid  
CT0308 Lukas Tinting Color Light Ochre Yellow 9510, Dried  
CT0309 Lukas Tinting Color Umbra Greenish 9522, Dried

CT0310 Lukas Tinting Color Umbra Greenish 9522, Liquid  
CT0311 Maize Starch  
CT0312 Methyl-Ethyl Keton Peroxide  
CT0313 Metolat LA 571  
CT0314 Metolat LA 571 (dried)  
CT0315 Microcrystalline cellulose  
CT0316 Monoethanolamine  
CT0317 Neofix R-800  
CT0318 Nitrocellulose  
CT0319 Pentaerytritol tetra mercapto propionate  
CT0320 Perlit SI/SW  
CT0321 Perlit SI/SW  
CT0322 Persoftal LU  
CT0323 Polyoxyethylene Polymethylene Siloxane DC470A (dried)  
CT0324 Polyquart H81  
CT0325 Polyquart H81 (dried)  
CT0326 Protorez FFO  
CT0327 Protorez FFO  
CT0328 Silquest A 1100  
CT0329 Silquest A 1100 (dried)  
CT0330 Silquest A 1120  
CT0331 Silquest A 171  
CT0332 Silquest A 171 (dried)  
CT0333 Silquest A 187  
CT0334 Silquest A 187 (dried)  
CT0335 Silquest A 189  
CT0336 Silquest A 189 (dried)  
CT0337 Syntegal V-7  
CT0338 Tanede Nova  
CT0339 Tanede RC  
CT0340 Tinting Oil Color Brillant Ochre Light Nr. 1507, Liquid  
CT0341 Tinting Oil Color Burned Sienna Nr. 1510, Dried  
CT0342 Tinting Oil Color Burned Sienna Nr. 1510, Liquid  
CT0343 Tinting Oil Color Chrom Yellow Light Nr. 1504, Dried  
CT0344 Tinting Oil Color Chrom Yellow Light Nr. 1504, Liquid  
CT0345 Tinting Oil Color Ivory Extra Burning, Dried  
CT0346 Tinting Oil Color Ivory Extra Burning, Liquid  
CT0347 Tinting Oil Color Real Orange Nr. 1600, Dried  
CT0348 Tinting Oil Color Real Yellow Light Nr. 1598, Dried  
CT0349 Tinting Oil Color Real Yellow Light Nr. 1598, Liquid  
CT0350 Tinting Oil Color Umbra Brown Nr. 1516, Dried  
CT0351 Tinting Oil Color Umbra Brown Nr. 1516, Liquid  
CT0352 Tinting Oil Color Umbra Greenish Nr. 1648, Dried  
CT0353 Tinting Oil Color Umbra Greenish Nr. 1648, Liquid  
CT0354 Tinting Oil Color Zinc Green Light Nr.1533, Dried  
CT0355 Tinting Oil Color Zinc Green Light Nr.1533, Liquid  
CT0356 Tinting Oil Color Zinc Green Middle Nr. 1532, Dried  
CT0357 Triethanolamine  
CT0358 Tungophen B  
CT0359 Zoldine MS Plus  
CT0360 Dow Corning Z6137  
CT0361 Dow Corning Z6137  
CT0362 Lubrizol 2061  
CT0363 Lubrizol 2061 (dried)  
CT0364 Lubrizol 2062  
CT0365 Lubrizol 2062 (dried)  
CT0366 Lubrizol 2063

CT0367 Lubrizol 2063 (dried)  
CT0368 Solthix 250  
CT0369 Solthix 250 (dried)  
CT0370 Ancamide 2353  
CT0371 Ancamide 2353 (dried)  
CT0372 Ancamine 2280  
CT0373 Ancamine 2280 (dried)  
CT0374 Ancamine 2432  
CT0375 Ancamine 2432 (dried)  
CT0376 Vestamin TMD (branched aliphatic diamine)  
CT0377 Etoxon AF-5  
CT0378 Etoxon AF-5 (dried)  
CT0379 Statexan HA  
CT0380 Statexan HA (dried)  
CT0381 Statexan PES  
CT0382 Statexan PES  
CT0383 Tanastat P  
CT0384 Tanastat P  
CT0385 Tanastat PH  
CT0386 Tanastat PH  
CT0387 1,4-Diazabicyclo [2.2.2]octan  
CT0388 Dabco T 12 CL (Dibutyltin Dilaurate)  
CT0389 DDBSA Acid catalyst (Dodecylbenzene Sulfonic Acid)  
CT0390 DDBSA Acid catalyst (Dodecylbenzene Sulfonic Acid) (dried)  
CT0391 DEEA (Diethyl Ethanolamine)  
CT0392 DMDA (Dimethyldodecyl Amine)  
CT0393 DMP 30, 2,4,6-tris (dimethylaminomethyl) phenol  
CT0394 DMP 30, 2,4,6-tris (dimethylaminomethyl) phenol (dried)  
CT0395 Isononanoic acid  
CT0396 K Kat XC 5218  
CT0397 K Kat XC 5218 (dried)  
CT0398 K-Kat 4205 (pentane 2,4-dione)  
CT0399 Nacure 5076  
CT0400 Nacure 5076 (dried)  
CT0401 Tegokat 226  
CT0402 Tegokat 226 (dried)  
CT0403 Versamine EH 50  
CT0404 Versamine EH 50 (dried)  
CT0405 Lusolvan FBH (Di-isobutyl ester of a mixture of dicarboxylic acids) (dried)  
CT0406 Optifilm Enhancer 300 (dried)  
CT0407 Plastilit 3060 (Poly(propylene glycol alkylphenyl ether) (dried)  
CT0408 Nofome AF  
CT0409 Nofome AF  
CT0410 Tego Airex 1435  
CT0411 Tego Airex 1435 (dried)  
CT0412 Tego Airex 902 W  
CT0413 Tego Airex 902 W (dried)  
CT0414 Tego Airex 907 W  
CT0415 Tego Airex 907 W (dried)  
CT0416 Tego Airex 935  
CT0417 Tego Airex 935 (dried)  
CT0418 Tego Airex 980  
CT0419 Tego Airex 985  
CT0420 Tego Airex 985 (dried)  
CT0421 Tego Airex 986  
CT0422 Tego Airex 986 (dried)  
CT0423 Agitan 100

CT0424 Agitan 100 (dried)  
CT0425 Agitan 120  
CT0426 Agitan 120 (dried)  
CT0427 Agitan 150 (dried)  
CT0428 Agitan 151  
CT0429 Agitan 151 (dried)  
CT0430 Agitan 217  
CT0431 Agitan 218  
CT0432 Agitan 232  
CT0433 Agitan 255  
CT0434 Agitan 255 (dried)  
CT0435 Agitan 260  
CT0436 Agitan 260 (dried)  
CT0437 Agitan 280  
CT0438 Agitan 281  
CT0439 Agitan 281 (dried)  
CT0440 Agitan 285  
CT0441 Agitan 285 (dried)  
CT0442 Agitan 288  
CT0443 Agitan 288 (dried)  
CT0444 Agitan 290  
CT0445 Agitan 295  
CT0446 Agitan 295 (dried)  
CT0447 Agitan 296  
CT0448 Agitan 296 (dried)  
CT0449 Agitan 305  
CT0450 Agitan 315  
CT0451 Agitan 315 (dried)  
CT0452 Agitan 380  
CT0453 Agitan 380 (dried)  
CT0454 Agitan 633  
CT0455 Agitan 650  
CT0456 Agitan 700  
CT0457 Agitan 702  
CT0458 Agitan 703 N  
CT0459 Agitan 770  
CT0460 Agitan E 257 (dried)  
CT0461 Antimussol WLN  
CT0462 Antimussol WLN  
CT0463 Dapro DF 1181  
CT0464 Dapro DF 1492  
CT0465 Dapro DF 2162  
CT0466 Dapro DF 880  
CT0467 Dapro DF 881 S  
CT0468 Dapro DF 881 S (dried)  
CT0469 Dapro DF 900  
CT0470 Dapro DF 944  
CT0471 Dapro DF 975  
CT0472 Dapro NA 1621  
CT0473 Dapro NA 1621 (dried)  
CT0474 Dehydran 1208  
CT0475 Dehydran 1208 (dried)  
CT0476 Perenol E 1  
CT0477 Tego Foamex 1495  
CT0478 Tego Foamex 1495 (dried)  
CT0479 Tego Foamex 7447 (dried)  
CT0480 Tego Foamex 8030

CT0481 Tego Foamex 8030 (dried)  
CT0482 Tego Foamex 805  
CT0483 Tego Foamex 805 (dried)  
CT0484 Tego Foamex 815 (dried)  
CT0485 Tego Foamex 825  
CT0486 Tego Foamex 825 (dried)  
CT0487 Tego Foamex 830  
CT0488 Dapro W 95 HS  
CT0489 Disperse AYD W 22  
CT0490 Disperse AYD W 22 (dried)  
CT0491 Dispex N 40  
CT0492 Dispex N 40 (dried)  
CT0493 Kortamol NNO  
CT0494 Lopon 890 (dried)  
CT0495 Lopon PO  
CT0496 Lopon PO (dried)  
CT0497 Metolat FC 514  
CT0498 Metolat FC 514 (dried)  
CT0499 Metolat FC 515  
CT0500 Metolat FC 515 (dried)  
CT0501 Orotan 681  
CT0502 Orotan 681 (dried)  
CT0503 Orotan 850  
CT0504 Sol Plus D 540  
CT0505 Sol Plus D 540 (dried)  
CT0506 Solsperse 27000  
CT0507 Solsperse 41090  
CT0508 Zetasperse 1200  
CT0509 Zetasperse 1200 (dried)  
CT0510 Zetasperse 1600  
CT0511 Zetasperse 1600 (dried)  
CT0512 Zetasperse 2100  
CT0513 Zetasperse 2300  
CT0514 Levegal DLP  
CT0515 Levegal DLP  
CT0516 Tanasperse OLG  
CT0517 Tanasperse OLG (dried)  
CT0518 Additol VXW 4940  
CT0519 Additol VXW 4940 (dried)  
CT0520 Additol VXW 6240  
CT0521 Additol VXW 6240 (dried)  
CT0522 Astragal PAN  
CT0523 Astragal PAN (dried)  
CT0524 Levalin VKU-F  
CT0525 Levalin VKU-F  
CT0526 Levogen BF 01  
CT0527 Levogen BF 01 (dried)  
CT0528 Mesitol AS  
CT0529 Mesitol AS  
CT0530 Mesitol BFN  
CT0531 Mesitol BFN  
CT0532 Persoftal L  
CT0533 Persoftal L  
CT0534 Aluminium stearate  
CT0535 Brij 35, Polyoxyethylene (35) Lauryl Ether  
CT0536 Diethyl adipate  
CT0537 Diethylene Glycol Monostearate

CT0538 Disponil AAP 436  
CT0539 Disponil AAP 436 (dried)  
CT0540 Disponil SUS 875 IC (Anionic di-2-ethylhexyl sulphosuccinate, sodium salt)  
CT0541 Disponil SUS 875 IC (Anionic di-2-ethylhexyl sulphosuccinate, sodium salt)  
(dried)  
CT0542 Glycerol  
CT0543 Glyceryl monostearate  
CT0544 Linoleic Acid  
CT0545 Magnesium stearate  
CT0546 Methyl esters of fatty acid  
CT0547 Methyl Linoleate  
CT0548 Methyl Palmitate  
CT0549 Oleic acid  
CT0550 Oleyl alcohol  
CT0551 Palmitic acid  
CT0552 PEG 400 Monooleate  
CT0553 PEG 400 Monostearate  
CT0554 Pentaerythritol Tetraoleate  
CT0555 Polyglycerol esters of fatty acids  
CT0556 Polyoxyethene (40) stearate  
CT0557 Polysorbate 20  
CT0558 PPG Dioleate  
CT0559 Silwet L-77  
CT0560 Slovamine SK-7  
CT0561 Slovamix PC (dried)  
CT0562 Slovamix TZ  
CT0563 Slovamix TZ (dried)  
CT0564 Slovamix ZF1  
CT0565 Sodium oleate  
CT0566 Sodium potassium and calcium salts of fatty acids  
CT0567 Sorbitan monolaurate  
CT0568 Sorbitan monooleate  
CT0569 Sorbitan monopalmitate  
CT0570 Sorbitan monostearate  
CT0571 Stearic acid  
CT0572 Tween 60  
CT0573 Plexene CA  
CT0574 Plexene FED  
CT0575 Plexene FED  
CT0576 Aluminium potassium sulfate  
CT0577 Aluminium Sulfate  
CT0578 Ammonium carbonate  
CT0579 Ammonium chloride  
CT0580 Ammonium Phosphate  
CT0581 Bentonite  
CT0582 Bentonite  
CT0583 Bicarbonate of soda  
CT0584 Calcium aluminium silicate  
CT0585 Calcium chloride  
CT0586 Calcium oxide  
CT0587 Calcium silicate  
CT0588 Calcium sulfate  
CT0589 Concrete Screed  
CT0590 Copper Sulfate  
CT0591 Glass wool  
CT0592 Granite Grey, dried  
CT0593 Gypsum white

CT0594 Chamotte  
CT0595 Kaolin  
CT0596 Kaolin, Aluminium silicate  
CT0597 Lime Extra White  
CT0598 Magnesium carbonate  
CT0599 Magnesium oxide  
CT0600 Monopotassium Phosphate  
CT0601 Monosodium phosphate  
CT0602 Oxide Aluminium  
CT0603 Plaster of Paris  
CT0604 Potassium carbonate  
CT0605 Potassium hydrogen carbonate  
CT0606 Potassium silicate  
CT0607 Sand, Silicon dioxide  
CT0608 Silica Gel  
CT0609 Sodium aluminium silicate  
CT0610 Sodium carbonate  
CT0611 Sodium metasilicate  
CT0612 Sodium nitrate  
CT0613 Sodium orthosilicate  
CT0614 Sodium tetraborate (Borax)  
CT0615 Talc  
CT0616 Tricalcium phosphate  
CT0617 Bayhydur 2336 (Hydrophylic Isocyanate)  
CT0618 Bayhydur 3100, based on hexamethylene diisocyanate (dried)  
CT0619 Bayhydur VP LS 2150, based on Isophorone Diisocyanate (IPDI)  
CT0620 Bayhydur VP LS 2150, based on Isophorone Diisocyanate (IPDI) (dried)  
CT0621 Cyclohexyl isocyanate  
CT0622 Desmodur VPLS 2098  
CT0623 Desmodur VPLS 2098 (dried)  
CT0624 Desmodur W (cycloaliphatic diisocyanate) (dried)  
CT0625 Desmodur XP 2410 (aliphatic polyisocyanate)  
CT0626 Desmodur XP 2410 (aliphatic polyisocyanate) (dried)  
CT0627 HDI - monomer Hexamethylene diisocyanate (aged) (dried)  
CT0628 HDI - monomer Hexamethylene diisocyanate (dried)  
CT0629 Tego Flow 300  
CT0630 Tego Flow 300 (dried)  
CT0631 Tego Flow 425  
CT0632 Acetyl tributyl citrate  
CT0633 Dibutylphthalate  
CT0634 Diisodecyl phthalate  
CT0635 Dioctyl adipate  
CT0636 Kodaflex TXIB  
CT0637 Persoftal NOBLE  
CT0638 Persoftal NOBLE (dried)  
CT0639 Persoftal OE  
CT0640 Persoftal OE (dried)  
CT0641 Persoftal PEK  
CT0642 Persoftal PEK (dried)  
CT0643 Persoftal PLUS  
CT0644 Acetic acid  
CT0645 Acticide MKC  
CT0646 Acticide MKC (dried)  
CT0647 Acticide SR 2278  
CT0648 Acticide SR 2278 (dried)  
CT0649 Benzoic acid  
CT0650 Eulan SPA 01

CT0651 Eulan SPA 01 (dried)  
CT0652 Formaldehyde  
CT0653 Naphthalene  
CT0654 Parmetol DF 18  
CT0655 Parmetol DF 18 (dried)  
CT0656 Parmetol DF 35  
CT0657 Parmetol DF 35 (dried)  
CT0658 Parmetol HF 25  
CT0659 Parmetol HF 25 (dried)  
CT0660 Potassium acetate  
CT0661 Preventol CI-7  
CT0662 Preventol CI-8  
CT0663 Preventol D7 (Isothiazolinones) (dried)  
CT0664 Preventol MP 260  
CT0665 Preventol MP 260 (dried)  
CT0666 Preventol® ON extra  
CT0667 Triclosan  
CT0668 Cardolite NC 513  
CT0669 Cardura E5 (glycidyl neopentanoate)  
CT0670 Erisys GE-22S (dried)  
CT0671 ChS Epoxy RR 700 (dried)  
CT0672 K Flex M B 301  
CT0673 K Flex M B 301 (dried)  
CT0674 Polypox 480 (Modified cycloaliphatic polyamine)  
CT0675 Polypox 480 (Modified cycloaliphatic polyamine) (dried)  
CT0676 Polypox H 016 (Phenol-free polyamine)  
CT0677 Polypox H 016 (Phenol-free polyamine) (dried)  
CT0678 Polypox H 354 L (modified cycloaliphatic polyamine)  
CT0679 Polypox H 354 L (modified cycloaliphatic polyamine) (dried)  
CT0680 Polypox R 6 (o-cresyl glycidyl ether)  
CT0681 Polypox R 6 (o-cresyl glycidyl ether) (dried)  
CT0682 t-Butyl Glycidyl Ether  
CT0683 Zoldine RD 20  
CT0684 Acrysol RM 1020  
CT0685 Acrysol RM 1020 (dried)  
CT0686 Acrysol RM 12W  
CT0687 Acrysol RM 12W (dried)  
CT0688 Acrysol RM 2020  
CT0689 Acrysol RM 2020 (dried)  
CT0690 Acrysol RM 825  
CT0691 Acrysol RM 825 (dried)  
CT0692 Collacral PU 85  
CT0693 Collacral PU 85 (dried)  
CT0694 Rheolate 266 (dried)  
CT0695 Rheolate 278  
CT0696 Rheolate 310  
CT0697 Rheolate 310 (dried)  
CT0698 Rohagit® SD 15  
CT0699 Tafigel AP 1/1  
CT0700 Tafigel AP 1/1 (dried)  
CT0701 Tafigel AP10  
CT0702 Tafigel AP10 (dried)  
CT0703 Tafigel PUR 40  
CT0704 Tafigel PUR 45 (dried)  
CT0705 Tafigel PUR 50  
CT0706 Tafigel PUR 50 (dried)  
CT0707 Tafigel PUR 60 (dried)

CT0708 Tego Visco Plus 3000  
CT0709 Tego Visco Plus 3000 (dried)  
CT0710 Tego Visco Plus 3010  
CT0711 Tego Visco Plus 3010 (dried)  
CT0712 Tego Visco Plus 3060 (dried)  
CT0713 Diadavin ANE  
CT0714 Diadavin ANE (dried)  
CT0715 Diadavin DSP  
CT0716 Diadavin DSP (dried)  
CT0717 Diadavin EPF 01  
CT0718 Diadavin EPF 01 (dried)  
CT0719 Diadavin UN  
CT0720 Diadavin UN (dried)  
CT0721 Tanaterge SD  
CT0722 Tanaterge SD (dried)  
CT0723 Edaplan LA 402  
CT0724 Edaplan LA 402 (dried)  
CT0725 Edaplan LA 410  
CT0726 Edaplan LA 411  
CT0727 Edaplan LA 412  
CT0728 Edaplan LA 431  
CT0729 Edaplan LA 431 (dried)  
CT0730 Tego Glide 406  
CT0731 Tego Glide 406 (dried)  
CT0732 Tego Glide 410  
CT0733 Tego Glide 420  
CT0734 Tego SMO  
CT0735 2,6-Di-tert-Butyl-p-Cresol  
CT0736 4,4-thiobis(t-butyl-m-cresol)  
CT0737 Antioxidant 2246  
CT0738 Cyasorb UV 3346  
CT0739 Didodecyl 3,3'-thiodipropionate, DLTPD  
CT0740 Dioctadecyl 3,3'-thiodipropionate  
CT0741 Hostanox 03  
CT0742 Hostavin N30 (HALS)  
CT0743 Chimasorb 81, 2-hydroxy-4-n-octyloxybenzophenone (UV-531)  
CT0744 Chimasorb 944  
CT0745 Irgafos 168  
CT0746 Irganox 1010  
CT0747 Irganox 1010 + Irganox 168  
CT0748 Irganox 1076  
CT0749 Irganox 1330  
CT0750 Irganox 3114  
CT0751 Leucopur EGM  
CT0752 Millad 3988  
CT0753 Polyvinylpyrrolidone  
CT0754 Tinuvin 622  
CT0755 Uvinul 4050 H  
CT0756 Byk - LD N 20048  
CT0757 Byk - LD N 20048 (dried)  
CT0758 Dapro S 65  
CT0759 Dapro S 65 (dried)  
CT0760 Dapro U 99  
CT0761 Dapro U 99 (dried)  
CT0762 2-(2-Hydroxy-5-methyl phenyl) Benzotriazole 97%  
CT0763 2,2,4,4 -Tetrahydroxy-Benzophone 97%  
CT0764 2-[3-(2H-Benzotriazol-2-yl)-4-hydroxyphenyl]ethyl methacrylate

CT0765 3-(2H-Benzotriazol-2-yl)-4-Hydroxy Phenetyl Alcohol  
CT0766 Cyasorb UV-416 Light Absorber  
CT0767 Tinuvin 1130 (dried)  
CT0768 Tinuvin 123-DW, HALS (dried)  
CT0769 Tinuvin 123-DW, HALS (water dispersion)  
CT0770 Tinuvin 292 (dried)  
CT0771 Tinuvin 477-DW (dried)  
CT0772 Tinuvin 477-DW (water dispersion)  
CT0773 Tinuvin 5155  
CT0774 Abeson S, dodecylbenzenesulfonic acid based  
CT0775 Alresat KM-40, modified maleic resin  
CT0776 Dow Corning 67  
CT0777 Edaplan 472  
CT0778 Edaplan 472 (dried)  
CT0779 Edaplan 480  
CT0780 Edaplan 480 (dried)  
CT0781 Edaplan 482  
CT0782 Edaplan VP 470 (dried)  
CT0783 Erkantol NR  
CT0784 Erkantol NR (dried)  
CT0785 Lutensol A8, C12-14 fatty alcohol 8EO  
CT0786 Metolat FC 355  
CT0787 Metolat FC 388  
CT0788 Metolat LA 524  
CT0789 Metolat LA 524 (dried)  
CT0790 Nofome BLF 01  
CT0791 Nofome BLF 01 (dried)  
CT0792 Tanaterge CSU (dried)  
CT0793 Tanaterge CSU Liquid  
CT0794 Tanawet RC-N  
CT0795 Tanawet RC-N (dried)  
CT0796 Tego Dispers 710  
CT0797 Tego Dispers 710 (dried)  
CT0798 Tego Dispers 725 W  
CT0799 Tego Dispers 725 W (dried)  
CT0800 Tego Dispers 735 W  
CT0801 Tego Dispers 735 W (dried)  
CT0802 Tego Dispers 745 W  
CT0803 Tego Dispers 745 W (dried)  
CT0804 Tego Dispers 750 W  
CT0805 Tego Dispers 750 W (dried)  
CT0806 Tego Wet 265  
CT0807 Tego Wet 265 (dried)  
CT0808 Tego Wet 270  
CT0809 Tego Wet 500  
CT0810 Tego Wet 505  
CT0811 Tego Wet 510  
CT0812 Acramin BS-N  
CT0813 Acramin BS-N (dried)  
CT0814 Acronal 290D  
CT0815 Acronal 500 D  
CT0816 Acronal A 603 (dried)  
CT0817 Acronal DS 6219  
CT0818 Acronal DS 6219 (dried)  
CT0819 Acronal LR 8977 (dried)  
CT0820 Acronal LR 8988  
CT0821 Acronal LR 8998

CT0822 Acronal S 533  
CT0823 Acronal S 533 (dried)  
CT0824 Acronal S 726  
CT0825 Acronal S 726 (dried)  
CT0826 Alberdingk AC 2514 (dried)  
CT0827 Alberdingk AC 2518 (dried)  
CT0828 Alberdingk AC 2522 (dried)  
CT0829 Alberdingk AC 2523  
CT0830 Alberdingk AC 2523 (dried)  
CT0831 Alberdingk AC 31  
CT0832 Alberdingk AC 31 (dried)  
CT0833 Alberdingk AC 5530  
CT0834 Alberdingk AC 5530 (dried)  
CT0835 Alberdingk ACU 2000  
CT0836 Alberdingk ACU 2000 (dried)  
CT0837 Alberdingk AS 2688  
CT0838 Alberdingk AS 2688 (dried)  
CT0839 Alberdingk CUR 69  
CT0840 Alberdingk CUR 69 (dried)  
CT0841 Alberdingk CUR 991 VP  
CT0842 Alberdingk CUR 991 VP (dried)  
CT0843 Alberdingk MAC 34 (dried)  
CT0844 Alberdingk U 520  
CT0845 Alberdingk U 520 (dried)  
CT0846 Alberdingk U 800  
CT0847 Alberdingk U 800 (dried)  
CT0848 Alberdingk U 9800 VP  
CT0849 Alberdingk U 9800 VP (dried)  
CT0850 Alberdingk UC 90  
CT0851 Alberdingk UC 90 (dried)  
CT0852 Albidur EP 2240 A  
CT0853 Albidur EP 2240 A (dried)  
CT0854 Albiflex 348  
CT0855 Alkydal F 26  
CT0856 Alkydal F 26 (dried)  
CT0857 Alnovol PN-320, Phenol Novolak  
CT0858 Alpex CK-450, Cyclized rubber  
CT0859 Ancamide 2445  
CT0860 Ancamide 2445 (dried)  
CT0861 Ancamide 2573  
CT0862 Ancamide 2573 (dried)  
CT0863 Asphalt  
CT0864 Bakelite  
CT0865 Bay 2427  
CT0866 Bay 2427 (dried)  
CT0867 Bay 2470  
CT0868 Bay 2470 (dried)  
CT0869 Bay 2500  
CT0870 Bay 2500 (dried)  
CT0871 Bay XP 2593/1  
CT0872 Bayderm Finish 90 UD  
CT0873 Bayderm Grund  
CT0874 Bayderm Grund (dried)  
CT0875 Baypret USV  
CT0876 Baypret USV  
CT0877 Carboset CR 717  
CT0878 Carboset CR 717 (dried)

CT0879 Cardolite LX 5206  
CT0880 Cardolite LX 5206 (dried)  
CT0881 Cardolite NC 556X80  
CT0882 Cardolite NC 556X80 (dried)  
CT0883 Cardolite NC 562  
CT0884 Cardolite NC 562 (dried)  
CT0885 Cellulose Acetobutyrate, Celidor  
CT0886 Colophony  
CT0887 Colophony, Amber  
CT0888 Colophony, Rosin  
CT0889 Copolymer ABS type  
CT0890 Copolymer EPDM Type  
CT0891 Copolymer EVA type  
CT0892 Copolymer SAN type  
CT0893 Crodasil SC 650W  
CT0894 Crodasil SC 650W (dried)  
CT0895 Crosslinker CX 100 (polyfunctional aziridine) (dried)  
CT0896 Cyclohexanon Formaldehyde resin  
CT0897 Cymel 301  
CT0898 Cymel 301 (dried)  
CT0899 Cymel 327 (dried)  
CT0900 Dammar  
CT0901 Dammar Aged  
CT0902 Daotan VTW 1265/36WA  
CT0903 Daotan VTW 1265/36WA (dried)  
CT0904 DER 660 B-80 (dried)  
CT0905 Desmodur N 3200, hexamethylene diisocyanate based resin  
CT0906 Desmodur N 3200, hexamethylene diisocyanate based resin (dried)  
CT0907 Desmodur N 3600, hexamethylene diisocyanate based resin  
CT0908 Desmodur N 3600, hexamethylene diisocyanate based resin (dried)  
CT0909 Desmophen 1150  
CT0910 Desmophen 2009/1  
CT0911 Desmophen 2249/1  
CT0912 Desmophen 2249/1 (dried)  
CT0913 Desmophen 670 (dried)  
CT0914 Desmophen 680/1  
CT0915 Desmophen 680/1 (dried)  
CT0916 Desmophen 800  
CT0917 Desmophen 870 (dried)  
CT0918 Desmophen A 160 SN  
CT0919 Desmophen A 160 SN (dried)  
CT0920 Desmophen A 365  
CT0921 Desmophen A 365 (dried)  
CT0922 Desmophen A VP LS 2350  
CT0923 Desmophen A VP LS 2350 (dried)  
CT0924 Desmophen A VPLS 2158  
CT0925 Desmophen A VPLS 2158 (dried)  
CT0926 Desmophen A VPLS 2177  
CT0927 Desmophen A VPLS 2177 (dried)  
CT0928 Desmophen C1200  
CT0929 Desmophen VP LS 2008/1  
CT0930 Desmophen VP LS 2008/1 (dried)  
CT0931 Dinomin MI-11-1  
CT0932 Dinomin MI-11-1 (dried)  
CT0933 Dispercoll U 54  
CT0934 Dispercoll U 54 (dried)  
CT0935 Duroftal VPI 2801/78 BAC

CT0936 Duroftal VPI 2801/78 BAC (dried)  
CT0937 Duroftal VPI 2803/78 BAC  
CT0938 Duroftal VPI 2803/78 BAC (dried)  
CT0939 Duroxyn VAX 6127W  
CT0940 Duroxyn VAX 6127W (dried)  
CT0941 Efka 1500  
CT0942 Elvacite 2046, Copolymer Butylmetacrylate/Isobutylmethacrylate  
CT0943 Enorex M 1630AV  
CT0944 Enorex M 1630AV (dried)  
CT0945 Enorex Print 2030  
CT0946 Enorex Print 2030  
CT0947 Enorex TE 221  
CT0948 Enorex TE 221 (dried)  
CT0949 Epalloy 5200 (dried)  
CT0950 Epon 828  
CT0951 Epon 828  
CT0952 Fluorosilicone Rubber  
CT0953 Fluorotelomer Vidax 1000  
CT0954 Glyptal Resin  
CT0955 Halweftal B 32 75 BAC  
CT0956 Halweftal B 32 75 BAC (dried)  
CT0957 Hostaflex CN-131, Polyvinyl Acetate/Polyvinyl Chloride  
CT0958 Hydrocarbon resin  
CT0959 Hydroxyfunctional Alkyd from talow oil  
CT0960 Hydroxyfunctional Alkyd from talow oil (dried)  
CT0961 Chs Alkyd L56 (dried)  
CT0962 CHS Epoxy 210X75  
CT0963 ChS Epoxy Resin 590  
CT0964 Incorez W 2200  
CT0965 Incorez W 2200 (dried)  
CT0966 Incorez W 2600  
CT0967 Incorez W 2600 (dried)  
CT0968 Incorez W 830/140  
CT0969 Incorez W 830/140 (dried)  
CT0970 Incorez W 830/256  
CT0971 Incorez W 830/256 (dried)  
CT0972 Joncryl 1915  
CT0973 Joncryl 1915 (dried)  
CT0974 Joncryl 8224  
CT0975 Joncryl 8224 (dried)  
CT0976 Joncryl 8226  
CT0977 Joncryl 8226 (dried)  
CT0978 Joncryl 8284  
CT0979 Joncryl 8284 (dried)  
CT0980 Joncryl 8330  
CT0981 Joncryl 8330 (dried)  
CT0982 Joncryl HPD 71 E  
CT0983 Joncryl HPD 71 E (dried)  
CT0984 K Flex 188  
CT0985 Keton Resin N  
CT0986 Laropal A-101, Aldehyde based Resin  
CT0987 Laropal A-81, Ketone based Resin  
CT0988 Levaflex EP-390, polypropylene based  
CT0989 Liocryl AS 680  
CT0990 Liocryl AS 680 (dried)  
CT0991 Liopur 2004-119  
CT0992 Liopur 2004-119 (dried)

CT0993 Long Soya Bean Alkyd  
CT0994 Long Soya Bean Alkyd (dried)  
CT0995 Long Urethanized Alkyd  
CT0996 Long Urethanized Alkyd (dried)  
CT0997 Luhydran A 848S  
CT0998 Luhydran A 848S (dried)  
CT0999 Macrynal SM 51ON 60 LG  
CT1000 Macrynal SM 51ON 60 LG (dried)  
CT1001 Maincote HG-86ER  
CT1002 Maincote HG-86ER (dried)  
CT1003 Maincote PR 71 (dried)  
CT1004 Manila Copal  
CT1005 Melamine Resin  
CT1006 Melamine-formaldehyde Resin  
CT1007 Melamine-urea-formaldehyde resin  
CT1008 Melform 200 I  
CT1009 Methacrylic Acid  
CT1010 Modocoll EK-20, Cellulose modified  
CT1011 Morez 300, Styrene/Acrylic Resin  
CT1012 Mowilith 611  
CT1013 Mowilith 611 (dried)  
CT1014 Mowilith 7410 (dried)  
CT1015 Mowilith 7416  
CT1016 Mowilith 7416 (dried)  
CT1017 Mowilith 7460  
CT1018 Mowilith 7460 (dried)  
CT1019 Mowilith 760  
CT1020 Mowilith 760 (dried)  
CT1021 Mowilith 7714 (dried)  
CT1022 Mowilith 7717  
CT1023 Mowilith 7729  
CT1024 Mowilith 7729 (dried)  
CT1025 Mowilith DM 2452  
CT1026 Mowilith® 20, PVAc  
CT1027 Mowilith® DM 778  
CT1028 Mowilith® SDM 5  
CT1029 Mowilith® SDM 5 (dried)  
CT1030 Mowiol® 4/98  
CT1031 Mowiol® 4-88, PVA  
CT1032 MS2A-Resin  
CT1033 Natural Rubber  
CT1034 NeoCryl A-1049  
CT1035 NeoCryl A-1049 (dried)  
CT1036 NeoCryl A-1091  
CT1037 NeoCryl A-1091 (dried)  
CT1038 NeoCryl A-1120  
CT1039 NeoCryl A-1120 (dried)  
CT1040 NeoCryl BT-24  
CT1041 NeoCryl BT-26 (dried)  
CT1042 NeoCryl SC-160  
CT1043 NeoCryl SC-160 (dried)  
CT1044 NeoCryl XK-12  
CT1045 NeoCryl XK-15  
CT1046 NeoCryl XK-205  
CT1047 NeoCryl XK-205 (dried)  
CT1048 NeoCryl XK-37  
CT1049 NeoCryl XK-52

CT1050 NeoCryl XK-52 (dried)  
CT1051 NeoCryl XK-62  
CT1052 NeoCryl XK-62 (dried)  
CT1053 NeoCryl XK-75  
CT1054 NeoCryl XK-75 (dried)  
CT1055 NeoCryl XK-78  
CT1056 NeoCryl XK-87  
CT1057 NeoCryl XK-87 (dried)  
CT1058 NeoCryl XK-90  
CT1059 NeoCryl XK-95 (dried)  
CT1060 NeoCryl XK-98  
CT1061 NeoPac E-106  
CT1062 NeoPac E-106 (dried)  
CT1063 NeoPac E-111  
CT1064 NeoPac E-111 (dried)  
CT1065 NeoPac E-114  
CT1066 NeoPac E-114 (dried)  
CT1067 NeoPac E-121  
CT1068 NeoPac E-121 (dried)  
CT1069 NeoPac E-125  
CT1070 NeoPac E-125 (dried)  
CT1071 NeoPac E-130  
CT1072 NeoPac E-130 (dried)  
CT1073 NeoPac E-180  
CT1074 NeoRez R-2100  
CT1075 NeoRez R-2100 (dried)  
CT1076 NeoRez R-9403  
CT1077 NeoRez R-9407  
CT1078 NeoRez R-9407 (dried)  
CT1079 Paraformaldehyde  
CT1080 Paraldehyde  
CT1081 Paraloid B66  
CT1082 Paraloid B66 (dried)  
CT1083 Paraloid™ B 48 N  
CT1084 Paraloid™ B 66  
CT1085 Paraloid™ B 72  
CT1086 PEG 7 Glyceryl Cocoate  
CT1087 PEG-200 Glyceryl Tallowate  
CT1088 Phenol resin  
CT1089 Phenolic Resin  
CT1090 Phthalopal LR8525, Phthalate based Resin  
CT1091 Plexigum® PQ 611  
CT1092 Plexisol® P 550-40  
CT1093 Poloxamer 105 (PEG-PPG copolymer)  
CT1094 Poly(acrylonitrile:butadiene)  
CT1095 Poly(butadiene::styrene)  
CT1096 Poly(butyl acrylate)  
CT1097 Poly(butyl methacrylate)  
CT1098 Poly(ethyl acrylate)  
CT1099 Poly(ethyl methacrylate)  
CT1100 Poly(ethylene adipate)  
CT1101 Poly(ethylene dimethacrylate)  
CT1102 Poly(ethylene glycol)  
CT1103 Poly(ethylene glycol) 1000  
CT1104 Poly(ethylene terephthalate)  
CT1105 Poly(ethylene:propylene)  
CT1106 Poly(isodecyl methacrylate)

CT1107 Poly(methacrylic acid)  
CT1108 Poly(methyl methacrylate)  
CT1109 Poly(methyl methacrylate):: (butyl methacrylate)  
CT1110 Poly(propylene glycol) 2000  
CT1111 Poly(styrene:ethylacrylate)  
CT1112 Poly(styrene-ethylene-butylene)  
CT1113 Poly(vinyl butyral)  
CT1114 Poly(vinyl isobutyl ether)  
CT1115 Poly(vinyl methyl ketone)  
CT1116 Polyacrylamide  
CT1117 Polyacrylonitrile  
CT1118 Polyamide resin  
CT1119 Polybutadiene  
CT1120 Polydimethylsiloxane  
CT1121 Polyethylene  
CT1122 Polychloroprene  
CT1123 Polyisobutylene  
CT1124 Polyphenylene Oxide, Noryl  
CT1125 Polypropylene  
CT1126 Polystyrene  
CT1127 Polyurethane  
CT1128 Polyvinylbutyral  
CT1129 Polyvinylchloride - hard  
CT1130 Polyvinylchloride - soft  
CT1131 Primal CL-3371 (dried)  
CT1132 Primal E 2573  
CT1133 Primal E 2573 (dried)  
CT1134 Primal E-3188  
CT1135 Primal E-3188 (dried)  
CT1136 Primal EP 6060  
CT1137 Primal EP 6060 (dried)  
CT1138 Primal HG 800  
CT1139 Primal WS 24E  
CT1140 Primal® AC 35 (dried)  
CT1141 Primal® WS 24  
CT1142 Primal® WS 24 (dried)  
CT1143 PUR Desmophan 460  
CT1144 Regalrez® 1126, Hydrocarbon based Resin  
CT1145 Resin AFS  
CT1146 Resin AW-2, Polycyclohexanone based Resin  
CT1147 Resin CA  
CT1148 Resin MH  
CT1149 Resin MH (dried)  
CT1150 Resin SK, Modified Acetophenone/Formaldehyde Resin  
CT1151 Resydrol AX 237W  
CT1152 Resydrol AX 237W (dried)  
CT1153 Resydrol AY 241W/40WA  
CT1154 Resydrol AY 241W/40WA (dried)  
CT1155 Resydrol AY 430W  
CT1156 Resydrol AY 430W (dried)  
CT1157 Resydrol AY 586W/38WA  
CT1158 Resydrol AY 586W/38WA (dried)  
CT1159 Resydrol UEA 209  
CT1160 Resydrol UEA 209 (dried)  
CT1161 Rheovis CR2, Polyalkyl acrylic acid based Resin  
CT1162 Rhodopas DS 850  
CT1163 Rhodopas DS 850 (dried)

CT1164 Sancure 20041  
CT1165 Sancure 20041 (dried)  
CT1166 SCX 910  
CT1167 SCX 910 (dried)  
CT1168 SCX 920  
CT1169 SCX 920 (dried)  
CT1170 Setal 168 SS 80  
CT1171 Setal 168 SS 80 (dried)  
CT1172 Setalux 6741 AQ 43  
CT1173 Setalux 6741 AQ 43 (dried)  
CT1174 Setalux 6742 AQ 40  
CT1175 Setalux 6742 AQ 40 (dried)  
CT1176 Setalux C 1151 XX 51  
CT1177 Setalux C 1151 XX 51 (dried)  
CT1178 Setalux C 1152 XX 51  
CT1179 Setalux C 1156 XS 60  
CT1180 Setalux C 1156 XS 60 (dried)  
CT1181 Setalux C 1182 SS 55  
CT1182 Setalux C 1182 SS 55 (dried)  
CT1183 Setaqua 6520  
CT1184 Setaqua 6520 (dried)  
CT1185 Setaqua 6755 (dried)  
CT1186 Setaqua 6756  
CT1187 Setaqua 6763  
CT1188 Setaqua 6763 (dried)  
CT1189 Setaqua 6771  
CT1190 Setaqua 6771 (dried)  
CT1191 Setaqua 6774  
CT1192 Setaqua 6774 (dried)  
CT1193 Setaqua 6776 (dried)  
CT1194 Setyrene 78 XS 55  
CT1195 Setyrene 78 XS 55 (dried)  
CT1196 Shellack Aged  
CT1197 Shellack Angelo 4621 Superblonde  
CT1198 Short Soya Bean Alkyd  
CT1199 Short Soya Bean Alkyd (dried)  
CT1200 Silikophen N P 50/300  
CT1201 Silikopon EF, Silicone-epoxy Resin  
CT1202 Silikopur 8001 F  
CT1203 Silikopur 8001 F (dried)  
CT1204 Sokrat 412  
CT1205 Styrofan D 422  
CT1206 Styrofan D 422 (dried)  
CT1207 Sulryn random copolymer poly(ethylene-co-methacrylic acid), EMAA  
CT1208 Synocure 220 HS  
CT1209 Synocure 220 HS (dried)  
CT1210 Synolac 1333  
CT1211 Synolac 1333 (dried)  
CT1212 Synthalan LS 768  
CT1213 Synthalan LS 768 (dried)  
CT1214 Teflon  
CT1215 Umacryl KI, Methylmethacrylate - styrene copolymer  
CT1216 Uracr ON CY 474 E (dried)  
CT1217 Uracron CY 474 E  
CT1218 Uralac SN 805 X-70  
CT1219 Uralac SN 805 X-70 (dried)  
CT1220 Uramex MF 821 B (dried)

CT1221 Viacryl SC 166  
CT1222 Viacryl VSC 6279W  
CT1223 Viacryl VSC 6279W (dried)  
CT1224 Viacryl VSC 6295W  
CT1225 Viacryl VSC 6295W (dried)  
CT1226 Vialkyd AC 406/75X  
CT1227 Vialkyd AC 406/75X (dried)  
CT1228 Vialkyd VAF 7200/70 BAC (dried)  
CT1229 Vinnapas CEP 15  
CT1230 Vinnapas CEP 15 (dried)  
CT1231 Vinnapas LL 6484  
CT1232 Vinyl Short Alkyd  
CT1233 Vinyl Short Alkyd (dried)  
CT1234 WorléeCryl 7930  
CT1235 WorléeCryl 7930 (dried)  
CT1236 Worlesin MS-240, Maleic resin  
CT1237 Zytel E103 HSLC NC-10  
CT1238 Abimal 125D, Glycerol ester + colophony + maleic  
CT1239 Alkaprene, Chloroprene based Adhesive  
CT1240 Refill Glue ER-S  
CT1241 Castor Oil  
CT1242 Coleseed Oil  
CT1243 Dehydrogenated Castor Oil (dried)  
CT1244 Dow Corning 41  
CT1245 Fluorolube  
CT1246 Linseed Oil  
CT1247 Mineral Oil, light and heavy  
CT1248 Oil of Pine  
CT1249 Oil of Thyme  
CT1250 Bergamot Oil  
CT1251 Caraway oil  
CT1252 Cedar Wood Oil  
CT1253 Cod-Liver Oil  
CT1254 Eucalyptus Oil  
CT1255 Lavender Oil  
CT1256 Marine fish oil  
CT1257 Sesame Oil  
CT1258 Yoyoba Oil  
CT1259 Owatrol Oil  
CT1260 Paraffin Oil  
CT1261 Perilla Oil  
CT1262 Polishing Oil  
CT1263 Poppy Oil  
CT1264 Soya Bean Oil  
CT1265 Tagat CH 60, PEG-60 Hydrogenated Castor Oil  
CT1266 Tallow Oil  
CT1267 Tung Oil  
CT1268 Walnut Oil  
CT1269 1,1,2,2-Tetrachloroethane  
CT1270 1,1,2-Trichlorotrifluoroethane  
CT1271 1,2-Dichloroethane  
CT1272 1,6-Hexandiol  
CT1273 1-Propanol  
CT1274 Acetone  
CT1275 Acetonitrile  
CT1276 Benzene  
CT1277 Benzyl acetate

CT1278 Benzyl alcohol  
CT1279 Benzyl benzoate  
CT1280 Butyl acetate  
CT1281 Carbon Tetrachloride  
CT1282 Cetyl alcohol  
CT1283 Cyclododecan  
CT1284 Cyclohexane  
CT1285 Diethyl ether  
CT1286 Diethylene Glycol Dibutyl Ether  
CT1287 Dichloromethane  
CT1288 Dimethylsulfoxide  
CT1289 Dodecane  
CT1290 Dowanol DPM, Dipropylene Glycol Methyl Ether  
CT1291 Dowanol EPH, 2-phenoxyethanol  
CT1292 Dowanol PNP, Propylene glycol n-propyl ether  
CT1293 Dowanol PPH, Propylene Glycol Phenyl Ether  
CT1294 Ethanol  
CT1295 Ethyl acetate  
CT1296 Ethylbenzene  
CT1297 Ethylene glycol  
CT1298 Ethylene glycol dimethyl ether  
CT1299 Ethylene glycol monoethyl ether  
CT1300 Heptane  
CT1301 Hexane  
CT1302 Hexanol  
CT1303 Chloroform  
CT1304 Isobutyl alcohol  
CT1305 m-Cresol  
CT1306 N-methylpyrrolidone  
CT1307 o-Cresol  
CT1308 o-Xylene  
CT1309 p-Cresol  
CT1310 Pentane  
CT1311 Petroleum ether  
CT1312 Phenol  
CT1313 Propylene glycol  
CT1314 Tetrahydrofuran  
CT1315 Tetrachloroethylene  
CT1316 Toluene  
CT1317 TPM, Tripropylene Glycol (Mono) Methyl Ether  
CT1318 TPNB, Tripropylene Glycol n-Butyl Ether  
CT1319 Triethylene Glycol  
CT1320 Triethylene Glycol Monobutylether  
CT1321 Water  
CT1322 Xylenes Mixed  
CT1323 2,4,7,9-Tetramethyl-5-Decyne-4,7-Diol  
CT1324 Alcopol O-70 PG, Aqueous/PPG contains dioctyl sulfosuccinate sodium  
CT1325 Altaran S-8  
CT1326 Altaran S-8 (dried)  
CT1327 Disponil A 3065 (modified ethoxylated fatty alcohols POE 30)  
CT1328 Disponil A 3065 (modified ethoxylated fatty alcohols POE 30) (dried)  
CT1329 Disponil A 4065 EXP (modified ethoxylated fatty alcohols POE 40)  
CT1330 Disponil A 4065 EXP (modified ethoxylated fatty alcohols POE 40) (dried)  
CT1331 Disponil A 5560 EXP (modified ethoxylated fatty alcohols) (dried)  
CT1332 Dow Corning 57  
CT1333 Dynasytan HS 2926  
CT1334 Dynasytan HS 2926 (dried)

CT1335 Empicol 0216, Oleoyl etherphosphoric acid  
CT1336 Empilan AMT 2, Alkylamine Ethoxylate  
CT1337 Empilan BQ 100, Fatty acid ethoxylate  
CT1338 Empilan BQ4/S, Ethoxylated Propoxylated fatty alcohol  
CT1339 Empilan KB-1, Ethoxylated C12-C14 1-EO  
CT1340 Empilan KB2/S, Lauryl ethoxylate (2EO)  
CT1341 Empilan KM11, Alcohol ethoxylate EO11  
CT1342 Empilan, Ethoxylated Propoxylated fatty alcohol  
CT1343 Empiphos Alkylphosphoric acid C16-C18/5632  
CT1344 Emulgator WNS, Arylpolyglykoether (dried)  
CT1345 Fluorad FC-129  
CT1346 Fluowet OTN  
CT1347 Fluowet OTN (dried)  
CT1348 Glycerol Trioleate  
CT1349 Marlophen 87 P  
CT1350 Nicca Sunsolt LM-7  
CT1351 Nicca Sunsolt RF-557  
CT1352 Nicca Sunsolt RF-557 (dried)  
CT1353 Sodium Dodecyl Sulfate  
CT1354 Spanscour Spark  
CT1355 Spanscour Spark (dried)  
CT1356 Surfynol 440 (dried)  
CT1357 Tanaterge DRC  
CT1358 Tanaterge DRC  
CT1359 TANEDE CL  
CT1360 Tannex BCA-N  
CT1361 Tannex BCA-N (dried)  
CT1362 Tego Twin 4000  
CT1363 Triton X-100 (dried)  
CT1364 Zonyl 9361  
CT1365 Zonyl 9361 (dried)  
CT1366 Zonyl FS 62  
CT1367 Zonyl FS 62 (dried)  
CT1368 Zonyl FSA  
CT1369 Zonyl FSA (dried)  
CT1370 Zonyl FSH  
CT1371 Zonyl FSH (dried)  
CT1372 Zonyl FSJ  
CT1373 Zonyl FSJ (dried)  
CT1374 Zonyl FSN  
CT1375 Zonyl FSN (dried)  
CT1376 Zonyl FSO 100  
CT1377 Zonyl FSO 100 (dried)  
CT1378 Zonyl FSP  
CT1379 Zonyl FSP (dried)  
CT1380 Byk Cera Aquacer 526  
CT1381 Byk Cera Aquacer 526 (dried)  
CT1382 Fluowax Pink  
CT1383 Fluowax Pink (dried)  
CT1384 Glycerol Tristearate  
CT1385 Chlorinated paraffin  
CT1386 Microcrystalline wax  
CT1387 Persoftal NPF  
CT1388 Persoftal NPF (dried)  
CT1389 Triacetin  
CT1390 Trilaurin  
CT1391 Trimyristin

CT1392 Wax V  
CT1393 1K Clearcoat  
CT1394 1K Clearcoat (dried)  
CT1395 Acryl White Universal Paint  
CT1396 Activ Wood Stain Colorless  
CT1397 Alkyd Lacquere Transparent  
CT1398 Amarol Triol Paint for Wood Aldyde Resin based  
CT1399 Anticorrosion Primer Paint, Binder part, Dried  
CT1400 Anticorrosion Primer Paint, Mixed, Dried  
CT1401 Anticorrosive Priming Paint  
CT1402 Antikor 200 Paint  
CT1403 Antirezin Black Paint  
CT1404 Asphaltum Laquer  
CT1405 Asphaltum Laquer (dried)  
CT1406 Balakryl Uni Mat Paint 0750 Orange  
CT1407 Balakryl Uni Mat Paint 0100 White  
CT1408 Balakryl Uni Mat Paint 0815 Red  
CT1409 Balakryl Uni Mat Paint 0510 Light Green  
CT1410 Balakryl Uni Mat Paint 0660 Ochre  
CT1411 Balakryl Uni Mat Paint 0220 Brown Light  
CT1412 Balakryl Uni Mat Paint 0250 Palisander  
CT1413 Balatix Paint, Tixotropic Enamel  
CT1414 Base Lissante Sheer Ivory Base Coat (for nails)  
CT1415 Base Lumere Protective Basecoat (for nails)  
CT1416 Base Paint for Metall, Glossy part, dried, from 1960  
CT1417 Base Paint for Metall, Matt Part, Dried, from 1960  
CT1418 Base Paint, from 1960  
CT1419 Bisil Light Blue Paint  
CT1420 Black Matt Paint for Tennis Tables, from 1960  
CT1421 Black Protective and Isolation Paint, Dried  
CT1422 Black Protective and Isolation Paint, Liquid  
CT1423 Black Toner MT 105B  
CT1424 Bondex Profi Wood Stain  
CT1425 Bondex Wood Care Product Natural Mahagony  
CT1426 Bondex Wood Care Product, Dried  
CT1427 Bondex Wood Care Product, Liquid  
CT1428 Bootlack, Paint for Boats  
CT1429 Brass Clear Varnish for Metal, Dried  
CT1430 Brass Clear Varnish for Metal, Liquid  
CT1431 Brillux Lacryl Glossy Finish Paint 275 RAL 1015  
CT1432 Bronze Goddess AJ (for nails)  
CT1433 C1036 Paint 0000, Nitrocellulose based  
CT1434 C1037 Celolesk 0000, Paint for wood, nitrocellulose based  
CT1435 C1038 Celomat Paint 0000, Wood Stain, nitrocellulose based  
CT1436 C2000 Paint  
CT1437 C2000 Paint (dried)  
CT1438 C5000 Sealant 110, Nitrocellulose based  
CT1439 Celox Paint  
CT1440 Cewesol Enamel Binder part, Liquid  
CT1441 Cewesol Enamel Settlement part, Dried  
CT1442 Cewesol Enamel Settlement part, Liquid  
CT1443 Cleaner 205 Bonding Agent  
CT1444 Clipper I Penguin Yachting Paint  
CT1445 Color-Matic Spray Paint Red  
CT1446 Consolan Wood Stain  
CT1447 Copper Rose Nail Lacquere  
CT1448 Cova Deep, Wall Paint Interior

CT1449 Cova Medium, Wall Paint Interior  
CT1450 Cova Pastel, Wall Paint Interior  
CT1451 Cova White, Wall Paint Interior  
CT1452 Cova Yellow, Wall Paint Interior  
CT1453 Cowacryl Acrylic Paint Waterthinnable  
CT1454 Cowidur Glossy Enamel, PU based, Component A+B, Aged A  
CT1455 Cowidur Glossy Enamel, PU based, Component A+B, Aged B  
CT1456 Cowidur Glossy Enamel, PU based, Component A+B, Liquid  
CT1457 Cowidur Hardener  
CT1458 DD Primer Paint for Yachts Grey  
CT1459 DD-Hartlack, PU based Lacquer for Yachts, Component A, Dried  
CT1460 DD-Hartlack, PU based Lacquer for Yachts, Component A, Liquid  
CT1461 DD-Hartlack, PUR based Lacquer for Yachts, Component B, Dried  
CT1462 DD-Hartlack, PUR based Lacquer for Yachts, Component B, Liquid  
CT1463 DD-Hartlack, PUR based Lacquer for Yachts, Components A+B mixed, Dried  
CT1464 DD-Hartlack, PUR based Lacquer for Yachts, Components A+B mixed, Liquid  
CT1465 Degrolux Synthetic Resin Based Enamel Paint, from 1960  
CT1466 Degronit 767 Floor Paint, from 1960  
CT1467 Deltalac Paint with Epicote, Binder part, Dried  
CT1468 Deremaqua Paint (dried)  
CT1469 Dermaline Nail Lacquere  
CT1470 Deyellowing Serum (for nails)  
CT1471 Dibenol (binder), from 1960  
CT1472 Dinocryl Dispersion Paint  
CT1473 Ducolux Enamel  
CT1474 Ducolux Metal Paint, Binder part  
CT1475 Ducolux Metal Paint, Settlement  
CT1476 Dulux Enamel Plus Red, Paint Silicone based  
CT1477 Duralit 1516 - 029, Industrial Paint for Metal  
CT1478 Duralit 1516 - 029, Industrial Paint for Metal (dried)  
CT1479 Durodur 3057-025, Industrial Paint for Metal  
CT1480 Durodur 3057-025, Industrial Paint for Metal (dried)  
CT1481 Eco Super White Enamel Inside, from 1960  
CT1482 Eggshell Paint  
CT1483 Ekohost Colorless Matt 0001, Paint urethane-acrylate based, for wood inside  
CT1484 Ekohost Colorless Semigloss 0000, Paint urethane-acrylate based, for wood inside  
CT1485 Enamel Spray for Cars Acrylate Combination Blue  
CT1486 Enamel Spray for Cars Acrylate Combination Blue (dried)  
CT1487 Enamel Spray for Cars Acrylate Combination Yellow (dried)  
CT1488 Enamel Spray for Cars Acrylate Green  
CT1489 Enamel Spray for Cars Acrylate Red  
CT1490 Enamel Spray for Cars Acrylate Red (dried)  
CT1491 Epolex Colorless Glossy Paint 0000  
CT1492 Epoxide Based Lacquere Two componets  
CT1493 Erbedol Semigloss Paint A + I  
CT1494 Erbedol Synthetic Resin Based, Impact-Proof Paint  
CT1495 Eternal Mat Revital Paint 217/201/218 Orange  
CT1496 Eternal Mat Revital Paint 217/216 Green  
CT1497 Eternal Mat Revital Paint 206/216 Dark Green  
CT1498 Eternal Mat Revital Paint 205/216 Blue-Green  
CT1499 Eternal Mat Revital Paint 205/202/217 Brown-Yellow  
CT1500 Eternal Mat Revital Paint 201/216/218 Purple  
CT1501 Eternal Mat Revital Paint 201/213 Grey  
CT1502 Eternal Paint 01 White  
CT1503 Eternal Paint 05 Yellow  
CT1504 Eternal Paint 07 Red-Brown

CT1505 Eternal Paint 08 Brick Red  
CT1506 Eternal Paint 14 Ivory  
CT1507 Eternal Paint 17 Light yellow  
CT1508 Eternal Paint 18 Strawberry Red  
CT1509 Eternal Paint 22 Dark Green  
CT1510 Euro Alkyd Based Lacquere Glossy  
CT1511 Express Finish 260 Exotic nail Lacquere  
CT1512  
Exterior Car Car Paint Chip

CT1513 Facade Paint TD  
CT1514 Fasadin, Facade Paint  
CT1515 Finadur 781 Stoving Varnish (dried)  
CT1516 Finadur 781, Stoving Varnish  
CT1517 Finishing Paint Alkyd Based, Matt, P199M Black  
CT1518 Finishing Paint Alkyd Based, Gloss, P6700 Ochre  
CT1519 Finishing Paint Alkyd Based, Matt, P100M White  
CT1520 Finishing Paint Alkyd Based, Gloss, P5321 Green  
CT1521 Finishing Paint Alkyd Based, Gloss, P6006 Ivory  
CT1522 Finishing Paint Alkyd based, wood and metal  
CT1523 Formex Extra Palisander 0280, Anticorrosion Single Layer Paint, for zinc  
CT1524 Formex Extra Sandy 0670, Anticorrosion Single Layer Paint, for zinc  
CT1525 Formex Extra White-Grey 0101, Anticorrosion Single Layer Paint, for zinc  
CT1526 Fortissimo, Wall Paint Interior  
CT1527 Glasurit Oil Color with Synthetic Resin  
CT1528 Glasurit Oil Color with Synthetic Resin  
CT1529 Gloss Paint  
CT1530 Glossy Colour Enamel Acrylic based, Dried  
CT1531 Glossy Colour Enamel Acrylic based, Liquid  
CT1532 Glossy Colour Enamel Alkyd based, Binder part, Dried  
CT1533 Glossy Colour Enamel Alkyd based, Mixed, Dried  
CT1534 Goldlack Synthetic Resin based Lacquere for effect, Binder part, Dried  
CT1535 Goldlack Synthetic Resin based Lacquere for effect, Binder part, Liquid  
CT1536 Goldlack Synthetic Resin based Lacquere for effect, Mixed, Dried  
CT1537 Goldlack Synthetic Resin based Lacquere for effect, Mixed, Liquid  
CT1538 Graphitti Spray Black - White Mate  
CT1539 Graphitti Spray Gold - White Creme  
CT1540 Graphitti Spray Molotow Premium - White  
CT1541 Graphitti Spray Pocket - White  
CT1542 GS Extra Mat Nail Lacquere  
CT1543 Gussofarb 66 Boat Paint, Dried  
CT1544 Gussofarb 66 Boat Paint, Liquid  
CT1545 H2203 Paint  
CT1546 Hairstyling Spray Ultralong Lasting  
CT1547 Hammerite Paint for Metal Green  
CT1548 Hammerite Paint for Metal Grey  
CT1549 Hammerite Paint for Metal Silver  
CT1550 Hammerite Paint for Metal Red  
CT1551 Hammerite Paint for Metal Blue  
CT1552 Hammerite Paint for Metal Dark Green  
CT1553 Herbol - Offenporig, Glazing Paint  
CT1554 Hetline Paint  
CT1555 High Build Stain  
CT1556 High Build Wood Stain  
CT1557 High Glossy Lacquere White  
CT1558 Hostagrund Aluminium 0910, Anticorrosion Single Layer Paint  
CT1559 Hostagrund Black 0199, Anticorrosion Single Layer Paint

CT1560 Hostagrund Brown 0243, Anticorrosion Single Layer Paint  
CT1561 Hostagrund Dark Blue 0455, Anticorrosion Single Layer Paint  
CT1562 Hostagrund Dark Green 0540, Anticorrosion Single Layer Paint  
CT1563 Hostagrund Light Blue 0440, Anticorrosion Single Layer Paint  
CT1564 Hostagrund Red-Brown 0840, Anticorrosion Single Layer Paint  
CT1565 Hostagrund White 0100, Anticorrosion Single Layer Paint  
CT1566 Hostagrund Yellow 6200, Anticorrosion Single Layer Paint  
CT1567 Hyundai Bumper ZW Dark Grey Effect Car Paint Chip  
CT1568 Icosit Dickschicht 5530 RAL 1014, Anticorrosion Paint  
CT1569 Icosit Dickschicht 5530 RAL 1014, Anticorrosion Paint  
CT1570 Impra Color Colorless Wood Stain  
CT1571 Impredur Glossy Lacquer 840 Alkyd Resin based, Dried  
CT1572 Impredur Glossy Lacquer 840 Alkyd Resin based, Liquid  
CT1573 Impredur Silk Glossy Lacquer 880 Alkyd Resin based, Dried  
CT1574 Impredur Silk Glossy Lacquer 880 Alkyd Resin based, Liquid  
CT1575 Industral Synthetic Paint for Outside - Brown  
CT1576 Industrol Black 1999, Paint Alkyd based  
CT1577 Industrol Blue 4265, Paint Alkyd based  
CT1578 Industrol Dark Red 8191, Paint Alkyd based  
CT1579 Industrol Green 5300, Paint Alkyd based  
CT1580 Industrol Ivory 6003, Paint Alkyd based  
CT1581 Industrol Light Red 8140, Paint Alkyd based  
CT1582 Industrol Orange 7551, Paint Alkyd based  
CT1583 Industrol Red-Brown 8440, Paint Alkyd based  
CT1584 Industrol Synthetic Paint for Outside - White  
CT1585 Industrol Yellow 6201, Paint Alkyd based  
CT1586 Inertol Paint for steel and concrete  
CT1587 Interior LK Rich Black Car Paint Chip  
CT1588 Interior Primer Paint, Dried, from 1960  
CT1589 Interior Primer Paint, Liquid, from 1960  
CT1590 Interior TP Paint  
CT1591 Interior ZK Beige Car Paint Chip  
CT1592 Izoban, Paint for swimming pools  
CT1593 Izoprim, Wall Paint Interior  
CT1594 Kronalux Street Marking Paint, Liquid PVC  
CT1595 Kunstschmiede Mattlack, PU-Alkyd Resin based Lacquer, Settlement part, Dried  
CT1596 Laksil Semiglossy Paint  
CT1597 Lazuro Colorless 0000, Alkyd Based Wood Stain with Oil and Wax  
CT1598 Lazuro Ebony 0099, Alkyd Based Wood Stain with Oil and Wax  
CT1599 Lazuro Green 0051, Alkyd Based Wood Stain with Oil and Wax  
CT1600 Lazuro Chestnut 0020, Alkyd Based Wood Stain with Oil and Wax  
CT1601 Lazuro Mahagony 0080, Alkyd Based Wood Stain with Oil and Wax  
CT1602 Lazuro Nut 0021, Alkyd Based Wood Stain with Oil and Wax  
CT1603 Lazuro Palisander 0022, Alkyd Based Wood Stain with Oil and Wax  
CT1604 Lazuro Pine Dark 0060, Alkyd Based Wood Stain with Oil and Wax  
CT1605 Lazuro Pine Light 0062, Alkyd Based Wood Stain with Oil and Wax  
CT1606 Lazuro Sipo 0025, Alkyd Based Wood Stain with Oil and Wax  
CT1607 Lazuro White 0010, Alkyd Based Wood Stain with Oil and Wax  
CT1608 Lazuro, Wood Stain, Light Oak  
CT1609 Le Vernis Multi Eclat 222 (for nails)  
CT1610 Le Vernis Nail Colour 461 Blue Satin  
CT1611 Le Vernis Nail Colour 465 Azur  
CT1612 Long Lasting Nail Lacquere 30  
CT1613 Luxol Extra Paint S1023, Wood Stain  
CT1614 Malak Primer Paint for Windows, from 1960  
CT1615 Malvena, Facade Paint  
CT1616 Matt, Synthetic Resin Based Paint, from 1960

CT1617 Meisterpreis Clear Varnish Seidenmatt with Polyurethane, Gel part, Dried  
CT1618 Meisterpreis Euro Ventilations, White Paint for Wood, Aldyde Resin Based  
CT1619 Meisterpreis Super Silky Polyurethane based Paint, Dried  
CT1620 Meisterpreis Super Silky Polyurethane based Paint, Liquid  
CT1621 Meisterpreis Tixoton High Build Wood Stain, Dried  
CT1622 Meisterpreis Tixoton High Build Wood Stain, Liquid  
CT1623 Meisterpreis Tixoton, High Build Stain  
CT1624 Meisterpreis White Enamel  
CT1625 ML AC 300 Stoving Varnish (dried)  
CT1626 ML AC 300, Stoving Varnish  
CT1627 Nail Art Multiglitter Tip Painter  
CT1628 Nail Enamel White  
CT1629 Nail Hardener  
CT1630 Nail Lacquere for Childern Strawberry  
CT1631 Nail Polish No. 48  
CT1632 Nail Stimulating Serum  
CT1633 Nissan Infinity Bumper AG05 Medium Grey Effect Car Paint Chip  
CT1634 O2004 Paint Undercoat for metals, industrial, outside, oil based  
CT1635 Oil Ground - Preparation for Oil Painting  
CT1636 Oknobal Palisander 2800, Paint Urethanised Alkyd Glossy  
CT1637 Oknobal White 1000, Paint Urethanised Alkyd Glossy  
CT1638 Oldorit Anticorrosive Paint Redbrown, Binder part, Dried  
CT1639 Oldorit Anticorrosive Paint Redbrown, Mixed, Dried  
CT1640 Oldorit Anticorrosive Paint Redbrown, Mixed, Liquid  
CT1641 Original French manicure for Nails  
CT1642 P05 Paint Remover  
CT1643 Pinotex Mahagony, Wood Care Product, Dried  
CT1644 Pinotex Mahagony, Wood Care Product, Liquid  
CT1645 Plus, Wall Paint Interior  
CT1646 Polar, Wall Paint Interior  
CT1647 Pragoprimer, Anticorrosion Paint Primer  
CT1648 Praktik Gloss Paint 9110  
CT1649 Praktik Primer Paint 100  
CT1650 Primalex Procolor S-PL, Paint for wood, base without pigments  
CT1651 Primer Paint 609 White, Dried, from 1960  
CT1652 Primer Paint 609 White, Liquid, from 1960  
CT1653 Primer Paint TP  
CT1654 Priming Paint For Wood Alkyd Based  
CT1655 Protective Concrete Paint, Dried  
CT1656 Protective Concrete Paint, Liquid  
CT1657 Quick Dry Top Coat For Nails  
CT1658 Quick Shine For Nails 510 Alin  
CT1659 Radbal White 1000, Paint Acrylate Based for Radiators  
CT1660 Rampolux Synthetic Resin Based Paint  
CT1661 Rampolux Synthetic Resin Based Paint, Dried  
CT1662 Rampolux Synthetic Resin Based Paint, Liquid  
CT1663 Rampotex Chlorinated Rubber based Paint  
CT1664 Rampozet Liquid Plastics  
CT1665 Ridgefiller (for nails)  
CT1666 Ripolin Ripoderm Primer Paint for Wood, Dried  
CT1667 Ripolin Ripoderm Primer Paint for Wood, Settlement part, Liquid  
CT1668 Ritual Color 13 Fatale (for nails)  
CT1669 Rostosam Anticorrosion Primer Paint, Dried  
CT1670 Rostosam Anticorrosion Primer Paint, Liquid  
CT1671 Rustic Satin High Build Stain, Binder part, Liquid  
CT1672 Rustic Wood Stain, Liquid, from 1960  
CT1673 S2028 Paint

CT1674 S2070 - Syntecol Primer 0100, Paint for wood, Alkyd based  
CT1675 S2222 Profi Radiator 1000, Alkyd/urethane based paint, for radiators  
CT1676 S2819 Paint  
CT1677 S2852 Krastenol 0100, Paint Polystyrene based for swimming pools  
CT1678 Sadosol, Wood Care Product  
CT1679 Sadotec, Wood Care Product  
CT1680 Seidenmattlack Matte Lacquere for Outside Black  
CT1681 Semigloss Enamel, Binder part, Liquid  
CT1682 Semigloss Enamel, Settlement part, Liquid  
CT1683 Semigloss Enamel, Settlement part, Liquid  
CT1684 SemiGloss Paint 530  
CT1685 S-he Stylezone Nail Polish  
CT1686 Sheer Shimmer Nail Colour 105  
CT1687 Sheer Strength Nail Foundation  
CT1688 Sigmastar Alkyde Resin based Paint  
CT1689 Specolor Priming Undercoat Interior  
CT1690 Specolux Titan Glossy Enamel  
CT1691 Specosal Glossy Color Enamel RAL 7035  
CT1692 Stove Pipe Paint Black Glossy for high Temperatures, Dried  
CT1693 Stove Pipe Paint Black Glossy for high Temperatures, Liquid  
CT1694 Stove Pipe Paint Black Glossy, Dried  
CT1695 Stove Pipe Paint Black Glossy, Liquid  
CT1696 Strahlenweiss Alkyd Resin Anticorrosion Paint  
CT1697 Supra Lux Synthetic Resin based Lacquere, Dried  
CT1698 Supra Lux Synthetic Resin based Lacquere, Liquid  
CT1699 Supra Oil Colors Glossy Enamel, Dried  
CT1700 Supra Oil Colors Glossy Enamel, Liquid  
CT1701 Swing Color, Anticorrosion Paint  
CT1702 Sylabal Colorless Matt 0001, Paint Alkyd based, transparent, for metal and wood inside  
CT1703 Teolin High Glossy Lacquere  
CT1704 Thick Layer Glazing Paint for Wood  
CT1705 Thin Wood Stain  
CT1706 Top Coat for Nails  
CT1707 Topcoat Paint 100 Synthetic  
CT1708 Topcoat Paint 100 Water Based  
CT1709 Transparent Sealer for Wood, from 1960  
CT1710 Trinity Nail Lacquere  
CT1711 Twinkle Shine Nail Lacquere  
CT1712 Undercoat Paint  
CT1713 Underhood Car Paint Chip 6517 Vermillon  
CT1714 Underhood Car Paint Chip 6519 Black  
CT1715 Underhood Car Paint Chip 6534 Titanium Frost  
CT1716 Underhood Car Paint Chip 6579 Grey  
CT1717 Underhood Car Paint Chip 6726 Red  
CT1718 Underhood Car Paint Chip 6901 Spruce Green  
CT1719 Underhood Car Paint Chip 6980 Chestnut  
CT1720 Underhood Car Paint Chip 7058 Dark SVT Blue  
CT1721 Underhood Car Paint Chip CS30 Silky White  
CT1722 Underhood Car Paint Chip Z26 Red  
CT1723 Underhood Car Paint Chip ZOC Light Grey  
CT1724 Unibunt Synthetic Resin Based Paint RAL 9001, from 1960  
CT1725 Unimat Black Paint 1999, Topcoat Matt  
CT1726 Unimat White Paint 1000, Topcoat Matt  
CT1727 Universal Brown Paint Urethane/Alkyd based  
CT1728 Universal Dispersion for Paint to enhance Elasticity  
CT1729 V 1314 Stoving Varnish

CT1730 V 1314 Stoving Varnish (dried)  
CT1731 V1600 Paint  
CT1732 V1605 Polyurex Mat 0000, Paint for wood, PU/acrylate based, flat  
CT1733 V1640 Paint  
CT1734 V2018 Latex Universal 1000, Wall Paint Interior  
CT1735 V2020 Paint  
CT1736 V2025 Okapol 842, Paint for Eaves  
CT1737 V2026 Paint  
CT1738 V2045 Balakryl Uni Mat 220, Paint for wood and walls, acrylate based, flat  
CT1739 V2068 Balakryl Gloss 1000, Paint for wood and walls, acrylate based  
CT1740 Velvet Orchid Nail Lacquere  
CT1741 Venticolor Window Enamel  
CT1742 Ventilux Paint for Windows, Dried  
CT1743 Ventilux Paint for Windows, Liquid  
CT1744 Vernis 001 (for nails)  
CT1745 Vernis 019 (for nails)  
CT1746 Viscolor Paint, from 1960  
CT1747 Viscolor Synthetic Resin based Enamel, Binder part, Dried  
CT1748 VP 169 Stoving Varnish  
CT1749 VP 169 Stoving Varnish (dried)  
CT1750 Waterproof Paint for Tiles, Isobutyl methacrylate based  
CT1751 Weissglut Paint 650  
CT1752 Wellaflex Hair Spray  
CT1753 Wilkenit Synthetic Resin Lacquere, Binder part, Dried  
CT1754 Wilkenit Synthetic Resin Lacquere, Mixed, Dried  
CT1755 Window Enamel White, from 1960  
CT1756 Wood Care Product Colorless  
CT1757 Wood Primer Paint  
CT1758 Wood Stain TC  
CT1759 Xylabrillant Wood Stain  
CT1760 Xyladecor 200 Wood Stain  
CT1761 Xyladecor Blue Wood Stain  
CT1762 Xylamon Colorless Wood Worm Killer  
CT1763 Yacht Systeme, Paint for Yachts  
CT1764 Powder Coating 19L, Epoxypolyester, glossy, RAL 1021 (powder)  
CT1765 Powder Coating 19L, Epoxypolyester, glossy, RAL 3020 (powder)  
CT1766 Powder Coating 3303, Epoxypolyester, matt, RAL 9005 (powder)  
CT1767 Powder Coating 4007, Acrylate, silk glossy, RAL 6005 (cured)  
CT1768 Powder Coating 4007, Acrylate, silk glossy, RAL 6005 (powder)  
CT1769 Powder Coating 500, Polyester, silk glossy, transparent (powder)  
CT1770 Powder Coating 504, Polyester, matt, RAL 9005 (powder)  
CT1771 Powder Coating 58L, Polyester, silk glossy, RAL 1023 (powder)  
CT1772 Powder Coating 58L, Polyester, silk glossy, RAL 3002 (powder)  
CT1773 Powder Coating 58L, Polyester, silk glossy, RAL 6005 (powder)  
CT1774 Powder Coating 58L, Polyester, silk glossy, RAL 9006 (powder)  
CT1775 Powder Coating 8009, Polyurethane, glossy, transparent (powder)  
CT1776 Powder Coating B23, Epoxide, matt, RAL 9005 (cured)  
CT1777 Powder Coating B23, Epoxide, matt, RAL 9005 (powder)  
CT1778 Powder Coating B44, Epoxide, matt, RAL 9005 (cured)  
CT1779 Powder Coating B44, Epoxide, matt, RAL 9005 (powder)  
CT1780 Powder Coating H50, Epoxypolyester, glossy, transparent (powder)  
CT1781 Matte Overprint Varnish V5080 (dried), Offset Print Color  
CT1782 Matte Overprint Varnish V5080, Offset Print Color  
CT1783 Novaspot Pantone Black (dried), Offset Print Color  
CT1784 Novavit Bio HKS44K (dried), Offset Print Color  
CT1785 Novavit Bio HKS44K, Offset Print Color  
CT1786 Pantone 186GA Red (dried), Offset Print Color

CT1787 Pantone 186GA Red, Offset Print Color  
CT1788 Pantone 8383 Gold (dried), Offset Print Color  
CT1789 Pantone Black Natural 10009 (dried), Offset Print Color  
CT1790 Pantone Green Natural 10008 (dried), Offset Print Color  
CT1791 Pantone Green Natural 10008, Offset Print Color  
CT1792 Pantone Orange 10021 (dried), Offset Print Color  
CT1793 Pantone Orange 10021, Offset Print Color  
CT1794 Pantone Process Blue 10007 (dried), Offset Print Color  
CT1795 Pantone Process Blue 10007, Offset Print Color  
CT1796 Pantone Purple Oxy (dried), Offset Print Color  
CT1797 Pantone Purple Oxy, Offset Print Color  
CT1798 Pantone Red 10032 (dried), Offset Print Color  
CT1799 Pantone Red 10032, Offset Print Color  
CT1800 Pantone Reflex Blue 10006 (dried), Offset Print Color  
CT1801 Pantone Reflex Blue 10006, Offset Print Color  
CT1802 Pantone Transparent White (dried), Offset Print Color  
CT1803 Pantone Transparent White, Offset Print Color  
CT1804 Pantone Violet Natural 10005 (dried), Offset Print Color  
CT1805 Pantone Violet Natural 10005, Offset Print Color  
CT1806 Pantone Warm Red 10001 (dried), Offset Print Color  
CT1807 Pantone Warm Red 10001, Offset Print Color  
CT1808 Pantone Yellow Natural 10000 (dried), Offset Print Color  
CT1809 Pantone Yellow Natural 10000, Offset Print Color  
CT1810 Print Varnish Glossy 10LW9500 (dried), Offset Print Color  
CT1811 Print Varnish Glossy 10LW9500, Offset Print Color  
CT1812 Turbo Board Plus Cyan 80B32 (dried), Offset Print Color  
CT1813 Turbo Board Plus Magenta 80B31 (dried), Offset Print Color  
CT1814 Turbo Board Plus Magenta 80B31, Offset Print Color  
CT1815 Turbo Board Plus Yellow 80B30 (dried), Offset Print Color  
CT1816 Turbo Board Plus Yellow 80B30, Offset Print Color

Nicodom IR Polymers Plus [www.ir-spectra.com](http://www.ir-spectra.com)