



NICODOM IR Fibers – List of Spectra

The infrared spectral library „NICODOM IR Fibers“ is a unique collection of 457 FTIR (ATR) digital spectra of fibers.

Examples of fiber types:

Acetate, Acrylic, Animal, Natural, Nylon, Olefin, Polyester, Rayon, Mixed Fibers, and most other fiber types.

The spectra were collected using the FTIR Spectrometer Nexus 670™ (Thermo) and the Miracle™ single bounce ATR accessory (Pike Technologies) equipped by Si crystal. The spectral library features resolution of 4cm⁻¹ in the spectral range 600-4000 cm⁻¹, measurement time 1 minute, the instrument was purged by dried air.

This library is available as digital searchable library compatible with your spectroscopic software. Also it can be delivered with SEARCH software or in *.pdf format.

This spectral library is available as a part of the package “Nicodom IR Polymers Extended” or “Nicodom IR Polymers Advanced”

To download the free Nicodom IR Demo Library (50 spectra) compatible with your software, to check the format compatibility as well as for pricelist, ordering info and information about all NICODOM IR/NIR Libraries and other NICODOM products please visit our webpage.

List of spectra

FB001, Comfortrel, Polyester, Wellman
FB002, Dry Step, Nylon 6, Royal-American
FB003, Arvil, Rayon, Avtex
FB004, Avisco, Acetate, Avisco
FB005, Wool Sheep Columbia, Wool, Animal
FB006, Dream, 13% Angora, 18% Rayon, 50% Acrylic, 10% Nylon, Almedahls
FB007, Feather Pigeon, Animal, NA
FB008, Silk India Assam Eri, Silk, Ethno
FB009, Lycra, Spandex, Invista
FB010, Spider Net (Araneus diadematus), Animal, NA
FB011, Fur Leopard, Animal, NA
FB012, Minnie, 23% Kid Mohair, 16% Merino, 38% Nylon, 23% Polyester, Classic Elite
FB013, Lin, 25% Flax, 25% Cotton, 50% Acrylic, Almedahls
FB014, Brighton , 30% Cotton, 28% Nylon, 24% Rayon, 18% Nylon, Trendsetter yarns
FB015, Fur Chinchilla, Animal, NA
FB016, Vlisseline, Mixed Yarn, Alysse
FB017, Kevlar 29, Aramid , Dupont
FB018, Chirlac Braided (surgery), Biodegradable, PLLA, Chirmac
FB019, Perlon, Nylon, Bayer
FB020, Sisal, Natural, Dorner Helmer
FB021, Hair Human Male, Human, NA
FB022, Barricaut, Polyester, Honeywell
FB023, Pelouche, 30% Nylon, 35% Acrylic, 35% wool, Bollicine
FB024, Firenze, 30% Wool, 30% Acrylic, 40% Nylon, Plymouth Brand
FB025, Sultana, 31% Acrylic, 29% Cotton, 40% Nylon, Berger Du Nord
FB026, Savvy Onyx, 34% Polyester, 56% Acrylic, 10% Cotton, Douglas Furniture
FB027, Yukon, 35% Mohair, 35% Wool, 30% Acrylic, Plymouth Brand
FB028, Pueblo Sky, 35% Polyester, 32% Olefin, 26% Cotton, 6% Acrylic, Douglas Furniture
FB029, Moda Color , 35% Rayon, 35% Nylon, 30% Acrylic, Austermann
FB030, For Keeps Sand, 37% Cotton, 63% Olefin, Douglas Furniture
FB031, Bamba , 38% Wool, 36% Nylon, 26% Polyester, Plymouth Brand
FB032, Caramelle, 40% Cotton, 34% Polyester, 26% Acrylic, Almedahls
FB033, Lavish , 40% Nylon, 32% Wool, 15% Polyester, 13% Acrylic, Berroco
FB034, Mystery , 40% Wool, 25% Cotton, 20% Acrylic, 15% Nylon, Artful Yarns
FB035, Bushlark Pecan, 41% Acrylic, 24% Polyester, 35% Polypropylene, Douglas Furniture
FB036, Clearview Beige, 41% Polyester, 20% Acrylic, 39% Polypropylene, Douglas Furniture
FB037, All Terrain Moss, 41% Polypropylene, 33% Acrylic, 26% Polyester, Douglas Furniture
FB038, Softwist, 41% Wool, 59% Rayon, Berroco
FB039, Olympic , 42% Rayon, 29% Cotton, 29% Acrylic, Artful Yarns
FB040, Chameleon , 42% Wool, 55% Acrylic, 3% Polyester, Wendy
FB041, Fox, 42% Wool, 58% Nylon, Aurora Yarns
FB042, Haitian Natural, 43% Cotton, 26% Rayon, 25% Polyester, 6% Acrylic, Douglas Furniture
FB043, Zaballa Multi, 43% Polyester, 41% Acrylic, 16% olefin, Douglas Furniture

NICODOM IR Fibers, 457 FTIR spectra of fibers

Copyright © NICODOM 2007

NICODOM Ltd., Hlavni 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

ftir@ftir.cz

FB044, Tapestry , 43% Rayon, 41% Nylon, 16% Polyester, Trendsetter yarns
FB045, Snap , 43% Wool, 37% Acrylic, 20% Nylon, Plymouth Brand
FB046, Spillo , 43% Wool, 55% Acrylic, 2% Rayon, Plymouth Brand
FB047, Prentice Froth, 44% Acrylic, 31% Polypropylene, 21% Polyester, 4% Cotton, Douglas Furniture
FB048, Caravan Multi, 44% Acrylic, 35% Polyester, 21% Cotton, Douglas Furniture
FB049, Tinker Confetti, 44% Polyester, 39% Acrylic, 17% Polypropylene, Douglas Furniture
FB050, Natasha , 45% Acrylic, 24% Nylon, 17% Rayon, 14% Polyester, Trendsetter yarns
FB051, Clown, 45% Cotton, 40% Superwash Wool, Almedahls
FB052, Sockotta Sock Yarn, 45% Cotton, 40% Superwash Wool, 15% Nylon, Plymouth Brand
FB053, Commando Blue, 45% Herculon, 40% Polyester, 15% Acrylic, Douglas Furniture
FB054, Frou Frou, 45% Nylon, 28% Cotton, 27% Acrylic, Plymouth Brand
FB055, Wishmaster Sage, 45% Olefin, 41% Polyester, 14% Cotton, Douglas Furniture
FB056, Wishmaster Beige, 45% Polyester, 41% Olefin, 14% Cotton, Douglas Furniture
FB057, Bettina , 45% Rayon, 35% Nylon, 20% Flax, Adrienne Vittadini
FB058, Chiara, 45% Rayon, 37% Nylon, 9% Cotton, 9% Acrylic, Adrienne Vittadini
FB059, Conifer Pine, 47% Acrylic, 25% Polypropylene, 20% Polyester, Douglas Furniture
FB060, Tuscany Bronze, 47% Polyester, 33% Acrylic, 17% Polypropylene, Douglas Furniture
FB061, Bonton Clay, 47% Polyester, 37% Acrylic, 16% Olefin, Douglas Furniture
FB062, Summer Stock Multi, 47% Polyester, 41% Acrylic, 12% Cotton, Douglas Furniture
FB063, Batika Color , 48% Cotton, 47% Acrylic, 5% Nylon, Austermann
FB064, Filbert Camel, 48% Polyester, 31% Acrylic, 21% Polypropylene, Douglas Furniture
FB065, Lauderdale Ebony, 48% Polyester, 32% Acrylic, 17% Polypropylene, 3% Nylon, Douglas Furniture
FB066, Palace , 48% Wool, 45% Acrylic, 7% Nylon, Artful Yarns
FB067, Fluffy, 49% Acrylic, 29% Wool, 19% Nylon, 3% Polyester, Almedahls
FB068, Speedway Springtime, 49% Polyester, 34% Acrylic, 17% Polypropylene, Douglas Furniture
FB069, Vabene , 5% Polyester, 48% Nylon, 26% Wool, 21% Acrylic, Austermann
FB070, Dayton Sage, 50% Acrylic, 22% Polyester, 28% Polypropylene, Douglas Furniture
FB071, Fireworks, 50% Acrylic, 38% Polyester, 12% Nylon, Karabella
FB072, Flax Isle, 50% Cotton, 30% Rayon, 20% Flax , Plymouth Brand
FB073, Madrid, 50% Cotton, 50% Rayon, Filatura Cervinia
FB074, Bom-Ull, 50% Cotton, 50% Wool, Almedahls
FB075, Tiffany , 50% Nylon, 31% Mohair, 19% Acrylic, Trendsetter yarns
FB076, Moisselle , 51% Acrylic, 30% Nylon, 9% Wool, Wendy
FB077, Bamboo Garden , 51% Bamboo, 47% Cotton, 2% Nylon, Plymouth Brand
FB078, Athena, 51% Cotton, 44% Nylon, 5% Rayon, Plymouth Brand
FB079, Massimo, 51% Wool, 27% Acrylic, 15% Nylon, 7% Polyester , Plymouth Brand
FB080, Antonia , 51% Wool, 49% Nylon, Adrienne Vittadini
FB081, Coventry, 52% Nylon, 18% Wool, 18% Acrylic, 12% Rayon, Trendsetter yarns
FB082, Mohican Beige, 52% Polyester, 43% Olefin, 5% Cotton, Douglas Furniture
FB083, Arlington Berry, 52% Polyester, 48% Cotton, Douglas Furniture
FB084, Foliage, 53% New Wool, 47% Acrylic, Berroco
FB085, Monarch , 53% Nylon, 47% Acrylic, Trendsetter yarns
FB086, Runway, 53% Nylon, 47% Polyester, Berroco
FB087, Lucia, 53% Wool, 47% Cotton, Adrienne Vittadini

FB088, Flame Navy, 54% Polyester, 46% Acrylic, Douglas Furniture
FB089, Max red, 55% Polyester, 35% Acrylic, 10% Polypropylene, Douglas Furniture
FB090, Toreador Multi, 55% Polyester, 45% Rayon, Lane Borgosesia
FB091, Iceland, 55% Wool, 35% Acrylic, 10% Nylon, Almedahls
FB092, Brillant, 56% Cotton, 44% Nylon, Berger Du Nord
FB093, Gordon Harvest, 56% Cotton, 44% Polyester, Douglas Furniture
FB094, Scarlet, 57% Acrylic, 29% Nylon, 14% Polyester, Austermann
FB095, Donahue Sand, 57% Polyester, 27% Olefin, 11% Acrylic, 5% Cotton, Douglas Furniture
FB096, Look Slo Mocha, 57% Polyester, 43% Olefin, Douglas Furniture
FB097, Mosaic , 58% Polyester, 34% Cotton, 8% Nylon, Trendsetter yarns
FB098, Hadley Loden, 58% Polyester, 37% Acrylic, 16% Olefin, Douglas Furniture
FB099, Chorus Azure, 58% Polyester, 39% Acrylic, 3% Nylon, Douglas Furniture
FB100, Zucca, 58% Tactel, 42% Polyester, Trendsetter yarns
FB101, Perle, 59% Nylon, 39% Acrylic, 2% Polyester, Almedahls
FB102, Vixen, 59% Rayon, 37% Cotton, 4% Polyester, Trendsetter yarns
FB103, Bubble Print , 60% Cotton, 40% Nylon, Austermann
FB104, Magique , 60% Nylon, 16% Cotton, 15% Acrylic, 9% Polyester, Trendsetter yarns
FB105, Montage , 60% Nylon, 35% Rayon, 5% Polyester, Trendsetter yarns
FB106, Odyssey Glitz , 60% Nylon, 37% Wool, 3% Lamé, Plymouth Brand
FB107, Heroes, 60% Nylon, 40% Acrylic, Artful Yarns
FB108, Odyssey, 60% Nylon, 40% Wool , Plymouth Brand
FB109, Accolade Dusk, 60% Olefin, 38% Acrylic, 2% Polyester, Douglas Furniture
FB110, Soulmate Multi, 60% Olefin, 40% Cotton, Douglas Furniture
FB111, Talley Carmel, 60% Polyester, 40% Cotton, Douglas Furniture
FB112, Stardust, 60% Polyester, 40% Nylon, Karabella
FB113, Roses , 60% Rayon, 40% Polyester, Karabella
FB114, Tyrol, 60% Wool, 30% Alpaca, 10% Acrylic, Grignasco
FB115, Bouclerone , 60% Wool, 34% Nylon, 6% Polyester, Austermann
FB116, Saville Garden Multi, 63% Olefin, 37% Polyester, Douglas Furniture
FB117, Puddel, 64% Mohair, 20% Wool, 16% Rayon, Aurora Yarns
FB118, Ischia Silk , 64% Silk, 32% Cotton, 4% Nylon, Bollicine
FB119, Aventura , 65% Cotton, 25% Nylon, 10% Rayon, Trendsetter yarns
FB120, Dune Beige, 65% Cotton, 35% Polyester, Douglas Furniture
FB121, Hot Hot Hot, 65% Nylon, 30% Polyester, 5% Acrylic , Plymouth Brand
FB122, P-zazz, 65% Nylon, 35% Rayon , Plymouth Brand
FB123, Topton Barley, 65% Olefin, 35% Dacron, Douglas Furniture
FB124, Charming , 65% Polyester, 35% Tactel, Trendsetter yarns
FB125, Combolo, 66% Nylon, 30% Tactel , Plymouth Brand
FB126, Persuade Brown, 66% Olefin, 21% Polyester, 13% Acrylic, Douglas Furniture
FB127, Tic Tac Multi, 68% Cotton, 30 % Polyester, 2% Olefin, Douglas Furniture
FB128, Horizon Blue, 69% Cotton, 19% Olefin, 12% Polyester, Douglas Furniture
FB129, Santana Blue, 69% Polyester, 17% Acrylic, 13% Cotton, 1% Nylon, Douglas Furniture
FB130, Camille , 70% Cotton, 13% Rayon, 11% Nylon, 6% Acrylic, Adrienne Vittadini
FB131, Teddy, 70% Courtelle Acrylic, 17% Wool, 13% Polyester, Almedahls
FB132, Martina , 70% Extrafine Merino Wool, 30% Silk, Adrienne Vittadini
FB133, Etoile, 70% Kid Mohair, 30% Acrylic, Bollicine

FB134, Maxi, 70% Merino, 30% Acrylic, Bollicine
FB135, Portrait , 70% Mohair, 25% Rayon, 5% Polyester, Artful Yarns
FB136, Baby Lash , 70% Nylon, 30% Acrylic, Bernat
FB137, La Scala , 70% Nylon, 30% Cotton, Trendsetter yarns
FB138, Elegante, 70% Rayon, 13% Cotton, 13% Acrylic, 4% Misc., Almedahls
FB139, Amsterdam, 73% Nylon, 27% Acrylic, Lana gatto
FB140, Frolic Fawn, 73% Olefin, 17% Polyester, 9% Cotton, Douglas Furniture
FB141, Encore Mega, 75% Acrylic, 25 %Wool, Plymouth Brand
FB142, Party , 75% Polyester, 25% Nylon, Wendy
FB143, Sunshine Print , 75% Rayon, 25% Polyester, Trendsetter yarns
FB144, Gabrielle, 75% Wool, 15% Rayon, 10% Nylon , Plymouth Brand
FB145, La Gran Mohair, 76% Mohair, 18% Wool, 6% Nylon, Classic Elite
FB146, Flora, 76% Rayon, 24% Polyester, Trendsetter yarns
FB147, Wyndham Blue, 78% Acrylic, 14% polyester, 8% Nylon, Douglas Furniture
FB148, Margarita, 78% Nylon, 22% Microtactel , Plymouth Brand
FB149, Baby Alpaca Brush, 80% Baby Alpaca, 20% Acrylic, Alpaca - Indiecita Brand
FB150, Rapunzel , 80% Mohair, 15% Wool, 5% Nylon, Trendsetter yarns
FB151, Trevira, Polyester, Hoechst
FB152, Polar, 100% Wool, Almedahls
FB153, Natura, 30% Alpaca, 30% Wool, 37% Acrylic, 3% Rayon, Almedahls
FB154, Imperiale Super Kid Mohair, 80% Mohair, 20% Nylon, Lane Cervinia Le Fibre Nobili
FB155, Ranee, 80% Nylon, 20% Polyester, Kertzer
FB156, Poco, 80% Rayon, 20% Nylon, Kertzer
FB157, Diamante , 80% Rayon, 20% Polyester, Karabella
FB158, Sinsation, 80% Rayon, 20% Wool, Plymouth Brand
FB159, 24K , 82% Nylon, 18% Lame, Plymouth Brand
FB160, Carnivala , 82% Rayon, 12% Nylon, 6% Polyester, Plymouth Brand
FB161, Bonneville Blue, 84% Herculon, 11% Polyester, 5% Cotton, Douglas Furniture
FB162, Naker Eros, 85% Acrylic, 15% Rayon, Plymouth Brand
FB163, Baker Beige, 85% Cotton, 15% Polyester, Douglas Furniture
FB164, Palazzi , 85% Nylon, 15% Wool, Plymouth Brand
FB165, Crisantemo , 85% Polyester, 15% Rayon, Trendsetter yarns
FB166, Country 8Ply Naturals, 85% Pure Wool, 10% Acrylic, 5% Rayon, Machine Washable, Cleckheaton
FB167, Iris , 86% Nylon, 14% Acrylic, Trendsetter yarns
FB168, Frizzultra , 90% Nylon, 10% Acrylic, Plymouth Brand
FB169, Diva, 95% Nylon, 5% Metallic Polyester, Almedahls
FB170, Acrylic >80%, Acrylic, Promod
FB171, Acrylan, Modacrylic, Monsanto
FB172, Leather - Pig - Velour, frontside, Animal, NA
FB173, Leather - Pig - Velour, underside, Animal, NA
FB174, Leather - Cow, Animal, NA
FB175, Leather - Lamb - Napa, Animal, NA
FB176, Leather - Goat - Velour, Animal, NA
FB177, Leather - Goat - Napa, Animal, NA
FB178, Fur Ocelot, Animal, NA
FB179, Fur Persian Breitschwanz, Animal, NA
FB180, Fur Persian Swakara, Animal, NA

FB181, Fur Persian 1 month, Animal, NA
FB182, Fur Persian 2 months, Animal, NA
FB183, Fur Antelope, Animal, NA
FB184, Fur Persian 3 months, Animal, NA
FB185, Skin Snake, frontside, Animal, NA
FB186, Skin Snake, underside, Animal, NA
FB187, Fur Ostrich, Animal, NA
FB188, Fur Raccoon, Animal, NA
FB189, Fur Red Fox Silver (Blue Colored), Animal, NA
FB190, Fur Red Fox White, Animal, NA
FB191, Fur Badger, Animal, NA
FB192, Fur Guereza Monkey, Animal, NA
FB193, Fur Rabbit, Animal, NA
FB194, Fur Goat frontside, Animal, NA
FB195, Fur Goat backside, Animal, NA
FB196, Fur Nutria frontside, Animal, NA
FB197, Fur Nutria backside, Animal, NA
FB198, Fur Muskrat, Animal, NA
FB199, Fur Mink, Animal, NA
FB200, Fur Marmot, Animal, NA
FB201, Fur Opossum, Animal, NA
FB202, Fur Stoat, Animal, NA
FB203, Fur Coon, Animal, NA
FB204, Fur Skunk frontside, Animal, NA
FB205, Fur Skunk backside, Animal, NA
FB206, Fur Marten, Animal, NA
FB207, Fur Bear, Animal, NA
FB208, Fur Wolf frontside, Animal, NA
FB209, Fur Wolf backside, Animal, NA
FB210, Lenka, Acrylic 15%, Cotton 64%, Flax 21%, Ariadne
FB211, Leacril, Acrylic, Monte Fibre
FB212, Acrylic/Polyamide #1, Mixed Yarn, C&A
FB213, Nylon 2%, Rayon 98%, Mixed Yarn, H&M
FB214, Elian Helen, Acrylic 33%, Nylon 34%, Wool 33%, VSV
FB215, Acrylic 80%, Wool 20%, Acrylic/Wool, Mango
FB216, Chemlon, Nylon , Chemko
FB217, Cotton/Rayon/Acrylic, Mixed Yarn, Zara
FB218, Acrylic 90%, Mohair 10%, Acrylic/Wool, Mango
FB219, Elaspán, Spandex, Invista
FB220, Glospan, Spandex, Globe
FB221, Hoechst, Polyester, Hoechst
FB222, Nylon > 85% #2, Nylon , Benetton
FB223, Nylon > 80% #3, Nylon, Tesco
FB224, Nylon > 80% #1, Nylon, Beltissimo
FB225, Acrylic 65%, Metalic 35%, Mixed Yarn, Mango
FB226, Rubber fiber #2, Mixed Yarn, Svadlena
FB227, Polarguard 3D, Polyester, KoSa
FB228, Elian Sport, Acrylic 54%, Cotton 16%, Rayon 4%, Wool 26%, VSV

FB229, Xylon, Olefin PP, Alladin
FB230, Innova, Olefin PP, American Fibers
FB231, Isover, Glass, Saint-Gobain
FB232, Rotaflex, Glass, Union
FB233, Wadding, Natural, Tesco
FB234, Cotton 80%, Natural, C&A
FB235, Nylon/Rayon, Mixed Yarn, Max Mara
FB236, Kevlar 49, Aramid , Dupont
FB237, Kevlar, Aramid , Dupont
FB238, Nomex #1, Aramid , Dupont
FB239, Nomex #2, Aramid , Dupont
FB240, Straw treated, Natural, NA
FB241, Sensura, Polyester, Wellman
FB242, PVC based fiber #2, PVC, Pepe
FB243, PVC based fiber #1, PVC, Pepe
FB244, Rayon 50%, Mixed Yarn, Max Mara
FB245, Cotton 90%, Natural, C&A
FB246, Amoco Cf Yarn, Olefin PP, Amoco
FB247, Synera, Olefin PP, Amoco
FB248, Vicryl (surgery), Biodegradable, PLLA, Chirmac
FB249, Dexon (surgery), Biodegradable, PLLA, Dexon
FB250, Spider Net (Tegenaria domestica), Animal, NA
FB251, Bamboo 80%, Cotton 20%, Mixed Yarn, Bambus Elegang
FB252, Bamboo fiber, Natural, Bambus Elegang
FB253, Fantasia, Acrylic 70%, Cotton 23%, Effect 7%, SDS
FB254, Tiffany forte, Acrylic 40%, Cotton 60%, SDS
FB255, Cotton 95%, Lycra 5%, Mixed Yarn, H&M
FB256, Rayon 92%, Elastan 8%, Mixed Yarn, Beltissimo
FB257, Polyester 86%, Elastan 14%, Mixed Yarn, Benetton
FB258, Nylon 92%, Elastan 8%, Mixed Yarn, Benetton
FB259, Acrylic/Polyamide #2, Mixed Yarn, Zara
FB260, Modal 92%, Spandex 8%, Mixed Yarn, H&M
FB261, Rayon 95%, Elastan 5%, Mixed Yarn, Tesco
FB262, Rayon 85%, Flax 15%, Mixed Yarn, Beltissimo
FB263, Ramie 42%, Flax 58%, Mixed Yarn, Beltissimo
FB264, Nylon 83%, Elastan 17%, Mixed Yarn, Benetton
FB265, Flax, Natural, Ethno
FB266, Nylon 30%, Silk 70%, Mixed Yarn, Benetton
FB267, Cotton 33%, Rayon 49%, Nylon 16%, Elastan 2%, Cotton Nylon, C&A
FB268, Hair of Horse #1, Animal, NA
FB269, Hair of Horse #2, Animal, NA
FB270, Feather Hen, Animal, NA
FB271, Feather Pheasant, Animal, NA
FB272, Hair of Pig, Animal, NA
FB273, Hair of Dog, Animal, NA
FB274, Hair of Rabbit, Animal, NA
FB275, Hair of Cat, Animal, NA
FB276, Silk India Assam Muga, Silk, Ethno

FB277, Straw, Natural, NA
FB278, Ramie treated, Natural, NA
FB279, Acrylic 50%, Polyester 37%, Wool 13%, Mixed Yarn, Mango
FB280, Hot Socks Cotton Stretch, Cotton 41%, Nylon 13%, Wool 39%, Grundl
FB281, Rayon 95%, Elastan 5%, Mixed Yarn, Tesco
FB282, Wool - Lamb, Animal, NA
FB283, Hair Human Female Colored 2, Human, NA
FB284, Hair Human Female Colored 1, Human, NA
FB285, Nylon Rion Taslan PU milky, Mixed Yarn, Promod
FB286, Fibro, Rayon, Courtaulds
FB287, Biofresh, Acrylic, Sterling Fibers
FB288, Crimplene, Polyester, Regal
FB289, Chromspun, Acetate, Voridian
FB290, Piuma Piccolo, Polyester 40%, Rayon 60%, Ariadne
FB291, Silk, Animal, Ethno
FB292, Rayon 80%, Rayon, H&M
FB293, Acrylic 70%, Wool 30%, Mixed Yarn, Clockhouse
FB294, Silkworm Cocoon, Animal, Africa
FB295, Amica, Acetate 5%, Acrylic 60%, Cotton 25%, Nylon 10%, Grundl
FB296, Beluga, Acrylic 30%, Nylon 30%, Polyester 20%, Wool 20%, Grundl
FB297, Sarah, Acrylic 62%, Cotton 36%, Nylon 2%, Toptex
FB298, Atlantis, Acrylic 41%, polyester 41%, Nylon 18%, SDS
FB299, Nefertiti, Acrylic 18%, Nylon 7%, Rayon 75%, SDS
FB300, Elian Betty, Acrylic 43%, Nylon 29%, Wool 28%, VSV
FB301, Evita, Acrylic 21%, Cotton 79%, Vlnap
FB302, Frote, Acrylic 60%, Nylon 40%, Ariadne
FB303, Elian Ester, Acrylic 65%, Polyester 34%, VSV
FB304, Anna, Acrylic 50%, Rayon 50%, Toptex
FB305, Sunshine, Cotton 35%, Rayon 65%, Grundl
FB306, Edita Glamour, Cotton 55%, Rayon 45%, Toptex
FB307, Laguna, Cotton 67%, Silk 33%, Ariadne
FB308, Atol, Cotton 70%, Wool 30%, Ariadne
FB309, Hot Socks Uni , Nylon 25%, Wool 75%, Grundl
FB310, Vinyon, PVC, Rhovyl
FB311, Gore-Tex, Fluorocarbon, Gore
FB312, Kuraray PVA, PVA, Kuraray
FB313, Dyneema, UHMW PE, DSM
FB314, Kuralon, PVA, Kuraray
FB315, Polyimide Fiber, Polyimide Fiber, Germany
FB316, Basofil, Melamine Formaldehyde, BASF
FB317, Rhovyl, PVC, Rhone-Poulenc
FB318, Ardil, Azlon, ICI
FB319, Fibrolane , Azlon, Courtaulds
FB320, Vicara, Azlon, Virginian
FB321, Tencel, Lyocell, Tencel
FB322, Nylon/Cotton/Acrylic, Mixed Yarn, Zara
FB323, Wellene, Polyester, Wellman
FB324, Nylon > 90%, Cotton, Nylon/Cotton, Promod

FB325, Silk Spider, Silk, Ethno
FB326, Rubber fiber, Rubber, Tesco
FB327, Wool Sheep Iceland, Wool, Animal
FB328, Estron, Acetate, Eastman Chemical
FB329, Nylon > 90% #2, Nylon, C&A
FB330, Microsafe, Acetate, Celanese
FB331, Wool 80%, Acrylic, Mixed Yarn, Promod
FB332, Regular, Acetate, Celanese
FB333, Wolpryla, Acrylic, Markische Faser
FB334, Wool Sheep Corriedale, Wool, Animal
FB335, Hwm Rayon, Rayon, Avtex
FB336, Angora, Wool, Wonoco
FB337, Pompadour Zephyr, Wool, Wonoco
FB338, Wool Sheep Cheviot, Wool, Animal
FB339, Wool English #1, Wool, Animal
FB340, Wool Ilama, Wool, Animal
FB341, Scotch Tweed, Wool, Wonoco
FB342, Wool Sheep Churro, Wool, Animal
FB343, Knitting Worsted, Wool, Wonoco
FB344, Zephyr Paradise, Wool, Wonoco
FB345, Zephyr Paradise green, Wool, Wonoco
FB346, Heather Tone, Wool, Wonoco
FB347, Silken Wool, Mixed Yarn, Iran
FB348, Camel Hair, Wool, Wonoco
FB349, Crepe Welnette, Mixed Yarn, Wonoco
FB350, Pebble Crepe, Natural, Wonoco
FB351, Suedetex, Cotton , Wonoco
FB352, Snow Flake, Cotton , Wonoco
FB353, Twistex, Cotton/Rayon, Wonoco
FB354, Mystic Tweed Crepe, Mixed Yarn, Wonoco
FB355, Bouncle, Mixed Yarn, Wonoco
FB356, Heathertone Sport, Wool, Wonoco
FB357, Saxony, Wool, Wonoco
FB358, German Town, Wool, Wonoco
FB359, Shetland Floss, Wool, Wonoco
FB360, Feather-Fleece, Wool, Wonoco
FB361, 2 Ply Silver Gloss, Wool, Wonoco
FB362, Sport Yarn, Wool, Wonoco
FB363, Sporting Angora, Wool, Wonoco
FB364, Speckle Crepe, Wool, Wonoco
FB365, Wool English #2, Wool, Animal
FB366, Cable Yarn, Wool, Wonoco
FB367, Silken Shetland, Mixed Yarn, Wonoco
FB368, Zephyr Tweed De Luxe, Wool, Wonoco
FB369, Wool Alpaca, Wool, Animal
FB370, Wool Sheep Border Leicester, Wool, Animal
FB371, Country Club, Wool, Wonoco
FB372, French Zephyr, Wool, Wonoco

FB373, Wool Flake, Wool, Wonoco
FB374, Bermuda Tweed, Mixed Yarn, Wonoco
FB375, Paisley Tweed, Mixed Yarn, Wonoco
FB376, Silk India Kanchi, Silk, Ethno
FB377, Silk Mud-Mee Thailand, Silk, Ethno
FB378, Dexon Plus (surgery), Biodegradable, PLLA, Dexon
FB379, Hemp, Natural, Ethno
FB380, Hemp Fine, Natural, Ethno
FB381, Agave, Natural, NA
FB382, Kapok, Natural, NA
FB383, Musa, Natural, NA
FB384, Pina, Natural, NA
FB385, Raffia, Natural, NA
FB386, Ramie, Natural, NA
FB387, Kapok treated, Natural, NA
FB388, Musa treated, Natural, NA
FB389, Palmtree fiber (treated), Natural, Spring
FB390, Bamboo - piece (treated), Natural, Spring
FB391, Jute, Natural, Spring
FB392, Avora FR, Polyester, KoSa
FB393, Rayon >80%, Nylon, Polyester, Mixed Yarn, Beltissimo
FB394, Rayon > 70%, Acetate Rayon, Retro
FB395, Quattro, Acrylic 25%, Nylon 15%, Wool 60%, Vlnap
FB396, Absorbit, Rayon, Enka
FB397, Comforel, Polyester, Invista
FB398, Dunova, Acrylic, Bayer
FB399, Arnel, Acetate, Celanese
FB400, Avril, Rayon, Avtex
FB401, Durvil Rayon, Rayon, Avtex
FB402, Beau Grip, Rayon, North American Rayon
FB403, Rustico, Acrylic 48%, Cotton 48%, Nylon 4%, SDS
FB404, Cotton , Cotton , China
FB405, Cotton , Cotton , Russia
FB406, Acrylic 90%, Nylon 10%, Acrylic/nylon, C&A
FB407, Nylon > 90% #3, Nylon, H&M
FB408, Cotton , Cotton , Korea
FB409, Nylon > 90% #1, Nylon, C&A
FB410, Nylon/Wool/Acrylic, Mixed Yarn, Promod
FB411, Flocking Tow, Rayon, North American Rayon
FB412, SXD Rayon, Rayon, Avtex
FB413, T-45, Rayon, North American Rayon
FB414, Cotton , Cotton , USA
FB415, Wool/Acrylic, Acrylic/Wool, C&A
FB416, Cotton , Cotton , India
FB417, Anso, Nylon 6, Honeywell
FB418, Antron, Nylon 6,6, Invista
FB419, Cordura, Nylon 6,6, Invista
FB420, Caprolan, Nylon 6, Honeywell

FB421, Seraceta, Acetate, Seraceta
FB422, Meryl, Nylon 6,6, Nylstar
FB423, Dyenamix, Nylon 6,6, Solutia
FB424, Shimmereen, Nylon 6, Berkley
FB425, Type 3000, Rayon, North American Rayon
FB426, Regular Rayon, Rayon, Avtex
FB427, Enkalon, Nylon 6, Enka
FB428, Powersilk, Nylon 6,6, BASF
FB429, Coloray, Rayon, Courtaulds
FB430, Samanta, Acrylic 20%, Cotton 60%, Rayon 20%, Vlnap
FB431, Lanese, Acetate Polyester, Celanese
FB432, Avtex, Acetate, Avtex
FB433, Acetate > 70%, Acetate Polyester, Retro
FB434, Duraspun, Acrylic, Solutia
FB435, Fortrel Microspun, Polyester, Wellman
FB436, Qiana, Nylon 6, Dupont
FB437, Nylon > 85% #1, Nylon, Benetton
FB438, Cotton , Cotton , Uzbekistan
FB439, Fillwell, Polyester, Wellman
FB440, Microdenier Sensura, Polyester, Wellman
FB441, Nylon > 80% #2, Nylon, Tesco
FB442, Steripur, Polyester, DAK Americas
FB443, Avlin, Polyester, Avtex
FB444, Synthetic Leather #4, Poly(ester urethane) based, Jena Furniture
FB445, Microlux, Polyester, KoSa
FB446, YKK Polyester, Polyester, YKK
FB447, Colorguard, Polyester, Martin Color
FB448, Synthetic Leather #1, Poly(ester urethane) based, Jena Furniture
FB449, Synthetic Leather #3, Poly(ester urethane) based, Jena Furniture
FB450, Dacron, Polyester, Dupont
FB451, Diolen, Polyester, Diolen Industrial Fibers
FB452, Synthetic Leather #2, Poly(ester urethane) based, Jena Furniture
FB453, Enkron, Polyester, Enka
FB454, Polysax-X, Olefin PP, Chase Bag
FB455, Vectra, Olefin PP, Chevron
FB456, Pina treated, Natural, NA
FB457, Agave treated, Natural, NA