



NICODOM NIR Pharmaceutical - List of spectra

The near infrared spectral library "NICODOM NIR Pharmaceuticals" contains 696 NIR spectra (4.200- 11.000 cm⁻¹) of active substances and excipients used in pharmaceutical industry. The source of the samples were local pharmaceutical companies.

The library is available as searchable digital library compatible with your spectroscopic software.

The sample information includes compound name (INN or commonly used), type (excipient, drug...), purpose (emulsifier, laxative...) characteristics (color, crystallinity...).

The spectra were collected by ThermoNicolet FTNIR spectrometer Antaris™, the powders were collected in glass vials on Integrating sphere module, transparent liquids were collected in transmittance or with a transreflectance mirror.

Available formats: searchable digital library compatible with your spectroscopic software or digital book in *.pdf format.

This library is available as digital searchable library compatible with your spectroscopic software. Also it can be delivered with SEARCH software or in *.pdf format.

To download the list of spectra from this library, to download a free demo library compatible with your software, to check the format compatibility as well as for pricelist, ordering info and information about all NICODOM IR/NIR Libraries and other NICODOM products please visit our webpage.

NICODOM NIR Pharmaceuticals, 696 NIR spectra of pharmaceuticals
Copyright © NICODOM 2005

NICODOM Ltd., Hlavni 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

ftir@ftir.cz

List of Spectra

B001 Vanilla odour
B002 Acesulfam K
B003 Acetazolamide
B004 Acetylcysteine
B005 Acetylsalicylic Acid
B006 Acetyltyrosine
B007 Linoleic Acid
B008 Acriflavine
B009 Hard Fat
B010 Adipic Acid
B011 Aesculine
B012 Agar
B013 Agar with maltose extract
B014 Alanine
B015 D-alpha-Alanine
B016 Opaspray
B017 Alcohol Cetylstearyl
B018 Allantoin (5- Ureidohydantoin)
B019 Allura red AC
B020 Alps Herbs Flower Aroma 740025 H
B021 Alprazolam
B022 Aluminium Chlorohydrate
B023 Aluminium Monostearate
B024 Aluminium Oxide 90
B025 Amaranthe
B026 Ambroxolic Chloride
B027 Aminoacetic Acid
B028 DL-2-Aminoadipic acid
B029 DL-2-Aminobutyric acid
B030 2-Aminoisobutyric Acid
B031 Aminophenazone
B032 Aminophylline
B033 Amitriptyline Chloride
B034 Ammonium Acetate
B035 Ammonium Bromide
B036 Ammonium Chloride
B037 Amylopectine
B038 Anise oil
B039 Antazoline mesylate
B040 Silver Albumin Acetyl Tannate
B041 Silver Albumin Acetyl Tannate with Borax
B042 Ascorbic Acid (Vitamin C)
B043 Ascorbyl palmitate
B044 Aspartate
B045 Aspartic Acid
B046 D-Aspartic acid

NICODOM NIR Pharmaceuticals, 696 NIR spectra of pharmaceuticals
Copyright © NICODOM 2005

NICODOM Ltd., Hlavni 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

ftir@ftir.cz

B047 Atenolol
B048 Atropine Sulphate
B049 Orange Oil
B050 Avocado Oil 3800010
B051 D-Azetidine-2-carboxylic acid
B052 Azidamfenicol
B053 Azorubine E 122
B054 Bacitracine Zinc
B055 Barbitol
B056 Barbitol Sodium
B057 Bees wax
B058 Bentonite
B059 Bentonite 750
B060 Benzoic Acid
B061 Bergamotte Oil
B062 Betacare
B063 BHT (butylated hydroxytoluene)
B064 Bilan
B065 Bisabolol
B066 Bismuth Subnitrate
B067 Bismuth Tribromophenol
B068 Bismuth Tribromophenol Bas.
B069 Bismuth Gallate Basic
B070 Bolus Alba (Kaolin)
B071 Boric Acid
B072 Bovexan-P (colagen)
B073 Brij 35
B074 Brij 58 P
B075 Brilliant blue FCF
B076 Bromisovale
B077 Bronopol
B078 Butamirate citrate
B079 Butylated hydroxyanisole
B080 Calcium Dobesilate
B081 Calcium Chloride
B082 Calcium Chloride Crystalline
B083 Calcium Lactate
B084 Calcium Acetate
B085 Calcium Carbonate
B086 Calcium Formate
B087 Calcium Gluconate
B088 Calcium Hydrogenphosphate
B089 Calcium Chloride Dihydrate
B090 Calcium Pantothenate (Vitamin B2 complex)
B091 Calcium Phosphate
B092 Calcium Stearate
B093 Camphor
B094 Candellila wax No.7455
B095 Capsula Gelatinous

NICODOM NIR Pharmaceuticals, 696 NIR spectra of pharmaceuticals
Copyright © NICODOM 2005

NICODOM Ltd., Hlavni 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

ftir@ftir.cz

B096 Caramel
B097 Caramel 15461
B098 Caraway oil
B099 Carbaethopendeciniate Bromide
B100 Carbomerum 940
B101 Carbopol 934
B102 Carboxycellulose Calcium (Traumacel)
B103 Carboxymethylamylopectine Sodium
B104 S-Carboxymethyl-L-cysteine
B105 Carmellose (Carboxymethylcellulose)
B106 Carmellose Sodium (Carboxymethylcellulose Sodium)
B107 Caryophyllic Oil
B108 Casein hydrolysate
B109 Cellacefate
B110 Cellactose
B111 Cellulose (Perlosa SF)
B112 Cellulose Microcrystallic (Avicel)
B113 White Bees Wax
B114 Wool Wax
B115 Cerasynt IP
B116 Ceresin white 6062
B117 Cetaceum
B118 Cetiol 1414 E
B119 Cetiol MM
B120 Cetiol OE
B121 Cetiol V
B122 Cetomacrogol 1200
B123 Cetyl Alcohol
B124 Cetylpyridinium Chloride
B125 2-Chloro-4-nitrophenol
B126 Cinchocaine Chloride
B127 Cinnamon Oil
B128 Ciprofloxacin
B129 Citric Oil
B130 Citric Acid
B131 Citric Acid Monohydrate
B132 Citronella oil
B133 Clopamide
B134 Clove oil
B135 Cocoa Butter
B136 Cocoamidopropyl Hydroxysultaine
B137 Codeine Phosphate Hemihydrate
B138 Caffeine
B139 Caffeine Citrate
B140 Caffeine in Sodium Benzoate
B141 Collodium
B142 Composition Maryna 70
B143 Composition Sonia 135/28
B144 Coriander Oil

NICODOM NIR Pharmaceuticals, 696 NIR spectra of pharmaceuticals
Copyright © NICODOM 2005

NICODOM Ltd., Hlavní 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

ftir@ftir.cz

B145 Cosmowax D
B146 Crodolan
B147 Crolactil SSL
B148 Croscarmellose sodium
B149 Crosilk powder
B150 Curacao flavour
B151 Cyanocobalamine (Vitamin B12)
B152 Cyanuric Acid
B153 DL-Cystathionine
B154 Cysteine
B155 DC Antifoam AF Emulsion
B156 DC Fluid 200/5000
B157 Decanol - Alfol
B158 Dexamethasone Acetate
B159 Dexpanthenol
B160 Dextrane 41
B161 Dextrane 70
B162 Diaveridine
B163 Dibutyl Phthalate
B164 Diclofenate Sodium
B165 Digoxine
B166 Dihydroxypropyladenin (DHPA)
B167 3;5-Diiodo- L-tyrosine
B168 Diltiazem Hydrochloride
B169 Dimeticone
B170 Dithiadene
B171 Dithiadene Chloride
B172 Dithranole
B173 DL-Ornithine hydrochloride
B174 DL-Phenylserine
B175 DL-DOPA = 3-(3;4-dihydroxyphenyl)- DL-alanine
B176 Dopamine
B177 Dosulepine
B178 Dow Corning Fluid 200/200 CS
B179 Dow Corning Fluid 345
B180 Doxycycline
B181 Edeta BD
B182 Elastin 30P-soluble
B183 Elfan AT 84 G
B184 Emetine dihydrochloride heptahydrate
B185 Emulgator C
B186 Emulgator D
B187 Emulsifier E 2149
B188 Enocyanine E163
B189 Ephedrine Chloride
B190 Epinephrine Ditartrate
B191 epsilon-Aminocaproic acid = 6-Aminohexanoic acid
B192 Ergotamine Tartrate
B193 Erythrosine

NICODOM NIR Pharmaceuticals, 696 NIR spectra of pharmaceuticals
Copyright © NICODOM 2005

NICODOM Ltd., Hlavní 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

ftir@ftir.cz

B194 Esculine Sesquihydrate
B195 Flavour citric and limonade C.K.51705
B196 Ethacridine Lactate
B197 Ethoxylated Lanoline 4000
B198 Ethyl Biscoumacetate
B199 Etofilline
B200 Eucalyptic Oil
B201 Eudragit L
B202 Eusolex 9020-Butyl methoxydibenzoylmethane
B203 Explotab
B204 Extract of Agrimonia eupatoria L. (10% in Lactose)
B205 Extract of Archangelica officinalis Hoffm. (10% in Lactose)
B206 Extract of Capsella bursa-pastoris L. (10% in Lactose)
B207 Extract of Crataegus oxyacantha L. (10% in Lactose)
B208 Extract of Eleuterococcus senticosus (10% in Lactose)
B209 Extract of Filipendula ulmaria (L.) Max. (10% in Lactose)
B210 Extract of Humulus lupulus L. (10% in Lactose)
B211 Extract of Lavandula officinalis L. (10% in Lactose)
B212 Extract of Melissa officinalis L. (10% in Lactose)
B213 Extract of Mentha piperita L. (10% in Lactose)
B214 Extract of Sambucus nigra L. (10% in Lactose)
B215 Extract of Urtica dioica L. (10% in Lactose)
B216 Extrapone camille 2/367290 Spec
B217 Famotidine
B218 Fenofibrate
B219 Fenpiverine bromide
B220 Iron(II) fumarate
B221 Floralix 1801014
B222 Fluid 3225 C Formulation AID
B223 Flunitrazepam
B224 Fluorescein
B225 Oil of Fennel
B226 Folic Acid (Vitamin B2 complex)
B227 Fuchsine
B228 Fumaric acid
B229 gamma-Aminobutyric acid
B230 Essential Oil Gearnii
B231 Gelatine
B232 Gelatine (180)
B233 Gelatine (Animal)
B234 Gelatine 150 blom, varietas B
B235 Germall 115
B236 Glucose
B237 Glucose Anhydrate
B238 Glucose Monohydrate
B239 Glutamic Acid
B240 Glycerin Ph. Eur. 86.5%
B241 Glycerol 85% W/W
B242 Glycerol 85% W/W /BOEH/

NICODOM NIR Pharmaceuticals, 696 NIR spectra of pharmaceuticals
Copyright © NICODOM 2005

NICODOM Ltd., Hlavni 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

ftir@ftir.cz

B243 Glycerol Monostearate
B244 Glycerol Stearate
B245 Glycine
B246 Glycylglycylglycine = Triglycine
B247 Glycyl-L-Tyrosine Dihydrate
B248 Grapefruit super base (KONTESA)
B249 Guaifenesine
B250 Haarkorbase 2384663
B251 Hamamelis Extract
B252 Hamamelis Extrapon
B253 Hartolan
B254 Sunflower oil
B255 Highly dispersed silicon dioxide
B256 Homatropine Bromide
B257 DL-meso-Homocystine
B258 Horse chestnut extract
B259 Hostacerin WO
B260 Mercury Aminochloride
B261 Hydrocortisone Acetate
B262 Hydrogenoricinomacrogol 1250
B263 Hydromyristenol-Neu 2/014260
B264 Hydroxyamphetamine hydrobromide
B265 Hydroxyethylcellulose
B266 Hydroxymethylcellulose
B267 Hymecromone
B268 Hymono 9004
B269 Hypromellose
B270 Hypromellose (4000)
B271 Chenodeoxycholic Acid
B272 Cherry brandy flavour
B273 Cherry flavour
B274 Quinidine Chloride
B275 Quinidine Sulphate
B276 Quinidine
B277 Chloral Hydrate
B278 Chloramphenicol
B279 Chloramphenicol Palmitate
B280 Chlorhexidine Digluconate, 20% solution
B281 Chlorhexidine Chloride
B282 Chloroxine
B283 Chlorpromazine Chloride
B284 Chlorprothixene Chloride
B285 Chlortalidone
B286 Chocolate brown 4004
B287 Cholesterol
B288 Cholesterol USP
B289 Choline Tartrate
B290 Ibuprofene
B291 Iminodiacetic acid = Diglycine

NICODOM NIR Pharmaceuticals, 696 NIR spectra of pharmaceuticals
Copyright © NICODOM 2005

NICODOM Ltd., Hlavni 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

ftir@ftir.cz

B292 Indapamide Hemihydrate
B293 Indigocarmine
B294 Indomethacine
B295 Inositol
B296 Iodoform
B297 Ioxitalamic Acid
B298 Extract of Ipecac
B299 Iriodin 100 Silver Pearl
B300 Iron Oxide Yellow
B301 Isosorbide 5-Mononitrate
B302 Isostearyl isostearate
B303 Cod-Liver Oil
B304 Yoyoba Oil
B305 Potassium Sorbate
B306 Kebuzone
B307 Kofola Flavour
B308 Lactic Acid
B309 Lactil
B310 Lactose Monohydrate
B311 L-Alanine
B312 Lanette 16 (Lorol C16)
B313 Lanoline Corona LA1439
B314 L-Arginine
B315 L-Arginine Basic
B316 L-Arginine Chloride
B317 Lavender Oil
B318 L-beta-Phenylalanine
B319 L-Cysteine Hydrochloride Hydrate
B320 L-Cystine
B321 L-Djenkolic Acid = L-s,s'-Methylenebis(cysteine)
B322 L-DOPA = 3-(3;4-Dihydroxyphenyl)- L-alanine
B323 Lemon Super Base (KONTESA)
B324 Leucine
B325 Levomenthol
B326 L-Glutamic acid hydrochloride
B327 L-Histidine
B328 5-Hydroxy-DL-lysine hydrochloride
B329 L-Hydroxyproline
B330 Light Liquid Paraffin
B331 Linseed Oil
B332 Liquid Paraffin
B333 L-Isoleucine
B334 L-Leucine
B335 L-Lysine
B336 L-meso-Lanthionine
B337 L-Methionine
B338 Loppa radix extract 2/033041
B339 L-Ornithine Chloride
B340 Betain K 30

NICODOM NIR Pharmaceuticals, 696 NIR spectra of pharmaceuticals
Copyright © NICODOM 2005

NICODOM Ltd., Hlavni 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

ftir@ftir.cz

B341 L-Phenylalanine
B342 L-Proline
B343 L-Serine
B344 L-Threonine
B345 L-Tryptophan
B346 L-Tyrosine
B347 Ludipress
B348 L-Valine
B349 Lycatab PGS
B350 Lysine Acetate
B351 Lysine Chloride
B352 Macrogol Glycerol Trioleate 1100
B353 Macrogol Sodium Sulphosuccinate 30%
B354 Macrogol 1500
B355 Macrogol 300
B356 Macrogol 4000
B357 Macrogol 6000
B358 Magnesium Citrate
B359 Magnesium Lactate
B360 Magnesium Stearate
B361 Magnesium Carbonate
B362 Magnesium Formate
B363 Magnesium Chloride
B364 Magnesium Chloride Hexahydrate
B365 Magnesium Lactate Dihydrate
B366 Magnesium Oxide
B367 Magnesium Sulphate
B368 Maize oil
B369 Maize starch extract
B370 Malic Acid
B371 Maltitol Syrup
B372 Maltodextrin KMS X-50
B373 Maltose
B374 Manganese Sulphate
B375 Mannitol
B376 Marine fish oil
B377 Corn Starch Pregelatinised
B378 Meglumine
B379 Oil of Peppermint
B380 Menthol
B381 Mephenoxalone
B382 Mercaptamine Chloride
B383 Merkur 674 Vaseline Yellow
B384 Mesalazine
B385 Metamizole sodium monohydrate
B386 Methacrylic acid copolymer type NE
B387 Methenamine
B388 Methionine
B389 Methylcellulose (450)

NICODOM NIR Pharmaceuticals, 696 NIR spectra of pharmaceuticals
Copyright © NICODOM 2005

NICODOM Ltd., Hlavní 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

ftir@ftir.cz

B390 Methylcellulose (Methocel E15)
B391 Methylcellulose (Methocel E3)
B392 Methylcellulose (Methocel K100)
B393 Methylcellulose (Methocel K15)
B394 Methylcellulose (Methocel K4)
B395 Methylcellulose(Methocel E5)
B396 Methylparaben
B397 Methylrosaniline Chloride
B398 Methylthionine Chloride
B399 Metronidazole
B400 Microcrystalline Cellulose
B401 Montane 80
B402 Montane 83
B403 Montanov 68
B404 Multifruit Flavour
B405 MYRJ 52 S
B406 MYRJ 53 P
B407 MYRJ 59
B408 Myrrh
B409 Tincture of Myrrh
B410 N-Acetyl-3;5-diiodo-L-tyrosine methyl ester
B411 Nacrelon S-30
B412 Naftidrofuryli Hydrogenooxalate
B413 Sodium Caseine
B414 Sodium Edetas
B415 Sodium ethyl parahydroxybenzoate
B416 Neo Dragocid Liquid
B417 Neomycine Sulphate
B418 Neoserine Methylsulfonate
B419 Nettle Extract
B420 Nicotinamide (Vitamin B2 complex)
B421 Nicotinic Acid
B422 Nipagin M
B423 Nipasol M
B424 Nitrofurantoin
B425 Nonivamide
B426 N-Phenylglycine
B427 Ofloxacin
B428 Oleoyl Macrogolglycerides
B429 Olive Oil
B430 Olive Oil for injections
B431 Opadry OY-B 22932 yellow
B432 Opadry OY-S-8720 orange
B433 Orange Flavour
B434 Orange Oil
B435 Orange Super Base (KONTESA)
B436 Ornidazole
B437 Rice Starch
B438 Oxiconazole Nitrate

NICODOM NIR Pharmaceuticals, 696 NIR spectra of pharmaceuticals
Copyright © NICODOM 2005

NICODOM Ltd., Hlavní 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

ftir@ftir.cz

B439 D-Panthenol USP
B440 Papaverine Chloride
B441 Paracetamol
B442 Paracetamol (DC-90)
B443 Paraffin
B444 Paraffin Solid
B445 Pareth-15-6
B446 Parsol 1789
B447 Parsol MCX
B448 Pecan-Nut Oil
B449 Pectine (Citrus 106)
B450 Perfume 102 556
B451 Perfume 202 695
B452 Perfume 37358
B453 Perfume AP 0662
B454 Perfume Apple 393
B455 Perfume Baby UN 166.326
B456 Perfume Brilliance UN150444
B457 Perfume BS 4186
B458 Perfume BS 4281
B459 Perfume Capricorn 027
B460 Perfume Cuddle 455
B461 Perfume Floral FEMME SMN
B462 Perfume FRAG 44451
B463 Perfume Greenol 507826
B464 Perfume Kiwi UN156.232
B465 Perfume L87-6452
B466 Perfume Laguna B30 550
B467 Perfume Lancia UN 81234
B468 Perfume LS3Q
B469 Perfume LS7
B470 Perfume Mereta UN156.179
B471 Perfume Phoenix 654
B472 Perfume PPL SC 2367
B473 Perfume PPL SC 2367
B474 Perfume Provide UN 437.134B
B475 Perfume Sangria UN134.886
B476 Perfume Serenity P/BR SM
B477 Perfume Shangrila UN134.638
B478 Perfume STELLA 238/01
B479 Perfume SUB 1881 NMF
B480 Perfume Toucana 830
B481 Perfume TS-3112B
B482 Perfume URSA 009
B483 Perfume WR07724 Baby Blue
B484 Phenacetine
B485 Phenazone
B486 Phenazone with Caffeine and Acid Citric
B487 Phenobarbital

NICODOM NIR Pharmaceuticals, 696 NIR spectra of pharmaceuticals
Copyright © NICODOM 2005

NICODOM Ltd., Hlavni 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

ftir@ftir.cz

B488 Phenobarbital Sodium
B489 Phenol
B490 Phenyl dimethicon modif.
B491 Phenylephrine hydrochloride
B492 Phenyl Mercury Borate
B493 Phenyl Salicylate
B494 Phthalylsulfathiazol
B495 Pilocarpine Chloride
B496 Oil of Pine
B497 Piroxicame
B498 Pitofenone hydrobromide
B499 Pitofenone Chloride
B500 Plantacare 2000 UP
B501 Resin of Podophyllum
B502 Polawax GP 200 (non-animal)
B503 Poloxamer 105
B504 Polyethylene 617A
B505 Polymer JR-400 UCARE
B506 Polynolum A
B507 Polysorbate
B508 Polysorbate 20
B509 Polysorbate 80
B510 Povidon 25 (Kollidon 25)
B511 Polyvidone 30 (Kollidonum 30)
B512 Povidon 90 (Kollidon 90)
B513 Povidon Insoluble (Kollidon CL)
B514 Ponceau 4R
B515 Portuguese Gum Turpentine
B516 Potassium Acetate
B517 Potassium Carbonate
B518 Potassium Citrate
B519 Potassium Citrate Monohydrate
B520 Potassium Dihydrogenfosfate
B521 Potassium Nitrate
B522 Prazosine Chloride
B523 Procaine Chloride
B524 Propafenone
B525 1;2-Propylene Glycol
B526 Propylene Glycol
B527 Propyl Gallate
B528 Propylparaben
B529 Propyphenazone
B530 Protachem IPS
B531 Pyridoxine Chloride (Vitamin B6)
B532 Quinine Chloride
B533 Quinine Sulfate
B534 Quinoline Yellow, E 104, CI 47005
B535 Ranitidine Chloride
B536 Raspberry Flavour 716K-5

NICODOM NIR Pharmaceuticals, 696 NIR spectra of pharmaceuticals
Copyright © NICODOM 2005

NICODOM Ltd., Hlavni 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

ftir@ftir.cz

B537 Raspberry Flavour D 9599
B538 Raspberry Flavour DU 04125
B539 Refined Maize Oil
B540 Reserpine
B541 Resorcinol
B542 Revitalin - BT
B543 Rewocid UTM 185
B544 Rewoderm ES 90
B545 Rewoderm LI 420-70
B546 Rewominox B 204
B547 Riboflavin (Vitamin B2 complex)
B548 Castor Oil
B549 Ricinomacrogol 1750
B550 Rose Perfume 0065071
B551 Rosmarine Oil
B552 Roxitromycine
B553 Rubor 10 B33, CI 17200
B554 Rubor 3 C.I.14720
B555 Rubor Congoensis
B556 Rutin Sodium Sulphate
B557 Rutoside
B558 Saccharin Sodium (Diaspon)
B559 Saccharin Sodium monohydrate
B560 Lactose
B561 Lactose 200M
B562 Lactose 100
B563 Lactose 200
B564 Salicylic Acid
B565 Sassafras Oil
B566 Sepicide LD
B567 Silica Gel, 70-230 mesh
B568 Siliceous Earth
B569 Siloxide
B570 Simethicone
B571 Slovapol N 6820 (ethoxylated alcohols)
B572 Sodium Salicylate
B573 Sodium Acetate Trihydrate
B574 Sodium Alkylpolyglycolethersulphate
B575 Sodium Benzoate
B576 Sodium Calcium Edetate
B577 Sodium Carbonate Anhydrate
B578 Sodium Cetylsulphate
B579 Sodium Citrate
B580 Sodium Citrate Dihydrate
B581 Sodium Cyclamate
B582 Sodium Dihydrogenphosphate
B583 Sodium Dichlorisocyanurate (Dikonit)
B584 Sodium Hydrogencitrate
B585 Sodium Hydrogenphosphate Dihydrate

NICODOM NIR Pharmaceuticals, 696 NIR spectra of pharmaceuticals
Copyright © NICODOM 2005

NICODOM Ltd., Hlavni 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

ftir@ftir.cz

B586 Sodium Hydrogensulphate
B587 Sodium Chloride
B588 Sodium Iodide
B589 Sodium Laurylsulphate
B590 Sodium Metabisulphite
B591 Sodium Perborate
B592 Sodium Phosphate dodecahydrate
B593 Sodium Picosulphate
B594 Sodium Sulphate
B595 Sodium Sulphate Anhydrate
B596 Sodium Sulphate Anhydrate
B597 Sodium Tetraborate
B598 Sodium Tetraborate Decahydrate
B599 Sodium Tetrahydroborate
B600 Sodium Thiosulphate
B601 Sodium-L-Glutamate
B602 Potato Starch Acetyl Oxidized
B603 Sorbic Acid
B604 Sorbitan Oleate
B605 Sorbitan Stearate (Arlacel 60)
B606 Sorbitol
B607 Sorbitol 70% (non crystallising)
B608 Spirit Caramel
B609 Sport Perfume 0062740
B610 Srodalan LA
B611 Stearic Acid
B612 Stearine
B613 Stearyl Fumarate Sodium
B614 Streptomycine Sulphate
B615 Struktol SB 2020
B616 Sucrose
B617 Sucrose
B618 Sucrose
B619 Sucrose (powder)
B620 Sulfadimidine Sodium
B621 Sulfanilic Acid
B622 Sulfathiazol
B623 Sulfoguaiacolum
B624 Sulphosalicylic Acid
B625 Sulphur Colloid
B626 Sulpiride
B627 Synthalen M
B628 Syntopon A 100
B629 Syrup Fructose
B630 Talc
B631 Tannin
B632 Tartaric Acid
B633 Tartrazine C-301
B634 Tartrazine E 102, CI 19140

NICODOM NIR Pharmaceuticals, 696 NIR spectra of pharmaceuticals
Copyright © NICODOM 2005

NICODOM Ltd., Hlavni 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

ftir@ftir.cz

B635 Tegin 4100
B636 Tego Alcanol 1618
B637 Tego Betain CK KB 5
B638 Tegosoft CT (cetearyl ethylhexanoate)
B639 Tegosoft EE
B640 Oil of Turpentine
B641 Oil of Turpentine Rectified
B642 Terfenadine
B643 Tetracaine Chloride
B644 Tetracycline Chloride
B645 Thalaspheeres (natural collagen)
B646 Theophylline
B647 Thiamine Chloride (Vitamin B1)
B648 Thioridazine
B649 Thiourea
B650 Oil of Thyme
B651 Thymol
B652 Tiaprofenic Acid
B653 Titanium Dioxide
B654 Tocopherol Acetate (50%) CWS/F (Vitamin E)
B655 Tragacantha
B656 Tramadole Chloride
B657 Triamcinolone Acetonide
B658 Tridecanamine Undecylenate
B659 Triethyl Citrate
B660 Triethylene Glycol
B661 Trimecaine
B662 Trimecaine Chloride
B663 Trimethoprim
B664 Tropical Extrakt
B665 Tropical Seaplant Extract
B666 Troxerutine
B667 Urea
B668 Ursodeoxycholic Acid
B669 Uvinul M 40
B670 Uvinul MS-40
B671 Valerianae Dry Extract
B672 Tincture of Valerian
B673 D-Valine
B674 Vaseline Oil 100 BDF
B675 Vaseline Pionier 0030-40 SYN
B676 Vaseline White Vara 4042 DAB10
B677 Veegum (coll. magnesium aluminum silicate)
B678 Verapamil Chloride
B679 Vinylpyrrolidon-vinylacetate Copolym. (VA64)
B680 Viprynum Embonate
B681 Brilliant Green
B682 Vitazyme C
B683 Wax

NICODOM NIR Pharmaceuticals, 696 NIR spectra of pharmaceuticals
Copyright © NICODOM 2005

NICODOM Ltd., Hlavni 2727, CZ-14100 Praha 4, Czech Republic, EU

Tel: +420-281914970, Fax: +420-281914971

<http://www.ir-spectra.com/>

ftir@ftir.cz

B684 WMK (Watermelon Keton)
B685 Xanthan Gum
B686 Xantiniol Nicotinate
B687 Xylitol
B688 Xylose
B689 Yellow orange S SY 75%
B690 Yellow orange S SY, E 110, CI 15958
B691 Yohimbine Chloride
B692 Zinc Stearate
B693 Zinc Sulfate
B694 Zinc Undecylenate
B695 Zinc Lactate
B696 Zinc Carbonate Basic

Last update: 18.09.07